



INVACARE®

# Matrix® Genera

Back



MODEL NO.	A		B		C		D		E		F		TRUNK WIDTH WITH CUSHION		CUSHION DEPTH		OVERALL WIDTH WITH CUSHION		OVERALL HEIGHT WITH CUSHION		OVERALL DEPTH WITH CUSHION		BACK WEIGHT (W/O HDWR)	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lbs	kgs
							N/A				N/A													
GB1616	15	38	15	38	4	10	N/A		11	28	N/A		13	33	3½	9	15½	39½	16	41	5½	14	3.8	1.7
GB1620	15	38	19	48	4	10	N/A		12½	32	N/A		13	33	3½	9	15½	39½	20	51	5½	14	4.3	1.9
GB1816	17	43	15	38	4	10	N/A		11	28	N/A		15	38	3½	9	17½	44½	16	41	5½	14	4.2	1.9
GB1820	17	43	19	48	4	10	N/A		12½	32	N/A		15	38	3½	9	17½	44½	20	51	5½	14	4.7	2.1
GB2016	19	48	15	38	4	10	N/A		11	28	N/A		17	43	3½	9	19½	49½	16	41	5½	14	4.5	2.0
GB2020	19	48	19	48	4	10	N/A		12½	32	N/A		17	43	3½	9	19½	49½	20	51	5½	14	5.0	2.3