





INVACARE®
Matrix® MX1
 Back



MODEL NO.	A		B		C		D		E		F		TRUNK WIDTH WITH CUSHION		CUSHION DEPTH		OVERALL WIDTH WITH CUSHION		OVERALL HEIGHT WITH CUSHION		OVERALL DEPTH WITH CUSHION		BACK WEIGHT (W/O HDWR)	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lbs	kgs
	MX1409	13	33	8	20	3½	9	N/A	5½	14	1½	4	11½	29	3	7½	14	36	10	26	6	15½	1.1	0.5
MX1412	13	33	11	28	3½	9	6½		16½	2½	6½	11½	29	3	7½	14	36	13	33	6	15½	1.2	0.5	
MX1416	13	33	15	38	3½	9	6½		16½	3	7½	11½	29	3	7½	14	36	17	43	6	15½	1.3	0.6	
MX1509	14	36	8	20	3½	9	N/A	5½	14	1½	4	12½	32	3	7½	15	38	10	26	6	15½	1.1	0.5	
MX1512	14	36	11	28	3½	9		5½	14	3½	9	12½	32	3	7½	15	38	13	33	6	15½	1.2	0.5	
MX1516	14	36	15	38	3½	9		7	18	3½	9	12½	32	3	7½	15	38	17	43	6	15½	1.4	0.6	
MX1609	15	38	8	20	3½	9	N/A	5½	14	1½	4	13½	34½	3	7½	16	41	10	26	6	15½	1.2	0.5	
MX1612	15	38	11	28	3½	9		5½	14	3½	9	13½	34½	3	7½	16	41	13	33	6	15½	1.4	0.6	
MX1616	15	38	15	38	3½	9		7	18	3½	9	13½	34½	3	7½	16	41	17	43	6	15½	1.6	0.7	
MX1709	16	41	8	20	3½	9	N/A	5½	14	1½	4	14½	37	3	7½	17	43	10	26	6	15½	1.2	0.5	
MX1712	16	41	11	28	3½	9		5½	14	3½	9	14½	37	3	7½	17	43	13	33	6	15½	1.4	0.6	
MX1716	16	41	15	38	3½	9		7	18	3½	9	14½	37	3	7½	17	43	17	43	6	15½	1.7	0.8	
MX1809	17	43	8	20	3½	9	N/A	5½	14	1½	4	15½	39½	3	7½	18	46	10	26	6	15½	1.3	0.6	
MX1812	17	43	11	28	3½	9		5½	14	3½	9	15½	39½	3	7½	18	46	13	33	6	15½	1.6	0.7	
MX1816	17	43	15	38	3½	9		7	18	3½	9	15½	39½	3	7½	18	46	17	43	6	15½	1.9	0.9	