


INVACARE®

# Matrix® PB Deep Heavy-Duty Back



MODEL NO.	A		B		C		D		E		F		TRUNK WIDTH WITH CUSHION		CUSHION DEPTH		OVERALL WIDTH WITH CUSHION		OVERALL HEIGHT WITH CUSHION		OVERALL DEPTH WITH CUSHION		BACK WEIGHT (W/O HDWR)	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lbs	kgs
	DHD2016	19	48	15½	39½	8	20	15	38	10	26	1	2½	18½	47	7	18	21	53	16	41	9	23	6.2
DHD2020	19	48	19½	48½	8	20	15	38	11	28	1	2½	18½	47	7	18	21	53	20	51	9	23	7.2	3.3
DHD2216	21	53	15½	39½	8	20	17	43	10	26	1	2½	20½	52	7	18	23	58	16	41	9	23	6.8	3.1
DHD2220	21	53	19½	48½	8	20	17	43	12	31	1	2½	20½	52	7	18	23	58	20	51	9	23	7.8	3.5
DHD2416	23	58	15½	39½	8	20	19	48	10	26	1	2½	22½	57	7	18	25	63	16	41	9	23	7.4	3.4
DHD2420	23	58	19½	48½	8	20	19	48	13	33	1	2½	22½	57	7	18	25	63	20	51	9	23	8.4	3.8
DHD2616	25	63	15½	39½	8	20	21	53	10	26	1	2½	24½	62	7	18	27	69	16	41	9	23	8.0	3.6
DHD2620	25	63	19½	48½	8	20	21	53	13	33	1	2½	24½	62	7	18	27	69	20	51	9	23	9.0	4.1
DHD2816	27	69	15½	39½	8	20	23	58	10	26	1	2½	26½	67	7	18	29	74	16	41	9	23	8.6	3.9
DHD2820	27	69	19½	48½	8	20	23	58	13	33	1	2½	26½	67	7	18	29	74	20	51	9	23	9.6	4.4
DHD3016	29	74	15½	39½	8	20	25	63	10	26	1	2½	28½	72	7	18	31	79	16	41	9	23	9.2	4.2
DHD3020	29	74	19½	48½	8	20	25	63	13	33	1	2½	28½	72	7	18	31	79	20	51	9	23	10.2	4.6