

AVAILABLE
NOW!

407 *maxx*

The
MATTER of
the **HEART**
EXCLUSIVELY ON

ROVI

A3

*7° of Anterior Assist allows for
repositioning and safe forward
transfer capability*



The all-NEW 407 Maxx - Single Pivot Tilt Power Positioning System from Motion Concepts was exclusively designed for the **NEW A3 Power Base from ROVI**. The simplicity of this system optimizes the geometry of the A3 base for maximum stability while offering up to 40° of Posterior Tilt with 7° of Anterior Assist as a standard configuration. With the addition of a Recline (170°) module, and LNX power articulating center mount foot platform, the **407 Maxx** system allows the client to achieve a unique position where the head is below and feet above the heart to assist with respiratory and circulatory conditions. The simple but effective design of the **407 Maxx** system maintains all the features, flexibility and good looks you have come to expect from Motion Concepts products.

Call to order

888.433.6818
motionconcepts.com



ROVI A3



407maxx



BASE VITALS

ROVI A3

HCPCS Codes	K0856, K0861
Suspension System	Active Ride Control (ARC)
Base Width	25.5"
Base Length	36.5"
Drive Wheel Size ⁽¹⁾	14"
Caster Wheel Size (front/rear)	6"
Base only Turning Radius	20.5"
Drivetrain	4-pole EAD Motors
Batteries	2 x M34 3 x M34* 2 x M34 and 1 x M50* (*option for Vent application)
Battery Charger	8A, Off Board
Batteries Weight (each)	38.5 lbs (M34)
Range ⁽²⁾	17.7 miles
Weight Capacity	300 lbs
Base Weight (w/o batteries)	159 lbs
Ground Clearance	3"
Maximum Speed	6.2 mph

⁽¹⁾ All A3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

SYSTEM VITALS

- Featuring up to 40° of Posterior Tilt with 7° of Anterior Assist
- Up to 170° of Recline
- The unique design of Single Pivot Tilt coupled with Recline module allows up to 185° of combined travel and enables a client to adopt a Trendelenburg position where feet are elevated above their head
- 16" - 22" W x 16" - 22"D, fully adjustable
- Optional LNX Power Center Mount Foot Platform or Fixed Center Mount Legrest
- Available with the full range of Motion Concepts' accessories and Invacare® Matrx® Seating products, mounted via our Ultra Rail and Ultra Track slotted interfaces
- Seat-to-Floor heights as low as 17"
- Up to 250 lb. weight capacity



Motion Concepts USA
700 Ensminger Rd., Suite 112
Tonawanda, New York 14150
888.433.6818

Motion Concepts CANADA
84 Citation Dr., Unit 1
Concord, ON L4K 1C3
866.748.7943

www.motionconcepts.com
info@motionconcepts.com

© 2019 Motion Concepts. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Motion Concepts unless otherwise noted. Specifications are subject to change without notification.

CONTACT YOUR SALES PROFESSIONAL

