


MATRIX® E2™ DEEP BACK BACK SPECIFICATIONS



E2D BACK SIZE	A BACK SHELL HEIGHT	B OVERAL BACK HEIGHT	C BACK SHELL WIDTH	D OVERAL INSIDE WIDTH	E BACK CANE WIDTH	F APEX OF LATERAL	G SLOT POSITION	H SLOT LENGTH
1406	6"/15cm	7.5"/191mm	12.75"/324mm	10.5"/267mm	13"-16"/33-41cm	3"/76mm	1.40"/36mm	3.25"/83mm
1408	8"/20cm	9.5"/241mm	12.75"/324mm	10.5"/267mm	13"-16"/33-41cm	4"/102mm	1.80"/46mm	3.50"/89mm
1410	10"/25cm	11.5"/292mm	12.75"/324mm	10.5"/267mm	13"-16"/33-41cm	6"/152mm	3.30"/84mm	4.00"/102mm
1412	12"/30cm	13.5"/343mm	12.75"/324mm	10.5"/267mm	13"-16"/33-41cm	7"/178mm	3.30"/84mm	4.39"/112mm
1414	14"/36cm	15.5"/393mm	12.75"/324mm	10.5"/267mm	13"-16"/33-41cm	7.5"/191mm	3.50"/89mm	5.57"/141mm
1416	16"/41cm	17.5"/446mm	12.75"/324mm	10.5"/267mm	13"-16"/33-41cm	8.25"/210mm	4.50"/114mm	5.57"/141mm
1418	18"/46cm	19.5"/495mm	12.75"/324mm	11.5"/292mm	13"-16"/33-41cm	8.25"/210mm	4.50"/114mm	6.50"/165mm
1420	20"/51cm	21.5"/546mm	12.75"/324mm	11.5"/292mm	13"-16"/33-41cm	8.25"/210mm	4.50"/114mm	6.50"/165mm
1506	6"/15cm	7.5"/191mm	13.75"/349mm	11.5"/292mm	14"-17"/36-43cm	3"/76mm	1.40"/36mm	3.25"/83mm
1508	8"/20cm	9.5"/241mm	13.75"/349mm	11.5"/292mm	14"-17"/36-43cm	4"/102mm	1.80"/46mm	3.50"/89mm
1510	10"/25cm	11.5"/292mm	13.75"/349mm	11.5"/292mm	14"-17"/36-43cm	6"/152mm	3.30"/84mm	4.00"/102mm
1512	12"/30cm	13.5"/343mm	13.75"/349mm	11.5"/292mm	14"-17"/36-43cm	7"/178mm	3.30"/84mm	4.39"/112mm
1514	14"/36cm	15.5"/393mm	13.75"/349mm	11.5"/292mm	14"-17"/36-43cm	7.5"/191mm	3.50"/89mm	5.57"/141mm
1516	16"/41cm	17.5"/446mm	13.75"/349mm	11.5"/292mm	14"-17"/36-43cm	8.25"/210mm	4.50"/114mm	5.57"/141mm
1518	18"/46cm	19.5"/495mm	13.75"/349mm	11.5"/292mm	14"-17"/36-43cm	8.25"/210mm	4.50"/114mm	6.50"/165mm
1520	20"/51cm	21.5"/546mm	13.75"/349mm	11.5"/292mm	14"-17"/36-43cm	8.25"/210mm	4.50"/114mm	6.50"/165mm
1606	6"/15cm	7.5"/191mm	14.75"/375mm	12.5"/318mm	15"-18"/38-46cm	3"/76mm	1.40"/36mm	3.25"/83mm
1608	8"/20cm	9.5"/241mm	14.75"/375mm	12.5"/318mm	15"-18"/38-46cm	4"/102mm	1.80"/46mm	3.50"/89mm
1610	10"/25cm	11.5"/292mm	14.75"/375mm	12.5"/318mm	15"-18"/38-46cm	6"/152mm	3.30"/84mm	4.00"/102mm
1612	12"/30cm	13.5"/343mm	14.75"/375mm	12.5"/318mm	15"-18"/38-46cm	7"/178mm	3.30"/84mm	4.39"/112mm
1614	14"/36cm	15.5"/393mm	14.75"/375mm	12.5"/318mm	15"-18"/38-46cm	7.5"/191mm	3.50"/89mm	5.57"/141mm
1616	16"/41cm	17.5"/446mm	14.75"/375mm	12.5"/318mm	15"-18"/38-46cm	8.25"/210mm	4.50"/114mm	5.57"/141mm
1618	18"/46cm	19.5"/495mm	14.75"/375mm	12.5"/318mm	15"-18"/38-46cm	8.25"/210mm	4.50"/114mm	6.50"/165mm
1620	20"/51cm	21.5"/546mm	14.75"/375mm	12.5"/318mm	15"-18"/38-46cm	8.25"/210mm	4.50"/114mm	6.50"/165mm

MORE SIZE SPECIFICATIONS ON THE BACK



MATRIX®
E2™ DEEP BACK
BACK SPECIFICATIONS

E2D BACK SIZE	A BACK SHELL HEIGHT	B OVERAL BACK HEIGHT	C BACK SHELL WIDTH	D OVERAL INSIDE WIDTH	E BACK CANE WIDTH	F APEX OF LATERAL	G SLOT POSITION	H SLOT LENGTH
1706	6"/15cm	7.5"/191mm	15.75"/400mm	13.5"/343mm	16-19"/41-48cm	3"/76mm	1.4"/36mm	3.25"/83mm
1708	8"/20cm	9.5"/241mm	15.75"/400mm	13.5"/343mm	16-19"/41-48cm	4"/102mm	1.8"/46mm	3.50"/89mm
1710	10"/25cm	11.5"/292mm	15.75"/400mm	13.5"/343mm	16-19"/41-48cm	6"/152mm	3.3"/84mm	4"/102mm
1712	12"/30cm	13.5"/343mm	15.75"/400mm	13.5"/343mm	16-19"/41-48cm	7"/178mm	3.3"/84mm	4.39"/112mm
1714	14"/36cm	15.5"/393mm	15.75"/400mm	13.5"/343mm	16-19"/41-48cm	7.5"/191mm	3.5"/89mm	5.57"/141mm
1716	16"/41cm	17.5"/446mm	15.75"/400mm	13.5"/343mm	16-19"/41-48cm	8.25"/210mm	4.5"/114mm	5.57"/141mm
1718	18"/46cm	19.5"/495mm	15.75"/400mm	13.5"/343mm	16-19"/41-48cm	8.25"/210mm	4.5"/114mm	6.5"/165mm
1720	20"/51cm	21.5"/546mm	15.75"/400mm	13.5"/343mm	16-19"/41-48cm	8.25"/210mm	4.50"/114mm	6.5"/165mm
1806	6"/15cm	7.5"/191mm	16.75"/425mm	14.5"/368mm	17-20"/43-51cm	3"/76mm	1.4"/36mm	3.25"/83mm
1808	8"/20cm	9.5"/241mm	16.75"/425mm	14.5"/368mm	17-20"/43-51cm	4"/102mm	1.8"/46mm	3.5"/89mm
1810	10"/25cm	11.5"/292mm	16.75"/425mm	14.5"/368mm	17-20"/43-51cm	6"/152mm	3.3"/84mm	4"/102mm
1812	12"/30cm	13.5"/343mm	16.75"/425mm	14.5"/368mm	17-20"/43-51cm	7"/178mm	3.3"/84mm	4.39"/112mm
1814	14"/36cm	15.5"/393mm	16.75"/425mm	14.5"/368mm	17-20"/43-51cm	7.5"/191mm	3.5"/89mm	5.57"/141mm
1816	16"/41cm	17.5"/446mm	16.75"/425mm	14.5"/368mm	17-20"/43-51cm	8.25"/210mm	4.5"/114mm	5.57"/141mm
1818	18"/46cm	19.5"/495mm	16.75"/425mm	14.5"/368mm	17-20"/43-51cm	8.25"/210mm	4.5"/114mm	6.5"/165mm
1820	20"/51cm	21.5"/546mm	16.75"/425mm	14.5"/368mm	17-20"/43-51cm	8.25"/210mm	4.5"/114mm	6.5"/165mm
1906	6"/15cm	7.5"/191mm	17.75"/451mm	15.5"/393mm	18-21"/46-53cm	3"/76mm	1.4"/36mm	3.25"/83mm
1908	8"/20cm	9.5"/241mm	17.75"/451mm	15.5"/393mm	18-21"/46-53cm	4"/102mm	1.8"/46mm	3.50"/89mm
1910	10"/25cm	11.5"/292mm	17.75"/451mm	15.5"/393mm	18-21"/46-53cm	6"/152mm	3.3"/84mm	4"/102mm
1912	12"/30cm	13.5"/343mm	17.75"/451mm	15.5"/393mm	18-21"/46-53cm	7"/178mm	3.3"/84mm	4.39"/112mm
1914	14"/36cm	15.5"/393mm	17.75"/451mm	15.5"/393mm	18-21"/46-53cm	7.5"/191mm	3.5"/89mm	5.57"/141mm
1916	16"/41cm	17.5"/446mm	17.75"/451mm	15.5"/393mm	18-21"/46-53cm	8.25"/210mm	4.5"/114mm	5.57"/141mm
1918	18"/46cm	19.5"/495mm	17.75"/451mm	15.5"/393mm	18-21"/46-53cm	8.25"/210mm	4.5"/114mm	6.5"/165mm
1920	20"/51cm	21.5"/546mm	17.75"/451mm	15.5"/393mm	18-21"/46-53cm	8.25"/210mm	4.5"/114mm	6.5"/165mm
2006	6"/15cm	7.5"/191mm	18.75"/476mm	16.5"/419mm	19-22"/48-56cm	3"/76mm	1.4"/36mm	3.25"/83mm
2008	8"/20cm	9.5"/241mm	18.75"/476mm	16.5"/419mm	19-22"/48-56cm	4"/102mm	1.8"/46mm	3.5"/89mm
2010	10"/25cm	11.5"/292mm	18.75"/476mm	16.5"/419mm	19-22"/48-56cm	6"/152mm	3.3"/84mm	4"/102mm
2012	12"/30cm	13.5"/343mm	18.75"/476mm	16.5"/419mm	19-22"/48-56cm	7"/178mm	3.3"/84mm	4.39"/112mm
2014	14"/36cm	15.5"/393mm	18.75"/476mm	16.5"/419mm	19-22"/48-56cm	7.5"/191mm	3.5"/89mm	5.57"/141mm
2016	16"/41cm	17.5"/446mm	18.75"/476mm	16.5"/419mm	19-22"/48-56cm	8.25"/210mm	4.5"/114mm	5.57"/141mm
2018	18"/46cm	19.5"/495mm	18.75"/476mm	16.5"/419mm	19-22"/48-56cm	8.25"/210mm	4.5"/114mm	6.5"/165mm
2020	20"/51cm	21.5"/546mm	18.75"/476mm	16.5"/419mm	19-22"/48-56cm	8.25"/210mm	4.5"/114mm	6.5"/165mm



Motion Concepts
www.motionconcepts.com

USA
 700 Ensminger Rd. Suite 112
 Tonawanda, NY 14150
 (888) 433-6818
 Fax: (888) 433-6834

Canada
 84 Citation Dr.,
 Unit 1-12
 Concord, ON L4K 3C1
 (866) 748-7943
 (905) 695-0134
 Fax: (905) 695-0138

©2022 Motion Concepts.
 All rights reserved. Trademarks are
 identified by the symbols ™ and ®. All
 trademarks are owned by or licensed to
 Motion Concepts unless otherwise noted.
 Specifications are subject to change
 without notification.
 Form No. TRD0603 REV B