


MATRIX®

MX2 BACK

BACK SPECIFICATIONS



MX2 BACK SIZE	A BACK SHELL HEIGHT	B OVERAL BACK HEIGHT	C BACK SHELL WIDTH	D OVERAL INSIDE WIDTH	E BACK CANE WIDTH	F APEX OF LATERAL	G SLOT POSITION	H SLOT LENGTH
1406	5"/12.5cm	7"/178mm	12.5"/318mm	10.5"/267mm	13"-16"/33-41cm	2.5"/64mm	1.38"/35mm	2.25"/57mm
1409	8"/20cm	10"/254mm	12.5"/318mm	10.5"/267mm	13"-16"/33-41cm	4"/102mm	2.13"/64mm	3.5"/89mm
1412	11"/28cm	13"/330mm	12.5"/318mm	10.5"/267mm	13"-16"/33-41cm	6.25"/159mm	3"/76mm	4"/102mm
1416	15"/38cm	17"/432mm	12.5"/318mm	10.5"/267mm	13"-16"/33-41cm	7.25"/184mm	3.25"/83mm	5.5"/140mm
1506	5"/12.5cm	7"/178mm	13.5"/343mm	11.5"/292mm	14"-17"/36-43cm	2.5"/64mm	1.38"/35mm	2.25"/57mm
1509	8"/20cm	10"/254mm	13.5"/343mm	11.5"/292mm	14"-17"/36-43cm	4"/102mm	2.13"/64mm	3.5"/89mm
1512	11"/28cm	13"/330mm	13.5"/343mm	11.5"/292mm	14"-17"/36-43cm	6.25"/159mm	3"/76mm	4"/102mm
1516	15"/38cm	17"/432mm	13.5"/343mm	11.5"/292mm	14"-17"/36-43cm	7.25"/184mm	3.25"/83mm	5.5"/140mm
1606	16"/41cm	7"/178mm	14.5"/368mm	12.5"/318mm	15"-18"/38-46cm	2.5"/64mm	1.38"/35mm	2.25"/57mm
1609	5"/12.5cm	10"/254mm	14.5"/368mm	12.5"/318mm	15"-18"/38-46cm	4"/102mm	2.13"/64mm	3.5"/89mm
1612	11"/28cm	13"/330mm	14.5"/368mm	12.5"/318mm	15"-18"/38-46cm	6.25"/159mm	3"/76mm	4"/102mm
1616	15"/38cm	17"/432mm	14.5"/368mm	12.5"/318mm	15"-18"/38-46cm	7.25"/184mm	3.25"/83mm	5.5"/140mm
1706	5"/13cm	7"/178mm	15.5"/393mm	13.5"/343mm	16-19"/41-48cm	2.5"/64mm	1.38"/35mm	2.25"/57mm
1709	8"/20cm	10"/254mm	15.5"/393mm	13.5"/343mm	16-19"/41-48cm	4"/102mm	2.13"/64mm	3.5"/89mm
1712	11"/28cm	13"/330mm	15.5"/393mm	13.5"/343mm	16-19"/41-48cm	6.25"/159mm	3"/76mm	4"/102mm
1716	15"/38cm	17"/432mm	15.5"/393mm	13.5"/343mm	16-19"/41-48cm	7.25"/184mm	3.25"/83mm	5.5"/140mm
1806	5"/12.5cm	7"/178mm	16.5"/419mm	14.5"/368mm	17-20"/43-51cm	2.5"/64mm	1.38"/35mm	2.25"/57mm
1809	8"/20cm	10"/254mm	16.5"/419mm	14.5"/368mm	17-20"/43-51cm	4"/102mm	2.13"/64mm	3.5"/89mm
1812	11"/28cm	13"/330mm	16.5"/419mm	14.5"/368mm	17-20"/43-51cm	6.25"/159mm	3"/76mm	4"/102mm
1816	15"/38cm	17"/432mm	16.5"/419mm	14.5"/368mm	17-20"/43-51cm	7.25"/184mm	3.25"/83mm	5.5"/140mm



Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr.,
Unit 1-12
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: (905) 695-0138

©2022 Motion Concepts.
All rights reserved. Trademarks are identified by the symbols ™ and ©. All trademarks are owned by or licensed to Motion Concepts unless otherwise noted. Specifications are subject to change without notification.
TRD0681 Rev A