


MATRIX® MAC BACK SPECIFICATIONS



MAC BACK SIZE	A* OVERALL BACK LENGTH	B* LOADED INNER CONTOUR WIDTH RANGE	C ADJUSTMENT RANGE TO APEX OF LATERAL	D MOUNTING SLOT LENGTH	E MOUNTING SLOT POSITION
1414	15.5" /356mm	10.5-11.5" /267-292mm	5.5-7.5" /140-191mm	5.5" /140mm	3.5" /89mm
1416	17.5" /406mm	10.5-11.5" /267-292mm	5.5-7.5" /140-191mm	5.5" /140mm	3.5" /89mm
1418	19.5" /457mm	10.5-11.5" /267-292mm	5.5-9.75" /140-248mm	6.75" /171mm	4.5" /114mm
1420	21.5" /508mm	10.5-11.5" /267-292mm	5.5-9.75" /140-248mm	6.75" /171mm	4.5" /114mm
1514	15.5" /356mm	11.5-12.5" /292-318mm	5.5-7.5" /140-191mm	5.5" /140mm	3.5" /89mm
1516	17.5" /406mm	11.5-12.5" /292-318mm	5.5-7.5" /140-191mm	5.5" /140mm	3.5" /89mm
1518	19.5" /457mm	11.5-12.5" /292-318mm	5.5-9.75" /140-248mm	6.75" /171mm	4.5" /114mm
1520	21.5" /508mm	11.5-12.5" /292-318mm	5.5-9.75" /140-248mm	6.75" /171mm	4.5" /114mm
1616	17.5" /406mm	12.5-13.5" /318-343mm	5.5-7.5" /140-191mm	5.5" /140mm	3.5" /89mm
1618	19.5" /457mm	12.5-13.5" /318-343mm	5.5-8.75" /140-222mm	6.5" /165mm	4.5" /114mm
1620	21.5" /508mm	12.5-13.5" /318-343mm	5.5-8.75" /140-222mm	6.5" /165mm	4.5" /114mm
1622	23.5" /559mm	12.5-13.5" /318-343mm	5.5-9.75" /140-248mm	6.5" /165mm	4.5" /114mm
1624	25.5" /610mm	12.5-13.5" /318-343mm	5.5-9.75" /140-248mm	6.5" /165mm	4.5" /114mm
1716	17.5" /406mm	13.5-14.5" /343-368mm	5.5-7.5" /140-191mm	5.5" /140mm	3.5" /89mm
1718	19.5" /457mm	13.5-14.5" /343-368mm	5.5-8.75" /140-222mm	6.5" /165mm	4.5" /114mm
1720	21.5" /508mm	13.5-14.5" /343-368mm	5.5-8.75" /140-222mm	6.5" /165mm	4.5" /114mm
1722	23.5" /559mm	13.5-14.5" /343-368mm	5.5-9.75" /140-248mm	6.5" /165mm	4.5" /114mm
1724	25.5" /610mm	13.5-14.5" /343-368mm	5.5-9.75" /140-248mm	6.5" /165mm	4.5" /114mm

*Please note measurements for soft components are approximate.




Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
505 Cityview Blvd.,
Unit 1
Vaughan, ON L4H 0L8
Canada
(866) 748-7943
(905) 695-0134
Fax: (905) 695-0138

©2026 Motion Concepts.
All rights reserved. Trademarks are identified by the symbols ™ and ©. All trademarks are owned by or licensed to Motion Concepts unless otherwise noted. Specifications are subject to change without notification.
Form No. TRD0741 Rev B

MATRIX® MAC BACK SPECIFICATIONS



MAC BACK SIZE	A* OVERALL BACK LENGTH	B* LOADED INNER CONTOUR WIDTH RANGE	C ADJUSTMENT RANGE TO APEX OF LATERAL	D MOUNTING SLOT LENGTH	E MOUNTING SLOT POSITION
1818	19.5" / 457mm	14.5-15.5" / 368-394mm	5.5-8.75" / 140-222mm	6.5" / 165mm	4.5" / 114mm
1820	21.5" / 508mm	14.5-15.5" / 368-394mm	5.5-8.75" / 140-222mm	6.5" / 165mm	4.5" / 114mm
1822	23.5" / 559mm	14.5-15.5" / 368-394mm	5.5-9.75" / 140-248mm	6.5" / 165mm	4.5" / 114mm
1824	25.5" / 610mm	14.5-15.5" / 368-394mm	5.5-9.75" / 140-248mm	6.5" / 165mm	4.5" / 114mm
1916	17.5" / 406mm	15.5-16.5" / 394-419mm	5.5-7.5" / 140-191mm	5.5" / 140mm	3.5" / 89mm
1918	19.5" / 457mm	15.5-16.5" / 394-419mm	5.5-8.75" / 140-222mm	6.5" / 165mm	4.5" / 114mm
1920	21.5" / 508mm	15.5-16.5" / 394-419mm	5.5-8.75" / 140-222mm	6.5" / 165mm	4.5" / 114mm
1922	23.5" / 559mm	15.5-16.5" / 394-419mm	5.5-9.75" / 140-248mm	6.5" / 165mm	4.5" / 114mm
1924	25.5" / 610mm	15.5-16.5" / 394-419mm	5.5-9.75" / 140-248mm	6.5" / 165mm	4.5" / 114mm
2016	17.5" / 406mm	16.5-17.5" / 419-445mm	5.5-7.5" / 140-191mm	5.5" / 140mm	3.5" / 89mm
2018	19.5" / 457mm	16.5-17.5" / 419-445mm	5.5-8.75" / 140-222mm	6.5" / 165mm	4.5" / 114mm
2020	21.5" / 508mm	16.5-17.5" / 419-445mm	5.5-8.75" / 140-222mm	6.5" / 165mm	4.5" / 114mm
2022	23.5" / 559mm	16.5-17.5" / 419-445mm	5.5-9.75" / 140-248mm	6.5" / 165mm	4.5" / 114mm
2024	25.5" / 610mm	16.5-17.5" / 419-445mm	5.5-9.75" / 140-248mm	6.5" / 165mm	4.5" / 114mm

*Please note measurements for soft components are approximate.



Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
505 Cityview Blvd.,
Unit 1
Vaughan, ON L4H 0L8
Canada
(866) 748-7943
(905) 695-0134
Fax: (905) 695-0138

©2026 Motion Concepts.
All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Motion Concepts unless otherwise noted. Specifications are subject to change without notification.
Form No. TRD0741 Rev B