

REQUIRED INFORMATION

Quote ☐ Order ☐

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured
please fill in all required

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



Customer Service: 1.866.748.7943 | Fax: 1.905.695.0583 | orderscdn@motionconcepts.com

ROVI BASE STATISTICS

SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH A ₃	25.25"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"
CASTER WHEEL SIZE (front/rear)	6"
TURNING RADIUS A ₃	20.5"
DRIVETRAIN	4-pole EAD Motors

BATTERIES	2 x M34
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs (M34) / 17.5 kg
RANGE ⁽²⁾	17.7 miles / 28.5 km/h
WEIGHT CAPACITY	300 lbs / 136 kg
BASE WEIGHT A ₃ (w/o batteries)	159 lbs / 72 kg
GROUND CLEARANCE	3"
MAXIMUM SPEED	5.8 mph

A3



*Shown with Motion Concepts Single Post Back for Contour Backrest. See page 7

⁽¹⁾ All X3 and A3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

BASE OPTIONS

☐ 888WNLEA3 ROVI A3⁽¹⁾⁽²⁾ \$6,800

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt. Weight capacity is 300 lbs.
2. When selecting an increased speed motor package there is a reduction in torque output by roughly 15%.

FENDER COLOUR

☐ 800914-00200-03	☐ Blizzard White
☐ 800914-00200-69	☐ Solar Yellow
☐ 800914-00200-23	☐ Original Orange
☐ 800914-00201-21	☐ Red Delicious
☐ 800914-00201-2C	☐ Party Pink
☐ BCF Black Carbon (Mat)	☐ Black Panther (Gloss)
☐ BPF	☐ Purple Reign
☐ PRF	☐ Indigo Blue
☐ IBF	☐ Motion Green
☐ MGF	

WHEELCHAIR TRANSPORT BRACKETS

☒ XBKTS Wheelchair Transport Brackets(1) STD
1. For unoccupied use only

TIRE OPTIONS

☐ STIRES Standard Solid Tires STD
☐ PTIRES Pneumatic Tires N/C

BATTERIES

☐ M34 SLD G MK Group Battery x 2 \$650

JOYSTICK SELECTION

☐ D51122 R-Net LED Joystick N/C

JOYSTICK MOUNT

☐ 800914-001-16-ML Fixed Joystick Mounting Left N/C
☐ 800914-001-16-MR Fixed Joystick Mounting Right N/C
☐ SAQL Height Adjustable Swing-Away Quad Link Left \$244
☐ SAQR Height Adjustable Swing-Away Quad Link Right \$244

ADDITIONAL ACCESSORIES

☐ MUSB Motion Dual USB Charger N/C
☐ UCH Universal Cup Holder⁽¹⁾ \$50

1. Not available on the joystick side.

SEAT-TO-FLOOR HEIGHT⁽¹⁾

☐ STF17 17" Seat-To-Floor N/C
☐ STF18 18" Seat-To-Floor N/C
☐ STF19 19" Seat-To-Floor N/C
☐ STF20 20" Seat-To-Floor N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"

REHAB SEAT

☐ RSMAXX Ultra Low Maxx Rehab Seat \$995

ADD: MODULES

☐ PRETILT 5° Fixed Posterior Tilt⁽¹⁾ N/C

1. This Fixed posterior tilt option adds approximately 2" to the front seat to floor height

ADULT SEAT WIDTH SETTINGS

☐ SW116 16" Wide (Adjustable from 16" to 20") N/C
☐ SW117 17" Wide (Adjustable from 16" to 20") N/C
☐ SW118 18" Wide (Adjustable from 16" to 20") N/C
☐ SW119 19" Wide (Adjustable from 16" to 20") N/C
☐ SW120 20" Wide (Adjustable from 16" to 20") N/C

1. This option is only available on the Rovi A3 Base.

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

☐ SD116 16" Deep (Adjustable from 16" to 21") N/C
☐ SD117 17" Deep (Adjustable from 16" to 21") N/C
☐ SD118 18" Deep (Adjustable from 16" to 21") N/C
☐ SD119 19" Deep (Adjustable from 16" to 21") N/C
☐ SD120 20" Deep (Adjustable from 16" to 21") N/C

LAP BELTS

☐ BELT48 Push Button Style Lap Belt, Length 48" STD
☐ BELT60 Push Button Style Lap Belt, Length 60" N/C
☐ PLB2 Bodypoint Padded Lap Belt 2 Point⁽¹⁾ \$250
☐ PLB4 Bodypoint Padded Lap Belt 4 Point⁽¹⁾ \$315

1. 67" length only option

MATRX SEAT CUSHION OPTIONS⁽¹⁾

☐ PSWWDD Matrx PS Seat Cushion \$569
☐ PSPWWDD Matrx PSP Seat Cushion \$645

1. Cushion size will match seat width and depth

FINISHED BACK HEIGHT ⁽¹⁾⁽²⁾

☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"	☐ 25"
FBT18	FBT19	FBT20	FBT21	FBT22	FBT23	FBT24	FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK CANES

☐ SBC18	18" Straight Back Canes	STD
☐ SBC20	20" Straight Back Canes	STD
☐ SBC22	22" Straight Back Canes	STD
☐ SBC24	24" Straight Back Canes	STD

Select Cane Angle

☐ 75°	☐ 79°	☐ 85°	☐ 90°	☐ 95°	☐ 101°	☐ 106°	☐ 112°	☐ 116°
BA75	BA79	BA85	BA90	BA95	BA101	BA106	BA112	BA116

BACK OPTIONS

☐ SRB	Standard Rehab Back & Cushion	\$125
☐ ORBC	Omit Rigid Back & Cushion	N/C

UPGRADABLE BACK OPTIONS⁽¹⁾

☐ E2SWW14	Matrix Elite E2 (3" Contour) 14" H	\$965
☐ E2SWW16	Matrix Elite E2 (3" Contour) 16" H	\$965
☐ E2SWW18	Matrix Elite E2 (3" Contour) 18" H	\$995
☐ E2SWW20	Matrix Elite E2 (3" Contour) 20" H	\$995
☐ E2DWW14	Matrix Elite E2D (6" Contour) 14" H	\$1,125
☐ E2DWW16	Matrix Elite E2D (6" Contour) 16" H	\$1,125
☐ E2DWW18	Matrix Elite E2D (6" Contour) 18" H	\$1,155
☐ E2DWW20	Matrix Elite E2D (6" Contour) 20" H	\$1,155

1. Back width will match seat width.

ARMRESTS

☐ DP1	Adult Dual Post 9.75" - 12.75"	STD
☐ DP2	Adult Dual Post 12.75" - 15.75"	STD
☐ DP3	Adult Dual Post 15.75" - 18.75"	STD

☐ RMCAN	Ultra Rail Seat Mnt Flip Back Cantilever MAXX Tilt Arm 9" - 12.5" ⁽¹⁾	\$488
--------------------------------	---	-------

1. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.

ARMREST TUBE LENGTH - PLEASE SPECIFY

☐ ATFL	Left Full	STD
☐ ATDL	Left Desk	STD
☐ ATFR	Right Full	STD
☐ ATDR	Right Desk	STD

ARMPADS**Standard**

☐ SAPFL	Left Full (2.25" x 14")	STD
☐ SAPDL	Left Desk (2.25" x 10")	STD
☐ SAPFR	Right Full (2.25" x 14")	STD
☐ SAPDR	Right Desk (2.25" x 10")	STD

Waterfall Arm Pads⁽¹⁾

☐ WAPFL	Left Full (13")	N/C
☐ WAPDL	Left Desk (9")	N/C
☐ WAPFR	Right Full (13")	N/C
☐ WAPDR	Right Desk (9")	N/C

Gel Pads

☐ GAPNL	Narrow Left (2" x 12")	\$88
☐ GAPNR	Narrow Right (2" x 12")	\$88
☐ GAPWL	Wide Left (3.5" x 12")	\$96
☐ GAPWR	Wide Right (3.5" x 12")	\$96
☐ GAPLL	Long Left (3.5" x 14")	\$103
☐ GAPLR	Long Right (3.5" x 14")	\$103

1. Full length pad requires full length arm tube if using Ultra Rail mounted Tilt Arms.

MANUAL CENTER MOUNT FRONT RIGGING



☐ MSFCMFP	Maxx Style Fixed Center Mount Foot Platform ⁽¹⁾	\$275
☐ IFP	Medium Rubber Coated Footplate (11.5" W x 10")	N/C

☐ ICP	Individual calf pads	\$100
Select seat pan to footplate:		
☐ SPTF1	9" - 13"	N/C
☐ SPTF14	13" - 19"	N/C

Set at:

☐ SA70	Set at 70°	N/C
☐ SA90	Set at 90°	N/C
☐ SA97	Set at 97°	N/C

1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.

**FOOT PLATFORM OPTIONS FOR MSFCMFP⁽¹⁾**

☐ SFFP	Small Platform 12"w x 8.5"d	\$180
☐ MFFP	Medium Platform 12"w x 10.75"d	\$180
☐ LFFP	Large Platform 14.5"w x 13"d	\$250
☐ XLFFP	Extra Large Platform 17"w x 13"d	\$250
☐ XXLFFP	Extra Extra Large Platform 21"w x 13"d	\$250

1. May need to adjust leg position to avoid interference with front castors

MANUAL ELEVATING CENTER MOUNT LEGREST

- ☐ **MECML** Manual Elevating Center Mount Legrest
With Plantar Flex. \$595
- Select legrest set up position:
- ☐ **90SU** Set at 90° (elevating range 90-40°) N/C
- ☐ **110TU** Tucked at 110° (elevating range 110-60°)⁽¹⁾⁽²⁾ N/C
- Rubber Coated Footplate:
- ☐ **MRCP** Rubber Coated Footplate (11.5" W x 10") N/C
- Select seat pan to footplate for Rubber Coated Footplate
- ☐ **CSPTF1** 11.5" - 15" N/C
- ☐ **CSPTF2** 14" - 19" N/C
- Optional:
- ☐ **CICP** Individual calf pads \$200

1. Not available on seat depths 16", 17".

SWING AWAY LEGREST

- ☐ **HD70** Heavy Duty 70° Swing Away Footrest⁽¹⁾⁽²⁾ \$495
- Select Receivers:
- ☐ **RS** Straight Receivers N/C
- ☐ **RF** Flared 1" Outward Receivers N/C
- Select Seat Pan to Footplate:
- ☐ **SPTF11** Length: 6"-10" N/C
- ☐ **SPTF8** Length: 11"-15" N/C
- ☐ **SPTF2** Length: 13"-17" N/C
- ☐ **SPTF7** Length: 16"-20" N/C
- Options:
- ☐ **HDFCP** Flat Calf Pad Option \$275
- ☐ **HD70CCP** Curved Calf Pad Option⁽³⁾ \$275
- ☐ **I70SAF** Invacare 70° Swing-away Footrests \$325
- ☐ **I70-CF** Comes with composite foot plates. N/C
- ☐ **70HDAIH** Invacare 70° Non-Tapered HD Swing-away
Footrests with AT5543 Adjustable Angle
Footplates, Heel Loops and Impact Guards \$457

1. Must choose footplate on page 6

2. Heavy duty durable construction (4.5lbs ea.).

3. For 16" and 17" seat widths you must select flared receivers.

SUBMIT

PRINT FORM

SAVE

CLEAR FORM

FOOTPLATE FOR SWING AWAY/ELEVATING LEGS**Adjustable Angle Footplates**

- ☐ **AFSL** Left Small Footplate (4.25"W x 8"D) \$165
- ☐ **AFSR** Right Small Footplate (4.25"W x 8"D) \$165
- ☐ **AFML** Left Medium Footplate (5.25"W x 8"D) \$165
- ☐ **AFMR** Right Medium Footplate (5.25"W x 8"D) \$165
- ☐ **AFL** Left Large Footplate (6.25"W x 8"D) \$165
- ☐ **AFLR** Right Large Footplate (6.25"W x 8"D) \$165
- ☐ **AFXL** Left Extra Large Footplate (9"W x 11"D) \$165
- ☐ **AFXR** Right Extra Large Footplate (9"W x 11"D) \$165

Multi-Axis Adjustable Angle Footplates

- ☐ **MFSL** Left Small Footplate (4.25"W x 8"D) \$225
- ☐ **MFSR** Right Small Footplate (4.25"W x 8"D) \$225
- ☐ **MFML** Left Medium Footplate (5.25"W x 8"D) \$225
- ☐ **MFMR** Right Medium Footplate (5.25"W x 8"D) \$225
- ☐ **MFL** Left Large Footplate (6.25"W x 8"D) \$225
- ☐ **MFLR** Right Large Footplate (6.25"W x 8"D) \$225
- ☐ **MFXL** Left Extra Large Footplate (9"W x 11"D) \$225
- ☐ **MFXR** Right Extra Large Footplate (9"W x 11"D) \$225

HEEL LOOPS, CUPS AND TOE STRAPS

- ☐ **HLL** Left Heel Loop (each) \$20
- ☐ **HLR** Right Heel Loop (each) \$20
- ☐ **HCL** Left Cup Loop (each) \$20
- ☐ **HCR** Right Cup Loop (each) \$20
- ☐ **TSL** Left Toe Strap (each) \$20
- ☐ **TSR** Right Toe Strap (each) \$20

FIXED ADUCTORS**Left Support**

- ☐ **HSPXLH** Size XS (5.5" W x 4" H) \$120
- ☐ **HSPSMLH** Size S/M (7" W x 4" H) \$120
- ☐ **HSPSLH** Size S (7"W x 6"H) \$120
- ☐ **HSPMLH** Size M (9" W x 4"H) \$120
- ☐ **HSPLLH** Size L (13" W x 4"H) \$120

Right Support

- ☐ **HSPXRH** Size XS (5.5"W x 4"H) \$120
- ☐ **HSPSMRH** Size S/M (7"W x 4"H) \$120
- ☐ **HSPSRH** Size S (7"W x 6"H) \$120
- ☐ **HSPMRH** Size M (9"W x 4"H) \$120
- ☐ **HSPLRH** Size L (13"W x 4"H) \$120

INSTALLATION

- ☒ **MCI** Motion Concepts Installation STD

Total Retail Value \$ _____