



Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

SAVE

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

SUBMIT

CLEAR FORM

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

888-433-6818

PRINT FORM

ULTRA LOW MAXX CG POWER POSITIONING SYSTEM ON TDX-SP2 HD BASE

US PRICE LIST AND ORDER FORM - Price Effective July 15, 2026



PRINT FORM

SAVE

CLEAR FORM

Motion Concepts Customer Service: 1-888-433-6818

quotes@motionconcepts.com

REQUIRED INFORMATION

☐ Quote ☐ Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag: _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

ATP/Therapist: _____

Ship to Address: _____

City: _____ State: _____


Zip: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all **REQUIRED*** measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness*	



POWER BASE OPTIONS

- ☐ TDXSP2HD-CG TDX-SP2 HD Base for Single Actuator System
HCPCS code K0858.....\$10,290
- ☐ TDXSP2HD-MCG TDX-SP2 HD Base for Multiple Actuator System
HCPCS code K0862.....\$10,290

BASE WIDTH

- ☐ STDFRM Standard Base 25.5" (MK 24 NF Batteries).....N/C

SEAT-TO-FLOOR HEIGHT⁽¹⁾

- ☐ SLOWH⁽²⁾ 17.25" Seat-To-Floor\$363
- ☐ LOWH 17.75" Seat-To-FloorN/C
- ☐ MEDH 18.75" Seat-To-FloorN/C
- ☐ TALLH 19.75" Seat-To-FloorN/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"

2. This STF needs a base modification to be performed by Motion Concepts Canada -MCI (see page 16)

USER WEIGHT LIMIT

- ☐ U350 Weight Capacity up to 350lbs..... STD
- ☐ U400 Weight Capacity up to 400lbs.....N/C

WHEELCHAIR TRANSPORT BRACKETS

- ☐ TRBKTS Wheelchair Transport Brackets⁽¹⁾..... STD
1. For unoccupied use only

TIRE OPTIONS

- ☐ B1431-3 14" x 3" Black Tire with Gel Foam Inserts..... STD
- ☐ B1430-3 14" x 3" Black Tire - Pneumatic Tire.....N/C

FORK OPTIONS

- ☐ DSFK Double-Sided Fork Package STD

BATTERY CHARGER OPTION

- ☐ 110CHARGER 110 Volt Battery Charger.....N/C
- ☐ L8900 Omit Battery Charger (\$84)

ENHANCED DRIVING

- ☐ GTRAC-LX LINX G-TRAC Module⁽¹⁾
HCPCS code K0108.....\$1,045

1. Recommended for use with Specialty Drive Controls

LIGHTS AND INDICATORS

- ☐ LIGHTS LED Lights and Indicators⁽¹⁾.....\$735

1. Must order REM215, REM216, REM400 or REM500

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FRAME FINISH

☐ 170P		Obsidian Black.....N/C
☐ 171P		Ruby Red.....N/C
☐ 172P		Peacock Blue.....N/C
☐ 174P		Purple Passion.....N/C
☐ 175P		Denim Navy.....N/C
☐ 176P		Lime Twist.....N/C
☐ 177P		Arctic Blue.....N/C

RIM INSERT COLORS - MUST PICK ONE

☐ 170PR		Obsidian Black.....N/C
☐ 171PR		Ruby Red.....N/C
☐ 172PR		Peacock Blue.....N/C
☐ 173PR		Pink Flamingo.....N/C
☐ 174PR		Purple Passion.....N/C
☐ 175PR		Denim Navy.....N/C
☐ 176PR		Lime Twist.....N/C
☐ 177PR		Arctic Blue.....N/C
☐ 178PR		Banana Split.....N/C
☐ 179PR		Orange Marmalade.....N/C
☐ 180PR		Morning Mist.....N/C

NON-EXPANDABLE SINGLE/DUAL JOYSTICKS

☐ REM110	LED Drive Only Remote/Joystick ⁽¹⁾N/C <i>Single drive profile</i>
☐ REM210	LED Non-Expandable Remote/Joystick ⁽¹⁾N/C <i>Three drive profiles, control of two seat functions</i>
☐ REM215	LED Non-Expandable Remote with Lights ⁽²⁾ ...\$368 <i>Three drive profiles, control of two seat functions and light operation</i>

1. Not available with lights and indicators.
2. Must choose lighting package separately

EXPANDABLE JOYSTICKS⁽¹⁾

☐ REM211	LED Expandable Remote/Joystick ⁽²⁾N/C <i>Three drive profiles, control of up to 6 actuators Not for use with specialty controls</i>
☐ REM216	LED Expandable Remote/Joystick with Lights ⁽³⁾\$368 <i>Three drive profiles, control of up to 6 actuators and light operation. Not for use with specialty controls</i>
☐ REM400	Color 3.5" Touch Screen Remote/Joystick ⁽⁴⁾\$1,050 <i>Four drive profiles, control of up to 6 actuators through the drive control and light operation.</i>
☐ SW400TGL	Toggle Kit for REM400 ⁽⁵⁾\$415

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Not available with lights and indicators
3. Must choose lighting package separately.
4. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.
5. Adds 2 toggle switches to REM400.

EXPANDABLE CONTROLLER AND HARNESS

☐ EXPC	Expandable Controller ⁽¹⁾ HCPCS code E2377.....\$775
☐ PWH	Harness for Expandable Controller HCPCS code E2313.....\$595

1. This item is required when ordering power tilt, power recline, and power elevating legrest together OR when ordering an Alternative Drive Control.

EXPANDABLE SPECIALTY DISPLAY

☐ REM500	LiNX Color 3.5" Touch Screen Display ⁽¹⁾⁽²⁾⁽³⁾\$1,575
---------------------------------	---

Select mounting position for REM500

☐ R5MR	Right.....N/C
☐ R5ML	Left.....N/C
☐ R5RPB	Rear Push Bar.....\$203

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Driver Control not included.
3. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.

SPECIALTY PROPORTIONAL CONTROLS⁽¹⁾⁽²⁾

☐ CREM	LiNX Compact Remote HCPCS code E2373...\$1,860
☐ CREM-LF	LiNX Compact Remote Low Force HCPCS code E2373.....\$1,860

1. Requires REM400 or REM500
2. Select Mount in Mounting For Driver Control section

ATTENDANT DRIVE CONTROL CHOICES

☐ ACU	Proportional Attendant Control ⁽¹⁾ HCPCS code E2331.....\$945
------------------------------	---

Select mounting position for ACU

☐ RPBMR	Right
☐ RPBML	Left
☐ RPBMC	Center

1. Requires Expandable Controller.

THIRD PARTY SPECIALTY CONTROLS

Motion Concepts Install recommended

Installed at Motion Concepts: ☐ YES ☐ NO

☐ Type: _____

☐ Specialty Control Mount:	☐ Right	☐ Left
☐ Display Mount:	☐ Right	☐ Left
☐ Mode Switch Mount:	☐ Right	☐ Left
☐ On/Off Mount:	☐ Right	☐ Left

Notes

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

JOYSTICK MOUNT

Arm Mounted

- ☐ FR Fixed Joystick Mounting - Right STD
- ☐ FL Fixed Joystick Mounting - Left STD
- ☐ MRR Maxx Resolve Swing-Away Joystick Mount
(Right Hand) HCPCS code E1032 \$573
- ☐ MRL Maxx Resolve Swing-Away Joystick Mount
(Left Hand) HCPCS code E1032 \$573
- ☐ SAQR Motion Height Adjustable Swing-Away
Quad Link (Right Hand)
HCPCS code E1032 \$341
- ☐ SAQL Motion Height Adjustable Swing-Away
Quad Link (Left Hand)
HCPCS code E1032 \$341

Rear Push Bar Mounted

- ☐ RPBMR Right \$203
- ☐ RPBML Left \$203
- ☐ RPBMC Center \$203

JOYSTICK KNOB - OPTIONAL

HCPCS code E2323

- ☐ PC101A BodyPoint 3" U-Shaped \$126
- ☐ PC102A BodyPoint 4" U-Shaped \$126
- ☐ PC104A BodyPoint Straight Handle \$126
- ☐ PC107A BodyPoint Domed Rubber \$126
- ☐ PC110A BodyPoint Grooved Mushroom \$126

1. To remove the joystick knob from the REM400, pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

ELECTRONICS - EXPANDABLE MODULES

- ☐ OUT500 Output Module \$525
- ☐ IN200 LiNX Switch Input Module with
Power ON/OFF for REM200 Series⁽¹⁾⁽²⁾ \$525
- ☐ IN500H LiNX 9 Pin Input Module for
Alternative Driver Controls⁽²⁾⁽³⁾⁽⁴⁾⁽⁶⁾ \$1,995
- ☐ IN500 LiNX 9 Pin Input Module with
Sip-N-Puff Nozzle⁽²⁾⁽³⁾⁽⁴⁾⁽⁶⁾ HCPCS code E2325. \$1,995
- ☐ PKG32666 Therafin Sip-N-Puff Breath Tube Kit \$509
HCPCS code E2326
- ☐ COVER-L Large Recline Shroud Cover⁽⁵⁾ N/C

1. Not for use with Alternative Drive Controls. (Uses switch ports).

2. Only one input module per chair.

3. Input module also required for any 9-pin specialty control.

4. 9 pin specialty control and Sip-N-Puff cannot both be used together.

5. May be required when choosing OUT500, IN200 or IN500 options.

6. Requires REM400 or REM500 for compatibility.

LINX ELECTRONIC ACCESSORIES

- ☐ LAK LINX Access Key⁽¹⁾ \$105

1. Required for programming

NOTES

Ultra Low Maxx CG Power Positioning System TDX-SP2 HD

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

POWER POSITIONING SYSTEMS⁽¹⁾

- ☐ MBT Ultra Low Maxx Bariatric CG Tilt (50°) Only System⁽¹⁾
Up to 350 lbs Capacity \$8,178
- ☐ MBT400 Ultra Low Maxx Bariatric CG Tilt (50°) Only System⁽¹⁾
350 lbs - 400 lbs Capacity \$9,108
- ☐ MBTRE Ultra Low Maxx Bariatric CG Tilt (50°)
and Power Recline (150°)
w/Extended Shear Reduction(ESR)⁽¹⁾
Up to 350 lbs Capacity \$12,630
- ☐ MBTRE400 Ultra Low Maxx Bariatric CG Tilt (50°)
and Power Recline (150°)
w/Extended Shear Reduction(ESR)⁽¹⁾
350 lbs to 400 lbs Capacity \$13,508
- ☐ MBRE Ultra Low Maxx Bariatric Power Recline (150°)
w/Extended Shear Reduction(ESR)
Up to 350 lbs Capacity \$6,894
- ☐ MBRE400 Ultra Low Maxx Bariatric Power Recline (150°)
w/Extended Shear Reduction(ESR)
350 lbs to 400 lbs Capacity \$7,757
- ☐ MBE Ultra Low Maxx Bariatric 10" Elevating
Lift Seat Only⁽²⁾ HCPCS code E2298 \$8,693
Up to 350 lbs Capacity

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module
2. Not available on seat sizes greater than 22" wide by 22" deep.

ADD: POWER ADJUSTABLE SEAT HEIGHT MODULE

- ☐ ESM10 10" HD Power Adjustable Seat Height Module⁽¹⁾⁽²⁾
HCPCS code E2298 \$4,193
Weight capacity up to 350lbs

1. Only available with tilt, uses an integrated Tilt/Lift mechanism, tilt is 45°.
2. Not available on seat sizes greater than 22" wide by 22" deep.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MODULES - MUST SELECT SYSTEM TO ADD MODULES

- ☐ PRM8 8° Precline Module⁽¹⁾ N/C
☐ PRM30 30° Precline Module⁽¹⁾ N/C
HCPCS code K0108 \$715
- ☐ ERSSP Programable Electronic Recline Function
Start & Stop⁽²⁾ \$580
- Specify:
Start: _____ degrees Stop: _____ degrees

1. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available. Must order Recline System
2. Requires Multi-Function M610 through Drive Control (MFCBD). Set range between 90° and 170°. Will be set to factory settings unless specified.

ADULT SEAT WIDTH SETTINGS

- ☐ BW120 20" Wide (Adjustable from 20" to 24") N/C
☐ BW121 21" Wide (Adjustable from 20" to 24") N/C
☐ BW122 22" Wide (Adjustable from 20" to 24") N/C
☐ BW123 23" Wide (Adjustable from 20" to 24") N/C
☐ BW124 24" Wide (Adjustable from 20" to 24") N/C
☐ BW222 22" Wide (Adjustable from 22" to 26") N/C
☐ BW223 23" Wide (Adjustable from 22" to 26") N/C
☐ BW224 24" Wide (Adjustable from 22" to 26") N/C
☐ BW225 25" Wide (Adjustable from 22" to 26") N/C
☐ BW226 26" Wide (Adjustable from 22" to 26") N/C

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

- ☐ SD117 17" Deep (Adjustable from 17" to 21") N/C
☐ SD118 18" Deep (Adjustable from 17" to 21") N/C
☐ SD119 19" Deep (Adjustable from 17" to 21") N/C
☐ SD120 20" Deep (Adjustable from 17" to 21") N/C
☐ SD121 21" Deep (Adjustable from 17" to 21") N/C
☐ SD219 19" Deep (Adjustable from 19" to 23") N/C
☐ SD220 20" Deep (Adjustable from 19" to 23") N/C
☐ SD221 21" Deep (Adjustable from 19" to 23") N/C
☐ SD222 22" Deep (Adjustable from 19" to 23") N/C
☐ SD223 23" Deep (Adjustable from 19" to 23") N/C

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX SEAT CUSHION OPTIONS

Matrix PS Seat Cushion 20" wide/reversable cover

HCPCS code E2605..... \$383

☐ PS2016 ☐ PS2018 ☐ PS2020

Matrix PS Seat Cushion 22" - 26" wide

HCPCS code E2606..... \$605

☐ PS2218 ☐ PS2220 ☐ PS2222 ☐ PS2418 ☐ PS2420
☐ PS2422 ☐ PS2618 ☐ PS2620 ☐ PS2622

Matrix PSP Seat Cushion 20" wide/reversable cover

HCPCS code E2607..... \$438

☐ PSP2016 ☐ PSP2018 ☐ PSP2020

Matrix PSP Seat Cushion 21"-22" wide

HCPCS code E2608..... \$660

☐ PSP2218 ☐ PSP2220 ☐ PSP2222

Matrix Vi Seat Cushion 20" wide

HCPCS code E2607..... \$497

☐ MA2016-VI ☐ MA2018-VI ☐ MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide

HCPCS code E2608..... \$556

☐ MA2218-VI ☐ MA2220-VI

Matrix Vi HD Seat Cushion 22" -24" wide

HCPCS code E2608..... \$731

☐ MA2222-VI ☐ MA2418-VI ☐ MA2420-VI

Matrix Vi HD Seat Cushion up to 26" wide

HCPCS code E2608..... \$731

☐ MWWDD-VI

Specify Size: _____

Matrix Libra Seat Cushion 20" wide

HCPCS code E2624..... \$567

☐ LC2016 ☐ LC2018 ☐ LC2020

Matrix Libra HD Seat Cushion 21" wide

HCPCS code E2624..... \$614

☐ LC2118 ☐ LC2120

Matrix Libra Seat Cushion 22" deep

HCPCS code E2624..... \$696

☐ LC2022

Matrix Libra HD Seat Cushion 22" wide

HCPCS code E2625..... \$696

☐ LC2122 ☐ LC2218 ☐ LC2220 ☐ LC2222 ☐ LC2318
☐ LC2320 ☐ LC2322 ☐ LC2418 ☐ LC2420 ☐ LC2422

Matrix Libra Adjust Seat Cushion 20" wide

HCPCS code E2624..... \$731

☐ LC2016-ADJ ☐ LC2018-ADJ ☐ LC2020-ADJ

Matrix Libra Adjust Seat Cushion 21" wide

HCPCS code E2624..... \$778

☐ LC2118-ADJ ☐ LC2120-ADJ

Matrix Libra Adjust Seat Cushion 22" deep

HCPCS code E2624..... \$860

☐ LC2022-ADJ

Matrix Libra Adjust Seat Cushion 21" - 24" wide

HCPCS code E2624 for 21" W / E2625 for 22 - 24" W \$860

☐ LC2122-ADJ ☐ LC2218-ADJ ☐ LC2220-ADJ ☐ LC2222-ADJ
☐ LC2318-ADJ ☐ LC2320-ADJ ☐ LC2322-ADJ ☐ LC2418-ADJ
☐ LC2420-ADJ ☐ LC2422-ADJ

Matrix Multi Cushion 20" wide w/reversable cover

HCPCS code E2607..... \$521

☐ MM2016 ☐ MM2018 ☐ MM2020

Matrix Multi Cushion 22" wide w/reversable cover

HCPCS code E2608..... \$696

☐ MM2216 ☐ MM2218 ☐ MM2220 ☐ MM2222

☐ OCCC "ON-CHAIR" Cushion Credit⁽¹⁾ (\$105)

1. Not available on Matrix PS or PSP Cushions.

FABRIC OPTIONS FOR CUSHIONS

Cushion	Meshtex	Startex	Spacetex	Onyx	Fusion
Libra/ Adjust	N/A	STD	N/A	N/A	N/A
PS/PSP	N/A	STD Reversible	N/A	N/A	N/A
Vi	STD	\$59	N/A	\$117	\$176

- ☐ LC-SP Spacetex Fabric Option for
Libra Cushion N/C
- ☐ LC-IC Infection Control Fabric Upgrade for
Libra Cushion \$59
- ☐ FSU-MA Fusion Fabric Upgrade for Matrix Vi Cushion... \$176
- ☐ OCU-MA Onyx Fabric Upgrade for Matrix Vi Cushion \$117
- ☐ STARTEX-S Smooth side Startex Fabric Upgrade for
Matrix Vi Cushion \$59
- ☐ STARTEX-F Fabric side Startex Fabric Upgrade for
Matrix Vi Cushion \$59

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

LAP BELTS

☐ BELT48	Push Button Style Lap Belt, Length 48" ⁽¹⁾ STD
☐ BELT60	Push Button Style Lap Belt, Length 60" ⁽¹⁾ STD
☐ BELT71	Push Button Style Lap Belt, Length 71" ⁽¹⁾ STD
☐ PLB2	Bodypoint Padded Lap Belt 2 Point ⁽¹⁾
	HCPCS code E0978. \$222
☐ PLB4	Bodypoint Padded Lap Belt 4 Point ⁽¹⁾
	HCPCS code E0978. \$275

1. 67" length only option

MATRIX HEADREST PADS⁽¹⁾

HCPCS code E0955

Standard Pads

☐ ESP6	Matrx Elan Standard Pad 6" \$275
☐ ESP10	Matrx Elan Standard Pad 10" \$275
☐ ESP14	Matrx Elan Standard Pad 14" \$275
☐ ESP6-HUG	Matrx Elan Hug Standard Pad 6" \$275
☐ ESP10-HUG	Matrx Elan Hug Standard Pad 10" \$275
☐ ESP14-HUG	Matrx Elan Hug Standard Pad 14" \$275
☐ SNUG14	Matrx SNUG Standard Pad 14" \$486
☐ SNUG14-A	Matrx SNUG Asymmetrical Pad 14" \$521
☐ SNUG18	Matrx SNUG Standard Pad 18" \$486
☐ SNUG18-A	Matrx SNUG Asymmetrical Pad 18" \$521

Standard Pads with Infection Control Cover

☐ ESP6-IC	Matrx Elan Infection Control Pad 6". \$275
☐ ESP10-IC	Matrx Elan Infection Control Pad 10". \$275
☐ ESP14-IC	Matrx Elan Infection Control Pad 14". \$275
☐ ESP6-IC-HUG	Matrx Elan Hug Infection Control Pad 6" \$275
☐ ESP10-IC-HUG	Matrx Elan Hug Infection Control Pad 10". \$275
☐ ESP14-IC-HUG	Matrx Elan Hug Infection Control Pad 14". \$275
☐ SNUG14-IC	Matrx SNUG Standard
	Infection Control Pad 14". \$486
☐ SNUG14-A-IC	Matrx SNUG Asymmetrical
	Infection Control Pad 14". \$521
☐ SNUG18-IC	Matrx SNUG Standard
	Infection Control Pad 18". \$486
☐ SNUG18-A-IC	Matrx SNUG Asymmetrical
	Infection Control Pad 18". \$521

Occipital /4-Point Pads

☐ EOP9	Matrx Elan Occipital Pad 9" \$275
☐ EOP12	Matrx Elan Occipital Pad 12" \$275
☐ E4PTSM	Matrx Elan 4-Point Pad 9"W x 7"H. \$380
☐ E4POINT	Matrx Elan 4-Point Pad 11"W x 10"H. \$380

Occipital / 4-Point Pads with Infection Control Cover

☐ EOP9-IC	Matrx Elan Occipital Pad 9" with Infection Control. \$275
☐ EOP12-IC	Matrx Elan Occipital Pad 12" with Infection Control. \$275
☐ E4PTSM-IC	Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control. \$380
☐ E4POINT-IC	Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control \$380

1. If mounting on PB Elite must order removable hardware

MOTION CONCEPTS HEADREST PAD

HCPCS code E0955

☐ MCSH	Motion Concepts Standard \$428
☐ MCOH	Motion Concepts Auto Style ⁽¹⁾ \$508

1. Comes Standard with Fixed Mounting hardware; for use on Matrx PB or Elite Backs you must select the removable hardware upgrade below.

HEADREST MOUNTING HARDWARE

HCPCS code E1033


☐ FM	Fixed Mounting ⁽¹⁾ N/C
☐ LXSHW	Loxx Standard Hardware 15" \$252
☐ LXMMH	Loxx Mini Hardware 11" \$252
☐ MEHW	Elan Hardware 15" \$263
☐ MEMH	Elan Mini Hardware 11" \$263
☐ MEHW-EXT	Elan Extended Hardware \$357
☐ FPLXSHW	FLIPP Loxx Standard Hardware 15" \$603
☐ FPLXSHW-EX	FLIPP Loxx Extended Hardware \$614
☐ FPLXMMH	FLIPP Loxx Mini Hardware 11" \$603
☐ FPLXMMH-EX	FLIPP Loxx Mini Extended Hardware \$614
☐ FPMEHW	FLIPP Elan Hardware 15" \$614
☐ FPMEHW	FLIPP Elan Mini Hardware 11" \$614
☐ FPMEHW-EXT	FLIPP Elan Extended Hardware \$708

☐ OHR	OMIT: Headrest N/C
------------------------------	------------------------------

1. HCPCS code does not apply.

TILT ONLY BACK CANE/SPREADER BAR OPTION⁽¹⁾

☐ MSBSB	Maxx Style Back Adjustable Spreader Bar STD
	Adjust height easily on back canes
☐ CSBSB	Classic Style Back Adjustable Spreader Bar(2) . . \$145
☐ OMR	Optional Multi-Rail. \$218
	Standard Rehab Back with Classic Style Spreader Bar Only.



1. Not available on the Single Post Backrest

2. Not available with Outback Cantilever Arms or the Freedom Interface Hardware

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

STANDARD BACK CANES - FOR TILT SYSTEMS

- | | | |
|--------------------------------|-------------------------------|-----|
| ☐ SBC18 | 18" Straight Back Canes | STD |
| ☐ SBC20 | 20" Straight Back Canes | STD |
| ☐ SBC22 | 22" Straight Back Canes | STD |
| ☐ SBC24 | 24" Straight Back Canes | STD |

Select Cane Angle

- | | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-----------------------------------|
| ☐ BA75 75° | ☐ BA79 79° | ☐ BA85 85° | ☐ BA90 90° | ☐ BA95 95° |
| ☐ BA101 101° | ☐ BA106 106° | ☐ BA112 112° | ☐ BA116 116° | |

- | | | |
|--------------------------------|----------------------------------|-----|
| ☐ ABC20 | 20" Angled (5°) Back Canes | N/C |
| ☐ ABC22 | 22" Angled (5°) Back Canes | N/C |
| ☐ ABC24 | 24" Angled (5°) Back Canes | N/C |

Select Cane Angle

- | | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ BA81 81° | ☐ BA84 84° | ☐ BA91 91° | ☐ BA96 96° | ☐ BA101 101° |
| ☐ BA106 106° | ☐ BA111 111° | ☐ BA117 117° | ☐ BA121 121° | |

- | | | |
|--------------------------------|---|-------|
| ☐ EBC22 | 22" Endomorph Back Canes ⁽¹⁾ | \$508 |
| ☐ EBC24 | 24" Endomorph Back Canes ⁽¹⁾ | \$508 |
| Please Select | | |
| ☐ EBCN1 | 1" Narrower | N/C |
| ☐ EBCN2 | 2" Narrower | N/C |
| ☐ EBCN4 | 4" Narrower | N/C |

- | | | |
|--------------------------------|--|-------|
| ☐ TBC18 | 18" 10° Mid-Angle Tilt Back Canes ⁽²⁾ | \$363 |
| ☐ TBC20 | 20" 10° Mid-Angle Tilt Back Canes ⁽²⁾ | \$363 |
| ☐ TBC22 | 22" 10° Mid-Angle Tilt Back Canes ⁽²⁾ | \$363 |
| ☐ TBC24 | 24" 10° Mid-Angle Tilt Back Canes ⁽²⁾ | \$363 |

1. Comes standard with Cantilever Armrest Cane Mounted.
2. Quick Height Adjustment on Backrest and Maxx Spreader Bars not available.
3. Minimum finished back widths of 16".
4. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
5. 4" of total vertical mounting adjustment.

CANE OPTIONS - FOR RECLINE SYSTEMS

- | | | |
|-------------------------------|--|---------|
| ☐ BFS2 | Recline Body Fit 2" Narrower ⁽¹⁾ | \$1,443 |
| ☐ BFS4 | Recline Body Fit 4" Narrower ⁽¹⁾ | \$1,443 |
| ☐ AHH | Adjustable Hooking Handles ⁽²⁾ (pair) | \$653 |
| ☐ RPHS | Push Handles for Recline - Straight ⁽³⁾ (pair) | \$508 |
| ☐ RPHA | Push Handles for Recline - Angled ⁽³⁾ (pair) | \$508 |

1. Minimum finished back widths of 16".
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ - FOR TILT ONLY/ELEVATE ONLY

- | | | | |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| ☐ FBT18 18" | ☐ FBT19 19" | ☐ FBT20 20" | ☐ FBT21 21" |
| ☐ FBT22 22" | ☐ FBT23 23" | ☐ FBT24 24" | ☐ FBT25 25" |

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ - FOR RECLINE OR TILT/RECLINE

- | | | | |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| ☐ FBR20 20" | ☐ FBR21 21" | ☐ FBR22 22" | ☐ FBR23 23" |
| ☐ FBR24 24" | ☐ FBR25 25" | ☐ FBR26 26" | ☐ FBR27 27" |

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

STANDARD BACK OPTIONS

- | | | |
|-------------------------------|--|---------|
| ☐ SRB | Standard Rehab Back & Cushion | N/C |
| ☐ SS | Additional Super Soft/HR Foam ⁽¹⁾ | \$653 |
| ☐ MFU | Matrx Fabric Upcharge | \$145 |
| ☐ OBC | OMIT: Back Cushion | (\$36) |
| ☐ PIPF | Flat Planar Interface/I-Back Plate ⁽²⁾⁽³⁾ | \$181 |
| ☐ PIPR | Recessed Planar/I-Back Interface Plate ⁽²⁾⁽³⁾⁽⁵⁾ | \$181 |
| ☐ ORBC | OMIT: Rigid Back & Cushion ⁽⁴⁾ | (\$109) |

1. Additional foam added may impact seat depth
2. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
3. Required with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
4. Not available with Recline/ESR or Recline/PSB.
5. Fits Comfort Company, Freedom, Pindot and Moulded Backs.

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

UPGRADABLE BACK OPTIONS

High Back 20" - 24" W \$696

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ HB2016 | ☐ HB2018 | ☐ HB2020 | ☐ HB2116 |
| ☐ HB2118 | ☐ HB2120 | ☐ HB2216 | ☐ HB2218 |
| ☐ HB2220 | ☐ HB2316 | ☐ HB2318 | ☐ HB2320 |
| ☐ HB2416 | ☐ HB2418 | ☐ HB2420 | |

Matrx Elite Back 20" W (3" Contour) \$731

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBE2014 | ☐ PBE2016 | ☐ PBE2018 | ☐ PBE2020 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrx Elite HD Back 21" - 24" W (3" Contour) \$930

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE2116-HD | ☐ PBE2120-HD | ☐ PBE2216-HD | ☐ PBE2220-HD |
| ☐ PBE2316-HD | ☐ PBE2320-HD | ☐ PBE2416-HD | ☐ PBE2420-HD |

Matrx Elite HD Back 26" W (3" Contour) \$1,632

- | | |
|-------------------------------------|-------------------------------------|
| ☐ PBE2616-HD | ☐ PBE2620-HD |
|-------------------------------------|-------------------------------------|

Matrx Elite TR Back 20" W (3" Contour) \$ 761

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE2016-TR | ☐ PBE2018-TR | ☐ PBE2020-TR |
|-------------------------------------|-------------------------------------|-------------------------------------|

Matrx Elite TR HD Back 21" - 24" W (3" Contour) \$989

- | | | | |
|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| ☐ PBE2120-TRHD | ☐ PBE2220-TRHD | ☐ PBE2320-TRHD | ☐ PBE2420-TRHD |
|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|

Matrx Elite Deep Back 14" - 20" W (6" Contour) \$889

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ EDB2014 | ☐ EDB2016 | ☐ EDB2018 | ☐ EDB2020 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrx Elite Deep HD Back 21" - 24" W (6" Contour) \$1,088

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ EDB2116-HD | ☐ EDB2118-HD | ☐ EDB2120-HD | ☐ EDB2216-HD |
| ☐ EDB2218-HD | ☐ EDB2220-HD | ☐ EDB2316-HD | ☐ EDB2318-HD |
| ☐ EDB2320-HD | ☐ EDB2416-HD | ☐ EDB2418-HD | ☐ EDB2420-HD |

Matrx PB Posture Back 20" W (5" Contour) \$731

- | | | |
|----------------------------------|----------------------------------|----------------------------------|
| ☐ MPB2012 | ☐ MPB2016 | ☐ MPB2020 |
|----------------------------------|----------------------------------|----------------------------------|

Matrx PB Heavy Duty Back⁽¹⁾ 20" - 22" W (5" Contour) \$977

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ HDB2016 | ☐ HDB2020 | ☐ HDB2216 | ☐ HDB2220 |
| ☐ HDB2416 | ☐ HDB2420 | | |

Matrx PB Posture Back Deep⁽¹⁾ 20" W (7" Contour) \$755

- | | | |
|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBD2012 | ☐ PBD2016 | ☐ PBD2020 |
|----------------------------------|----------------------------------|----------------------------------|

Matrx PB Heavy Duty Deep Back⁽¹⁾
20"-22" W (5" Contour) \$1,094

- | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ DHDB2016 | ☐ DHDB2020 | ☐ DHDB2216 | ☐ DHDB2220 |
| ☐ DHDB2416 | ☐ DHDB2420 | | |

☐ OCBC "ON-CHAIR" Matrx Back Credit (\$146)

1. The Heavy Duty PB Back and PB Deep Back are not available on Recline systems with seat widths of 21" and 22"

OTHER UPGRADABLE BACK OPTIONS

Cushion	Meshtex	4-Way Stretch Polyester	Onyx	Fusion
High	STD	N/A	\$117	\$145
E2, Elite, Elite HD, Elite Deep	STD	N/A	N/A	N/A
E2, Elite TR, Elite TR HD	STD	N/A	N/A	N/A

- | | | |
|--------------------------------|--|-------|
| ☐ EBSS | Super Soft Foam for Elite Back Upgrade..... | \$117 |
| ☐ FFUHB | Fusion Fabric Upgrade for High Backs..... | \$145 |
| ☐ OCUHB | Onyx Fabric Upgrade for High Backs..... | \$117 |
| ☐ CHIEB | Chest Harness Interface Kit for Elite Back ⁽¹⁾⁽²⁾ | \$111 |
| ☐ CHIPB | Chest Harness Interface Kit for PB Back ⁽¹⁾⁽²⁾ | \$53 |

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harness

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

TILT ARMRESTS

Adult Dual Post

☐ DP1	9.75" - 12.75"	STD
☐ DP2	12.75" - 15.75"	STD
☐ DP3	15.75" - 18.75"	STD
☐ FP	Flat Pouch	STD
☐ GB	Glove Boxes	N/C
☐ OAL	Outback Arm Option Left	\$283
☐ OAR	Outback Arm Option Right	\$283

Ultra Rail Seat Mount Flip Back Cantilever MAXX Tilt Arm⁽¹⁾

☐ RMCAN	9" - 12.5"	\$363
☐ RMCAN2	11" - 14"	\$363
☐ DPM	Dual Post Module to RMCAN (pair)	\$363

Outback Cantilever Tilt Armrest

☐ OCTA1	9"-12"	\$508
☐ OCTA2	12"-15"	\$508
☐ OCTA3	6"-12"	\$508

Cantilver Arms for Endomorph Back Canes⁽¹⁾

☐ CAFEB	See page 4	STD
--------------------------------	------------	-----

RECLINE ARMRESTS

Reclining Adjustable Height

☐ RC1	9.5"-13"	STD
☐ RC2	13"-16"	STD
☐ OAL	Outback Arm Option Left	\$283
☐ OAR	Outback Arm Option Right	\$283

Maxx Style Cane Mounted Height, Width, Angle Adj.

Cantilever Flip Back Armrests⁽¹⁾

☐ MCAN1	9.5"-13"	\$363
☐ MCAN2	12.5"-16"	\$363

Outback Cantilever Recline Armrest

☐ OCRA1	9"-12"	\$508
☐ OCRA2	12"-15"	\$508
☐ OCRA3	7" - 9"	\$508

Cane Mounted Height and Angle Adj.

Cantilever Flip Back Locking Armrest⁽¹⁾

☐ CANQ1	9.5"-12"	\$508
☐ CANQ2	10"-13"	\$508
☐ CANQ3	13"-16"	\$508

Adult Dual Post On Recline

☐ DPOR1	9.75" - 12.75"	\$363
☐ DPOR2	12.75" - 15.75"	\$363
☐ DPOR3	15.75" - 18.75"	\$363
☐ FP	Flat Pouch	STD
☐ GB	Glove Boxes	N/C

ARMREST TUBE LENGTH - PLEASE SPECIFY

☐ ATFL	Left Full	STD
☐ ATDL	Left Desk	STD
☐ ATFR	Right Full	STD
☐ ATDR	Right Desk	STD

ARMPADS⁽⁵⁾

Standard

☐ SAPFL	Left Full (2.25" x 14")	STD
☐ SAPDL	Left Desk (2.25" x 10")	STD
☐ SAPML	Left Mini (2.25" x 7")	STD
☐ SAPFR	Right Full (2.25" x 14")	STD
☐ SAPDR	Right Desk (2.25" x 10")	STD
☐ SAPMR	Right Mini (2.25" x 7")	STD

Waterfall Arm Pads⁽¹⁾

☐ WAPFL	Left Full (13")	\$73
☐ WAPDL	Left Desk (9")	\$73
☐ WAPFR	Right Full (13")	\$73
☐ WAPDR	Right Desk (9")	\$73

Modular Arm Pad - Left⁽²⁾

☐ MAPTFL	Left Tray Full (14")	\$80
☐ MAPTDL	Left Tray Desk (10")	\$80

Select Pad Insert below:

☐ ISSFL	Iskin - Integrated soft skin foam
☐ SCVFL	Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

☐ MAPTFR	Right Tray Full (14")	\$80
☐ MAPTDR	Right Tray Desk (10")	\$80

Select Pad Insert below:

☐ ISSFR	Iskin - Integrated soft skin foam
☐ SCVFR	Startex Covered Visco foam

Flat Multi Position Arm Pads

☐ MPAFL	Left Full (14")	\$64
☐ MPADL	Left Desk (10")	\$64
☐ MPAML	Left Mini (7")	\$64
☐ MPAFR	Right Full (14")	\$64
☐ MPADR	Right Desk (10")	\$64
☐ MPAMR	Right Mini (7")	\$64

Gel Pads

☐ GAPNL	Narrow Left Desk (2" x 12")	\$88
☐ GAPNR	Narrow Right Desk (2" x 12")	\$88
☐ GAPWL	Wide Left Desk (3.5" x 12")	\$88
☐ GAPWR	Wide Right Desk (3.5" x 12")	\$88
☐ GAPLL	Long Left Full (3.5" x 14")	\$105
☐ GAPLR	Long Right Full (3.5" x 14")	\$105

Arm Troughs

☐ EATL	Ergonomic Arm Trough - Left HCPCS code E2209	\$99
☐ EATR	Ergonomic Arm Trough - Right HCPCS code E2209	\$99
☐ FHPL	Flat Hand Pad - Left ⁽³⁾ HCPCS code K0108	\$64
☐ FHPR	Flat Hand Pad - Right ⁽³⁾ HCPCS code K0108	\$64

Arm Pad Storage/Cup Holder⁽⁴⁾

☐ APSL	Left	\$428
☐ APSR	Right	\$428
☐ OAP	Omit Armpads(7/8" diameter tube)	(\$36)

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side, or with Outback Cantilever arms.
5. Full length arm pads used on the same side as the Joystick result in a "Back to Joystick" measurement of 20" minimum (24" including switch mounting)

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ARM ADDUCTOR PAD⁽¹⁾⁽²⁾

☐ AADSL	Left Small (6" x 2.25")	\$145
☐ AADSR	Right Small (6" x 2.25")	\$145
☐ AADLL	Left Large (8.5" x 2.25")	\$145
☐ AADLR	Right Large (8.5" x 2.25")	\$145

1. Only available with Modular or Flat Multi Position Arm Pads.
2. 1" thick pad, unlimited fore/aft adjustment.

ELBOW BLOCKS

Left Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

☐ EBXL	Size XS (3.25" x 4.25")	\$283
☐ EBSL	Size S (3.5" x 5.25")	\$283
☐ EBML	Size M (3.75" x 5.75")	\$283
☐ EBL	Size L (4.25" x 6.5")	\$283

☐ LV3	Left Vertical Extension 3"
☐ LV5	Left Vertical Extension 5"
☐ LH5	Left Horizontal Extension 5"
☐ LH7	Left Horizontal Extension 7"

Right Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

☐ EBXR	Size XS (3.25" x 4.25")	\$283
☐ EBSR	Size S (3.5" x 5.25")	\$283
☐ EBMR	Size M (3.75" x 5.75")	\$283
☐ EBLR	Size L (4.25" x 6.5")	\$283

☐ RV3	Right Vertical Extension 3"
☐ RV5	Right Vertical Extension 5"
☐ RH5	Right Horizontal Extension 5"
☐ RH7	Right Horizontal Extension 7"

1. Prevents arm from slipping.

ADDITIONAL ARM ACCESSORIES

☐ MACESL	Multi Axis Upper Extremity Support ⁽¹⁾ (each)	
	Left Side	\$263
☐ MACESR	Multi Axis Upper Extremity Support ⁽¹⁾ (each)	
	Right Side	\$263
☐ OAS2L	1" Offset Armrest Spacer Left(each)	\$181
☐ OAS2R	1" Offset Armrest Spacer Right(each)	\$181
☐ OAS1L	1/2" Offset Armrest Spacer Left(each)	\$181
☐ OAS1R	1/2" Offset Armrest Spacer Right(each)	\$181
☐ MPAPA	Adapter Block ⁽²⁾ (each) Quantity ()	\$51

1. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
2. Not available with Standard Arm Pads or Waterfall Arm Pads

MAXX LATERAL TRUNK SUPPORTS⁽¹⁾

Lateral Trunk Support Left HCPCS code E0956

☐ LATXL	Size XS (3.25" x 4.25")	\$145
☐ LATSL	Size S (3.5" x 5.25")	\$145
☐ LATML	Size M (3.75" x 5.75")	\$145
☐ LATLL	Size L (4.25" x 6.5")	\$145
Add:		
☐ SALATL	MAXX Style Swing Away Hardware - Left ⁽²⁾	
	HCPCS code E1034	\$1283
☐ ATLL	Additional Telescoping Link - Left ⁽³⁾	\$109

Lateral Trunk Support Right HCPCS code E0956

☐ LATXR	Size XS (3.25" x 4.25")	\$145
☐ LATSR	Size S (3.5" x 5.25")	\$145
☐ LATMR	Size M (3.75" x 5.75")	\$145
☐ LATLR	Size L (4.25" x 6.5")	\$145
Add:		
☐ SALATR	MAXX Style Swing Away Hardware - Right ⁽²⁾	
	HCPCS code E1034	\$283
☐ ATLR	Additional Telescoping Link - Right ⁽³⁾	\$109

1. Comes standard with Fixed Mount Hardware. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Comes standard with 1" telescoping link which gives 0-2.75" medial offset.
3. The amount of offset available may be impacted by back thickness.

MATRIX FIXED LATERAL TRUNK SUPPORTS⁽¹⁾

HCPCS code E1034

☐ PBSSL	Matrix Small Lateral Fixed-Left (3x5.5")	\$117
☐ PBSSR	Matrix Small Lateral Fixed-Right (3x5.5")	\$117
☐ PBLSL	Matrix Standard Lateral Fixed-Left (7.5x5")	\$117
☐ PBLSR	Matrix Standard Lateral Fixed-Right (7.5x5")	\$117
☐ PBOLL	Matrix Offset Lateral Fixed Left (7.5x5")	\$129
☐ PBOLR	Matrix Offset Lateral Fixed Right (7.5x5")	\$129

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

MATRIX SWING AWAY LATERAL TRUNK SUPPORTS⁽¹⁾

Matrix Lateral Trunk Support Left HCPCS code E0956

☐ EPADXL	Size XS (3.25" x 4.25")	\$117
☐ EPADSL	Size S (3.5" x 5.25")	\$117
☐ EPADML	Size M (3.75" x 5.75")	\$117
☐ EPADLL	Size L (4.25" x 6.5")	\$117
Must select one option:		
☐ SAHEBL	Matrix Elite Swing Away Hardware Left	
	HCPCS code E1034	\$228
☐ OSAHL	Matrix Offset Elite Swing Away Hardware Left	
	HCPCS code E1034	\$228

Matrix Lateral Trunk Support Right HCPCS code E0956

☐ EPADXR	Size XS (3.25" x 4.25")	\$117
☐ EPADSR	Size S (3.5" x 5.25")	\$117
☐ EPADMR	Size M (3.75" x 5.75")	\$117
☐ EPADLR	Size L (4.25" x 6.5")	\$117
Must select one option:		
☐ SAHEBR	Matrix Elite Swing Away Hardware Right	
	HCPCS code E1034	\$228
☐ OSAHR	Matrix Offset Elite Swing Away Hardware Right	
	HCPCS code E1034	\$228

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HIP AND KNEE SUPPORTS

Lateral Hip Support Left HCPCS code E0953

☐ HSPXLH	Size XS (5.5" W x 4" H)	\$174
☐ HSPSMLH	Size S/M (7" W x 4" H)	\$174
☐ HSPSLH	Size S (7" W x 6" H)	\$174
☐ HSPMLH	Size M (9" W x 4" H)	\$174
☐ HSPLLH	Size L (13" W x 4" H)	\$174

Lateral Hip Support Right⁽¹⁾ HCPCS code E0953

☐ HSPXRH	Size XS (5.5" W x 4" H)	\$174
☐ HSPSMRH	Size S/M (7" W x 4" H)	\$174
☐ HSPSRH	Size S (7" W x 6" H)	\$174
☐ HSPMRH	Size M (9" W x 4" H)	\$174
☐ HSPLRH	Size L (13" W x 4" H)	\$174

Lateral Knee Support Left⁽¹⁾ HCPCS code E0953

☐ HSPXLK	Size XS (5.5" W x 4" H)	\$174
☐ HSPSMLK	Size S/M (7" W x 4" H)	\$174
☐ HSPSLK	Size S (7" W x 6" H)	\$174
☐ HSPMLK	Size M (9" W x 4" H)	\$174
☐ HSPLLK	Size L (13" W x 4" H)	\$174

Lateral Knee Support Right⁽¹⁾ HCPCS code E0953

☐ HSPXRK	Size XS (5.5" W x 4" H)	\$174
☐ HSPSMRK	Size S/M (7" W x 4" H)	\$174
☐ HSPSRK	Size S (7" W x 6" H)	\$174
☐ HSPMRK	Size M (9" W x 4" H)	\$174
☐ HSPLRK	Size L (13" W x 4" H)	\$174

Add Removable Hardware:

MAXX Style Quick Release, Multi-Axis Mounting

Hardware for Hip Supports HCPCS code E1034

☐ MHSHLH	Left mounted at Hip	\$312
☐ MHSHLK	Left mounted at Knee	\$312
☐ MHSHRH	Right mounted at Hip	\$312
☐ MHSHRK	Right mounted at Knee	\$312

Swing Away Hardware for Hip Supports HCPCS code E1034

☐ SHSHLH	Left mounted at Hip	\$312
☐ SHSHLK	Left mounted at Knee	\$312
☐ SHSHRH	Right mounted at Hip	\$312
☐ SHSHRK	Right mounted at Knee	\$312

Extended Maxx Style Quick Release, Multi-Axis Mounting

Hardware for Hip Supports⁽²⁾ HCPCS code E1034

☐ EMQHLH	Left mounted at Hip	\$312
☐ EMQHLK	Left mounted at Knee	\$312
☐ EMQHRH	Right mounted at Hip	\$312
☐ EMQHRK	Right mounted at Knee	\$312

Lift Off Hardware for Hip Support

(Hip Position Only) HCPCS code E1034

☐ LOHL	Left	\$283
☐ LOHR	Right	\$283

1. Comes standard with Fixed Mount Hardware.

2. Recommended when finished seat cushion thickness will be greater than 3.5" high.

POWER CENTER MOUNT FOOT PLATFORM

LNK Power Center Mount Foot Platform HCPCS code E1012

☐ LNX	9" Articulation	\$3,594
☐ LSPTF1	10" - 14" Seat pan to footplate	
☐ LSPTF2	13" - 19" Seat pan to footplate	
☐ LNXA	6" Articulation ⁽³⁾	\$3,594
☒ LASPTF2	10" - 14" Seat pan to footplate	

Options:

☐ LIFP	Rubber Coated Footplate 11.5" W x 10" D	STD
☐ LICP	Individual calf pads	STD
☐ OPCP	One Piece Malleable Calf Panel ⁽¹⁾	\$363
☐ LNXPF	Plantar Flex Option for LNX Foot platform ⁽²⁾ ...	\$363

1. Not available with Plantar Flex Option for LNX Foot platform, or with LNXA.

2. Only available with Foot Platform Options.

3. For clients under 4' only.

MANUAL ELEVATING CENTER MOUNT LEGREST

☐ MECML	Manual Elevating Center Mount Legrest	
	With Plantar Flex	\$863

Select legrest set up position:

☐ 90SU	Set at 90° (elevating range 90-40°)	STD
☐ 110TU	Tucked at 110° (elevating range 110-60°) ⁽¹⁾	N/C

Rubber Coated Footplate:

☐ MRCP	Rubber coated footplate 11.5" W x 10" D	\$239
	Select seat pan to footplate for Rubber Coated Footplate	
☐ CSPTF1	11.5" - 15"	N/C
☐ CSPTF2	14" - 19"	N/C

Individual Footplates:

☐ IFSL	Small - (5" x 7.5" deep) - Left	\$239
☐ IFSR	Small - (5" x 7.5" deep) - Right	\$239
☐ IFML	Medium - (5.5" x 9.5" deep) - Left	\$239
☐ IFMR	Medium - (5.5" x 9.5" deep) - Right	\$239
☐ IFL	Large - (6" x 11.5" deep) - Left	\$239
☐ IFLR	Large - (6" x 11.5" deep) - Right	\$239


Select seat pan to footplate for Individual Footplates

☐ CSPTF3	9.5" - 14.5"	N/C
☐ CSPTF4	12" - 19"	N/C

Optional:

☐ CICP	Individual calf pads	\$290
-------------------------------	----------------------------	-------

1. Not available on 17" seat depth.



Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

STANDARD FRONT RIGGING


- ☐ FCMP Fixed Center Mount Platform.....N/C
☐ FBIFP Comes with a standard footplate
 17" W x 13" D.....N/C

Select seat pan to footplate:

- ☐ FSPTF1 9" - 13".....N/C
☐ FSPTF14 13" -19".....N/C

Set at:

- ☐ FSA70 Set at 70°N/C
☐ FSA83 Set at 83°N/C
☐ FSA90 Set at 90°N/C
☐ FSA97 Set at 97°N/C



MANUAL CENTER MOUNT FRONT RIGGING

- ☐ MSFCMFP Maxx Style Fixed Center Mount
 Foot Platform⁽¹⁾.....\$573

- ☐ MBIFP Comes with a rubber coated footplate
 17" W x 13" D.....N/C


- ☐ MICP Comes standard with individual calf pads.....N/C

Select seat pan to footplate:

- ☐ MSPTF1 9" - 13".....N/C
☐ MSPTF14 13" -19".....N/C


Set at:

- ☐ MSA70 Set at 70°N/C
☐ MSA83 Set at 83°N/C
☐ MSA90 Set at 90°N/C
☐ MSA97 Set at 97°N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.





FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP



- ☐ MFFP Medium Platform (12"w x10.75"d)\$261
☐ LFFP Large Platform (14.5"w x13"d)\$363
☐ XLFFP Extra Large Platform (17"w x13"d)\$363
☐ XXLFFP Extra Extra Large Platform (21"w x13"d)\$363

1. May need to adjust leg position to avoid interference with front casters

INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP



- ☐ IFSL Small (5" W x 7.5" D) - Left\$239
☐ IFSR Small (5" W x 7.5" D) - Right.....\$239
☐ IFML Medium (5.5" W x 9.5" D) - Left\$239
☐ IFMR Medium (5.5" W x 9.5" D) - Right.....\$239
☐ IFLL Large (6" W x 11.5" D) - Left.....\$239
☐ IFLR Large (6" W x 11.5" D) - Right\$239

Add:

- ☐ TSSL Left Toe Strap Small 6.5"\$181
☐ TSSR Right Toe Strap Small 6.5".....\$181
☐ TSML Left Toe Strap Medium 8.5"\$181
☐ TSMR Right Toe Strap Medium 8.5"\$181
☐ TSLL Left Toe Strap Large 10.5"\$181
☐ TSLR Right Toe Strap Large 10.5".....\$181

- ☐ MFOL Multi-Axis Footplate (6"W x 10"D) - Left\$326
☐ MFOR Multi-Axis Footplate (6"W x 10"D) - Right\$326

1. Foot plates may interfere with casters in certain configurations

MAXX INDEPENDENT ADJUSTABLE POSITIONING



- 5.5" - 6" WIDTH ADJUSTMENT
- 9.5" - 11.5" DEPTH ADJUSTMENT
- UP TO 9° ADDUCTION
- UP TO 45° ABDUCTION

- ☐ MPFPL Maxx Independent Adjustable Positioning
 Foot plate⁽¹⁾ - Left\$290
☐ MPFPR Maxx Independent Adjustable Positioning
 Foot plate⁽¹⁾ - Right\$290
☐ MTSL Left Toe Strap (each) HCPCS code E0952\$51
☐ MTSR Right Toe Strap (each) HCPCS code E0952\$51

1. Foot plates may interfere with casters in certain configurations.

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MAXX STYLE POWER ELEVATING LEGRESTS

- ☐ MSPPL MAXX Swing Away
Power Elevating⁽¹⁾ HCPCS code E1010 \$3,611

Select Power Legrest Function:

- ☐ INDA Individual Movement N/C
☐ COMA Combined Movement N/C
☐ SYNA Synchronized/Recline Movement N/C

Select Receivers:

- ☐ MRS Straight Receivers N/C
☐ MRF Flared 1" Outward Receivers N/C
☐ MRF2 Flared 2" Outward Receivers \$363

Select Seat Pan to Footplate:

- ☐ PSPTF5 Length: 10"-14" N/C
☐ PSPTF2 Length: 13"-17" N/C
☐ PSPTF7 Length: 16"-20" N/C

Optional:

- ☐ PMCP Moulded I-Skin Calf Pads STD
☐ PFCPMS Flat Calf Pad Option (pair) \$399
☐ PCCPMS Curved Calf Pad Option⁽²⁾ (pair) \$399

1. Must choose footplate on page 15.

2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

MAXX STYLE MANUAL ELEVATING LEGRESTS

- ☐ MMPPSA MAXX Swing Away
Manual Elevating⁽¹⁾ HCPCS code E0990 \$986

Select Receivers:

- ☐ ERS Straight Receivers N/C
☐ ERF Flared 1" Outward Receivers N/C
☐ ERF2 Flared 2" Outward Receivers \$363

Select Seat Pan to Footplate:

- ☐ ESPTF5 Length: 10"-14" N/C
☐ ESPTF2 Length: 13"-17" N/C
☐ ESPTF7 Length: 16"-20" N/C

Optional:

- ☐ EMCP Moulded I-Skin Calf Pads STD
☐ EFCPMS Flat Calf Pad Option (pair) \$399
☐ ECCPMS Curved Calf Pad Option⁽²⁾ (pair) \$399

1. Must choose footplate on page 15.

2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

SWING AWAY LEGREST

- ☐ HD70 Heavy Duty 70° Swing Away Footrest(pair)⁽¹⁾⁽²⁾ .. \$718

Select Receivers:

- ☐ RS Straight Receivers N/C
☐ RF Flared 1" Outward Receivers N/C
☐ RF2 Flared 2" Outward Receivers \$363

Select Seat Pan to Footplate:

- ☐ HSPTF11 Length: 6"-10" N/C
☐ HSPTF8 Length: 10"-14" N/C
☐ HSPTF2 Length: 12"-16" N/C
☐ HSPTF7 Length: 15"-19" N/C

Options:

- ☐ HD70FCP Flat Calf Pad Option (pair) \$399
☐ HD70CCP Curved Calf Pad Option (pair) \$399

- ☐ RLRL Residual Limb Support - Left (8.5" x 6")⁽³⁾
HCPCS code E1020 \$718
☐ RLRS Residual Limb Support - Right (8.5" x 6")⁽³⁾
HCPCS code E1020 \$718

1. Must choose footplate on page 14..

2. Heavy duty durable construction (4.5lbs ea.).

3. Swing away and removable, 50° of angle adjustment, 4.5" height adjustment, 3" horizontal and 2" of for/aft adjustment.

DUAL FRONT RIGGING ADAPTOR⁽¹⁾

- ☐ DFRA Dual Front Rigging Receivers. \$653

1. Must purchase each set of front rigging separately. Not available on no charge front rigging options.

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____





Dealer Account Number: _____

PO #: _____

Tag: _____

FOOTPLATES - SWING AWAY/ELEVATING LEGRESTS⁽¹⁾

HCPCS code K0040






Angle Adjustable Footplates

☐ AFSL	Left Small Footplate (4.25"W x 8"D)	\$239
☐ AFSR	Right Small Footplate (4.25"W x 8"D)	\$239
☐ AFML	Left Medium Footplate (5.25"W x 8"D)	\$239
☐ AFMR	Right Medium Footplate (5.25"W x 8"D)	\$239
☐ AFLL	Left Large Footplate (6.25"W x 8"D)	\$239
☐ AFLR	Right Large Footplate (6.25"W x 8"D)	\$239
☐ AFXL	Left Extra Large Footplate (9"W x 11"D)	\$239
☐ AFXR	Right Extra Large Footplate (9"W x 11"D)	\$239

Multi-Axis Adjustable Angle Footplates

☐ MFSL	Left Small Footplate (4.25"W x 8"D)	\$326
☐ MFSR	Right Small Footplate (4.25"W x 8"D)	\$326
☐ MFML	Left Medium Footplate (5.25"W x 8"D)	\$326
☐ MFMR	Right Medium Footplate (5.25"W x 8"D)	\$326
☐ MFLL	Left Large Footplate (6.25"W x 8"D)	\$326
☐ MFLR	Right Large Footplate (6.25"W x 8"D)	\$326
☐ MFXL	Left Extra Large Footplate (9"W x 11"D)	\$326
☐ MFXR	Right Extra Large Footplate (9"W x 11"D)	\$326



Cut Out Angle Adjustable Footplates

☐ AFCSL	Left Small Footplate (7.25"W x 9"D)	\$239
☐ AFCSR	Right Small Footplate (7.25"W x 9"D)	\$239
☐ AFCML	Left Medium Footplate (7.5"W x 11"D)	\$239
☐ AFCMR	Right Medium Footplate (7.5"W x 11"D)	\$239
☐ AFCLL	Left Large Footplate (10"W x 14"D)	\$239
☐ AFCLR	Right Large Footplate (10"W x 14"D)	\$239

Cut Out Multi-Axis Adjustable Angle Footplates

☐ MFCSL	Left Small Footplate (7.25"W x 9"D)	\$326
☐ MFCSR	Right Small Footplate (7.25"W x 9"D)	\$326
☐ MFCML	Left Medium Footplate (7.5"W x 11"D)	\$326
☐ MFCMR	Right Medium Footplate (7.5"W x 11"D)	\$326
☐ MFCLL	Left Large Footplate (10"W x 14"D)	\$326
☐ MFCLR	Right Large Footplate (10"W x 14"D)	\$326

1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

HEEL LOOPS, CUPS AND TOE STRAPS

☐ HLL	Left Heel Loop (each) HCPCS code E0951	\$29
☐ HLR	Right Heel Loop (each) HCPCS code E0951	\$29
☐ HCL	Left Cup Loop (each) HCPCS code E0951	\$29
☐ HCR	Right Cup Loop (each) HCPCS code E0951	\$29
☐ TSL	Left Toe Strap (each) HCPCS code E0952	\$29
☐ TSR	Right Toe Strap (each) HCPCS code E0952	\$29

CALF STRAPS AND CALF PANELS⁽¹⁾

☐ CSM	Body Point Airmesh Calf Strap for Chairs 15"-18" wide with 3.5" H Pad ⁽¹⁾	\$218
☐ CSL	Body Point Airmesh Calf Strap for Chairs 18"-22" wide with 3.5" H Pad ⁽¹⁾	\$218
☐ CPM	Body Point Airmesh Calf Panel for Chairs 15"-18" wide with 35.5" W x 9" H Pad ⁽¹⁾	\$319
☐ CPL	Body Point Airmesh Calf Panel for Chairs 18"-22" wide with 38.5" W x 9" H Pad ⁽¹⁾	\$341

1. Not available for Center Mounted Front Rigging.

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SINGLE-DUAL FUNCTION ELECTRONICS

- ☐ SFCB Single/Dual Function Electronics⁽¹⁾ STD
- Select Switch:
- ☐ SPBSO Single Push Button
- ☐ STO Single Toggle
- ☐ DPBO Dual Push Button
- ☐ OSO Omit Switch
- Optional:
- ☐ RDE Reduced Drive for Elevating Seat⁽²⁾ \$580
-
- ☐ SFCBA Single/Dual Function Drive Control⁽¹⁾
- HCPSC code E2310 \$2,190
- Optional:
- ☐ RDES Reduced Drive for Elevating Seat⁽²⁾ \$580
-
- ☐ DFCBD Dual Function Through Drive Control⁽¹⁾
- HCPSC code E2311 \$2,919
- Optional:
- ☐ RDED Reduced Drive for Elevating Seat⁽²⁾ \$580

1. These electronics include drive lock-out and tilt limit.

2. If reduced drive is not selected the system will be in drive lock-out when elevated.

MULTI- FUNCTION ELECTRONICS

- ☐ MFCBS Multi-Function through Switch⁽¹⁾ N/C
- Select Switch:
- ☐ QPBO Quad Push Button N/C
- ☐ FWT Four-Way Toggle Switch \$218
- ☐ EWR 8-Way Rocker \$297
- Optional:
- ☐ RDEM Reduced Drive for Elevating Seat⁽²⁾ \$480
-
- ☐ MFCBD Multi-Function M610 Through Drive Control⁽¹⁾ HCPSC code E2311 \$3,504
- Up to 6 Seating Functions
 - On Board Programming Seating Functions
 - Specialty Products (i.e. PSB Power Sliding Back)
 - Includes reduced drive for Elevating Seat

1. Requires Expandable Controller (EXPC) & Harness for Expandable System (PWH), select on page 2.

2. If reduced drive is not selected system will be in drive lock-out when elevated.

3. These electronics include seating attendant control, drive lock-out, tilt/recline limit and actuator speed control.

ELECTRONIC SWITCH MOUNTING

- ☐ ESML Left Side (Armpad)
- ☐ ESMR Right Side (Armpad)
- ☐ ESMJL Joystick Left (Mounts to Joystick Bracket)
- ☐ ESMJR Joystick Right (Mounts to Joystick Bracket)
- ☐ ESMD Do not mount

ELECTRONICS OPTIONS

Additional Switches

- ☐ SPBSA Single Push Button \$218
- ☐ STA Single Toggle Switch⁽⁴⁾ \$218
- ☐ DPBA Dual Push Button⁽⁴⁾ \$218
- ☐ QPBA Quad Push Button⁽⁴⁾ \$218
- ☐ 10WA 10-Way Push Button⁽¹⁾⁽⁴⁾ \$370
- ☐ ASL320R ASL Pulse Mechanical Switch - Red \$181
- ☐ ASL320G ASL Pulse Mechanical Switch - Green \$181
- ☐ ASL320BK ASL Pulse Mechanical Switch - Black \$181
- ☐ SESG Stealth Egg Switch - Green \$196
- ☐ SESB Stealth Egg Switch - Black \$196
- ☐ FTSQ Feather Touch - Single Quarter Size \$295
- ☐ FTMD Feather Touch - Dual Dime Size⁽¹⁾⁽³⁾ \$428
- ☐ FTMQ Feather Touch - Quad Dime Size⁽¹⁾⁽³⁾ \$428

Additional Switches (not available with ACT4 Electronics)

- ☐ FWTM Four-Way Toggle Switch⁽¹⁾⁽³⁾⁽⁶⁾ \$218
- ☐ EWRA 8-Way Rocker⁽¹⁾⁽³⁾⁽⁶⁾ \$297
- ☐ TWRA 10-Way Rocker⁽¹⁾⁽³⁾⁽⁶⁾ \$297
- ☐ TWTA 10-Way Toggle⁽¹⁾⁽³⁾⁽⁶⁾ \$297
- ☐ EWPA 8-Way Push Button⁽¹⁾⁽³⁾⁽⁶⁾ \$370

Specify Use for Additional Switch:

- ☐ STG Seating
- ☐ JM Joystick Mode⁽²⁾⁽⁵⁾
- ☐ JOO Joystick On/Off⁽²⁾⁽⁵⁾

Specify Mounting for Additional Switch:

- ☐ ESMLA Left Side (Armpad)
- ☐ ESMRA Right Side (Armpad)
- ☐ ESMJLA Joystick Left (Mounts to Joystick Bracket)
- ☐ ESMJRA Joystick Right (Mounts to Joystick Bracket)
- ☐ ESMDA Do not mount

Miscellaneous

- ☐ MJBS Motion Junction Box for Switch Control⁽¹⁾⁽⁵⁾
For Mode and/or Emergency Stop \$254
- ☐ S46L Single to Dual Switch Adapter \$384
Allows for dual use of a single switch
- ☐ M616S M616 Converter⁽¹⁾ HCPSC code E2399 \$181
Allows up to 4 mono phono switches into Multi-Function box
- ☐ LJS Lazarus Jump Starter HCPSC code A9270 \$218
- ☐ APPS Accessory Port Power Supply⁽¹⁾
3 amp capacity \$283
- ☐ TPSMH 12 Piece Switch Mounting Hardware
Kit (additional kits) \$57

1. Not available with Single-Dual Function Electronics.

2. IN200 is required with REM210, REM211, REM216 joysticks.

3. Not available with REM110 joystick.

4. Not available for Joystick Mode or Joystick On/Off. Must Select Seating (STG)

5. Cannot be used for power On/Off.

6. Not Available on Multi Function ACT4-LiNX Electronics.

Ultra Low Maxx CG Power Positioning System

TDX-SP2 HD

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ARTICULATING TRAYS⁽¹⁾

- ☐ USAT Universal Style Articulating Tray \$1,958
14"W x Adjustable 6"-10"D
- Add:
- ☐ TMA Trilogy Mounting Adapter \$73
- ☐ AVT Articulating Vent (LTV) Tray \$1,958
- ☐ AVBH Articulating Vent Bag Holder⁽²⁾ \$1,958
- Add:
- ☐ ABMC Accessory Bag Mounting and Clips \$145

1. The user weight limit is 225 lbs. The maximum weight capacity of the vent tray is 35 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
2. Motion Concepts does not supply the Vent Bag.

GIMBAL VENT TRAYS⁽¹⁾

- ☐ GUT Gimbal Universal Tray \$1,958
14"W x Adjustable 6"-10"D
- Add:
- ☐ GTMA Trilogy Mounting Adapter \$73
- ☐ GVT Gimbal Vent (LTV) Tray \$1,958
- ☐ GVBH Gimbal Vent Bag Holder⁽²⁾ \$1,958
- Add:
- ☐ GABMC Accessory Bag Mounting and Clips \$145

1. The Gimbal Vent Trays are not available with Tilt Only and/or Elevate Only Systems. The user weight limit is 200 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. The maximum weight capacity of the vent tray is 20 lbs.
2. Motion Concepts does not supply the Vent Bag.

O2 HOLDERS⁽¹⁾

- ☐ NAO2H Non-Articulating O2 Gas Holder \$573
- ☐ AOHG Articulating O2 Gas Holder Cylinder Style⁽²⁾⁽³⁾ ... \$1,813
- ☐ AOHL Articulating Liquid O2 Holder (6.5" x 4.5")⁽²⁾⁽³⁾ ... \$1,813
- ☐ NAIP L Non-Articulating IV Pole User's Left⁽⁴⁾ \$573
- ☐ NAIP R Non-Articulating IV Pole User's Right⁽⁴⁾ \$573
- ☐ AIVP Articulating IV Pole⁽²⁾⁽³⁾⁽⁴⁾ \$1,813

1. The user weight limit is 225 lbs. The maximum weight capacity of the O2 holder is 20 lbs.
2. Required for use on Recline Systems
3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. The maximum weight capacity of the IV pole is 5 lbs.

VENTILATOR INFORMATION Required when ordering a Vent Tray

- ☐ BIB Built in battery
- ☐ DTSVPS Dealer to supply vent power source

REAR STORAGE BASKET⁽¹⁾⁽²⁾

- ☐ RSB Removable Rear Storage Basket \$1,878
Exterior dimensions: 18"W x 6"H x 7"D
Interior dimensions: 17.8"W x 6"H x 6.8"D
Interior capacity: 3 gallon
Max Load: 20 lbs

1. Tall items extending beyond the basket's height may be damaged upon tilting.
2. Not available with XRO, XROE brackets

ACCESSORIES

- ☐ MUSB Motion Dual USB Charger N/C
- ☐ UCH Universal Cup Holder⁽³⁾ \$51
- ☐ PH X-Grip Phone Holder \$109
- ☐ AMC Accessory Mount with Clip (QTY ____) \$73
- ☐ CRVMR Rear View Mirror Right Mounted \$65
- ☐ CRVML Rear View Mirror Left Mounted \$65
- ☐ TAG571850 Wheelchair Cover \$152
- ☐ HK Hydration Kit⁽¹⁾ \$363
- ☐ MUB Medical Utility Bag w/ Mounting Hardware. \$363
- ☐ THS Transfer Handles 4" - 6.5" \$508
- ☐ THL Transfer Handles 9" - 11.5" \$508
- ☐ FFLTL Flip and Fold Away 1/2 Lap Tray
Left Mounting(14"W x 10"D)⁽¹⁾⁽²⁾ \$1,153
Flip and Fold Away 1/2 Lap Tray
Right Mounting(14"W x 10"D)⁽¹⁾⁽²⁾ \$1,153
- ☐ FFLTR
- ☐ CFH Flag Holder \$94
- ☐ Left Mount
- ☐ Right Mount
- ☐ CSCC Cane and Crutch Holder \$573
- ☐ Left Mount
- ☐ Right Mount

1. Not available on joystick side and only available with a flat arm pad or modular arm pad. Not available with Outback Cantilever Arms.
2. The maximum weight capacity of the tray is 10 lbs.
3. Not available on the joystick side

BATTERIES

- ☐ M24 SLD G MK Group 24 NF Gel Battery (ea) Quantity()
HCPCS code E2363 \$529

INSTALLATION

- ☐ IVCI Motion Concepts Installation⁽¹⁾ N/C
- ☐ MCI Motion Concepts Installation⁽¹⁾ N/C

1. Batteries may be required if items on Accessories order form are selected at time of installation, please call to confirm.

NOTES

RETAIL \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM

Ultra Low Maxx CG Power Positioning System

Dealer Name: _____

Dealer Account Number:_____

PO #: _____

Tag: _____

ADDITIONAL MEASUREMENTS FOR BARIATRIC CLIENTS

ADDITIONAL NOTES

To ensure system is accurately configured please fill in all required measurements.

Current weight:	
Weight history:	
1. Back of knee /calf to back of buttocks (seat pan depth)	
2. Back of knee/calf to thoracic-lumber trunk (for seat depth)	
3. Seat pan to under forearm (armrest height)	
4. Seat pan to top of gluteal tissue (lower aspect of back support height)	
5. Width at toes (lateral aspect)	
6. Width from lateral calf to lateral calf (at widest aspect)	
7. Overall hip width	

