



ULTRA LOW MAXX REHAB SEAT ON TDX-SP2 AND TDX-SP2 HD BASE CDN PRICE LIST AND ORDER FORM

Price Effective July 15, 2026



PRINT FORM

SAVE

CLEAR FORM

REQUIRED INFORMATION

Quote ☐ Order ☐

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured
please fill in all REQUIRED* measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness	



Customer Service: 1.866.748-7943 | Fax: 1.905.695.0583 | www.motionconcepts.com

INVACARE POWER BASES

- ☐ TDXSP2 TDX-SP2 Base..... \$7,575
- Select STF(1): ☐ SLOW 16.75^{”(2)}..... \$250
- ☐ LOW 17.25"
- ☐ MED 18.25"
- ☐ TALL 19.25"

- ☐ U300 Weight Capacity up to 300 lbs..... STD

- ☐ 22TRY Narrow Base 24"..... N/C

- ☐ 24TRY Standard Base 25.5"..... N/C

- ☐ TDXSP2HD TDX-SP2 HD Base..... \$8,275
- Select STF(1): ☐ HD-SLOW 17.25^{”(2)}..... \$250
- ☐ HD-LOW 17.75"
- ☐ HD-MED 18.75"
- ☐ HD-TALL 19.75"

- ☐ U400 Weight Capacity up to 400 lbs..... STD

- ☐ 24TRY Group 24 Style Tray - Standard Base 25.5"..... N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size.
2. Must order 24TRY battery tray type and Transport Brackets TRBKTS.
Only available on 17" seat depths and greater.

PERFORMANCE MOTORS

- ☐ HSMTR Performance Motors - 7.5 mph⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾..... \$995

1. Must order 24TRY Battery Tray Type
2. Extended high speed use on rugged terrain is not recommended.
3. Performance/High Speed Motors may have less torque and may have less range than standard speed motors.
4. Not available on TDX-SP2 HD base

TRANSPORT TIE DOWN

- ☐ TRBKTS Wheelchair Transport Brackets⁽¹⁾..... STD

1. For unoccupied use only

TIRE OPTIONS

- ☐ B1431-3 14" x 3" Black Tire with Gel Foam Inserts STD
- ☐ B1430-3T 14" x 3" Black Tire - Pneumatic Tire N/C

FORK OPTIONS

- ☐ DSFK Double-Sided Fork Package STD
- ☐ FKPKG Single-Sided Fork Package \$375

INSTALLATION OF BATTERIES

- ☐ 22NFBATTERY 22 NF Battery (On Chair)..... \$700
- ☐ 24BATTERY 24 Gel Battery (On Chair) \$840

BATTERY CHARGER OPTION

- ☐ 110CHARGER 110 Volt Battery Charger N/C
- ☐ L8900 Omit Battery Charger N/C

LIGHTS AND INDICATORS

- ☐ LIGHTS LED Lights and Indicators⁽¹⁾ \$995

1. Must order REM216, REM400 or REM500

LINXS ELECTRONICS ACCESSORIE

- ☐ LAK LINX Access Key⁽¹⁾ \$50
- ☐ MUSB Motion Dual USB Charger N/C

1. Required for programming

ADDITIONAL ACCESSORIES

- ☐ UCH Universal Cup Holder⁽¹⁾ \$50

1. Not available on the joystick side.

FRAME FINISH

☐ 170P		Obsidian Black.....	N/C
☐ 171P		Ruby Red.....	N/C
☐ 172P		Peacock Blue.....	N/C
☐ 174P		Purple Passion.....	N/C
☐ 175P		Denim Navy.....	N/C
☐ 176P		Lime Twist.....	N/C
☐ 177P		Arctic Blue.....	N/C

RIM INSERT COLORS - MUST PICK ONE

☐ 170PR		Obsidian Black.....	N/C
☐ 171PR		Ruby Red.....	N/C
☐ 172PR		Peacock Blue.....	N/C
☐ 173PR		Pink Flamingo.....	N/C
☐ 174PR		Purple Passion.....	N/C
☐ 175PR		Denim Navy.....	N/C
☐ 176PR		Lime Twist.....	N/C
☐ 177PR		Arctic Blue.....	N/C
☐ 178PR		Banana Split.....	N/C
☐ 179PR		Orange Marmalade.....	N/C
☐ 180PR		Morning Mist.....	N/C

CONTROLLERS

☐ PWH	Harness Required for Expandable System.....	STD
☐ GTRAC-LX	LiNX G-TRAC Module.....	\$1,200

NON-EXPANDABLE DRIVE CONTROL

☐ REM110	LED Drive Only Remote/Joystick ⁽¹⁾	N/C
	Single drive profile	
☐ REM210	LED Non-Expandable Remote/Joystick ⁽¹⁾	N/C
	Three drive profiles, control of two seat functions	

1. Not available with lights and indicators.
2. Must choose lighting package separately

EXPANDABLE JOYSTICKS⁽¹⁾

☐ REM211	LED Expandable Remote/Joystick ⁽²⁾	\$875
	Three drive profiles, control of up to 6 actuators. Not for use with specialty controls.	
☐ REM216	LED Expandable Remote/Joystick with Lights ⁽³⁾	\$875
	Three drive profiles, control of up to 6 actuators and light operation. Not for use with specialty controls.	
☐ REM400	Color 3.5" Touch Screen Remote/Joystick ⁽⁴⁾	\$1,000
	Four drive profiles, control of up to 6 actuators through the drive control and light operation.	
☐ SW400TGL	Toggle Kit for REM400 ⁽⁵⁾	\$200

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Not available with lights and indicators
3. Must choose lighting package separately.
4. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.
5. Adds 2 toggle switches to REM400.

ATTENDANT DRIVE CONTROL CHOICES

☐ ACU	Proportional Attendant Control ⁽¹⁾	\$1,024
	Select mounting position for ACU	
☐ RPBMR	Right	
☐ RPBML	Left	
☐ RPBMC	Center	

MOUNTING FOR DRIVER CONTROLS

☐ FM-ML	Fixed Mounting Left.....	N/C
☐ FM-MR	Fixed Mounting Right.....	N/C
☐ SAQR	Motion Height Adjustable Swing-Away Quad Link (Right Hand).....	\$344
☐ SAQL	Motion Height Adjustable Swing-Away Quad Link (Left Hand).....	\$344

REHAB SEAT

☐ D0195	Ultra Low Maxx Rehab Seat ⁽¹⁾	\$2,995
--------------------------------	--	---------

1. Must select TDX-SP2 Base

BARIATRIC REHAB SEAT

☐ MBRS	Ultra Low Maxx Bariatric Rehab Seat ¹⁽²⁾	\$3,995
-------------------------------	---	---------

1. Not available with 12" Elevating Seat Lift Module.
2. Must order TDX-SP2 HD base.

ELEVATING SEAT MODULE

☐ MD12L 12" Elevating	Lift for Rehab Seat ⁽¹⁾	\$1,995
	Weight capacity up to 250lbs	

☐ SFCB	Single Function Electronics ⁽²⁾
-------------------------------	--

Select Switch:

☐ SPBO	Single Push Button
☐ STO	Single Toggle
☐ DPBO	Dual Push Button
☐ SOS	Omit Switch
☐ JOY	Through Joystick

Select Switch Mounting:

☐ ESML	Left Armpad
☐ ESMR	Right Armpad
☐ ESMJL	Left Joystick (Mounts to Joystick Bracket)
☐ ESMJR	Right Joystick (Mounts to Joystick Bracket)
☐ ESMD	Do not mount

Optional:

☐ RDE	Reduced Drive for Elevating Seat ⁽³⁾	\$400
------------------------------	---	-------

1. Not available on seat sizes greater than 20" wide by 20" deep.
2. Required when selecting 12" Elevating Lift. These electronics include drive lock-out.
3. If reduced drive is not selected the system will be in drive lock-out when elevated.

ADD: MODULES

☐ PRETILT 5° Fixed Pre-Tilt ⁽¹⁾	N/C
---	-----

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height.

Maxx Rehab Seat
Invacare TDX-SP2 / TDX-SP2 HD Base

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

ADULT SEAT WIDTH FOR TDX SP2⁽¹⁾

☐ SW116	16" Wide (Adjustable from 16" to 20")	N/C
☐ SW117	17" Wide (Adjustable from 16" to 20")	N/C
☐ SW118	18" Wide (Adjustable from 16" to 20")	N/C
☐ SW119	19" Wide (Adjustable from 16" to 20")	N/C
☐ SW120	20" Wide (Adjustable from 16" to 20")	N/C
☐ W219	19" Wide (Adjustable from 19" to 22")	\$460
☐ W220	20" Wide (Adjustable from 19" to 22")	\$460
☐ W221	21" Wide (Adjustable from 19" to 22")	\$460
☐ W222	22" Wide (Adjustable from 19" to 22")	\$460

1. Not available on TDX-SP2 HD base. You must select both Adult Seat Width and Depth together.

ADULT ULTRA RAIL SEAT DEPTH FOR TDX SP2⁽¹⁾

☐ SD116	16" Deep (Adjustable from 16" to 21")	N/C
☐ SD117	17" Deep (Adjustable from 16" to 21")	N/C
☐ SD118	18" Deep (Adjustable from 16" to 21")	N/C
☐ SD119	19" Deep (Adjustable from 16" to 21")	N/C
☐ SD120	20" Deep (Adjustable from 16" to 21")	N/C
☐ D219	19" Deep (Adjustable from 19" to 22")	\$460
☐ D220	20" Deep (Adjustable from 19" to 22")	\$460
☐ D221	21" Deep (Adjustable from 19" to 22")	\$460
☐ D222	22" Deep (Adjustable from 19" to 22")	\$460
☐ D223	23" Deep (Adjustable from 19" to 22")	\$460

1. Not available on TDX-SP2 HD base. You must select both Adult Seat Width and Depth together.

BARIATRIC SEAT WIDTH FOR TDX SP2 HD⁽¹⁾

☐ BW120	20" Wide (Adjustable from 20" to 24")	N/C
☐ BW121	21" Wide (Adjustable from 20" to 24")	N/C
☐ BW122	22" Wide (Adjustable from 20" to 24")	N/C
☐ BW123	23" Wide (Adjustable from 20" to 24")	N/C
☐ BW124	24" Wide (Adjustable from 20" to 24")	N/C
☐ BW222	22" Wide (Adjustable from 22" to 26")	N/C
☐ BW223	23" Wide (Adjustable from 22" to 26")	N/C
☐ BW224	24" Wide (Adjustable from 22" to 26")	N/C
☐ BW225	25" Wide (Adjustable from 22" to 26")	N/C
☐ BW226	26" Wide (Adjustable from 22" to 26")	N/C

1. Not available on TDX-SP2 base. You must select both Bariatric Seat Width and Depth together.

BARIATRIC ULTRA RAIL SEAT DEPTH FOR TDX SP2 HD⁽¹⁾

☐ SD117	17" Deep (Adjustable from 17" to 21")	N/C
☐ SD118	18" Deep (Adjustable from 17" to 21")	N/C
☐ SD119	19" Deep (Adjustable from 17" to 21")	N/C
☐ SD120	20" Deep (Adjustable from 17" to 21")	N/C
☐ SD219	19" Deep (Adjustable from 19" to 23")	N/C
☐ SD220	20" Deep (Adjustable from 19" to 23")	N/C
☐ SD221	21" Deep (Adjustable from 19" to 23")	N/C
☐ SD222	22" Deep (Adjustable from 19" to 23")	N/C
☐ SD223	23" Deep (Adjustable from 19" to 23")	N/C

1. Not available on TDX-SP2 base. You must select both Bariatric Seat Width and Depth together.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

LAP BELTS

☐ BELT48	Push Button Style Lap Belt, Length 48"	STD
☐ BELT60	Push Button Style Lap Belt, Length 60"	N/C
☐ BELT71	Push Button Style Lap Belt, Length 71"	N/C
☐ PLB2	Bodypoint Padded Lap Belt 2 Point ⁽¹⁾	\$250
☐ PLB4	Bodypoint Padded Lap Belt 4 Point ⁽¹⁾	\$315

1. 67" length only option

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾

☐ 18" FBT18	☐ 19" FBT19	☐ 20" FBT20	☐ 21" FBT21	☐ 22" FBT22	☐ 23" FBT23	☐ 24" FBT24	☐ 25" FBT25
------------------------------------	------------------------------------	------------------------------------	------------------------------------	------------------------------------	------------------------------------	------------------------------------	------------------------------------

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK CANES - NON RECLINE SYSTEMS

☐ SBC18	18" Straight Back Canes	STD
☐ SBC20	20" Straight Back Canes	STD
☐ SBC22	22" Straight Back Canes	STD
☐ SBC24	24" Straight Back Canes	STD

Select Cane Angle

☐ 75° BA75	☐ 79° BA79	☐ 85° BA85	☐ 90° BA90	☐ 95° BA95	☐ 101° BA101	☐ 106° BA106	☐ 112° BA112	☐ 116° BA116
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-------------------------------------	-------------------------------------	-------------------------------------	-------------------------------------

☐ ABC20	20" Angled 5° Back Canes	N/C
☐ ABC22	22" Angled 5° Back Canes	N/C
☐ ABC24	24" Angled 5° Back Canes	N/C

Select Cane Angle

☐ 81° BA81	☐ 84° BA84	☐ 91° BA91	☐ 96° BA96	☐ 101° BA101	☐ 106° BA106	☐ 111° BA111	☐ 117° BA117	☐ 121° BA121
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-------------------------------------	-------------------------------------	-------------------------------------	-------------------------------------	-------------------------------------

BACK OPTIONS

☐ SRB	Standard Rehab Back & Cushion	\$125
☐ ORBC	Omit Rigid Back & Cushion ⁽³⁾	N/C

ARMRESTS - NON RECLINE SYSTEMS

☐ DP1	Adult Dual Post 9.75" - 12.75"	STD
☐ DP2	Adult Dual Post 12.75" - 15.75"	STD
☐ DP3	Adult Dual Post 15.75" - 18.75"	STD

Select:

☐ FP	Flat Pouch	STD
☐ GB	Glove Boxes	N/C

Optional:

☐ OAL	Outback Arm Option Left ⁽¹⁾ each	\$195
☐ OAR	Outback Arm Option Right ⁽¹⁾ each	\$195

☐ RMCAN	Ultra Rail Seat Mnt Flip Back Cantilever MAXX Tilt Arm 9" - 12.5" ⁽²⁾	\$488
--------------------------------	---	-------

Optional:

☐ DPM	Dual Post Module to RMCAN (pair)	\$250
------------------------------	----------------------------------	-------

☐ OCTA1	Outback Cantilever Tilt Armrest 9"-12" ⁽³⁾	\$488
☐ OCTA2	Outback Cantilever Titl Armrest 12"-15" ⁽³⁾	\$488

1. Only available with Dual Post Armrests.

2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.

ARMREST TUBE LENGTH - PLEASE SPECIFY

☐ ATFL	Left Full	STD
☐ ATDL	Left Desk	STD
☐ ATFR	Right Full	STD
☐ ATDR	Right Desk	STD

ARMPADS

Standard

☐ SAPFL	Left Full (2.25" x 14")	STD
☐ SAPDL	Left Desk (2.25" x 10")	STD
☐ SAPML	Left Mini (2.25" x 7")	STD
☐ SAPFR	Right Full (2.25" x 14")	STD
☐ SAPDR	Right Desk (2.25" x 10")	STD
☐ SAPMR	Right Mini (2.25" x 7")	STD

Waterfall Arm Pads⁽¹⁾

☐ WAPFL	Left Full (13")	N/C
☐ WAPDL	Left Desk (9")	N/C
☐ WAPFR	Right Full (13")	N/C
☐ WAPDR	Right Desk (9")	N/C

Gel Pads

☐ GAPNL	Narrow Left (2" x 12")	\$88
☐ GAPNR	Narrow Right (2" x 12")	\$88
☐ GAPWL	Wide Left (3.5" x 12")	\$96
☐ GAPWR	Wide Right (3.5" x 12")	\$96
☐ GAPLL	Long Left (3.5" x 14")	\$103
☐ GAPLR	Long Right (3.5" x 14")	\$103

1. Full length pad requires full length arm tube if using Ultra Rail mounted Tilt Arms.

MANUAL CENTER MOUNT FRONT RIGGING

☐ MSFCMFP	Maxx Style Fixed Center Mount Foot Platform ⁽¹⁾	\$395
----------------------------------	---	-------

☐ IFP	Medium Rubber Coated Footplate 11.5" W x 10"	N/C
------------------------------	---	-----

☐ IFP	Large Rubber Coated Footplate 17" W x 13"	N/C
------------------------------	--	-----


☐ ICP	Individual calf pads	\$200
------------------------------	----------------------------	-------

Select seat pan to footplate:

☐ SPTF1	9" - 13"	N/C
☐ SPTF14	13" - 19"	N/C

Set at:

☐ SA70	Set at 70°	N/C
☐ SA90	Set at 90°	N/C
☐ SA97	Set at 97°	N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.

FOOT PLATFORM OPTIONS FOR CENTER MOUNTS⁽¹⁾



☐ SFFP	Small Foot Platform (12"w x 8.5"d)	\$180
☐ MFFP	Medium Platform (12"w x 10.75"d)	\$180
☐ LFFP	Large Platform (14.5"w x 13"d)	\$250
☐ XLFFP	Extra Large Platform (17"w x 13"d)	\$250
☐ XXLFFP	Extra Extra Large Platform (21"w x 13"d)	\$250\

1. May need to adjust leg position to avoid interference with front casters

MANUAL ELEVATING CENTER MOUNT LEGREST

☐ MECML	Manual Elevating Center Mount Legrest With Plantar Flex.	\$595
--------------------------------	--	-------

Select legrest set up position:

☐ 90SU	Set at 90° (elevating range 90-40°)	STD
☐ 110TU	Tucked at 110° (elevating range 110-60°) ⁽¹⁾⁽²⁾	N/C

Rubber Coated Footplate:

☐ MRCP	Rubber coated footplate 11.5" W x 10" D	N/C
-------------------------------	---	-----

Select seat pan to footplate for Rubber Coated Footplate

☐ CSPTF1	11.5" - 15"	N/C
☐ CSPTF2	14" - 19"	N/C

Individual Footplates:

☐ IFSL	Small - (5" x 7.5" deep) - Left	\$165
☐ IFSR	Small - (5" x 7.5" deep) - Right	\$165
☐ IFML	Medium - (5.5" x 9.5" deep) - Left	\$165
☐ IFMR	Medium - (5.5" x 9.5" deep) - Right	\$165
☐ IFLL	Large - (6" x 11.5" deep) - Left	\$165
☐ IFLR	Large - (6" x 11.5" deep) - Right	\$165

Select seat pan to footplate for Individual Footplates

☐ CSPTF3	9.5" - 14.5"	N/C
☐ CSPTF4	12" - 19"	N/C

Optional:

☐ CICP	Individual calf pads	\$200
-------------------------------	----------------------------	-------

1. Not available on SP2 Narrow base.

2. Not available on seat depths 16", 17".

MAXX STYLE MANUAL ELEVATING LEGRESTS

- ☐ MMPPSA MAXX Swing Away Manual Elevating⁽¹⁾⁽²⁾ \$680
- Select Receivers:
- ☐ RS Straight Receivers N/C
- ☐ RF Flared 1" Outward Receivers N/C
- Select Seat Pan to Footplate:
- ☐ SPTF5 Length: 10"-14" N/C
- ☐ SPTF2 Length: 13"-17" N/C
- ☐ SPTF7 Length: 16"-20" N/C
- Options:
- ☐ FCPMS Flat Calf Pad Option \$275
- ☐ CCPMS Curved Calf Pad Option⁽³⁾ \$275





1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

SWING AWAY LEGREST

- ☐ HD70 Heavy Duty 70° Swing Away Footrest⁽¹⁾⁽²⁾ \$495
- Select Receivers:
- ☐ RS Straight Receivers N/C
- ☐ RF Flared 1" Outward Receivers N/C
- Select Seat Pan to Footplate:
- ☐ SPTF11 Length: 6"-10" N/C
- ☐ SPTF8 Length: 11"-15" N/C
- ☐ SPTF2 Length: 13"-17" N/C
- ☐ SPTF7 Length: 16"-20" N/C
- Options:
- ☐ HDFCP Flat Calf Pad Option \$275
- ☐ HD70CCP Curved Calf Pad Option⁽³⁾ \$275
- ☐ DCS Dynamic Coil Springs \$190
- ☐ I70SAF Invacare 70° Swing-away Footrests \$325
- ☐ I70-CF Comes with composite foot plates N/C
- ☐ 70HDAIH Invacare 70° Non-Tapered HD Swing-away Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards \$457
- ☐ 70TAPAIH Invacare Tapered 70° Swing-away Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards \$497
- ☐ ALPT Longer Pivot and Slide Tube for 70° Swingaways (Pair) \$132
- ☐ RLSL Residual Limb Support - Left⁽⁴⁾
6.75"W x 9.75"L x 2"H (each) \$495
- ☐ RLSR Residual Limb Support - Right⁽⁴⁾
6.75"W x 9.75"L x 2"H (each) \$495

1. Must choose footplate on page 8
2. Heavy duty durable construction (4.5lbs ea.).
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.
4. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.

FOOTPLATES FOR SWING AWAY/ELEVATING LEGREST⁽¹⁾



Angle Adjustable Footplates

- ☐ AFSL Left Small Footplate (4.25"W x 8"D) \$165
- ☐ AFSR Right Small Footplate (4.25"W x 8"D) \$165
- ☐ AFML Left Medium Footplate (5.25"W x 8"D) \$165
- ☐ AFMR Right Medium Footplate (5.25"W x 8"D) \$165
- ☐ AFLL Left Large Footplate (6.25"W x 8"D) \$165
- ☐ AFLR Right Large Footplate (6.25"W x 8"D) \$165
- ☐ AFXL Left Extra Large Footplate (9"W x 11"D) \$165
- ☐ AFXR Right Extra Large Footplate (9"W x 11"D) \$165

Multi-Axis Adjustable Angle Footplates

- ☐ MFSL Left Small Footplate (4.25"W x 8"D) \$225
- ☐ MFSR Right Small Footplate (4.25"W x 8"D) \$225
- ☐ MFML Left Medium Footplate (5.25"W x 8"D) \$225
- ☐ MFMR Right Medium Footplate (5.25"W x 8"D) \$225
- ☐ MFLL Left Large Footplate (6.25"W x 8"D) \$225
- ☐ MFLR Right Large Footplate (6.25"W x 8"D) \$225
- ☐ MFXL Left Extra Large Footplate (9"W x 11"D) \$225
- ☐ MFXR Right Extra Large Footplate (9"W x 11"D) \$225

HEEL LOOPS, CUPS AND TOE STRAPS

- ☐ HLL Left Heel Loop (each) \$20
- ☐ HLR Right Heel Loop (each) \$20
- ☐ HCL Left Cup Loop (each) \$20
- ☐ HCR Right Cup Loop (each) \$20
- ☐ TSL Left Toe Strap (each) \$20
- ☐ TSR Right Toe Strap (each) \$20

INSTALLATION

- ☒ MCI Motion Concepts Installation STD

PLEASE SPECIFY PAGES INCLUDED WITH ORDER

- ☐ Pages 1-5 Excludes Matrix
- ☐ Pages 1-8 Complete Order Form

Total Retail Value \$

SUBMIT

PRINT FORM

SAVE

CLEAR FORM

MATRIX SEAT CUSHIONS

Matrx Libra Seat Cushion \$745

☐ LC1616	☐ LC1618	☐ LC1620	☐ LC1716	☐ LC1718
☐ LC1720	☐ LC1816	☐ LC1818	☐ LC1820	☐ LC1916
☐ LC1918	☐ LC1920	☐ LC2016	☐ LC2018	☐ LC2020

Matrx Libra Seat Cushion \$1,025

☐ LC1622	☐ LC1722	☐ LC1822	☐ LC1922	☐ LC2022
---------------------------------	---------------------------------	---------------------------------	---------------------------------	---------------------------------

Matrx Libra HD \$1,025

☐ LC2118	☐ LC2120	☐ LC2122	☐ LC2218	☐ LC2220
☐ LC2320	☐ LC2322	☐ LC2418	☐ LC2420	☐ LC2422

Matrx PS Seat Cushion \$569

☐ PS1616	☐ PS1618	☐ PS1620	☐ PS1716	☐ PS1718
☐ PS1720	☐ PS1816	☐ PS1818	☐ PS1820	☐ PS2016
☐ PS2018	☐ PS2020			

Matrx PS HD Seat Cushion \$975

☐ PS2218	☐ PS2220	☐ PS2222	☐ PS2418	☐ PS2420
☐ PS2422	☐ PS2618	☐ PS2620	☐ PS2622	

Matrx PSP Seat Cushion \$645

☐ PSP1616	☐ PSP1618	☐ PSP1620	☐ PSP1716	☐ PSP1718
☐ PSP1720	☐ PSP1816	☐ PSP1818	☐ PSP1820	☐ PSP2016
☐ PSP2018	☐ PSP2020			

Matrx PSP HD Seat Cushion \$995

☐ PSP2218	☐ PSP2220	☐ PSP2222	☐ PSP2418	☐ PSP2420
☐ PSP2422	☐ PSP2618	☐ PSP2620	☐ PSP2622	

Matrx Multi Cushion with reversable cover \$745

☐ MM1616	☐ MM1618	☐ MM1620	☐ MM1716	☐ MM1718
☐ MM1720	☐ MM1816	☐ MM1818	☐ MM1820	☐ MM2016
☐ MM2018	☐ MM2020			

Matrx Multi Cushion with reversable cover \$995

☐ MM2216	☐ MM2218	☐ MM2220	☐ MM2222	
---------------------------------	---------------------------------	---------------------------------	---------------------------------	--

UPGRADABLE MATRIX BACK OPTIONS

Matrx Elite Back (3" Contour) \$766

☐ PBE1614	☐ PBE1616	☐ PBE1618	☐ PBE1620	☐ PBE1714
☐ PBE1716	☐ PBE1718	☐ PBE1720	☐ PBE1814	☐ PBE1816
☐ PBE1818	☐ PBE1820	☐ PBE1916	☐ PBE1918	☐ PBE1920
☐ PBE2014	☐ PBE2016	☐ PBE2018	☐ PBE2020	

Matrx Elite HD Back (3" Contour) \$1,095

☐ PBE2116-HD	☐ PBE2120-HD	☐ PBE2216-HD		
☐ PBE2220-HD		☐ PBE2316-HD	☐ PBE2320-HD	
☐ PBE2416-HD	☐ PBE2420-HD			

Matrx Elite HD Back (3" Contour) \$1,666

☐ EDB2616-HD	☐ EDB2620-HD			
-------------------------------------	-------------------------------------	--	--	--

Matrx Elite Deep Back (6" Contour) \$918

☐ EDB1614	☐ EDB1616	☐ EDB1618	☐ EDB1620	☐ EDB1714
☐ EDB1716	☐ EDB1718	☐ EDB1720	☐ EDB1814	☐ EDB1816
☐ EDB1818	☐ EDB1820	☐ EDB1914	☐ EDB1916	☐ EDB1918
☐ EDB1920	☐ EDB2014	☐ EDB2016	☐ EDB2018	☐ EDB2020

UPGRADABLE MATRIX BACK OPTIONS CONT.

Matrx Elite Deep HD Back (6" Contour) \$1,095

☐ EDB2116-HD	☐ EDB2120-HD	☐ EDB2216-HD
☐ EDB2220-HD	☐ EDB2316-HD	☐ EDB2320-HD
☐ EDB2416-HD	☐ EDB2420-HD	

Matrx Elite E2 Back (3" Contour) \$965

☐ E2S1614	☐ E2S1616	☐ E2S1714	☐ E2S1716	☐ E2S1814
☐ E2S1816	☐ E2S1914	☐ E2S1916	☐ E2S2014	☐ E2S2016

Matrx Elite E2 Back (3" Contour) \$995

☐ E2S1618	☐ E2S1620	☐ E2S1718	☐ E2S1720	☐ E2S1818
☐ E2S1820	☐ E2S1918	☐ E2S1920	☐ E2S2018	☐ E2S2020

Matrx Elite E2 Deep Back (6" Contour) \$1,125

☐ E2D1614	☐ E2D1616	☐ E2D1714	☐ E2D1716	☐ E2D1814
☐ E2D1816	☐ E2D1914	☐ E2D1916	☐ E2D2014	☐ E2D2016

Matrx Elite E2 Deep Back (6" Contour) \$1,155

☐ E2D1618	☐ E2D1620	☐ E2D1718	☐ E2D1720	☐ E2D1818
☐ E2D1820	☐ E2D1918	☐ E2D1920	☐ E2D2018	☐ E2D2020

Matrx Elite TR Back (3" Contour) \$866

☐ PBE1616-TR	☐ PBE1618-TR	☐ PBE1620-TR
☐ PBE1716-TR	☐ PBE1718-TR	☐ PBE1720-TR
☐ PBE1816-TR	☐ PBE1818-TR	☐ PBE1820-TR
☐ PBE1916-TR	☐ PBE1918-TR	☐ PBE1920-TR
☐ PBE2016-TR	☐ PBE2018-TR	☐ PBE2020-TR

Matrx Elite TR HD Back (3" Contour) \$1,151

☐ PBE2120-TRHD	☐ PBE2220-TRHD
☐ PBE2320-TRHD	☐ PBE2420-TRHD

Matrx MAC Back⁽¹⁾ \$1,495

☐ MAC1616	☐ MAC1618	☐ MAC1620	☐ MAC1622	☐ MAC1716
☐ MAC1718	☐ MAC1720	☐ MAC1722	☐ MAC1816	☐ MAC1818
☐ MAC1820	☐ MAC1822	☐ MAC1916	☐ MAC1918	☐ MAC1920
☐ MAC1922	☐ MAC2016	☐ MAC2018	☐ MAC2020	☐ MAC2022

1. Must select Pad Contour in Matrx MAC Back Pad Contour section below.

Matrx MAC Back w/ HUG⁽¹⁾ \$1,695

☐ MAC1616-HUG	☐ MAC1618-HUG	☐ MAC1620-HUG
☐ MAC1622-HUG	☐ MAC1716-HUG	☐ MAC1718-HUG
☐ MAC1720-HUG	☐ MAC1722-HUG	☐ MAC1816-HUG
☐ MAC1818-HUG	☐ MAC1820-HUG	☐ MAC1822-HUG
☐ MAC1916-HUG	☐ MAC1918-HUG	☐ MAC1920-HUG
☐ MAC1922-HUG	☐ MAC2016-HUG	☐ MAC2018-HUG
☐ MAC2020-HUG	☐ MAC2022-HUG	

1. Must select Pad Contour in Matrx MAC Back Pad Contour section below.

Matrx MAC Back Pad Contour

☐ MLCSL	Left 3.5"	☐ MLCSR	Right 3.5"
☐ MLCDL	Left 6"	☐ MLCXR	Right 6"
☐ MLCXDL	Left 9"	☐ MLCXDR	Right 9"

MATRX ELAN HEADREST PAD⁽¹⁾

☐ ESP6	Matrx Elan Standard Pad 6"	\$310
☐ ESP10	Matrx Elan Standard Pad 10"	\$310
☐ ESP14	Matrx Elan Standard Pad 14"	\$310
☐ ESP6-HUG	Matrx Elan Hug Standard Pad 6"	\$310
☐ ESP10-HUG	Matrx Elan Hug Standard Pad 10"	\$310
☐ ESP14-HUG	Matrx Elan Hug Standard Pad 14"	\$310

1. If mounting on PB Elite must order removable hardware

HEADREST MOUNTING HARDWARE OPTIONS

☐ LXSHW	Loxx Standard Hardware 15"	\$275
☐ MEHWRE	Elan Hardware 15"	\$285
☐ FPLXSHW	FLIPP Loxx Standard Hardware 15"	\$705
☐ FPLXSHW-EX	FLIPP Loxx Extended Hardware.....	\$705
☐ FPMEHW	FLIPP Elan Hardware 15"	\$715

MATRX FIXED LATERALS - FOR ELITE / ELITE TR BACKS

☐ PBSSL	Matrx Small Lateral Fixed-Left(3x5.5").....	\$113
☐ PBSSR	Matrx Small Lateral Fixed-Right(3x5.5")	\$113
☐ PBLSL	Matrx Standard Lateral Fixed-Left(7.5x5")	\$116
☐ PBLSR	Matrx Standard Lateral Fixed-Right (7.5x5")	\$116
☐ PBOLL	Matrx Offset Lateral Fixed Left	\$125
☐ PBOLR	Matrx Offset Lateral Fixed Right	\$125

1. Not available on the MAC backrests.

MATRX SWING AWAY LATERAL TRUNK SUPPORTS⁽¹⁾

Matrx Lateral Trunk Support Left

☐ EPADXL	Size XS (3.25" x 4.25")	\$115
☐ EPADSL	Size S (3.5" x 5.25")	\$115
☐ EPADML	Size M (3.75" x 5.75")	\$115
☐ EPADLL	Size L (4.25" x 6.5")	\$115

Must select one option:

☐ SAHEBL	Matrx Elite Swing Away Hardware Left	\$250
☐ OSAHL	Matrx Offset Elite Swing Away Hardware Left.....	\$250

Matrx Lateral Trunk Support Right

☐ EPADXR	Size XS (3.25" x 4.25")	\$115
☐ EPADSR	Size S (3.5" x 5.25")	\$115
☐ EPADMR	Size M (3.75" x 5.75")	\$115
☐ EPADLR	Size L (4.25" x 6.5")	\$115

Must select one option:

☐ SAHEBR	Matrx Elite Swing Away Hardware Right	\$250
☐ OSAHR	Matrx Offset Elite Swing Away Hardware Right	\$250

1. Not available on the MAC backrests.

HIP AND KNEE SUPPORTS

Lateral Hip Support Left⁽¹⁾

☐ HSPXLH	Size XS (5.5" W x 4" H)	\$120
☐ HSPSMLH	Size S/M (7" W x 4" H)	\$120
☐ HSPSLH	Size S (7"W x 6"H)	\$120
☐ HSPMLH	Size M (9" W x 4"H)	\$120
☐ HSPLLH	Size L (13" W x 4"H)	\$120

Lateral Hip Support Right⁽¹⁾

☐ HSPXRH	Size XS (5.5"W x 4"H)	\$120
☐ HSPSMRH	Size S/M (7"W x 4"H)	\$120
☐ HSPSRH	Size S (7"W x 6"H)	\$120
☐ HSPMRH	Size M (9"W x 4"H)	\$120
☐ HSPLRH	Size L (13"W x 4"H)	\$120

Lateral Knee Support Left⁽¹⁾

☐ HSPXLK	Size XS (5.5"W x 4" H)	\$120
☐ HSPSMLK	Size S/M (7"W x 4" H)	\$120
☐ HSPSLK	Size S (7"W x 6" H).....	\$120
☐ HSPMLK	Size M (9"W x 4"H)	\$120
☐ HSPLLK	Size L (13"W x 4"H)	\$120

Lateral Knee Support Right⁽¹⁾

☐ HSPXRK	Size XS (5.5"W x 4"H)	\$120
☐ HSPSMRK	Size S/M (7"W x 4"H)	\$120
☐ HSPSRK	Size S (7"W x 6"H)	\$120
☐ HSPMRK	Size M (9"W x 4"H)	\$120
☐ HSPLRK	Size L (13"W x 4"H)	\$120

Add Removable Hardware:

MAXX Style Quick Release, Multi-Axis Mounting Hardware for Hip Supports

☐ MHSHLH	Left mounted at Hip	\$215
☐ MHSHLK	Left mounted at Knee	\$215
☐ MHSHRH	Right mounted at Hip	\$215
☐ MHSHRK	Right mounted at Knee	\$215

Extended Maxx Style Quick Release, Multi-Axis Mounting Hardware for Hip Supports⁽²⁾

☐ EMQHLH	Left mounted at Hip	\$215
☐ EMQHLK	Left mounted at Knee	\$215
☐ EMQHRH	Right mounted at Hip	\$215
☐ EMQHRK	Right mounted at Knee	\$215

1. Comes standard with Fixed Mount Hardware.

2. Recomend when finished seat cushion thickness will be greater then 3.5" high.

NOTES

Select Pad Size and Both Extensions:

- ### Right Elbow Block⁽¹⁾

- ### Select Pad Size and Both Extensions:

- | | | | |
|---|------|-------------------------------|-------|
|  | EBX | Size XS (3.25" x 4.25") | \$195 |
|  | EBSR | Size S (3.5" x 5.25") | \$195 |
|  | EBMR | Size M (3.75" x 5.75") | \$195 |
|  | EBLR | Size L (4.25" x 6.5") | \$195 |
|  | V3 | Right Verticle Extension 3" | |
|  | V5 | Right Verticle Extension 5" | |
|  | H5 | Right Horizontal Extension 5" | |
|  | H7 | Right Horizontal Extension 7" | |

1. Prevents arm from slipping.

Total Retail Value \$

SUBMIT

PRINT FORM

SAVE

CLEAR FORM