



## Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



### Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

**SAVE**

### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

**SUBMIT**

**CLEAR FORM**

### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**888-433-6818**

**PRINT FORM**



**TURNING,  
DRIVING  
& RIDE  
COMFORT!**

# ULTRA LOW MAXX CG POWER POSITIONING SYSTEM ON STORM TORQUE RX20 BASE

CDN PRICE LIST AND ORDER FORM - Price Effective May 15, 2026



**PRINT FORM**

**SAVE**

**CLEAR FORM**

Motion Concepts Customer Service: 1-866-748-7943 quotescdn@motionconcepts.com

## REQUIRED INFORMATION

☐ Quote ☐ Order

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_ Tag \_\_\_\_\_

Purchasing Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_

RTS/Therapist: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Province: \_\_\_\_\_


Postal Code: \_\_\_\_\_ Client Gender: ☐ M ☐ F

Special Client Conditions: \_\_\_\_\_

## REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all **REQUIRED\*** measurements

|                               |  |
|-------------------------------|--|
| <b>Client Height*</b>         |  |
| <b>Client Weight (lbs)*</b>   |  |
| <b>A. Seat to Shoulder</b>    |  |
| <b>B. Trunk Depth</b>         |  |
| <b>C. Chest Width</b>         |  |
| <b>D. Knee to Back</b>        |  |
| <b>E. Seat to Top of Head</b> |  |
| <b>F. Back to Joystick</b>    |  |
| <b>G. Seat to Elbow</b>       |  |
| <b>H. Hip Width</b>           |  |
| <b>I. Knee to Heel*</b>       |  |
| <b>Cushion Thickness*</b>     |  |



## BASE OPTIONS

**NEW!**

☐ SRX-20 Storm **Torque** RX20 Base<sup>(1)</sup> ..... \$7,575

- Enhanced Front Caster Housing and Bearings
- New Front Caster Design
- Improved double-sided Front Forks
- Upgraded Rear Suspension
- Optimized Linx Parameters

1. Comes Standard with Lighting Package

## BASE WIDTH

☐ NARROWFRM Narrow Base 22.4" (MK 50 Batteries) ..... N/C

☐ STDFRM Standard Base 24.4" (MK 24 FT Batteries) ..... N/C

## SEAT-TO-FLOOR HEIGHT<sup>(1)</sup>

☐ SLOW 16.75" Seat-To-Floor<sup>(2)</sup> ..... \$250

☐ LOW 17.25" Seat-To-Floor ..... N/C

☐ MED 18.25" Seat-To-Floor ..... N/C

☐ TALL 19.25" Seat-To-Floor ..... N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"

2. This seat to floor height is only available with the Group 24 Battery Tray (24TRY) and on 17" seat depths and greater.

## USER WEIGHT LIMIT

☐ UW250 Weight Capacity up to 250lbs<sup>(1)</sup> ..... STD

☐ U300 Weight Capacity up to 300lbs ..... STD

1. Required for Elevate Systems and/or NARROWFRM

## WHEELCHAIR TRANSPORT BRACKETS

☐ TRBKTS Wheelchair Transport Brackets<sup>(1)</sup> ..... STD

☐ TRRO Transport Ready Option<sup>(2)(4)</sup> ..... \$499

☐ TRRO-E Transport Ready Option for Elevate<sup>(2)(3)(4)</sup> ..... \$499

1. For unoccupied use only

2. TRRO meets ISO-7176 Part 19

## WHEELCHAIR TRANSPORT BRACKETS

3. Required for Elevate Systems to meet ISO-7176 Part 19, includes two additional tie-down loops on rail.

4. Must select a Motion Concepts or Matrix Elan headrest with hardware and Motion Concepts Standard Rehab, High Back or Matrix Upgradeable back. Comes with pelvic safety belt standard, please select TROBELT on page 5.

## TIRE OPTIONS

☐ AB1431-3 14" x 3" Black Tire with Gel Foam Inserts ..... N/C

☐ AB1430-3 14" x 3" Black Tire - Pneumatic Tire ..... N/C

## CASTER OPTION

**NEW!**

☐ AB1228-3 8" x 2.5" Black Ribbed Casters with Solid Inlay .. STD

## BATTERY CHARGER OPTION

☐ 110CHARGER 110 Volt Battery Charger ..... N/C

☐ OMIT Omit Battery Charger ..... N/C

## BATTERIES

☐ M24 SLD G FT MK Group 24 FT Gel Battery (Pair) ..... \$840

☐ M50-12 SLD M MK Group 50 Battery (Pair) ..... \$650

## HIGHSPEED MOTORS AND ENHANCED DRIVING

☐ HSMTR Performance Motors 7.5 mph<sup>(1)(2)(3)</sup> ..... \$995

☐ GTRAC-LX LINX G-TRAC Module<sup>(4)</sup> ..... \$1,200

1. Not available on the NARROWFRM.

2. Extended high speed use on rugged terrain is not recommended.

3. Performance/High Speed Motors may have less torque and may have less range than standard speed motors.

4. Recommended for use with Specialty Drive Controls

# Ultra Low Maxx CG Power Positioning System

## Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### FRAME FINISH

|                               |                                                                                   |                     |     |
|-------------------------------|-----------------------------------------------------------------------------------|---------------------|-----|
| <input type="checkbox"/> 170P |  | Obsidian Black..... | N/C |
| <input type="checkbox"/> 171P |  | Ruby Red .....      | N/C |
| <input type="checkbox"/> 172P |  | Peacock Blue .....  | N/C |
| <input type="checkbox"/> 173P |  | Pink Flamingo ..... | N/C |
| <input type="checkbox"/> 174P |  | Purple Passion..... | N/C |
| <input type="checkbox"/> 175P |  | Denim Navy.....     | N/C |
| <input type="checkbox"/> 176P |  | Lime Twist .....    | N/C |
| <input type="checkbox"/> 177P |  | Arctic Blue .....   | N/C |

### RIM INSERT COLORS - MUST PICK ONE

|                                |                                                                                     |                       |     |
|--------------------------------|-------------------------------------------------------------------------------------|-----------------------|-----|
| <input type="checkbox"/> 170PR |    | Obsidian Black.....   | N/C |
| <input type="checkbox"/> 171PR |    | Ruby Red .....        | N/C |
| <input type="checkbox"/> 172PR |    | Peacock Blue .....    | N/C |
| <input type="checkbox"/> 173PR |    | Pink Flamingo .....   | N/C |
| <input type="checkbox"/> 174PR |    | Purple Passion.....   | N/C |
| <input type="checkbox"/> 175PR |    | Denim Navy.....       | N/C |
| <input type="checkbox"/> 176PR |    | Lime Twist .....      | N/C |
| <input type="checkbox"/> 177PR |    | Arctic Blue .....     | N/C |
| <input type="checkbox"/> 178PR |    | Banana Split .....    | N/C |
| <input type="checkbox"/> 179PR |   | Orange Marmalade..... | N/C |
| <input type="checkbox"/> 180PR |  | Morning Mist .....    | N/C |

### NON-EXPANDABLE SINGLE/DUAL JOYSTICK

|                                                                                |                                             |     |
|--------------------------------------------------------------------------------|---------------------------------------------|-----|
| <input type="checkbox"/> REM215                                                | LED Non-Expandable Remote with Lights ..... | N/C |
| <i>Three drive profiles, control of two seat functions and light operation</i> |                                             |     |

### EXPANDABLE JOYSTICKS

|                                                                                                                        |                                                                    |         |
|------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------|---------|
| <input type="checkbox"/> REM216                                                                                        | LED Expandable Remote/Joystick<br>with Lights <sup>(2)</sup> ..... | \$875   |
| <i>Three drive profiles, control of up to 6 actuators<br/>and light operation. Not for use with specialty controls</i> |                                                                    |         |
| <input type="checkbox"/> REM400                                                                                        | Color 3.5" Touch Screen<br>Remote/Joystick <sup>(1)</sup> .....    | \$1,100 |
| <i>Four drive profiles, control of up to 6 actuators through the<br/>drive control and light operation.</i>            |                                                                    |         |
| <input type="checkbox"/> SW400TGL                                                                                      | Toggle Kit for REM400 <sup>(2)</sup> .....                         | \$200   |

1. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.  
2. Adds 2 toggle switches to REM400.

### EXPANDABLE SPECIALTY DISPLAY

|                                 |                                                                 |         |
|---------------------------------|-----------------------------------------------------------------|---------|
| <input type="checkbox"/> REM500 | LiNX Color 3.5" Touch Screen<br>Display <sup>(1)(2)</sup> ..... | \$1,100 |
|---------------------------------|-----------------------------------------------------------------|---------|

Select mounting position for REM500

|                                |                     |       |
|--------------------------------|---------------------|-------|
| <input type="checkbox"/> R5MR  | Right.....          | N/C   |
| <input type="checkbox"/> R5ML  | Left .....          | N/C   |
| <input type="checkbox"/> R5RPB | Rear Push Bar ..... | \$140 |

1. Driver Control not included.  
2. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.

### SPECIALTY PROPORTIONAL CONTROLS<sup>(1)(2)</sup>

|                                  |                                     |         |
|----------------------------------|-------------------------------------|---------|
| <input type="checkbox"/> CREM    | LiNX Compact Remote .....           | \$1,097 |
| <input type="checkbox"/> CREM-LF | LiNX Compact Remote Low Force ..... | \$1,097 |

1. Requires REM400 or REM500

2. Select Mount in Mounting For Driver Control section  
controls purchased separately.

### ATTENDANT DRIVE CONTROL CHOICES

|                              |                                                     |         |
|------------------------------|-----------------------------------------------------|---------|
| <input type="checkbox"/> ACU | Proportional Attendant Control <sup>(1)</sup> ..... | \$1,024 |
|------------------------------|-----------------------------------------------------|---------|

Select mounting position for ACU

|                                |        |
|--------------------------------|--------|
| <input type="checkbox"/> RPBMR | Right  |
| <input type="checkbox"/> RPBML | Left   |
| <input type="checkbox"/> RPBMC | Center |

### JOYSTICK MOUNT

#### Arm Mounted

|                               |                                                                     |       |
|-------------------------------|---------------------------------------------------------------------|-------|
| <input type="checkbox"/> FR   | Fixed Joystick Mounting - Right .....                               | STD   |
| <input type="checkbox"/> FL   | Fixed Joystick Mounting - Left .....                                | STD   |
| <input type="checkbox"/> MRR  | Maxx Resolve Swing-Away Joystick Mount<br>(Right Hand) .....        | \$450 |
| <input type="checkbox"/> MRL  | Maxx Resolve Swing-Away Joystick Mount<br>(Left Hand) .....         | \$450 |
| <input type="checkbox"/> SAQR | Motion Height Adjustable Swing-Away<br>Quad Link (Right Hand) ..... | \$244 |
| <input type="checkbox"/> SAQL | Motion Height Adjustable Swing-Away<br>Quad Link (Left Hand) .....  | \$244 |

#### Rear Push Bar Mounted

|                                |              |       |
|--------------------------------|--------------|-------|
| <input type="checkbox"/> RPBMR | Right .....  | \$140 |
| <input type="checkbox"/> RPBML | Left .....   | \$140 |
| <input type="checkbox"/> RPBMC | Center ..... | \$140 |

### JOYSTICK KNOB - OPTIONAL

|                                 |                                 |       |
|---------------------------------|---------------------------------|-------|
| <input type="checkbox"/> PC101A | BodyPoint 3" U-Shaped.....      | \$120 |
| <input type="checkbox"/> PC102A | BodyPoint 4" U-Shaped .....     | \$120 |
| <input type="checkbox"/> PC104A | BodyPoint Straight Handle ..... | \$120 |
| <input type="checkbox"/> PC107A | BodyPoint Domed Rubber .....    | \$120 |
| <input type="checkbox"/> PC110A | BodyPoint Grooved Mushroom..... | \$120 |

1. To remove the joystick knob from the REM400, pull straight up on the joystick  
(DO NOT TWIST) otherwise damage may occur and may void the warranty.

# Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number:\_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

## ELECTRONICS - EXPANDABLE MODULES

- |                                                                                          |                                                                                         |         |
|------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|---------|
|  OUT500   | Output Module .....                                                                     | \$645   |
|  IN200    | LiNX Switch Input Module with<br>Power ON/OFF for REM200 Series <sup>(1)(2)</sup> ..... | \$645   |
|  IN500    | LiNX 9 Pin Input Module with<br>Sip-N-Puff Nozzle <sup>(2)(3)(4)(6)</sup> .....         | \$1,800 |
|  PKG32666 | Therafin Sip-N-Puff Breath Tube Kit .....                                               | \$440   |
|  COVER-L  | Large Recline Shroud Cover <sup>(5)</sup> .....                                         | N/C     |

1. Not for use with Alternative Drive Controls. (Uses switch ports).
2. Only one input module per chair.
3. Input module also required for any 9-pin specialty control.
4. 9 pin specialty control and Sip-N-Puff cannot both be used together.
5. May be required when choosing OUT500, IN200 or IN500 options.
6. Requires REM400 or REM500 for compatibility.

## LINX ELECTRONIC ACCESSORIES

- |                                                                                       |                                      |      |
|---------------------------------------------------------------------------------------|--------------------------------------|------|
|  LAK | LiNX Access Key <sup>(1)</sup> ..... | \$50 |
|---------------------------------------------------------------------------------------|--------------------------------------|------|

1. Required for programming

## NOTES

[illegible]

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

#### POWER POSITIONING SYSTEMS<sup>(1)</sup>

|                                 |                                                                                                    |          |
|---------------------------------|----------------------------------------------------------------------------------------------------|----------|
| <input type="checkbox"/> ULMT   | Ultra Low Maxx CG Tilt (50°) Only System <sup>(1)</sup>                                            | \$6,969  |
| <input type="checkbox"/> ULMTRE | Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Extended Shear Reduction(ESR) <sup>(1)</sup> | \$11,569 |
| <input type="checkbox"/> ULMTRP | Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Power Sliding Back(PSB) <sup>(1)(2)</sup>    | \$12,069 |
| <input type="checkbox"/> ULMRE  | Ultra Low Maxx Power Recline (170°) w/Extended Shear Reduction(ESR)                                | \$10,350 |
| <input type="checkbox"/> ULMRP  | Ultra Low Maxx Power Recline (170°) w/Power Sliding Back(PSB) <sup>(2)</sup>                       | \$10,850 |
| <input type="checkbox"/> ULME   | Ultra Low Maxx 12 " Elevating Lift Seat Only <sup>(3)</sup>                                        | \$4,995  |
|                                 | Weight capacity up to 250lbs                                                                       |          |

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module
2. Requires Multi-Function M610 through Drive Control (MFCBD).
3. Not available on seat sizes greater than 20" wide by 20" deep.

#### ADD: POWER ADJUSTABLE SEAT HEIGHT MODULE

|                                |                                                              |         |
|--------------------------------|--------------------------------------------------------------|---------|
| <input type="checkbox"/> ESM12 | 12" Power Adjustable Seat Height Module <sup>(1)(2)(3)</sup> | \$2,495 |
|                                | Weight capacity up to 250lbs                                 |         |

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.
2. Not available on seat sizes greater than 20" wide by 20" deep
3. For weights greater than 175lbs tilt while elevated is not available.

#### MODULES - MUST SELECT SYSTEM TO ADD MODULES

|                                   |                                                                     |       |
|-----------------------------------|---------------------------------------------------------------------|-------|
| <input type="checkbox"/> EPRETILT | 5° Electronic Pre-Tilt <sup>(1)</sup>                               | N/C   |
| <input type="checkbox"/> PRM8     | 8° Precline Module <sup>(2)</sup>                                   | N/C   |
| <input type="checkbox"/> PRM30    | 30° Precline Module <sup>(2)</sup>                                  | \$500 |
| <input type="checkbox"/> ERSSP    | Programable Electronic Recline Function Start & Stop <sup>(3)</sup> | \$400 |

Specify:

Start: \_\_\_\_\_degrees Stop: \_\_\_\_\_degrees

1. This Pre-tilt option can add approximately 2" (5°) to the front seat to floor height. Multi-Function M610 through Drive Control (MFCBD).
2. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available. Must order Recline System
3. Requires Multi-Function M610 through Drive Control (MFCBD). The range is dependent on the Recline module selected. Will be set to factory settings unless specified.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### ADULT SEAT WIDTH SETTINGS

|                                |                                                         |       |
|--------------------------------|---------------------------------------------------------|-------|
| <input type="checkbox"/> SW116 | 16" Wide (Adjustable from 16" to 20")                   | N/C   |
| <input type="checkbox"/> SW117 | 17" Wide (Adjustable from 16" to 20")                   | N/C   |
| <input type="checkbox"/> SW118 | 18" Wide (Adjustable from 16" to 20")                   | N/C   |
| <input type="checkbox"/> SW119 | 19" Wide (Adjustable from 16" to 20")                   | N/C   |
| <input type="checkbox"/> SW120 | 20" Wide (Adjustable from 16" to 20")                   | N/C   |
| <input type="checkbox"/> SW219 | 19" Wide (Adjustable from 19" to 22") <sup>(1)(2)</sup> | \$460 |
| <input type="checkbox"/> SW220 | 20" Wide (Adjustable from 19" to 22") <sup>(1)(2)</sup> | \$460 |
| <input type="checkbox"/> SW221 | 21" Wide (Adjustable from 19" to 22") <sup>(1)(2)</sup> | \$460 |
| <input type="checkbox"/> SW222 | 22" Wide (Adjustable from 19" to 22") <sup>(1)(2)</sup> | \$460 |

1. This seat size is not available on the NARROWFRM.
2. This seat range is not available with ESM12

#### ADULT ULTRA RAIL SEAT DEPTH SETTINGS

|                                |                                                         |       |
|--------------------------------|---------------------------------------------------------|-------|
| <input type="checkbox"/> SD116 | 16" Deep (Adjustable from 16" to 21")                   | N/C   |
| <input type="checkbox"/> SD117 | 17" Deep (Adjustable from 16" to 21")                   | N/C   |
| <input type="checkbox"/> SD118 | 18" Deep (Adjustable from 16" to 21")                   | N/C   |
| <input type="checkbox"/> SD119 | 19" Deep (Adjustable from 16" to 21")                   | N/C   |
| <input type="checkbox"/> SD120 | 20" Deep (Adjustable from 16" to 21")                   | N/C   |
| <input type="checkbox"/> SD121 | 21" Deep (Adjustable from 16" to 21")                   | N/C   |
| <input type="checkbox"/> SD219 | 19" Deep (Adjustable from 19" to 23") <sup>(1)(2)</sup> | \$460 |
| <input type="checkbox"/> SD221 | 21" Deep (Adjustable from 19" to 23") <sup>(1)(2)</sup> | \$460 |
| <input type="checkbox"/> SD222 | 22" Deep (Adjustable from 19" to 23") <sup>(1)(2)</sup> | \$460 |
| <input type="checkbox"/> SD223 | 23" Deep (Adjustable from 19" to 23") <sup>(1)(2)</sup> | \$460 |

1. This seat size is not available on the NARROWFRM.
2. This seat range is not available with ESM12

#### SMALL ADULT SEAT WIDTH SETTINGS<sup>(1)</sup>

|                                |                                       |       |
|--------------------------------|---------------------------------------|-------|
| <input type="checkbox"/> SW314 | 14" Wide (Adjustable from 14" to 15") | \$250 |
| <input type="checkbox"/> SW315 | 15" Wide (Adjustable from 14" to 15") | \$250 |

1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

#### SMALL ADULT SEAT DEPTH SETTINGS<sup>(1)</sup>

|                                |                                       |       |
|--------------------------------|---------------------------------------|-------|
| <input type="checkbox"/> SD314 | 14" Deep (Adjustable from 14" to 18") | \$250 |
| <input type="checkbox"/> SD315 | 15" Deep (Adjustable from 14" to 18") | \$250 |
| <input type="checkbox"/> SD316 | 16" Deep (Adjustable from 14" to 18") | \$250 |
| <input type="checkbox"/> SD317 | 17" Deep (Adjustable from 14" to 18") | \$250 |
| <input type="checkbox"/> SD318 | 18" Deep (Adjustable from 14" to 18") | \$250 |

1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### MATRIX SEAT CUSHION OPTIONS

##### Matrix PS Seat Cushion 14" - 20" wide/reversable cover

..... \$569

|                                 |                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> PS1414 | <input type="checkbox"/> PS1416 | <input type="checkbox"/> PS1616 | <input type="checkbox"/> PS1618 | <input type="checkbox"/> PS1620 |
| <input type="checkbox"/> PS1716 | <input type="checkbox"/> PS1718 | <input type="checkbox"/> PS1720 | <input type="checkbox"/> PS1816 | <input type="checkbox"/> PS1818 |
| <input type="checkbox"/> PS1820 | <input type="checkbox"/> PS2016 | <input type="checkbox"/> PS2018 | <input type="checkbox"/> PS2020 |                                 |

##### Matrix PS HD Seat Cushion 21"-22" wide

..... \$975

|                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> PS2218 | <input type="checkbox"/> PS2220 | <input type="checkbox"/> PS2222 |
|---------------------------------|---------------------------------|---------------------------------|

##### Matrix PSP Seat Cushion 14" - 20" wide/reversable cover

..... \$645

|                                  |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PSP1414 | <input type="checkbox"/> PSP1416 | <input type="checkbox"/> PSP1616 | <input type="checkbox"/> PSP1618 | <input type="checkbox"/> PSP1620 |
| <input type="checkbox"/> PSP1716 | <input type="checkbox"/> PSP1718 | <input type="checkbox"/> PSP1720 | <input type="checkbox"/> PSP1816 | <input type="checkbox"/> PSP1818 |
| <input type="checkbox"/> PSP1820 | <input type="checkbox"/> PSP2016 | <input type="checkbox"/> PSP2018 | <input type="checkbox"/> PSP2020 |                                  |

##### Matrix PSP HD Seat Cushion 21"-22" wide

..... \$995

|                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PSP2218 | <input type="checkbox"/> PSP2220 | <input type="checkbox"/> PSP2222 |
|----------------------------------|----------------------------------|----------------------------------|

##### Matrix Vi Seat Cushion 14" - 20" wide

..... \$624

|                                    |                                    |                                    |                                    |                                    |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| <input type="checkbox"/> KA1414-VI | <input type="checkbox"/> MA1415-VI | <input type="checkbox"/> MA1417-VI | <input type="checkbox"/> MA1616-VI | <input type="checkbox"/> MA1618-VI |
| <input type="checkbox"/> MA1620-VI | <input type="checkbox"/> MA1716-VI | <input type="checkbox"/> MA1718-VI | <input type="checkbox"/> MA1720-VI | <input type="checkbox"/> MA1816-VI |
| <input type="checkbox"/> MA1818-VI | <input type="checkbox"/> MA1820-VI | <input type="checkbox"/> MA2016-VI | <input type="checkbox"/> MA2018-VI |                                    |
| <input type="checkbox"/> MA2020-VI |                                    |                                    |                                    |                                    |

##### Matrix Vi Seat Cushion 21" - 22" wide

..... \$681

|                                    |                                    |
|------------------------------------|------------------------------------|
| <input type="checkbox"/> MA2218-VI | <input type="checkbox"/> MA2220-VI |
|------------------------------------|------------------------------------|

##### Matrix Vi HD Seat Cushion 22" wide

..... \$840

|                                    |
|------------------------------------|
| <input type="checkbox"/> MA2222-VI |
|------------------------------------|

##### Matrix Libra Seat Cushion 14" - 20" wide

..... \$745

|                                 |                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> LC1414 | <input type="checkbox"/> LC1416 | <input type="checkbox"/> LC1516 | <input type="checkbox"/> LC1518 | <input type="checkbox"/> LC1520 |
| <input type="checkbox"/> LC1616 | <input type="checkbox"/> LC1618 | <input type="checkbox"/> LC1620 | <input type="checkbox"/> LC1716 | <input type="checkbox"/> LC1718 |
| <input type="checkbox"/> LC1720 | <input type="checkbox"/> LC1816 | <input type="checkbox"/> LC1818 | <input type="checkbox"/> LC1820 | <input type="checkbox"/> LC1916 |
| <input type="checkbox"/> LC1918 | <input type="checkbox"/> LC1920 | <input type="checkbox"/> LC2016 | <input type="checkbox"/> LC2018 | <input type="checkbox"/> LC2020 |

##### Matrix Libra HD Seat Cushion 21" wide

..... \$1,025

|                                 |                                 |
|---------------------------------|---------------------------------|
| <input type="checkbox"/> LC2118 | <input type="checkbox"/> LC2120 |
|---------------------------------|---------------------------------|

##### Matrix Libra Seat Cushion 22" deep

..... \$1,025

|                                 |                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> LC1622 | <input type="checkbox"/> LC1722 | <input type="checkbox"/> LC1822 | <input type="checkbox"/> LC1922 | <input type="checkbox"/> LC2022 |
| <input type="checkbox"/> LC2122 |                                 |                                 |                                 |                                 |

##### Matrix Libra HD Seat Cushion 22" wide

..... \$1,025

|                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> LC2218 | <input type="checkbox"/> LC2220 | <input type="checkbox"/> LC2222 |
|---------------------------------|---------------------------------|---------------------------------|

##### Matrix Libra Adjust Seat Cushion 14" -20" wide

..... \$985

|                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> LC1414-ADJ | <input type="checkbox"/> LC1416-ADJ | <input type="checkbox"/> LC1516-ADJ | <input type="checkbox"/> LC1518-ADJ |
| <input type="checkbox"/> LC1520-ADJ | <input type="checkbox"/> LC1616-ADJ | <input type="checkbox"/> LC1618-ADJ | <input type="checkbox"/> LC1620-ADJ |
| <input type="checkbox"/> LC1716-ADJ | <input type="checkbox"/> LC1718-ADJ | <input type="checkbox"/> LC1720-ADJ | <input type="checkbox"/> LC1816-ADJ |
| <input type="checkbox"/> LC1818-ADJ | <input type="checkbox"/> LC1820-ADJ | <input type="checkbox"/> LC1916-ADJ | <input type="checkbox"/> LC1918-ADJ |
| <input type="checkbox"/> LC1920-ADJ | <input type="checkbox"/> LC2016-ADJ | <input type="checkbox"/> LC2018-ADJ | <input type="checkbox"/> LC2020-ADJ |

##### Matrix Libra Adjust Seat Cushion 22" deep

..... \$1,265

|                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> LC1622-ADJ | <input type="checkbox"/> LC1722-ADJ | <input type="checkbox"/> LC1822-ADJ | <input type="checkbox"/> LC1922-ADJ |
| <input type="checkbox"/> LC2022-ADJ | <input type="checkbox"/> LC2122-ADJ |                                     |                                     |

##### Matrix Libra Adjust HD Seat Cushion 22" wide

..... \$1,265

|                                     |                                     |
|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> LC2118-ADJ | <input type="checkbox"/> LC2120-ADJ |
|-------------------------------------|-------------------------------------|

##### Matrix Multi Cushion 14" - 20" wide w/reversable cover

..... \$745

|                                 |                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> MM1414 | <input type="checkbox"/> MM1416 | <input type="checkbox"/> MM1616 | <input type="checkbox"/> MM1618 | <input type="checkbox"/> MM1620 |
| <input type="checkbox"/> MM1716 | <input type="checkbox"/> MM1718 | <input type="checkbox"/> MM1720 | <input type="checkbox"/> MM1816 | <input type="checkbox"/> MM1818 |
| <input type="checkbox"/> MM1820 | <input type="checkbox"/> MM2016 | <input type="checkbox"/> MM2018 | <input type="checkbox"/> MM2020 |                                 |

##### Matrix Multi Cushion 22" wide w/reversable cover

..... \$995

|                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> MM2216 | <input type="checkbox"/> MM2218 | <input type="checkbox"/> MM2220 | <input type="checkbox"/> MM2222 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|

#### FABRIC OPTIONS FOR CUSHIONS

| Cushion      | Meshtex | Startex        | Infection Control | Spacetex | Onyx  | Fusion |
|--------------|---------|----------------|-------------------|----------|-------|--------|
| Libra/Adjust | N/A     | STD            | \$50              | N/A      | N/A   | N/A    |
| PS/PSP       | N/A     | STD Reversible | STD Reversible    | N/A      | N/A   | N/A    |
| Vi           | STD     | \$50           | \$50              | N/A      | \$100 | \$150  |

##### ☐ LC-SP

Spacetex Fabric Option for  
Libra Cushion ..... N/C

##### ☐ LC-IC

Infection Control Fabric Upgrade for  
Libra Cushion ..... \$50

##### ☐ FSU-MA

Fusion Fabric Upgrade for Matrix Vi Cushion... \$150

##### ☐ OCU-MA

Onyx Fabric Upgrade for Matrix Vi Cushion .... \$100

##### ☐ STARTEX-S

Smooth side Startex Fabric Upgrade for  
Matrix Vi Cushion ..... \$50

##### ☐ STARTEX-F

Fabric side Startex Fabric Upgrade for  
Matrix Vi Cushion ..... \$50

#### LAP BELTS

##### ☐ BELT48

Push Button Style Lap Belt, Length 48"<sup>(1)</sup> ..... STD

##### ☐ BELT60

Push Button Style Lap Belt, Length 60"<sup>(1)</sup> ..... STD

##### ☐ BELT71

Push Button Style Lap Belt, Length 71"<sup>(1)</sup> ..... STD

##### ☐ PLB2

Bodypoint Padded Lap Belt 2 Point<sup>(2)</sup> ..... \$250

##### ☐ PLB4

Bodypoint Padded Lap Belt 4 Point<sup>(2)</sup> ..... \$315

##### ☐ TROBELT

Safety Belt for use only with TRRO & TRRO-E... N/C

1. Not available with TRRO and TRRO-E.

2. 67" length only option

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### MATRIX HEADREST PADS<sup>(1)</sup>

##### Standard Pads

|                                    |                                 |       |
|------------------------------------|---------------------------------|-------|
| <input type="checkbox"/> ESP6      | Matrx Elan Standard Pad 6"      | \$310 |
| <input type="checkbox"/> ESP10     | Matrx Elan Standard Pad 10"     | \$310 |
| <input type="checkbox"/> ESP14     | Matrx Elan Standard Pad 14"     | \$310 |
| <input type="checkbox"/> ESP6-HUG  | Matrx Elan Hug Standard Pad 6"  | \$310 |
| <input type="checkbox"/> ESP10-HUG | Matrx Elan Hug Standard Pad 10" | \$310 |
| <input type="checkbox"/> ESP14-HUG | Matrx Elan Hug Standard Pad 14" | \$310 |
| <input type="checkbox"/> SNUG14    | Matrx SNUG Standard Pad 14"     | \$625 |
| <input type="checkbox"/> SNUG14-A  | Matrx SNUG Asymmetrical Pad 14" | \$645 |
| <input type="checkbox"/> SNUG18    | Matrx SNUG Standard Pad 18"     | \$625 |
| <input type="checkbox"/> SNUG18-A  | Matrx SNUG Asymmetrical Pad 18" | \$645 |

##### Standard Pads with Infection Control Cover

|                                       |                                                   |       |
|---------------------------------------|---------------------------------------------------|-------|
| <input type="checkbox"/> ESP6-IC      | Matrx Elan Infection Control Pad 6"               | \$335 |
| <input type="checkbox"/> ESP10-IC     | Matrx Elan Infection Control Pad 10"              | \$335 |
| <input type="checkbox"/> ESP14-IC     | Matrx Elan Infection Control Pad 14"              | \$335 |
| <input type="checkbox"/> ESP6-IC-HUG  | Matrx Elan Hug Infection Control Pad 6"           | \$335 |
| <input type="checkbox"/> ESP10-IC-HUG | Matrx Elan Hug Infection Control Pad 10"          | \$335 |
| <input type="checkbox"/> ESP14-IC-HUG | Matrx Elan Hug Infection Control Pad 14"          | \$335 |
| <input type="checkbox"/> SNUG14-IC    | Matrx SNUG Standard Infection Control Pad 14"     | \$665 |
| <input type="checkbox"/> SNUG14-A-IC  | Matrx SNUG Asymmetrical Infection Control Pad 14" | \$685 |
| <input type="checkbox"/> SNUG18-IC    | Matrx SNUG Standard Infection Control Pad 18"     | \$665 |
| <input type="checkbox"/> SNUG18-A-IC  | Matrx SNUG Asymmetrical Infection Control Pad 18" | \$685 |

##### Occipital /4-Point Pads

|                                  |                                    |       |
|----------------------------------|------------------------------------|-------|
| <input type="checkbox"/> EOP9    | Matrx Elan Occipital Pad 9"        | \$341 |
| <input type="checkbox"/> EOP12   | Matrx Elan Occipital Pad 12"       | \$341 |
| <input type="checkbox"/> E4PTSM  | Matrx Elan 4-Point Pad 9"W x 7"H   | \$448 |
| <input type="checkbox"/> E4POINT | Matrx Elan 4-Point Pad 11"W x 10"H | \$448 |

##### Occipital / 4-Point Pads with Infection Control Cover

|                                     |                                                           |       |
|-------------------------------------|-----------------------------------------------------------|-------|
| <input type="checkbox"/> EOP9-IC    | Matrx Elan Occipital Pad 9" with Infection Control        | \$366 |
| <input type="checkbox"/> EOP12-IC   | Matrx Elan Occipital Pad 12" with Infection Control       | \$366 |
| <input type="checkbox"/> E4PTSM-IC  | Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control   | \$473 |
| <input type="checkbox"/> E4POINT-IC | Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control | \$473 |

1. If mounting on PB Elite must order removable hardware

#### MOTION CONCEPTS HEADREST PAD

|                               |                                           |       |
|-------------------------------|-------------------------------------------|-------|
| <input type="checkbox"/> MCSH | Motion Concepts Standard                  | N/C   |
| <input type="checkbox"/> MCOH | Motion Concepts Auto Style <sup>(1)</sup> | \$235 |

1. Comes Standard with Fixed Mounting hardware; for use on Matrx PB or Elite Backs you must select the removable hardware upgrade below.

#### HEADREST MOUNTING HARDWARE

|                                     |                                   |       |
|-------------------------------------|-----------------------------------|-------|
| <input type="checkbox"/> FM         | Fixed Mounting                    | N/C   |
| <input type="checkbox"/> LXSHW      | Loxx Standard Hardware 15"        | \$275 |
| <input type="checkbox"/> LXMMH      | Loxx Mini Hardware 11"            | \$275 |
| <input type="checkbox"/> MEHWRE     | Elan Hardware 15"                 | \$285 |
| <input type="checkbox"/> MEMH       | Elan Mini Hardware 11"            | \$285 |
| <input type="checkbox"/> MEHW-EXT   | Elan Extended Hardware            | \$320 |
| <input type="checkbox"/> FPLXSHW    | FLIPP Loxx Standard Hardware 15"  | \$705 |
| <input type="checkbox"/> FPLXSHW-EX | FLIPP Loxx Extended Hardware      | \$705 |
| <input type="checkbox"/> FPLXMMH    | FLIPP Loxx Mini Hardware 11"      | \$705 |
| <input type="checkbox"/> FPLXMMH-EX | FLIPP Loxx Mini Extended Hardware | \$705 |
| <input type="checkbox"/> FPMEHW     | FLIPP Elan Hardware 15"           | \$715 |
| <input type="checkbox"/> FPMEH      | FLIPP Elan Mini Hardware 11"      | \$715 |
| <input type="checkbox"/> FPMEHW-EXT | FLIPP Elan Extended Hardware      | \$780 |

|                              |                |     |
|------------------------------|----------------|-----|
| <input type="checkbox"/> OHR | OMIT: Headrest | N/C |
|------------------------------|----------------|-----|

#### STANDARD BACK CANES - FOR TILT SYSTEMS

|                                |                         |     |
|--------------------------------|-------------------------|-----|
| <input type="checkbox"/> SBC18 | 18" Straight Back Canes | STD |
| <input type="checkbox"/> SBC20 | 20" Straight Back Canes | STD |
| <input type="checkbox"/> SBC22 | 22" Straight Back Canes | STD |
| <input type="checkbox"/> SBC24 | 24" Straight Back Canes | STD |

##### Select Cane Angle

|                                     |                                     |                                     |                                     |                                   |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-----------------------------------|
| <input type="checkbox"/> BA75 75°   | <input type="checkbox"/> BA79 79°   | <input type="checkbox"/> BA85 85°   | <input type="checkbox"/> BA90 90°   | <input type="checkbox"/> BA95 95° |
| <input type="checkbox"/> BA101 101° | <input type="checkbox"/> BA106 106° | <input type="checkbox"/> BA112 112° | <input type="checkbox"/> BA116 116° |                                   |

|                                |                            |     |
|--------------------------------|----------------------------|-----|
| <input type="checkbox"/> ABC20 | 20" Angled (5°) Back Canes | N/C |
| <input type="checkbox"/> ABC22 | 22" Angled (5°) Back Canes | N/C |
| <input type="checkbox"/> ABC24 | 24" Angled (5°) Back Canes | N/C |

##### Select Cane Angle

|                                     |                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> BA81 81°   | <input type="checkbox"/> BA84 84°   | <input type="checkbox"/> BA91 91°   | <input type="checkbox"/> BA96 96°   | <input type="checkbox"/> BA101 101° |
| <input type="checkbox"/> BA106 106° | <input type="checkbox"/> BA111 111° | <input type="checkbox"/> BA117 117° | <input type="checkbox"/> BA121 121° |                                     |

|                                |                                         |       |
|--------------------------------|-----------------------------------------|-------|
| <input type="checkbox"/> EBC22 | 22" Endomorph Back Canes <sup>(1)</sup> | \$250 |
| <input type="checkbox"/> EBC24 | 24" Endomorph Back Canes <sup>(1)</sup> | \$250 |

##### Please Select

|                                |             |     |
|--------------------------------|-------------|-----|
| <input type="checkbox"/> EBCN1 | 1" Narrower | N/C |
| <input type="checkbox"/> EBCN2 | 2" Narrower | N/C |
| <input type="checkbox"/> EBCN4 | 4" Narrower | N/C |

|                                |                                                  |       |
|--------------------------------|--------------------------------------------------|-------|
| <input type="checkbox"/> TBC18 | 18" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup> | \$250 |
| <input type="checkbox"/> TBC20 | 20" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup> | \$250 |
| <input type="checkbox"/> TBC22 | 22" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup> | \$250 |
| <input type="checkbox"/> TBC24 | 24" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup> | \$250 |

|                               |                                                 |       |
|-------------------------------|-------------------------------------------------|-------|
| <input type="checkbox"/> FDBC | Fold Down Back Canes w/ Pull Pin <sup>(3)</sup> | \$250 |
|-------------------------------|-------------------------------------------------|-------|

1. Comes standard with Cantilever Armrest Cane Mounted.  
2. Quick Height Adjustment on Backrest and Maxx Spreader Bars not available.  
3. 250 lbs weight capacity. Must select canes above.

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### CANE OPTIONS - FOR RECLINE SYSTEMS

|                               |                                                                 |       |
|-------------------------------|-----------------------------------------------------------------|-------|
| <input type="checkbox"/> BFS2 | Recline Body Fit 2" Narrower <sup>(1)</sup> .....               | \$995 |
| <input type="checkbox"/> BFS4 | Recline Body Fit 4" Narrower <sup>(1)</sup> .....               | \$995 |
| <input type="checkbox"/> AHH  | Adjustable Hooking Handles <sup>(2)</sup> (pair) .....          | \$450 |
| <input type="checkbox"/> RPHS | Push Handles for Recline - Straight <sup>(3)</sup> (pair) ..... | \$350 |
| <input type="checkbox"/> RPHA | Push Handles for Recline - Angled <sup>(3)</sup> (pair) .....   | \$350 |

1. Minimum finished back widths of 16".

2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.

3. 4" of total vertical mounting adjustment.

#### FINISHED BACK HEIGHT<sup>(1)(2)</sup> - FOR TILT ONLY/ELEVATE ONLY

|                                    |                                    |                                    |                                    |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| <input type="checkbox"/> FBT18 18" | <input type="checkbox"/> FBT19 19" | <input type="checkbox"/> FBT20 20" | <input type="checkbox"/> FBT21 21" |
| <input type="checkbox"/> FBT22 22" | <input type="checkbox"/> FBT23 23" | <input type="checkbox"/> FBT24 24" | <input type="checkbox"/> FBT25 25" |

1. 4" Gap from the seat pan to the bottom of back pan is recommended.

2. Finished back height is measured from seat pan to top of backrest.

#### FINISHED BACK HEIGHT<sup>(1)(2)</sup> - FOR RECLINE OR TILT/RECLINE

|                                    |                                    |                                    |                                    |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| <input type="checkbox"/> FBR20 20" | <input type="checkbox"/> FBR21 21" | <input type="checkbox"/> FBR22 22" | <input type="checkbox"/> FBR23 23" |
| <input type="checkbox"/> FBR24 24" | <input type="checkbox"/> FBR25 25" | <input type="checkbox"/> FBR26 26" | <input type="checkbox"/> FBR27 27" |

1. 6" Gap from the seat pan to the bottom of back pan is recommended.

2. Finished back height is measured from seat pan to top of backrest.

#### STANDARD BACK OPTIONS

|                               |                                                                      |       |
|-------------------------------|----------------------------------------------------------------------|-------|
| <input type="checkbox"/> SRB  | Standard Rehab Back & Cushion .....                                  | N/C   |
| <input type="checkbox"/> SS   | Additional Super Soft/HR Foam <sup>(1)</sup> .....                   | \$450 |
| <input type="checkbox"/> MFU  | Matrix Fabric Upcharge .....                                         | \$450 |
| <input type="checkbox"/> OBC  | OMIT: Back Cushion .....                                             | N/C   |
| <input type="checkbox"/> PIPF | Flat Planar Interface/I-Back Plate <sup>(2)(3)(5)</sup> .....        | \$125 |
| <input type="checkbox"/> PIPR | Recessed Planar/I-Back Interface Plate <sup>(2)(3)(5)(6)</sup> ..... | \$125 |
| <input type="checkbox"/> ORBC | OMIT: Rigid Back & Cushion <sup>(4)</sup> .....                      | N/C   |

1. Additional foam added may impact seat depth

2. Not available on tilt only or tilt/elevate systems, for use with Recline Module.

3. Required with all Recline Systems when no back is selected. May also be used as an interface for after market backs.

4. Not available with Recline/ESR or Recline/PSB.

5. Not available with TRRO and TRRO-E.

6. Fits Comfort Company, Freedom, Pindot and Moulded Backs.

#### MOTION CONCEPTS SINGLE POST BACKREST<sup>(1)(2)(3)</sup>

|                                |                                                                                           |                                 |
|--------------------------------|-------------------------------------------------------------------------------------------|---------------------------------|
| <input type="checkbox"/> SPBMI | Single Post Angle Adjustable Fold Down<br>Back interface-Matrix Back.....                 | \$750<br>(Folds to 25° forward) |
| <input type="checkbox"/> SPBCI | Single Post Angle Adjustable Fold Down<br>Back for Contour/High Back <sup>(4)</sup> ..... | \$750<br>(Folds to 25° forward) |

1. Single Post Back is an angle adjustable back between 73° and 109° in 6° increments from 17° forward to 19° back. The back folds down and is only available for 16"-20" seat widths &amp; 21"-24" back height. A 4" gap from bottom of back to seat pan is required.

2. This option is not available with Recline and/or a vent tray and has a 250lb weight capacity.

3. Cane mounted armrest are not available with this option.

4. Not available with TRRO/TRRO-E

#### UPGRADEABLE BACK OPTIONS

Contour Back 15" wide - 20" wide ..... \$495

|                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> CB1512 | <input type="checkbox"/> CB1515 | <input type="checkbox"/> CB1517 | <input type="checkbox"/> CB1612 |
| <input type="checkbox"/> CB1615 | <input type="checkbox"/> CB1617 | <input type="checkbox"/> CB1621 | <input type="checkbox"/> CB1712 |
| <input type="checkbox"/> CB1715 | <input type="checkbox"/> CB1717 | <input type="checkbox"/> CB1721 | <input type="checkbox"/> CB1812 |
| <input type="checkbox"/> CB1815 | <input type="checkbox"/> CB1817 | <input type="checkbox"/> CB1821 | <input type="checkbox"/> CB1912 |
| <input type="checkbox"/> CB1915 | <input type="checkbox"/> CB1917 | <input type="checkbox"/> CB1921 | <input type="checkbox"/> CB2012 |
| <input type="checkbox"/> CB2015 | <input type="checkbox"/> CB2017 | <input type="checkbox"/> CB2021 |                                 |

Contour Back 21"-22" wide ..... \$595

|                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> CB2115 | <input type="checkbox"/> CB2117 | <input type="checkbox"/> CB2121 | <input type="checkbox"/> CB2215 |
| <input type="checkbox"/> CB2217 | <input type="checkbox"/> CB2221 |                                 |                                 |

High Back 16" - 22" W ..... \$595

|                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> HB1616 | <input type="checkbox"/> HB1618 | <input type="checkbox"/> HB1620 | <input type="checkbox"/> HB1716 |
| <input type="checkbox"/> HB1718 | <input type="checkbox"/> HB1720 | <input type="checkbox"/> HB1816 | <input type="checkbox"/> HB1818 |
| <input type="checkbox"/> HB1820 | <input type="checkbox"/> HB1916 | <input type="checkbox"/> HB1918 | <input type="checkbox"/> HB1920 |
| <input type="checkbox"/> HB2016 | <input type="checkbox"/> HB2018 | <input type="checkbox"/> HB2020 | <input type="checkbox"/> HB2116 |
| <input type="checkbox"/> HB2118 | <input type="checkbox"/> HB2120 | <input type="checkbox"/> HB2216 | <input type="checkbox"/> HB2218 |
| <input type="checkbox"/> HB2220 |                                 |                                 |                                 |

Matrix Elite Back 14" - 20" W (3" Contour)..... \$766

|                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PBE1414 | <input type="checkbox"/> PBE1416 | <input type="checkbox"/> PBE1418 | <input type="checkbox"/> PBE1420 |
| <input type="checkbox"/> PBE1514 | <input type="checkbox"/> PBE1516 | <input type="checkbox"/> PBE1518 | <input type="checkbox"/> PBE1520 |
| <input type="checkbox"/> PBE1614 | <input type="checkbox"/> PBE1616 | <input type="checkbox"/> PBE1618 | <input type="checkbox"/> PBE1620 |
| <input type="checkbox"/> PBE1714 | <input type="checkbox"/> PBE1716 | <input type="checkbox"/> PBE1718 | <input type="checkbox"/> PBE1720 |
| <input type="checkbox"/> PBE1814 | <input type="checkbox"/> PBE1816 | <input type="checkbox"/> PBE1818 | <input type="checkbox"/> PBE1820 |
| <input type="checkbox"/> PBE1916 | <input type="checkbox"/> PBE1918 | <input type="checkbox"/> PBE1920 | <input type="checkbox"/> PBE2014 |
| <input type="checkbox"/> PBE2016 | <input type="checkbox"/> PBE2018 | <input type="checkbox"/> PBE2020 |                                  |

Matrix Elite HD Back 21" - 22" W (3" Contour)

..... \$1,095

|                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> PBE2116-HD | <input type="checkbox"/> PBE2120-HD | <input type="checkbox"/> PBE2216-HD | <input type="checkbox"/> PBE2220-HD |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|

Matrix Elite TR Back 15" - 20" W (3" Contour)

..... \$866

|                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> PBE1516-TR | <input type="checkbox"/> PBE1518-TR | <input type="checkbox"/> PBE1520-TR | <input type="checkbox"/> PBE1616-TR |
| <input type="checkbox"/> PBE1618-TR | <input type="checkbox"/> PBE1620-TR | <input type="checkbox"/> PBE1716-TR | <input type="checkbox"/> PBE1718-TR |
| <input type="checkbox"/> PBE1720-TR | <input type="checkbox"/> PBE1816-TR | <input type="checkbox"/> PBE1818-TR | <input type="checkbox"/> PBE1820-TR |
| <input type="checkbox"/> PBE1916-TR | <input type="checkbox"/> PBE1918-TR | <input type="checkbox"/> PBE1920-TR | <input type="checkbox"/> PBE2016-TR |
| <input type="checkbox"/> PBE2018-TR | <input type="checkbox"/> PBE2020-TR |                                     |                                     |

Matrix Elite TR HD Back 21" - 22" W (3" Contour)

..... \$1,151

|                                       |                                       |
|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> PBE2120-TRHD | <input type="checkbox"/> PBE2220-TRHD |
|---------------------------------------|---------------------------------------|

Matrix Elite Deep Back 14" - 20" W (6" Contour)

..... \$918

|                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> EDB1416 | <input type="checkbox"/> EDB1418 | <input type="checkbox"/> EDB1420 | <input type="checkbox"/> EDB1514 |
| <input type="checkbox"/> EDB1516 | <input type="checkbox"/> EDB1518 | <input type="checkbox"/> EDB1520 | <input type="checkbox"/> EDB1614 |
| <input type="checkbox"/> EDB1616 | <input type="checkbox"/> EDB1618 | <input type="checkbox"/> EDB1620 | <input type="checkbox"/> EDB1714 |
| <input type="checkbox"/> EDB1716 | <input type="checkbox"/> EDB1718 | <input type="checkbox"/> EDB1720 | <input type="checkbox"/> EDB1814 |
| <input type="checkbox"/> EDB1816 | <input type="checkbox"/> EDB1818 | <input type="checkbox"/> EDB1820 | <input type="checkbox"/> EDB1914 |
| <input type="checkbox"/> EDB1916 | <input type="checkbox"/> EDB1918 | <input type="checkbox"/> EDB1920 | <input type="checkbox"/> EDB2014 |
| <input type="checkbox"/> EDB2016 | <input type="checkbox"/> EDB2018 | <input type="checkbox"/> EDB2020 |                                  |

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### UPGRADEABLE BACK OPTIONS

Matrx Elite Deep HD Back 21" - 22" W (6" Contour) ..... \$1,095

☐ EDB2116-HD    ☐ EDB2118-HD    ☐ EDB2120-HD  
☐ EDB2216-HD    ☐ EDB2218-HD    ☐ EDB2220-HD

Matrx Elite E2S Back 14" - 20" W x 14"/16"H (3" Contour) ..... \$1,125

☐ E2S1414    ☐ E2S1416    ☐ E2S1514    ☐ E2S1516  
☐ E2S1614    ☐ E2S1616    ☐ E2S1714    ☐ E2S1716  
☐ E2S1814    ☐ E2S1816    ☐ E2S1914    ☐ E2S1916  
☐ E2S2014    ☐ E2S2016

Matrx Elite E2S Back 14" - 20" W x 18"/20"H (3" Contour) ..... \$995

☐ E2S1418    ☐ E2S1420    ☐ E2S1518    ☐ E2S1520  
☐ E2S1618    ☐ E2S1620    ☐ E2S1718    ☐ E2S1720  
☐ E2S1818    ☐ E2S1820    ☐ E2S1918    ☐ E2S1920  
☐ E2S2018    ☐ E2S2020

Matrx Elite E2 Deep Back  
14" - 20" W x 14"/16"H (6" Contour) ..... \$1,125

☐ E2D1414    ☐ E2D1416    ☐ E2D1514    ☐ E2D1516  
☐ E2D1614    ☐ E2D1616    ☐ E2D1714    ☐ E2D1716  
☐ E2D1814    ☐ E2D1816    ☐ E2D1914    ☐ E2D1916  
☐ E2D2014    ☐ E2D2016

Matrx Elite E2 Deep Back  
14" - 20" W x 18"/20"H (6" Contour) ..... \$1,155

☐ E2D1418    ☐ E2D1420    ☐ E2D1518    ☐ E2D1520  
☐ E2D1618    ☐ E2D1620    ☐ E2D1718    ☐ E2D1720  
☐ E2D1818    ☐ E2D1820    ☐ E2D1918    ☐ E2D1920  
☐ E2D2018    ☐ E2D2020

Matrx PB Posture Back  
14" - 20" W (5" Contour) ..... \$734

☐ MPB1412    ☐ MPB1612    ☐ MPB1616    ☐ MPB1620  
☐ MPB1712    ☐ MPB1716    ☐ MPB1720    ☐ MPB1812  
☐ MPB1816    ☐ MPB1820    ☐ MPB2012    ☐ MPB2016  
☐ MPB2020

Matrx PB Heavy Duty Back<sup>(1)</sup>  
20" - 22" W (5" Contour) ..... \$1,051

☐ HDB2016    ☐ HDB2020    ☐ HDB2216    ☐ HDB2220

Matrx PB Posture Back Deep<sup>(1)</sup>  
14" - 20" W (7" Contour) ..... \$892

☐ PBD1412    ☐ PBD1612    ☐ PBD1616    ☐ PBD1620  
☐ PBD1712    ☐ PBD1716    ☐ PBD1720    ☐ PBD1812  
☐ PBD1816    ☐ PBD1820    ☐ PBD2012    ☐ PBD2016  
☐ PBD2020

Matrx PB Heavy Duty Deep Back<sup>(1)</sup>  
20"-22" W (5" Contour) ..... \$1,209

☐ DHDB2016    ☐ DHDB2020    ☐ DHDB2216    ☐ DHDB2220

Matrx MAC Back<sup>(2)</sup> ..... \$1,495

☐ MAC1616    ☐ MAC1618    ☐ MAC1620    ☐ MAC1622  
☐ MAC1716    ☐ MAC1718    ☐ MAC1720    ☐ MAC1722  
☐ MAC1816    ☐ MAC1818    ☐ MAC1820    ☐ MAC1822  
☐ MAC1916    ☐ MAC1918    ☐ MAC1920    ☐ MAC1922  
☐ MAC2016    ☐ MAC2018    ☐ MAC2020    ☐ MAC2022

Matrx MAC Back w/ HUG<sup>(2)</sup> ..... \$1,645

☐ MAC1616-HUG    ☐ MAC1618-HUG    ☐ MAC1620-HUG  
☐ MAC1622-HUG    ☐ MAC1716-HUG    ☐ MAC1718-HUG  
☐ MAC1720-HUG    ☐ MAC1722-HUG    ☐ MAC1816-HUG  
☐ MAC1818-HUG    ☐ MAC1820-HUG    ☐ MAC1822-HUG  
☐ MAC1916-HUG    ☐ MAC1918-HUG    ☐ MAC1920-HUG  
☐ MAC1922-HUG    ☐ MAC2016-HUG    ☐ MAC2018-HUG  
☐ MAC2020-HUG    ☐ MAC2022-HUG

Matrx MAC Lateral Contour

☐ MLCSL    Left 3.5"    ☐ MLC SR    Right 3.5"  
☐ MLC DL    Left 6"    ☐ MLC DR    Right 6"  
☐ MLC XDL    Left 9"    ☐ MLC XDR    Right 9"

1. The Heavy Duty PB Back and PB Deep Back are not available on Recline systems with seat widths of 21" and 22"

2. Must select Pad Contour in Matrx MAC Lateral Contour section.

#### OTHER UPGRADABLE BACK OPTIONS

| Cushion                         | Meshtex | 4-Way Stretch Polyester | Infection Control | Onyx  | Fusion |
|---------------------------------|---------|-------------------------|-------------------|-------|--------|
| Contour                         | N/A     | STD                     | N/A               | \$100 | N/A    |
| High                            | STD     | N/A                     | N/A               | \$100 | \$100  |
| E2, Elite, Elite HD, Elite Deep | STD     | N/A                     | \$100             | N/A   | N/A    |
| E2, Elite TR, Elite TR HD       | STD     | N/A                     | \$100             | N/A   | N/A    |

#### OTHER UPGRADABLE BACK OPTIONS

☐ EBSS    Super Soft Foam for Elite Back Upgrade. .... \$150  
☐ FFUHB    Fusion Fabric Upgrade for High Backs. .... \$150  
☐ OCUCB    Onyx Fabric Upgrade for Contour Backs. .... \$150  
☐ OCUHB    Onyx Fabric Upgrade for High Backs. .... \$150  
☐ CHIEB    Chest Harness Interface Kit for Elite Back<sup>(1)(2)</sup> .... \$103  
☐ CHIPB    Chest Harness Interface Kit for PB Back<sup>(1)(2)</sup> .... \$52  
☐ E2CHI-S    Chest Harness Interface Kit for E2 Back for 14"-15" Backs<sup>(1)(2)</sup> ..... \$85  
☐ E2CHI-M    Chest Harness Interface Kit for E2 Back for 16"-20" Backs<sup>(1)(2)</sup> ..... \$85  
☐ E2BKB    Contoured anatomical PSIS/Sacrolumbar Support ..... \$75

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harness

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### CHEST STRAPS - ABDOMINAL SUPPORT

|                                  |                                                                                 |       |
|----------------------------------|---------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> MFUS    | BP Monoflex, Underarm RLS S(3" x 15.5")                                         | \$195 |
| <input type="checkbox"/> MFUM    | BP Monoflex, Underarm RLS M(4" x 18.25")                                        | \$195 |
| <input type="checkbox"/> MFUL    | BP Monoflex, Underarm RLS L(4.5" x 21.5")                                       | \$195 |
| <input type="checkbox"/> MFCS    | BP Monoflex, Center RLS S(2.5" x 15.5")                                         | \$195 |
| <input type="checkbox"/> MFCM    | BP Monoflex, Center RLS M(3.25" x 18.25")                                       | \$195 |
| <input type="checkbox"/> MFCL    | BP Monoflex, Center RLS L(3.5" x 21.5")                                         | \$195 |
| <input type="checkbox"/> PFLEX-M | PostureFlex Abdominal Support<br>Height 6.5", length 16" to 23" <sup>(1)</sup>  | \$245 |
| <input type="checkbox"/> PFLEX-L | PostureFlex Abdominal Support<br>Height 7", length 20" to 32" <sup>(1)</sup>    | \$245 |
| <input type="checkbox"/> QRTAB   | Quad release tab for easy independent<br>application and removal of PostureFlex | \$35  |

1. Only available on the E2 and E2 Deep Backs.

#### TILT ARMRESTS

##### Adult Dual Post

|                              |                          |       |
|------------------------------|--------------------------|-------|
| <input type="checkbox"/> DP1 | 9.75" - 12.75"           | STD   |
| <input type="checkbox"/> DP2 | 12.75" - 15.75"          | STD   |
| <input type="checkbox"/> DP3 | 15.75" - 18.75"          | STD   |
| <input type="checkbox"/> FP  | Flat Pouch               | STD   |
| <input type="checkbox"/> GB  | Glove Boxes              | N/C   |
| <input type="checkbox"/> OAL | Outback Arm Option Left  | \$195 |
| <input type="checkbox"/> OAR | Outback Arm Option Right | \$195 |

##### Ultra Rail Seat Mount Flip Back Cantilever MAXX Tilt Arm<sup>(1)</sup>

|                                 |                                  |       |
|---------------------------------|----------------------------------|-------|
| <input type="checkbox"/> RMCAN  | 9" - 12.5"                       | \$488 |
| <input type="checkbox"/> RMCAN2 | 11" - 14"                        | \$488 |
| <input type="checkbox"/> DPM    | Dual Post Module to RMCAN (pair) | \$250 |

##### Outback Cantilever Tilt Armrest

|                                |         |       |
|--------------------------------|---------|-------|
| <input type="checkbox"/> OCTA1 | 9"-12"  | \$488 |
| <input type="checkbox"/> OCTA2 | 12"-15" | \$488 |
| <input type="checkbox"/> OCTA2 | 6"-12"  | \$488 |

##### Cantilever Arms for Endomorph Back Canes<sup>(1)</sup>

|                                |            |     |
|--------------------------------|------------|-----|
| <input type="checkbox"/> CAFEB | See page 4 | STD |
|--------------------------------|------------|-----|

1. Not available with Small Adult Sizes.

#### RECLINE ARMRESTS

##### Reclining Adjustable Height

|                              |                          |       |
|------------------------------|--------------------------|-------|
| <input type="checkbox"/> RC1 | 9.5"-13"                 | STD   |
| <input type="checkbox"/> RC2 | 13"-16"                  | STD   |
| <input type="checkbox"/> OAL | Outback Arm Option Left  | \$195 |
| <input type="checkbox"/> OAR | Outback Arm Option Right | \$195 |

##### Maxx Style Cane Mounted Height, Width, Angle Adj.

##### Cantilever Flip Back Armrests<sup>(1)</sup>

|                                |           |       |
|--------------------------------|-----------|-------|
| <input type="checkbox"/> MCAN1 | 9.5"-13"  | \$250 |
| <input type="checkbox"/> MCAN2 | 12.5"-16" | \$250 |

##### Outback Cantilever Recline Armrest

|                                |         |       |
|--------------------------------|---------|-------|
| <input type="checkbox"/> OCRA1 | 9"-12"  | \$350 |
| <input type="checkbox"/> OCRA2 | 12"-15" | \$350 |
| <input type="checkbox"/> OCRA3 | 7" - 9" | \$350 |

##### Cane Mounted Height and Angle Adj.

##### Cantilever Flip Back Armrest<sup>(1)</sup>

|                                |          |       |
|--------------------------------|----------|-------|
| <input type="checkbox"/> CANQ1 | 9.5"-12" | \$350 |
| <input type="checkbox"/> CANQ2 | 10"-13"  | \$350 |
| <input type="checkbox"/> CANQ3 | 13"-16"  | \$350 |

##### Adult Dual Post On Recline

|                                |                 |       |
|--------------------------------|-----------------|-------|
| <input type="checkbox"/> DPOR1 | 9.75" - 12.75"  | \$250 |
| <input type="checkbox"/> DPOR2 | 12.75" - 15.75" | \$250 |
| <input type="checkbox"/> DPOR3 | 15.75" - 18.75" | \$250 |
| <input type="checkbox"/> FP    | Flat Pouch      | STD   |
| <input type="checkbox"/> GB    | Glove Boxes     | N/C   |

1. Not available with Small Adult Sizes.

#### ARMREST TUBE LENGTH - PLEASE SPECIFY

|                               |            |     |
|-------------------------------|------------|-----|
| <input type="checkbox"/> ATFL | Left Full  | STD |
| <input type="checkbox"/> ATDL | Left Desk  | STD |
| <input type="checkbox"/> ATFR | Right Full | STD |
| <input type="checkbox"/> ATDR | Right Desk | STD |

#### ARMPADS<sup>(5)</sup>

##### Standard

|                                |                          |     |
|--------------------------------|--------------------------|-----|
| <input type="checkbox"/> SAPFL | Left Full (2.25" x 14")  | STD |
| <input type="checkbox"/> SAPDL | Left Desk (2.25" x 10")  | STD |
| <input type="checkbox"/> SAPML | Left Mini (2.25" x 7")   | STD |
| <input type="checkbox"/> SAPFR | Right Full (2.25" x 14") | STD |
| <input type="checkbox"/> SAPDR | Right Desk (2.25" x 10") | STD |
| <input type="checkbox"/> SAPMR | Right Mini (2.25" x 7")  | STD |

##### Waterfall Arm Pads<sup>(1)</sup>

|                                |                  |     |
|--------------------------------|------------------|-----|
| <input type="checkbox"/> WAPFL | Left Full (13")  | N/C |
| <input type="checkbox"/> WAPDL | Left Desk (9")   | N/C |
| <input type="checkbox"/> WAPFR | Right Full (13") | N/C |
| <input type="checkbox"/> WAPDR | Right Desk (9")  | N/C |

##### Modular Arm Pad - Left<sup>(2)</sup>

|                                 |                      |      |
|---------------------------------|----------------------|------|
| <input type="checkbox"/> MAPTFL | Left Tray Full (14") | \$95 |
| <input type="checkbox"/> MAPTDL | Left Tray Desk (10") | \$95 |

Select Pad Insert below:

|                                |                                   |
|--------------------------------|-----------------------------------|
| <input type="checkbox"/> ISSFL | Iskin - Integrated soft skin foam |
| <input type="checkbox"/> SCVFL | Startex Covered Visco foam        |

##### Modular Arm Pad - Right<sup>(2)</sup>

|                                 |                       |      |
|---------------------------------|-----------------------|------|
| <input type="checkbox"/> MAPTFR | Right Tray Full (14") | \$95 |
| <input type="checkbox"/> MAPTDR | Right Tray Desk (10") | \$95 |

Select Pad Insert below:

|                                |                                   |
|--------------------------------|-----------------------------------|
| <input type="checkbox"/> ISSFR | Iskin - Integrated soft skin foam |
| <input type="checkbox"/> SCVFR | Startex Covered Visco foam        |

##### Flat Multi Position Arm Pads

|                                |                  |       |
|--------------------------------|------------------|-------|
| <input type="checkbox"/> MPAFL | Left Full (14")  | \$150 |
| <input type="checkbox"/> MPADL | Left Desk (10")  | \$150 |
| <input type="checkbox"/> MPAML | Left Mini (7")   | \$150 |
| <input type="checkbox"/> MPAFR | Right Full (14") | \$150 |
| <input type="checkbox"/> MPADR | Right Desk (10") | \$150 |
| <input type="checkbox"/> MPAMR | Right Mini (7")  | \$150 |

##### Gel Pads

|                                |                              |       |
|--------------------------------|------------------------------|-------|
| <input type="checkbox"/> GAPNL | Narrow Left Desk (2" x 12")  | \$88  |
| <input type="checkbox"/> GAPNR | Narrow Right Desk (2" x 12") | \$88  |
| <input type="checkbox"/> GAPWL | Wide Left Desk (3.5" x 12")  | \$96  |
| <input type="checkbox"/> GAPWR | Wide Right Desk (3.5" x 12") | \$96  |
| <input type="checkbox"/> GAPLL | Long Left Full (3.5" x 14")  | \$103 |
| <input type="checkbox"/> GAPLR | Long Right Full (3.5" x 14") | \$103 |

##### Arm Troughs

|                               |                                      |       |
|-------------------------------|--------------------------------------|-------|
| <input type="checkbox"/> EATL | Ergonomic Arm Trough - Left          | \$121 |
| <input type="checkbox"/> EATR | Ergonomic Arm Trough - Right         | \$121 |
| <input type="checkbox"/> FHPL | Flat Hand Pad - Left <sup>(3)</sup>  | \$91  |
| <input type="checkbox"/> FHPR | Flat Hand Pad - Right <sup>(3)</sup> | \$91  |

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### ARMPADS<sup>(5)</sup> - CONTINUED

##### Arm Pad Storage/Cup Holder<sup>(4)</sup>

- ☐ APSL Left .....\$295  
☐ APSR Right .....\$295  
☐ OAP Omit Armpads(7/8" diameter tube) .....N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side, or with Outback Cantilever arms.
5. Full length arm pads used on the same side as the Joystick result in a "Back to Joystick" measurement of 20" minimum (24" including switch mounting)

#### ARM ADDUCTOR PAD<sup>(1)(2)</sup>

- ☐ AADSL Left Small (6" x 2.25") .....\$140  
☐ AADSR Right Small (6" x 2.25") .....\$140  
☐ AADLL Left Large (8.5" x 2.25") .....\$140  
☐ AADLR Right Large (8.5" x 2.25") .....\$140

1. Only available with Modular or Flat Multi Position Arm Pads.
2. 1" thick pad, unlimited fore/aft adjustment.

#### ELBOW BLOCKS

##### Left Elbow Block<sup>(1)</sup>

Select Pad Size and Both Extensions:

- ☐ EBXL Size XS (3.25" x 4.25") .....\$195  
☐ EBSL Size S (3.5" x 5.25") .....\$195  
☐ EBML Size M (3.75" x 5.75") .....\$195  
☐ EBL L Size L (4.25" x 6.5") .....\$195

- ☐ LV3 Left Vertical Extension 3"  
☐ LV5 Left Vertical Extension 5"  
☐ LH5 Left Horizontal Extension 5"  
☐ LH7 Left Horizontal Extension 7"

##### Right Elbow Block<sup>(1)</sup>

Select Pad Size and Both Extensions:

- ☐ EBXR Size XS (3.25" x 4.25") .....\$195  
☐ EBSR Size S (3.5" x 5.25") .....\$195  
☐ EBMR Size M (3.75" x 5.75") .....\$195  
☐ EBLR Size L (4.25" x 6.5") .....\$195

- ☐ RV3 Right Vertical Extension 3"  
☐ RV5 Right Vertical Extension 5"  
☐ RH5 Right Horizontal Extension 5"  
☐ RH7 Right Horizontal Extension 7"

1. Prevents arm from slipping.

#### ADDITIONAL ARM ACCESSORIES

- ☐ MACESL Multi Axis Upper Extremity Support<sup>(1)</sup>(each)  
Left Side .....\$225  
☐ MACESR Multi Axis Upper Extremity Support<sup>(1)</sup> (each)  
Right Side .....\$225  
☐ OAS2L 1" Offset Armrest Spacer Left(each) .....\$125  
☐ OAS2R 1" Offset Armrest Spacer Right(each) .....\$125  
☐ OAS1L 1/2" Offset Armrest Spacer Left(each) .....\$125  
☐ OAS1R 1/2" Offset Armrest Spacer Right(each) .....\$125  
☐ MPAPA Adapter Block<sup>(2)</sup> (each) Quantity ( ) .....\$35

1. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
2. Not available with Standard Arm Pads or Waterfall Arm Pads

#### MAXX LATERAL TRUNK SUPPORTS<sup>(1)</sup>

##### Lateral Trunk Support Left

- ☐ LATXL Size XS (3.25" x 4.25") .....\$100  
☐ LATSL Size S (3.5" x 5.25") .....\$100  
☐ LATML Size M (3.75" x 5.75") .....\$100  
☐ LATLL Size L (4.25" x 6.5") .....\$100

Add:

- ☐ SALATL MAXX Style Swing Away Hardware - Left<sup>(2)</sup> .....\$195  
☐ ATLL Additional Telescoping Link - Left<sup>(3)</sup> .....\$75

##### Lateral Trunk Support Right

- ☐ LATXR Size XS (3.25" x 4.25") .....\$100  
☐ LATSR Size S (3.5" x 5.25") .....\$100  
☐ LATMR Size M (3.75" x 5.75") .....\$100  
☐ LATLR Size L (4.25" x 6.5") .....\$100

Add:

- ☐ SALATR MAXX Style Swing Away Hardware - Right<sup>(2)</sup> .....\$195  
☐ ATLR Additional Telescoping Link - Right<sup>(3)</sup> .....\$75

1. Comes standard with Fixed Mount Hardware. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Comes standard with 1" telescoping link which gives 0-2.75" medial offset.
3. The amount of offset available may be impacted by back thickness.

#### MATRIX FIXED LATERAL TRUNK SUPPORTS<sup>(1)</sup>

- ☐ E2FLL Matrix E2 Fixed Laterals-Left(6.5x5")<sup>(2)</sup> .....\$195  
☐ E2FLR Matrix E2 Fixed Laterals-Right(6.5x5")<sup>(2)</sup> .....\$195  
☐ PBSSL Matrix Small Lateral Fixed-Left (3x5.5") .....\$113  
☐ PBSSR Matrix Small Lateral Fixed-Right (3x5.5") .....\$113  
☐ PBLSL Matrix Standard Lateral Fixed-Left (7.5x5") .....\$116  
☐ PBLSR Matrix Standard Lateral Fixed-Right (7.5x5") .....\$116  
☐ PBOLL Matrix Offset Lateral Fixed Left (7.5x5") .....\$125  
☐ PBOLR Matrix Offset Lateral Fixed Right (7.5x5") .....\$125

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Only available with E2 backs.

#### MATRIX SWING AWAY LATERAL TRUNK SUPPORTS<sup>(1)</sup>

##### Matrix Lateral Trunk Support Left

- ☐ EPADXL Size XS (3.25" x 4.25") .....\$115  
☐ EPADSL Size S (3.5" x 5.25") .....\$115  
☐ EPADML Size M (3.75" x 5.75") .....\$115  
☐ EPADLL Size L (4.25" x 6.5") .....\$115

Must select one option:

- ☐ SAHEBL Matrix Elite Swing Away Hardware Left .....\$250  
☐ OSAHL Matrix Offset Elite Swing Away Hardware Left .....\$250

##### Matrix Lateral Trunk Support Right

- ☐ EPADXR Size XS (3.25" x 4.25") .....\$115  
☐ EPADSR Size S (3.5" x 5.25") .....\$115  
☐ EPADMR Size M (3.75" x 5.75") .....\$115  
☐ EPADLR Size L (4.25" x 6.5") .....\$115

Must select one option:

- ☐ SAHEBR Matrix Elite Swing Away Hardware Right .....\$250  
☐ OSAHR Matrix Offset Elite Swing Away Hardware Right .....\$250

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### HIP AND KNEE SUPPORTS

#### Lateral Hip Support Left

- ☐ HSPXLH Size XS (5.5" W x 4" H) ..... \$120
- ☐ HSPSMLH Size S/M (7" W x 4" H) ..... \$120
- ☐ HSPSLH Size S (7" W x 6" H) ..... \$120
- ☐ HSPMLH Size M (9" W x 4" H) ..... \$120
- ☐ HSPLLH Size L (13" W x 4" H) ..... \$120

#### Lateral Hip Support Right<sup>(1)</sup>

- ☐ HSPXRH Size XS (5.5" W x 4" H) ..... \$120
- ☐ HSPSMRH Size S/M (7" W x 4" H) ..... \$120
- ☐ HSPSRH Size S (7" W x 6" H) ..... \$120
- ☐ HSPMRH Size M (9" W x 4" H) ..... \$120
- ☐ HSPLRH Size L (13" W x 4" H) ..... \$120

#### Lateral Knee Support Left<sup>(1)</sup>

- ☐ HSPXLK Size XS (5.5" W x 4" H) ..... \$120
- ☐ HSPSMLK Size S/M (7" W x 4" H) ..... \$120
- ☐ HSPSLK Size S (7" W x 6" H) ..... \$120
- ☐ HSPMLK Size M (9" W x 4" H) ..... \$120
- ☐ HSPLLK Size L (13" W x 4" H) ..... \$120

#### Lateral Knee Support Right<sup>(1)</sup>

- ☐ HSPXRK Size XS (5.5" W x 4" H) ..... \$120
- ☐ HSPSMRK Size S/M (7" W x 4" H) ..... \$120
- ☐ HSPSRK Size S (7" W x 6" H) ..... \$120
- ☐ HSPMRK Size M (9" W x 4" H) ..... \$120
- ☐ HSPLRK Size L (13" W x 4" H) ..... \$120

Add Removable Hardware:

#### MAXX Style Quick Release, Multi-Axis Mounting

##### Hardware for Hip Supports

- ☐ MHSHLH Left mounted at Hip ..... \$215
- ☐ MHSHLK Left mounted at Knee ..... \$215
- ☐ MHSHRH Right mounted at Hip ..... \$215
- ☐ MHSHRK Right mounted at Knee ..... \$215

##### Swing Away Hardware for Hip Supports

- ☐ SHSHLH Left mounted at Hip ..... \$215
- ☐ SHSHLK Left mounted at Knee ..... \$215
- ☐ SHSHRH Right mounted at Hip ..... \$215
- ☐ SHSHRK Right mounted at Knee ..... \$215

#### Extended Maxx Style Quick Release, Multi-Axis Mounting

##### Hardware for Hip Supports<sup>(2)</sup>

- ☐ EMQHLH Left mounted at Hip ..... \$215
- ☐ EMQHLK Left mounted at Knee ..... \$215
- ☐ EMQHRH Right mounted at Hip ..... \$215
- ☐ EMQHRK Right mounted at Knee ..... \$215

#### Lift Off Hardware for Hip Support

##### (Hip Position Only)

- ☐ LOHL Left ..... \$195
- ☐ LOHR Right ..... \$195

1. Comes standard with Fixed Mount Hardware.

2. Recommended when finished seat cushion thickness will be greater than 3.5" high.

### POWER CENTER MOUNT FOOT PLATFORM

| <b>MUST BE COMPLETED</b>              | <i>Measurement in Inches</i> |
|---------------------------------------|------------------------------|
| <i>Knee to Heel:</i>                  |                              |
| <i>Cushion Thickness: (Subtract)</i>  |                              |
| <i>Cushion Compression: (Add)</i>     |                              |
| <i>Seat Pan to Footplate: (Total)</i> |                              |

#### LNK Power Center Mount Foot Platform

- ☐ LNX 9" Articulation<sup>(1)(2)</sup> ..... \$3,000
  - ☐ LSPTF1 10" - 14" Seat pan to footplate
  - ☐ LSPTF2 13" - 19" Seat pan to footplate
- ☐ LNXA 6" Articulation<sup>(1)(2)(6)</sup> ..... \$3,000
  - ☒ LASPTF2 10" - 14" Seat pan to footplate

##### Options:

- ☐ LIFPN Rubber coated footplate 10.5" W x 10" D ..... STD
- ☐ LICP Individual calf pads ..... STD
- ☐ OPCP One Piece Malleable Calf Panel<sup>(3)</sup> ..... \$250
- ☐ LNXPF Plantar Flex Option for LNX Foot platform<sup>(4)</sup> ..... \$250

#### Latitude Style Foot Platform

- ☐ LSFP Latitude Style Foot Platform<sup>(5)</sup> ..... \$2,995  
*For finished knee to heel measurements of 6.5"-12"*
- ☐ LSPFP Latitude Style Padded Foot Plate ..... \$250

1. 3" ground clearance required.

2. Only available on adult range seat systems 16" deep or greater.

For shorter seat depths call Customer Service for other options.

3. Not available with Plantar Flex Option for LNX Foot platform, or with LNXA.

4. Only available with Foot Platform Options.

5. Not available with TRRO and TRRO-E.

6. For clients under 4' only.

### MANUAL ELEVATING CENTER MOUNT LEGREST

- ☐ MECML Manual Elevating Center Mount Legrest  
With Plantar Flex ..... \$595

Select legrest set up position:

- ☐ 90SU Set at 90° (elevating range 90-40°) ..... STD

##### Rubber Coated Footplate:

- ☐ MRCP Rubber Coated Footplate 10.5" W x 10" D ..... N/C

Select seat pan to footplate for One Piece Foot Plate

- ☐ CSPTF1 11.5" - 15" ..... N/C

- ☐ CSPTF2 14" - 19" ..... N/C

##### Individual Footplates:

- ☐ IFSL Small - (5" x 7.5" deep) - Left ..... \$165
- ☐ IFSR Small - (5" x 7.5" deep) - Right ..... \$165
- ☐ IFML Medium - (5.5" x 9.5" deep) - Left ..... \$165
- ☐ IFMR Medium - (5.5" x 9.5" deep) - Right ..... \$165
- ☐ IFLR Large - (6" x 11.5" deep) - Left ..... \$165
- ☐ IFLR Large - (6" x 11.5" deep) - Right ..... \$165


Select seat pan to footplate for Individual Footplates

- ☐ CSPTF3 9.5" - 14.5" ..... N/C

- ☐ CSPTF4 12" - 19" ..... N/C

##### Optional:

- ☐ CICP Individual calf pads ..... \$200



# Ultra Low Maxx CG Power Positioning System

## Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_


Tag: \_\_\_\_\_

### MANUAL CENTER MOUNT FRONT RIGGING

- ☐ MSFCMFP Maxx Style Fixed Center Mount Foot Platform<sup>(1)(2)</sup> .....\$275
- ☐ MIFPN Comes with a rubber coated footplate 10.5" W x 10" D. ....N/C
- ☐ MICP Individual calf pads .....\$100


Select seat pan to footplate:

- ☐ MSPTF1 9" - 13" .....N/C
- ☐ MSPTF14 13" -19" .....N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
2. 3" ground clearance required.


### FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP



- ☐ SFFP Small Foot Platform (12"w x 8.5"d) .....\$180
- ☐ MFFP Medium Platform (12"w x 10.75"d) .....\$180
- ☐ LFFP Large Platform (14.5"w x 13"d) .....\$250
- ☐ XLFFP Extra Large Platform (17"w x 13"d) .....\$250
- ☐ XXLFFP Extra Extra Large Platform (21"w x 13"d) .....\$250

1. May need to adjust leg position to avoid interference with front casters

### INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP



- ☐ IFSL Small (5" W x 7.5" D) - Left .....\$165
- ☐ IFSR Small (5" W x 7.5" D) - Right .....\$165
- ☐ IFML Medium (5.5" W x 9.5" D) - Left .....\$165
- ☐ IFMR Medium (5.5" W x 9.5" D) - Right .....\$165
- ☐ IFLL Large (6" W x 11.5" D) - Left .....\$165
- ☐ IFLR Large (6" W x 11.5" D) - Right .....\$165


Add:

- ☐ TSSL Left Toe Strap Small 6.5" .....\$125
- ☐ TSSR Right Toe Strap Small 6.5" .....\$125
- ☐ TSML Left Toe Strap Medium 8.5" .....\$125
- ☐ TSMR Right Toe Strap Medium 8.5" .....\$125
- ☐ TSLL Left Toe Strap Large 10.5" .....\$125
- ☐ TSLR Right Toe Strap Large 10.5" .....\$125

- ☐ MFOL Multi-Axis Footplate (6"W x 10"D) - Left .....\$225
- ☐ MFOR Multi-Axis Footplate (6"W x 10"D) - Right .....\$225

1. Foot plates may interfere with casters in certain configurations

### MAXX INDEPENDENT ADJUSTABLE POSITIONING



- 5.5" - 6" WIDTH ADJUSTMENT
- 9.5" - 11.5" DEPTH ADJUSTMENT
- UP TO 9° ADDUCTION
- UP TO 45° ABDUCTION

- ☐ MPFPL Maxx Independent Adjustable Positioning Foot plate<sup>(1)</sup> - Left .....\$200
- ☐ MPFPR Maxx Independent Adjustable Positioning Foot plate<sup>(1)</sup> - Right .....\$200
- ☐ MTSL Left Toe Strap (each) .....\$35
- ☐ MTSR Right Toe Strap (each) .....\$35

1. Foot plates may interfere with casters in certain configurations.

### MAXX STYLE POWER ELEVATING LEGRESTS

- ☐ MSPPL MAXX Swing Away Power Elevating<sup>(1)</sup> .....\$2,490

Select Power Legrest Function:

- ☐ INDIA Individual Movement .....N/C
- ☐ COMA Combined Movement .....N/C
- ☐ SYNA Synchronized/Recline Movement .....N/C

Select Receivers:

- ☐ MRS Straight Receivers .....N/C
- ☐ MRF Flared 1" Outward Receivers .....N/C
- ☐ MRF2 Flared 2" Outward Receivers .....\$250

Select Seat Pan to Footplate:

- ☐ PSPTF5 Length: 10"-14" .....N/C
- ☐ PSPTF2 Length: 13"-17" .....N/C
- ☐ PSPTF7 Length: 16"-20" .....N/C

Optional:

- ☐ PMCP Moulded I-Skin Calf Pads .....STD
- ☐ PFCPMS Flat Calf Pad Option (pair) .....\$275
- ☐ PCCPMS Curved Calf Pad Option<sup>(2)</sup> (pair) .....\$275

1. Must choose footplate on page 13.

2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

### MAXX STYLE MANUAL ELEVATING LEGRESTS

- ☐ MMPPSA MAXX Swing Away Manual Elevating<sup>(1)</sup> .....\$680

Select Receivers:

- ☐ ERS Straight Receivers .....N/C
- ☐ ERF Flared 1" Outward Receivers .....N/C
- ☐ ERF2 Flared 2" Outward Receivers .....\$250

Select Seat Pan to Footplate:

- ☐ ESPTF5 Length: 10"-14" .....N/C
- ☐ ESPTF2 Length: 13"-17" .....N/C
- ☐ ESPTF7 Length: 16"-20" .....N/C

Optional:

- ☐ EMCP Moulded I-Skin Calf Pads .....STD
- ☐ EFCPMS Flat Calf Pad Option (pair) .....\$275
- ☐ ECCPMS Curved Calf Pad Option<sup>(2)</sup> (pair) .....\$275

1. Must choose footplate on page 13.

2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

# Ultra Low Maxx CG Power Positioning System

## Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### SWING AWAY LEGREST

☐ HD70 Heavy Duty 70° Swing Away Footrest(pair)<sup>(1)(2)</sup> . \$495

#### Select Receivers:

☐ RS Straight Receivers .....N/C  
☐ RF Flared 1" Outward Receivers .....N/C  
☐ RF2 Flared 2" Outward Receivers.....\$250

#### Select Seat Pan to Footplate:

☐ HSPTF11 Length: 6"-10" .....N/C  
☐ HSPTF8 Length: 10"-14" .....N/C  
☐ HSPTF2 Length: 12"-16" .....N/C  
☐ HSPTF7 Length: 15"-19" .....N/C

#### Options:

☐ HD70FCP Flat Calf Pad Option (pair) .....\$275  
☐ HD70CCP Curved Calf Pad Option<sup>(3)</sup> (pair) .....\$275

☐ I70SAFT Invacare 70° Swingaway Footrests.....\$325  
☐ I70-CFT Comes with composite foot plates.....N/C

☐ RLSL Residual Limb Support - Left (8.5" x 6")<sup>(4)</sup>  
.....\$495

☐ RLSR Residual Limb Support - Right (8.5" x 6")<sup>(4)</sup>  
.....\$495

1. Must choose footplate on page 13.

2. Heavy duty durable construction (4.5lbs ea.).

3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.





4. Swing away and removable, 50° of angle adjustment, 4.5" height adjustment, 3" horizontal and 2" of for/aft adjustment.

### DUAL FRONT RIGGING ADAPTOR<sup>(1)</sup>

☐ DFRA Dual Front Rigging Receivers.....\$450

1. Must purchase each set of front rigging separately. Not available on no charge front rigging options.

### FOOTPLATES - SWING AWAY/ELEVATING LEGRESTS<sup>(1)</sup>






#### Angle Adjustable Footplates

☐ AFSL Left Small Footplate (4.25"W x 8"D) .....\$165  
☐ AFSR Right Small Footplate (4.25"W x 8"D) .....\$165  
☐ AFML Left Medium Footplate (5.25"W x 8"D) .....\$165  
☐ AFMR Right Medium Footplate (5.25"W x 8"D) .....\$165  
☐ AFLL Left Large Footplate (6.25"W x 8"D) .....\$165  
☐ AFLR Right Large Footplate (6.25"W x 8"D) .....\$165  
☐ AFXL Left Extra Large Footplate (9"W x 11"D) .....\$165  
☐ AFXR Right Extra Large Footplate (9"W x 11"D) .....\$165

#### Multi-Axis Adjustable Angle Footplates

☐ MFSL Left Small Footplate (4.25"W x 8"D) .....\$225  
☐ MFSR Right Small Footplate (4.25"W x 8"D) .....\$225  
☐ MFML Left Medium Footplate (5.25"W x 8"D) .....\$225  
☐ MFMR Right Medium Footplate (5.25"W x 8"D) .....\$225  
☐ MFLL Left Large Footplate (6.25"W x 8"D) .....\$225  
☐ MFLR Right Large Footplate (6.25"W x 8"D) .....\$225  
☐ MFXL Left Extra Large Footplate (9"W x 11"D) .....\$225  
☐ MFXR Right Extra Large Footplate (9"W x 11"D) .....\$225



#### Cut Out Angle Adjustable Footplates

☐ AFCSL Left Small Footplate (7.25"W x 9"D) .....\$165  
☐ AFCSR Right Small Footplate (7.25"W x 9"D) .....\$165  
☐ AFCML Left Medium Footplate (7.5"W x 11"D) .....\$165  
☐ AFCMR Right Medium Footplate (7.5"W x 11"D) .....\$165  
☐ AFCLL Left Large Footplate (10"W x 14"D) .....\$165  
☐ AFCLR Right Large Footplate (10"W x 14"D) .....\$165

#### Cut Out Multi-Axis Adjustable Angle Footplates

☐ MFCSL Left Small Footplate (7.25"W x 9"D) .....\$225  
☐ MFCSR Right Small Footplate (7.25"W x 9"D) .....\$225  
☐ MFCML Left Medium Footplate (7.5"W x 11"D) .....\$225  
☐ MFCMR Right Medium Footplate (7.5"W x 11"D) .....\$225  
☐ MFCLL Left Large Footplate (10"W x 14"D) .....\$225  
☐ MFCLR Right Large Footplate (10"W x 14"D) .....\$225

1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### HEEL LOOPS, CUPS AND TOE STRAPS

|                              |                        |           |
|------------------------------|------------------------|-----------|
| <input type="checkbox"/> HLL | Left Heel Loop (each)  | .....\$20 |
| <input type="checkbox"/> HLR | Right Heel Loop (each) | .....\$20 |
| <input type="checkbox"/> HCL | Left Cup Loop (each)   | .....\$20 |
| <input type="checkbox"/> HCR | Right Cup Loop (each)  | .....\$20 |
| <input type="checkbox"/> TSL | Left Toe Strap (each)  | .....\$20 |
| <input type="checkbox"/> TSR | Right Toe Strap (each) | .....\$20 |

#### CALF STRAPS AND CALF PANELS<sup>(1)</sup>

|                              |                                                                                                 |             |
|------------------------------|-------------------------------------------------------------------------------------------------|-------------|
| <input type="checkbox"/> CSM | Body Point Airmesh Calf Strap for Chairs<br>15"-18" wide with 3.5" H Pad <sup>(1)</sup>         | ..... \$200 |
| <input type="checkbox"/> CSL | Body Point Airmesh Calf Strap for Chairs<br>18"-22" wide with 3.5" H Pad <sup>(1)</sup>         | ..... \$200 |
| <input type="checkbox"/> CPM | Body Point Airmesh Calf Panel for Chairs<br>15"-18" wide with 35.5" W x 9" H Pad <sup>(1)</sup> | .....\$295  |
| <input type="checkbox"/> CPL | Body Point Airmesh Calf Panel for Chairs<br>18"-22" wide with 38.5" W x 9" H Pad <sup>(1)</sup> | ..... \$315 |

1. Not available for Center Mounted Front Rigging.

#### STANDARD FRONT RECEIVERS

|                                 |                                          |          |
|---------------------------------|------------------------------------------|----------|
| <input type="checkbox"/> PHDR   | Pride HD Drop-in Receivers Only          | .....N/C |
| <input type="checkbox"/> PSAR   | Pride Style #8 Swing Away Receivers Only | .....N/C |
| <input type="checkbox"/> IAR    | Invacare 70° Swing Away Receivers Only   | .....N/C |
| <input type="checkbox"/> Q70R   | Quickie 70° Swing Away Receivers Only    | .....N/C |
| <input type="checkbox"/> Q70HDR | Quickie 70HD Swing Away Receivers Only   | .....N/C |

#### SINGLE-DUAL FUNCTION ELECTRONICS

|                                |                                                   |             |
|--------------------------------|---------------------------------------------------|-------------|
| <input type="checkbox"/> SFCB  | Single/Dual Function Electronics <sup>(1)</sup>   | ..... STD   |
| Select Switch:                 |                                                   |             |
| <input type="checkbox"/> SPBSO | Single Push Button                                |             |
| <input type="checkbox"/> STO   | Single Toggle                                     |             |
| <input type="checkbox"/> DPBO  | Dual Push Button                                  |             |
| <input type="checkbox"/> OSO   | Omit Switch                                       |             |
| <input type="checkbox"/> SFCBA | Single/Dual Function Drive Control <sup>(1)</sup> | .....\$750  |
| Select Control:                |                                                   |             |
| <input type="checkbox"/> JOY   | Through Joystick                                  |             |
| <input type="checkbox"/> SC    | Through Specialty Controls                        |             |
| <input type="checkbox"/> OSO   | Reduced Drive for Elevating Seat <sup>(2)</sup>   | ..... \$400 |

1. These electronics include drive lock-out and tilt limit.

2. If reduced drive is not selected the system will be in drive lock-out when elevated.

#### MULTI- FUNCTION ELECTRONICS

|                                           |                                                                |               |
|-------------------------------------------|----------------------------------------------------------------|---------------|
| <input type="checkbox"/> MFCB             | Multi-Function through Switch <sup>(1)</sup>                   | ..... \$1,950 |
| Select Switch:                            |                                                                |               |
| <input type="checkbox"/> QPBO             | Quad Push Button                                               | .....N/C      |
| <input type="checkbox"/> FWT              | Four-Way Toggle Switch                                         | .....\$150    |
| <input type="checkbox"/> EWR              | 8-Way Rocker <sup>(2)</sup>                                    | .....\$205    |
| <input type="checkbox"/> MFCBD            | Multi-Function M610 Through<br>Drive Control <sup>(1)</sup>    | ..... \$1,950 |
| Includes reduced drive for Elevating Seat |                                                                |               |
| Select Control:                           |                                                                |               |
| <input type="checkbox"/> JOY              | Through Joystick                                               |               |
| <input type="checkbox"/> SC               | Through Specialty Controls                                     |               |
| Optional:                                 |                                                                |               |
| <input type="checkbox"/> RDEY             | Include Reduced Drive for Elevating Seat <sup>(3)</sup>        |               |
| <input type="checkbox"/> RDEN             | Do Not Include Reduced Drive for Elevating Seat <sup>(3)</sup> |               |

1. These electronics include attendant control, drive lock-out, tilt/recline limit and actuator speed control.

2. Dual function is not available with 8-way Rocker.

3. If reduced drive is not selected system will be in drive lock-out when elevated.

#### ELECTRONIC SWITCH MOUNTING

|                                |                                             |
|--------------------------------|---------------------------------------------|
| <input type="checkbox"/> ESML  | Left Side (Armpad)                          |
| <input type="checkbox"/> ESMR  | Right Side (Armpad)                         |
| <input type="checkbox"/> ESMJL | Joystick Left (Mounts to Joystick Bracket)  |
| <input type="checkbox"/> ESMJR | Joystick Right (Mounts to Joystick Bracket) |
| <input type="checkbox"/> ESMD  | Do not mount                                |

## Ultra Low Maxx CG Power Positioning System

### Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### ELECTRONICS OPTIONS

##### Additional Switches<sup>(4)</sup>

|                                   |                                                        |       |
|-----------------------------------|--------------------------------------------------------|-------|
| <input type="checkbox"/> SPBSA    | Single Push Button .....                               | \$150 |
| <input type="checkbox"/> STA      | Single Toggle Switch <sup>(3)</sup> .....              | \$150 |
| <input type="checkbox"/> DPBA     | Dual Push Button <sup>(3)</sup> .....                  | \$150 |
| <input type="checkbox"/> QPBA     | Quad Push Button <sup>(1)(3)</sup> .....               | \$150 |
| <input type="checkbox"/> 10WA     | 10-Way Push Button <sup>(1)(3)</sup> .....             | \$255 |
| <input type="checkbox"/> ASL320R  | ASL Pulse Mechanical Switch - Red .....                | \$200 |
| <input type="checkbox"/> ASL320G  | ASL Pulse Mechanical Switch - Green .....              | \$200 |
| <input type="checkbox"/> ASL320BK | ASL Pulse Mechanical Switch - Black .....              | \$200 |
| <input type="checkbox"/> SESG     | Stealth Egg Switch - Green .....                       | \$135 |
| <input type="checkbox"/> SESB     | Stealth Egg Switch - Black .....                       | \$135 |
| <input type="checkbox"/> FTSQ     | Feather Touch - Single Quarter Size .....              | \$295 |
| <input type="checkbox"/> FTMD     | Feather Touch - Dual Dime Size <sup>(1)(3)</sup> ..... | \$295 |
| <input type="checkbox"/> FTMQ     | Feather Touch - Quad Dime Size <sup>(1)(3)</sup> ..... | \$295 |
| <input type="checkbox"/> FWTA     | Four-Way Toggle Switch <sup>(1)(3)</sup> .....         | \$150 |
| <input type="checkbox"/> EWRA     | 8-Way Rocker <sup>(1)(3)</sup> .....                   | \$255 |
| <input type="checkbox"/> TWRA     | 10-Way Rocker <sup>(1)</sup> .....                     | \$255 |
| <input type="checkbox"/> TWTA     | 10-Way Toggle <sup>(1)</sup> .....                     | \$255 |
| <input type="checkbox"/> EWPA     | 8-Way Push Button <sup>(1)</sup> .....                 | \$255 |

##### Specify Use for Additional Switch:

- ☐ STG Seating
- ☐ JM Joystick Mode<sup>(2)</sup>
- ☐ JOO Joystick On/Off<sup>(2)</sup>

##### Specify Mounting for Additional Switch:

- ☐ ESMLA Left Side (Armpad)
- ☐ ESMRA Right Side (Armpad)
- ☐ ESMJLA Joystick Left (Mounts to Joystick Bracket)
- ☐ ESMJRA Joystick Right (Mounts to Joystick Bracket)
- ☐ ESMDA Do not mount

##### Miscellaneous

|                                |                                                                                                  |       |
|--------------------------------|--------------------------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> MJBS  | Motion Junction Box for Switch Control <sup>(1)(5)</sup><br>For Mode and/or Emergency Stop ..... | \$175 |
| <input type="checkbox"/> S46L  | Single to Dual Switch Adapter .....                                                              | \$265 |
|                                | <i>Allows for dual use of a single switch</i>                                                    |       |
| <input type="checkbox"/> M616S | M616 Converter <sup>(1)</sup> .....                                                              | \$125 |
|                                | <i>Allows up to 4 mono phono switches into<br/>Multi-Function box</i>                            |       |
| <input type="checkbox"/> LJS   | Lazarus Jump Starter .....                                                                       | \$150 |
| <input type="checkbox"/> APPS  | Accessory Port Power Supply <sup>(1)</sup><br>3 amp capacity .....                               | \$195 |
| <input type="checkbox"/> TPSMH | 12 Piece Switch Mounting Hardware<br>Kit (additional kits) .....                                 | \$39  |

1. Not available with Single-Dual Function Electronics.

2. IN200 is required with REM210, REM211, REM216 joysticks.

3. Not available with REM110 joystick.

4. Not available for Joystick Mode or Joystick On/Off. Must Select Seating (STG)

5. Cannot be used for power On/Off.

#### ARTICULATING TRAYS<sup>(1)</sup>

|                               |                                                   |         |
|-------------------------------|---------------------------------------------------|---------|
| <input type="checkbox"/> USAT | Universal Style Articulating Tray .....           | \$1,350 |
|                               | <i>14"W x Adjustable 6"-10"D</i>                  |         |
| Add:                          |                                                   |         |
| <input type="checkbox"/> TMA  | Trilogy Mounting Adapter .....                    | \$50    |
| <input type="checkbox"/> AVT  | Articulating Vent (LTV) Tray .....                | \$1,350 |
| <input type="checkbox"/> AVBH | Articulating Vent Bag Holder <sup>(2)</sup> ..... | \$1,350 |
| Add:                          |                                                   |         |
| <input type="checkbox"/> ABMC | Accessory Bag Mounting and Clips .....            | \$100   |

1. The user weight limit is 225 lbs. The maximum weight capacity of the vent tray is 35 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. Not available with TRRO and TRRO-E.

2. Motion Concepts does not supply the Vent Bag.

#### GIMBAL VENT TRAYS<sup>(1)</sup>

|                                |                                             |         |
|--------------------------------|---------------------------------------------|---------|
| <input type="checkbox"/> GUT   | Gimbal Universal Tray .....                 | \$1,350 |
|                                | <i>14"W x Adjustable 6"-10"D</i>            |         |
| Add:                           |                                             |         |
| <input type="checkbox"/> GTMA  | Trilogy Mounting Adapter .....              | \$50    |
| <input type="checkbox"/> GVT   | Gimbal Vent (LTV) Tray .....                | \$1,350 |
| <input type="checkbox"/> GVBH  | Gimbal Vent Bag Holder <sup>(2)</sup> ..... | \$1,350 |
| Add:                           |                                             |         |
| <input type="checkbox"/> GABMC | Accessory Bag Mounting and Clips .....      | \$100   |

1. The Gimbal Vent Trays are not available with Tilt Only and/or Elevate Only Systems. The user weight limit is 200 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. The maximum weight capacity of the vent tray is 20 lbs. Not available with TRRO and TRRO-E.

2. Motion Concepts does not supply the Vent Bag.

#### O2 HOLDERS<sup>(1)</sup>

|                                 |                                                                  |         |
|---------------------------------|------------------------------------------------------------------|---------|
| <input type="checkbox"/> NAO2H  | Non-Articulating O2 Gas Holder .....                             | \$395   |
| <input type="checkbox"/> AOHG   | Articulating O2 Gas Holder Cylinder Style <sup>(2)(3)</sup> ..   | \$1,250 |
| <input type="checkbox"/> AOHL   | Articulating Liquid O2 Holder (6.5" x 4.5") <sup>(2)(3)</sup> .. | \$1,250 |
| <input type="checkbox"/> NAIP L | Non-Articulating IV Pole User's Left <sup>(4)</sup> .....        | \$395   |
| <input type="checkbox"/> NAIP R | Non-Articulating IV Pole User's Right <sup>(4)</sup> .....       | \$395   |
| <input type="checkbox"/> AIVP   | Articulating IV Pole <sup>(2)(3)(4)</sup> .....                  | \$1,250 |

1. The user weight limit is 225 lbs. The maximum weight capacity of the O2 holder is 20 lbs. Not available with TRRO and TRRO-E.

2. Required for use on Recline Systems

3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems or with Single Post Backs.

4. The maximum weight capacity of the IV pole is 5 lbs.

#### VENTILATOR INFORMATION

Required when ordering a Vent Tray

- ☐ BIB Built in battery
- ☐ DTSVPS Dealer to supply vent power source

# Storm Torque RX20

Tag: \_\_\_\_\_

## NOTES

1. Not available with TRRO, TRRO-E brackets
2. Not available on joystick side and only available with a flat arm pad or modular arm pad. Not available with Outback Cantilever Arms.
3. The maximum weight capacity of the tray is 10 lbs.
4. Not available on the joystick side.

|                                         |                                                      |
|-----------------------------------------|------------------------------------------------------|
| <input checked="" type="checkbox"/> MCI | Motion Concepts Installation <sup>(1)</sup> .....N/C |
|-----------------------------------------|------------------------------------------------------|

1. Batteries may be required if items on Accessories order form are selected at time of installation, please call to confirm.

**RETAIL \$** \_\_\_\_\_

CLEAR FORM