



Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

SAVE

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

SUBMIT

CLEAR FORM

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

905-695-0583

PRINT FORM



**TURNING,
DRIVING
& RIDE
COMFORT!**

STORM TORQUE RX20 BASE WITH ULTRA LOW MINI POWER POSITIONING SYSTEM

CDN PRICE LIST AND ORDER FORM

Price Effective July 15, 2026



PRINT FORM

SAVE

CLEAR FORM

REQUIRED INFORMATION

Quote ☐ Order ☐

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured
please fill in all REQUIRED* measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness	



Customer Service: 1.800.748.7943 | Fax: 1.905.695.0583 | www.motionconcepts.com

INVACARE POWER BASES

NEW!

☐ SRX-20 Storm Torque RX20\$7,575

Select STF⁽¹⁾: ☐ LOWP 18.125" ☐ MEDP 18.625" ☐ TALLP 19.625"
(Base comes standard with LED Lights)

☐ U175 Weight Capacity up to 175 lbsSTD

☐ NARROWFRM Narrow Base 22.4" (uses 2 M50 Batteries)N/C

☐ STDFRM Standard Base 24.4" (uses 2 Group 24 Batteries)N/C



-Enhanced Front Caster Housing and Bearings
-New Front Caster Design
-Improved double-sided Front Forks
-Upgraded Rear Suspension
-Optimized Linx Parameters

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size.

PERFORMANCE MOTORS

☐ STDMTR Standard Motor - 5.8mphSTD

☐ HSMTR High Speed 4-Pole Motor - 7.5 mph⁽¹⁾⁽²⁾⁽³⁾\$995

- Not available on the NARROWFRM.
- Extended high speed use on rugged terrain is not recommended.
- Performance/High Speed Motors may have less torque and may have less range than standard speed motors.

TRANSPORT TIE DOWN

☐ TRBKTS Wheelchair Transport Brackets⁽¹⁾STD

1. For unoccupied use only

CASTER OPTION

NEW!

☐ AB1228-3 8" x 2.5" Black Casters w/Solid InlayN/C

TIRE OPTION

☐ AB1430-3 14" x 3" Black Tires PneumaticN/C

☐ AB1431-3 14" x 3" Black Tires w/Gel Foam InsertsN/C

BATTERY CHARGER OPTION

☐ 110CHARGER 110 Volt Battery ChargerN/C

☐ OMIT Omit Battery ChargerN/C

INSTALLATION OF BATTERIES

☐ M24 SLD G FT MK Group 24 Battery (On Chair)\$840

☐ M50-12 SLD M MK Group 50 Battery (On Chair)\$650

FRAME FINISH

☐ 170P	Obsidian Black	N/C
☐ 171P	Ruby Red	N/C
☐ 172P	Peacock Blue	N/C
☐ 173P	Pink Flamingo	N/C
☐ 174P	Purple Passion	N/C
☐ 175P	Denim Navy	N/C
☐ 176P	Lime Twist	N/C
☐ 177P	Arctic Blue	N/C

RIM INSERT COLORS - MUST PICK ONE

☐ 170PR	Obsidian Black	N/C
☐ 171PR	Ruby Red	N/C
☐ 172PR	Peacock Blue	N/C
☐ 173PR	Pink Flamingo	N/C
☐ 174PR	Purple Passion	N/C
☐ 175PR	Denim Navy	N/C
☐ 176PR	Lime Twist	N/C
☐ 177PR	Arctic Blue	N/C
☐ 178PR	Banana Split	N/C
☐ 179PR	Orange Marmalade	N/C
☐ 180PR	Morning Mist	N/C

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CONTROLLERS

- ☐ EXPC Expandable Controller STD
- ☐ PWH Harness Required for Expandable System N/C
- ☐ GTRAC-LX LINX G-TRAC Module \$1,200

DRIVER CONTROL

- ☐ REM215 LED Non-Expandable Remote/Drive Only/Joystick⁽¹⁾ STD

1. Controls up to 2 actuators. Must choose Single/Dual function through drive control and light operation.

- ☐ REM216 LED Expandable Remote/Joystick⁽¹⁾ \$875
Three drive functions, control of up to 6 actuators through the drive control and light operation.

- ☐ REM400 Color 3.5" Touch Screen
Remote/Joystick⁽¹⁾⁽²⁾ \$1,100

- ☐ SW400TGL Toggle Kit for REM400⁽³⁾ \$200

- ☐ REM500 LINX Color 3.5" Touch Screen
Display Only - No Driver Control⁽¹⁾⁽²⁾ \$1,100

Select mounting position

- ☐ R5MR Right
☐ R5ML Left

1. Must choose EXPC Expandable Controller and PWH Harness.
2. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.
3. Adds 2 toggle switches to REM400.

ELECTRONICS - EXPANDABLE MODULES

- ☐ OUT500 Output Module \$645
- ☐ IN200 LINX Switch Input Module with
Power ON/OFF for REM200 Series⁽¹⁾ \$645
- ☐ IN500 LINX 9 Pin Input Module
with Sip-N-Puff Nozzle⁽²⁾⁽³⁾⁽⁴⁾ \$1,800
- ☐ PKG32666 Therafin Sip-N-Puff Breath Tube Kit \$440

1. Not for use with Alternative Drive Controls.
2. Only one input module per chair, requires REM400 or REM500.
3. Input module also required for any 9-pin specialty control.
4. 9 pin specialty control and Sip-N-Puff cannot both be used together.

MOUNTING FOR DRIVER CONTROLS

- ☐ FR Fixed Joystick Mounting - Right STD
- ☐ FL Fixed Joystick Mounting - Left STD
- ☐ SAQR Motion Height Adjustable Swing-Away
Quad Link (Right Hand)⁽¹⁾ \$244
- ☐ SAQL Motion Height Adjustable Swing-Away
Quad Link (Left Hand)⁽¹⁾ \$244
- ☐ MRR Maxx Resolve Swing-Away Joystick Mount
(Right Hand)⁽¹⁾ \$450
- ☐ MRL Maxx Resolve Swing-Away Joystick Mount
(Left Hand)⁽¹⁾ \$450

1. Not available with the REM500.

ATTENDANT DRIVE CONTROL CHOICES

- ☐ ACU Proportional Attendant Control⁽¹⁾ \$1,024

Select mounting position

- ☐ RIGHT-ACU Right
☐ LEFT-ACU Left

1. Requires Expandable Controller.

SPECIALTY PROPORTIONAL CONTROLS REQUIRES REM400 OR REM500

- ☐ CREM LINX Compact Remote \$1,097
- ☐ CREM-LF LINX Compact Remote Low Force \$1,097

Select mounting position

- ☐ MR Right
☐ ML Left

LINX ELECTRONIC ACCESSORIES

- ☐ LAK LINX Access Key⁽¹⁾ \$50
- ☐ MUSB Motion Dual USB Charger N/C

1. Required for programming

JOYSTICK TOPS

- ☐ PC101A Bodypoint U Shaped Handle 3"⁽¹⁾ \$120
- ☐ PC102A Bodypoint U Shaped Handle 4"⁽¹⁾ \$120
- ☐ PC107A Bodypoint Rubber Dome⁽¹⁾ \$120
- ☐ PC110A Grooved Mushroom Joystick Handle
3/16" Stem⁽¹⁾ \$120
- ☐ 1560 T Handle Flexible Joystick
Extension⁽¹⁾ \$120
- ☐ 1561 Straight Handle Flexible Joystick
Extension⁽¹⁾ \$120
- ☐ 1826 Chin Cup⁽¹⁾ \$120

1. To remove/disengage the joystick knob from the REM400 Remote pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

Ultra Low Mini Power Positioning System on
Invacare Storm Torque RX20 Base - CDN Price List

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

REHAB SEAT

☐ PED-RS Mini Pediatric Rehab Seat \$3,995

1. Not available with 12" Elevating Seat Lift Module, Order Ultra Low Maxx Elevating Lift Seat Only.

POWER POSITIONING SYSTEMS

☐ MDCGT-50 Ultra Low Mini Maxx CG Tilt (50°) Only \$6,969

☐ MDCGT-45 Ultra Low Mini Maxx CG Tilt (45°) Only⁽¹⁾ \$6,969

☐ MDCGTRESR-50 Ultra Low Mini Maxx CG Tilt (50°) & Power Recline^(163°) w/Extended Shear Reduction(ESR) \$11,569

☐ MDCGTRESR-45 Ultra Low Mini Maxx CG Tilt (45°) & Power Recline^(163°) w/Extended Shear Reduction(ESR)⁽¹⁾ \$11,569

☐ MDRESR Ultra Low Mini Maxx Power Recline (163°) w/Extended Shear Reduction(ESR) \$10,350

☐ MDES12 Ultra Low Mini Maxx Elevate Lift Seat Only \$4,995
Weight capacity up to 175lbs

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module

ADD: ELEVATING LIFT SEAT MODULE

☐ ESM12 12" Elevating Lift Seat Module⁽¹⁾ \$2,495
Weight capacity up to 175lbs

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.

MINI MODULES MUST SELECT SYSTEM TO ADD MODULES

☐ EPRETILT 5° Electronic Pre-Tilt⁽¹⁾⁽²⁾ N/C

☐ ERSSP Programable Electronic Recline Function Start & Stop⁽²⁾ Specify Start ____degrees & Stop ____degrees \$500

☐ PRM8 8° Precline Module⁽³⁾ only available with Recline N/C

☐ PRM30 30° Precline Module⁽³⁾ only available with Recline \$500

1. This Electronic Pre-tilt option adds approximately 2" to the front seat to floor height
2. Requires Multi-Function Electronics through Drive Control (MFCBD). Set range between 90° and 163°. Will be set to factory settings unless specified.
3. Total degrees of recline is 163° ie. if 30° of precline is used, then 133° of recline will be available.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SEAT WIDTH SETTINGS

☐ MSW11 11" Wide (Adjustable from 11" to 15") N/C

☐ MSW12 12" Wide (Adjustable from 11" to 15") N/C

☐ MSW13 13" Wide (Adjustable from 11" to 15") N/C

☐ MSW14 14" Wide (Adjustable from 11" to 15") N/C

☐ MSW15 15" Wide (Adjustable from 11" to 15") N/C

ULTRA RAIL SEAT DEPTH SETTINGS

☐ MSD11 11" Deep⁽¹⁾ (Adjustable from 11" to 15") N/C

☐ MSD12 12" Deep (Adjustable from 11" to 15") N/C

☐ MSD13 13" Deep (Adjustable from 11" to 15") N/C

☐ MSD14 14" Deep (Adjustable from 11" to 15") N/C

☐ MSD15 15" Deep (Adjustable from 11" to 15") N/C

1. For seat depths of 11" weight capacity is 150lbs.

MATRIX SEATING OPTIONS

☐ Matrix Cushions available on page 8

☐ Matrix Backs available on page 8

☐ Matrix Headrests available on page 9

matrix®
SEATING SERIES

LAP BELTS

☐ BELT35 Push Button Style Lap Belt, Length 35" STD

☐ BELT48 Push Button Style Lap Belt, Length 48" N/C

☐ MBPHB Mini BP Padded Hip Belt 2 point⁽¹⁾ \$175

1. Adjustable from 14" to 41" in length, pad is 13" x 1.5" in length

CHEST STRAP

☐ MFUXS BP Monoflex, Underarm RLS - XS^(2.75"x12.25") \$195

☐ MFUS BP Monoflex, Underarm RLS - S^(3"x15.5") \$195

☐ MFCXS BP Monoflex, Centre RLS - XS^(2.75"x15.25") \$195

☐ MFCSBP Monoflex, Centre RLS - S^(3"x15.5") \$195

☐ MINIFLEX-XS Miniflex Abdominal Support Band 4.5" H for Matrix Mini backs 11"-14" wide \$245

☐ MINIFLEX-S Miniflex Abdominal Support Band 6" H for Matrix Mini backs 11"-14" wide \$245

MOTION CONCEPTS HEADREST OPTIONS

☐ MCSH Motion Concepts Standard⁽¹⁾ N/C

☐ MCOH Motion Concepts Auto Style⁽¹⁾ \$235

☐ FM Fixed Mounting (non-removable)⁽²⁾ N/C

☐ OHR Omit: Headrest N/C

1. Must select Fixed Mounting (FM) or Matrix Elan Headrest & Hardware on page 9.
2. Only compatible with Standard Rehab, Contour or High Back Models. All other upgradeable backs require removeable hardware.

MATRIX HEADRESTS OPTIONS

☐ Matrix Elan Headrests and Hardware available on page 9

matrix®
SEATING SERIES

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

BACK CANES (NON RECLINE SYSTEMS)⁽¹⁾

☐ SBC14	14" Straight Back Canes	STD
☐ SBC16	16" Straight Back Canes	STD
☐ SBC18	18" Straight Back Canes	STD
☐ SBC20	20" Straight Back Canes	STD

Select Cane Angle

☐ 75°	☐ 79°	☐ 85°	☐ 90°	☐ 95°	☐ 101°	☐ 106°	☐ 112°	☐ 116°
BA75	BA79	BA85	BA90	BA95	BA101	BA106	BA112	BA116

1. Canes will be set at 90° angle and are adjustable from 75° to 116°

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ FOR REHAB, TILT AND /ELEVATE ONLY

☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"
FBT14	FBT15	FBT16	FBT17	FBT18	FBT19	FBT20
☐ 21"	☐ 23"	☐ 24"	☐ 25"	☐ 22"		
FBT21	FBT22	FBT23	FBT24	FBT25		

- 4" Gap from the seat pan to the bottom of back pan is recommended.
- Finished back height is measured from seat pan to top of backrest.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ FOR RECLINE OR TILT/RECLINE

☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"
FBR16	FBR17	FBR18	FBR19	FBR20	FBR21	FBR22
☐ 23"	☐ 24"	☐ 25"				
FBR23	FBR24	FBR25				

- 6" Gap from the seat pan to the bottom of back pan is recommended.
- Finished back height is measured from seat pan to top of backrest.

BACK OPTIONS

☐ SRB	Standard Rehab Back & Cushion	\$125
☐ PIPF	Flat Planar/I-Back Interface Plate ⁽¹⁾⁽²⁾⁽⁴⁾	\$125
☐ PIPR	Recessed Planar/I-Back Interface Plate ⁽¹⁾⁽²⁾⁽⁴⁾	\$125
☐ ORBC	Omit Rigid Back & Cushion ⁽³⁾	N/C

- Not available on tilt only or tilt/elevate systems, for use with Recline Module.
- Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
- Not available with Recline/ESR.
- Fits Comfort Company, Freedom, Pindot and Moulded Backs.

MATRIX BACKREST OPTIONS

- ☐ Matrix Backs available on page 6/7

matrix
SEATING SERIES

ARMRESTS (NON RECLINE SYSTEMS)

☐ PDP1	Pediatric Dual Post Adj.Height 6.5"-9.5"	STD
☐ PDP2	Pediatric Dual Post Adj.Height 9.75"-12.75"	STD
☐ PDP3	Pediatric Dual Post Adj. Height 12.75"-15.75"	STD

RECLINING ARMRESTS

☐ PDR1	Pediatric Reclining Adj.Height 7"-8.25"	STD
☐ PDR2	Pediatric Reclining Adj.Height 9.5"-13"	STD
☐ PDR3	Pediatric Reclining Adj. Height 13"-16"	STD

ARMREST TUBE LENGTH - PLEASE SPECIFY

☐ ATFL	Left Full	STD
☐ ATDL	Left Desk	STD
☐ ATFR	Right Full	STD
☐ ATDR	Right Desk	STD

ARM PADS

Standard

☐ SAPFL	Left Full (2.25" x 14")	STD
☐ SAPDL	Left Desk (2.25" x 10")	STD
☐ SAPFR	Right Full (2.25" x 14")	STD
☐ SAPDR	Right Desk (2.25" x 10")	STD

Waterfall Arm Pads

☐ WAPFL	Left Full (13")	N/C
☐ WAPDL	Left Desk (9")	N/C
☐ WAPFR	Right Full (13")	N/C
☐ WAPDR	Right Desk (9")	N/C

Modular Arm Pad - Left⁽¹⁾

☐ MAPTFL	Left Tray Full (14")	\$95
☐ MAPTDL	Left Tray Desk (10")	\$95

Select Pad Insert below:

☐ ISSFL	Iskin - Integrated soft skin foam
☐ SCVFL	Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

☐ MAPTFR	Right Tray Full (14")	\$95
☐ MAPTDR	Right Tray Desk (10")	\$95

Select Pad Insert below:

☐ ISSFR	Iskin - Integrated soft skin foam
☐ SCVFR	Startex Covered Visco foam

Flat Multi Position Arm Pads

☐ MPAFL	Left Full (14")	\$150
☐ MPADL	Left Desk (10")	\$150
☐ MPAFR	Right Full (14")	\$150
☐ MPADR	Right Desk (10")	\$150

Gel Pads

☐ GAPNL	Narrow Left (2" x 12")	\$88
☐ GAPNR	Narrow Right (2" x 12")	\$88
☐ GAPWL	Wide Left (3.5" x 12")	\$96
☐ GAPWR	Wide Right (3.5" x 12")	\$96
☐ GAPLL	Long Left (3.5" x 14")	\$103
☐ GAPLR	Long Right (3.5" x 14")	\$103

☐ OMIT	Omit Armpads(7/8" diameter tube)	(\$25)
-------------------------------	----------------------------------	--------

1. Requires interchangeable pad insert.

ARM ADDUCTOR PAD⁽¹⁾

☐ AADSL	Left Small (6" x 2.25")	\$135
☐ AADSR	Right Small (6" x 2.25")	\$135
☐ AADLL	Left Large (8.5" x 2.25")	\$135
☐ AADLR	Right Large (8.5" x 2.25")	\$135

1. Only available with Modular or Flat Multi Position Arm Pads.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

POWER CENTRE MOUNT FOOT PLATFORM

LNx Power Center Mount Foot Platform⁽¹⁾ HCPCS code E1012


☐ LNxA	10" - 14" Seat pan to footplate	\$3,000
Options		
☐ LIFPN	Comes with a rubber coated footplate (10.5" W x 10")	STD
☐ LICP	Comes standard with individual calf pads.	STD
☐ OPCP	One Piece Calf Panel ⁽²⁾	\$250
☐ LNXPF	Plantar Flex Option for LNx Foot platform ⁽¹⁾	\$325
☐ LSFP	Latitude Style Foot Platform.....	\$2,995
	For Finished Knee to heel measurements of 6.5"-12"	
☐ LSPFP	Latitude Style Padded Foot platform	\$250

1. Only available with Foot Platform Options.

2. Not available with TRRO and TRRO-E.


STANDARD FRONT RIGGING

☐ FCMP	Fixed Centre Mount Platform	\$275
☐ IFP	Comes with a rubber coated footplate (10.5" W x 10")	N/C
Select seat pan to footplate:		
☐ SPTF12	6" - 13"	N/C



MANUAL CENTRE MOUNT FRONT RIGGING

☐ MSFCMFP	Maxx Style Fixed Centre Mount Foot Platform ⁽¹⁾	\$275
	With Quick Release Height Adjustment	
☐ IFP	Comes with a rubber coated footplate (10.5" W x 10")	N/C
☐ ICP	Individual calf pads (pair)	\$100
Select seat pan to footplate:		
☐ SPTF13	8" - 13"	N/C
☐ SPTF14	13" -19"	N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.

FOOT PLATFORM OPTIONS FOR MAXX CMFP



☐ SFFP	Small Foot Platform ⁽¹⁾ (12"w x 8.5"d)	\$180
☐ MFFP	Medium Platform ⁽¹⁾ (12"w x 10.75"d)	\$180

1. May need to adjust leg position to avoid interference with front castors

SWING AWAY FOOTRESTS

☐ HD70	Heavy Duty 70° Swing Away Footrest(pair) ⁽¹⁾⁽²⁾	\$495
Select Receivers:		
☐ RS	Straight Receivers	N/C
☐ RF	Flared 1" Outward Receivers.....	N/C
Select Seat Pan to Footplate:		
☐ SPTF11	Seat pan to footplate 6"-10"	N/C
☐ SPTF8	Seat pan to footplate 11"-15"	N/C
☐ SPTF2	Seat pan to footplate 13"-17"	N/C
☐ HDFCP	Flat Calf Pad Option (pair)	\$275

1. Must choose footplate on page 6


2. Heavy duty durable construction (4.5lbs ea.).

SWING AWAY FOOTRESTS

☐ I70SAF	Invacare 70° Swingaway Footrests (Pair)	\$325
☐ I70-CF	Comes with composite foot plates	N/C
☐ 90SWING	Invacare 90° Swingaway Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards ⁽¹⁾ (Range 4.5"-7.5")(Pair)	\$450
Select Seat Pan to Footplate⁽¹⁾:		
☐ SPTF11	Seat pan to footplate 6"-10"	N/C
☐ SPTF8	Seat pan to footplate 11"-15"	N/C
☐ SPTF2	Seat pan to footplate 13"-17"	N/C

1. For optional Seat pan to footplate extensions on 90SWING you must order with Motion Concepts Foot plates

FOOTPLATE SIZE SELECTION



1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ADJUSTABLE ANGLE FOOTPLATES

☐ AFSL	Left Small Footplate (4.25"W x 8"D)	\$165
☐ AFSR	Right Small Footplate (4.25"W x 8"D)	\$165
☐ AFML	Left Medium Footplate (5.25"W x 8"D)	\$165
☐ AFMR	Right Medium Footplate (5.25"W x 8"D)	\$165
☐ AFLL	Left Large Footplate (6.25"W x 8"D)	\$165
☐ AFLR	Right Large Footplate (6.25"W x 8"D)	\$165

MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES

☐ MFSL	Left Small Footplate (4.25"W x 8"D)	\$225
☐ MFSR	Right Small Footplate (4.25"W x 8"D)	\$225
☐ MFML	Left Medium Footplate (5.25"W x 8"D)	\$225
☐ MFMR	Right Medium Footplate (5.25"W x 8"D)	\$225
☐ MFL	Left Large Footplate (6.25"W x 8"D)	\$225
☐ MFLR	Right Large Footplate (6.25"W x 8"D)	\$225

HEEL LOOPS, CUPS AND TOE STRAPS

☐ HLL	Left Heel Loop (each)	\$20
☐ HLR	Right Heel Loop (each)	\$20
☐ HCL	Left Cup Loop (each)	\$20
☐ HCR	Right Cup Loop (each)	\$20
☐ TSL	Left Toe Strap (each)	\$20
☐ TSR	Right Toe Strap (each)	\$20

SINGLE-DUAL FUNCTION ELECTRONICS

☐ SFCB	Single/Dual Function Electronics ⁽¹⁾	STD
-------------------------------	---	-----

Select Switch:

- ☐ SPBO Single Push Button
☐ STO Single Toggle
☐ DPBO Dual Push Button
☐ SOS Omit Switch
☐ JOY Through Joystick
☐ SC Through Specialty Control

Optional:

☐ RDE	Reduced Drive for Elevating Seat ⁽²⁾	\$400
☐ SPBSA	Additional Single Push Button	\$150
☐ STA	Additional Single Toggle Switch	\$150
☐ DPBA	Additional Dual Push Button	\$150

1. These electronics include drive lock-out and tilt limit.

2. If reduced drive is not selected the system will be in drive lock-out when elevated.

MULTI-FUNCTION ELECTRONICS

☐ MFCB	Multi-Function through Switch/Toggle ⁽¹⁾	\$1,950
Select Switch:		
☐ QPBO	Quad Push Button	
☐ FWTO	Four-Way Toggle Switch	
☐ EWRO	8-Way Rocker ⁽²⁾	

☐ MFCBD	Multi-Function through Drive Control ⁽¹⁾	\$1,950
Select Control:		
☐ JOY	Through Joystick	
☐ SC	Through Specialty Control	

Additional Switches

☐ QPBA	Additional Quad Button	\$150
☐ 10WA	Additional 10-Way Push Button	\$255

Additional Switches ⁽⁴⁾ (not available with ACT4 Electronics)

☐ FWTA	Four-Way Toggle Switch ^{(1)(3) (4)}	\$150
☐ EWRA	8-Way Rocker ^{(1)(2) (4)}	\$255
☐ TWRA	10-Way Rocker ^{(1)(2) (4)}	\$255
☐ TWTA	10-Way Toggle ^{(1)(2) (4)}	\$255
☐ EWPA	8-Way Push Button ^{(1)(2) (4)}	\$255

1. These electronics include attendant control, drive lock-out, tilt/recline limit and actuator speed control

2. Dual function is not available.

3. If reduced drive is not selected system will be in drive lock-out when elevated.

4. Not available on Multi Function ACT4-LiNX Electronics

ELECTRONIC SWITCH MOUNTING

☐ ESML	Left Side (Armpad)
☐ ESMR	Right Side (Armpad)
☐ ESMJL	Joystick Left (Mounts to Joystick Bracket)
☐ ESMJR	Joystick Right (Mount to Joystick Bracket)
☐ ESMD	Do not mount

Ultra Low Mini Power Positioning System on Invacare Storm Torque RX20 Base - CDN Price List

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ELECTRONICS OPTIONS

- | | | | |
|---|----------|--|-------|
|  | M616S | M616 Converter ⁽¹⁾ | \$125 |
|  | FTSQ | Feather Touch Switch Single Function
Quarter Size | \$295 |
|  | FTMD | Feather Touch Switch Multi-Function
Dual Dime Size | \$295 |
|  | FTMQ | Feather Touch Switch Multi-Function
Quad Dime Size | \$295 |
|  | ASL320R | ASL Pulse Mechanical Switch Red | \$200 |
|  | ASL320G | ASL Pulse Mechanical Switch Green | \$200 |
|  | ASL320BK | ASL Pulse Mechanical Switch Black | \$200 |
|  | PSOS | Power Shut Off Switch - Toggle ⁽¹⁾ | \$175 |
|  | LJS | Lazarus Jump Starter ⁽¹⁾ | \$150 |
|  | APPS | Accessory Port Power Supply ⁽¹⁾
(3 amp capacity) | \$195 |
|  | TPSMH | 12 Piece Switch Mounting Hardware
Kit (additional kits) | \$39 |

1. Not available with Single-Dual Function Electronics.

GIMBAL VENT TRAYS⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾

- | | | |
|-------------------------------------|---|---------|
| ☐ GUT | Gimbal Universal Tray..... | \$1,350 |
| | (14"W x Adjustable 6"-10"D) | |
| Add: ☐ GTMA | Trilogy Mounting Adapter..... | \$50 |
| ☐ GVT | Gimbal Vent (LTV)Tray | \$1,350 |
| ☐ GVBH | Gimbal Vent Bag Holder ⁽⁵⁾ | \$1,350 |
| Add: ☐ GABMC | Accessory Bag Mounting and Clips | \$100 |

1. The Gimbal Vent Trays are not available with Tilt Only Systems.
2. The user weight limit is 150 lbs.
3. Vents are not standard configuration on Elevating or Precline Systems. Please call Customer Service for RFQ Quotation.
4. The maximum weight capacity of the vent tray is 20 lbs.
5. Motion Concepts does not supply the Vent Bag.

ADDITIONAL OPTIONS & ACCESSORIES

- | | | |
|---|--|-------|
|  THS | Transfer Handles 4" - 6.5" | \$350 |
|  THL | Transfer Handles 9" - 11.5" | \$350 |
|  FFLTL | Flip and Fold Away 1/2 Lap Tray
Left Mounting(14"W x 10"D) ⁽¹⁾⁽²⁾ | \$795 |
|  FFLTR | Flip and Fold Away 1/2 Lap Tray
Right Mounting(14"W x 10"D) ⁽¹⁾⁽²⁾ | \$795 |

1. Not available on joystick side and only available with a flat arm pad or modular arm pad. For seat widths 14" and greater.
2. The maximum weight capacity of the tray is 10 lbs.

INSTALLATION

- ☒ MCI Installation by Motion Concepts N/C

PLEASE SPECIFY PAGES INCLUDED

- | | |
|--|--|
|  Pages 1-7 | Excludes Matrix and Maxx Seating Accessories |
|  Pages 1-10 | Complete Order Form |

NOTES

[illegible]

Total Retail Value \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CUSHION COVER REFERENCE GUIDE

Cushion	Kid*ab*ra	Meshtex	Startex Spacetex	Pink Lollipop Blue Blizzard	Infection Control
PS/PSP	N/C	N/A	STD	STD	Reversible
Kid*ab*ra	STD	\$103	N/A	\$103	N/A
Vi	\$103	STD	\$103	\$103	N/A
Libra/Adjust	N/A	N/A	STD	N/A	\$50

MATRIX SEAT CUSHION OPTIONS

Matrix Libra Adjust Seat Cushion 14" - 15" wide \$985

☐ LC1414-ADJ ☐ LC1416-ADJ ☐ LC1516-ADJ

Matrix Libra Seat Cushion 14" - 15" wide \$745

☐ LC1414 ☐ LC1416 ☐ LC1516

Fabric Options for Libra/Libra Adjust Cushions

☐ SFLF Spacetex N/C

☐ ICFL Infection Control N/C

Matrix PS Seat Cushion 11" - 14" wide \$569

☐ PS1112 ☐ PS1114 ☐ PS1212 ☐ PS1214 ☐ PS1312 ☐ PS1314..
☐ PS1412 ☐ PS1414

Matrix PSP Seat Cushion 11" - 14" wide \$645

☐ PSP1112 ☐ PSP1114 ☐ PSP1212 ☐ PSP1214 ☐ PSP1312
☐ PSP1314 ☐ PSP1412 ☐ PSP1414

Fabric Options for PS/PSP Cushions

☐ CV-K Kid*ab*ra (Non Reversible) N/C

☐ CV-P Pink Lollipop (Non Reversible) N/C

☐ CV-B Blue Blizzard (Non Reversible) N/C

Matrix Kid*ab*ra Cushion 12" - 14" wide \$576
(Comes standard with Kid*ab*bra fabric)

☐ KA1212-VI ☐ KA1214-VI ☐ KA1412-VI ☐ KA1414-VI

Fabric Options for Kidabra Cushions

☐ CV-M Meshtex (Non Reversible) \$103

☐ CV-P Pink Lollipop (Non Reversible) \$103

☐ CV-B Blue Blizzard (Non Reversible) \$103

Matrix Vi Seat Cushion 14" wide \$624

☐ MA1415-VI

BACK COVER REFERENCE GUIDE

Back	Kidabra	Meshtex	Startex	Pink Lollipop	Blue Blizzard	Infection Control
Mini	STD	N/A	N/C	N/C	N/C	N/A
Elite / Elite TR	N/A	N/A	N/A	N/A	N/A	N/A
Elite Deep	N/A	N/A	STD	N/A	N/A	N/A
PB/PB Deep	N/A	N/A	STD	N/A	N/A	N/A

UPGRADABLE BACK OPTIONS

Matrix Mini Back 10" - 14" wide \$845
(Comes standard with Kid*ab*bra fabric)

☐ MINI1010 ☐ MINI1012 ☐ MINI1110 ☐ MINI1112
☐ MINI1210 ☐ MINI1212 ☐ MINI1214 ☐ MINI1310
☐ MINI1312 ☐ MINI1314 ☐ MINI1410 ☐ MINI1412
☐ MINI1414

Fabric Options for Mini Back

☐ MINI-S Startex (Non Reversible) N/C

☐ MINI-P Pink Lollipop (Non Reversible) N/C

☐ MINI-B Blue Blizzard (Non Reversible) N/C

☐ HIMINI Chest Harness Interface Kit for PB Back⁽¹⁾⁽²⁾ \$110

Matrix Elite Back 14" - 15" wide (3" Contour) \$766

☐ PBE1410 ☐ PBE1412 ☐ PBE1414 ☐ PBE1416 ☐ PBE1418
☐ PBE1420 ☐ PBE1510 ☐ PBE1512 ☐ PBE1514 ☐ PBE1516
☐ PBE1518 ☐ PBE1520

Matrix Elite TR Back 15" wide \$866

☐ PBE1516-TR ☐ PBE1518-TR ☐ PBE1520-TR

Matrix Elite Deep Back 14" - 15" wide \$918

☐ EDB1412 ☐ EDB1416 ☐ EDB1418 ☐ EDB1420 ☐ EDB1510
☐ EDB1512 ☐ EDB1514 ☐ EDB1516 ☐ EDB1518 ☐ EDB1520

Matrix Elite E2 Back 14" - 15" wide/10"-16"H(3" Contour) \$965

☐ E2S1410 ☐ E2S1412 ☐ E2S1414 ☐ E2S1416
☐ E2S1514 ☐ E2S1516

Matrix Elite E2 Back 14" - 15" wide/18"/20"H(3" Contour) \$995

☐ E2S1418 ☐ E2S1420 ☐ E2S1518 ☐ E2S1520

Matrix Elite E2 Deep Back 14" - 15" wide/10"-16"H(6" Contour)
..... \$1,125

☐ E2D1410 ☐ E2D1412 ☐ E2D1414 ☐ E2D1416
☐ E2D1514 ☐ E2D1516

Matrix Elite E2 Deep Back 14" - 15" wide/18"/20"H(6" Contour)
..... \$1,155

☐ E2D1418 ☐ E2D1420 ☐ E2D1518 ☐ E2D1520

Matrix PB Back 14"(5" Contour) \$734

☐ MPB1412 ☐ MPB1416

Matrix PB Deep Back 14" (7" Contour) \$892

☐ MPBD1412 ☐ MPBD1416

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harness



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

OTHER UPGRADABLE BACK OPTIONS

☐ EBSS	Super Soft Foam for Elite Back Upgrade.	\$103
☐ WSCU	Infection Control Cover for Elite/PB Backs.	\$103
	purchased with back	
☐ WSC	Infection Control Cover for Elite/PB Backs.	\$150
	purchased with without back	
☐ CHIEB	Chest Harness Interface Kit for Elite Back ⁽¹⁾⁽²⁾	\$103
☐ CHIEB	Chest Harness Interface Kit for PB Back ⁽¹⁾⁽²⁾	\$52

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harnesses.

HEADREST MOUNTING HARDWARE OPTIONS

☐ LXSHW	Loxx Standard Hardware 15"	\$275
☐ LXMMH	Loxx Mini Hardware 11"	\$275
☐ MEHWRE	Elan Hardware 15"	\$285
☐ MEMH	Elan Mini Hardware 11"	\$285
☐ MEHW-EXT	Elan Extended Hardware	\$320
☐ FPLXSHW	FLIPP Loxx Standard Hardware 15"	\$705 NEW
☐ FPLXSHW-EX	FLIPP Loxx Extended Hardware.	\$705 NEW
☐ FPLXMMH	FLIPP Loxx Mini Hardware 11"	\$705
☐ FPLXMMH-EX	FLIPP Loxx Mini Extended Hardware.	\$705
☐ FPMEHW	FLIPP Elan Hardware 15"	\$715
☐ FPMEMH	FLIPP Elan Mini Hardware 11"	\$715
☐ FPMEHW-EXT	FLIPP Elan Extended Hardware	\$780

MATRIX ELAN HEADREST PAD⁽¹⁾

Standard Pads

☐ ESP6	Matrx Elan Standard Pad 6"	\$310
☐ ESP10	Matrx Elan Standard Pad 10"	\$310
☐ ESP14	Matrx Elan Standard Pad 14"	\$310
☐ ESP6-HUG	Matrx Elan Hug Standard Pad 6"	\$310
☐ ESP10-HUG	Matrx Elan Hug Standard Pad 10"	\$310
☐ ESP14-HUG	Matrx Elan Hug Standard Pad 14"	\$310
☐ SNUG14	Matrx SNUG Standard Pad 14"	\$625
☐ SNUG14-A	Matrx SNUG Asymmetrical Pad 14"	\$645
☐ SNUG18	Matrx SNUG Standard Pad 18"	\$625
☐ SNUG18-A	Matrx SNUG Asymmetrical Pad 18"	\$645

Standard Pads with Infection Control Cover

☐ ESP6-IC	Matrx Elan Infection Control Pad 6"	\$335
☐ ESP10-IC	Matrx Elan Infection Control Pad 10"	\$335
☐ ESP14-IC	Matrx Elan Infection Control Pad 14"	\$335
☐ ESP6-IC-HUG	Matrx Elan Hug Infection Control Pad 6"	\$335
☐ ESP10-IC-HUG	Matrx Elan Hug Infection Control Pad 10"	\$335
☐ ESP14-IC-HUG	Matrx Elan Hug Infection Control Pad 14"	\$335
☐ SNUG14-IC	Matrx SNUG Standard Infection Control Pad 14"	\$665
☐ SNUG14-A-IC	Matrx SNUG Asymmetrical Infection Control Pad 14"	\$685
☐ SNUG18-IC	Matrx SNUG Standard Infection Control Pad 18"	\$665
☐ SNUG18-A-IC	Matrx SNUG Asymmetrical Infection Control Pad 18"	\$685

Occipital /4-Point Pads

☐ EOP9	Matrx Elan Occipital Pad 9"	\$341
☐ EOP12	Matrx Elan Occipital Pad 12"	\$341
☐ E4PTSM	Matrx Elan 4-Point Pad 9"W x 7"H	\$448
☐ E4POINT	Matrx Elan 4-Point Pad 11"W x 10"H	\$448

Occipital / 4-Point Pads with Infection Control Cover

☐ EOP9-IC	Matrx Elan Occipital Pad 9" with Infection Control.	\$366
☐ EOP12-IC	Matrx Elan Occipital Pad 12" with Infection Control.	\$366
☐ E4PTSM-IC	Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control.	\$473
☐ E4POINT-IC	Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control	\$473

1. If mounting on PB Elite must order removable hardware



MATRIX MINI FIXED LATERALS - FORMINI BACKS

Matrix MINI Fixed Lateral with Pad

☐ MINILATL	LEFT (4.5" x 3")	\$130
☐ MINILATR	RIGHT(4.5" x 3")	\$130

MATRIX MINI SWING AWAY LATERALS - FORMINI BACKS

MINI Swing Away Lateral Mounting Hardware

☐ SAHMINIL	LEFT SIDE	\$250
☐ SAHMINIR	RIGHT SIDE	\$250

Mini Lateral Trunk Support RIGHT(each)

Select Pad Size:

☐ EPADXSR	XS(3.25"x4.25")	\$115
☐ EPADSR	S(3.5"x5.25")	\$115

Mini Lateral Trunk Support LEFT(each)

Select Pad Size:

☐ EPADXSL	XS(3.25"x4.25")	\$115
☐ EPADSL	S(3.5"x5.25")	\$115

Fabric Options for Mini Lateral Trunk Supports

- ☐ MINILAT-K Startex (None Reversible) N/C
☐ MINILAT-P Pink Lollipop (None Reversible) N/C
☐ MINILAT-B Blue Blizzard (None Reversible) N/C

MATRIX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS

☐ PBSSL	Matrix Small Lateral Fixed-Left(3x5.5")	\$113
☐ PBSSR	Matrix Small Lateral Fixed-Right(3x5.5")	\$113
☐ PBLSL	Matrix Standard Lateral Fixed-Left(7.5x5")	\$116
☐ PBLSR	Matrix Standard Lateral Fixed-Right (7.5x5")	\$116
☐ PBOLL	Matrix Offset Lateral Fixed Left	\$125
☐ PBOLR	Matrix Offset Lateral Fixed Right	\$125

MATRIX SWING AWAY LATERALS - FORELITE AND ELITE TR BACKS

Matrix Lateral Trunk Support Left⁽¹⁾

☐ EPADXL	Size XS (3.25x4.25")	\$115
☐ EPADSL	Size S (3.5x5.25")	\$115
☐ EPADML	Size M (3.75x5.75")	\$115
☐ EPADL L	Size L (4.25x6.5")	\$115

Matrix Lateral Trunk Support Right⁽¹⁾

☐ EPADXR	Size XS (3.25x4.25")	\$115
☐ EPADSR	Size S (3.5x5.25")	\$115
☐ EPADMR	Size M (3.75x5.75")	\$115
☐ EPADLR	Size L (4.25x6.5")	\$115

1. Must choose swing away hardware on next page

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX SWING AWAY HARDWARE - FORELITE AND ELITE TR BACKS

☐ SAEHBL	Matrix Elite Swing Away Hardware	
	Left	\$250
☐ SAHEBR	Matrix Elite Swing Away Hardware	
	Right	\$250
☐ OSAHL	Matrix Offset Elite Swing Away	
	Hardware Left	\$250
☐ OSAHR	Matrix Offset Elite Swing Away	
	Hardware Right	\$250

OTHER MATRIX ACCESSORIES

☐ MACESL	Multi Axis Arm Support - Left ⁽¹⁾	\$225
☐ MACESR	Multi Axis Arm Support - Right ⁽¹⁾	\$225

1. Allows rotation of the armpad, only compatible with Gel pad, Flat pad or Ergonomic Arm Troughs

MAXX LATERALS TRUNK SUPPORTS⁽¹⁾

Left Fixed Mount Lateral Trunk Support⁽¹⁾

Select Left Trunk Support Size

☐ LATXL	Size XS (3.25" x 4.25")	\$100
☐ LATSL	Size S (3.5" x 5.25")	\$100
☐ LATML	Size M (3.75" x 5.75")	\$100
☐ LATLL	Size L (4.25" x 6.5")	\$100

Right Fixed Mount Lateral Trunk Support⁽¹⁾

Select Right Trunk Support Size

☐ LATXR	Size XS (3.25" x 4.25")	\$100
☐ LATSR	Size S (3.5" x 5.25")	\$100
☐ LATMR	Size M (3.75" x 5.75")	\$100
☐ LATLR	Size L (4.25" x 6.5")	\$100

1. Not compatible with "deep" matrix backs or sling back.

Swing Away Trunk Support Hardware⁽⁴⁾

☐ SALATL	Left MAXX Style Swing Away HDW ⁽¹⁾	\$195
---------------------------------	---	-------

Select Left Trunk Support Size

☐ LATXL	Size XS (3.25" x 4.25")	\$100
☐ LATSL	Size S (3.5" x 5.25")	\$100
☐ LATML	Size M (3.75" x 5.75")	\$100
☐ LATLL	Size L (4.25" x 6.5")	\$100

☐ SALATR	Right MAXX Style Swing Away HDW ⁽¹⁾	\$195
---------------------------------	--	-------

Select Right Trunk Support Size

☐ LATXR	Size XS (3.25" x 4.25")	\$100
☐ LATSR	Size S (3.5" x 5.25")	\$100
☐ LATMR	Size M (3.75" x 5.75")	\$100
☐ LATLR	Size L (4.25" x 6.5")	\$100

☐ ATLL	Additional Telescoping Link - Left	\$75
☐ ATLR	Additional Telescoping Link - Right	\$75

1. Up to 2.75" medial offset
2. Recommended for all backs except rehab back, PB and elite deep backs
3. Medial offset depends on back thickness
4. Must choose fixed mount and pads

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HIP SUPPORTS

Lateral Hip Support with Fixed Mounting Hardware Left⁽¹⁾

Select Pad Size and mounting position at Hip:

- ☐ HSPXLH Size XS (5.5" W x 4" H) mounted at Hip\$120
☐ HSPSMLH Size S/M (7" W x 4" H) mounted at Hip\$120
☐ HSPSLH Size S (7" W x 6" H) mounted at Hip\$120
☐ HSPMLH Size M (9" W x 4" H) mounted at Hip\$120
☐ HSPLLH Size L (13" W x 4" H) mounted at Hip\$120

Select Pad Size and mounting position at Knee:

- ☐ HSPXLK Size XS (5.5" W x 4" H) mounted at Knee \$120
☐ HSPSMLK Size S/M (7" W x 4" H) mounted at Knee \$120
☐ HSPSLK Size S (7" W x 6" H) mounted at Knee \$120
☐ HSPMLK Size M (9" W x 4" H) mounted at Knee \$120
☐ HSPLLK Size L (13" W x 4" H) mounted at Knee \$120

Lateral Hip Support with Fixed Mounting Hardware Right⁽¹⁾

Select Pad Size and mounting position:

- ☐ HSPXRH Size XS (5.5" W x 4" H) mounted at Hip \$120
☐ HSPSMRH Size S/M (7" W x 4" H) mounted at Hip \$120
☐ HSPSRH Size S (7" W x 6" H) mounted at Hip \$120
☐ HSPMRH Size M (9" W x 4" H) mounted at Hip \$120
☐ HSPLRH Size L (13" W x 4" H) mounted at Hip \$120

Select Pad Size and mounting position at Knee:

- ☐ HSPXRK Size XS (5.5" W x 4" H) mounted at Knee \$120
☐ HSPSMRK Size S/M (7" W x 4" H) mounted at Knee \$120
☐ HSPSRK Size S (7" W x 6" H) mounted at Knee \$120
☐ HSPMRK Size M (9" W x 4" H) mounted at Knee \$120
☐ HSPLRK Size L (13" W x 4" H) mounted at Knee \$120

1. Must choose swing away hardware below if required.

REMOVABLE HIP SUPPORT HARDWARE

Lift Off Removable Hardware for Hip Support

- ☐ LOHL Left \$195
☐ LOHR Right \$195

MAXX Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports

Select mounting position:

- ☐ MHSHLH Left mounted at Hip \$215
☐ MHSHLK Left mounted at Knee \$215
☐ MHSHRH Right mounted at Hip \$215
☐ MHSHRK Right mounted at Knee \$215

Swing Away Removable Hardware for Hip Supports

Select mounting position:

- ☐ SHSHLH Left mounted at Hip \$215
☐ SHSHLK Left mounted at Knee \$215
☐ SHSHRH Right mounted at Hip \$215
☐ SHSHRK Right mounted at Knee \$215

Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports⁽¹⁾

Select mounting position:

- ☐ EMQHLH Left mounted at Hip \$215
☐ EMQHLLK Left mounted at Knee \$215
☐ EMQHRH Right mounted at Hip \$215
☐ EMQHRK Right mounted at Knee \$215

1. Recommended when finished seat cushion thickness will be greater than 3.5" high

NOTES

Total Retail Value \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM