



Date of Order _____
 Dealer Name: _____
 Dealer Account #: _____
 PO #: _____
 Tag: _____
 Purchasing Contact: _____

Phone: _____
 Fax: _____
 Ship to Address: _____
 City: _____
 State: _____ Zip Code: _____
 Comments: _____

Order Form Content

PAGE

Matrix® MAC Back Dual Layer Foam Back Model	2
Matrix® MAC Back Dual Layer Foam Back Model with Infection Control Cover	3
Matrix® MAC Back HUG StaminaFibre® Back Model	4
Matrix® MAC Back HUG StaminaFibre® Back Model with Infection Control Cover	5
Mounting Hardware for Matrix® MAC Backs	6
Optional Accessories for Matrix® MAC Backs	7-8

Center Panel

- Dual Layer Foam or HUG StaminaFibre® Model

Headrest Interface

- Heavy-duty interface included

- Fixed and Swing Away Options available

EasySet Hardware

- Height, depth, angle and rotation adjustable





MAC Back

Dual Layer Foam Back Model

MAC BACK: Complete with fixed laterals (specify below) and dual-layer HR foam

Weight capacity: 300 lbs.

Standard hardware:

- EasySet fits 3/4", 7/8" and 1" diameter back canes (Other hardware options available on page 4)
- EasySet is angle (40°), depth (2" in front or 2" behind back canes) & rotation adjustable

Heavy-duty headrest interface included

Standard cover: Meshtex fabric



BACK WIDTH									MSRP	HCPCS
SIZE	14"	15"	16"	17"	18"	19"	20"			
FITS W/C	13" - 16"	14" - 17"	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"			
14" (1) (2)	☐ MAC1414	☐ MAC1514	-	-	-	-	-		\$1117	E2620
16"	☐ MAC1416	☐ MAC1516	☐ MAC1616	☐ MAC1716	☐ MAC1816	☐ MAC1916	☐ MAC2016			
18"	☐ MAC1418	☐ MAC1518	☐ MAC1618	☐ MAC1718	☐ MAC1818	☐ MAC1918	☐ MAC2018			
20"	☐ MAC1420	☐ MAC1520	☐ MAC1620	☐ MAC1720	☐ MAC1820	☐ MAC1920	☐ MAC2020			
22"	-	-	☐ MAC1622	☐ MAC1722	☐ MAC1822	☐ MAC1922	☐ MAC2022			

1. Backs 14" Length can only be ordered with 8" Length Laterals

2. 14" and 15" Back lengths include 1/2" RBR Foam

LATERAL CONTOUR EXTENSIONS			
LENGTH	DEPTH	LEFT LATERAL	RIGHT LATERAL
10"	3.5"	☐ MLC-S	☐ MLC-S
10"	6"	☐ MLC-D	☐ MLC-D
10"	9"	☐ MLC-XD	☐ MLC-XD
8"	3.5"	☐ MLS-S	☐ MLS-S
8"	6"	☐ MLS-D	☐ MLS-D
8"	9"	☐ MLS-XD	☐ MLS-XD



Extra Cover for MAC Back Dual Layer Foam			MSRP
☐	Meshtex Fabric Cover with Dual Layer Foam	MACC(wwwll)* WIDTH _____ LENGTH _____	\$170

*Refer to back sizes above for proper fit



m MAC Back

Dual Layer Foam with Infection Control Back Model

MAC BACK: Complete with fixed laterals (specify below) and dual-layer HR foam
Weight capacity: 300 lbs.

Standard hardware:

- EasySet fits 3/4", 7/8" and 1" diameter back canes (Other hardware options available on page 4)
- EasySet is angle (40°), depth (2" in front or 2" behind back canes) & rotation adjustable

Heavy-duty headrest interface included

Standard cover: Front panel in contact with the user is made of wipeable Startex (PU coated). Rear panel is Meshtex.



BACK WIDTH									
SIZE	14"	15"	16"	17"	18"	19"	20"	MSRP	HCPCS
FITS W/C	13" - 16"	14" - 17"	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"	\$1205	E2620
14" (1)	☐ MAC1414-IC	☐ MAC1514-IC	-	-	-	-	-		
16"	☐ MAC1416-IC	☐ MAC1516-IC	☐ MAC1616-IC	☐ MAC1716-IC	☐ MAC1816-IC	☐ MAC1916-IC	☐ MAC2016-IC		
18"	☐ MAC1418-IC	☐ MAC1518-IC	☐ MAC1618-IC	☐ MAC1718-IC	☐ MAC1818-IC	☐ MAC1918-IC	☐ MAC2018-IC		
20"	☐ MAC1420-IC	☐ MAC1520-IC	☐ MAC1620-IC	☐ MAC1720-IC	☐ MAC1820-IC	☐ MAC1920-IC	☐ MAC2020-IC		
22"	-	-	☐ MAC1622-IC	☐ MAC1722-IC	☐ MAC1822-IC	☐ MAC1922-IC	☐ MAC2022-IC		

1. Backs 14" Length can only be ordered with 8" Length Laterals

LATERAL CONTOUR EXTENSIONS			
LENGTH	DEPTH	LEFT LATERAL	RIGHT LATERAL
10"	3.5"	☐ MLC-S-IC	☐ MLC-S-IC
10"	6"	☐ MLC-D-IC	☐ MLC-D-IC
10"	9"	☐ MLC-XD-IC	☐ MLC-XD-IC
8"	3.5"	☐ MLS-S-IC	☐ MLS-S-IC
8"	6"	☐ MLS-D-IC	☐ MLS-D-IC
8"	9"	☐ MLS-XD-IC	☐ MLS-XD-IC



Extra Infection Control Cover for MAC Back Dual Layer Foam			MSRP
☐	Startex Wipeable Seating Surface with Dual Layer Foam	MACC(wwwll)-IC* WIDTH _____ LENGTH _____	\$225

*Refer to back sizes above for proper fit



MAC Back

HUG StaminaFibre® Back Model

MAC BACK: Complete with fixed laterals (specify below) and conformable StaminaFibre® insert
Weight capacity: 300 lbs.

Standard hardware:

- EasySet fits 3/4", 7/8" and 1" diameter back canes (Other hardware options available on page 4)
- EasySet is angle (40°), depth (2" in front or 2" behind back canes) & rotation adjustable

Heavy-duty headrest interface included


Standard cover: Meshtex fabric



BACK WIDTH									MSRP	HCPCS
SIZE	14"	15"	16"	17"	18"	19"	20"			
FITS W/C	13" - 16"	14" - 17"	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"			
14" ⁽¹⁾	☐ MAC1414-HUG	☐ MAC1415-HUG	-	-	-	-	-		\$1234	E2620
16"	☐ MAC1416-HUG	☐ MAC1516-HUG	☐ MAC1616-HUG	☐ MAC1716-HUG	☐ MAC1816-HUG	☐ MAC1916-HUG	☐ MAC2016-HUG			
18"	☐ MAC1418-HUG	☐ MAC1518-HUG	☐ MAC1618-HUG	☐ MAC1718-HUG	☐ MAC1818-HUG	☐ MAC1918-HUG	☐ MAC2018-HUG			
20"	☐ MAC1420-HUG	☐ MAC1520-HUG	☐ MAC1620-HUG	☐ MAC1720-HUG	☐ MAC1820-HUG	☐ MAC1920-HUG	☐ MAC2020-HUG			
22"	-	-	☐ MAC1622-HUG	☐ MAC1722-HUG	☐ MAC1822-HUG	☐ MAC1922-HUG	☐ MAC2022-HUG			

1. Backs 14" Length can only be ordered with 8" Length Laterals

LATERAL CONTOUR EXTENSIONS			
LENGTH	DEPTH	LEFT LATERAL	RIGHT LATERAL
10"	3.5"	☐ MLC-S	☐ MLC-S
10"	6"	☐ MLC-D	☐ MLC-D
10"	9"	☐ MLC-XD	☐ MLC-XD
8"	3.5"	☐ MLS-S	☐ MLS-S
8"	6"	☐ MLS-D	☐ MLS-D
8"	9"	☐ MLS-XD	☐ MLS-XD



Extra Cover for MAC Back with HUG StaminaFibre			MSRP
☐	Meshtex Fabric Cover with HUG pad	MACC-H(wwll)* WIDTH _____ LENGTH _____	\$199

*Refer to back sizes above for proper fit



m MAC Back

HUG StaminaFibre[®] with Infection Control Back Model

MAC BACK: Complete with fixed laterals (specify below) and conformable StaminaFibre[®] insert
Weight capacity: 300 lbs.

Standard hardware:

- EasySet fits 3/4", 7/8" and 1" diameter back canes (Other hardware options available on page 4)
- EasySet is angle (40°), depth (2" in front or 2" behind back canes) & rotation adjustable

Heavy-duty headrest interface included


Standard cover: Front panel in contact with the user is made of wipeable Startex (PU coated). Rear panel is Meshtex



BACK WIDTH									MSRP	HCPCS
SIZE	14"	15"	16"	17"	18"	19"	20"			
FITS W/C	13" - 16"	14" - 17"	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"			
BACK LENGTH	14"	☐ MAC1414-HUG-IC	☐ MAC1415-HUG-IC	-	-	-	-	-	\$1322	E2620
	16"	☐ MAC1416-HUG-IC	☐ MAC1516-HUG-IC	☐ MAC1616-HUG-IC	☐ MAC1716-HUG-IC	☐ MAC1816-HUG-IC	☐ MAC1916-HUG-IC	☐ MAC2016-HUG-IC		
	18"	☐ MAC1418-HUG-IC	☐ MAC1518-HUG-IC	☐ MAC1618-HUG-IC	☐ MAC1718-HUG-IC	☐ MAC1818-HUG-IC	☐ MAC1918-HUG-IC	☐ MAC2018-HUG-IC		
	20"	☐ MAC1420-HUG-IC	☐ MAC1520-HUG-IC	☐ MAC1620-HUG-IC	☐ MAC1720-HUG-IC	☐ MAC1820-HUG-IC	☐ MAC1920-HUG-IC	☐ MAC2020-HUG-IC		
	22"	-	-	☐ MAC1622-HUG-IC	☐ MAC1722-HUG-IC	☐ MAC1822-HUG-IC	☐ MAC1922-HUG-IC	☐ MAC2022-HUG-IC		

1. Backs 14" Length can only be ordered with 8" Length Laterals

LATERAL CONTOUR EXTENSIONS			
LENGTH	DEPTH	LEFT LATERAL	RIGHT LATERAL
10"	3.5"	☐ MLC-S-IC	☐ MLC-S-IC
10"	6"	☐ MLC-D-IC	☐ MLC-D-IC
10"	9"	☐ MLC-XD-IC	☐ MLC-XD-IC
8"	3.5"	☐ MLS-S-IC	☐ MLS-S-IC
8"	6"	☐ MLS-D-IC	☐ MLS-D-IC
8"	9"	☐ MLS-XD-IC	☐ MLS-XD-IC



Extra Infection Control Cover for MAC Back with HUG StaminaFibre			MSRP
☐	Startex Wipeable Seating Surface with Hug pad	MACC-H(wwll)-IC* WIDTH _____ LENGTH _____	\$265


*Refer to back sizes above for proper fit



Optional Back Cane Inserts

MODEL

MSRP



For 16mm back Canes (Pair)



BCI16

\$5.00

For 20mm back canes (Pair)



BCI20

\$5.00

For Quickie fold down back canes (Pair)



BCIQFD

\$5.00

Standard Hardware

AT TIME OF ORDER

SOLD SEPARATELY



EasySet Hardware Fits 3/4", 7/8", & 1" back canes
Angle (40°), Depth (2" in front or 2" behind back canes) & rotation adjustable
Hardware will include Standard EasySet, unless upgrades are selected below:

(Included)



PBMH

\$263

Optional Hardware

AT TIME OF ORDER

SOLD SEPARATELY



Easy Set Hardware Fits 1 1/8" back canes



PBMH-QU

n/c



PBMH-Q

\$263

Easy Set Hardware Fits 1 1/4" back canes



PBMH-XU

n/c



PBMH-X

\$263



Heavy-Duty EasySet

Recommended for users with high tone or heavy intended use



HDMHU

n/c



HDMH

\$380



Extended Range Mounting

Provides up to 4" depth adjustment in front or behind canes
Recommended for HUG & HUG with Shape Insert



ERMHU

n/c



ERMH

\$322



Fixed EasySet

Recommended for Tilt/Recline applications, requires tools for removal
Fits chairs 1" narrower to 2" wider than back model width



ESFXU

\$59



ESFX

\$263



Split Mounting Hardware

For avoidance of obstacles on frame



SMHU

n/c



SMH

\$263

MAC Back Interface Kit

MODEL

MSRP



Permobil Mounting Hardware, E2, MAC Back



PERM-E2MAC

\$545

Quantum Mounting Hardware, E2, MAC Back, Fixed



QUANMH-E2 FIXED

\$545

Quantum Mounting Hardware, E2 MAC, Removable



QUAN-E2MAC

\$545

Quickie Mounting Hardware, E2, Mac Back



QUIKMH-E2

\$545



Optional Accessories


MSRP

	Matrix SHAPE Insert <i>Customizable shape interface for accommodation of orthopedic asymmetries.</i> <i>Recommended for use with HUG models.</i> <i>Extended Range Mounting Hardware recommended.</i>	☐ SHAPE		\$111
	MAC HUG Insert <i>Does not include cover</i> SPECIFY: WIDTH____ LENGTH____ <i>*Refer to Back Sizes on Page 2</i>	☐ HUGI		\$111
	PostureFlex Abdominal Support with Tether straps (Mounting Kit) Length 6.5", Width 16" to 23"	☐ PFLEX-MAC-M		\$305
	PostureFlex Abdominal Support with Tether straps (Mounting Kit) Length 7", Width 20" 32"	☐ PFLEX-MAC-L		\$305
	Quad release tab <i>For easy independent application and removal of PostureFlex</i>	☐ QRTAB		\$29
	MAC Chest Harness Interface <i>Allows precise placement of Shoulder Strap Guides;</i> <i>Fits MAC backs 14" - 15" wide</i>	☐ E2CHI-S		\$88
	MAC Chest Harness Interface <i>Allows precise placement of Shoulder Strap Guides;</i> <i>Fits MAC backs 16" - 20" wide</i>	☐ E2CHI-M		\$88
	MAC BackBone PSIS Support <i>Can be trimmed to desired size</i>	☐ E2BKB		\$64
	Small Privacy Flap 13" W x 16" L - fits 15"-16" wide backs	☐ PFLAP-S		\$47
	Medium Privacy Flap 15" W- 16" L - fits 17"-18" wide backs	☐ PFLAP-M		\$47
	Large Privacy Flap 17"W - 16" L - fits 19"-20" wide backs	☐ PFLAP-L		\$47
	Additional Lateral 3.5" depth x 10" length (single)	☐ MLC-3.5	☐ MLC-3.5-IC	\$170
	Additional Lateral 6" depth x 10" length (single)	☐ MLC-6	☐ MLC-6-IC	\$170
	Additional Lateral 9" depth x 10" length (single)	☐ MLC-9	☐ MLC-9-IC	\$170
	Additional Lateral 3.5" depth x 8" length (single)	☐ MLS-3.5	☐ MLS-3.5-IC	\$170
	Additional Lateral 6" depth x 8" length (single)	☐ MLS-6	☐ MLS-6-IC	\$170
	Additional Lateral 9" depth x 8" length (single)	☐ MLS-9	☐ MLS-9-IC	\$170
*Default fabric is Meshtex				
	Protect - Infection Control Slip Cover <i>Stands up to repeated use of hospital-grade disinfectants</i> <i>Machine-washable at temperatures required to kill bacteria and viruses</i>	☐ PROTMACWWLL		\$95

*Default fabric is Meshtex



Optional Swing Away Lateral Pads
(Price per each)



10"L Swing-Away Standard Lateral with Pad ⁽¹⁾⁽²⁾								HCPCS E0956/E0134	
Users Left					Users Right				
Pad Length	Pad Depth	MODEL #			Pad Length	Pad Depth	MODEL #		
3.25"	4.25"		SAMLL-XS	\$462	3.25"	4.25"		SAMLR-XS	\$462
3.5"	5.25"		SAMLL-S	\$462	3.5"	5.25"		SAMLR-S	\$462
3.75"	5.75"		SAMLL-M	\$462	3.75"	5.75"		SAMLR-M	\$462
4.25"	6.5"		SAMLL-L	\$462	4.25"	6.5"		SAMLR-L	\$462

8"L Swing-Away Short Lateral with Pad ⁽¹⁾							HCPCS E0956/E0134		
Users Left					Users Right				
Pad Length	Pad Depth	MODEL #			Pad Length	Pad Depth	MODEL #		
3.25"	4.25"		SAMSL-XS	\$462	3.25"	4.25"		SAMSR-XS	\$462
3.5"	5.25"		SAMSL-S	\$462	3.5"	5.25"		SAMSR-S	\$462
3.75"	5.75"		SAMSL-M	\$462	3.75"	5.75"		SAMSR-M	\$462
4.25"	6.5"		SAMSL-L	\$462	4.25"	6.5"		SAMSR-L	\$462

1. Wipeable Startex not available with swing-away laterals.
2. 14" back lengths can only be ordered with 8" length laterals.