

CONTACT MOTION CONCEPTS CUSTOMER SERVICE: 1.888.433.6818 / QUOTES@MOTIONCONCEPTS.COM

REQUIRED INFORMATION

☐ Quote ☐ Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag: _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

ATP/Therapist: _____

Ship to Address: _____

City: _____ State: _____

Zip: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

Client Height:		
Client Weight:		
Use the table below to calculate your Seat Pan to Footplate Measurement. Then select the appropriate range in the legrest section.	Solid Seat	Captain Seat
Knee to Heel:		
Cushion Thickness: (Subtract)		4"
Cushion Compression:(Add)		
Seat Pan to Footplate:(Total)		

Notes:

BASE	X3	A3
HCPCS CODES	K0856, K0861	K0848, K0849
BASE WIDTH	23.25"	25.5"
BASE LENGTH	36.5"	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"	14"
CASTER WHEEL SIZE	6" (front/rear)	6" (front/rear)
BASE ONLY TURNING RADIUS	19.5"	20.5"
DRIVETRAIN	4-pole Motors	4-pole Motors
BATTERIES	2 x M34	2 x M34

BASE	X3	A3
BATTERY CHARGER	8A, 110V	8A, 110V
BATTERIES WEIGHT	38.5 lbs (each)	38.5 lbs (each)
RANGE ⁽²⁾	17.7 miles	17.7 miles
SUSPENSION SYSTEM	Active Ride Control (ARC)	
WEIGHT CAPACITY	300 lbs	300lbs
BASE WEIGHT	157 lbs	159 lbs
GROUND CLEARANCE	3"	3"
MAXIMUM SPEED	6.2 mph	6.2 mph

⁽¹⁾ All A3/X3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

SOLID SEAT



CAPTAIN SEAT



ROVI POWER BASE☐ 888WNLEX3-SEMIRROVI X3 Base for Captain Seat
HCPCS code K0849\$7,746☐ 888WNLEX3-SEMIRBROVI X3 Base for Solid Seat
HCPCS code K0848\$11,161☐ 888WNLEA3-SEMIRROVI A3 Base for Captain Seat
HCPCS code K0849\$7,746☐ 888WNLEA3-SEMIRBROVI A3 Base for Solid Seat
HCPCS code K0848\$11,161

Select seat to floor height: (inches)		STF Height	Captain's Seat		Solid Seat	
			No Elevate	With Elevate	No Elevate	With Elevate
☐ TLOW	☐ TLOW-S	Low	20 ³ / ₄	22 ¹ / ₂	17	18 ³ / ₄
☐ TMED	☐ TMED-S	Medium	21 ³ / ₄	23 ¹ / ₂	18	19 ³ / ₄
☐ TTALL	☐ TTALL-S	Tall	22 ³ / ₄	24 ¹ / ₂	19	20 ³ / ₄
☐ TXTALL	☐ TXTALL-S	X-Tall	23 ³ / ₄	N/A	20	N/A

SEAT STYLE☐ SEMIR

Semi-Recline (35°) Captain's Seat STD

☐ SEMIRB

Semi-Recline (35°) with Solid Seat \$145

SEAT COVER☐ BSC

Black Slip Cover \$573

1. Includes seat, back, and headrest cover

PERFORMANCE MOTORS☐ HSM62Optional 6.2 mph High Speed Motors \$998
Must order OnTraxx Enhanced Driving Module☐ OEDMOnTraxx Enhanced Driving Module
HCPCS code K0108\$788**TRANSPORT TIE DOWN**☒ XBKTSWheelchair Transport Brackets⁽¹⁾ \$263

1. For unoccupied use only

TIRE OPTIONS☐ STIRES

Standard Solid Tires STD

☐ PTIRES

Pneumatic Tires \$209

INSTALLATION OF BATTERIES

HCPCS code E2359

☐ M34 SLD G

MK Group 34 Battery (ea) Quantity () \$442

LIGHTED FENDER

Lights operate through single push button

☐ BWFL

Blizzard White (Metallic)

☐ OCFL

Orange Crush

☐ MGFL

Motion Green

☐ BCFL

Black Carbon (Matte)

☐ BPFL

Black Panther (Gloss)

☐ CARFL

Candy Apple Red (Metallic)

☐ GGFL

Graphite Grey (Metallic)

☐ IBFL

Indigo Blue

☐ MYFL

Mellow Yellow

☐ NBFL

Nordic Blue

☐ RDPFL

Razzle Dazzle Pink (Metallic)

☐ PRFL

Purple Reign

MATCHING BATTERY DOORS☐ MBDMatching Battery Doors \$495
(Select fender color first)**CUSTOM DESIGNS**Please see form number TRD0607 - ROVI Custom Designs Price Sheet,
for Custom/Dipped Fender colors, battery door and other personalized options.**ELEVATING SEAT**☐ ESM12CElevate lift 12" HCPCS code E2298 \$3,504
Weight capacity up to 300lbs☐ RDEReduced Drive for Elevating Seat⁽¹⁾ \$580

1. If reduced drive is not selected the system will be in drive lock-out when elevated.

ELECTRONICS FOR ELEVATING SEAT☐ SFCB

Single Function Electronics STD

Select Switch:

☐ SPBSO

Single Push Button

☐ STO

Single Toggle Switch

☐ SFCBASingle/Dual Function Drive Control⁽¹⁾
HCPCS code E2310 \$2,190

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ADDITIONAL SWITCH FOR POWERED SEATING

- ☐ ASL320R ASL Pulse Mechanical Switch - Red\$181
- ☐ ASL320G ASL Pulse Mechanical Switch - Green\$181
- ☐ ASL320BK ASL Pulse Mechanical Switch - Black\$181
- ☐ SESG Stealth Egg Switch - Green\$196
- ☐ SESB Stealth Egg Switch - Black\$196

MOUNTING FOR POWERED SEATING SWITCH

- ☐ ESML Left Side (Armpad)
- ☐ ESMR Right Side (Armpad)
- ☐ ESMJL Joystick Left (Mounts to Joystick Bracket)
- ☐ ESMJR Joystick Right (Mounts to Joystick Bracket)
- ☐ ESMD Do not mount

SEAT WIDTH SETTING

- ☐ W18 18" Wide (Actual Width 17.7" Solid / 19" Captain)..... N/C
- ☐ W20 20" Wide (Actual Width 18.5" Solid / 20" Captain)..... N/C
- ☐ W22 22" Wide (Actual Width 20.5" Solid / 22" Captain)..... N/C

SEAT DEPTH SETTING

- ☐ D16 16" Deep N/C
- ☐ D17 17" Deep N/C
- ☐ D18 18" Deep N/C
- ☐ D19 19" Deep N/C
- ☐ D20 20" Deep N/C
- ☐ D21 21" Deep N/C

Seat Width and Depth Compatibility						
	16D	17D	18D	19D	20D	21D
18W	✓	✓	✓	✓	N/A	N/A
20W	N/A	N/A	N/A	✓	✓	✓
22W	N/A	N/A	N/A	✓	✓	✓

JOYSTICK SELECTION

- ☐ D50677 VR2 Non Expandable Joystick N/C
- ☐ D50680 VR2 Non Expandable Joystick
w/ 2 Actuator Control N/C
- ☐ D51122 R-Net LED Joystick⁽¹⁾..... N/C
- ☐ D51165 R-Net Color Remote Joystick
with Monojacks⁽¹⁾ \$630
- ☐ D51711 R-Net Advanced Joystick⁽¹⁾⁽²⁾ \$1,045

1. Requires Expandable Controller and Harness (D51109 and REH)
2. 3.5" high-resolution VGA colour display, IR Control as standard, includes integrated Bluetooth Module.

EXPANDABLE CONTROLLER⁽¹⁾

- ☐ D51109 R-Net Expandable Controller
HCPCS code E2377\$775
- ☐ REH R-Net Harness for Expandable Controller
HCPCS code E2313\$595

1. Required with R-Net Joysticks

MOUNTING FOR DRIVER CONTROL

- ☐ FR Fixed Joystick Mounting - Right STD
- ☐ FL Fixed Joystick Mounting - Left STD
- ☐ MRR Maxx Resolve Swing-Away Joystick Mount
(Right Hand) HCPCS code E1032 \$573
- ☐ MRL Maxx Resolve Swing-Away Joystick Mount
(Left Hand) HCPCS code E1032 \$573
- ☐ SAQR Motion Height Adjustable Swing-Away
Quad Link (Right Hand) HCPCS code E1032 \$341
- ☐ SAQL Motion Height Adjustable Swing-Away
Quad Link (Left Hand) HCPCS code E1032 \$341

ELECTRONIC OPTIONS

- ☐ D50883 R-Net Output Module (EMC)⁽¹⁾ \$1,045
- ☐ D51580 R-Net Bluetooth Module IOS Apple⁽¹⁾ \$940
- ☐ D51111 R-Net Bluetooth Module Android/PC⁽¹⁾ \$940
- ☐ SPH Smart Phone Holder \$53
- ☐ MUSB Motion Dual USB Charger N/C

1. Requires D51165 or D51711 Joystick.

OPTIONAL JOYSTICK KNOB

HCPCS code E2323

- ☐ SA80941 "T" Handle for Joystick Knob \$121
- ☐ SA72441 Ball Handle for Joystick Knob \$121
- ☐ PC101B-20 BodyPoint Joystick Handle 3" U-Shaped \$126
- ☐ PC102B-20 BodyPoint Joystick Handle 4" U-Shaped \$126
- ☐ PC104B-20 BodyPoint Straight Handle \$126
- ☐ PC107B-20 BodyPoint Joystick Handle Domed Rubber \$126

LAP BELTS

- ☐ BELT48 Push Button Style Lap Belt, Length 48" STD
- ☐ BELT60 Push Button Style Lap Belt, Length 60" N/C
- ☐ BELT71 Push Button Style Lap Belt, Length 71" N/C
- ☐ PLB2 Bodypoint Padded Lap Belt 2 Point⁽¹⁾ \$222
- ☐ PLB4 Bodypoint Padded Lap Belt 4 Point⁽¹⁾ \$275

1. 67" length only option

ARMRESTS

HCPCS code E0973

- ☒ LVAHF Left Full Length, Adjustable Height\$167
- ☒ RVAHF Right Full Length, Adjustable Height \$167

MATRIX SEAT CUSHIONS - FOR SOLID SEAT SYSTEMS

Matrix PS Seat Cushion 14" - 20" wide/reversible cover
HCPCS code E2605 \$383

☐ PS1816 ☐ PS1818 ☐ PS2018 ☐ PS2020

Matrix PS Seat Cushion 21"-22" wide
HCPCS code E2606 \$605

☐ PS2218 ☐ PS2220

Matrix PSP Seat Cushion 14" - 20" wide/reversible cover
HCPCS code E2607 \$438

☐ PSP1816 ☐ PSP1818 ☐ PSP2018 ☐ PSP2020

Matrix PSP Seat Cushion 21"-22" wide
HCPCS code E2608 \$660

☐ PSP2218

Matrix Vi Seat Cushion 14" - 20" wide
HCPCS code E2607 \$497

☐ MA1816-VI ☐ MA1818-VI ☐ MA2018-VI ☐ MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide
HCPCS code E2608 \$556

☐ MA2218-VI ☐ MA2220-VI

Matrix Libra Seat Cushion 14" - 20" wide
HCPCS code E2624 \$567

☐ LC1816 ☐ LC1818 ☐ LC2018 ☐ LC2020

Matrix Libra HD Seat Cushion 22" wide
HCPCS code E2625 \$696

☐ LC2218 ☐ LC2220

Matrix Libra Adjust Seat Cushion 14" - 20" wide
HCPCS code E2624 \$731

☐ LC1816-ADJ ☐ LC1818-ADJ ☐ LC2018-ADJ ☐ LC2020-ADJ

Matrix Libra Adjust HD Seat Cushion 22" wide
HCPCS code E2625 \$860

☐ LC2218-ADJ ☐ LC2220-ADJ

Matrix Multi Cushion 14" - 20" wide/reversible cover
HCPCS code E2607 \$521

☐ MM1816 ☐ MM1818 ☐ MM2018 ☐ MM2020

Matrix Multi Cushion 22" wide/reversible cover
HCPCS code E2608 \$696

☐ MM2218 ☐ MM2220 ☐ MM2222

☐ OCCC "ON-CHAIR" Cushion Credit⁽¹⁾ (\$105)

1. Not available on Matrix PS or PSP Cushions.

STANDARD FRONT RIGGING

☐ FCMP Fixed Center Mount Platform⁽¹⁾ N/C

☐ FIFP Rubber Coated Footplate (11.5" W x 10") N/C

☐ FBIFP Large Rubber Coated Footplate (17" W x 13" d) \$250

Select seat pan to footplate:

☐ FSPTF3 10.25" - 14.25" N/C


☐ FSPTF4 14.25" - 20.25" N/C

Set at:

☐ FSA83 Set at 83° N/C

☐ FSA90 Set at 90°⁽²⁾ N/C

☐ FSA97 Set at 97°⁽²⁾ N/C



1. Not available on 18" W seats.

MANUAL ELEVATING CENTER MOUNT LEG REST

☐ MECML Manual Elevating Center Mount Legrest
With Plantar Flex. \$863

Select legrest set up position:

☐ 90SU Set at 90° (elevating range 90-40°) N/C

☐ 110TU Tucked at 110° (elevating range 110-60°)⁽¹⁾ N/C

Rubber Coated Footplate:

☐ MRCP Rubber Coated Footplate (11.5" W x 10") \$239

Select seat pan to footplate for Rubber Coated Footplate

☐ CSPTF1 12.75" - 16.25" N/C

☐ CSPTF2 15.25" - 20.25" N/C

Individual Footplates:

☐ IFSL Small - (5" x 7.5" deep) - Left \$239

☐ IFSR Small - (5" x 7.5" deep) - Right \$239

☐ IFML Medium - (5.5" x 9.5" deep) - Left \$239

☐ IFMR Medium - (5.5" x 9.5" deep) - Right \$239

☐ IFLL Large - (6" x 11.5" deep) - Left \$239

☐ IFLR Large - (6" x 11.5" deep) - Right \$239


Select seat pan to footplate for Individual Footplates

☐ CSPTF3 10" - 15" N/C

☐ CSPTF4 12.5" - 19.5" N/C

Optional:

☐ CICP Individual calf pads \$290



1. Not available on 18" W seats.

Dealer Name: _____

Dealer Account Number:

PO #:

Tag: _____

MAXX STYLE FRONT RIGGING


☐ MSFCMFP	Maxx Style Fixed Center Mount Platform ⁽¹⁾⁽²⁾	\$573
☐ MIFP	Rubber Coated Footplate (11.5" W x 10")	N/C
☐ MBIFP	Large Rubber Coated Footplate (17" w x 13"d).....	\$250
☐ MICP	Individual calf pads	N/C

Select seat pan to footplate:

☐ MSPTF3	10.25" - 14.25"	N/C
☐ MSPTF4	14.25" - 20.25"	N/C

Set at:

☐ MSA83	Set at 83°	N/C
☐ MSA90	Set at 90° ⁽³⁾	N/C
☐ MSA97	Set at 97° ⁽³⁾	N/C



The diagram illustrates the Maxx Style Fixed Center Mount Platform. It features a central vertical support structure. At the top, there is a horizontal bar with two circular components, likely for securing the seat pan. Below this, a footplate is attached to the support. The footplate has a rectangular base with a central recessed area and two side flaps. Above the footplate, there are two rectangular calf pads mounted on the support structure. The entire assembly is designed for a fixed center mount.

1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
2. 3" ground clearance required.
3. Not available on 18" W seats.

FOOT PLATFORM OPTIONS ⁽¹⁾



☐ SFFP	Small Foot Platform (12" x 8.5" deep)	\$261
☐ MFFP	Medium Platform (12" x10.75" deep).....	\$261
☐ LFFP	Large Platform (14.5" x13" deep).....	\$363
☐ XLFFP	Extra Large Platform (17" x13" deep).....	\$363
☐ XXLFFP	Extra Extra Large Platform (21" x13" deep)	\$363

1. May need to adjust leg position to avoid interference with front casters

ADDITIONAL ACCESSORIES

☐ PH	X-Grip Phone Holder	\$109
☐ UCH	Universal Cup Holder ⁽¹⁾	\$51
☐ CSCC	Cane and Crutch Holder	\$573
☐ CSO2	Non-Articulating O ² Gas Holder	\$363
☐ CFH	Flag Holder	\$94
☐ CRVMR	Rear View Mirror Right Mounted	\$65
☐ CRVML	Rear View Mirror Left Mounted	\$65
☐ TAG571850	Captain Seat Wheelchair Cover	\$152
☐ AMC	Accessory Mount with Clip	\$73
☐ CSWH	Captain Seat Walker Holder ⁽²⁾	\$145

1. Not available on the joystick side.
2. Maximum load is 20 lbs

MILITARY DECAL⁽¹⁾

☐ DECALA	Army Decal 3" Diameter.....	\$0
☐ DECALF	Air Force Decal 3" Diameter	\$0
☐ DECALN	Navy Decal 3" Diameter	\$0
☐ DECALM	Marines Decal 3" Diameter	\$0

1. Comes uninstalled off chair.

INSTALLATION

MCI	Motion Concepts Installation ⁽¹⁾	N/C
---	---	-----

1. Batteries may be required if items on Accessories order form are selected at time of installation, please call to confirm.

NOTES

[illegible]

Total Retail Value \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM