

IN V A C A R E®

POWER

At Invacare America we make Power Wheelchairs with Seating that help improve the lives of our customers. By listening to our customers and their care teams, we focus on meaningful design and purposeful technology that enhances performance.

Chair Base/Drive Type	STORM TORQUE™ RX 20 REAR Wheel Drive	ROVI X3 MID Wheel Drive	ROVI A3 MID Wheel Drive	TDX® SP2 MID Wheel Drive	AVIVA® FX FRONT Wheel Drive
Design. Technology. Performance.					
Base Weight Capacity	300 lbs. - Standard 250 lbs. - Narrow 175 lbs. - Ultra Low Mini	300 lbs.	300 lbs.	300lbs. 400lbs. HD w/Ultra Low Maxx, 450lbs w/captain seat	300 lbs.
HCPCS Codes	K0848, K0849, K0856, K0861,	K0848, K0849, K0856, K0861	K0848, K0849, K0856, K0861	K0848, K0849, K0850, K0851, K0856, K0858, K0861, K0862	K0848, K0849, K0856, K0861
Speed	Standard - 5.8 mph Performance - 7.5 mph (Std width only)	Standard - 5.0 & 5.8 mph Performance - 6.2 mph	Standard - 5.8 & mph Performance - 6.2 mph	Standard - 5.8 mph Performance - 7.5 mph	Standard - 5.8 mph Performance - 7.5 mph
Drive Wheels	14"	14"	14"	14"	14"
Caster Size	8" x 2.5" Black Casters w/Solid Inlay	6" Solid (Front/Rear)	6" Solid (Front/Rear)	6" Solid (Front/Rear)	8" x 2.5" Black Casters w/Solid Inlay
Ground Clearance	>3.2" (under battery box)	3"	3"	2.5"	3.15"
Minimum Base Turning Radius	24.8" - Narrow 25" - Standard	19.5"	20.5"	20"	26.8"
Driving Range	Up to 19.7 miles - Standard Up to 13.8 miles - Narrow	Up to 17.7 miles	Up to 17.7 miles	Up to 20 miles	Up to 20 miles
Base Width	22.4" - Narrow 24.4" - Standard	23.25"	25.5"	24" & 25.5"	24.25"
Base Length	33.5"	36.5"	36.6"	35.6"	33.5" (includes footrest)
Suspension	SureGrip™ Suspension System	Active Ride Control (ARC)	Active Ride Control (ARC)	Surestep® and Stability Lock Suspension	4Sure™ Suspension
Tracking Technology	G-Trac™ - Digital Gyroscopic Technology (Optional)	OnTraxx - Enhanced Driving Module	OnTraxx - Enhanced Driving Module	G-Trac - Digital Gyroscopic Technology (Optional)	G-Trac - Digital Gyroscopic Technology
Motor Type	DuraWatt™ 4-Pole Sealed Housing	4-Pole EAD Motors	4-Pole EAD Motors	DuraWatt 4-Pole Sealed Housing	DuraWatt 4-Pole Sealed Housing
Battery Type	Group 24 - Standard M50 - Narrow	Group 34 & Group 50	Group 34 & Group 50	22NF & Group 24	Group 24
Battery Isolation Switch	YES	NO	NO	NO	YES
Incline Climbing	9 degrees	8 degrees	8 degrees	9 degrees	9 degrees
EZ Lock Compatible	YES	YES	YES	YES	YES
Lights	LED integrated into the front caster and a rear LED strip built into the rear shroud. (standard)	Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. (optional)	Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. (optional)	LED Light Package Available for front and rear lights. (optional)	LED Light package available for front and rear lights. (optional)

Pelvic positioning strap not shown on-chair for clarity purposes.

SEATING					
Captain Seat	YES	YES	YES	YES	YES
Rehab Seating	YES	YES	YES	YES	YES
Ultra Low Mini Maxx Power Positioning System	YES	YES	NO	YES	YES
Ultra Low Maxx Power Positioning System	YES	YES	YES	YES	YES
Ultra Low Bariatric Power Positioning System	NO	NO	NO	YES/HD	NO
MPS Mini Maxx - Power Standing System	NO	YES	YES	NO	NO
MPS Maxx - Modular Power Standing System	NO	YES	YES	NO	YES
407 Maxx - Single Pivot Tilt Power Positioning System	NO	YES	YES	NO	NO

ELECTRONICS					
Technology	LiNX®	R-Net	R-Net	LiNX®	LiNX®
Driver Control Options	▶ LED Joystick ▶ 3.5" Color Touch Screen ▶ Color Touchscreen Display	▶ LED Joystick ▶ Color Remote Joystick ▶ Advanced Joystick ▶ Omni Display	▶ LED Joystick ▶ Color Remote Joystick ▶ Advanced Joystick ▶ Omni Display	▶ LED Joystick ▶ 3.5" Color Touch Screen ▶ Color Touchscreen Display	▶ LED Joystick ▶ 3.5" Color Touch Screen ▶ Color Touchscreen Display
Remote Diagnostics	MyLiNX	N/A	N/A	MyLiNX	MyLiNX
Bluetooth	YES	YES	YES	YES	YES

OTHER FEATURES					
Limited Warranty	▶ 1 Year Electronics/Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries	▶ 1 Year Electronics/Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries	▶ 1 Year Electronics/Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries	▶ 1 Year Electronics/Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries	▶ 1 Year Electronics/Motors ▶ 2 Years Seating Electronics ▶ 1 Year Base Frame ▶ 5 Years Seat Frame ▶ 13 Months Batteries

FRAME COLORS																	
Obsidian Black	Ruby Red	Denim Navy	Black Carbon (Matte)	Mellow Yellow	Motion Green	Black Carbon (Matte)	Mellow Yellow	Motion Green	Obsidian Black	Ruby Red	Denim Navy	Obsidian Black	Ruby Red	Denim Navy	Obsidian Black	Ruby Red	Denim Navy
Morning Mist (Rim Insert Color)	Orange Marmalade (Rim Insert Color)	Pink Flamingo	Orange Crush	Indigo Blue	Razzle Dazzle Pink (Metallic)	Orange Crush	Indigo Blue	Razzle Dazzle Pink (Metallic)	Morning Mist (Rim Insert Color)	Orange Marmalade (Rim Insert Color)	Pink Flamingo (Rim Insert Color)	Morning Mist (Rim Insert Color)	Orange Marmalade (Rim Insert Color)	Pink Flamingo (Rim Insert Color)	Morning Mist (Rim Insert Color)	Orange Marmalade (Rim Insert Color)	Pink Flamingo (Rim Insert Color)
Lime Twist	Purple Passion	Arctic Blue	Nordic Blue	Purple Passion	Blizzard White (Metallic)	Nordic Blue	Purple Passion	Blizzard White (Metallic)	Lime Twist	Purple Passion	Arctic Blue	Lime Twist	Purple Passion	Arctic Blue	Lime Twist	Purple Passion	Arctic Blue
Blue Peacock	Banana Split (Rim Insert Color)		Candy Apple Red (Metallic)	Graphite Grey (Metallic)	Black Panther (Gloss)	Candy Apple Red (Metallic)	Graphite Grey (Metallic)	Black Panther (Gloss)	Blue Peacock	Banana Split (Rim Insert Color)		Blue Peacock	Banana Split (Rim Insert Color)		Blue Peacock	Banana Split (Rim Insert Color)	