



Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

SAVE

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

SUBMIT

CLEAR FORM

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

905-695-0583

PRINT FORM



TDX-SP2 HD BASE WITH
MAXX BARIATRIC POWER POSITIONING SYSTEM
CDN PRICE LIST AND ORDER FORM -AADL W801HD/W990
Price Effective July 15, 2026



PRINT FORM

SAVE

CLEAR FORM

REQUIRED INFORMATION

Quote ☐ Order ☐

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: ☐M ☐F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured
please fill in all REQUIRED* measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness	



Customer Service: 1.866.748-7943 | Fax: 1.905.695.0583 | www.motionconcepts.com

INVACARE POWER BASES

☐ TDXSP2HD TDX-SP2 HD Base\$8,275

Select STF⁽¹⁾:
HD-SLOW ☐ 17.25"⁽²⁾\$250
HD-LOW ☐ 17.75"
HD-MED ☐ 18.75"
HD-TALL ☐ 19.75"

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size.
2. Must order Transport Brackets (TRBKTS). Only available on 17" seat depths and greater.

USER WEIGHT LIMITS

☐ U400 Weight Capacity up to 400 lbsSTD

TRANSPORT TIE DOWN

CLIENT COST

☐ TRBKTS Wheelchair Transport Brackets⁽¹⁾ STD
☐ TRRO Transport Ready Option⁽²⁾⁽⁴⁾ \$499
☐ TRRO-E Transport Ready Option for Elevate⁽²⁾⁽³⁾⁽⁴⁾ \$499

1. For unoccupied use only
2. TRRO meets ISO-7176 Part 19
3. Required for Elevate Systems to meet ISO-7176 Part 19, includes two additional tie-down loops on rail.
4. Must select a Motion Concepts or Matrix Elan headrest with hardware and Motion Concepts Standard Rehab, High Back or Matrix Upgradeable back. Comes with pelvic safety belt standard, please select TROBELT on page 4.

TIRE OPTIONS

☐ B1431-3 14"x3" Black Tire with Gel Foam InsertsSTD
☐ B1430-3T 14"x3" Black Tire - Pneumatic TireN/C

FORK OPTIONS

☐ DSFK Double-Sided Fork Package STD
☐ FPKPG Single-Sided Fork Package \$375

BATTERY TRAY TYPE

☐ 24TRY Group 24 Style Tray - Standard Base 25.5"N/C

FRAME FINISH

☐ 170P		Obsidian Black	AADL COLOUR.....N/C
☐ 173P		Peacock Blue	AADL COLOUR.....N/C
☐ 177P		Artic Blue	N/C
☐ 175P		Lime Twist	N/C
☐ 171P		Ruby Red	N/C
☐ 174P		Purple Passion	N/C
☐ 172P		Denim Navy	AADL COLOUR.....N/C

RIM INSERT COLORS - MUST PICK ONE

☐ 170PR		Obsidian Black	N/C
☐ 172PR		Denim Navy	N/C
☐ 173PR		Peacock Blue	N/C
☐ 177PR		Artic Blue	N/C
☐ 175PR		Lime Twist	N/C
☐ 178PR		Banana Split	N/C
☐ 179PR		Orange Marmalade	N/C
☐ 171PR		Ruby Red	N/C
☐ 176PR		Pink Flamingo	N/C
☐ 174PR		Purple Passion	N/C
☐ 180PR		Morning Mist	N/C



LIGHTS AND INDICATORS

CLIENT COST

☐ LIGHTS LED Lights and Indicators⁽¹⁾\$995

1. Must order REM216, REM400 or REM500

CONTROLLERS

☐ EXPC Expandable ControllerSTD
☐ PWH Harness Required for Expandable System N/C
☐ GTRAC-LX LINX G-TRAC Module. \$1,200

EXPANDABLE DRIVE CONTROL

☐ REM211 LED Expandable Remote/Joystick⁽¹⁾⁽²⁾ \$875
Three drive functions, control of up to 6 actuators through the drive control. Not for use with specialty controls.

☐ REM216 LED Expandable Remote/Joystick
with Lights⁽²⁾⁽³⁾\$875
Three drive functions, control of up to 6 actuators through the drive control and light operation. Not for use with specialty controls.

☐ REM400 Colour 3.5" Touch Screen
Remote/Joystick⁽²⁾⁽⁴⁾⁽⁵⁾ \$1,100
☐ SW400TGL Toggle Kit for REM400⁽⁶⁾\$200

☐ REM500 LINX Colour 3.5" Touch Screen
Display Only - No Driver Control⁽²⁾⁽³⁾⁽⁴⁾ \$1,100
Select mounting position
☐ R5MR Right
☐ R5ML Left

1. Not available with lights and indicators
2. Must choose EXPC Expandable Controller and PWH Harness.
3. Must choose lighting package separately.
4. Can control lighting, must choose lighting package separately.
5. May also be used as display for specialty controls, chose specialty controls separately.
6. Adds 2 toggle switches to REM400.

ATTENDANT DRIVE CONTROL CHOICES

PRIOR APPROVAL REQUIRED

☐ ACU Proportional Attendant Control \$1,024
Select seating style for mounting type
☐ TAM Tilt & Elevate Mount
☐ RAM Recline Mount

Select mounting position
☐ RIGHT-ACU Right
☐ LEFT-ACU Left

SPECIALTY PROPORTIONAL CONTROLS

AADL CODE - W100

☐ CREM LINX Compact Remote\$1,097
☐ CREM-LF LINX Compact Remote Low Force\$1,097

Select mounting position
☐ MR Right
☐ ML Left

1. Requires REM400 or REM500

EXPANDABLE MODULES

☐ OUT500 Output Module \$645
☐ IN200 LiNX Switch Input Module with
Power ON/OFF for REM200 Series⁽¹⁾ \$645
☐ IN500 LiNX 9 Pin Input Module
with Sip-N-Puff Nozzle⁽²⁾⁽³⁾⁽⁴⁾\$1,800
AADL CODE - W100.....
☐ PKG32666 Therafin Sip-N-Puff Breath Tube Kit
AADL CODE - W100..... \$440

1. Not for use with Alternative Drive Controls.
2. Only one input module per chair, requires REM400 or REM500.
3. Input module also required for any 9-pin specialty control.
4. 9 pin specialty control and Sip-N-Puff cannot both be used together.

MOUNTING FOR DRIVER CONTROLS

☐ FM-ML Fixed Mounting LeftN/C
☐ FM-MR Fixed Mounting RightN/C
☐ SAQR Motion Height Adjustable Swing-Away
Quad Link (Right Hand)⁽¹⁾ \$344
☐ SAQL Motion Height Adjustable Swing-Away
Quad Link (Left Hand)⁽¹⁾ \$344
☐ MRR Maxx Resolve Swing-Away Joystick Mount
(Right Hand)⁽¹⁾ \$450
☐ MRL Maxx Resolve Swing-Away Joystick Mount
(Left Hand)⁽¹⁾ \$450

1. Not available with the REM500.

LINX ELECTRONIC ACCESSORIES

☐ LAK LINX Access Key⁽¹⁾
NOT COVERED UNDER AADL\$50
☐ MUSB Motion Dual USB ChargerN/C

1.Required for programming

JOYSTICK TOPS

☐ PC101A Bodypoint U Shaped Handle 3" \$120
☐ PC102A Bodypoint U Shaped Handle 4" \$120
☐ PC107A Bodypoint Rubber Dome..... \$120
☐ 1560 T Handle Flexible Joystick Extension \$120
☐ 1561 Straight Handle Flexible Joystick Extension \$120
☐ 1826 Chin Cup \$120

1. To remove/disengage the joystick knob from the REM400 Remote pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

INSTALLATION OF BATTERIES

☐ 22NFBATTERY 22 NF Battery (On Chair) AADL PRICING
☐ 24BATTERY 24 Gel Battery (On Chair) AADL PRICING

BATTERY CHARGER OPTION

☐ 110CHARGER 110 Volt Battery Charger.....N/C
☐ L8900 Omit Battery ChargerN/C

Maxx Bariatric Power Positioning
on Invacare TDX-SP2 HD Base

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

BARIATRIC REHAB SEAT

☐ MBRS Maxx Bariatric Rehab Seat \$4,995

BARIATRIC POWER POSITIONING SYSTEMS

☐ MBT Maxx Bariatric CG Tilt \$7,995
☐ MBTRE⁽¹⁾ Maxx Bariatric CG Tilt (45°) & Power Recline
(150°) w/Extended Shear Reduction(ESR) \$12,595
☐ MBRESR Modular Maxx Power Recline (150°)
w/Extended Shear Reduction(ESR) \$10,350

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module

ADD: HD ELEVATING LIFT SEAT MODULE

CLIENT COST

☐ ESM10 10" Elevating Lift Seat Module⁽¹⁾⁽²⁾ \$2,995
Weight capacity up to 300lbs

1. Only Available with Tilt. Uses an integrated Tilt/Lift mechanism, tilt is 45°
2. Not available on seat sizes greater than 22" wide by 22" deep

MODULES MUST SELECT SYSTEM TO ADD MODULES

☐ PRETILT 5° Fixed Pre-Tilt⁽¹⁾ N/C
☐ 8PR 8° Precline Module⁽¹⁾ only available with Recline N/C
☐ 30PR 30° Precline Module⁽¹⁾ only available with Recline \$500

1. Total degrees of recline is 150° ie. if 30° of precline is used, then 120° of recline will be available.

ADULT SEAT WIDTH SETTINGS

☐ BW120 20" Wide (Adjustable from 20" to 24") N/C
☐ BW121 21" Wide (Adjustable from 20" to 24") N/C
☐ BW122 22" Wide (Adjustable from 20" to 24") N/C
☐ BW123 23" Wide (Adjustable from 20" to 24") N/C
☐ BW124 24" Wide (Adjustable from 20" to 24") N/C
☐ BW222 22" Wide (Adjustable from 22" to 26") N/C
☐ BW223 23" Wide (Adjustable from 22" to 26") N/C
☐ BW224 24" Wide (Adjustable from 22" to 26") N/C
☐ BW225 25" Wide (Adjustable from 22" to 26") N/C
☐ BW226 26" Wide (Adjustable from 22" to 26") N/C

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

☐ SD117 17" Deep (Adjustable from 17" to 21") N/C
☐ SD118 18" Deep (Adjustable from 17" to 21") N/C
☐ SD119 19" Deep (Adjustable from 17" to 21") N/C
☐ SD120 20" Deep (Adjustable from 17" to 21") N/C
☐ SD121 21" Deep (Adjustable from 17" to 21") N/C
☐ SD219 19" Deep (Adjustable from 19" to 23") N/C
☐ SD220 20" Deep (Adjustable from 19" to 23") N/C
☐ SD221 21" Deep (Adjustable from 19" to 23") N/C
☐ SD222 22" Deep (Adjustable from 19" to 23") N/C
☐ SD223 23" Deep (Adjustable from 19" to 23") N/C

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX SEATING OPTIONS

☐ Matrix Cushions available on page 10

AADL SEATING BENEFIT



CHEST STRAPS

AADL SEATING BENEFIT

☐ MFUM BP Monoflex, Underarm RLS M(4" x 18.25") \$195
☐ MFUL BP Monoflex, Underarm RLS L(4.5" x 21.5") \$195
☐ MFCM BP Monoflex, Center RLS M(3.25" x 18.25") \$195
☐ MFCL BP Monoflex, Center RLS L(3.5" x 21.5") \$195

LAP BELTS

☐ BELT60 Push Button Style Lap Belt, Length 60" N/C
☐ BELT71 Push Button Style Lap Belt, Length 71" N/C
☐ PLB2 Bodypoint Padded Lap Belt 2 Point⁽¹⁾ \$250
☐ PLB4 Bodypoint Padded Lap Belt 4 Point⁽¹⁾ \$315

1. 67" length only option

MOTION CONCEPTS HEADREST OPTIONS

☐ MCSH Motion Concepts Standard⁽¹⁾ N/C
☐ MCOH Motion Concepts Auto Style⁽¹⁾ \$235
☐ FM Fixed Mounting (non-removable)⁽²⁾ N/C
☐ OHR Omit: Headrest N/C

1. Must select Fixed Mounting (FM) or Matrix Headrest Hardware on page 12.
2. Only compatible with Standard Rehab, Contour or High Back Models. All other upgradeable backs require removeable hardware.

MATRIX HEADRESTS OPTIONS

☐ Matrix Elan Headrests and Hardware
available on page 12



FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ FOR TILT ONLY/ELEVATING SEATING SERIES

☐ 18" ☐ 19" ☐ 20" ☐ 21" ☐ 22" ☐ 23" ☐ 24" ☐ 25"
FBT18 FBT19 FBT20 FBT21 FBT22 FBT23 FBT24 FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ FOR RECLINE OR TILT/RECLINE

☐ 20" ☐ 21" ☐ 22" ☐ 23" ☐ 24" ☐ 25" ☐ 26" ☐ 27"
FBR20 FBR21 FBR22 FBR23 FBR24 FBR25 FBR26 FBR27

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK CANES (NON RECLINE SYSTEMS)

☐ SBC18	18" Straight Back Canes	STD
☐ SBC20	20" Straight Back Canes	STD
☐ SBC22	22" Straight Back Canes	STD
☐ SBC24	24" Straight Back Canes	STD
Select Cane Angle ☐ 75° ☐ 79° ☐ 85° ☐ 90° ☐ 95° ☐ 101° ☐ 106° ☐ 112° ☐ 116° BA75 BA79 BA85 BA90 BA95 BA101 BA106 BA112 BA116		
☐ ABC20	20" Angled (5°) Back Canes	N/C
☐ ABC22	22" Angled (5°) Back Canes	N/C
☐ ABC24	24" Angled (5°) Back Canes	N/C
Select Cane Angle ☐ 81° ☐ 84° ☐ 91° ☐ 96° ☐ 101° ☐ 106° ☐ 111° ☐ 117° ☐ 121° BA81 BA84 BA91 BA96 BA101 BA106 BA111 BA117 BA121		
☐ EBC22	22" Endomorph Back Canes ⁽¹⁾	\$350
☐ EBC24	24" Endomorph Back Canes ⁽¹⁾	\$350
☐ EBCN1	1" Narrower	N/C
☐ EBCN2	2" Narrower	N/C
☐ EBCN4	4" Narrower	N/C

1. Comes standard with Cantilever Armrest Cane Mounted.

CANE OPTIONS FOR RECLINE SYSTEMS

☐ BFS2	Recline Body Fit 2" Narrower ⁽¹⁾	\$995
☐ BFS4	Recline Body Fit 4" Narrower ⁽¹⁾	\$995
☐ AHH	Adjustable Hooking Handles ⁽²⁾ (pair)	\$450
☐ RPHS	Push Handles for Recline - Straight ⁽³⁾ (pair)	\$350
☐ RPHA	Push Handles for Recline - Angled ⁽³⁾ (pair)	\$350

1. Minimum finished back widths of 16".
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

BACK OPTIONS

☐ SRB	Standard Rehab Back & Cushion	N/C
☐ PIPF	Flat Planar Interface Plate ⁽¹⁾	N/C
☐ PIPR	Recessed Planar Interface Plate ⁽¹⁾⁽²⁾	N/C
☐ ORBC	Omit Rigid Back & Cushion ⁽³⁾	N/C

1. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
2. Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
3. Not available with Recline/ESR.

MATRIX BACKREST OPTIONS

AADL SEATING BENEFIT

- ☐ Matrix Backs available on page 11



ARMRESTS (NON RECLINE SYSTEMS)

☐ DP1	Adult Dual Post Adj. 9.75"-12.75"	STD
☐ DP2	Adult Dual Post Adj. 12.75"-15.75"	STD
☐ DP3	Adult Dual Post Adj. 15.75"-18.75"	STD
Select:		
☐ OAL	ADD: OutBack Arm Option Left ⁽¹⁾ each	\$195
☐ OAR	ADD: OutBack Arm Option Right ⁽¹⁾ each	\$195
☐ OCTA1	Outback Cantilever Tilt Armrest 9"-12" ⁽³⁾	\$488
☐ OCTA2	Outback Cantilever Tilt Armrest 12"-15" ⁽³⁾	\$488
☐ RMCAN	Ultra Rail Seat Mnt Flip Back Cantilever MAXX Tilt Arm (9-12.5") ⁽²⁾	\$488
☐ DPM	ADD: Dual Post Module to RMCAN (pair)	\$250
☐ CAFEB	Cantilever Arms for Endomorph Back Canes (pair) (see page 4)	STD

1. Only available with Dual Post Armrests.
2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.
3. Not available with the Universal Cup Holder or Flip and Fold Lap Tray.

RECLINE ARMRESTS

☐ RC1	Reclining Adjustable Height. 9.5"-13"	STD
☐ RC2	Reclining Adjustable Height. 13"-16"	STD
Select:		
☐ OAL	ADD: OutBack Arm Option Left ⁽¹⁾ each	\$195
☐ OAR	ADD: OutBack Arm Option Right ⁽¹⁾ each	\$195
☐ OCRA1	Outback Cantilever Recline Armrest 9"-12" ⁽²⁾	\$488
☐ OCRA2	Outback Cantilever Recline Armrest 12"-15" ⁽²⁾	\$488
☐ MCAN1	Maxx Style Cane Mounted Height, Width and Angle Adjustable Cantilever Flip Back Armrests for Recline Systems 9.5"-13" ⁽²⁾	\$488
☐ MCAN2	Maxx Style Cane Mounted Height, Width and Angle Adjustable Cantilever Flip Back Armrests for Recline Systems 12.5"-16" ⁽²⁾	\$488
☐ CANQ1	Cane Mounted Height & Angle Adjustable Cantilever Flip Back Armrest 9.5"-12" ⁽³⁾	\$488
☐ CANQ2	Cane Mounted Height & Angle Adjustable Cantilever Flip Back Armrest 10"-13" ⁽³⁾	\$488
☐ CANQ3	Cane Mounted Height & Angle Adjustable Cantilever Flip Back Armrest 13"-16" ⁽³⁾	\$488

1. Only available with Reclining Armrests.
2. Not available with the Universal Cup Holder or Flip and Fold Lap Tray.
3. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.

ARMREST TUBE LENGTH - PLEASE SPECIFY

☐ ATFL	Left Full	STD
☐ ATDL	Left Desk	STD
☐ ATFR	Right Full	STD
☐ ATDR	Right Desk	STD

ARMPADS

Standard

☐ SAPFL	Left Full (2.25" x 14")	STD
☐ SAPDL	Left Desk (2.25" x 10")	STD
☐ SAPFR	Right Full (2.25" x 14")	STD
☐ SAPDR	Right Desk (2.25" x 10")	STD

Waterfall Arm Pads⁽¹⁾

☐ WAPFL	Left Full (13")	N/C
☐ WAPDL	Left Desk (9")	N/C
☐ WAPFR	Right Full (13")	N/C
☐ WAPDR	Right Desk (9")	N/C

Modular Arm Pad - Left⁽²⁾

AADL SEATING BENEFIT

☐ MAPTFL	Left Tray Full (14")	\$95
☐ MAPTDL	Left Tray Desk (10")	\$95

Select Pad Insert below:

☐ ISSFL	Iskin - Integrated soft skin foam
☐ SCVFL	Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

AADL SEATING BENEFIT

☐ MAPTFR	Right Tray Full (14")	\$95
☐ MAPTDR	Right Tray Desk (10")	\$95

Select Pad Insert below:

☐ ISSFR	Iskin - Integrated soft skin foam
☐ SCVFR	Startex Covered Visco foam

Flat Multi Position Arm Pads

CLIENT JUSTIFICATION

☐ MPAFL	Left Full (14")	\$150
☐ MPADL	Left Desk (10")	\$150
☐ MPAFR	Right Full (14")	\$150
☐ MPADR	Right Desk (10")	\$150

Gel Pads

AADL SEATING BENEFIT

☐ GAPNL	Narrow Left (2" x 12")	\$88
☐ GAPNR	Narrow Right (2" x 12")	\$88
☐ GAPWL	Wide Left (3.5" x 12")	\$96
☐ GAPWR	Wide Right (3.5" x 12")	\$96
☐ GAPLL	Long Left (3.5" x 14")	\$103
☐ GAPLR	Long Right (3.5" x 14")	\$103

Arm Troughs

AADL SEATING BENEFIT

☐ EATL	Ergonomic Arm Trough - Left	\$121
☐ EATR	Ergonomic Arm Trough - Right	\$121
☐ FHPL	Flat Hand Pad - Left ⁽³⁾	\$91
☐ FHPR	Flat Hand Pad - Right ⁽³⁾	\$91

Arm Pad Storage/Cup Holder⁽⁴⁾

CLIENT COST

☐ APSL	Left	\$295
☐ APSR	Right	\$295

ARMPADS - CONTINUED

☐ OMIT	Omit Armpads(7/8" diameter tube)	N/C
-------------------------------	--	-----

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side.

ADDITIONAL ARM ACCESSORIES

☐ UCH	Universal Cup Holder ⁽¹⁾	\$35
☐ MACESL	Multi Axis Upper Extremity Support ⁽²⁾ (each) Left Side	\$225
☐ MACESR	Multi Axis Upper Extremity Support ⁽²⁾ (each) Right Side	\$225
☐ OAS2L	1" Offset Armrest Spacer Left(each)	\$125
☐ OAS2R	1" Offset Armrest Spacer Right(each)	\$125
☐ OAS1L	1/2" Offset Armrest Spacer Left(each)	\$125
☐ OAS1R	1/2" Offset Armrest Spacer Right(each)	\$125
☐ MPAPA	Multi Purpose Arm Pad Adapter Block ⁽³⁾ (each) Quantity ()	\$35

1. Not available on the joystick side.
2. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
3. Not available with Standard Arm Pads or Waterfall Arm Pads.

ELBOW BLOCKS

☐ Elbow Blocks available on page 14
--

POWER CENTER MOUNT FOOT PLATFORM

LNX Power Center Mount Foot Platform⁽¹⁾


CLIENT COST

☐ LNX	13" - 19" Seat pan to footplate.	\$3,200
	With Quick Release Height Adjustment	
☐ IFP	Comes with a rubber coated footplate	
	(11.5" W x 10")	N/C
☐ ICP	Comes standard with individual calf pads.	N/C
☐ OPCP	One Piece Malleable Calf Panel	\$250

1. Comes standard with moulded i-skin calf pads.

MANUAL CENTER MOUNT FRONT RIGGING

☐ MSFCMFP	Maxx Style Fixed Center Mount Foot Platform ⁽¹⁾	\$395
	With Quick Release Height Adjustment	
☐ FBIFP	Comes with a standard footplate with	
	(17" w x 13"d)	N/C
☐ ICP	Individual calf pads	\$200
Select seat pan to footplate:		
☐ SPTF1	9" - 13"	N/C
☐ SPTF14	13" -19"	N/C
Set at:		
☐ SA70	Set at 70°	N/C
☐ SA90	Set at 90°	N/C
☐ SA97	Set at 97°	N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.


FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP



☐ MFFP	Medium Platform ⁽¹⁾ (12"w x 10.75"d)	\$180
☐ LFFP	Large Platform ⁽¹⁾ (14.5"w x 13"d)	\$250
☐ XLFFP	Extra Large Platform ⁽¹⁾ (17"w x 13"d)	\$250
☐ XXLFFP	Extra Extra Large Platform ⁽¹⁾ (21"w x 13"d)	\$250

1. May need to adjust leg position to avoid interference with front castors

INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP



☐ IFSL	Small ⁽¹⁾ - (5" x 7.5" deep) - Left	\$165
☐ IFSR	Small ⁽¹⁾ - (5" x 7.5" deep) - Right	\$165
☐ IFML	Medium ⁽¹⁾ - (5.5" x 9.5" deep) - Left	\$165
☐ IFMR	Medium ⁽¹⁾ - (5.5" x 9.5" deep) - Right	\$165
☐ IFLL	Large ⁽¹⁾ - (6" x 11.5" deep) - Left	\$165
☐ IFLR	Large ⁽¹⁾ - (6" x 11.5" deep) - Right	\$165
☐ MFOL	Multi-Axis Optional Footplate Left (6"W x 10"D)	\$225
☐ MFOR	Multi-Axis Optional Footplate Right (6"W x 10"D)	\$225

1. Foot plates may interfere with casters in certain configurations

MAXX STYLE POWER ELEVATING LEGRESTS

CLIENT COST

☐ MSPPL MAXX Swing Away
Power Elevating⁽¹⁾⁽²⁾\$2490

Select Power Legrest Function:

☐ INDIA Individual Movement.....N/C
☐ COMA Combined Movement.....N/C
☐ SYNA Synchronized/Recline Movement.....N/C

Select Receivers:

☐ RS Straight Receivers.....N/C
☐ RF Flared 1" Outward Receivers.....N/C

Select Seat Pan to Footplate:

☐ SPTF5 Length: 10"-14".....N/C
☐ SPTF2 Length: 13"-17".....N/C
☐ SPTF7 Length: 16"-20".....N/C

☐ FCPMS Flat Calf Pad Option (pair) \$275

☐ CCPMS Curved Calf Pad Option (inside with 6") (pair) \$275

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads.

SWING AWAY LEGREST

☐ HD70 Heavy Duty 70° Swing Away Footrest(pair)⁽¹⁾⁽²⁾\$495

Select Receivers:

☐ RS Straight Receivers.....N/C
☐ RF Flared 1" Outward Receivers.....N/C

Select Seat Pan to Footplate:

☐ SPTF11 Length: 6"-10".....N/C
☐ SPTF8 Length: 11"-15".....N/C
☐ SPTF2 Length: 13"-17".....N/C
☐ SPTF7 Length: 16"-20".....N/C

☐ HDFCP Flat Calf Pad Option (pair) \$275

☐ HD70CCP Curved Calf Pad Option (inside with 6") (pair) \$275

☐ DCS Dynamic Coil Springs (Pair) for Heavy Duty 70°
Swingaway Footrests \$190

1. Must choose footplate on page 8
2. Heavy duty durable construction (4.5lbs ea.).

MAXX STYLE MANUAL ELEVATING LEGRESTS

☐ MMPPSA MAXX Swing Away
Manual Elevating⁽¹⁾⁽²⁾ \$680

Select Receivers:

☐ RS Straight Receivers.....N/C
☐ RF Flared 1" Outward Receivers.....N/C

Select Seat Pan to Footplate:


☐ SPTF5 Length: 10"-14".....N/C
☐ SPTF2 Length: 13"-17".....N/C
☐ SPTF7 Length: 16"-20".....N/C

☐ FCPMS Flat Calf Pad Option (pair) \$275

☐ CCPMS Curved Calf Pad Option (inside with 6") (pair) \$275

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads

FOOTPLATE FOR SWING AWAY/ELEVATING LEGS



1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

ADJUSTABLE ANGLE FOOTPLATES

PRIOR APPROVAL

☐ AFSL	Left Small Footplate (4.25"W x 8"D)	\$165
☐ AFSR	Right Small Footplate (4.25"W x 8"D)	\$165
☐ AFML	Left Medium Footplate (5.25"W x 8"D)	\$165
☐ AFMR	Right Medium Footplate (5.25"W x 8"D)	\$165
☐ AFLL	Left Large Footplate (6.25"W x 8"D)	\$165
☐ AFLR	Right Large Footplate (6.25"W x 8"D)	\$165
☐ AFXL	Left Extra Large Footplate (9"W x 11"D)	\$165
☐ AFXR	Right Extra Large Footplate (9"W x 11"D)	\$165

MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES


PRIOR APPROVAL

☐ MFSL	Left Small Footplate (4.25"W x 8"D)	\$225
☐ MFSR	Right Small Footplate (4.25"W x 8"D)	\$225
☐ MFML	Left Medium Footplate (5.25"W x 8"D)	\$225
☐ MFMR	Right Medium Footplate (5.25"W x 8"D)	\$225
☐ MFLL	Left Large Footplate (6.25"W x 8"D)	\$225
☐ MFRL	Right Large Footplate (6.25"W x 8"D)	\$225
☐ MFXL	Left Extra Large Footplate (9"W x 11"D)	\$225
☐ MFXR	Right Extra Large Footplate (9"W x 11"D)	\$225

HEEL LOOPS, CUPS AND TOE STRAPS

☐ HLL	Left Heel Loop (each)	\$20
☐ HLR	Right Heel Loop (each)	\$20
☐ HCL	Left Cup Loop (each)	\$20
☐ HCR	Right Cup Loop (each)	\$20
☐ TSL	Left Toe Strap (each)	\$20
☐ TSR	Right Toe Strap (each)	\$20

FOOTPLATE/CUT OUTS SIZE SELECTION



1. Certain footplate sizes may not fit all seat widths depending on legrests selected

ADJUSTABLE ANGLE FOOTPLATES/CUT OUTS

PRIOR APPROVAL

☐ AFCSL	Left Small Footplate (7.25"W x 9"D)	\$165
☐ AFCSR	Right Small Footplate (7.25"W x 9"D)	\$165
☐ AFCML	Left Medium Footplate (7.5"W x 11"D)	\$165
☐ AFCMR	Right Medium Footplate (7.5"W x 11"D)	\$165
☐ AFCLL	Left Large Footplate (10"W x 14"D)	\$165
☐ AFCLR	Right Large Footplate (10"W x 14"D)	\$165

MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES/ CUT OUTS

PRIOR APPROVAL

☐ MFCSL	Left Small Footplate (7.25"W x 9"D)	\$225
☐ MFCSR	Right Small Footplate (7.25"W x 9"D)	\$225
☐ MFCML	Left Medium Footplate (7.5"W x 11"D)	\$225
☐ MFCMR	Right Medium Footplate (7.5"W x 11"D)	\$225
☐ MFCLL	Left Large Footplate (10"W x 14"D)	\$225
☐ MFCLR	Right Large Footplate (10"W x 14"D)	\$225

CALF STRAPS AND CALF PANELS

☐ CSL	Body Point Airmesh Calf Strap for Chairs 18"-22" wide with 3.5" H Pad ⁽¹⁾	\$200
☐ CPL	Body Point Airmesh Calf Panel for Chairs 18"-22" wide with 38.5" W x 9" H Pad ⁽¹⁾	\$315

THORACIC LATERALS AND HIP/KNEE SUPPORTS

☐ Matrx Laterals please see page 12

☐ Maxx Laterals and Hip/Knee supports please see page 13

matrx
SEATING SERIES

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ELECTRONIC SWITCH MOUNTING

- ☐ ESML Left Side (Armpad)
☐ ESMR Right Side (Armpad)
☐ ESMJ Joystick (Mounts to Joystick Bracket)
☐ ESMD Do not mount

SINGLE-DUAL FUNCTION ELECTRONICS

- ☐ SFCB Single/Dual Function Electronics⁽¹⁾ STD

Select Switch: ☐ SPBO Single Push Button

☐ STO Single Toggle

☐ DPBO Dual Push Button

☐ SOS Omit Switch

☐ JOY Through Joystick

☐ SC Through Specialty Control

Optional:

☐ RDE Reduced Drive
for Elevating Seat⁽²⁾ **CLIENT COST** \$400

☐ SPBSA Additional Single Push Button \$97

☐ STA Additional Single Toggle Switch \$97

☐ DPBA Additional Dual Push Button \$97

1. These electronics include drive lock-out and tilt limit.

2. If reduced drive is not selected the system will be in drive lock-out when elevated.

MULTI- FUNCTION ELECTRONICS

- ☐ MFCB Multi-Function through
Switch/Toggle⁽¹⁾ \$1,950

Select Switch:

☐ QPBO Quad Push Button

☐ FWTO Four-Way Toggle

☐ EWRO 8-Way Rocker⁽²⁾

- ☐ MFCBD Multi-Function through
Drive Control⁽¹⁾ \$1,950

Select Control:

☐ JOY Through Joystick

☐ SC Through Specialty Control

Optional:

☐ RDY Include Reduced Drive for Elevating Seat⁽³⁾
☐ RDN Do not Induce Reduced Drive for Elevating
Seat⁽³⁾

☐ QPBA Additional Quad Button \$97

☐ FWTA Additional Four-Way Toggle Switch \$97

☐ EWRA Additional 8-Way Rocker \$150

1. These electronics include seating attendant control, drive lock-out, tilt/recline limit
and actuator speed control.

2. Dual function is not available with 8-Way Rocker.

3. If reduced drive is not selected system will be in drive lock-out when elevated.

ELECTRONICS OPTIONS

- ☐ M616S M616 Converter⁽¹⁾ \$125
- ☐ FTSQ Feather Touch Switch Single Function
Quarter Size \$295
- ☐ FTMD Feather Touch Switch Multi-Function
Dual Dime Size \$295
- ☐ FTMQ Feather Touch Switch Multi-Function
Quad Dime Size \$295
- ☐ ASL320R ASL Pulse Mechanical Switch Red \$125
- ☐ ASL320G ASL Pulse Mechanical Switch Green \$125
- ☐ ASL320BK ASL Pulse Mechanical Switch Black \$125
- ☐ SESG Stealth Egg Switch - Single Function Green \$135
- ☐ SESB Stealth Egg Switch - Single Function Black \$135
- ☐ PSOS Power Shut Off Switch - Toggle⁽¹⁾ \$175
- ☐ LJS Lazarus Jump Starter⁽¹⁾ \$150
- ☐ APPS Accessory Port Power Supply⁽¹⁾
(3 amp capacity) \$195
- ☐ TPSMH 12 Piece Switch Mounting Hardware
Kit (additional kits) \$39

1. Not available with Single-Dual Function Electronics.

VENT TRAYS, O2 HOLDERS, ACCESSORIES

- ☐ For O2 Holders and Additional Accessories
Please see page 14

INSTALLATION

- ☒ MCI Motion Concepts Installation Included

PLEASE SPECIFY PAGES INCLUDED WITH ORDER

- ☐ Pages 1-9 Excludes Matrix and Maxx Seating Accessories
Thoracic Supports, Hip Supports, Vent Trays and
Accessories
- ☐ Pages 1-14 Complete Order Form

Total Retail Value \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM



AADL SEATING BENEFIT

CUSHION COVER REFERENCE GUIDE

Meshtex	Startex	Infection Control			
Spacetex	Onyx				
Cushion	Meshtex	Startex	Infection	Spacetex	Onyx
Libra/Adjust	N/A	STD	\$50	N/C	N/A
PS/PSP	N/A	STD-Reversible	STD-Reversible	N/A	N/A
Vi	STD	\$103	\$103	N/A	\$103

MATRIX SEAT CUSHION OPTIONS

Matrix Libra Adjust Seat Cushion - 20" wide \$985

☐ LC2016-ADJ ☐ LC2018-ADJ ☐ LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" depth \$1,265

☐ LC2022-ADJ

Matrix Libra Adjust Seat Cushion 21" - 24" wide \$1,265

☐ LC2118-ADJ ☐ LC2120-ADJ ☐ LC2122-ADJ
☐ LC2218-ADJ ☐ LC2220-ADJ ☐ LC2222-ADJ
☐ LC2318-ADJ ☐ LC2320-ADJ ☐ LC2322-ADJ
☐ LC2418-ADJ ☐ LC2420-ADJ ☐ LC2422-ADJ

Matrix Libra Seat Cushion - 20" wide \$745

☐ LC2016 ☐ LC2018 ☐ LC2020

Matrix Libra Seat Cushion 22" depth \$1,025

☐ LC2022

Matrix Libra Seat Cushion 21" - 24" wide \$1,025

☐ LC2118 ☐ LC2120 ☐ LC2122 ☐ LC2218 ☐ LC2220
☐ LC2320 ☐ LC2322 ☐ LC2418 ☐ LC2420
☐ LC2422

Matrix PS Seat Cushion 20" wide \$569

☐ PS2016 ☐ PS2018 ☐ PS2020

Matrix PS HD Seat Cushion 22" - 26" wide \$975

☐ PS2218 ☐ PS2220 ☐ PS2222 ☐ PS2418 ☐ PS2420
☐ PS2422 ☐ PS2618 ☐ PS2620 ☐ PS2622

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX SEAT CUSHION OPTIONS CONTINUED

Matrix PSP Seat Cushion - 20" wide \$645

☐ PSP2016 ☐ PSP2018 ☐ PSP2020

Matrix PSP HD Seat Cushion 22" - 22" wide \$995

☐ PSP2218 ☐ PSP2220 ☐ PSP2222 ☐ PSP2418
☐ PSP2420 ☐ PSP2422 ☐ PSP2618 ☐ PSP2620
☐ PSP2622

Matrix Vi Seat Cushion 20" - 22" wide \$681

☐ MA2016-VI ☐ MA2018-VI ☐ MA2020-VI ☐ MA2218-VI
☐ MA2220-VI

Matrix Vi HD Seat Cushion up to 26" wide \$840

(For cushions 21" width and depth up to 26")

Specify size: _____

Matrix Multi Cushion 14" - 20" wide/reversible cover

HCPCS code E2607 \$445

☐ MM2016 ☐ MM2018 ☐ MM2020

FABRIC OPTIONS FOR CUSHIONS

☐ SFLF Spacetex Fabric Option for
Libra Cushion N/C

☐ ICFL Infection Control Fabric Upgrade for
Libra Cushion \$50

☐ OCUMA Onyx Fabric Upgrade for
Matrix Vi Cushion \$103

☐ SFVS Smooth side Startex Fabric Upgrade for
Matrix Vi Cushion \$103

☐ SFVF Fabric side Startex Fabric Upgrade for
Matrix Vi Cushion \$103



AADL SEATING BENEFIT

BACK COVER REFERENCE GUIDE

	Meshtex	Infection Control	Onyx
Back	Meshtex	Infection	Onyx
High	STD	N/A	\$103
Elite/ Elite HD/Elite Deep	STD	\$103	N/A
Elite TR/ Elite TR HD	STD	\$103	N/A
PB/PB HD/PB Deep/PB HD Deep	STD	\$103	N/A

UPGRADABLE BACK OPTIONS

High Back 20" - 24" wide\$595

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ HB2016 | ☐ HB2018 | ☐ HB2020 | ☐ HB2116 |
| ☐ HB2118 | ☐ HB2120 | ☐ HB2216 | ☐ HB2218 |
| ☐ HB2220 | ☐ HB2316 | ☐ HB2318 | ☐ HB2320 |
| ☐ HB2416 | ☐ HB2418 | ☐ HB2420 | |

Matrix Elite Back 20" wide (3" Contour)\$766

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBE2014 | ☐ PBE2016 | ☐ PBE2018 | ☐ PBE2020 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrix Elite HD Back 21" - 24" wide (3" Contour).....\$1,095

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE2116-HD | ☐ PBE2118-HD | ☐ PBE2120-HD | ☐ PBE2216-HD |
| ☐ PBE2218-HD | ☐ PBE2220-HD | ☐ PBE2316-HD | ☐ PBE2318-HD |
| ☐ PBE2320-HD | ☐ PBE2416-HD | ☐ PBE2418-HD | ☐ PBE2420-HD |

Matrix Elite HD Back 26" wide (3" Contour).....\$1,666

- | | |
|-------------------------------------|-------------------------------------|
| ☐ PBE2616-HD | ☐ PBE2620-HD |
|-------------------------------------|-------------------------------------|

Matrix Elite Deep Back 20" wide (6" Contour)\$918

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ EDB2014 | ☐ EDB2016 | ☐ EDB2018 | ☐ EDB2020 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrix Elite Deep HD Back 21" - 24" wide (6" Contour).....\$1,095

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|---------------------------------------|
| ☐ EDB2116-HD | ☐ EDB2118-HD | ☐ EDB2120-HD | ☐ EDB2216-HD |
| ☐ EDB2218-HD | ☐ EDB2220-HD | ☐ EDB2316-HD | ☐ EDB2318-HD |
| ☐ EDB2320-HD | ☐ EDB2416-HD | ☐ EDB2418-HD | ☐ EDB2420-HD.. |

Matrix Elite TR Back 20" wide (3" Contour)\$866

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE2016-TR | ☐ PBE2018-TR | ☐ PBE2020-TR |
|-------------------------------------|-------------------------------------|-------------------------------------|

Matrix Elite TR HD Back 21" - 24" wide (3" Contour).....\$1,151

- | | | |
|---------------------------------------|---------------------------------------|---------------------------------------|
| ☐ PBE2120-TRHD | ☐ PBE2220-TRHD | ☐ PBE2320-TRHD |
| ☐ PBE2420-TRHD | | |

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

UPGRADABLE BACK OPTIONS

Matrix PB Posture Back - 20" wide (5" Contour)\$734

- | | | |
|----------------------------------|----------------------------------|----------------------------------|
| ☐ MPB2012 | ☐ MPB2016 | ☐ MPB2020 |
|----------------------------------|----------------------------------|----------------------------------|

Matrix PB Posture Back Deep - 20" wide (7" Contour).....\$892

- | | | |
|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBD2012 | ☐ PBD2016 | ☐ PBD2020 |
|----------------------------------|----------------------------------|----------------------------------|

Matrix PB Heavy Duty Back 20" - 24" wide (5" Contour).....\$1,051

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ HDB2016 | ☐ HDB2020 | ☐ HDB2216 | ☐ HDB2220 |
| ☐ HDB2416 | ☐ HDB2420 | | |

Matrix PB Heavy Duty Deep Back 20"-22" wide (5" Contour).....\$1,209

- | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ DHDB2016 | ☐ DHDB2020 | ☐ DHDB2216 | ☐ DHDB2220 |
| ☐ DHDB2416 | ☐ DHDB2420 | | |

OTHER UPGRADABLE BACK OPTIONS

- | | |
|--------------------------------|--|
| ☐ EBSS | Super Soft Foam for Elite Back Upgrade. \$103 |
| ☐ OCUHB | Onyx Fabric Upgrade for High Backs..... \$103 |
| ☐ WSCU | Infection Control Cover for Elite/PB Backs..... \$103
purchased with back |
| ☐ WSC | Infection Control Cover for Elite/PB Backs..... \$150
purchased with without back |
| ☐ CHIEB | Chest Harness Interface Kit for Elite Back ⁽¹⁾⁽²⁾ \$103 |
| ☐ CHIEB | Chest Harness Interface Kit for PB Back ⁽¹⁾⁽²⁾ \$52 |

1. Does not include the chest strap.
2. This Interface Kit is only Required for 4 point harness



MATRIX ELAN HEADREST PAD⁽¹⁾⁽²⁾

☐ ESP6	Matrx Elan Standard Pad 6"	\$310
☐ ESP10	Matrx Elan Standard Pad 10"	\$310
☐ ESP14	Matrx Elan Standard Pad 14"	\$310
☐ EOP9	Matrx Elan Occipital Pad 9"	\$341
☐ EOP12	Matrx Elan Occipital Pad 12"	\$341
☐ E4PTSM	Matrx Elan 4-Point Pad 9"W x 7"H.....	\$448
☐ E4POINT	Matrx Elan 4-Point Pad 11"W x 10"H	\$448
☐ ESP6-IC	Matrx Elan Standard Pad 6" with Infection Control	\$335
☐ ESP10-IC	Matrx Elan Standard Pad 10" with Infection Control	\$335
☐ ESP14-IC	Matrx Elan Standard Pad 14" with Infection InfectionControl	\$335
☐ EOP9-IC	Matrx Elan Occipital Pad 9" with Infection Infection Control	\$366
☐ EOP12-IC	Matrx Elan Occipital Pad 12"with Infection Control	\$366
☐ E4PTSM-IC	Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control	\$473
☐ E4POINT-IC	Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control	\$473

1. Not compatible with sling back (ASBA-2)

2. If mounting on PB Elite must order removable hardware

HEADREST MOUNTING HARDWARE OPTIONS

☐ MEHWRE	Elan Multi-Axis Removable HDW	\$285
---------------------------------	-------------------------------------	-------

1. Only compatible w/Standard Rehab, Contour & High Back Models. All other head rests REQUIRE an upgrade to Removable Hardware.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS

AADL SEATING BENEFIT

☐ PBSSL	Matrx Small Lateral Fixed-Left(3x5.5")	\$113
☐ PBSSR	Matrx Small Lateral Fixed-Right(3x5.5")	\$113
☐ PBLSL	Matrx Standard Lateral Fixed-Left(7.5x5")	\$116
☐ PBLSR	Matrx Standard Lateral Fixed-Right (7.5x5")	\$116
☐ PBOLL	Matrx Offset Lateral Fixed Left	\$125
☐ PBOLR	Matrx Offset Lateral Fixed Right	\$125

MATRIX SWING AWAY LATERALS - FORELITE AND ELITE TR BACKS

AADL SEATING BENEFIT

Matrx Lateral Trunk Support Left⁽¹⁾

☐ EPADXL	Size XS (3.25x4.25").....	\$115
☐ EPADSL	Size S (3.5x5.25")	\$115
☐ EPADML	Size M (3.75x5.75")	\$115
☐ EPADL L	Size L (4.25x6.5")	\$115

Matrx Lateral Trunk Support Right⁽¹⁾

☐ EPADXR	Size XS (3.25x4.25").....	\$115
☐ EPADSR	Size S (3.5x5.25")	\$115
☐ EPADMR	Size M (3.75x5.75")	\$115
☐ EPADLR	Size L (4.25x6.5")	\$115

1. Must choose swing away hardware below

MATRIX SWING AWAY HARDWARE - FORELITE AND ELITE TR BACKS

AADL SEATING BENEFIT

☐ SAEHBL	Matrx Elite Swing Away Hardware Left	\$250
☐ SAHEBR	Matrx Elite Swing Away Hardware Right	\$250
☐ OSAHL	Matrx Offset Elite Swing Away Hardware Left	\$250
☐ OSAHR	Matrx Offset Elite Swing Away HardwareRight	\$250

MAXX LATERALS TRUNK SUPPORTS⁽¹⁾

AADL SEATING BENEFIT

Left Fixed Mount Lateral Trunk Support⁽¹⁾

Select Left Trunk Support Size

☐ LATXL	Size XS (3.25" x 4.25")	\$100
☐ LATSL	Size S (3.5" x 5.25")	\$100
☐ LATML	Size M (3.75" x 5.75")	\$100
☐ LATLL	Size L (4.25" x 6.5")	\$100

Right Fixed Mount Lateral Trunk Support⁽¹⁾

Select Right Trunk Support Size

☐ LATXR	Size XS (3.25" x 4.25")	\$100
☐ LATSR	Size S (3.5" x 5.25")	\$100
☐ LATMR	Size M (3.75" x 5.75")	\$100
☐ LATLR	Size L (4.25" x 6.5")	\$100

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs.

MAXX STYLE SWING-AWAY TRUNK SUPPORT HARDWARE⁽¹⁾

AADL SEATING BENEFIT

Left MAXX Style Swing Away HDW

Select Left Trunk Support Size

☐ LATXL	Size XS (3.25" x 4.25")	\$100
☐ LATSL	Size S (3.5" x 5.25")	\$100
☐ LATML	Size M (3.75" x 5.75")	\$100
☐ LATLL	Size L (4.25" x 6.5")	\$100

Right MAXX Style Swing Away HDW⁽¹⁾

Select Right Trunk Support Size

☐ LATXR	Size XS (3.25" x 4.25")	\$100
☐ LATSR	Size S (3.5" x 5.25")	\$100
☐ LATMR	Size M (3.75" x 5.75")	\$100
☐ LATLR	Size L (4.25" x 6.5")	\$100

☐ ATLL Additional Telescoping Link - Left⁽²⁾ \$75

☐ ATLR Additional Telescoping Link - Right⁽²⁾ \$75

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs
2. Up to 2.75" medial offset

FIXED HIP/KNEE SUPPORTS

AADL SEATING BENEFIT

Lateral Hip Support with Fixed Mounting Hardware Left⁽¹⁾

Select Pad Size and mounting position at Hip:

☐ HSPXLH	Size XS (5.5" W x 4" H) mounted at Hip	\$120
☐ HSPSMLH	Size S/M (7"W x 4"H) mounted at Hip	\$120
☐ HSPSLH	Size S (7"W x 6"H) mounted at Hip	\$120
☐ HSPMLH	Size M (9" W x 4"H) mounted at Hip	\$120
☐ HSPLLH	Size L (13" W x 4"H) mounted at Hip	\$120

Select Pad Size and mounting position at Knee:

☐ HSPXLK	Size XS (5.5"W x 4" H) mounted at Knee	\$120
☐ HSPSMLK	Size S/M (7"W x 4" H) mounted at Knee	\$120
☐ HSPSLK	Size S (7"W x 6" H) mounted at Knee	\$120
☐ HSPMLK	Size M (9"W x 4"H) mounted at Knee	\$120
☐ HSPLLK	Size L (13"W x 4"H) mounted at Knee	\$120

FIXED HIP/KNEE SUPPORTS- CONTINUED

AADL SEATING BENEFIT

Lateral Hip Support with Fixed Mounting Hardware Right⁽¹⁾

Select Pad Size and mounting position:

☐ HSPXRH	Size XS (5.5"W x 4"H) mounted at Hip	\$120
☐ HSPSMRH	Size S/M (7"W x 4"H) mounted at Hip	\$120
☐ HSPSRH	Size S (7"W x 6"H) mounted at Hip	\$120
☐ HSPMRH	Size M (9"W x 4"H) mounted at Hip	\$120
☐ HSPLRH	Size L (13"W x 4"H) mounted at Hip	\$120

Select Pad Size and mounting position at Knee:

☐ HSPXRK	Size XS(5.5"W x 4"H) mounted at Knee	\$120
☐ HSPSMRK	Size S/M (7"W x 4"H) mounted at Knee	\$120
☐ HSPSRK	Size S (7"W x 6"H) mounted at Knee	\$120
☐ HSPMRK	Size M (9"W x 4"H) mounted at Knee	\$120
☐ HSPLRK	Size L (13"W x 4"H) mounted at Knee	\$120

1. Must choose swing away hardware below if required.

REMOVABLE HIP/KNEE SUPPORT HARDWARE

AADL SEATING BENEFIT

Lift Off Removable Hardware for Hip Support

☐ LOHL	Left	\$195
☐ LOHR	Right	\$195

MAXX Style Quick Release, Removable, Multi-Axis Mounting

Hardware for Hip Supports

Select mounting position:

☐ MHSHLH	Left mounted at Hip	\$215
☐ MHSHLK	Left mounted at Knee	\$215
☐ MHSHRH	Right mounted at Hip	\$215
☐ MHSHRK	Right mounted at Knee	\$215

Swing Away Removable Hardware for Hip Supports

Select mounting position:

☐ SHSHLH	Left mounted at Hip	\$215
☐ SHSHLK	Left mounted at Knee	\$215
☐ SHSHRH	Right mounted at Hip	\$215
☐ SHSHRK	Right mounted at Knee	\$215

Extended Maxx Style Quick Release, Removable, Multi-Axis

Mounting Hardware for Hip Supports⁽¹⁾

Select mounting position:

☐ EMQHLH	Left mounted at Hip	\$215
☐ EMQHLK	Left mounted at Knee	\$215
☐ EMQHRH	Right mounted at Hip	\$215
☐ EMQHRK	Right mounted at Knee	\$215

1. Recommended when finished seat cushion thickness will be greater than 3.5" high.

ELBOW BLOCKS

AADL SEATING BENEFIT

Left Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

- ☐ EBSL Size XS (3.25" x 4.25") \$195
☐ EBSL Size S (3.5" x 5.25") \$195
☐ EBSL Size M (3.75" x 5.75") \$195
☐ EBSL Size L (4.25" x 6.5") \$195

- ☐ V3 Left Verticle Extension3"
☐ V5 Left Verticle Extension 5"
☐ H5 Left Horizontal Extension 5"
☐ H7 Left Horizontal Extension 7"

Right Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

- ☐ EBSL Size XS (3.25" x 4.25") \$195
☐ EBSL Size S (3.5" x 5.25") \$195
☐ EBSL Size M (3.75" x 5.75") \$195
☐ EBSL Size L (4.25" x 6.5") \$195

- ☐ V3 Right Verticle Extension3"
☐ V5 Right Verticle Extension 5"
☐ H5 Right Horizontal Extension 5"
☐ H7 Right Horizontal Extension 7"

1. Prevents arm from slipping.

ADDITIONAL ARM ACCESSORIES

- ☐ UCH Universal Cup Holder⁽¹⁾ CLIENT COST \$50
☐ MACESL Multi Axis Upper Extremity Support⁽²⁾ (each)
Left Side \$225
☐ MACESR Multi Axis Upper Extremity Support⁽²⁾ (each)
Right Side \$225
☐ OAS2L 1" Offset Armrest Spacer Left(each) \$125
☐ OAS2R 1" Offset Armrest Spacer Right(each) \$125
☐ OAS1L 1/2" Offset Armrest Spacer Left(each) \$125
☐ OAS1R 1/2" Offset Armrest Spacer Right(each) \$125
☐ MPAPA Multi Purpose Arm Pad
Adapter Block⁽³⁾ (each) Quantity () \$35

1. Not available on the joystick side.
2. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
3. Not available with Standard Arm Pads or Waterfall Arm Pads.

ARM ADDUCTOR PAD⁽¹⁾

- ☐ AADSL Left Small (6" x 2.25") \$135
☐ AADSR Right Small (6" x 2.25") \$135
☐ AADLL Left Large (8.5" x 2.25") \$135
☐ AADLR Right Large (8.5" x 2.25") \$135

1. Only available with Modular or Flat Multi Position Arm Pads.

O2 HOLDERS⁽¹⁾

- ☐ NAO2H Non-Articulating O² Gas Holder \$395
☐ AOHG Articulating O² Gas Holder
Cylinder Style⁽²⁾⁽³⁾ \$1,250
☐ AOHL Articulating Liquid O² Holder
(6.5" x 4.5")⁽²⁾⁽³⁾ \$1,250

1. The user weight limit is 325 lbs. The maximum weight capacity of the O2 holder is 20 lbs.
2. Required for use on Recline Systems
3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems.

ADDITIONAL OPTIONS & ACCESSORIES

CLIENT COST

- ☐ THS Transfer Handles 4" - 6.5" \$350
☐ THL Transfer Handles 9" - 11.5" \$350
☐ FFLTL Flip and Fold Away 1/2 Lap Tray
Left Mounting(14"W x 10"D)⁽¹⁾⁽²⁾ \$795
☐ FFLTR Flip and Fold Away 1/2 Lap Tray
Right Mounting(14"W x 10"D)⁽¹⁾⁽²⁾ \$795

1. Not available on joystick side and only available with a flat arm pad or modular arm pad.
2. The maximum weight capacity of the tray is 10 lbs.

SUBMIT

PRINT FORM

SAVE

CLEAR FORM