

PRODUCTS LISTED BELOW ARE NOW EXCLUSIVELY AVAILABLE FROM MOTION CONCEPTS USA

Matrx[®] Elite Back Options



1. SuperSoft Foam is now standard on all Elite Full Backs
Optional upgrade for all other Elite Backs & Elite Deep Backs
 - Laminate of 1" soft over 1" firm HR foam for extra immersion and comfort
2. New 14" back heights - now available in Elite backs 12-18" wide
HCPCS Code: E2620

Matrx[®] Elite Heavy-duty Back



- 500 lb. weight capacity
 - 21"-24" wide x 16", 20" high
 - Extra Wide Heavy-duty sizes available up to 30" wide
- HCPCS Code: Pending

Matrx[®] MX1 SuperLite Back



- Lightest weight ever - 1.5lbs and up including hardware!
 - Minimalist hardware when no width/angle adjustment required
 - New style cover with attractive leatherette panels
- HCPCS Code: Pending

Matrx[®] Elite Back Headrest Interface (EHRA)*



- A no charge option when headrest is ordered together with the back
- Elite backs now come standard with pre-drilled holes for headrest attachment (Exceptions include 15"-18" wide x 10", 12" high)

* Also available when ordered on Invacare[®] Solara[®] 3G Wheelchair with Stealth headrest

Matrx[®] Chest Harness Interface



- Now available for Elite and PB backs
- Allows easy attachment and positioning of various chest harnesses with no drilling of back shell

Matrx[®] Genera Back Hardware*



- New single point hardware for easiest installation and adjustment ever
 - The best choice for support, adjustability and comfort at a great price
 - 16", 18", 20" side x 16", 20" high
- HCPCS Code: E2611 General Use

* Also available through Invacare