

... ANOTHER SOLUTION

IN THESE ARMS ■

Multi-Axis Upper Extremity Arm Support


- ✓ Single locking bolt or ratchet handle allows infinite adjustment of angle and rotation of armpad or arm trough
- ✓ Ability to capture complex positions of client's arm/hand
- ✓ Universal 7/8" clamp allows it to fit on any wheelchair armrest, manual or power
- ✓ Low profile bracket easily interfaces with the Motion Concept's Gel Armpad, Arm Trough or Multi-Position Flat Armpad



Multi-Axis Hardware



Upper Extremity Arm Support



OutBack Arm Flipped Out – 30° Angle



OutBack Arm —

Available for the Motion Concept's dual post tilt armrest or the recline armrest – the entire armrest can be flipped out at a 30° angle to eliminate interference of the arm when using the Multi-Axis Upper Extremity support or when there is a deep back or laterals that may interfere with the flip back of a traditional armrest.

www.motionconcepts.com

usa

700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
Toll Free Tel: 1.888.433.6818
Toll Free Fax: 1.888.433.6834

canada

84 Citation Dr.,
Concord, ON L4K 3C1
Toll Free Tel: 1.866.748.7943
Tel: 905.695.0134
Fax: 905.695.0138

Motion Concepts

matrix
SEATING SERIES