

REQUIRED INFORMATION

Quote ☐ Order ☐

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



Customer Service: 1.866.748.7943 | Fax: 1.905.695.0583 | orderscdn@motionconcepts.com

ROVI

X3

BASE STATISTICS

ADP CODES	WP0001539
SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH	23.25"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"
CASTER WHEEL SIZE (front/rear)	6"
TURNING RADIUS	20"
DRIVETRAIN	4-pole EAD Motors

⁽¹⁾ All X3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

BATTERIES	2 x M34
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs (M34) / 17.5 kg
RANGE ⁽²⁾	17.7 miles / 28.5 km/h
WEIGHT CAPACITY	300 lbs / 136 kg
BASE WEIGHT (w/o batteries)	157 lbs / 71 kg
GROUND CLEARANCE	3"
MAXIMUM SPEED	5.0 mph



Shown with Motion Concepts Single Post
Back for Contour Backrest. See page 7



Maxx Rehab Seat on ROVI VR2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

X3 BASE OPTIONS

☒ 888WNLEX3-VR2 ROVI X3 5.0 MPH MOTORS⁽¹⁾\$6,800

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt. Weight capacity is 300 lbs.

FENDER COLOUR

☐ 800914-00200-03	Blizzard White		N/C
☐ 800914-00200-50	Matte Black		N/C
☐ 800914-00200-69	Solar Yellow		N/C
☐ 800914-00200-36	Motion Green		N/C
☐ 800914-00200-23	Original Orange		N/C
☐ 800914-00200-23	Bright Blue		N/C
☐ 800914-00200-21	Red Delicious		N/C
☐ 800914-00200-2C	Party Pink		N/C
☐ 800914-00200-80	Purple Reign		N/C

WHEELCHAIR TRANSPORT BRACKETS

☒ XBKTS Wheelchair Transport Brackets⁽¹⁾ STD

1. For unoccupied use only

TIRE OPTIONS

☐ STIRES	Standard Solid Tires	STD
☐ PTIRES	Pneumatic Tires	N/C

BATTERIES

☐ M34 SLD G MK Group Battery x 2 \$650

JOYSTICK SELECTION

☒ D50677 VR2 Non-Expandable Joystick STD

JOYSTICK MOUNT

☐ 800914-001-16-ML	Fixed Joystick Mounting Left	N/C
☐ 800914-001-16-MR	Fixed Joystick Mounting Right	N/C
☐ SAQL	Height Adjustable Swing-Away Quad Link Left	\$344
☐ SAQR	Height Adjustable Swing-Away Quad Link Right	\$344

JOYSTICK KNOB - OPTIONAL

☐ SA80941	'T' Handle for Joystick Knob	\$115
☐ SA72441	Ball Handle for Joystick Knob	\$115
☐ PC101B-20	BodyPoint Joystick Handle 3" U-Shaped	\$120
☐ PC102B-20	BodyPoint Joystick Handle 4" U-Shaped	\$120
☐ PC107B-20	BodyPoint Joystick Handle Domed Rubber	\$120

ELECTRONIC PROGRAMMER

☐ D49510	VR2 Programmer (PPI)	\$603
☐ D50996	Diagnostic Test Tool (DTT) with VR2 & RNet Cables	\$1,876
☐ SPH	Smart Phone Holder	\$50
☐ MUSB	Motion Dual USB Charger	N/C

SEAT-TO-FLOOR HEIGHT⁽¹⁾

☐ STF17	17" Seat-To-Floor	N/C
☐ STF18	18" Seat-To-Floor	N/C
☐ STF19	19" Seat-To-Floor	\$250
☐ STF20	20" Seat-To-Floor	\$250

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"

NOTES

Maxx Rehab Seat on ROVI VR2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

REHAB SEATING SYSTEM

- ☒ D0195 Maxx Adjustable Rehab Seat \$1,200
- ☐ PRETILT 5° Fixed Static Tilt Angle (Pre-Tilt)⁽¹⁾ \$250
- ☐ ADJUST Adjustable Static-Tilt⁽¹⁾ \$250
(in 1.25° increments from 0° to 5°)
Specify set at _____°

1. The front seat to floor height will be approximately 2" greater than the STF selected above. It will also be dependant on the seat depth selected.

ADULT SEAT WIDTH SETTINGS

- ☐ SW116 16" Wide (Adjustable from 16" to 20") N/C
- ☐ SW117 17" Wide (Adjustable from 16" to 20") N/C
- ☐ SW118 18" Wide (Adjustable from 16" to 20") N/C
- ☐ SW119 19" Wide (Adjustable from 16" to 20") N/C
- ☐ SW120 20" Wide (Adjustable from 16" to 20") N/C
- ☐ W219 19" Wide (Adjustable from 19" to 22") \$460
- ☐ W220 20" Wide (Adjustable from 19" to 22") \$460
- ☐ W221 21" Wide (Adjustable from 19" to 22") \$460
- ☐ W222 22" Wide (Adjustable from 19" to 22") \$460

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

- ☐ SD116 16" Deep (Adjustable from 16" to 21") N/C
- ☐ SD117 17" Deep (Adjustable from 16" to 21") N/C
- ☐ SD118 18" Deep (Adjustable from 16" to 21") N/C
- ☐ SD119 19" Deep (Adjustable from 16" to 21") N/C
- ☐ SD120 20" Deep (Adjustable from 16" to 21") N/C
- ☐ D219 19" Deep (Adjustable from 19" to 23") \$460
- ☐ D220 20" Deep (Adjustable from 19" to 23") \$460
- ☐ D221 21" Deep (Adjustable from 19" to 23") \$460
- ☐ D222 22" Deep (Adjustable from 19" to 23") \$460
- ☐ D223 23" Deep (Adjustable from 19" to 23") \$460

MATRIX SEATING OPTIONS

- ☐ Matrix Cushions available on page 8
- ☐ Matrix Backs available on page 8
- ☐ Matrix Headrests available on page 9



CHEST STRAPS

- ☐ MFUS BP Monoflex, Underarm RLS S(3" x 15.5") . \$195
- ☐ MFUM BP Monoflex, Underarm RLS M(4" x 18.25") \$195
- ☐ MFUL BP Monoflex, Underarm RLS L(4.5" x 21.5") . \$195
- ☐ MFCS BP Monoflex, Center RLS S(2.5" x 15.5") . . \$195
- ☐ MFCM BP Monoflex, Center RLS M(3.25" x 18.25") \$195
- ☐ MFCL BP Monoflex, Center RLS L(3.5" x 21.5") . . \$195

LAP BELTS

- ☐ BELT48 Push Button Style Lap Belt, Length 48" STD
- ☐ BELT60 Push Button Style Lap Belt, Length 60" N/C
- ☐ BELT71 Push Button Style Lap Belt, Length 71" N/C
- ☐ PLB2 Bodypoint Padded Lap Belt 2 Point⁽¹⁾ \$250
- ☐ PLB4 Bodypoint Padded Lap Belt 4 Point⁽¹⁾ \$315
1. 67" length only option

MOTION CONCEPTS HEADREST OPTIONS

- ☐ MCSH Motion Concepts Standard N/C
- ☐ MCOH Motion Concepts Auto Style \$235
- ☐ FM Fixed Mounting (non-removable) N/C
- ☐ MAHU Motion Concepts Multi-Axis Rem.HDW \$195
- ☐ OHR Omit: Headrest N/C

MATRIX HEADRESTS OPTIONS

- ☐ Matrix Elan Headrests and Hardware
available on page 9



†STANDARD BACK CANES⁽¹⁾

- ☐ SBC18 18" Straight Back Canes STD
- ☐ SBC20 20" Straight Back Canes STD
- ☐ SBC22 22" Straight Back Canes STD
- ☐ SBC24 24" Straight Back Canes STD
- ☐ TBC18 18" 10° Mid-Angle Canes for Sling Back . . . N/C
- ☐ TBC20 20" 10° Mid-Angle Canes for Sling Back . . . N/C
- ☐ TBC22 22" 10° Mid-Angle Canes for Sling Back . . . N/C
- ☐ TBC24 24" 10° Mid-Angle Canes for Sling Back . . . N/C
- ☐ ABC20 20" Angled (5°) Back Canes N/C
- ☐ ABC22 22" Angled (5°) Back Canes N/C
- ☐ ABC24 24" Angled (5°) Back Canes N/C
- ☐ SBCF18 18" Straight Fold Down Back Canes \$250
- ☐ SBCF20 20" Straight Fold Down Back Canes \$250
- ☐ SBCF22 22" Straight Fold Down Back Canes \$250
- ☐ SBCF24 24" Straight Fold Down Back Canes \$250

1. Canes will be set at 90° angle and are adjustable from 75° to 116°

Maxx Rehab Seat on ROVI VR2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾

☐18" ☐19" ☐20" ☐21" ☐22" ☐23" ☐24" ☐25"
FBT18 FBT19 FBT20 FBT21 FBT22 FBT23 FBT24 FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK OPTIONS

- ☐SRB Standard Rehab Back & CushionN/C
☐SLING Sling Back Fabric \$125
☐ORBC Omit Rigid Back & CushionN/C

MATRIX BACKREST OPTIONS

- ☐ Matrix Backs available on page 9



ARMRESTS

- ☐DP1 Adult Dual Post Adj. 9.75"-12.75"STD
☐DP2 Adult Dual Post Adj. 12.75"-15.75"STD
☐DP3 Adult Dual Post Adj. 15.75"-18.75"STD
Select: ☐FP Flat PouchSTD
☐GB Glove BoxesN/C
☐OAL ADD: OutBack Arm Option Left⁽¹⁾ each \$195
☐OAR ADD: OutBack Arm Option Right⁽¹⁾ each \$195
☐OCTA1 Outback Cantilever Tilt Armrest 9"-12"⁽³⁾ . . . \$488
☐OCTA2 Outback Cantilever Tilt Armrest 12"-15"⁽³⁾ . . \$488
☐RMCAN Ultra Rail Seat Mnt Flip Back Cantilever
MAXX Tilt Arm (9-12.5")⁽²⁾ \$488
☐RMCAN2 Ultra Rail Seat Mnt Flip Back Cantilever
MAXX Tilt Arm (11-14")⁽²⁾ \$488
☐DPM ADD: Dual Post Module to RMCAN (pair) . . . \$250

1. Only available with Dual Post Armrests.

2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.

3. Not available with the Universal Cup Holder or Flip and Fold Lap Tray.

ARMREST TUBE LENGTH - PLEASE SPECIFY

- ☐ATFL Left FullSTD
☐ATDL Left DeskSTD

☐ATFR Right FullSTD
☐ATDR Right DeskSTD

ARMPADS

Standard

- ☐SAPFL Left Full (2.25" x 14")STD
☐SAPDL Left Desk (2.25" x 10")STD
☐SAPFR Right Full (2.25" x 14")STD
☐SAPDR Right Desk (2.25" x 10")STD

ARMPADS

Waterfall Arm Pads⁽¹⁾

- ☐WAPFL Left Full (13") N/C
☐WAPDL Left Desk (9") N/C
☐WAPFR Right Full (13") N/C
☐WAPDR Right Desk (9") N/C

Modular Arm Pad - Left⁽²⁾

- ☐MAPTFL Left Tray Full (14") \$95
☐MAPTDL Left Tray Desk (10") \$95

Select Pad Insert below:

- ☐ISSFL Iskin - Integrated soft skin foam
☐SCVFL Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

- ☐MAPTFR Right Tray Full (14") \$95
☐MAPTDR Right Tray Desk (10") \$95

Select Pad Insert below:

- ☐ISSFR Iskin - Integrated soft skin foam
☐SCVFR Startex Covered Visco foam

Flat Multi Position Arm Pads

- ☐MPAFL Left Full (14") \$150
☐MPADL Left Desk (10") \$150
☐MPAFR Right Full (14") \$150
☐MPADR Right Desk (10") \$150

Gel Pads

- ☐GAPNL Narrow Left (2" x 12") \$88
☐GAPNR Narrow Right (2" x 12") \$88
☐GAPWL Wide Left (3.5" x 12") \$96
☐GAPWR Wide Right (3.5" x 12") \$96
☐GAPLL Long Left (3.5" x 14") \$103
☐GAPLR Long Right (3.5" x 14") \$103

Arm Troughs

- ☐EATL Ergonomic Arm Trough - Left \$121
☐EATR Ergonomic Arm Trough - Right \$121
☐FHPL Flat Hand Pad - Left⁽³⁾ \$91
☐FHPR Flat Hand Pad - Right⁽³⁾ \$91

Arm Pad Storage/Cup Holder⁽⁴⁾

- ☐APSL Left \$295
☐APSR Right \$295
☐OMIT Omit Armpads(7/8" diameter tube) N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.

2. Requires interchangeable pad insert.

3. Must order with ergonomic arm trough.

4. Not available on the joystick side.

Maxx Rehab Seat on ROVI VR2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

STANDARD FRONT RIGGING


- ☐ FCMP Fixed Center Mount Platform⁽¹⁾ N/C
- ☐ IFP Comes with a rubber coated footplate
(11.5" W x 10") N/C

Select seat pan to footplate:

- ☐ SPTF1 9" - 13" N/C
- ☐ SPTF14 13" -19" N/C

Set at:

- ☐ SA70 Set at 70° N/C
- ☐ SA90 Set at 90° N/C
- ☐ SA97 Set at 97° N/C



1. 3" ground clearance required.

MANUAL CENTER MOUNT FRONT RIGGING

- ☐ MSFCMFP Maxx Style Fixed Center Mount
Foot Platform⁽¹⁾⁽²⁾ \$395
- With Quick Release Height Adjustment


- ☐ IFP Comes with a rubber coated footplate
(11.5" W x 10") N/C
- ☐ ICP Comes standard with individual calf pads. N/C

Select seat pan to footplate:

- ☐ SPTF1 9" - 13" N/C
- ☐ SPTF14 13" -19" N/C

Set at:

- ☐ SA70 Set at 70° N/C
- ☐ SA90 Set at 90° N/C
- ☐ SA97 Set at 97° N/C



1. Depth and height adjustable.
2. 3" ground clearance required.

INVACARE CENTER MOUNT

- ☐ IAKACM Invacare Center⁽¹⁾ "CTMT" \$595

Select seat pan to footplate:

- ☐ SPTF3 9.25" - 12" N/C
- ☐ SPTF14 12.25" -15" N/C

Select foot plates/footboard

- ☐ SM Individual Foot Plates (5.25"W x 7.5"D) N/C
- ☐ MDM Individual Foot Plates (5.5"W x 9.5"D) N/C
- ☐ LRG Individual Foot Plates (5.9"W x 11.5"D) N/C
- ☐ FB Footboard. STD

Optional:

- ☐ CALF Invacare Calf Pads for use with CTMT⁽²⁾ \$268

1. 3" ground clearance required.
2. Not available with the SPTF3 of 9.25"-12"

MAXX STYLE ELEVATING LEGRESTS

- ☐ MMPPSA MAXX Swing Away
Manual Elevating - Straight⁽¹⁾⁽²⁾ \$680

Select Receivers:

- ☐ RS Straight Receivers N/C
- ☐ RF Flared 1" Outward Receivers N/C


Select Seat Pan to Footplate:

- ☐ SPTF5 Length: 10"-14" N/C
- ☐ SPTF2 Length: 13"-17" N/C
- ☐ SPTF7 Length: 16"-20" N/C

- ☐ FCPMS Flat Calf Pad Option (pair) \$275
- ☐ CCPMS Curved Calf Pad Option (inside with 6")⁽³⁾ (pair) \$275

1. Must choose footplate on page 6
2. Comes standard with moulded i-skin calf pads.
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.


FOOT PLATFORM OPTIONS FOR FIXED/MAXX CMFP



- ☐ SFFP Small Foot Platform⁽¹⁾ (12"w x 8.5"d) \$180
- ☐ MFFP Medium Platform⁽¹⁾ (12"w x 10.75"d) \$180
- ☐ LFFP Large Platform⁽¹⁾ (14.5"w x 13"d) \$250
- ☐ XLFFP Extra Large Platform⁽¹⁾ (17"w x 13"d) \$250
- ☐ XXLFFP Extra Extra Large Platform⁽¹⁾ (21"w x 13"d) \$250

1. May need to adjust leg position to avoid interference with front castors

INDIVIDUAL FOOT PLATES FOR FIXED/MAXX CMFP



- ☐ IFSL Small⁽¹⁾ - (5" x 7.5" deep) - Left \$165
- ☐ IFSR Small⁽¹⁾ - (5" x 7.5" deep) - Right \$165
- ☐ IFML Medium⁽¹⁾ - (5.5" x 9.5" deep) - Left \$165
- ☐ IFMR Medium⁽¹⁾ - (5.5" x 9.5" deep) - Right \$165
- ☐ IFLR Large⁽¹⁾ - (6" x 11.5" deep) - Left \$165
- ☐ IFLR Large⁽¹⁾ - (6" x 11.5" deep) - Right \$165

1. Foot plates may interfere with casters in certain configurations

SWING AWAY LEGREST


- ☐ HD70 Heavy Duty 70° Swing Away Footrest(pair)⁽¹⁾⁽²⁾ . \$495
- Select Receivers:
- ☐ RS Straight Receivers N/C
- ☐ RF Flared 1" Outward Receivers N/C
- Select Seat Pan to Footplate:
- ☐ SPTF11 Length: 6"-10" N/C
- ☐ SPTF8 Length: 11"-15" N/C
- ☐ SPTF2 Length: 13"-17" N/C
- ☐ SPTF7 Length: 16"-20" N/C
- ☐ HDFCP Flat Calf Pad Option (pair) \$275
- ☐ HD70CCP Curved Calf Pad Option⁽³⁾ (inside with 6") (pair) \$275
- ☐ DCS Dynamic Coil Springs (Pair) for Heavy Duty 70° Swingaway Footrests \$190
- ☐ I70SAF Invacare 70° Swingaway Footrests (Pair) \$325
- ☐ I70-CF Comes with composite foot plates N/C
- ☐ 70HDAIH Invacare 70° Non-Tapered HD Swingaway Footrests with AT5543 Adjustable Angle Foot plates, Heel Loops and Impact Guards⁽¹⁾ (Pair) \$457
- ☐ 70TAPAIH Invacare Tapered 70° Swingaway Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards(Pair) \$497
- ☐ ALPT Longer Pivot and Slide Tube for 70° Swingaways (Pair) \$132
- ☐ RLSL Residual Limb Support - Left⁽⁴⁾ (6.75"W x 9.75"L x 2"H) (each) \$495
- ☐ RLSR Residual Limb Support - Right⁽⁴⁾ (6.75"W x 9.75"L x 2"H) (each) \$495

1. Must choose footplate on below.

2. Heavy duty durable construction (4.5lbs ea.).

3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

4. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.

FOOTPLATE SIZE SELECTION



1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

ADJUSTABLE ANGLE FOOTPLATES

- ☐ AFSL Left Small Footplate (4.25"W x 8"D) \$165
- ☐ AFSR Right Small Footplate (4.25"W x 8"D) \$165
- ☐ AFML Left Medium Footplate (5.25"W x 8"D) \$165
- ☐ AFMR Right Medium Footplate (5.25"W x 8"D) \$165
- ☐ AFLL Left Large Footplate (6.25"W x 8"D) \$165
- ☐ AFLR Right Large Footplate (6.25"W x 8"D) \$165
- ☐ AFXL Left Extra Large Footplate (9"W x 11"D) \$165
- ☐ AFXR Right Extra Large Footplate (9"W x 11"D) \$165
- ☐ HLL Left Heel Loop (each) \$20
- ☐ HLR Right Heel Loop (each) \$20

MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES

- ☐ MFSL Left Small Footplate (4.25"W x 8"D) \$225
- ☐ MFSR Right Small Footplate (4.25"W x 8"D) \$225
- ☐ MFML Left Medium Footplate (5.25"W x 8"D) \$225
- ☐ MFMR Right Medium Footplate (5.25"W x 8"D) \$225
- ☐ MFLR Left Large Footplate (6.25"W x 8"D) \$225
- ☐ MFLR Right Large Footplate (6.25"W x 8"D) \$225
- ☐ MFXL Left Extra Large Footplate (9"W x 11"D) \$225
- ☐ MFXR Right Extra Large Footplate (9"W x 11"D) \$225
- ☐ HLL Left Heel Loop (each) \$20
- ☐ HLR Right Heel Loop (each) \$20

FOOTPLATE/CUT OUTS SIZE SELECTION


1. Certain footplate sizes may not fit all seat widths depending on legrests selected

ADJUSTABLE ANGLE FOOTPLATES/CUT OUTS

- ☐ AFCSL Left Small Footplate (7.25"W x 9"D) \$165
- ☐ AFCSR Right Small Footplate (7.25"W x 9"D) \$165
- ☐ AFCML Left Medium Footplate (7.5"W x 11"D) \$165
- ☐ AFCMR Right Medium Footplate (7.5"W x 11"D) \$165
- ☐ AFCLL Left Large Footplate (10"W x 14"D) \$165
- ☐ AFCLR Right Large Footplate (10"W x 14"D) \$165
- ☐ HLL Left Heel Loop (each) \$20
- ☐ HLR Right Heel Loop (each) \$20

MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES/ CUT OUTS

- ☐ MFCSL Left Small Footplate (7.25"W x 9"D) \$225
- ☐ MFCSR Right Small Footplate (7.25"W x 9"D) \$225
- ☐ MFCML Left Medium Footplate (7.5"W x 11"D) \$225
- ☐ MFCMR Right Medium Footplate (7.5"W x 11"D) \$225
- ☐ MFCLL Left Large Footplate (10"W x 14"D) \$225
- ☐ MFCLR Right Large Footplate (10"W x 14"D) \$225
- ☐ HLL Left Heel Loop (each) \$20
- ☐ HLR Right Heel Loop (each) \$20

THORACIC LATERALS AND HIP/KNEE SUPPORTS

- ☐ Matrx Laterals please see page 10
- ☐ Maxx Laterals and Hip/Knee supports please see page 11




Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CUSHION COVER REFERENCE GUIDE

					
					
Cushion	Meshtex	Startex	Infection	Spacetex	Onyx
Libra/Adjust	N/A	STD	\$50	N/C	N/A
PS/PSP	N/A	STD -Reversible	STD -Reversible	N/A	N/A
Vi	STD	\$103	\$103	N/A	\$103

MATRIX SEAT CUSHION OPTIONS
Matrx Libra Adjust Seat Cushion 16" - 20" wide \$985

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ LC1616-ADJ | ☐ LC1618-ADJ | ☐ LC1620-ADJ | ☐ LC1716-ADJ |
| ☐ LC1718-ADJ | ☐ LC1720-ADJ | ☐ LC1816-ADJ | ☐ LC1818-ADJ |
| ☐ LC1820-ADJ | ☐ LC1916-ADJ | ☐ LC1918-ADJ | ☐ LC1920-ADJ |
| ☐ LC2016-ADJ | ☐ LC2018-ADJ | ☐ LC2020-ADJ | |

Matrx Libra Adjust Seat Cushion 22" depth \$1,265

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ LC1622-ADJ | ☐ LC1722-ADJ | ☐ LC1822-ADJ |
| ☐ LC1922-ADJ | ☐ LC2022-ADJ | |

Matrx Libra Adjust Seat Cushion 21" - 22" wide \$1,265

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ LC2118-ADJ | ☐ LC2120-ADJ | ☐ LC2122-ADJ |
| ☐ LC2218-ADJ | ☐ LC2220-ADJ | ☐ LC2222-ADJ |

Matrx Libra Seat Cushion 16" - 20" wide \$745

- | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ LC1616 | ☐ LC1618 | ☐ LC1620 | ☐ LC1716 | ☐ LC1718 |
| ☐ LC1720 | ☐ LC1816 | ☐ LC1818 | ☐ LC1820 | ☐ LC1916 |
| ☐ LC1918 | ☐ LC1920 | ☐ LC2016 | ☐ LC2018 | ☐ LC2020 |

Matrx Libra Seat Cushion 22" deep \$1,025

- | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ LC1622 | ☐ LC1722 | ☐ LC1822 | ☐ LC1922 | ☐ LC2022 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|

Matrx Libra HD Seat Cushion 21" - 22" wide \$1,025

- | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ LC2118 | ☐ LC2120 | ☐ LC2122 | ☐ LC2218 | ☐ LC2220 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|

Matrx PS Seat Cushion 16" - 20" wide. \$569

- | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ PS1616 | ☐ PS1618 | ☐ PS1620 | ☐ PS1716 | ☐ PS1718 |
| ☐ PS1720 | ☐ PS1816 | ☐ PS1818 | ☐ PS1820 | ☐ PS2016 |
| ☐ PS2018 | ☐ PS2020 | | | |

Matrx PS HD Seat Cushion 22" - 22" wide \$975

- | | | |
|---------------------------------|---------------------------------|---------------------------------|
| ☐ PS2218 | ☐ PS2220 | ☐ PS2222 |
|---------------------------------|---------------------------------|---------------------------------|

Matrx PSVF Seat Cushion 16" - 20" wide \$589

- | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ PSVF1616 | ☐ PSVF1618 | ☐ PSVF1620 |
| ☐ PSVF1816 | ☐ PSVF1818 | ☐ PSVF1820 |
| ☐ PSVF2016 | ☐ PSVF2018 | ☐ PSVF2020 |

MATRIX SEAT CUSHION OPTIONS
Matrx PSP Seat Cushion 16" - 20" wide \$645

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PSP1616 | ☐ PSP1618 | ☐ PSP1620 | ☐ PSP1716 |
| ☐ PSP1718 | ☐ PSP1720 | ☐ PSP1816 | ☐ PSP1818 |
| ☐ PSP1820 | ☐ PSP2016 | ☐ PSP2018 | ☐ PSP2020 |

Matrx Vi Seat Cushion 14" - 20" wide \$624

- | | | | |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| ☐ MA1415-VI | ☐ MA1417-VI | ☐ MA1616-VI | ☐ MA1618-VI |
| ☐ MA1620-VI | ☐ MA1716-VI | ☐ MA1718-VI | ☐ MA1720-VI |
| ☐ MA1816-VI | ☐ MA1818-VI | ☐ MA1820-VI | ☐ MA2016-VI |
| ☐ MA2018-VI | ☐ MA2020-VI | | |

Matrx Vi Seat Cushion 21" - 22" wide \$681

- | | | | |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| ☐ MA2016-VI | ☐ MA2018-VI | ☐ MA2020-VI | ☐ MA2218-VI |
| ☐ MA2220-VI | | | |


Matrx Vi HD Seat Cushion 22" wide \$840

- | |
|------------------------------------|
| ☐ MA2222-VI |
|------------------------------------|

FABRIC OPTIONS FOR CUSHIONS

- | | | |
|--------------------------------|--|-------|
| ☐ SFLF | Spacetex Fabric Option for
Libra Cushion | N/C |
| ☐ ICFL | Infection Control Fabric Upgrade for
Libra Cushion | \$50 |
| ☐ OCUMA | Onyx Fabric Upgrade for
Matrx Vi Cushion. | \$100 |
| ☐ SFVS | Smooth side Startex Fabric Upgrade for
Matrx Vi Cushion | \$50 |
| ☐ SFVF | Fabric side Startex Fabric Upgrade for
Matrx Vi Cushion | \$50 |

BACK COVER REFERENCE GUIDE

			
Back	Meshtex	Infection	Onyx
High	STD	N/A	\$103
Elite/ Elite HD/Elite Deep	STD	\$103	N/A
Elite TR/ Elite TR HD	STD	\$103	N/A
PB/PB HD/PB Deep/PB HD Deep	STD	\$103	N/A

OTHER UPGRADABLE BACK OPTIONS

- | | | |
|--------------------------------|---|-------|
| ☐ EBSS | Super Soft Foam for Elite Back Upgrade . . . | \$103 |
| ☐ OCUHG | Onyx Fabric Upgrade for High Backs. | \$103 |
| ☐ WSCU | Infection Control Cover for Elite/PB Backs . . . | \$103 |
| | <i>purchased with back</i> | |
| ☐ WSC | Infection Control Cover for Elite/PB Backs . . . | \$150 |
| | <i>purchased with without back</i> | |
| ☐ CHIEB | Chest Harness Interface Kit for Elite Back ⁽¹⁾⁽²⁾ . | \$103 |
| ☐ CHIEB | Chest Harness Interface Kit for PB Back ⁽¹⁾⁽²⁾ . . . | \$52 |

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harnesses.

UPGRADABLE BACK OPTIONS

High Back 16" - 22" wide \$595

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ HB1616 | ☐ HB1618 | ☐ HB1620 | ☐ HB1716 |
| ☐ HB1718 | ☐ HB1720 | ☐ HB1816 | ☐ HB1818 |
| ☐ HB1820 | ☐ HB1916 | ☐ HB1918 | ☐ HB1920 |
| ☐ HB2016 | ☐ HB2018 | ☐ HB2020 | ☐ HB2116 |
| ☐ HB2118 | ☐ HB2120 | ☐ HB2216 | ☐ HB2218 |
| ☐ HB2220 | | | |

Matrx Elite Back 16" - 20" wide (3" Contour) \$766

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBE1614 | ☐ PBE1616 | ☐ PBE1618 | ☐ PBE1620 |
| ☐ PBE1714 | ☐ PBE1716 | ☐ PBE1718 | ☐ PBE1720 |
| ☐ PBE1814 | ☐ PBE1816 | ☐ PBE1818 | ☐ PBE1820 |
| ☐ PBE1916 | ☐ PBE1918 | ☐ PBE1920 | ☐ PBE2014 |
| ☐ PBE2016 | ☐ PBE2018 | ☐ PBE2020 | |

Matrx Elite HD Back 21" - 22" wide (3" Contour) . . . \$1,095

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE2116-HD | ☐ PBE2120-HD | ☐ PBE2216-HD |
| ☐ PBE2220-HD | | |

Matrx Elite Deep Back 16" - 20" wide (3" Contour) . . . \$918

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ EDB1614 | ☐ EDB1616 | ☐ EDB1618 | ☐ EDB1620 |
| ☐ EDB1714 | ☐ EDB1716 | ☐ EDB1718 | ☐ EDB1720 |
| ☐ EDB1814 | ☐ EDB1816 | ☐ EDB1818 | ☐ EDB1820 |
| ☐ EDB1914 | ☐ EDB1916 | ☐ EDB1918 | ☐ EDB1920 |
| ☐ EDB2014 | ☐ EDB2016 | ☐ EDB2018 | ☐ EDB2020 |

Matrx Elite Deep HD Back 21" - 22" wide (6" Contour). \$1,095

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ EDB2116-HD | ☐ EDB2120-HD | ☐ EDB2216-HD |
| ☐ EDB2220-HD | | |

NEW Matrx Elite E2 Back 16" - 20" wide/14"/16"H(3" Contour) . \$965

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2S1614 | ☐ E2S1616 | ☐ E2S1714 | ☐ E2S1716 |
| ☐ E2S1814 | ☐ E2S1816 | ☐ E2S1914 | ☐ E2S1916 |
| ☐ E2S2014 | ☐ E2S2016 | | |

NEW Matrx Elite E2 Back 16" - 20" wide/18"/20"H(3" Contour) \$995

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2S1618 | ☐ E2S1620 | ☐ E2S1718 | ☐ E2S1720 |
| ☐ E2S1818 | ☐ E2S1820 | ☐ E2S1918 | ☐ E2S1920 |
| ☐ E2S2018 | ☐ E2S2020 | | |

NEW Matrx Elite E2 Deep Back 16" - 20" wide/14"/16"H (6" Contour) \$1,125

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2D1614 | ☐ E2D1616 | ☐ E2D1714 | ☐ E2D1716 |
| ☐ E2D1814 | ☐ E2D1816 | ☐ E2D1914 | ☐ E2D1916 |
| ☐ E2D2014 | ☐ E2D2016 | | |

NEW Matrx Elite E2 Deep Back 14" - 20" wide/18"/20"H (6" Contour) \$1,155

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2D1418 | ☐ E2D1420 | ☐ E2D1518 | ☐ E2D1520 |
| ☐ E2D1618 | ☐ E2D1620 | ☐ E2D1718 | ☐ E2D1720 |
| ☐ E2D1818 | ☐ E2D1820 | ☐ E2D1918 | ☐ E2D1920 |
| ☐ E2D2018 | ☐ E2D2020 | | |

UPGRADABLE BACK OPTION

Matrx Elite TR Back 15" - 20" wide (3" Contour). \$866

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE1516-TR | ☐ PBE1518-TR | ☐ PBE1520-TR |
| ☐ PBE1616-TR | ☐ PBE1618-TR | ☐ PBE1620-TR |
| ☐ PBE1716-TR | ☐ PBE1718-TR | ☐ PBE1720-TR |
| ☐ PBE1816-TR | ☐ PBE1818-TR | ☐ PBE1820-TR |
| ☐ PBE1916-TR | ☐ PBE1918-TR | ☐ PBE1920-TR |
| ☐ PBE2016-TR | ☐ PBE2018-TR | ☐ PBE2020-TR |

Matrx Elite TR HD Back 21" - 22" wide (3" Contour). . \$1,151

- | | |
|---------------------------------------|---------------------------------------|
| ☐ PBE2120-TRHD | ☐ PBE2220-TRHD |
|---------------------------------------|---------------------------------------|

Matrx PB Posture Back 14" - 20" wide (5" Contour) \$734

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ MPB1412 | ☐ MPB1612 | ☐ MPB1616 | ☐ MPB1620 |
| ☐ MPB1712 | ☐ MPB1716 | ☐ MPB1720 | ☐ MPB1812 |
| ☐ MPB1816 | ☐ MPB1820 | ☐ MPB2012 | ☐ MPB2016 |
| ☐ MPB2020 | | | |

Matrx PB Posture Back Deep 14" - 20" wide (7" Contour) . \$892

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBD1412 | ☐ PBD1612 | ☐ PBD1616 | ☐ PBD1620 |
| ☐ PBD1712 | ☐ PBD1716 | ☐ PBD1720 | ☐ PBD1812 |
| ☐ PBD1816 | ☐ PBD1820 | ☐ PBD2012 | ☐ PBD2016 |
| ☐ PBD2020 | | | |

Matrx PB Heavy Duty Back 20" - 22" wide (5" Contour) . \$1,051

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ HDB2016 | ☐ HDB2020 | ☐ HDB2216 | ☐ HDB2220 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrx PB Heavy Duty Deep Back 20"-22" wide (5" Contour) \$1,209

- | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ DHDB2016 | ☐ DHDB2020 | ☐ DHDB2216 | ☐ DHDB2220 |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|

Maxx Rehab Seat on ROVI VR2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX ELAN HEADREST PAD⁽¹⁾⁽²⁾

☐ ESP6	Matrx Elan Standard Pad 6"	\$310
☐ ESP10	Matrx Elan Standard Pad 10"	\$310
☐ ESP14	Matrx Elan Standard Pad 14"	\$310
☐ EOP9	Matrx Elan Occipital Pad 9"	\$341
☐ EOP12	Matrx Elan Occipital Pad 12"	\$341
☐ E4POINT	Matrx Elan 4-Point Pad 11"W x 10"H	\$448
☐ ESP6-IC	Matrx Elan Standard Pad 6" with Infection Control	\$335
☐ ESP10-IC	Matrx Elan Standard Pad 10" with Infection Control	\$335
☐ ESP14-IC	Matrx Elan Standard Pad 14" with Infection Control	\$335
☐ EOP9-IC	Matrx Elan Occipital Pad 9" with Infection Control	\$366
☐ EOP12-IC	Matrx Elan Occipital Pad 12" with Infection Control	\$366
☐ E4PTSM-IC	Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control	\$473
☐ E4POINT-IC	Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control	\$473

1. Not compatible with sling back

2. If mounting on PB Elite must order removable hardware

HEADREST MOUNTING HARDWARE OPTIONS

 ☐ LXSHW	Loxx Standard Hardware (15" Vertical Post)	\$275
 ☐ LXMMH	Loxx Mini Hardware (11" Vertical Post)	\$275
☐ MEHWRE	Elan Hardware (15" Vertical Post)	\$285
☐ MEMH	Elan Mini Hardware (11" Vertical Post)	\$285
☐ MEHW-EXT	Elan Extended Hardware (Extra Link)	\$320

1. Only compatible w/Standard Rehab, Contour & High Back Models. All other head rests REQUIRE an upgrade to Removable Hardware.

MATRIX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS

☐ E2FLL	Matrx E2 Fixed Laterals-Left(6.5x5") ⁽¹⁾	\$195
☐ E2FLR	Matrx E2 Fixed Laterals-Right(6.5x5") ⁽¹⁾	\$195
☐ PBSSL	Matrx Small Lateral Fixed-Left(3x5.5")	\$113
☐ PBSSR	Matrx Small Lateral Fixed-Right(3x5.5")	\$113
☐ PBLSL	Matrx Standard Lateral Fixed-Left(7.5x5")	\$116
☐ PBLSR	Matrx Standard Lateral Fixed-Right (7.5x5")	\$116
☐ PBOLL	Matrx Offset Lateral Fixed Left	\$125
☐ PBOLR	Matrx Offset Lateral Fixed Right	\$125

1. Only available on the E2 backrests.

MATRIX SWING AWAY LATERALS - FOR ELITE AND ELITE TR BACKS

Matrx Lateral Trunk Support Left⁽¹⁾

☐ EPADXL	Size XS (3.25x4.25")	\$115
☐ EPADSL	Size S (3.5x5.25")	\$115
☐ EPADML	Size M (3.75x5.75")	\$115
☐ EPADLL	Size L (4.25x6.5")	\$115

Matrx Lateral Trunk Support Right⁽¹⁾

☐ EPADXR	Size XS (3.25x4.25")	\$115
☐ EPADSR	Size S (3.5x5.25")	\$115
☐ EPADMR	Size M (3.75x5.75")	\$115
☐ EPADLR	Size L (4.25x6.5")	\$115

1. Must choose swing away hardware below

MATRIX SWING AWAY HARDWARE - FOR ELITE AND ELITE TR BACKS

☐ SAEHBL	Matrx Elite Swing Away Hardware Left	\$250
☐ SAHEBR	Matrx Elite Swing Away Hardware Right	\$250
☐ OSAHL	Matrx Offset Elite Swing Away Hardware Left	\$250
☐ OSAHR	Matrx Offset Elite Swing Away Hardware Right	\$250

OTHER MATRIX ACCESSORIES

Multi Axis Upper Extreimity Support

☐ MACESL	Multi Axis Arm Support - Left ⁽¹⁾	\$225
☐ MACESR	Multi Axis Arm Support - Right ⁽¹⁾	\$225

☐ PFLEX-M	PostureFlex Abdominal Support Height 6.5", length 16" to 23"	\$295
----------------------------------	--	-------

☐ PFLEX-L	PostureFlex Abdominal Support Height 7", length 20" to 32"	\$295
----------------------------------	--	-------

☐ QRTAB	Quad release tab for easy independent application and removal of PostureFlex	\$35
--------------------------------	--	------

1. Allows rotation of the armpad, only compatible with Gel pad, Flat pad or Ergonomic Arm Troughs

Maxx Rehab Seat on ROVI VR2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MAXX LATERALS TRUNK SUPPORTS⁽¹⁾

Left Fixed Mount Lateral Trunk Support⁽¹⁾

Select Left Trunk Support Size

☐ LATXL	Size XS (3.25" x 4.25")		\$100
☐ LATSL	Size S (3.5" x 5.25")		\$100
☐ LATML	Size M (3.75" x 5.75")		\$100
☐ LATLL	Size L (4.25" x 6.5")		\$100

Right Fixed Mount Lateral Trunk Support⁽¹⁾

Select Right Trunk Support Size

☐ LATXR	Size XS (3.25" x 4.25")		\$100
☐ LATSR	Size S (3.5" x 5.25")		\$100
☐ LATMR	Size M (3.75" x 5.75")		\$100
☐ LATLR	Size L (4.25" x 6.5")		\$100

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs.

MAXX STYLE SWING-AWAY TRUNK SUPPORT HARDWARE⁽¹⁾

☐ SALATL Left MAXX Style Swing Away HDW \$195

Select Left Trunk Support Size

☐ LATXL	Size XS (3.25" x 4.25")		\$100
☐ LATSL	Size S (3.5" x 5.25")		\$100
☐ LATML	Size M (3.75" x 5.75")		\$100
☐ LATLL	Size L (4.25" x 6.5")		\$100

☐ SALATR Right MAXX Style Swing Away HDW⁽¹⁾ . . . \$195

Select Right Trunk Support Size

☐ LATXR	Size XS (3.25" x 4.25")		\$100
☐ LATSR	Size S (3.5" x 5.25")		\$100
☐ LATMR	Size M (3.75" x 5.75")		\$100
☐ LATLR	Size L (4.25" x 6.5")		\$100

☐ ATLL Additional Telescoping Link - Left⁽²⁾ . . . \$75

☐ ATLR Additional Telescoping Link - Right⁽²⁾ . . . \$75

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs

2. Up to 2.75" medial offset

FIXED HIP/ KNEE SUPPORTS

Lateral Hip Support with Fixed Mounting Hardware Left⁽¹⁾

Select Pad Size and mounting position at Hip:

☐ HSPXLH	Size XS (5.5" W x 4" H) mounted at Hip	..	\$120
☐ HSPSLH	Size S (7" W x 6" H) mounted at Hip	..	\$120
☐ HSPMLH	Size M (9" W x 4" H) mounted at Hip	..	\$120
☐ HSPLLH	Size L (13" W x 4" H) mounted at Hip	..	\$120

Select Pad Size and mounting position at Knee:

☐ HSPXLK	Size XS (5.5" W x 4" H) mounted at Knee	..	\$120
☐ HSPSLK	Size S (7" W x 6" H) mounted at Knee	..	\$120
☐ HSPMLK	Size M (9" W x 4" H) mounted at Knee	..	\$120
☐ HSPLLK	Size L (13" W x 4" H) mounted at Knee	..	\$120

FIXED HIP/ KNEE SUPPORTS - CONTINUED

Lateral Hip Support with Fixed Mounting Hardware Right⁽¹⁾

Select Pad Size and mounting position:

☐ HSPXRH	Size XS (5.5" W x 4" H) mounted at Hip	..	\$120
☐ HSPSRH	Size S (7" W x 6" H) mounted at Hip	..	\$120
☐ HSPMRH	Size M (9" W x 4" H) mounted at Hip	..	\$120
☐ HSPLRH	Size L (13" W x 4" H) mounted at Hip	..	\$120

Select Pad Size and mounting position at Knee:

☐ HSPXRK	Size XS (5.5" W x 4" H) mounted at Knee	..	\$120
☐ HSPSRK	Size S (7" W x 6" H) mounted at Knee	..	\$120
☐ HSPMRK	Size M (9" W x 4" H) mounted at Knee	..	\$120
☐ HSPLRK	Size L (13" W x 4" H) mounted at Knee	..	\$120

1. Must choose swing away hardware below if required.

REMOVABLE HIP/ KNEE SUPPORT HARDWARE

Lift Off Removable Hardware for Hip Support

☐ LOHL	Left		\$195
☐ LOHR	Right		\$195

MAXX Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports

Select mounting position:

☐ MHSHLH	Left mounted at Hip		\$215
☐ MHSHLK	Left mounted at Knee		\$215
☐ MHSHRH	Right mounted at Hip		\$215
☐ MHSHRK	Right mounted at Knee		\$215

Swing Away Removable Hardware for Hip Supports

Select mounting position:

☐ SHSHLH	Left mounted at Hip		\$215
☐ SHSHLK	Left mounted at Knee		\$215
☐ SHSHRH	Right mounted at Hip		\$215
☐ SHSHRK	Right mounted at Knee		\$215

Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports⁽¹⁾

Select mounting position:

☐ EMQHLH	Left mounted at Hip		\$215
☐ EMQHLK	Left mounted at Knee		\$215
☐ EMQHRH	Right mounted at Hip		\$215
☐ EMQHRK	Right mounted at Knee		\$215

1. Recommended when finished seat cushion thickness will be greater than 3.5" high.