



REQUIRED INFORMATION


Quote Order

Date of Order: _____ Dealer Account #: _____
 Dealer Name: _____
 PO #: _____ Tag _____
 Purchasing Contact: _____
 Phone: _____ Fax: _____
 E-mail: _____
 RTS/Therapist: _____
 Ship to Address: _____
 City: _____ Province: _____
 Postal Code: _____ Client Gender: M F
 Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

| | |
|-------------------------|------------------------|
| Client Height: | |
| Client Weight (lbs): | |
| A. Seat to Shoulder: | |
| B. Trunk Depth: | |
| C. Chest Width: | |
| D. Knee to Back: | LEFT SIDE RIGHT SIDE |
| E. Seat to Top of Head: | |
| F. Elbow to Hand: | |
| G. Seat to Elbow: | |
| H. Hip Width: | |
| I. Knee to Heel: | LEFT SIDE RIGHT SIDE |
| Cushion Thickness: | |



Weight Capacity is up to 200lbs unless specified otherwise

MULTI-POSITION STANDING MAXX SYSTEM



BASE VITALS

| | |
|--|---------------------------|
| ADP CODES | WP0001287 |
| SUSPENSION SYSTEM | Active Ride Control (ARC) |
| BASE WIDTH | 25.5" |
| BASE LENGTH | 36.5" |
| DRIVE WHEEL SIZE ⁽¹⁾ | 14" |
| CASTER WHEEL SIZE (front/rear) | 6" |
| TURNING RADIUS | 20.5" |
| DRIVETRAIN | 4-pole EAD Motors |
| BATTERIES | 2 x M34 |
| BATTERY CHARGER | 8A, Off Board |
| BATTERIES WEIGHT (each) | 38.5 lbs (M34) / 17.5 kg |
| RANGE ⁽²⁾ | 17.7 miles /28.5 km/h |
| WEIGHT CAPACITY | 200 lbs / 136 kg |
| BASE WEIGHT (w/o batteries) | 157 lbs /71 kg |
| GROUND CLEARANCE | 3" |
| MAXIMUM SPEED | 5.8 mph |

⁽¹⁾ All A3 wheels come with solid inserts
⁽²⁾ Driving range is affected by selection of battery, drive motor, and tire options. It may also vary due to other factors, such as user weight, terrain, battery condition, tire condition, as well as variations in drive and control systems.

MPS Mini Power Standing System Rovi A3

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

BASE FENDER COLOR

- Select Color: **STD** **With Lights⁽¹⁾ . \$995**
- Blizzard White Part # 800914-00200-03 Part # 800914-00201-03
 - Matte Black Part # 800914-00200-50 Part # 800914-00201-50
 - Solar Yellow Part # 800914-00200-69 Part # 800914-00201-69
 - Motion Green Part # 800914-00200-36 Part # 800914-00201-36
 - Original Orange Part # 800914-00200-23 Part # 800914-00201-23
 - Bright Blue Part # 800914-00200-40 Part # 800914-00201-40
 - Red Delicious Part # 800914-00200-21 Part # 800914-00201-21
 - Party Pink Part # 800914-00200-2C Part # 800914-00201-2C
 - Purple Reign Part # 800914-00200-80 Part # 800914-00201-80

1. Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. Must select Joystick for Lights feature or Omni for operation.



TIRE OPTIONS

- Standard Solid Tires** Part # STIRES **STD**
- Pneumatic Tires** Part # PTIRES **N/C**

BATTERIES

- 2 x MK GROUP 34 Battery** Part # M34 SLD G (Two) **\$650**

WHEELCHAIR TRANSPORT BRACKETS

- Wheelchair Transport Brackets⁽¹⁾** Part # XBKTS **N/C**

1. For unoccupied use only

ROVI R-NET BASE SELECTION

- ROVI A₃ 5.8 mph Quick Speed Motors⁽¹⁾**
Part # 888WNLEA3-RNET58. **\$7,595**
- OnTraxx Enhanced Driving Module⁽²⁾**
Part # OEMD. **\$1,195**
- R-Net Expandable Controller**
Part # D51109. **N/C**

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt, Weight capacity is 250 lbs.
2. OnTraxx must be selected with the Standing Module.

ELECTRONIC PROGRAMMER

- R-Net PC Programmer with Dongle** Part # D50612. **\$603**
- Diagnostic Test Tool (DTT) with RNet & VR2 Cables**
Part # D50996. **\$1,876**

OPTIONAL R-NET JOYSTICKS

- R-Net Color Remote Joystick with Monojacks⁽¹⁾⁽²⁾**
Part # D51165 **\$875**
- R-Net Color Remote Joystick with Monojacks for Lights⁽¹⁾** Part # D50982. **\$975**
- R-Net Advanced Joystick⁽¹⁾⁽²⁾** Part # D51711 **\$1,175**
3.5" high-resolution VGA color display IR Control as standard. Includes Integrated Bluetooth module.
- R-Net Omni Display for Specialty Drive Input⁽¹⁾**
Part # D50881 **\$1,395**
- R-Net Attendant Control⁽¹⁾⁽²⁾** Part # D50882. **\$995**

1. Requires an R-Net Controller
2. Light package not available with this joystick

JOYSTICK MOUNT

- Fixed** **N/C**
Select Mounting Position: Left Right
 Part # 800907-001-16-ML 800907-001-16-MR
- Height Adjustable Swing-Away Quad Link** **\$344**
Select Mounting Position: Left Right
 Part # SAQL SAQR
- Adjustable Swing-Away Mount for Omni** **\$339**
Select Mounting Position: Left Right
 Part # ARMS260-Rovi-ML ARMS260-Rovi-MR

JOYSTICK HANDLES

- "T" Handle for Joystick Knob** Part # SA80941 **\$115**
- Ball Handle for Joystick Knob** Part # SA72441 **\$115**
- BodyPoint Joystick Handle.** **\$120**
 3" U-Shaped Part # PC101B-20
 4" U-Shaped Part # PC102B-20
 Domed Rubber Part # PC107B-20

ELECTRONIC OPTIONS

- R-Net Output Module (EMC)⁽¹⁾⁽²⁾** Part # D50883 **\$995**
- R-Net Bluetooth Module⁽¹⁾⁽²⁾** **\$895**
Select OS: iOS Apple Part # D51580 Android/PC Part # D51111
- Smart Phone Holder** Part # SPH **\$50**
- Motion Dual USB Charger** Part # MUSB **N/C**

1. Requires an R-Net Controller
2. Must order R-Net Color Remote joystick with these options

MPS Mini Power Standing System Rovi A3

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

Sections in yellow boxes are mandatory. Power modules are not retro-fittable in the field. Add-ons in the future must be done at Motion Concepts.

SELECT SEAT-TO-FLOOR HEIGHT

Select STF⁽¹⁾: 18.5" 19.5"
Part # STF185 STF195

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on conditions of the power base and seat depth greater than 18".

NOTE: Systems Include Seat Pan, Back Pan, Standard Armrests and Chest Strap. They are designed for use with a 2"- 3" cushion.

SYSTEM

Modular Maxx CG Tilt (50°) and Power Recline (170°) System w/Extended Shear Reduction (ESR)
Part # MDCGTRESR-45 **\$11,569**

POWER STAND MODULE MUST BE ORDERED WITH SYSTEM ABOVE

MPS Maxx Power Stand Module for Tilt/Recline/ESR⁽¹⁾
Part # PSMTRE **\$7,995**

POWER ADJUSTABLE SEAT HEIGHT MODULE -
MUST BE ORDERED WITH SYSTEM ABOVE

7" Power Adjustable Seat Height Module⁽¹⁾
Part # ESM7 **\$1,995**

1. The legrest will extend slightly during elevation.

SPECIFICATIONS

Select Seat Width and Ultra Rail Width Adjustment

12" 13" 14" 15" 16" 17" (Adjustable from 12"- 17")

Part# MSSW12 MSSW13 MSSW14 MSSW15 MSSW16

Select Seat Depth and Ultra Rail Depth Adjustment

12" 13" 14" 15" 16" 17" (Adjustable from 12"- 17")

Part# MSSD12 MSSD13 MSSD14 MSSD15 MSSD16 MSSD17

CUSHION COVER REFERENCE GUIDE

| Cushion | Kid*ab*ra | Meshtex | Startex/Spacetex | Pink Lolipop/Blue Blizzard | Infection Control |
|--------------------|------------|------------|------------------|----------------------------|-------------------|
| PS/PSP | N/C | N/A | STD | STD | Reversible |
| Kid*ab*ra | STD | \$75 | N/A | \$75 | N/A |
| Vi | \$75 | STD | \$50 | \$50 | N/A |
| Libra/Libra Adjust | N/A | N/A | STD | N/A | \$50 |

MATRIX SEAT CUSHION OPTIONS

Matrix PS Cushion Comes Standard in Startex Fabric (Reversible) **.\$569**

PS1212 PS1214 PS1216 PS1312 PS1314
 PS1316 PS1412 PS1414 PS1416 PS1616
 PS1716

Matrix PSP Cushion Comes Standard in Startex Fabric (Reversible) **\$645**

PSP1212 PSP1214 PSP1216 PSP1312 PSP1314
 PSP1316 PSP1412 PSP1414 PSP1416 PSP1616
 PS1716

Fabric Options for PS and PSP Cushion N/C

Kid*ab*ra (None Reversible) Part # PSC-K
 Pink Lollipop (None Reversible) Part # PSC-P
 Blue Blizzard (None Reversible) Part # PSC-B

Matrix Kid*ab*ra Cushion Comes Standard in Kid*ab*ra Fabric **.\$576**

KA1212-VI KA1214-VI KA1412-VI KA1414-VI

Fabric Options for Kid*ab*ra Cushion⁽¹⁾

Meshtex Part # CV-M \$75
 Pink Lollipop Part # CV-P \$75
 Blue Blizzard Part # CV-B \$75

Matrix Vi Cushion Comes Standard in Meshtex Fabric **.\$624**

MA1415-VI MA1616-Vi MA1716-Vi

Fabric Options for Matrix Vi Cushion⁽¹⁾

Kid*ab*ra (Non Reversible) Part # MA-K \$75
 Startex Fabric Smooth Side up Part # SFVS \$50 (Non Reversible)
 Startex Fabric Fabric Side up Part # SFVF \$50 (Non Reversible)

Matrix Libra Adjust Cushion Comes Standard in Startex Fabric **.\$985**

LC1412-ADJ LC1414-ADJ LC1416-ADJ
 LC1516-ADJ LC1616-ADJ LC1716-ADJ

Matrix Libra Cushion Comes Standard in Startex Fabric **.\$745**

LC1412 LC1414 LC1416 LC1516 LC1616
 LC1716

Fabric Options for Libra and Libra Adjust Cushions⁽¹⁾

Spacetex Part # SFLF N/C
 Infection Control Part # ICFL \$50

1. Lead time up to 2 weeks

MPS Mini Power Standing System Rovi A3

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CHEST STRAPS

Mini BP Padded Chest Strap 2 point⁽¹⁾
Part # MCS..... **STD**

BP Monoflex, Underarm Release\$195

XS(2.75"x12.25") Part # MFUXS S(3"x15.5") Part # MFUS
 M(4"x18.25") Part # MFUM L(4.5"x21.5") Part # MFUL

BP Monoflex, Center Release\$195

XS(2.75"x12.25") Part # MFCXS S(3"x15.5") Part # MFCUS
 M(4"x18.25") Part # MFCM L(4.5"x21.5") Part # MFCL

Miniflex Abdominal support Band 4.5" H
for Matrx Mini backs 11"-14" wide Part # MINIFLEX-XS ..\$245

Miniflex Abdominal Support Band 6" H
for Matrx Mini backs 11"-14" wide Part # MINIFLEX-S ..\$245

1. Adjustable from 14" to 41" in length, pad is 13" x 1.5" in length.

LAP BELTS

Push Button Style Seat Belt **STD**
Part # 35" BELT35 48" BELT48

Mini BP Padded Hip Belt 2 point⁽¹⁾
Part # MBPHB.....**\$175**

1. Adjustable from 14" to 41" in length, pad is 13" x 1.5" in length.

CHEST BAR OPTION

Foam Wrapped Chest Bar Part # FWCB **STD**

Flat Pad Chest Bar Part # FPCB\$450

Select Swing Away Mounting Side: Left Right
Part # ML MR

LAP TRAY OPTION

Height Adjustable Lap Tray in Clear Polycarbonate Plastic (16"w x 12"d) for 12"-15" widths Part # MHALT **\$750**

Height Adjustable Lap Tray in Clear Polycarbonate Plastic (21"w x 12"d) for 16"-17" widths Part # HALT..... **\$750**

MATRIX ELAN HEADREST PAD

Matrx Elan Standard Pad**\$310**

Part # 6" Pad ESP6 10" Pad ESP10 14" Pad ESP14

Matrx Elan Standard Pad with Infection Control .. **\$335**

Part # 6" Pad ESP6-IC 10" Pad ESP10-IC 14" Pad ESP14-IC

Matrx Elan Occipital Pad**\$341**

Part # 9" Pad EOP9 12" Pad EOP12

Matrx Elan Occipital Pad with Infection Control**\$366**

Part # 9" Pad EOP9-IC 12" Pad EOP12-IC

MATRIX ELAN HEADREST PAD - CONTINUED

Matrx Elan 4-Point Pad (9"W x 7"H")

Part # Standard Cover E4PTSM**\$448**

Infection Control E4PTSM-IC**\$473**

Matrx Elan 4-Point Pad (11"W x 10"H")

Part # Standard Cover E4POINT **\$448**

Infection Control E4POINT-IC **\$473**

MOTION CONCEPTS HEADREST PAD

Motion Concepts Standard⁽¹⁾ Part # MGS **N/C**

1. Unless selecting a "Removable" hdw option, this Headrest will come standard w/ Fixed Straight Mounting Hardware (non-removable)

HEADREST MOUNTING HARDWARE OPTIONS

Fixed Mounting (non-removable) Part # FM..... **N/C**

Elan Mini Headrest Hardware - Multi-Axis Removable
Part # MEMH**\$285**

Elan Headrest Hardware - Multi-Axis Removable
Part # MEHW**\$285**

Motion Concepts Multi - Axis Removable Mounting
Part # MAHU**\$225**

OMIT: Headrest Part # OHR **N/C**

FINISHED BACK HEIGHTS⁽¹⁾⁽²⁾

Part# 16" FBR16 17" FBR17 18" FBR18 19" FBR19 20" FBR20 21" FBR21 22" FBR22

1. 5"-6" gap is required for Stander systems.

2. Recline push bar will be at 16" height from seat pan.

BACK OPTIONS

Standard Rehab Back and Back cushion Part # SRB
Finished height from 16" to 22" **N/C**

Additional Super Soft/HR Foam⁽¹⁾ Part # SS..... **\$100**

Matrx Fabric Upcharge Part # MFU..... **\$100**

Planar interface plate

Select one: Flat Part # PIPF **\$125**

1. Additional foam added may impact seat depth.

BACK COVER REFERENCE GUIDE

| Back | Kid*ab*ra | Meshtex | Startex | Pink Lolipop | Blue Blizzard | Infection Control |
|------------|-----------|---------|---------|--------------|---------------|-------------------|
| Mini | STD | N/A | N/C | N/C | N/C | N/A |
| Elite | N/A | N/A | N/A | N/A | N/A | N/A |
| Elite TR | N/A | N/A | N/A | N/A | N/A | N/A |
| Elite Deep | N/A | N/A | STD | N/A | N/A | N/A |
| PB | N/A | N/A | STD | N/A | N/A | N/A |
| PB Deep | N/A | N/A | STD | N/A | N/A | N/A |

UPGRADABLE BACK OPTIONS

Matrix Mini Back (3" Contour)\$845
Comes Standard in Kid*ab*ra Fabric

- MINI1207 MINI1210 MINI1212 MINI1214
- MINI1307 MINI1310 MINI1312 MINI1314
- MINI1407 MINI1410 MINI1412 MINI1414

Fabric Options for Mini BackN/C

- Startex (None Reversible) Part # MINI-S
- Pink Lollipop (None Reversible) Part # MINI-P
- Blue Blizzard (None Reversible) Part # MINI-B

High Back 16"-17" wide\$595

- HB1616 HB1618 HB1620 HB1716 HB1718 HB1720

Matrix Elite Back 14"-17" wide (3" Contour)\$766

Comes Standard in Meshtex Fabric

- PBE1410 PBE1412 PBE1414 PBE1416
- PBE1418 PBE1420 PBE1510 PBE1512
- PBE1514 PBE1516 PBE1518 PBE1520
- PBE1610 PBE1612 PBE1614 PBE1616
- PBE1618 PBE1620 PBE1710 PBE1714
- PBE1716 PBE1718 PBE1720

Matrix Elite TR Back 15"-17" wide (3" Contour)\$866

Comes Standard in Meshtex Fabric

- PBE1516-TR PBE1518-TR PBE1520-TR
- PBE1616-TR PBE1618-TR PBE1620-TR
- PBE1716-TR PBE1718-TR PBE1720-TR

Matrix Elite Deep Back 14"-17" wide (6" Contour)\$766

Comes Standard in Meshtex Fabric

- EDB1412 EDB1416 EDB1418 EDB1420
- EDB1510 EDB1512 EDB1514 EDB1516
- EDB1518 EDB1520 EDB1610 EDB1614
- EDB1616 EDB1618 EDB1620 EDB1710
- EDB1712 EDB1714 EDB1716 EDB1720

Matrix PB Back 14"-17" wide (5" Contour)\$734

Comes Standard in Meshtex Fabric

- MPB1412 MPB1416 MPB1612 MPB1616
- MPB1712 MPB1716 MPB1720

Matrix PB Deep Back 14"-17" wide (7" Contour)\$766

Comes Standard in Meshtex Fabric

- PBD1412 PBD1416 PBD1612 PBD1616
- PBD1712 PBD1716 PBD1720

OTHER UPGRADABLE BACK OPTIONS

- Super Soft Foam for Elite Back Upgrade Part # EBSS\$100
- Infection Control Cover for Elite/PB Backs Part # WSCU\$100
(purchased with back)
- Infection Control Cover for Elite/PB Backs Part # WSC\$150
(purchased without back)

ARMRESTS

Stander Mini Cantilever Cane Mounted Flip Back Armrest (Minimum armrest height is 7"-10")
Part # MSCMFBA **STD**

ARMREST TUBE LENGTH

- LEFT Armrest Tube** Select: Full Part # ATFL
 Desk Part # ATDL
- RIGHT Armrest Tube** Select: Full Part # ATRF
 Desk Part # ATDR

ARM PADS

Modular Arm Pads⁽¹⁾ (each)\$95

- Left** Full(14") Part # MAPTFL Desk(10") Part # MAPTDL

Must select Pad Insert:

- Iskin** - Integrated soft skin foam Part # ISSFL
- Startex Covered Visco foam** Part # SCVFL

- Right** Full(14") Part # MAPTRF Desk(10") Part # MAPTRD

Must select Pad Insert:

- Iskin** - Integrated soft skin foam Part # ISSFR
- Startex Covered Visco foam** Part # SCVFR

Standard Pads (each) **STD**

- Left** Full(14") Part # SAPFL Desk(10") Part # SAPDL
- Right** Full(14") Part # SAPFR Desk(10") Part # SAPDR

Flat Multi Position Pads (each)\$150

- Left** Full(14") Part # MPAFL Desk(10") Part # MPADL
- Right** Full(14") Part # MPAFR Desk(10") Part # MPADR

Gel Pads (each)

- Narrow (2"x12")**\$88
 Left Part # GAPNL **Right** Part # GAPNR
- Wide (3.5"x12")**\$96
 Left Part # GAPWL **Right** Part # GAPWR
- Long (3.5"x14")**\$103
 Left Part # GAPLL **Right** Part # GAPLR

1. Requires interchangeable pad insert.

MATRIX MINI FIXED LATERALS - FOR MINI BACKS

Matrix MINI Fixed Lateral with Pad each)
 Select pad: LEFT (4.5" x 3") Part # MINILATL \$85
 RIGHT(4.5" x 3") Part # MINILATR \$85

MATRIX MINI SWING AWAY LATERALS - FOR MINI BACKS

MINI Swing Away Lateral Mounting Hardware(each)
 LEFT Part # SAHMINIL \$225
 RIGHT Part # SAHMINIR \$225

Mini Lateral Trunk Support RIGHT(each) \$110
 Select Pad Size:
 XS(3.25" x 4.25") Part # EPADXR
 S(3.5" x 5.25") Part # EPADXR

Mini Lateral Trunk Support LEFT(each) \$110
 Select Pad Size:
 XS(3.25" x 4.25") Part # EPADXL
 S(3.5" x 5.25") Part # EPADSL

MATRIX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS

Matrix Standard Fixed LEFT \$116
 Select pad: Standard (7.5" x 5") Part # PBLSL
 Small (5.5" x 3") Part # PBSSL

Matrix Standard Fixed RIGHT \$116
 Select pad: Standard (7.5" x 5") Part # PBLSR
 Small (5.5" x 3") Part # PBSSR

Matrix Offset Fixed LEFT Part # PBOLL \$125

Matrix Offset Fixed RIGHT Part # PBOLR \$125

MATRIX SWING AWAY LATERALS - FOR ELITE AND ELITE TR BACKS

Matrix Lateral Trunk Support LEFT pad (each) \$115
 Select Pad Size:
 XS(3.25" x 4.25") Part # EPADXL S(3.5" x 5.25") Part # EPADSL
 M(3.75" x 5.75") Part # EPADML L(4.25" x 6.5") Part # EPADLL

Matrix Lateral Trunk Support RIGHT pad (each) \$115
 Select Pad Size:
 XS(3.25" x 4.25") Part # EPADXR S(3.5" x 5.25") Part # EPADSR
 M(3.75" x 5.75") Part # EPADMR L(4.25" x 6.5") Part # EPADLR

SELECT SWING AWAY HARDWARE TO CONVERT ABOVE:

Matrix Elite Back Swing Away Hardware
 LEFT Part # SAHEBL \$250
 RIGHT Part # SAHEBR \$250

Matrix Offset Elite Swing Away Hardware
 LEFT Part # OSAHL \$250
 RIGHT Part # OSAHR \$250

MAXX LATERALS PADS (1)

Fixed Mounted Lateral Trunk Support LEFT. \$100

Select Pad Size:
 XS(3.25" x 4.25") Part # LATXL S(3.5" x 5.25") Part # LATSL
 M(3.75" x 5.75") Part # LATML L(4.25" x 6.5") Part # LATLL

Fixed Mounted Lateral Trunk Support RIGHT \$100

Select Pad Size:
 XS(3.25" x 4.25") Part # LATXR S(3.5" x 5.25") Part # LATSR
 M(3.75" x 5.75") Part # LATMR L(4.25" x 6.5") Part # LATLR

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs.
2. The amount of offset available may be impacted by back thickness.

SELECT MAXX SWING-AWAY HARDWARE(NEXT) TO CONVERT ABOVE SUPPORTS:

MAXX SWING-AWAY HARDWARE(1)

Maxx Style Swing-away Multi-Adjustable Mounting Hardware LEFT Part # SALATL \$195
 (Standard with 1" telescoping link gives 0-2.75" medial offset (2))

Maxx Style Swing-away Multi-Adjustable Mounting Hardware RIGHT Part # SALATR \$195
 (Standard with 1" telescoping link gives 0-2.75" medial offset (2))

Additional Telescoping Link(ea) Quantity () \$75
 (Adding an additional 1" telescoping link gives up to 5.5" medial offset(2))
 LEFT Part # ATLL RIGHT Part # ATLR

1. These laterals are not available with Matrix PB, PB Elite Deep, PB Deep Backs and PB HD and HD Deep .
2. The amount of offset available may be impacted by back thickness.

FIXED HIP/KNEE SUPPORTS

Lateral Hip Support Pad with Fixed Mounting LEFT . \$120

Mounted at HIP Pad Size:

- XS(5.5"W x 4"H) Part # HSPXLH
- M(9"W x 4"H) Part # HSPMLH
- S(7"W x 6"H) Part # HSPSLH
- L(13"W x 4"H) Part # HSPLLH

Mounted at KNEE Pad Size:

- XS(5.5"W x 4"H) Part # HSPXLK
- M(9"W x 4"H) Part # HSPMLK
- S(7"W x 6"H) Part # HSPSLK
- L(13"W x 4"H) Part # HSPLLK

Lateral Hip Support Pad with Fixed Mounting RIGHT . \$120

Mounted at HIP Pad Size:

- XS(5.5"W x 4"H) Part # HSPXRH
- M(9"W x 4"H) Part # HSPMRH
- S(7"W x 6"H) Part # HSPSRH
- L(13"W x 4"H) Part # HSPLRH

Mounted at KNEE Pad Size:

- XS(5.5"W x 4"H) Part # HSPXRK
- M(9"W x 4"H) Part # HSPMRK
- S(7"W x 6"H) Part # HSPSRK
- L(13"W x 4"H) Part # HSPLRK

SELECT MAXX SWING-AWAY HARDWARE(NEXT) TO CONVERT ABOVE SUPPORTS

REMOVABLE & SWING AWAY HARDWARE UPGRADE

HARDWARE TO CONVERT ABOVE HIP SUPPORT TO UPGRADE TO REMOVEABLE/SWINGAWAY HIP OR KNEE SUPPORTS

Lift Off Removable Hardware for Hip Support

(Hip position only)

- LEFT Part # LOHL \$195
- RIGHT Part # LOHR \$195

Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports

- LEFT Mounted at Hip Part # MHSHLH \$215
- RIGHT Mounted at Hip Part # MHSHRH \$215
- LEFT Mounted at Knee Part # MHSHLK \$215
- RIGHT Mounted at Knee Part # MHSHRK \$215

Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports⁽¹⁾

- LEFT Mounted at Hip Part # EMQHLH \$215
- RIGHT Mounted at Hip Part # EMQHRH \$215
- LEFT Mounted at Knee Part # EMQHLK \$215
- RIGHT Mounted at Knee Part # EMQHRK \$215

1. Recommended when finished seat cushion thickness will be greater than 3.5" high.

POWER CENTER MOUNT FOOT PLATFORM OPTIONS

| MUST BE COMPLETED | Measurement in Inches |
|-------------------------------|-----------------------|
| Knee to Heel: | |
| Cushion Thickness: (Subtract) | |
| Cushion Compression:(Add) | |
| Seat Pan to Footplate:(Total) | |

POWER STANDER KNEE BLOCKS/SUPPORT

Modular Knee Blocks with Height, Depth and Angle Adjustable Mounting Hardware⁽¹⁾⁽²⁾ STD

Select Pad Size:

- (6") Part # MKB6
- (8") Part # MKB8

Eccentric Knee Bolster Mounting Hardware

(Allows individual adjustment of knee blocks to accommodate leg discrepancies)

Part# EKBMH \$250

1. 1.5" of height adjustment and 3" of depth adjustment and 2" of width independent adjustment.
2. Knee blocks are removable and stored on the side of the chair when not in use.

POWER CENTER MOUNT FOOT PLATFORM OPTIONS

Mini Power Elevating BDC Belt Drive Center Mount Articulating Foot Platform Part # MBDC \$3,295


Seat pan to footplate: 8"-12" Part # SPTF22
Comes standard with individual calf pads (5.5" H x 4" W)

Seat pan to footplate: 11"-15.5"/16.5⁽¹⁾ Part # SPTF18
Comes standard with individual calf pads (7.5" H x 6" W)

1. With 18.5" seat-to-floor selected the maximum seat pan to footplate is 15.5", with 19.5" seat-to-floor height selected the seat pan to footplate is 16.5".

MUST SELECT INDIVIDUAL FOOT PLATE SIZE

Select: Footplate Option:



- Left: S(5"W x 7.5"D) M(5.5"Wx9.5"D) L(6"Wx11.5"D) . . . \$165
- Part # IFSL IFML IFLL
- Right: S(5"W x 7.5"D) M(5.5"Wx9.5"D) L(6"Wx11.5"D) . . . \$165
- Part # IFSR IFMR IFLR

