

## REQUIRED INFORMATION

Quote  Order

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_ Tag \_\_\_\_\_

Purchasing Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_

RTS/Therapist: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Province: \_\_\_\_\_


Postal Code: \_\_\_\_\_ Client Gender:  M  F

Special Client Conditions: \_\_\_\_\_

## REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



**Customer Service: 1.866.748.7943 | Fax: 1.905.695.0583 | orderscdn@motionconcepts.com**

### ROVI BASE STATISTICS

ADP CODES	WP0001287
SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH X <sub>3</sub>	23.25"
BASE WIDTH A <sub>3</sub>	25.25"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE <sup>(1)</sup>	14"
CASTER WHEEL SIZE (front/rear)	6"
TURNING RADIUS X <sub>3</sub>	20"
TURNING RADIUS A <sub>3</sub>	20.5"
DRIVETRAIN	4-pole EAD Motors

BATTERIES	2 x M34
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs (M34) / 17.5 kg
RANGE <sup>(2)</sup>	17.7 miles /28.5 km/h
WEIGHT CAPACITY	300 lbs / 136 kg
BASE WEIGHT X <sub>3</sub> (w/o batteries)	157 lbs /71 kg
BASE WEIGHT A <sub>3</sub> (w/o batteries)	159 lbs /72 kg
GROUND CLEARANCE	3"
MAXIMUM SPEED	6.2 mph




\*Shown with Motion Concepts Single Post Back for Contour Backrest. See page 7

<sup>(1)</sup> All X3 and A3 wheels come with solid inserts

<sup>(2)</sup> Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

# Maxx Rehab Seat on ROVI X3/A3 Base

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

## X3/A3 BASE OPTIONS

- 888WNLEX3-RNET58  
ROVI X<sub>3</sub> 5.8 MPH QUICK SPEED MOTORS<sup>(1)(2)</sup> . \$7,395
- 888WNLEX3-RNET62  
ROVI X<sub>3</sub> 6.2 MPH HIGH SPEED MOTORS<sup>(1)(2)</sup> . \$7,395  
*must order OnTraxx Enhanced Driving Module*

---

- 888WNLEA3-RNET58  
ROVI A<sub>3</sub> 5.8 MPH QUICK SPEED MOTORS<sup>(1)(2)</sup> \$7,595
- 888WNLEA3-RNET62  
ROVI A<sub>3</sub> 6.2 MPH HIGH SPEED MOTORS<sup>(1)(2)</sup> . \$7,595  
*must order OnTraxx Enhanced Driving Module*

---

- OEDM OnTraxx Enhanced Driving Module . . . . . \$1,175
- D51109 R-Net Expandable Controller . . . . . N/C

- Price of the powerbase includes: Active Ride Control suspension, tie-down hard ware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt. Weight capacity is 300 lbs.
- When selecting an increased speed motor package there is a reduction in torque output by roughly 15%.

## FENDER COLOUR



- 800914-00200-03 Blizzard White. . . . . N/C
- 800914-00200-50 Matte Black . . . . . N/C
- 800914-00200-69 Solar Yellow . . . . . N/C
- 800914-00200-36 Motion Green. . . . . N/C
- 800914-00200-23 Original Orange. . . . . N/C
- 800914-00200-40 Bright Blue . . . . . N/C
- 800914-00200-21 Red Delicious . . . . . N/C
- 800914-00200-2C Party Pink . . . . . N/C
- 800914-00200-80 Purple Reign . . . . . N/C

## FENDER COLOUR WITH LIGHTS<sup>(1)</sup>

- 800914-00201-03 Blizzard White. . . . . \$995
- 800914-00201-50 Matte Black . . . . . \$995
- 800914-00201-69 Solar Yellow . . . . . \$995
- 800914-00201-36 Motion Green. . . . . \$995
- 800914-00201-23 Original Orange. . . . . \$995
- 800914-00201-40 Bright Blue . . . . . \$995
- 800914-00201-21 Red Delicious . . . . . \$995
- 800914-00201-2C Party Pink . . . . . \$995
- 800914-00201-80 Purple Reign . . . . . \$995

- Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. Must select Joystick for Lights feature or Omni for operation.

## WHEELCHAIR TRANSPORT BRACKETS

- XBKTS Wheelchair Transport Brackets<sup>(1)</sup> . . . . . STD

- For unoccupied use only

## TIRE OPTIONS

- STIRES Standard Solid Tires . . . . . STD
- PTIRES Pneumatic Tires. . . . . N/C

## BATTERIES

- M34 SLD G MK Group Battery x 2. . . . . \$650

## JOYSTICK SELECTION

- D51122 R-Net LED Joystick<sup>(1)(2)</sup> . . . . . N/C
- D51108 R-Net LED Joystick for Lights<sup>(1)</sup> . . . . . \$500
- D51165 R-Net Color Remote Joystick with Monojacks<sup>(1)(2)</sup> . . . . . \$875
- D50982 R-Net Color Remote Joystick / Monojacksfor Lights<sup>(1)</sup> . . . \$975
- D51711 R-Net Advanced Joystick<sup>(1)(3)</sup> . . . . . \$1,595
- D50881 R-Net Omni Display for Specialty Drive Input<sup>(1)</sup> . . . . \$1,395
- D50882 R-Net Attendant Control . . . . . \$995

- Requires an R-Net Controller
- Light package not available with this joystick
- 3.5" high-resolution VGA colour display, IR Control as standard, includes integrated Bluetooth Module.

## ELECTRONIC OPTIONS

- D50883 R-Net Output Module (EMC)<sup>(1)(2)</sup> . . . . . \$995
- D51580 R-Net Bluetooth Module IOS Apple<sup>(1)(2)</sup> . . . . . \$895
- D51111 R-Net Bluetooth Module Android/PC<sup>(1)(2)</sup> . . . . . \$895
- SPH Smart Phone Holder . . . . . \$50
- MUSB Motion Dual USB Charger . . . . . N/C

- Requires an R-Net Controller
- Must order R-Net Colour Remote joystick with these options

## JOYSTICK MOUNT

- 800914-001-16-ML Fixed Joystick Mounting Left . . . . . N/C
- 800914-001-16-MR Fixed Joystick Mounting Right . . . . . N/C
- SAQL Height Adjustable Swing-Away Quad Link Left. . . . . \$344
- SAQR Height Adjustable Swing-Away Quad Link Right . . . . . \$344
- ARMS260-ROVI-ML Adjustable Swing-Away Mount - Omni Left . . \$339
- ARMS260-ROVI-MR Adjustable Swing-Away Mount - Omni Right . \$339

## JOYSTICK KNOB - OPTIONAL

- SA80941 "T" Handle for Joystick Knob . . . . . \$115
- SA72441 Ball Handle for Joystick Knob . . . . . \$115
- PC101B-20 BodyPoint Joystick Handle 3" U-Shaped . . . . . \$120
- PC102B-20 BodyPoint Joystick Handle 4" U-Shaped . . . . . \$120
- PC107B-20 BodyPoint Joystick Handle Domed Rubber . . . . . \$120

## ELECTRONIC PROGRAMMER - OPTIONS

- D50612 R-Net PC Programmer with Dongle. . . . . \$603
- D50996 Diagnostic Test Tool (DTT) with VR2 & RNet Cables . \$1,876

## SEAT-TO-FLOOR HEIGHT<sup>(1)</sup>

- STF17 17" Seat-To-Floor . . . . . N/C
- STF18 18" Seat-To-Floor . . . . . N/C
- STF19 19" Seat-To-Floor . . . . . N/C
- STF20 20" Seat-To-Floor . . . . . N/C

- Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"

**Maxx Rehab Seat on ROVI X3/A3 Base**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**REHAB SEATING SYSTEM**

- D0195 Maxx Adjustable Rehab Seat . . . . . \$1,200
- PRETILT 5° Fixed Static Tilt Angle (Pre-Tilt)<sup>(1)</sup> . . . . . N/C
- ADJUST Adjustable Static-Tilt<sup>(1)</sup> . . . . . N/C  
(in 1.25° increments from 0° to 5°)  
Specify set at \_\_\_\_\_°

1. The front seat to floor height will be approximately 2" greater than the STF selected above. It will also be dependant on the seat depth selected.

**ADULT SEAT WIDTH SETTINGS**

- SW116 16" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW117 17" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW118 18" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW119 19" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW120 20" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW219 19" Wide (Adjustable from 19" to 22") . . . . . \$460
- SW220 20" Wide (Adjustable from 19" to 22") . . . . . \$460
- SW221 21" Wide (Adjustable from 19" to 22") . . . . . \$460
- SW222 22" Wide (Adjustable from 19" to 22") . . . . . \$460
- BW123 23" Wide (Adjustable from 20" to 24")<sup>(1)</sup> . . . . . \$460
- BW124 24" Wide (Adjustable from 20" to 24")<sup>(1)</sup> . . . . . \$460

1. This option is only available on the Rovi A3 Base.

**ADULT ULTRA RAIL SEAT DEPTH SETTINGS**

- SD116 16" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD117 17" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD118 18" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD119 19" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD120 20" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD219 19" Deep (Adjustable from 19" to 23") . . . . . \$460
- SD220 20" Deep (Adjustable from 19" to 23") . . . . . \$460
- SD221 21" Deep (Adjustable from 19" to 23") . . . . . \$460
- SD222 22" Deep (Adjustable from 19" to 23") . . . . . \$460
- SD223 23" Deep (Adjustable from 19" to 23") . . . . . \$460

**SMALL ADULT SEAT WIDTH SETTINGS<sup>(1)</sup>**

- SW314 14" Wide (Adjustable from 14" to 15") . . . . . N/C
- SW315 15" Wide (Adjustable from 14" to 15") . . . . . N/C

**SMALL ADULT SEAT DEPTH SETTINGS<sup>(1)</sup>**

- SD314 14" Deep (Adjustable from 14" to 18") . . . . . N/C
- SD315 15" Deep (Adjustable from 14" to 18") . . . . . N/C
- SD316 16" Deep (Adjustable from 14" to 18") . . . . . N/C
- SD317 17" Deep (Adjustable from 14" to 18") . . . . . N/C
- SD318 18" Deep (Adjustable from 14" to 18") . . . . . N/C

1. Small Adult Seat Size - weight limit is 175 lbs on these seat sizes.

**MATRIX SEATING OPTIONS**

- Matrix Cushions available on page 8
- Matrix Backs available on page 8
- Matrix Headrests available on page 9



**CHEST STRAPS**

- MFUS BP Monoflex, Underarm RLS S(3" x 15.5") . \$195
- MFUM BP Monoflex, Underarm RLS M(4" x 18.25") \$195
- MFUL BP Monoflex, Underarm RLS L(4.5" x 21.5") . \$195
- MFCS BP Monoflex, Center RLS S(2.5" x 15.5") . . \$195
- MFCM BP Monoflex, Center RLS M(3.25" x 18.25") \$195
- MFCL BP Monoflex, Center RLS L(3.5" x 21.5") . . . \$195

**LAP BELTS**

- BELT48 Push Button Style Lap Belt, Length 48" . . . . . STD
- BELT60 Push Button Style Lap Belt, Length 60" . . . . . N/C
- BELT71 Push Button Style Lap Belt, Length 71" . . . . . N/C
- PLB2 Bodypoint Padded Lap Belt 2 Point<sup>(1)</sup> . . . . . \$250
- PLB4 Bodypoint Padded Lap Belt 4 Point<sup>(1)</sup> . . . . . \$315

1. 67" length only option

**MOTION CONCEPTS HEADREST OPTIONS**

- MCSH Motion Concepts Standard . . . . . N/C
- MCOH Motion Concepts Auto Style . . . . . \$235
- FM Fixed Mounting (non-removable) . . . . . N/C
- MAHU Motion Concepts Multi-Axis Rem.HDW . . . . . \$195
- OHR Omit: Headrest . . . . . N/C

**MATRIX HEADRESTS OPTIONS**

- Matrix Elan Headrests and Hardware available on page 9

**STANDARD BACK CANES<sup>(1)</sup>**

- SBC18 18" Straight Back Canes . . . . . STD
- SBC20 20" Straight Back Canes . . . . . STD
- SBC22 22" Straight Back Canes . . . . . STD
- SBC24 24" Straight Back Canes . . . . . STD
- TBC18 18" 10° Mid-Angle Canes for Sling Back. . . . . N/C
- TBC20 20" 10° Mid-Angle Canes for Sling Back. . . . . N/C
- TBC22 22" 10° Mid-Angle Canes for Sling Back. . . . . N/C
- TBC24 24" 10° Mid-Angle Canes for Sling Back. . . . . N/C
- ABC20 20" Angled (5°) Back Canes . . . . . N/C
- ABC22 22" Angled (5°) Back Canes . . . . . N/C
- ABC24 24" Angled (5°) Back Canes . . . . . N/C
- SBCF18 18" Straight Fold Down Back Canes . . . . . \$250
- SBCF20 20" Straight Fold Down Back Canes . . . . . \$250
- SBCF22 22" Straight Fold Down Back Canes . . . . . \$250
- SBCF24 24" Straight Fold Down Back Canes . . . . . \$250



1. Canes will be set at 90° angle and are adjustable from 75° to 116°

**Maxx Rehab Seat on ROVI X3/A3 Base**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**FINISHED BACK HEIGHT<sup>(1)(2)</sup>**

- 18" FBT18    19" FBT19    20" FBT20    21" FBT21    22" FBT22    23" FBT23    24" FBT24    25" FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.  
2. Finished back height is measured from seat pan to top of backrest.

**BACK OPTIONS**

- SRB Standard Rehab Back & Cushion . . . . .N/C  
 SLING Sling Back Fabric . . . . . \$125  
 ORBC Omit Rigid Back & Cushion . . . . .N/C

**MATRIX BACKREST OPTIONS**

Matrix Backs available on page 8



**ARMRESTS**

- DP1 Adult Dual Post Adj. 9.75"-12.75" . . . . . STD  
 DP2 Adult Dual Post Adj. 12.75"-15.75" . . . . . STD  
 DP3 Adult Dual Post Adj. 15.75"-18.75" . . . . . STD  
 Select:  FP Flat Pouch . . . . . STD  
            GB Glove Boxes . . . . . N/C  
 OAL ADD: OutBack Arm Option Left<sup>(1)</sup> each . . . . . \$195  
 OAR ADD: OutBack Arm Option Right<sup>(1)</sup> each . . . . . \$195  
 OCTA1 Outback Cantilever Tilt Armrest 9"-12"<sup>(3)</sup> . . . \$488  
 OCTA2 Outback Cantilever Titl Armrest 12"-15"<sup>(3)</sup> . . \$488  
 RMCAN Ultra Rail Seat Mnt Flip Back Cantilever  
 MAXX Tilt Arm (9-12.5")<sup>(2)</sup> . . . . . \$488  
 DPM ADD: Dual Post Module to RMCAN (pair) . . . . \$250

1. Only available with Dual Post Armrests.  
2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.  
3. Not available with the Universal Cup Holder or Flip and Fold Lap Tray.

**ARMREST TUBE LENGTH - PLEASE SPECIFY**

- ATFL Left Full . . . . . STD  
 ATDL Left Desk . . . . . STD  
 ATRF Right Full . . . . . STD  
 ATDR Right Desk . . . . . STD

**ARMPADS**

- Standard**  
 SAPFL Left Full (2.25" x 14") . . . . . STD  
 SAPDL Left Desk (2.25" x 10") . . . . . STD  
 SAPFR Right Full (2.25" x 14") . . . . . STD  
 SAPDR Right Desk (2.25" x 10") . . . . . STD

**ARMPADS**

- Waterfall Arm Pads<sup>(1)</sup>**  
 WAPFL Left Full (13") . . . . . N/C  
 WAPDL Left Desk (9") . . . . . N/C  
 WAPFR Right Full (13") . . . . . N/C  
 WAPDR Right Desk (9") . . . . . N/C

- Modular Arm Pad - Left<sup>(2)</sup>**  
 MAPTFL Left Tray Full (14") . . . . . \$95  
 MAPTDL Left Tray Desk (10") . . . . . \$95

Select Pad Insert below:

- ISSFL Iskin - Integrated soft skin foam  
 SCVFL Startex Covered Visco foam

- Modular Arm Pad - Right<sup>(2)</sup>**  
 MAPTFR Right Tray Full (14") . . . . . \$95  
 MAPTDR Right Tray Desk (10") . . . . . \$95

Select Pad Insert below:

- ISSFR Iskin - Integrated soft skin foam  
 SCVFR Startex Covered Visco foam

- Flat Multi Position Arm Pads**  
 MPAFL Left Full (14") . . . . . \$150  
 MPADL Left Desk (10") . . . . . \$150  
 MPAFR Right Full (14") . . . . . \$150  
 MPADR Right Desk (10") . . . . . \$150

- Gel Pads**  
 GAPNL Narrow Left (2" x 12") . . . . . \$88  
 GAPNR Narrow Right (2" x 12") . . . . . \$88  
 GAPWL Wide Left (3.5" x 12") . . . . . \$96  
 GAPWR Wide Right (3.5" x 12") . . . . . \$96  
 GAPLL Long Left (3.5" x 14") . . . . . \$103  
 GAPLR Long Right (3.5" x 14") . . . . . \$103

- Arm Troughs**  
 EATL Ergonomic Arm Trough - Left . . . . . \$121  
 EATR Ergonomic Arm Trough - Right . . . . . \$121  
 FHPL Flat Hand Pad - Left<sup>(3)</sup> . . . . . \$91  
 FHPR Flat Hand Pad - Right<sup>(3)</sup> . . . . . \$91

- Arm Pad Storage/Cup Holder<sup>(4)</sup>**  
 APSL Left . . . . . \$295  
 APSR Right . . . . . \$295  
 OMIT Omit Armpads(7/8" diameter tube) . . . . . N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.  
2. Requires interchangeable pad insert.  
3. Must order with ergonomic arm trough.  
4. Not available on the joystick side.



# Maxx Rehab Seat on ROVI X3/A3 Base

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

## STANDARD FRONT RIGGING


- FCMP Fixed Center Mount Platform<sup>(1)</sup> . . . . . N/C
- IFP Comes with a rubber coated footplate (11.5" W x 10") . . . . . N/C

Select seat pan to footplate:

- SPTF1 9" - 13" . . . . . N/C
- SPTF14 13" - 19" . . . . . N/C

Set at:

- SA70 Set at 70° . . . . . N/C
- SA90 Set at 90° . . . . . N/C
- SA97 Set at 97° . . . . . N/C



1. 3" ground clearance required.

## MANUAL CENTER MOUNT FRONT RIGGING

- MSFCMFP Maxx Style Fixed Center Mount Foot Platform<sup>(1)(2)</sup> . . . . . \$395  
With Quick Release Height Adjustment


- IFP Comes with a rubber coated footplate (11.5" W x 10") . . . . . N/C
- ICP Comes standard with individual calf pads. . . . . N/C

Select seat pan to footplate:

- SPTF1 9" - 13" . . . . . N/C
- SPTF14 13" - 19" . . . . . N/C

Set at:

- SA70 Set at 70° . . . . . N/C
- SA90 Set at 90° . . . . . N/C
- SA97 Set at 97° . . . . . N/C



1. Depth and height adjustable.  
2. 3" ground clearance required.

## INVACARE CENTER MOUNT

- IAKACM Invacare Center<sup>(1)</sup> "CTMT" . . . . . \$595

Select seat pan to footplate:

- SPTF3 9.25" - 12" . . . . . N/C
- SPTF14 12.25" - 15" . . . . . N/C

Select foot plates/footboard

- SM Individual Foot Plates (5.25"W x 7.5"D) . . . . . N/C
- MDM Individual Foot Plates (5.5"W x 9.5"D) . . . . . N/C
- LRG Individual Foot Plates (5.9"W x 11.5"D) . . . . . N/C
- FB Footboard. . . . . STD

Optional:

- CALF Invacare Calf Pads for use with CTMT<sup>(2)</sup> . . . . . \$268

1. 3" ground clearance required.  
2. Not available with the SPTF3 of 9.25"-12"

## MAXX STYLE ELEVATING LEGRESTS

- MMPPSA MAXX Swing Away Manual Elevating - Straight<sup>(1)(2)</sup> . . . . . \$680

Select Receivers:

- RS Straight Receivers . . . . . N/C
- RF Flared 1" Outward Receivers . . . . . N/C


Select Seat Pan to Footplate:

- SPTF5 Length: 10"-14" . . . . . N/C
- SPTF2 Length: 13"-17" . . . . . N/C
- SPTF7 Length: 16"-20" . . . . . N/C

- FCPMS Flat Calf Pad Option (pair) . . . . . \$275
- CCPMS Curved Calf Pad Option (inside with 6")<sup>(3)</sup> (pair) \$275

1. Must choose footplate on page 6  
2. Comes standard with moulded i-skin calf pads.  
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.


## FOOT PLATFORM OPTIONS FOR FIXED/MAXX CMFP



- SFFP Small Foot Platform<sup>(1)</sup> (12"w x 8.5"d) . . . . \$180
- MFFP Medium Platform<sup>(1)</sup> (12"w x 10.75"d) . . . . \$180
- LFFP Large Platform<sup>(1)</sup> (14.5"w x 13"d) . . . . \$250
- XLFFP Extra Large Platform<sup>(1)</sup> (17"w x 13"d) . . . . \$250
- XXLFFP Extra Extra Large Platform<sup>(1)</sup> (21"w x 13"d) \$250

1. May need to adjust leg position to avoid interference with front castors

## INDIVIDUAL FOOT PLATES FOR FIXED/MAXX CMFP



- IFSL Small<sup>(1)</sup> - (5" x 7.5" deep) - Left . . . . . \$165
- IFSR Small<sup>(1)</sup> - (5" x 7.5" deep) - Right . . . . . \$165
- IFML Medium<sup>(1)</sup> - (5.5" x 9.5" deep) - Left . . . . . \$165
- IFMR Medium<sup>(1)</sup> - (5.5" x 9.5" deep) - Right . . . . . \$165
- IFLL Large<sup>(1)</sup> - (6" x 11.5" deep) - Left . . . . . \$165
- IFLR Large<sup>(1)</sup> - (6" x 11.5" deep) - Right . . . . . \$165
- MFOL Multi-Axis Optional Footplate Left (6"W x 10"D) . \$225
- MFOR Multi-Axis Optional Footplate Right (6"W x 10"D) \$225

1. Foot plates may interfere with casters in certain configurations

**Maxx Rehab Seat on ROVI X3/A3 Base**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_


PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**SWING AWAY LEGREST**

- HD70 Heavy Duty 70° Swing Away Footrest(pair)<sup>(1)(2)</sup> . \$495
  - Select Receivers:
    - RS Straight Receivers ..... N/C
    - RF Flared 1" Outward Receivers ..... N/C
  - Select Seat Pan to Footplate:
    - SPTF11 Length: 6"-10" .....N/C
    - SPTF8 Length: 11"-15" .....N/C
    - SPTF2 Length: 13"-17" .....N/C
    - SPTF7 Length: 16"-20" .....N/C
  - HDFCP Flat Calf Pad Option (pair) ..... \$275
  - HD70CCP Curved Calf Pad Option<sup>(3)</sup> (inside with 6") (pair) \$275
  - DCS Dynamic Coil Springs (Pair) for Heavy Duty 70° Swingaway Footrests ..... \$190
  - I70SAF Invacare 70° Swingaway Footrests (Pair) ..... \$325
  - I70-CF Comes with composite foot plates ..... N/C
  - 70HDAIH Invacare 70° Non-Tapered HD Swingaway Footrests with AT5543 Adjustable Angle Foot plates, Heel Loops and Impact Guards<sup>(1)</sup> (Pair) ..... \$457
  - 70TAPAIH Invacare Tapered 70° Swingaway Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards(Pair) ..... \$497
  - ALPT Longer Pivot and Slide Tube for 70° Swingaways (Pair) ..... \$132
  - RLSL Residual Limb Support - Left<sup>(4)</sup> (6.75"W x 9.75"L x 2"H) (each) ..... \$495
  - RLSR Residual Limb Support - Right<sup>(4)</sup> (6.75"W x 9.75"L x 2"H) (each) ..... \$495
1. Must choose footplate on below.  
 2. Heavy duty durable construction (4.5lbs ea.).  
 3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.  
 4. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.

**FOOTPLATE SIZE SELECTION**



1. Certain footplate sizes may not fit all seat widths depending on legrests selected.


**ADJUSTABLE ANGLE FOOTPLATES**

- AFSL Left Small Footplate (4.25"W x 8"D) ..... \$165
- AFSR Right SmallFootplate (4.25"W x 8"D) ..... \$165
- AFML Left Medium Footplate (5.25"W x 8"D) ..... \$165
- AFMR Right Medium Footplate (5.25"W x 8"D) ..... \$165
- AFLL Left Large Footplate (6.25"W x 8"D) ..... \$165
- AFLR Right Large Footplate (6.25"W x 8"D) ..... \$165
- AFXL Left Extra Large Footplate (9"W x 11"D) ..... \$165
- AFXR Right Extra Large Footplate (9"W x 11"D) ..... \$165
- HLL Left Heel Loop (each) ..... \$20
- HLR Right Heel Loop (each) ..... \$20

**MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES**

- MFSL Left Small Footplate (4.25"W x 8"D) ..... \$225
- MFSR Right Small Footplate (4.25"W x 8"D) ..... \$225
- MFML Left Medium Footplate (5.25"W x 8"D) ..... \$225
- MFMR Right Medium Footplate (5.25"W x 8"D) ..... \$225
- MFLR Left Large Footplate (6.25"W x 8"D) ..... \$225
- MFLR Right Large Footplate (6.25"W x 8"D) ..... \$225
- MFXL Left Extra Large Footplate (9"W x 11"D) ..... \$225
- MFXR Right Extra Large Footplate (9"W x 11"D) ..... \$225
- HLL Left Heel Loop (each) ..... \$20
- HLR Right Heel Loop (each) ..... \$20

**FOOTPLATE/CUT OUTS SIZE SELECTION**



1. Certain footplate sizes may not fit all seat widths depending on legrests selected

**ADJUSTABLE ANGLE FOOTPLATES/CUT OUTS**

- AFCSL Left Small Footplate (7.25"W x 9"D) ..... \$165
- AFCSR Right Small Footplate (7.25"W x 9"D) ..... \$165
- AFCML Left Medium Footplate (7.5"W x 11"D) ..... \$165
- AFCMR Right Medium Footplate (7.5"W x 11"D) ..... \$165
- AFCLL Left Large Footplate (10"W x 14"D) ..... \$165
- AFCLR Right Large Footplate (10"W x 14"D) ..... \$165
- HLL Left Heel Loop (each) ..... \$20
- HLR Right Heel Loop (each) ..... \$20

**MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES/ CUT OUTS**

- MFCSL Left Small Footplate (7.25"W x 9"D) ..... \$225
- MFCSR Right Small Footplate (7.25"W x 9"D) ..... \$225
- MFCML Left Medium Footplate (7.5"W x 11"D) ..... \$225
- MFCMR Right Medium Footplate (7.5"W x 11"D) ..... \$225
- MFCLL Left Large Footplate (10"W x 14"D) ..... \$225
- MFCLR Right Large Footplate (10"W x 14"D) ..... \$225
- HLL Left Heel Loop (each) ..... \$20
- HLR Right Heel Loop (each) ..... \$20

**THORACIC LATERALS AND HIP/KNEE SUPPORTS**

- Matrx Laterals please see page 10
- Maxx Laterals and Hip/Knee supports please see page 11



## Maxx Rehab Seat on ROVI X3/A3 Base

Dealer Name: \_\_\_\_\_  
 Dealer Account Number: \_\_\_\_\_  
 PO #: \_\_\_\_\_  
 Tag: \_\_\_\_\_

### ADDITIONAL OPTIONS & ACCESSORIES

- NAO2H Non-Articulating O<sub>2</sub> Gas Holder<sup>(1)</sup> . . . . . \$395
- NA IPL Non-Articulating IV Pole User's Left. . . . . \$395
- NA IPR Non-Articulating IV Pole User's Right . . . . \$395
- THS Transfer Handles 4" - 6.5" . . . . . \$350
- THL Transfer Handles 9" - 11.5" . . . . . \$350
- FFLTL Flip and Fold Away 1/2 Lap Tray  
Left Mounting(14"W x 10"D)<sup>(2)(3)</sup> . . . . . \$795
- FFLTR Flip and Fold Away 1/2 Lap Tray  
Right Mounting(14"W x 10"D)<sup>(2)(3)</sup> . . . . . \$795

1. Client weight capacity is 225 lbs when using an O<sub>2</sub> holder. The maximum weight capacity of the O<sub>2</sub> holder is 20 lbs.
2. Not available on joystick side and only available with a flat arm pad or modular arm pad.
3. The maximum weight capacity of the tray is 10 lbs.

### MOTION CONCEPTS SINGLE POST BACKREST

- SPBMI Single Post Angle Adjustable Fold Down  
Back interface-Matrix Back<sup>(1)</sup> (Folds to 25° forward) **\$750**  
 Matrix Backs available on page 8
- SPBCI Single Post Angle Adjustable Fold Down  
Back for Contour Back<sup>(1)(2)</sup> (Folds to 25° forward) . . . **\$995**  
 CB1617  CB1717  CB1817  
 CB1917  CB2017
- OCUCB Onyx Fabric Upgrade for Contour Backs . . . \$100

1. Single Post Back is an angle adjustable back between 73° and 109° in 6° increments from 17° forward to 19° back. The back folds down and is only available for 16"-20" seat widths & 21"-24" back height. A 4" gap from bottom of back to seat pan is required.
2. Includes a Contour Back with fabric cover.

### ADDITIONAL ARM ACCESSORIES

- UCH Universal Cup Holder<sup>(1)</sup> . . . . . \$35
- MACESL Multi Axis Upper Extremity Support<sup>(2)</sup>(each)  
Left Side . . . . . \$225
- MACESR Multi Axis Upper Extremity Support<sup>(2)</sup> (each)  
Right Side . . . . . \$225
- OAS2L 1" Offset Armrest Spacer Left(each) . . . . . \$125
- OAS2R 1" Offset Armrest Spacer Right(each) . . . . \$125
- OAS1L 1/2" Offset Armrest Spacer Left(each) . . . \$125
- OAS1R 1/2" Offset Armrest Spacer Right(each) . . \$125
- MPAPA Multi Purpose Arm Pad  
Adapter Block<sup>(3)</sup> (each) Quantity (   ) . . . \$35

1. Not available on the joystick side.
2. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.

### INSTALLATION

- MCI Installation by Motion Concepts . . . . . N/C

**Total Retail Value \$** \_\_\_\_\_

### PLEASE SPECIFY PAGES INCLUDED

- Pages 1-7 Excludes Matrix and Maxx Seating Accessories
- Pages 1-11 Complete Order Form

### NOTES

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---



Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**CUSHION COVER REFERENCE GUIDE**

Cushion	Meshtex	Startex	Infection	Spacetex	Onyx
Libra/Adjust	N/A	STD	\$50	N/C	N/A
PS/PSP	N/A	STD -Reversible	STD -Reversible	N/A	N/A
Vi	STD	\$103	\$103	N/A	\$103

**MATRIX SEAT CUSHION OPTIONS**

Matrix Libra Adjust Seat Cushion 14" -20" wide . . . . . \$985

- LC1416-ADJ     LC1516-ADJ     LC1518-ADJ     LC1520-ADJ
- LC1616-ADJ     LC1618-ADJ     LC1620-ADJ     LC1716-ADJ
- LC1718-ADJ     LC1720-ADJ     LC1816-ADJ     LC1818-ADJ
- LC1820-ADJ     LC1916-ADJ     LC1918-ADJ     LC1920-ADJ
- LC2016-ADJ     LC2018-ADJ     LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" depth . . . . . \$1,265

- LC1622-ADJ     LC1722-ADJ     LC1822-ADJ
- LC1922-ADJ     LC2022-ADJ

Matrix Libra Adjust Seat Cushion 21" - 24" wide . . . . \$1,265

- LC2118-ADJ     LC2120-ADJ     LC2122-ADJ     LC2218-ADJ
- LC2220-ADJ     LC2222-ADJ     LC2318-ADJ     LC2320-ADJ
- LC2322-ADJ     LC2418-ADJ     LC2420-ADJ     LC2422-ADJ

Matrix Libra Seat Cushion 14" - 20" wide . . . . . \$745

- LC1416     LC1516     LC1518     LC1520     LC1616
- LC1618     LC1620     LC1716     LC1718     LC1720
- LC1816     LC1818     LC1820     LC1916     LC1918
- LC1920     LC2016     LC2018     LC2020

Matrix Libra Seat Cushion 22" depth . . . . . \$1,025

- LC1622     LC1722     LC1822     LC1922     LC2022
- LC2118

Matrix Libra Seat Cushion 21" - 24" wide . . . . . \$1,025

- LC2118     LC2120     LC2122     LC2218     LC2220
- LC2222     LC2318     LC2320     LC2322     LC2418
- LC2420     LC2422

Matrix PS Seat Cushion 14" - 20" wide . . . . . \$569

- PS1416     PS1616     PS1618     PS1620     PS1816
- PS1818     PS1820     PS2016     PS2018     PS2020

Matrix PS HD Seat Cushion 22" - 24" wide . . . . . \$975

- PS2218     PS2220     PS2222     PS2418     PS2420
- PS2422

Matrix PSVF Seat Cushion 16" - 20" wide . . . . . \$589

- PSVF1616     PSVF1618     PSVF1620
- PSVF1816     PSVF1818     PSVF1820
- PSVF2016     PSVF2018     PSVF2020

**MATRIX SEAT CUSHION OPTIONS**

Matrix PSP Seat Cushion 14" - 20" wide . . . . . \$645

- PSP1416     PSP1616     PSP1618     PSP1620
- PSP1816     PSP1818     PSP1820     PSP2016
- PSP2018     PSP2020

Matrix Vi Seat Cushion 14" - 20" wide . . . . . \$624

- MA1415-VI     MA1417-VI     MA1616-VI
- MA1618-VI     MA1620-VI     MA1816-VI     MA1818-VI
- MA1820-VI

Matrix Vi Seat Cushion 20" - 22" wide . . . . . \$681

- MA2016-VI     MA2018-VI     MA2020-VI     MA2218-VI
- MA2220-VI

Matrix Vi HD Seat Cushion 22" wide . . . . . \$840

- MA2222-VI

**FABRIC OPTIONS FOR CUSHIONS**

- SFLF    Spacetex Fabric Option for Libra Cushion . . . . . N/C
- ICFL    Infection Control Fabric Upgrade for Libra Cushion . . . . . \$50
- OCUMA    Onyx Fabric Upgrade for Matrix Vi Cushion . . . . . \$103
- SFVS    Smooth side Startex Fabric Upgrade for Matrix Vi Cushion . . . . . \$103
- SFVF    Fabric side Startex Fabric Upgrade for Matrix Vi Cushion . . . . . \$103

**BACK COVER REFERENCE GUIDE**

Back	Meshtex	Infection	Onyx
High	STD	N/A	\$103
Elite/ Elite HD/Elite Deep	STD	\$103	N/A
Elite TR/ Elite TR HD	STD	\$103	N/A
PB/PB HD/PB Deep/PB HD Deep	STD	\$103	N/A

**OTHER UPGRADABLE BACK OPTIONS**

- EBSS    Super Soft Foam for Elite Back Upgrade . . . \$103
- OCUHG    Onyx Fabric Upgrade for High Backs . . . . \$103
- WSCU    Infection Control Cover for Elite/PB Backs . . \$103  
*purchased with back*
- WSC    Infection Control Cover for Elite/PB Backs . . \$150  
*purchased with without back*
- CHIEB    Chest Harness Interface Kit for Elite Back<sup>(1)(2)</sup> . \$103
- CHIEB    Chest Harness Interface Kit for PB Back<sup>(1)(2)</sup> . . \$52

1. Does not include the chest strap.  
2. This Interface Kit is only Required for 4 point harnesses.





Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**UPGRADABLE BACK OPTIONS**

High Back 16" - 22" wide ..... \$595 <sup>(NEW)</sup>

- |                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> HB1616 | <input type="checkbox"/> HB1618 | <input type="checkbox"/> HB1620 | <input type="checkbox"/> HB1716 |
| <input type="checkbox"/> HB1718 | <input type="checkbox"/> HB1720 | <input type="checkbox"/> HB1816 | <input type="checkbox"/> HB1818 |
| <input type="checkbox"/> HB1820 | <input type="checkbox"/> HB1916 | <input type="checkbox"/> HB1918 | <input type="checkbox"/> HB1920 |
| <input type="checkbox"/> HB2016 | <input type="checkbox"/> HB2018 | <input type="checkbox"/> HB2020 | <input type="checkbox"/> HB2116 |
| <input type="checkbox"/> HB2118 | <input type="checkbox"/> HB2120 | <input type="checkbox"/> HB2216 | <input type="checkbox"/> HB2218 |
| <input type="checkbox"/> HB2220 | <input type="checkbox"/> HB2316 | <input type="checkbox"/> HB2318 | <input type="checkbox"/> HB2320 |
| <input type="checkbox"/> HB2416 | <input type="checkbox"/> HB2418 | <input type="checkbox"/> HB2420 |                                 |

Matrix Elite Back 14" - 20" wide (3" Contour) ..... \$766

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PBE1414 | <input type="checkbox"/> PBE1416 | <input type="checkbox"/> PBE1418 | <input type="checkbox"/> PBE1420 |
| <input type="checkbox"/> PBE1514 | <input type="checkbox"/> PBE1516 | <input type="checkbox"/> PBE1518 | <input type="checkbox"/> PBE1520 |
| <input type="checkbox"/> PBE1614 | <input type="checkbox"/> PBE1616 | <input type="checkbox"/> PBE1618 | <input type="checkbox"/> PBE1620 |
| <input type="checkbox"/> PBE1714 | <input type="checkbox"/> PBE1716 | <input type="checkbox"/> PBE1718 | <input type="checkbox"/> PBE1720 |
| <input type="checkbox"/> PBE1814 | <input type="checkbox"/> PBE1816 | <input type="checkbox"/> PBE1818 | <input type="checkbox"/> PBE1820 |
| <input type="checkbox"/> PBE1916 | <input type="checkbox"/> PBE1918 | <input type="checkbox"/> PBE1920 | <input type="checkbox"/> PBE2014 |
| <input type="checkbox"/> PBE2016 | <input type="checkbox"/> PBE2018 | <input type="checkbox"/> PBE2020 |                                  |

Matrix Elite HD Back 21" - 24" wide (3" Contour) . . . . \$1,095

- |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> PBE2116-HD | <input type="checkbox"/> PBE2120-HD | <input type="checkbox"/> PBE2216-HD |
| <input type="checkbox"/> PBE2220-HD | <input type="checkbox"/> PBE2316-HD | <input type="checkbox"/> PBE2320-HD |
| <input type="checkbox"/> PBE2416-HD | <input type="checkbox"/> PBE2420-HD |                                     |

Matrix Elite Deep Back 14" - 20" wide (6" Contour) . . . . \$918

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> EDB1416 | <input type="checkbox"/> EDB1418 | <input type="checkbox"/> EDB1420 | <input type="checkbox"/> EDB1514 |
| <input type="checkbox"/> EDB1516 | <input type="checkbox"/> EDB1518 | <input type="checkbox"/> EDB1520 | <input type="checkbox"/> EDB1614 |
| <input type="checkbox"/> EDB1616 | <input type="checkbox"/> EDB1618 | <input type="checkbox"/> EDB1620 | <input type="checkbox"/> EDB1714 |
| <input type="checkbox"/> EDB1716 | <input type="checkbox"/> EDB1718 | <input type="checkbox"/> EDB1720 | <input type="checkbox"/> EDB1814 |
| <input type="checkbox"/> EDB1816 | <input type="checkbox"/> EDB1818 | <input type="checkbox"/> EDB1820 | <input type="checkbox"/> EDB1914 |
| <input type="checkbox"/> EDB1916 | <input type="checkbox"/> EDB1918 | <input type="checkbox"/> EDB1920 | <input type="checkbox"/> EDB2014 |
| <input type="checkbox"/> EDB2016 | <input type="checkbox"/> EDB2018 | <input type="checkbox"/> EDB2020 |                                  |

Matrix Elite Deep HD Back 21" - 24" wide (6" Contour) . . . \$1,095

- |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> EDB2116-HD | <input type="checkbox"/> EDB2120-HD | <input type="checkbox"/> EDB2216-HD |
| <input type="checkbox"/> EDB2220-HD | <input type="checkbox"/> EDB2316-HD | <input type="checkbox"/> EDB2320-HD |
| <input type="checkbox"/> EDB2416-HD | <input type="checkbox"/> EDB2420-HD |                                     |

**UPGRADABLE BACK OPTIONS**

<sup>(NEW)</sup> Matrix Elite E2 Deep Back 14" - 20" wide/14"/16"H (6" Contour) ..... \$1,125

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2D1414 | <input type="checkbox"/> E2D1416 | <input type="checkbox"/> E2D1514 | <input type="checkbox"/> E2D1516 |
| <input type="checkbox"/> E2D1614 | <input type="checkbox"/> E2D1616 | <input type="checkbox"/> E2D1714 | <input type="checkbox"/> E2D1716 |
| <input type="checkbox"/> E2D1814 | <input type="checkbox"/> E2D1816 | <input type="checkbox"/> E2D1914 | <input type="checkbox"/> E2D1916 |
| <input type="checkbox"/> E2D2014 | <input type="checkbox"/> E2D2016 |                                  |                                  |

<sup>(NEW)</sup> Matrix Elite E2 Deep Back 14" - 20" wide/18"/20"H (6" Contour) ..... \$1,155

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2D1418 | <input type="checkbox"/> E2D1420 | <input type="checkbox"/> E2D1518 | <input type="checkbox"/> E2D1520 |
| <input type="checkbox"/> E2D1618 | <input type="checkbox"/> E2D1620 | <input type="checkbox"/> E2D1718 | <input type="checkbox"/> E2D1720 |
| <input type="checkbox"/> E2D1818 | <input type="checkbox"/> E2D1820 | <input type="checkbox"/> E2D1918 | <input type="checkbox"/> E2D1920 |
| <input type="checkbox"/> E2D2018 | <input type="checkbox"/> E2D2020 |                                  |                                  |

Matrix Elite TR Back 15" - 20" wide (3" Contour) . . . . \$866

- |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> PBE1516-TR | <input type="checkbox"/> PBE1518-TR | <input type="checkbox"/> PBE1520-TR |
| <input type="checkbox"/> PBE1616-TR | <input type="checkbox"/> PBE1618-TR | <input type="checkbox"/> PBE1620-TR |
| <input type="checkbox"/> PBE1716-TR | <input type="checkbox"/> PBE1718-TR | <input type="checkbox"/> PBE1720-TR |
| <input type="checkbox"/> PBE1816-TR | <input type="checkbox"/> PBE1818-TR | <input type="checkbox"/> PBE1820-TR |
| <input type="checkbox"/> PBE1916-TR | <input type="checkbox"/> PBE1918-TR | <input type="checkbox"/> PBE1920-TR |
| <input type="checkbox"/> PBE2016-TR | <input type="checkbox"/> PBE2018-TR | <input type="checkbox"/> PBE2020-TR |

Matrix Elite TR HD Back 21" - 24" wide (3" Contour) . . \$1,151

- |                                       |                                       |                                       |
|---------------------------------------|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> PBE2120-TRHD | <input type="checkbox"/> PBE2220-TRHD | <input type="checkbox"/> PBE2320-TRHD |
| <input type="checkbox"/> PBE2420-TRHD |                                       |                                       |

Matrix PB Posture Back 14" - 20" wide (5" Contour) . . . . \$734

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> MPB1412 | <input type="checkbox"/> MPB1612 | <input type="checkbox"/> MPB1616 | <input type="checkbox"/> MPB1620 |
| <input type="checkbox"/> MPB1712 | <input type="checkbox"/> MPB1716 | <input type="checkbox"/> MPB1720 | <input type="checkbox"/> MPB1812 |
| <input type="checkbox"/> MPB1816 | <input type="checkbox"/> MPB1820 | <input type="checkbox"/> MPB2012 | <input type="checkbox"/> MPB2016 |
| <input type="checkbox"/> MPB2020 |                                  |                                  |                                  |

Matrix PB Posture Back Deep 14" - 20" wide (7" Contour) . \$892

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PBD1412 | <input type="checkbox"/> PBD1612 | <input type="checkbox"/> PBD1616 | <input type="checkbox"/> PBD1620 |
| <input type="checkbox"/> PBD1712 | <input type="checkbox"/> PBD1716 | <input type="checkbox"/> PBD1720 | <input type="checkbox"/> PBD1812 |
| <input type="checkbox"/> PBD1816 | <input type="checkbox"/> PBD1820 | <input type="checkbox"/> PBD2012 | <input type="checkbox"/> PBD2016 |
| <input type="checkbox"/> PBD2020 |                                  |                                  |                                  |

Matrix PB Heavy Duty Back 20" - 22" wide (5" Contour) . \$1,051

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> HDB2016 | <input type="checkbox"/> HDB2020 | <input type="checkbox"/> HDB2216 | <input type="checkbox"/> HDB2220 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrix PB Heavy Duty Deep Back 20"-22" wide (5" Contour) \$1,209

- |                                   |                                   |                                   |                                   |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> DHDB2016 | <input type="checkbox"/> DHDB2020 | <input type="checkbox"/> DHDB2216 | <input type="checkbox"/> DHDB2220 |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|

<sup>(NEW)</sup> Matrix Elite E2 Back 14" - 20" wide/14"/16"H(3" Contour) . \$965

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2S1414 | <input type="checkbox"/> E2S1416 | <input type="checkbox"/> E2S1514 | <input type="checkbox"/> E2S1516 |
| <input type="checkbox"/> E2S1614 | <input type="checkbox"/> E2S1616 | <input type="checkbox"/> E2S1714 | <input type="checkbox"/> E2S1716 |
| <input type="checkbox"/> E2S1814 | <input type="checkbox"/> E2S1816 | <input type="checkbox"/> E2S1914 | <input type="checkbox"/> E2S1916 |
| <input type="checkbox"/> E2S2014 | <input type="checkbox"/> E2S2016 |                                  |                                  |

<sup>(NEW)</sup> Matrix Elite E2 Back 14" - 20" wide/18"/20"H(3" Contour) ..... \$995

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2S1418 | <input type="checkbox"/> E2S1420 | <input type="checkbox"/> E2S1518 | <input type="checkbox"/> E2S1520 |
| <input type="checkbox"/> E2S1618 | <input type="checkbox"/> E2S1620 | <input type="checkbox"/> E2S1718 | <input type="checkbox"/> E2S1720 |
| <input type="checkbox"/> E2S1818 | <input type="checkbox"/> E2S1820 | <input type="checkbox"/> E2S1918 | <input type="checkbox"/> E2S1920 |
| <input type="checkbox"/> E2S2018 | <input type="checkbox"/> E2S2020 |                                  |                                  |

**Maxx Rehab Seat on ROVI X3/A3 Base**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRX ELAN HEADREST PAD<sup>(1)(2)</sup>**

<input type="checkbox"/> ESP6	Matrx Elan Standard Pad 6" .....	\$310
<input type="checkbox"/> ESP10	Matrx Elan Standard Pad 10" .....	\$310
<input type="checkbox"/> ESP14	Matrx Elan Standard Pad 14" .....	\$310
<input type="checkbox"/> EOP9	Matrx Elan Occipital Pad 9" .....	\$341
<input type="checkbox"/> EOP12	Matrx Elan Occipital Pad 12" .....	\$341
<input type="checkbox"/> E4PTSM	Matrx Elan 4-Point Pad 9"W x 7"H .....	\$448
<input type="checkbox"/> E4POINT	Matrx Elan 4-Point Pad 11"W x 10"H .....	\$448
<input type="checkbox"/> ESP6-IC	Matrx Elan Standard Pad 6" with Infection Control .....	\$335
<input type="checkbox"/> ESP10-IC	Matrx Elan Standard Pad 10" with Infection Control .....	\$335
<input type="checkbox"/> ESP14-IC	Matrx Elan Standard Pad 14" with Infection Control .....	\$335
<input type="checkbox"/> EOP9-IC	Matrx Elan Occipital Pad 9" with Infection Control .....	\$366
<input type="checkbox"/> EOP12-IC	Matrx Elan Occipital Pad 12" with Infection Control .....	\$366
<input type="checkbox"/> E4PTSM-IC	Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control .....	\$473
<input type="checkbox"/> E4POINT-IC	Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control .....	\$473

1. Not compatible with sling back (ASBA-2)  
 2. If mounting on PB Elite must order removable hardware

**HEADREST MOUNTING HARDWARE OPTIONS**

<input type="checkbox"/> MEHWRE	Elan Multi-Axis Removable HDW.....	\$285
---------------------------------	------------------------------------	-------

1. Only compatible w/Standard Rehab, Contour & High Back Models. All other head rests REQUIRE an upgrade to Removable Hardware.

**MATRX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS**

<input type="checkbox"/> E2FLL	Matrx E2 Fixed Laterals-Left(6.5x5") <sup>(1)</sup> .....	\$195
<input type="checkbox"/> E2FLR	Matrx E2 Fixed Laterals-Right(6.5x5") <sup>(1)</sup> .....	\$195
<input type="checkbox"/> PBSSL	Matrx Small Lateral Fixed-Left(3x5.5") .....	\$113
<input type="checkbox"/> PBSSR	Matrx Small Lateral Fixed-Right(3x5.5") . . . . .	\$113
<input type="checkbox"/> PBLSL	Matrx Standard Lateral Fixed-Left(7.5x5") . . . . .	\$116
<input type="checkbox"/> PBLSR	Matrx Standard Lateral Fixed-Right (7.5x5")	\$116
<input type="checkbox"/> PBOLL	Matrx Offset Lateral Fixed Left . . . . .	\$125
<input type="checkbox"/> PBOLR	Matrx Offset Lateral Fixed Right . . . . .	\$125

1. Only available on the E2 backrests.

**MATRX SWING AWAY LATERALS - FOR ELITE AND ELITE TR BACKS**

Matrx Lateral Trunk Support Left <sup>(1)</sup>		
<input type="checkbox"/> EPADXL	Size XS (3.25x4.25") .....	\$115
<input type="checkbox"/> EPADSL	Size S (3.5x5.25") .....	\$115
<input type="checkbox"/> EPADM L	Size M (3.75x5.75") .....	\$115
<input type="checkbox"/> EPADLL	Size L (4.25x6.5") .....	\$115

Matrx Lateral Trunk Support Right <sup>(1)</sup>		
<input type="checkbox"/> EPADXR	Size XS (3.25x4.25") .....	\$115
<input type="checkbox"/> EPADSR	Size S (3.5x5.25") .....	\$115
<input type="checkbox"/> EPADMR	Size M (3.75x5.75") .....	\$115
<input type="checkbox"/> EPADLR	Size L (4.25x6.5") .....	\$115

1. Must choose swing away hardware below

**MATRX SWING AWAY HARDWARE - FOR ELITE AND ELITE TR BACKS**

<input type="checkbox"/> SAEHBL	Matrx Elite Swing Away Hardware Left . . . . .	\$250
<input type="checkbox"/> SAHEBR	Matrx Elite Swing Away Hardware Right . . . . .	\$250
<input type="checkbox"/> OSAHL	Matrx Offset Elite Swing Away Hardware Left . . . . .	\$250
<input type="checkbox"/> OSAHR	Matrx Offset Elite Swing Away HardwareRight . . . . .	\$250

1. Recomendend when finished seat cushion thickness will be greater then 3.5" high.

**Maxx Rehab Seat on ROVI X3/A3 Base**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MAXX LATERALS TRUNK SUPPORTS<sup>(1)</sup>**

**Left Fixed Mount Lateral Trunk Support<sup>(1)</sup>**

Select Left Trunk Support Size

- LATXL Size XS (3.25" x 4.25") . . . . . \$100
- LATSL Size S (3.5" x 5.25") . . . . . \$100
- LATML Size M (3.75" x 5.75") . . . . . \$100
- LATLL Size L (4.25" x 6.5") . . . . . \$100

**Right Fixed Mount Lateral Trunk Support<sup>(1)</sup>**

Select Right Trunk Support Size

- LATXR Size XS (3.25" x 4.25") . . . . . \$100
- LATSR Size S (3.5" x 5.25") . . . . . \$100
- LATMR Size M (3.75" x 5.75") . . . . . \$100
- LATLR Size L (4.25" x 6.5") . . . . . \$100

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs.

**MAXX STYLE SWING-AWAY TRUNK SUPPORT HARDWARE<sup>(1)</sup>**

- SALATL Left MAXX Style Swing Away HDW . . . . . \$195

Select Left Trunk Support Size

- LATXL Size XS (3.25" x 4.25") . . . . . \$100
- LATSL Size S (3.5" x 5.25") . . . . . \$100
- LATML Size M (3.75" x 5.75") . . . . . \$100
- LATLL Size L (4.25" x 6.5") . . . . . \$100

- SALATR Right MAXX Style Swing Away HDW<sup>(1)</sup> . . . \$195

Select Right Trunk Support Size

- LATXR Size XS (3.25" x 4.25") . . . . . \$100
- LATSR Size S (3.5" x 5.25") . . . . . \$100
- LATMR Size M (3.75" x 5.75") . . . . . \$100
- LATLR Size L (4.25" x 6.5") . . . . . \$100

- ATLL Additional Telescoping Link - Left<sup>(2)</sup> . . . \$75

- ATLR Additional Telescoping Link - Right<sup>(2)</sup> . . \$75

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs  
2. Up to 2.75" medial offset

**FIXED HIP/ KNEE SUPPORTS**

**Lateral Hip Support with Fixed Mounting Hardware Left<sup>(1)</sup>**

Select Pad Size and mounting position at Hip:

- HSPXLH Size XS (5.5" W x 4" H) mounted at Hip . \$120
- HSPSLH Size S (7" W x 6" H) mounted at Hip . . . \$120
- HSPMLH Size M (9" W x 4" H) mounted at Hip . . . \$120
- HSPLLH Size L (13" W x 4" H) mounted at Hip . . . \$120

Select Pad Size and mounting position at Knee:

- HSPXLK Size XS (5.5" W x 4" H) mounted at Knee . . . \$120
- HSPSLK Size S (7" W x 6" H) mounted at Knee . . . . \$120
- HSPMLK Size M (9" W x 4" H) mounted at Knee . . . . \$120
- HSPLLK Size L (13" W x 4" H) mounted at Knee . . . . \$120

**FIXED HIP/ KNEE SUPPORTS - CONTINUED**

**Lateral Hip Support with Fixed Mounting Hardware Right<sup>(1)</sup>**

Select Pad Size and mounting position:

- HSPXRH Size XS (5.5" W x 4" H) mounted at Hip . . . . \$120
- HSPSRH Size S (7" W x 6" H) mounted at Hip . . . . . \$120
- HSPMRH Size M (9" W x 4" H) mounted at Hip . . . . . \$120
- HSPLRH Size L (13" W x 4" H) mounted at Hip . . . . . \$120

Select Pad Size and mounting position at Knee:

- HSPXRK Size XS (5.5" W x 4" H) mounted at Knee . . . \$120
- HSPSRK Size S (7" W x 6" H) mounted at Knee . . . . \$120
- HSPMRK Size M (9" W x 4" H) mounted at Knee . . . . \$120
- HSPLRK Size L (13" W x 4" H) mounted at Knee . . . . \$120

1. Must choose swing away hardware below if required.

**REMOVABLE HIP/ KNEE SUPPORT HARDWARE**

**Lift Off Removable Hardware for Hip Support**

- LOHL Left . . . . . \$195
- LOHR Right . . . . . \$195

**MAXX Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports**

Select mounting position:

- MHSHLH Left mounted at Hip . . . . . \$215
- MHSHLK Left mounted at Knee . . . . . \$215
- MHSHRH Right mounted at Hip . . . . . \$215
- MHSHRK Right mounted at Knee . . . . . \$215

**Swing Away Removable Hardware for Hip Supports**

Select mounting position:

- SHSHLH Left mounted at Hip . . . . . \$215
- SHSHLK Left mounted at Knee . . . . . \$215
- SHSHRH Right mounted at Hip . . . . . \$215
- SHSHRK Right mounted at Knee . . . . . \$215

**Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports<sup>(1)</sup>**

Select mounting position:

- EMQHLH Left mounted at Hip . . . . . \$215
- EMQHLLK Left mounted at Knee . . . . . \$215
- EMQHRH Right mounted at Hip . . . . . \$215
- EMQHRK Right mounted at Knee . . . . . \$215

1. Recommended when finished seat cushion thickness will be greater than 3.5" high.