

REQUIRED INFORMATION


Quote Order

Date of Order: _____ Dealer Account #: _____
 Dealer Name: _____
 PO #: _____ Tag _____
 Purchasing Contact: _____
 Phone: _____ Fax: _____
 E-mail: _____
 RTS/Therapist: _____
 Ship to Address: _____
 City: _____ Province: _____
 Postal Code: _____ Client Gender: M F
 Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



Customer Service: 1.866.748.7943 | Fax: 1.905.695.0583 | orderscdn@motionconcepts.com

**ROVI
BASE STATISTICS**

SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH A ₃	25.25"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"
CASTER WHEEL SIZE (front/rear)	6"
TURNING RADIUS A ₃	20.5"
DRIVETRAIN	4-pole EAD Motors

BATTERIES	2 x M34
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs (M34) / 17.5 kg
RANGE ⁽²⁾	17.7 miles / 28.5 km/h
WEIGHT CAPACITY	300 lbs / 136 kg
BASE WEIGHT A ₃ (w/o batteries)	159 lbs / 72 kg
GROUND CLEARANCE	3"
MAXIMUM SPEED	5.8 mph

A₃



*Shown with Motion Concepts Single Post Back for Contour Backrest. See page 7

⁽¹⁾ All X3 and A3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

BASE OPTIONS

888WNLEA3 ROVI A3⁽¹⁾⁽²⁾ \$6,800

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt. Weight capacity is 300 lbs.
2. When selecting an increased speed motor package there is a reduction in torque output by roughly 15%.

FENDER COLOUR

- 800914-00200-03 Blizzard White
- 800914-00200-69 Solar Yellow
- 800914-00200-23 Original Orange
- 800914-00201-21 Red Delicious
- 800914-00201-2C Party Pink
- BCF Black Carbon (Mat) Black Panther (Gloss)
- BPF Purple Reign
- PRF Indigo Blue
- IBF Motion Green
- MGF Motion Green

WHEELCHAIR TRANSPORT BRACKETS

XBKTS Wheelchair Transport Brackets(1) STD
1. For unoccupied use only

TIRE OPTIONS

STIRES Standard Solid Tires STD
 PTIRES Pneumatic Tires N/C

BATTERIES

M34 SLD G MK Group Battery x 2 \$650

JOYSTICK SELECTION

D51122 R-Net LED Joystick N/C

JOYSTICK MOUNT

800914-001-16-ML Fixed Joystick Mounting Left N/C
 800914-001-16-MR Fixed Joystick Mounting Right N/C
 SAQL Height Adjustable Swing-Away Quad Link Left \$244
 SAQR Height Adjustable Swing-Away Quad Link Right \$244

ADDITIONAL ACCESSORIES

MUSB Motion Dual USB Charger N/C
 UCH Universal Cup Holder⁽¹⁾ \$50

1. Not available on the joystick side.

SEAT-TO-FLOOR HEIGHT⁽¹⁾

STF17 17" Seat-To-Floor N/C
 STF18 18" Seat-To-Floor N/C
 STF19 19" Seat-To-Floor N/C
 STF20 20" Seat-To-Floor N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"

REHAB SEAT

RSMAXX Ultra Low Maxx Rehab Seat \$995

ADD: MODULES

PRETILT 5° Fixed Posterior Tilt⁽¹⁾ N/C

1. This Fixed posterior tilt option adds approximately 2" to the front seat to floor height

ADULT SEAT WIDTH SETTINGS

SW116 16" Wide (Adjustable from 16" to 20") N/C
 SW117 17" Wide (Adjustable from 16" to 20") N/C
 SW118 18" Wide (Adjustable from 16" to 20") N/C
 SW119 19" Wide (Adjustable from 16" to 20") N/C
 SW120 20" Wide (Adjustable from 16" to 20") N/C

1. This option is only available on the Rovi A3 Base.

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

SD116 16" Deep (Adjustable from 16" to 21") N/C
 SD117 17" Deep (Adjustable from 16" to 21") N/C
 SD118 18" Deep (Adjustable from 16" to 21") N/C
 SD119 19" Deep (Adjustable from 16" to 21") N/C
 SD120 20" Deep (Adjustable from 16" to 21") N/C

LAP BELTS

BELT48 Push Button Style Lap Belt, Length 48" STD
 BELT60 Push Button Style Lap Belt, Length 60" N/C
 PLB2 Bodypoint Padded Lap Belt 2 Point⁽¹⁾ \$250
 PLB4 Bodypoint Padded Lap Belt 4 Point⁽¹⁾ \$315

1. 67" length only option

MATRIX SEAT CUSHION OPTIONS⁽¹⁾

PSWWDD Matrix PS Seat Cushion \$569
 PSPWWDD Matrix PSP Seat Cushion \$645

1. Cushion size will match seat width and depth

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FINISHED BACK HEIGHT ⁽¹⁾⁽²⁾

- 18" FBT18
- 19" FBT19
- 20" FBT20
- 21" FBT21
- 22" FBT22
- 23" FBT23
- 24" FBT24
- 25" FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK CANES

- SBC18 18" Straight Back Canes STD
- SBC20 20" Straight Back Canes STD
- SBC22 22" Straight Back Canes STD
- SBC24 24" Straight Back Canes STD

Select Cane Angle

- 75° BA75
- 79° BA79
- 85° BA85
- 90° BA90
- 95° BA95
- 101° BA101
- 106° BA106
- 112° BA112
- 116° BA116

BACK OPTIONS

- SRB Standard Rehab Back & Cushion \$125
- ORBC Omit Rigid Back & Cushion N/C

UPGRADABLE BACK OPTIONS⁽¹⁾

- E2SWW14 Matrix Elite E2 (3" Contour) 14" H \$965
- E2SWW16 Matrix Elite E2 (3" Contour) 16" H \$965
- E2SWW18 Matrix Elite E2 (3" Contour) 18" H \$995
- E2SWW20 Matrix Elite E2 (3" Contour) 20" H \$995

- E2DWW14 Matrix Elite E2D (6" Contour) 14" H \$1,125
- E2DWW16 Matrix Elite E2D (6" Contour) 16" H \$1,125
- E2DWW18 Matrix Elite E2D (6" Contour) 18" H \$1,155
- E2DWW20 Matrix Elite E2D (6" Contour) 20" H \$1,155

1. Back width will match seat width.

ARMRESTS

- DP1 Adult Dual Post 9.75" - 12.75" STD
- DP2 Adult Dual Post 12.75" - 15.75" STD
- DP3 Adult Dual Post 15.75" - 18.75" STD

- RMCAN Ultra Rail Seat Mnt Flip Back Cantilever
MAXX Tilt Arm 9" - 12.5"⁽¹⁾ \$488

1. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.

ARMREST TUBE LENGTH - PLEASE SPECIFY

- ATFL Left Full STD
- ATDL Left Desk STD
- ATFR Right Full STD
- ATDR Right Desk STD

ARMPADS

Standard

- SAPFL Left Full (2.25" x 14") STD
- SAPDL Left Desk (2.25" x 10") STD
- SAPFR Right Full (2.25" x 14") STD
- SAPDR Right Desk (2.25" x 10") STD

Waterfall Arm Pads⁽¹⁾

- WAPFL Left Full (13") N/C
- WAPDL Left Desk (9") N/C
- WAPFR Right Full (13") N/C
- WAPDR Right Desk (9") N/C

Gel Pads

- GAPNL Narrow Left (2" x 12") \$88
- GAPNR Narrow Right (2" x 12") \$88
- GAPWL Wide Left (3.5" x 12") \$96
- GAPWR Wide Right (3.5" x 12") \$96
- GAPLL Long Left (3.5" x 14") \$103
- GAPLR Long Right (3.5" x 14") \$103


1. Full length pad requires full length arm tube if using Ultra Rail mounted Tilt Arms.

MANUAL CENTER MOUNT FRONT RIGGING

- MSFCMFP Maxx Style Fixed Center Mount
Foot Platform⁽¹⁾⁽²⁾ \$275


 - IFP Medium Rubber Coated Footplate
(11.5" W x 10") N/C

 - ICP Individual calf pads \$100
- Select seat pan to footplate:
- SPTF1 9" - 13" N/C
 - SPTF14 13" - 19" N/C
- Set at:
- SA70 Set at 70° N/C
 - SA90 Set at 90° N/C
 - SA97 Set at 97° N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
2. 3" ground clearance required.

FOOT PLATFORM OPTIONS FOR MSFCMFP⁽¹⁾



- SFFP Small Platform 12"w x 8.5"d \$180
- MFFP Medium Platform 12"w x 10.75"d \$180
- LFFP Large Platform 14.5"w x 13"d \$250
- XLFFP Extra Large Platform 17"w x 13"d \$250
- XXLFFP Extra Extra Large Platform 21"w x 13"d \$250

1. May need to adjust leg position to avoid interference with front castors

MANUAL ELEVATING CENTER MOUNT LEGREST

- MECML Manual Elevating Center Mount Legrest With Plantar Flex. \$595
 - Select legrest set up position:
 - 90SU Set at 90° (elevating range 90-40°) N/C
 - 110TU Tucked at 110° (elevating range 110-60°)⁽¹⁾⁽²⁾ N/C
 - Rubber Coated Footplate:
 - MRCP Rubber Coated Footplate (11.5" W x 10") N/C
 - Select seat pan to footplate for Rubber Coated Footplate
 - CSPTF1 11.5" - 15" N/C
 - CSPTF2 14" - 19" N/C
 - Optional:
 - CICP Individual calf pads \$200

1. Not available on seat depths 16", 17".

SWING AWAY LEGREST

- HD70 Heavy Duty 70° Swing Away Footrest⁽¹⁾⁽²⁾ \$495
 - Select Receivers:
 - RS Straight Receivers. N/C
 - RF Flared 1" Outward Receivers N/C
 - Select Seat Pan to Footplate:
 - SPTF11 Length: 6"-10" N/C
 - SPTF8 Length: 11"-15" N/C
 - SPTF2 Length: 13"-17" N/C
 - SPTF7 Length: 16"-20" N/C
 - Options:
 - HDFCP Flat Calf Pad Option. \$275
 - HD70CCP Curved Calf Pad Option⁽³⁾. \$275
- I70SAF Invacare 70° Swing-away Footrests \$325
 - I70-CF Comes with composite foot plates. N/C
- 70HDAIH Invacare 70° Non-Tapered HD Swing-away Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards \$457

- 1. Must choose footplate on page 6
- 2. Heavy duty durable construction (4.5lbs ea.)
- 3. For 16" and 17" seat widths you must select flared receivers.

SUBMIT

PRINT FORM

SAVE

CLEAR FORM

FOOTPLATE FOR SWING AWAY/ELEVATING LEGS

Adjustable Angle Footplates

- AFSL Left Small Footplate (4.25"W x 8"D) \$165
- AFSR Right Small Footplate (4.25"W x 8"D) \$165
- AFML Left Medium Footplate (5.25"W x 8"D) \$165
- AFMR Right Medium Footplate (5.25"W x 8"D) \$165
- AFLL Left Large Footplate (6.25"W x 8"D) \$165
- AFLR Right Large Footplate (6.25"W x 8"D) \$165
- AFXL Left Extra Large Footplate (9"W x 11"D) \$165
- AFXR Right Extra Large Footplate (9"W x 11"D) \$165

Multi-Axis Adjustable Angle Footplates

- MFSL Left Small Footplate (4.25"W x 8"D) \$225
- MFSR Right Small Footplate (4.25"W x 8"D) \$225
- MFML Left Medium Footplate (5.25"W x 8"D) \$225
- MFMR Right Medium Footplate (5.25"W x 8"D) \$225
- MFLL Left Large Footplate (6.25"W x 8"D) \$225
- MFLR Right Large Footplate (6.25"W x 8"D) \$225
- MFXL Left Extra Large Footplate (9"W x 11"D) \$225
- MFXR Right Extra Large Footplate (9"W x 11"D) \$225

HEEL LOOPS, CUPS AND TOE STRAPS

- HLL Left Heel Loop (each) \$20
- HLR Right Heel Loop (each) \$20
- HCL Left Cup Loop (each) \$20
- HCR Right Cup Loop (each) \$20
- TSL Left Toe Strap (each) \$20
- TSR Right Toe Strap (each) \$20

FIXED ADUCTORS

Left Support

- HSPXLH Size XS (5.5" W x 4" H) \$120
- HSPSMLH Size S/M (7" W x 4" H) \$120
- HSPSLH Size S (7"W x 6"H) \$120
- HSPMLH Size M (9" W x 4"H) \$120
- HSPLLH Size L (13" W x 4"H) \$120

Right Support

- HSPXRH Size XS (5.5"W x 4"H) \$120
- HSPSMRH Size S/M (7" W x 4" H) \$120
- HSPSRH Size S (7"W x 6"H) \$120
- HSPMRH Size M (9"W x 4"H) \$120
- HSPLRH Size L (13"W x 4"H) \$120

INSTALLATION

- MCI Motion Concepts Installation. STD

Total Retail Value \$ _____