



Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

SAVE

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

SUBMIT

CLEAR FORM

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

905-695-0583

PRINT FORM



TDX-SP2 BASE WITH ULTRA LOW MAXX POWER POSITIONING SYSTEM

CDN PRICE LIST AND ORDER FORM

Price Effective November 26, 2024



ULTRA^{Low}
maxx

PRINT FORM

SAVE

CLEAR FORM

REQUIRED INFORMATION

Quote ☐ Order ☐

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: ☐M ☐F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured
please fill in all **REQUIRED*** measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness	



Customer Service: 1.800.748.7943 | Fax: 1.905.695.0583 | www.motionconcepts.com

INVACARE POWER BASES

- ☐ TDXSP2 TDX-SP2 Base \$7,575
- Select STF(1): ☐ SLOW 16.75" (2) \$250
- ☐ LOW 17.25"
- ☐ MED 18.25"
- ☐ TALL 19.25"
- ☐ 22TRY Narrow Base 24" N/C
- ☐ 24TRY Standard Base 25.5" N/C
- ☐ U300 Weight Capacity up to 300 lbs STD

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size.

2. Must order 24TRY battery tray type and Transport Brackets TRBKTS.
Only available on 17" seat depths and greater.

TRANSPORT TIE DOWN

- ☐ TRBKTS Wheelchair Transport Brackets(1) STD
- ☐ TRRO Transport Ready Option(2)(4) \$275
- ☐ TRRO-E Transport Ready Option for Elevate(2)(3)(4) . . . \$275

1. For unoccupied use only

2. TRRO meets ISO-7176 Part 19

3. Required for Elevate Systems to meet ISO-7176 Part 19, includes two additional tie-down loops on rail.

4. Must select a Motion Concepts or Matrix Elan headrest with hardware and Motion Concepts Standard Rehab, High Back or Matrix Upgradeable back. Comes with pelvic safety belt standard, please select TROBELT on page 5.

TIRE OPTIONS

- ☐ B1431-3 14" x 3" Black Tire with Gel Foam Inserts . . STD
- ☐ B1430-3T 14" x 3" Black Tire - Pneumatic Tire N/C

FORK OPTIONS

- ☐ DSFK Double-Sided Fork Package STD
- ☐ FKPKG Single-Sided Fork Package \$375

INSTALLATION OF BATTERIES

- ☐ 22NFBATTERY 22 NF Battery (On Chair) \$700
- ☐ 24BATTERY 24 Gel Battery (On Chair) \$840

BATTERY CHARGER OPTION

- ☐ 110CHARGER 110 Volt Battery Charger N/C
- ☐ L8900 Omit Battery Charger N/C

FRAME FINISH

- ☐ 170P Obsidian Black N/C
- ☐ 172P Peacock Blue N/C
- ☐ 177P Artic Blue N/C
- ☐ 176P Lime Twist N/C
- ☐ 171P Ruby Red N/C
- ☐ 174P Purple Passion N/C
- ☐ 175P Denim Navy N/C

RIM INSERT COLORS - MUST PICK ONE

- ☐ 170PR Obsidian Black N/C
- ☐ 175PR Denim Navy N/C
- ☐ 172PR Peacock Blue N/C
- ☐ 177PR Artic Blue N/C
- ☐ 176PR Lime Twist N/C
- ☐ 178PR Banana Split N/C
- ☐ 179PR Orange Marmalade N/C
- ☐ 171PR Ruby Red N/C
- ☐ 173PR Pink Flamingo N/C
- ☐ 174PR Purple Passion N/C
- ☐ 180PR Morning Mist N/C

Ultra Low Maxx Power Positioning
on Invacare TDX-SP2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

NON-EXPANDABLE SINGLE/DUAL ACTUATOR JOYSTICKS

- ☐ REM110 LED Drive Only Remote/Joystick⁽¹⁾ N/C
Single drive profile
- ☐ REM210 LED Non-Expandable Remote/Joystick⁽¹⁾ N/C
Three drive profiles, control of two seat functions

1. Not available with lights and indicators.
2. Must choose lighting package separately

EXPANDABLE JOYSTICKS

- ☐ REM211 LED Expandable Remote/Joystick⁽¹⁾ \$875
Three drive profiles, control of up to 6 actuators.
Not for use with specialty controls.
- ☐ REM216 LED Expandable Remote/Joystick
with Lights⁽²⁾ \$875
Three drive profiles, control of up to 6 actuators
and light operation. Not for use with specialty controls.
- ☐ REM400 Color 3.5" Touch Screen
Remote/Joystick⁽³⁾ \$1,100
Four drive profiles, control of up to 6 actuators through the
drive control and light operation.
- ☐ SW400TGL Toggle Kit for REM400⁽⁴⁾ \$200

1. Not available with lights and indicators
2. Must choose lighting package separately.
3. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.
4. Adds 2 toggle switches to REM400.

EXPANDABLE SPECIALTY DISPLAY

- ☐ REM500 LiNX Color 3.5" Touch Screen
Display⁽¹⁾⁽²⁾⁽³⁾ \$1,100
- Select mounting position for REM500
- ☐ R5MR Right N/C
- ☐ R5ML Left N/C
- ☐ R5RPB Rear Push Bar N/C

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Driver Control not included.
3. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.

PERFORMANCE MOTORS AND ENHANCED DRIVING

- ☐ HSMTR Performance Motors - 7.5 mph⁽¹⁾ \$995
- ☐ GTRAC-LX LiNX G-TRAC Module⁽²⁾ \$1,200

1. Must order 24TRY Battery Tray Type
2. Recommended for use with Specialty Drive Controls

EXPANDABLE MODULES

- ☐ OUT500 Output Module \$645
- ☐ IN200 LiNX Switch Input Module for
Power ON/OFF for REM200 Series⁽¹⁾ \$645
- ☐ IN500 LiNX 9 Pin Input Module
with Sip-N-Puff Nozzle⁽²⁾⁽³⁾⁽⁴⁾ \$1,800
- ☐ PKG32666 Therafin Sip-N-Puff Breath Tube Kit. \$440

1. Not for use with Alternative Drive Controls.
2. Only one input module per chair, requires REM400 or REM500.
3. Input module also required for any 9-pin specialty control.
4. 9 pin specialty control and Sip-N-Puff cannot both be used together.

SPECIALTY PROPORTIONAL CONTROLS⁽¹⁾⁽²⁾

- ☐ CREM LiNX Compact Remote \$1,097
- ☐ CREM-LF LiNX Compact Remote Low Force \$1,097

1. Requires REM400 or REM500
2. Select Mount in Mounting For Driver Control section

ATTENDANT DRIVE CONTROL

- ☐ ACU Proportional Attendant Control⁽¹⁾ \$1,024
- Select mounting position for ACU
- ☐ RPBMR Right
- ☐ RPBML Left
- ☐ RPBMC Center

1. Requires Expandable Controller.

LIGHTS AND INDICATORS

- ☐ LIGHTS LED Lights and Indicators⁽¹⁾ \$995

1. Must order REM216, REM400 or REM500

LiNX ELECTRONIC ACCESSORIES

- ☐ LAK LiNX Access Key⁽¹⁾ \$50
- ☐ MUSB Motion Dual USB Charger. N/C

1. Required for programming

MOUNTING FOR DRIVER CONTROLS

Arm Mounted

- ☐ FR Fixed Joystick Mounting - Right STD
- ☐ FL Fixed Joystick Mounting - Left STD
- ☐ MRR Maxx Resolve Swing-Away Joystick Mount
(Right Hand) \$450
- ☐ MRL Maxx Resolve Swing-Away Joystick Mount
(Left Hand) \$450
- ☐ SAQR Motion Height Adjustable Swing-Away
Quad Link (Right Hand) \$244
- ☐ SAQL Motion Height Adjustable Swing-Away
Quad Link (Left Hand). \$244

Rear Push Bar Mounted

- ☐ RPBMR Right \$140
- ☐ RPBML Left \$140
- ☐ RPBMC Center \$140

JOYSTICK TOPS

- ☐ PC101A Bodypoint U Shaped Handle 3"⁽¹⁾ \$120
- ☐ PC102A Bodypoint U Shaped Handle 4"⁽¹⁾ \$120
- ☐ PC104A BodyPoint Straight Handle \$120
- ☐ PC107A Bodypoint Rubber Dome⁽¹⁾. \$120
- ☐ PC110A Grooved Mushroom Joystick Handle
3/16" Stem⁽¹⁾ \$120

1. To remove the joystick knob from the REM400, pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

Ultra Low Maxx Power Positioning
on Invacare TDX-SP2 Base

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

POWER POSITIONING SYSTEMS

☐ MDCGT-50	Ultra Low Maxx CG Tilt Only 50°. . . . \$6,969
☐ MDCGT-45	Ultra Low Maxx CG Tilt Only 45° ⁽¹⁾ . . \$6,969
☐ MDCGTRESR-50	Ultra Low Maxx CG Tilt 50° & Power Recline 170° with Extended Shear Reduction. . . . \$11,569
☐ MDCGTRESR-45	Ultra Low Maxx CG Tilt 45° ⁽¹⁾ & Power Recline 170° with Extended Shear Reduction. . . . \$11,569
☐ MDCGTRPSB-50	Ultra Low Maxx CG Tilt 50° & Power Recline 170° with Power Sliding Back ⁽²⁾ \$12,069
☐ MDCGTRPSB-45	Ultra Low Maxx CG Tilt 45° ⁽¹⁾ & Power Recline 170° with Power Sliding Back ⁽²⁾ \$12,069
☐ MDRESR	Ultra Low Maxx Power Recline 170° with Extended Shear Reduction. . . . \$10,350
☐ MDRPSB	Ultra Low Maxx Power Recline 170° with Power Sliding Back ⁽²⁾ \$10,850

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module
2. Only available on seat sized 17" wide and greater

ADD: ELEVATING LIFT SEAT MODULE

☐ ESM12	12" Elevating Lift Seat Module ⁽¹⁾⁽²⁾⁽³⁾ \$2,495 Weight capacity up to 250lbs
--------------------------------	---

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.
2. Not available on seat sizes greater than 20" wide by 20" deep
3. For weights greater than 175lbs tilt while elevated is not available.

ADD: MODULES

☐ PRETILT	5° Fixed Pre-Tilt ⁽¹⁾ N/C
☐ ERSSP	Electronic Recline Start & Stop Program Function ⁽²⁾ \$500
☐ PRM8	8° Precline Module ⁽³⁾ N/C
☐ PRM30	30° Precline Module ⁽³⁾ \$500

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height
2. Requires Multi-Function M610 through Drive Control (MFCBD). Set range between 90° and 170°.
3. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available.

ADULT SEAT WIDTH SETTINGS

☐ SW116	16" Wide (Adjustable from 16" to 20") N/C
☐ SW117	17" Wide (Adjustable from 16" to 20") N/C
☐ SW118	18" Wide (Adjustable from 16" to 20") N/C
☐ SW119	19" Wide (Adjustable from 16" to 20") N/C
☐ SW120	20" Wide (Adjustable from 16" to 20") N/C
☐ W219	19" Wide (Adjustable from 19" to 22") \$460
☐ W220	20" Wide (Adjustable from 19" to 22") \$460
☐ W221	21" Wide (Adjustable from 19" to 22") \$460
☐ W222	22" Wide (Adjustable from 19" to 22") \$460

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

☐ SD116	16" Deep (Adjustable from 16" to 21") N/C
☐ SD117	17" Deep (Adjustable from 16" to 21") N/C
☐ SD118	18" Deep (Adjustable from 16" to 21") N/C
☐ SD119	19" Deep (Adjustable from 16" to 21") N/C
☐ SD120	20" Deep (Adjustable from 16" to 21") N/C
☐ D219	19" Deep (Adjustable from 19" to 22") \$460
☐ D220	20" Deep (Adjustable from 19" to 22") \$460
☐ D221	21" Deep (Adjustable from 19" to 22") \$460
☐ D222	22" Deep (Adjustable from 19" to 22") \$460
☐ D223	23" Deep (Adjustable from 19" to 22") \$460

SMALL ADULT SEAT WIDTH SETTINGS⁽¹⁾

☐ W314	14" Wide (Adjustable from 14" to 15") \$250
☐ W315	15" Wide (Adjustable from 14" to 15") \$250

SMALL ADULT SEAT DEPTH SETTINGS⁽¹⁾

☐ D314	14" Deep (Adjustable from 14" to 18") \$250
☐ D315	15" Deep (Adjustable from 14" to 18") \$250
☐ D316	16" Deep (Adjustable from 14" to 18") \$250
☐ D317	17" Deep (Adjustable from 14" to 18") \$250
☐ D318	18" Deep (Adjustable from 14" to 18") \$250

1. Small Adult Seat Size - weight limit is 175 lbs on these seat sizes and if recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

MATRIX SEATING OPTIONS

- ☐ Matrix Cushions available on page 10
☐ Matrix Backs available on page 11
☐ Matrix Headrests available on page 12

matrix
SEATING SERIES

CHEST STRAPS

☐ MFUS	BP Monoflex, Underarm RLS S(3" x 15.5") . \$195
☐ MFUM	BP Monoflex, Underarm RLS M(4" x 18.25") \$195
☐ MFUL	BP Monoflex, Underarm RLS L(4.5" x 21.5") . \$195
☐ MFCS	BP Monoflex, Center RLS S(2.5" x 15.5") . . \$195
☐ MFCM	BP Monoflex, Center RLS M(3.25" x 18.25") . \$195
☐ MFCL	BP Monoflex, Center RLS L(3.5" x 21.5") . . \$195
☐ PFLEX-M	PostureFlex Abdominal Support Height 6.5", length 16" to 23" \$295
☐ PFLEX-L	PostureFlex Abdominal Support Height 7", length 20" to 32" \$295
☐ QRTAB	Quad release tab for easy independent application and removal of PostureFlex \$35

LAP BELTS

☐ BELT48	Push Button Style Lap Belt, Length 48"	STD
☐ BELT60	Push Button Style Lap Belt, Length 60"	N/C
☐ BELT71	Push Button Style Lap Belt, Length 71"	N/C
☐ PLB2	Bodypoint Padded Lap Belt 2 Point ⁽¹⁾	\$250
☐ PLB4	Bodypoint Padded Lap Belt 4 Point ⁽¹⁾	\$315

1. 67" length only option

MOTION CONCEPTS HEADREST OPTIONS

☐ MCSH	Motion Concepts Standard ⁽¹⁾	N/C
☐ MCOH	Motion Concepts Auto Style ⁽¹⁾	\$235
☐ FM	Fixed Mounting (non-removable) ⁽²⁾	N/C
☐ OHR	Omit: Headrest.	N/C

1. Must select Fixed Mounting (FM) or Matrix Headrest Hardware on page 12.
2. Only compatible with Standard Rehab, Contour or High Back Models. All other upgradeable backs require removeable hardware.

MATRIX HEADRESTS OPTIONS

☐ Matrix Elan Headrests and Hardware available on page 12



FINISHED BACK HEIGHT - NON RECLINE SYSTEMS⁽¹⁾⁽²⁾

☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"	☐ 25"
FBT18	FBT19	FBT20	FBT21	FBT22	FBT23	FBT24	FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK CANES - NON RECLINE SYSTEMS

☐ SBC18	18" Straight Back Canes	STD
☐ SBC20	20" Straight Back Canes	STD
☐ SBC22	22" Straight Back Canes	STD
☐ SBC24	24" Straight Back Canes	STD

Select Canes Angle

☐ 75°	☐ 79°	☐ 85°	☐ 90°	☐ 95°	☐ 101°	☐ 106°	☐ 112°	☐ 116°
BA75	BA79	BA85	BA90	BA95	BA101	BA106	BA112	BA116

☐ ABC20	20" Angled 5° Back Canes	N/C
☐ ABC22	22" Angled 5° Back Canes	N/C
☐ ABC24	24" Angled 5° Back Canes	N/C

Select Canes Angle

☐ 81°	☐ 84°	☐ 91°	☐ 96°	☐ 101°	☐ 106°	☐ 111°	☐ 117°	☐ 121°
BA81	BA84	BA91	BA96	BA101	BA106	BA111	BA117	BA121

☐ EBC22	22" Endomorph Back Canes ⁽¹⁾	\$250
☐ EBC24	24" Endomorph Back Canes ⁽¹⁾	\$250
☐ EBCN1	1" Narrower	N/C
☐ EBCN2	2" Narrower	N/C
☐ EBCN4	4" Narrower	N/C

☐ FDBC Fold Down Back Canes w/ Pull Pin⁽²⁾ \$250

1. Comes standard with Cantilever Armrest Canes Mounted.
2. Must select back canes above. 250lb user weight capacity.

FINISHED BACK HEIGHT - RECLINE SYSTEMS⁽¹⁾⁽²⁾

☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"	☐ 25"	☐ 26"	☐ 27"
FBR20	FBR21	FBR22	FBR23	FBR24	FBR25	FBR26	FBR27

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK CANE OPTIONS - RECLINE SYSTEMS

☐ BFS2	Recline Body Fit 2" Narrower ⁽¹⁾	\$995
☐ BFS4	Recline Body Fit 4" Narrower ⁽¹⁾	\$995
☐ AHH	Adjustable Hooking Handles ⁽²⁾ (pair)	\$450
☐ RPHS	Push Handles for Recline - Straight ⁽³⁾ (pair) . .	\$350
☐ RPHA	Push Handles for Recline - Angled ⁽³⁾ (pair) . .	\$350

1. Minimum finished back widths of 16".
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

BACK OPTIONS

☐ SRB	Standard Rehab Back & Cushion	\$125
☐ PIPF	Flat Planar Interface Plate ⁽¹⁾	\$125
☐ PIPR	Recessed Planar Interface Plate ⁽¹⁾⁽²⁾	\$125
☐ ORBC	Omit Rigid Back & Cushion ⁽³⁾	N/C

1. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
2. Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
3. Not available with Recline/ESR or Recline/PSB

RIDE STYLE MULTI-AXIS INTERFACE

☐ RSMI	Ride Style Multi-Axis Interface for Recline Systems 14" to 20" wide	\$250
-------------------------------	---	-------

1. Not available with TRRO and TRRO-E.

SLING BACK FOR TILT SYSTEMS⁽¹⁾

☐ SLINGB-TA	TRX Sling Backrest with Adjustable Fabric . .	\$950
☐ SLINGB-TT	TRX Sling Backrest with Tension Fabric . .	\$1,150

1. Finished Sling Back height must be 2" shorter than back cane height.

MATRIX BACKREST OPTIONS

☐ Matrix Backs available on page 11



ARMRESTS - NON RECLINE SYSTEMS

Adult Dual Post

☐ DP1	9.75" - 12.75"	STD
☐ DP2	12.75" - 15.75"	STD
☐ DP3	15.75" - 18.75"	STD
☐ FP	Flat Pouch	STD
☐ GB	Glove Boxes	N/C
☐ OAL	Outback Arm Option Left	\$195
☐ OAR	Outback Arm Option Right	\$195

Ultra Rail Seat Mnt Flip Back Cantilever MAXX Tilt Arm

☐ RMCAN	9" - 12.5" ⁽¹⁾	\$488
☐ RMCAN2	11" - 14" ⁽¹⁾⁽²⁾	\$488
☐ DPM	Dual Post Module to RMCAN (pair)	\$250

Outback Cantilever Tilt Armrest

☐ OCTA1	9"-12" ⁽²⁾	\$488
☐ OCTA2	12"-15" ⁽²⁾	\$488
☐ OCTA2	6"-12" ⁽²⁾	\$488

Cantilver Arms for Endomorph Back Canes⁽³⁾

☐ CAFEB	See page 4	STD
--------------------------------	------------	-----

1. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.
2. Not available with Small Adult Sizes.

RECLINE ARMRESTS

Reclining Adjustable Height

☐ RC1	9.5"-13"	STD
☐ RC2	13"-16"	STD
☐ OAL	Outback Arm Option Left	\$195
☐ OAR	Outback Arm Option Right	\$195

Maxx Style Cane Mounted Height, Width, Angle Adj.

Cantilever Flip Back Armrests

☐ MCAN1	9.5"-13"	\$488
☐ MCAN2	12.5"-16"	\$488

Outback Cantilever Recline Armrest

☐ OCRA1	9"-12" ⁽¹⁾	\$488
☐ OCRA2	12"-15" ⁽¹⁾	\$488
☐ OCRA3	7" - 9" ⁽¹⁾	\$488

Cane Mounted Height and Angle Adj.

Cantilver Flip Back Armrest

☐ CANQ1	9.5"-12" ⁽¹⁾	\$488
☐ CANQ2	10"-13" ⁽¹⁾	\$488
☐ CANQ3	13"-16" ⁽¹⁾	\$488

Adult Dual Post On Recline

☐ DPOR1	9.75" - 12.75"	\$250
☐ DPOR2	12.75" - 15.75"	\$250
☐ DPOR3	15.75" - 18.75"	\$250
☐ FP	Flat Pouch	STD
☐ GB	Glove Boxes	N/C

1. Not available with Small Adult Sizes.

ARMREST TUBE LENGTH - PLEASE SPECIFY

☐ ATFL	Left Full	STD
☐ ATDL	Left Desk	STD
☐ ATFR	Right Full	STD
☐ ATDR	Right Desk	STD

ARMPADS

Standard

☐ SAPFL	Left Full (2.25" x 14")	STD
☐ SAPDL	Left Desk (2.25" x 10")	STD
☐ SAPML	Left Mini (2.25" x 7")	STD
☐ SAPFR	Right Full (2.25" x 14")	STD
☐ SAPDR	Right Desk (2.25" x 10")	STD
☐ SAPMR	Right Mini (2.25" x 7")	STD

Waterfall Arm Pads⁽¹⁾

☐ WAPFL	Left Full (13")	N/C
☐ WAPDL	Left Desk (9")	N/C
☐ WAPFR	Right Full (13")	N/C
☐ WAPDR	Right Desk (9")	N/C

Modular Arm Pad - Left⁽²⁾

☐ MAPTFL	Left Tray Full (14")	\$95
☐ MAPTDL	Left Tray Desk (10")	\$95

Select Pad Insert below:

☐ ISSFL	Iskin - Integrated soft skin foam
☐ SCVFL	Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

☐ MAPTFR	Right Tray Full (14")	\$95
☐ MAPTDR	Right Tray Desk (10")	\$95

Select Pad Insert below:

☐ ISSFR	Iskin - Integrated soft skin foam
☐ SCVFR	Startex Covered Visco foam

Flat Multi Position Arm Pads

☐ MPAFL	Left Full (14")	\$150
☐ MPADL	Left Desk (10")	\$150
☐ MPAML	Left Mini (7")	\$150
☐ MPAFR	Right Full (14")	\$150
☐ MPADR	Right Desk (10")	\$150
☐ MPAMR	Right Mini (7")	\$150

Gel Pads

☐ GAPNL	Narrow Left (2" x 12")	\$88
☐ GAPNR	Narrow Right (2" x 12")	\$88
☐ GAPWL	Wide Left (3.5" x 12")	\$96
☐ GAPWR	Wide Right (3.5" x 12")	\$96
☐ GAPLL	Long Left (3.5" x 14")	\$103
☐ GAPLR	Long Right (3.5" x 14")	\$103

Arm Troughs

☐ EATL	Ergonomic Arm Trough - Left	\$121
☐ EATR	Ergonomic Arm Trough - Right	\$121
☐ FHPL	Flat Hand Pad - Left ⁽³⁾	\$91
☐ FHPR	Flat Hand Pad - Right ⁽³⁾	\$91

Arm Pad Storage/Cup Holder⁽⁴⁾

☐ APSL	Left	\$295
☐ APSR	Right	\$295
☐ OMIT	Omit Armpads(7/8" diameter tube)	N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side, or with Outback Cantilever arms.
5. Not available with TRRO/TRRO-E
6. Full length arm pads used on the same side as the Joystick result in a "Back to Joystick" measurement of 20" minimum (24" including switch mounting)

Dealer Name: _____


Dealer Account Number: _____

PO #: _____

Tag: _____

MANUAL CENTER MOUNT FRONT RIGGING

- ☐ **MSFCMFP** Maxx Fixed Center Mount Platform⁽¹⁾ \$275
- ☐ **IFP** Comes with a rubber coated footplate (11.5" W x 10") N/C
- ☐ **ICP** Individual calf pads \$100
- Select seat pan to footplate:
- ☐ **SPTF1** 9" - 13" N/C
- ☐ **SPTF14** 13" - 19" N/C
- Set at:
- ☐ **SA70** Set at 70° N/C
- ☐ **SA90** Set at 90° N/C
- ☐ **SA97** Set at 97° N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
2. 3" ground clearance required.


MANUAL ELEVATING CENTER MOUNT LEGREST⁽¹⁾

- ☐ **MECML** Manual Elevating Center Mount Legrest With Plantar Flex \$595
- Select legrest set up position:
- ☐ **90SU** Set at 90° (elevating range 90-40°) STD
- ☐ **110TU** Tucked at 110° (elevating range 110-60°)⁽²⁾⁽³⁾ N/C
- Rubber Coated Footplate:
- ☐ **MRCP** Rubber coated footplate 11.5" W x 10" D \$165
- Select seat pan to footplate for Rubber Coated Footplate
- ☐ **CSPTF1** 11.5" - 15" N/C
- ☐ **CSPTF2** 14" - 19" N/C
- Individual Footplates:
- ☐ **IFSL** Small - (5" x 7.5" deep) - Left \$165
- ☐ **IFSR** Small - (5" x 7.5" deep) - Right \$165
- ☐ **IFML** Medium - (5.5" x 9.5" deep) - Left \$165
- ☐ **IFMR** Medium - (5.5" x 9.5" deep) - Right \$165
- ☐ **IFLL** Large - (6" x 11.5" deep) - Left \$165
- ☐ **IFLR** Large - (6" x 11.5" deep) - Right \$165
- Select seat pan to footplate for Individual Footplates
- ☐ **CSPTF3** 9.5" - 14.5" N/C
- ☐ **CSPTF4** 12" - 19" N/C
- Optional:
- ☐ **CICP** Individual calf pads \$200

1. 3" ground clearance required.
2. Not available on SP2 Narrow base.
3. Not available on seat depths 16", 17".

INDIVIDUAL FOOT PLATES FOR CENTER MOUNTS⁽¹⁾

- ☐ **IFSL** Small Left - 5" w x 7.5" d \$165
- ☐ **IFSR** Small Right - 5" w x 7.5" d \$165
- ☐ **IFML** Medium Left - 5.5" w x 9.5" d \$165
- ☐ **IFMR** Medium Right - 5.5" w x 9.5" d \$165
- ☐ **IFLL** Large Left - 6" w x 11.5" d \$165
- ☐ **IFLR** Large Right - 6" w x 11.5" d \$165
- ☐ **MFOL** Multi-Axis Footplate Left 6" w x 10" d \$225
- ☐ **MFOR** Multi-Axis Footplate Right 6" w x 10" d \$225



1. Foot plates may interfere with casters in certain configurations.

POWER CENTER MOUNT FOOT PLATFORM

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression: (Add)	
Seat Pan to Footplate: (Total)	

LNx Power Center Mount Foot Platform HCPCS code E1012

- ☐ **LNx** 9" Articulation⁽¹⁾⁽²⁾ \$3,000
- ☐ **LSPTF1** 10" - 14" Seat pan to footplate
- ☐ **LSPTF2** 13" - 19" Seat pan to footplate
- ☐ **LNxA** 6" Articulation⁽¹⁾⁽²⁾⁽⁶⁾ \$3,000
- ☒ **LASPTF2** 10" - 14" Seat pan to footplate

Options:

- ☐ **LIFP** Rubber coated footplate 11.5" W x 10" D STD
- ☐ **LICP** Individual calf pads STD
- ☐ **OPCP** One Piece Malleable Calf Panel⁽³⁾ \$250
- ☐ **LNXPf** Plantar Flex Option for LNx Foot platform⁽⁴⁾. \$250

Latitude Style Foot Platform

- ☐ **LSFP** Latitude Style Foot Platform⁽⁵⁾ \$2,995
- For finished knee to heel measurements of 6.5"-12"
- ☐ **LSPFP** Latitude Style Padded Foot Plate \$250

1. 3" ground clearance required.
2. Only available on adult range seat systems 16" deep or greater. For shorter seat depths call Customer Service for other options.
3. Not available with Plantar Flex Option for LNx Foot platform, or with LNxA.
4. Only available with Foot Platform Options.
5. Not available with TRRO and TRRO-E.
6. For clients under 4' only.

FOOT PLATFORM OPTIONS FOR CENTER MOUNTS⁽¹⁾



- ☐ **SFFP** Small Platform 12"w x 8.5"d \$180
- ☐ **MFFP** Medium Platform 12"w x 10.75"d \$180
- ☐ **LFPP** Large Platform 14.5"w x 13"d \$250
- ☐ **XLFFP** Extra Large Platform 17"w x 13"d \$250
- ☐ **XXLFFP** Extra Extra Large Platform 21"w x 13"d \$250

1. May need to adjust leg position to avoid interference with front castors. Not available with Invacare Center Mount

Ultra Low Maxx Power Positioning
on Invacare TDX-SP2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MAXX STYLE POWER ELEVATING LEGRESTS

☐ **MSPPL MAXX Swing Away Power Elevating⁽¹⁾⁽²⁾ . . \$2,490**

Select Power Legrest Function:

- ☐ **INDA Individual Movement N/C**
☐ **COMA Combined Movement N/C**
☐ **SYNA Synchronized/Recline Movement N/C**

Select Receivers:

- ☐ **RS Straight Receivers N/C**
☐ **RF Flared 1" Outward Receivers. N/C**

Select Seat Pan to Footplate:

- ☐ **SPTF5 Length: 10"-14" N/C**
☐ **SPTF2 Length: 13"-17" N/C**
☐ **SPTF7 Length: 16"-20" N/C**

Options:

- ☐ **FCPMS Flat Calf Pad Option \$275**
☐ **CCPMS Curved Calf Pad Option⁽³⁾ \$275**

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads.
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

MAXX STYLE MANUAL ELEVATING LEGRESTS

☐ **MMPPSA MAXX Swing Away Manual Elevating⁽¹⁾⁽²⁾ . . \$680**

Select Receivers:

- ☐ **RS Straight Receivers N/C**
☐ **RF Flared 1" Outward Receivers. N/C**

Select Seat Pan to Footplate:

- ☐ **SPTF5 Length: 10"-14" N/C**
☐ **SPTF2 Length: 13"-17" N/C**
☐ **SPTF7 Length: 16"-20" N/C**

Options:

- ☐ **FCPMS Flat Calf Pad Option \$275**
☐ **CCPMS Curved Calf Pad Option⁽³⁾ \$275**

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

SWING AWAY LEGREST

☐ **HD70 Heavy Duty 70° Swing Away Footrest⁽¹⁾⁽²⁾ . . . \$495**

Select Receivers:

- ☐ **RS Straight Receivers N/C**
☐ **RF Flared 1" Outward Receivers. N/C**

Select Seat Pan to Footplate:

- ☐ **SPTF11 Length: 6"-10" N/C**
☐ **SPTF8 Length: 11"-15" N/C**
☐ **SPTF2 Length: 13"-17" N/C**
☐ **SPTF7 Length: 16"-20" N/C**

Options:

- ☐ **HDFCP Flat Calf Pad Option \$275**
☐ **HD70CCP Curved Calf Pad Option⁽³⁾ \$275**
☐ **DCS Dynamic Coil Springs \$190**

☐ **I70SAF Invacare 70° Swing-away Footrests \$325**

- ☐ **I70-CF Comes with composite foot plates N/C**

☐ **70HDAIH Invacare 70° Non-Tapered HD Swing-away Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards . \$457**

☐ **70TAPAIH Invacare Tapered 70° Swing-away Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards . \$497**

- ☐ **ALPT Longer Pivot and Slide Tube for 70° Swingaways (Pair). \$132**

☐ **RLSL Residual Limb Support - Left⁽⁴⁾ 6.75"W x 9.75"L x 2"H (each) \$495**

☐ **RLSR Residual Limb Support - Right⁽⁴⁾ 6.75"W x 9.75"L x 2"H (each) \$495**

1. Must choose footplate on page 8
2. Heavy duty durable construction (4.5lbs ea.).
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.
4. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.





Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FOOTPLATES FOR SWING AWAY/ELEVATING LEGREST⁽¹⁾






Angle Adjustable Footplates HCPCS code K0040

☐ AFSL	Left Small Footplate (4.25"W x 8"D)	\$165
☐ AFSR	Right Small Footplate (4.25"W x 8"D)	\$165
☐ AFML	Left Medium Footplate (5.25"W x 8"D)	\$165
☐ AFMR	Right Medium Footplate (5.25"W x 8"D)	\$165
☐ AFLL	Left Large Footplate (6.25"W x 8"D)	\$165
☐ AFLR	Right Large Footplate (6.25"W x 8"D)	\$165
☐ AFXL	Left Extra Large Footplate (9"W x 11"D)	\$165
☐ AFXR	Right Extra Large Footplate (9"W x 11"D)	\$165

Multi-Axis Adjustable Angle Footplates HCPCS code K0040

☐ MFSL	Left Small Footplate (4.25"W x 8"D)	\$225
☐ MFSR	Right Small Footplate (4.25"W x 8"D)	\$225
☐ MFML	Left Medium Footplate (5.25"W x 8"D)	\$225
☐ MFMR	Right Medium Footplate (5.25"W x 8"D)	\$225
☐ MFL	Left Large Footplate (6.25"W x 8"D)	\$225
☐ MFLR	Right Large Footplate (6.25"W x 8"D)	\$225
☐ MFXL	Left Extra Large Footplate (9"W x 11"D)	\$225
☐ MFXR	Right Extra Large Footplate (9"W x 11"D)	\$225



Cut Out Angle Adjustable Footplates HCPCS code K0040

☐ AFCSL	Left Small Footplate (7.25"W x 9"D)	\$165
☐ AFCSR	Right Small Footplate (7.25"W x 9"D)	\$165
☐ AFCML	Left Medium Footplate (7.5"W x 11"D)	\$165
☐ AFCMR	Right Medium Footplate (7.5"W x 11"D)	\$165
☐ AFCLL	Left Large Footplate (10"W x 14"D)	\$165
☐ AFCLR	Right Large Footplate (10"W x 14"D)	\$165

Cut Out Multi-Axis Adjustable Angle Footplates HCPCS code K0040

☐ MFCSL	Left Small Footplate (7.25"W x 9"D)	\$225
☐ MFCSR	Right Small Footplate (7.25"W x 9"D)	\$225
☐ MFCML	Left Medium Footplate (7.5"W x 11"D)	\$225
☐ MFCMR	Right Medium Footplate (7.5"W x 11"D)	\$225
☐ MFCLL	Left Large Footplate (10"W x 14"D)	\$225
☐ MFCLR	Right Large Footplate (10"W x 14"D)	\$225

1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

HEEL LOOPS, CUPS AND TOE STRAPS

☐ HLL	Left Heel Loop (each)	\$20
☐ HLR	Right Heel Loop (each)	\$20
☐ HCL	Left Cup Loop (each)	\$20
☐ HCR	Right Cup Loop (each)	\$20
☐ TSL	Left Toe Strap (each)	\$20
☐ TSR	Right Toe Strap (each)	\$20

CALF STRAPS AND CALF PANELS

☐ CSM	Body Point Airmesh Calf Strap for Chairs 15"-18" wide with 3.5" H Pad ⁽¹⁾	\$200
☐ CSL	Body Point Airmesh Calf Strap for Chairs 18"-22" wide with 3.5" H Pad ⁽¹⁾	\$200
☐ CPM	Body Point Airmesh Calf Panel for Chairs 15"-18" wide with 35.5" W x 9" H Pad ⁽¹⁾	\$295
☐ CPL	Body Point Airmesh Calf Panel for Chairs 18"-22" wide with 38.5" W x 9" H Pad ⁽¹⁾	\$315

THORACIC LATERALS AND HIP/KNEE SUPPORTS

- ☐ Matrix Laterals please see page 12
- ☐ Maxx Laterals and Hip/Knee supports please see page 13

SINGLE-DUAL FUNCTION ELECTRONICS

☐ SFCB Single/Dual Function Electronics⁽¹⁾ STD

Select Switch:

- ☐ SPBO Single Push Button
☐ STO Single Toggle
☐ DPBO Dual Push Button
☐ SOS Omit Switch
☐ JOY Through Joystick
☐ SC Through Specialty Control

Optional:

- ☐ RDE Reduced Drive for Elevating Seat⁽²⁾ \$400

☐ SPBSA Additional Single Push Button \$150
☐ STA Additional Single Toggle Switch \$150
☐ DPBA Additional Dual Push Button \$150

1. These electronics include drive lock-out and tilt limit.

2. If reduced drive is not selected the system will be in drive lock-out when elevated.

MULTI- FUNCTION ELECTRONICS

☐ MFCB Multi-Function through
Switch/Toggle⁽¹⁾ \$1,950

Select Switch:

- ☐ QPBO Quad Push Button
☐ FWTO Four-Way Toggle
☐ EWRO 8-Way Rocker⁽²⁾

☐ MFCBD Multi-Function through
Drive Control⁽¹⁾ \$1,950

Select Control:

- ☐ JOY Through Joystick
☐ SC Through Specialty Control

Optional:

- ☐ RDY Reduced Drive for Elevating Seat⁽³⁾
☐ RDN No Reduced Drive for Elevating Seat⁽³⁾

☐ QPBA Additional Quad Button \$150
☐ FWTA Additional Four-Way Toggle Switch \$150
☐ EWRA Additional 8-Way Rocker⁽²⁾ \$255
☐ 10WA Additional 10-Way Switch \$255

1. These electronics include seating attendant control, drive lock-out, tilt/recline limit and actuator speed control

2. Dual function is not available with 8-Way Rocker.

3. If reduced drive is not selected system will be in drive lock-out when elevated.

ELECTRONIC SWITCH MOUNTING

- ☐ ESML Left Armpad
☐ ESMR Right Armpad
☐ ESMJL Left Joystick (Mounts to Joystick Bracket)
☐ ESMJR Right Joystick (Mounts to Joystick Bracket)
☐ ESMD Do not mount

ELECTRONIC OPTIONS

Switches

- ☐ SPBSA Single Push Button \$150
☐ ASL320R ASL Pulse Mechanical Switch - Red \$200
☐ ASL320G ASL Pulse Mechanical Switch - Green \$200
☐ ASL320BK ASL Pulse Mechanical Switch - Black \$200
☐ FTSQ Feather Touch - Single Quarter Size \$295
☐ FTMD Feather Touch - Dual Dime Size \$295
☐ FTMQ Feather Touch - Quad Dime Size⁽⁶⁾ \$295

Specify Use:

- ☐ Seating
☐ Joystick Mode⁽⁴⁾⁽⁵⁾
☐ Joystick On/Off⁽⁴⁾⁽⁵⁾

Miscellaneous

- ☐ S46L Single to Dual Switch Adapter \$265
Allows for dual use of a single switch
☐ M616S M616 Converter⁽¹⁾ \$125
Allows up to 4 mono phono switches into Multi-Function box
☐ LJS Lazarus Jump Starter \$150
☐ APPS Accessory Port Power Supply⁽¹⁾
(3 amp capacity) \$195
☐ TPSMH 12 Piece Switch Mounting Hardware
Kit (additional kits) \$39

1. Not available with Single-Dual Function Electronics.

2. Not available for Mode or On/Off.

3. Requires M616 Converter when using Multi-function electronics.

4. IN200 is required with REM210, REM211, REM216 joysticks.

5. Not available with REM110 joystick.

6. Only available for Seating.

VENT TRAYS, O2 HOLDERS, ACCESSORIES

- ☐ Vent Trays, O2 Holders, IV Pole and Additional Accessories available
on page 14

INSTALLATION

- ☒ MCI Motion Concepts Installation STD

PLEASE SPECIFY PAGES INCLUDED WITH ORDER

- ☐ Pages 1-9 Excludes Matrix and Maxx Seating Accessories
Thoracic Supports, Hip Supports, Vent Trays and
Accessories
☐ Pages 1-15 Complete Order Form

Total Retail Value \$

SUBMIT

PRINT FORM

SAVE

CLEAR FORM



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX SEAT CUSHION OPTIONS

Matrix Libra Adjust Seat Cushion 16" - 20" wide \$985

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ LC1416-ADJ | ☐ LC1516-ADJ | ☐ LC1518-ADJ | ☐ LC1520-ADJ |
| ☐ LC1616-ADJ | ☐ LC1618-ADJ | ☐ LC1620-ADJ | ☐ LC1716-ADJ |
| ☐ LC1718-ADJ | ☐ LC1720-ADJ | ☐ LC1816-ADJ | ☐ LC1818-ADJ |
| ☐ LC1820-ADJ | ☐ LC1916-ADJ | ☐ LC1918-ADJ | ☐ LC1920-ADJ |
| ☐ LC2016-ADJ | ☐ LC2018-ADJ | ☐ LC2020-ADJ | |

Matrix Libra Adjust Seat Cushion 22" depth. \$1,265 ^{NEW}

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ LC1622-ADJ | ☐ LC1722-ADJ | ☐ LC1822-ADJ | ☐ LC1922-ADJ |
| ☐ LC2022-ADJ | | | |

Matrix Libra Adjust Seat Cushion 21" - 22" wide \$1,265

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ LC2118-ADJ | ☐ LC2120-ADJ | ☐ LC2122-ADJ |
| ☐ LC2218-ADJ | ☐ LC2220-ADJ | ☐ LC2222-ADJ |

Matrix Libra Seat Cushion 16" - 20" wide \$745

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ LC1416 | ☐ LC1516 | ☐ LC1518 | ☐ LC1520 |
| ☐ LC1616 | ☐ LC1618 | ☐ LC1620 | ☐ LC1716 |
| ☐ LC1718 | ☐ LC1720 | ☐ LC1816 | ☐ LC1818 |
| ☐ LC1820 | ☐ LC1916 | ☐ LC1918 | ☐ LC1920 |
| ☐ LC2016 | ☐ LC2018 | ☐ LC2020 | |

Matrix Libra Seat Cushion 22" deep \$1,025

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ LC1622 | ☐ LC1722 | ☐ LC1822 | ☐ LC1922 |
| ☐ LC2022 | | | |

Matrix Libra HD Seat Cushion 21" - 22" wide \$1,025

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ LC2118 | ☐ LC2120 | ☐ LC2122 | ☐ LC2218 |
| ☐ LC2220 | | | |

Matrix PS Seat Cushion 16" - 20" wide. \$569

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ PS1416 | ☐ PS1616 | ☐ PS1618 | ☐ PS1620 |
| ☐ PS1716 | ☐ PS1718 | ☐ PS1720 | ☐ PS1816 |
| ☐ PS1818 | ☐ PS1820 | ☐ PS2016 | ☐ PS2018 |
| ☐ PS2020 | | | |

Matrix PS HD Seat Cushion 22" - 22" wide \$975

- | | | |
|---------------------------------|---------------------------------|---------------------------------|
| ☐ PS2218 | ☐ PS2220 | ☐ PS2222 |
|---------------------------------|---------------------------------|---------------------------------|

Matrix PSP Seat Cushion 14" - 20" wide \$645

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PSP1416 | ☐ PSP1616 | ☐ PSP1618 | ☐ PSP1620 |
| ☐ PSP1716 | ☐ PSP1718 | ☐ PSP1720 | ☐ PSP1816 |
| ☐ PSP1818 | ☐ PSP1820 | ☐ PSP2016 | ☐ PSP2018 |
| ☐ PSP2020 | | | |

Matrix PSP HD Seat Cushion 22" - 22" wide \$995

- | | | |
|----------------------------------|----------------------------------|----------------------------------|
| ☐ PSP2218 | ☐ PSP2220 | ☐ PSP2222 |
|----------------------------------|----------------------------------|----------------------------------|

Matrix Vi Seat Cushion 14" - 20" wide \$624

- | | | | |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| ☐ MA1415-VI | ☐ MA1417-VI | ☐ MA1616-VI | ☐ MA1618-VI |
| ☐ MA1620-VI | ☐ MA1716-VI | ☐ MA1718-VI | ☐ MA1720-VI |
| ☐ MA1816-VI | ☐ MA1818-VI | ☐ MA1820-VI | |

MATRIX SEAT CUSHION OPTIONS CONTINUED

Matrix Vi Seat Cushion 21" - 22" wide \$681

- | | | | |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| ☐ MA2016-VI | ☐ MA2018-VI | ☐ MA2020-VI | ☐ MA2218-VI |
| ☐ MA2220-VI | | | |

Matrix Vi HD Seat Cushion 22" wide \$840

- | |
|------------------------------------|
| ☐ MA2222-VI |
|------------------------------------|

Matrix Multi Cushion 14" - 20" wide

with reversable cover \$745

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ MM1414 | ☐ MM1416 | ☐ MM1616 | ☐ MM1616 |
| ☐ MM1618 | ☐ MM1620 | ☐ MM1716 | ☐ MM1718 |
| ☐ MM1720 | ☐ MM1816 | ☐ MM1818 | ☐ MM1820 |
| ☐ MM2016 | ☐ MM2018 | ☐ MM2020 | |

^{NEW} Matrix Multi Cushion 22" wide


with reversable cover \$995

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ MM2216 | ☐ MM2218 | ☐ MM2220 | ☐ MM2222 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|

FABRIC OPTIONS FOR CUSHIONS

- | | |
|-------------------------------|--|
| ☐ SFLF | Spacetex Fabric Option for
Libra Cushion N/C |
| ☐ ICFL | Infection Control Fabric Upgrade for
Libra Cushion \$50 |

CUSHION COVER REFERENCE GUIDE

					
					
Cushion	Meshtex	Startex	Infection	Spacetex	Onyx
Libra/Adjust	N/A	STD	\$50	N/C	N/A



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

UPGRADABLE BACK OPTIONS

Matrx Elite Back 14" - 20" W (3" Contour) \$766

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBE1414 | ☐ PBE1416 | ☐ PBE1418 | ☐ PBE1420 |
| ☐ PBE1514 | ☐ PBE1516 | ☐ PBE1518 | ☐ PBE1520 |
| ☐ PBE1614 | ☐ PBE1616 | ☐ PBE1618 | ☐ PBE1620 |
| ☐ PBE1714 | ☐ PBE1716 | ☐ PBE1718 | ☐ PBE1720 |
| ☐ PBE1814 | ☐ PBE1816 | ☐ PBE1818 | ☐ PBE1820 |
| ☐ PBE1916 | ☐ PBE1918 | ☐ PBE1920 | ☐ PBE2014 |
| ☐ PBE2016 | ☐ PBE2018 | ☐ PBE2020 | |

Matrx Elite HD Back 21" - 22" W (3" Contour) \$1,095

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE2116-HD | ☐ PBE2120-HD | ☐ PBE2216-HD | ☐ PBE2220-HD |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|

Matrx Elite Deep Back 14" - 20" W (6" Contour) \$918

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ EDB1416 | ☐ EDB1418 | ☐ EDB1420 | ☐ EDB1514 |
| ☐ EDB1516 | ☐ EDB1518 | ☐ EDB1520 | ☐ EDB1614 |
| ☐ EDB1616 | ☐ EDB1618 | ☐ EDB1620 | ☐ EDB1714 |
| ☐ EDB1716 | ☐ EDB1718 | ☐ EDB1720 | ☐ EDB1814 |
| ☐ EDB1816 | ☐ EDB1818 | ☐ EDB1820 | ☐ EDB1914 |
| ☐ EDB1916 | ☐ EDB1918 | ☐ EDB1920 | ☐ EDB2014 |
| ☐ EDB2016 | ☐ EDB2018 | ☐ EDB2020 | |

Matrx Elite Deep HD Back 21" - 22" W (6" Contour) . . . \$1,095

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ EDB2116-HD | ☐ EDB2120-HD | ☐ EDB2216-HD |
| ☐ EDB2220-HD | | |

Matrx Elite E2 Back

14" - 20" W x 14"/16"H (3" Contour) \$965

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2S1414 | ☐ E2S1416 | ☐ E2S1514 | ☐ E2S1516 |
| ☐ E2S1614 | ☐ E2S1616 | ☐ E2S1714 | ☐ E2S1716 |
| ☐ E2S1814 | ☐ E2S1816 | ☐ E2S1914 | ☐ E2S1916 |
| ☐ E2S2014 | ☐ E2S2016 | | |

Matrx Elite E2 Back

14" - 20" W x 18"/20"H (3" Contour) \$995

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2S1418 | ☐ E2S1420 | ☐ E2S1518 | ☐ E2S1520 |
| ☐ E2S1618 | ☐ E2S1620 | ☐ E2S1718 | ☐ E2S1720 |
| ☐ E2S1818 | ☐ E2S1820 | ☐ E2S1918 | ☐ E2S1920 |
| ☐ E2S2018 | ☐ E2S2020 | | |

Matrx Elite E2 Deep Back

14" - 20" W x 14"/16"H (6" Contour) \$1,125

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2D1414 | ☐ E2D1416 | ☐ E2D1514 | ☐ E2D1516 |
| ☐ E2D1614 | ☐ E2D1616 | ☐ E2D1714 | ☐ E2D1716 |
| ☐ E2D1814 | ☐ E2D1816 | ☐ E2D1914 | ☐ E2D1916 |
| ☐ E2D2014 | ☐ E2D2016 | | |

Matrx Elite E2 Deep Back

14" - 20" W x 18"/20"H (6" Contour) \$1,155

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ E2D1418 | ☐ E2D1420 | ☐ E2D1518 | ☐ E2D1520 |
| ☐ E2D1618 | ☐ E2D1620 | ☐ E2D1718 | ☐ E2D1720 |
| ☐ E2D1818 | ☐ E2D1820 | ☐ E2D1918 | ☐ E2D1920 |
| ☐ E2D2018 | ☐ E2D2020 | | |

UPGRADABLE BACK OPTION

Matrx Elite TR Back 15" - 20" W (3" Contour) \$866

- | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ PBE1516-TR | ☐ PBE1518-TR | ☐ PBE1520-TR | ☐ PBE1616-TR |
| ☐ PBE1618-TR | ☐ PBE1620-TR | ☐ PBE1716-TR | ☐ PBE1718-TR |
| ☐ PBE1720-TR | ☐ PBE1816-TR | ☐ PBE1818-TR | ☐ PBE1820-TR |
| ☐ PBE1916-TR | ☐ PBE1918-TR | ☐ PBE1920-TR | ☐ PBE2016-TR |
| ☐ PBE2018-TR | ☐ PBE2020-TR | | |

Matrx Elite TR HD Back 21" - 22" W (3" Contour) \$1,151

- | | |
|---------------------------------------|---------------------------------------|
| ☐ PBE2120-TRHD | ☐ PBE2220-TRHD |
|---------------------------------------|---------------------------------------|

Matrx PB Posture Back 14" - 20" W (5" Contour) \$734

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ MPB1412 | ☐ MPB1612 | ☐ MPB1616 | ☐ MPB1620 |
| ☐ MPB1712 | ☐ MPB1716 | ☐ MPB1720 | ☐ MPB1812 |
| ☐ MPB1816 | ☐ MPB1820 | ☐ MPB2012 | ☐ MPB2016 |
| ☐ MPB2020 | | | |

Matrx PB Posture Back Deep 14" - 20" W (7" Contour) . . . \$892

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PBD1412 | ☐ PBD1612 | ☐ PBD1616 | ☐ PBD1620 |
| ☐ PBD1712 | ☐ PBD1716 | ☐ PBD1720 | ☐ PBD1812 |
| ☐ PBD1816 | ☐ PBD1820 | ☐ PBD2012 | ☐ PBD2016 |
| ☐ PBD2020 | | | |

Matrx PB Heavy Duty Back

20" - 22" W (5" Contour) \$1,051

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ HDB2016 | ☐ HDB2020 | ☐ HDB2216 | ☐ HDB2220 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrx PB Heavy Duty Deep Back

20"-22" W (5" Contour) \$1,209

- | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ DHDB2016 | ☐ DHDB2020 | ☐ DHDB2216 | ☐ DHDB2220 |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|

NEW Matrx MAC Back⁽¹⁾ \$1,495

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ MAC1616 | ☐ MAC1618 | ☐ MAC1620 | ☐ MAC1622 |
| ☐ MAC1716 | ☐ MAC1718 | ☐ MAC1720 | ☐ MAC1722 |
| ☐ MAC1816 | ☐ MAC1818 | ☐ MAC1820 | ☐ MAC1822 |
| ☐ MAC1916 | ☐ MAC1918 | ☐ MAC1920 | ☐ MAC1922 |
| ☐ MAC2016 | ☐ MAC2018 | ☐ MAC2020 | ☐ MAC2022 |

1. Must select Pad Contour in Matrx MAC Back Pad Contour section below.

NEW Matrx MAC Back w/ HUG⁽¹⁾ \$1,695

- | | | |
|--------------------------------------|--------------------------------------|--------------------------------------|
| ☐ MAC1616-HUG | ☐ MAC1618-HUG | ☐ MAC1620-HUG |
| ☐ MAC1622-HUG | ☐ MAC1716-HUG | ☐ MAC1718-HUG |
| ☐ MAC1720-HUG | ☐ MAC1722-HUG | ☐ MAC1816-HUG |
| ☐ MAC1818-HUG | ☐ MAC1820-HUG | ☐ MAC1822-HUG |
| ☐ MAC1916-HUG | ☐ MAC1918-HUG | ☐ MAC1920-HUG |
| ☐ MAC1922-HUG | ☐ MAC2016-HUG | ☐ MAC2018-HUG |
| ☐ MAC2020-HUG | ☐ MAC2022-HUG | |

1. Must select Pad Contour in Matrx MAC Back Pad Contour section below.

Matrx MAC Back Pad Contour

- | | | | |
|---------------------------------|-----------|---------------------------------|------------|
| ☐ MLCSL | Left 3.5" | ☐ MLCSR | Right 3.5" |
| ☐ MLCDL | Left 6" | ☐ MLCDR | Right 6" |
| ☐ MLCXDL | Left 9" | ☐ MLCXDR | Right 9" |



BACK COVER REFERENCE GUIDE

	Meshtex	Infection Control	Onyx
Back	Meshtex	Infection	Onyx
High	STD	N/A	\$103
Elite/ Elite HD/Elite Deep	STD	\$103	N/A
Elite TR/ Elite TR HD	STD	\$103	N/A
PB/PB HD/PB Deep/PB HD Deep	STD	\$103	N/A

FABRIC OPTIONS FOR BACKS

- ☐ EBSS Super Soft Foam for Elite Back Upgrade . . . \$103
- ☐ WSCU Infection Control Cover for Elite/PB Backs . . . \$103
purchased with back
- ☐ WSC Infection Control Cover for Elite/PB Backs . . . \$150
purchased with without back
- ☐ CHIEB Chest Harness Interface Kit for Elite Back⁽¹⁾⁽²⁾ . \$103
- ☐ CHIPB Chest Harness Interface Kit for PB Back⁽¹⁾⁽²⁾ . . . \$52

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harness

MOTION CONCEPTS SINGLE POST BACKREST

- ☐ SPBBI Single Post Angle Adjustable
Fold Down Back - Matrix Back⁽¹⁾⁽³⁾⁽⁴⁾ \$750
Select Matrix Back seperally
 - ☐ SPBCI Single Post Angle Adjustable
Fold Down Back - Contour Back⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾ \$995
- Select Back: ☐ CB1617 ☐ CB1717 ☐ CB1817 ☐ CB1917 ☐ CB2017
- ☐ OCUCB Onyx Fabric Upgrade for Contour Backs . . . \$100

1. Single Post Back is an angle adjustable back between 73° and 109° in 6° incre. .
ments from 17° forward to 19° back. The back folds down and is only available . .
for 16"-20" seat widths & 21"-24" back height. A 4" gap from bottom of back to
seat pan is required.

2. Includes a Contour Back with fabric cover.

3. This option is not available with a vent tray and has a 250lb weight capacity.

4. Cane mounted armrest are not available with this option.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX ELAN HEADREST PAD⁽¹⁾

Standard Pads

- ☐ ESP6 Matrix Elan Standard Pad 6" \$310
- ☐ ESP10 Matrix Elan Standard Pad 10" \$310
- ☐ ESP14 Matrix Elan Standard Pad 14" \$310
- ☐ ESP6-HUG Matrix Elan Hug Standard Pad 6" \$310
- ☐ ESP10-HUG Matrix Elan Hug Standard Pad 10" \$310
- ☐ ESP14-HUG Matrix Elan Hug Standard Pad 14" \$310

Standard Pads with Infection Control Cover

- ☐ ESP6-IC Matrix Elan Standard Pad 6"
with Infection Control \$335
- ☐ ESP10-IC Matrix Elan Standard Pad 10"
with Infection Control \$335
- ☐ ESP14-IC Matrix Elan Standard Pad 14"
with Infection Control \$335
- ☐ ESP6-IC-HUG Matrix Elan Hug Standard Pad 6"
with Infection Control \$335
- ☐ ESP10-IC-HUG Matrix Elan Hug Standard Pad 10"
with Infection Control \$335
- ☐ ESP14-IC-HUG Matrix Elan Hug Standard Pad 14"
with Infection Control \$335

Occipital /4-Point Pads

- ☐ EOP9 Matrix Elan Occipital Pad 9" \$341
- ☐ EOP12 Matrix Elan Occipital Pad 12" \$341
- ☐ E4PTSM Matrix Elan 4-Point Pad 9"W x 7"H \$448
- ☐ E4POINT Matrix Elan 4-Point Pad 11"W x 10"H \$448

Occipital / 4-Point Pads with Infection Control Cover

- ☐ EOP9-IC Matrix Elan Occipital Pad 9"
with Infection Control \$366
- ☐ EOP12-IC Matrix Elan Occipital Pad 12"
with Infection Control \$366
- ☐ E4PTSM-IC Matrix Elan 4-Point Pad 9"W x 7"H
with Infection Control \$473
- ☐ E4POINT-IC Matrix Elan 4-Point Pad 11"W x 10"H
with Infection Control \$473

1. If mounting on PB Elite must order removable hardware

HEADREST MOUNTING HARDWARE OPTIONS

- ☐ LXSHW Loxx Standard Hardware 15" \$275
- ☐ LXMMH Loxx Mini Hardware 11" \$275
- ☐ MEHWRE Elan Hardware 15" \$285
- ☐ MEMH Elan Mini Hardware 11" \$285
- ☐ MEHW-EXT Elan Extended Hardware \$320
- ☐ FPLXSHW FLIPP Loxx Standard Hardware 15" \$705
- ☐ FPLXSHW-EX FLIPP Loxx Extended Hardware \$705
- ☐ FPLXMMH FLIPP Loxx Mini Hardware 11" \$705
- ☐ FPLXMMH-EX FLIPP Loxx Mini Extended Hardware \$705
- ☐ FPMEHW FLIPP Elan Hardware 15" \$715
- ☐ FPMEMH FLIPP Elan Mini Hardware 11" \$715
- ☐ FPMEHW-EXT FLIPP Elan Extended Hardware \$780

MAXX LATERAL TRUNK SUPPORTS⁽¹⁾

Lateral Trunk Support Left

☐ LATXL	Size XS (3.25" x 4.25")	\$100
☐ LATSL	Size S (3.5" x 5.25")	\$100
☐ LATML	Size M (3.75" x 5.75")	\$100
☐ LATLL	Size L (4.25" x 6.5")	\$100

Add:

☐ SALATL	MAXX Swing Away Hardware - Left ⁽²⁾	\$195
☐ ATLL	Additional Telescoping Link - Left ⁽³⁾	\$75

Lateral Trunk Support Right

☐ LATXR	Size XS (3.25" x 4.25")	\$100
☐ LATSR	Size S (3.5" x 5.25")	\$100
☐ LATMR	Size M (3.75" x 5.75")	\$100
☐ LATLR	Size L (4.25" x 6.5")	\$100

Add:

☐ SALATR	MAXX Swing Away Hardware - Right ⁽²⁾	\$195
☐ ATLR	Additional Telescoping Link - Right ⁽³⁾	\$75

1. Comes standard with Fixed Mount Hardware. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

2. Comes standard with 1" telescoping link which gives 0-2.75" medial offset.

3. The amount of offset available may be impacted by back thickness.

MATRIX FIXED LATERALS - FOR ELITE / ELITE TR BACKS

☐ E2FLL	Matrix E2 Fixed Laterals-Left(6.5x5") ⁽¹⁾	\$195
☐ E2FLR	Matrix E2 Fixed Laterals-Right(6.5x5") ⁽¹⁾	\$195
☐ PBSSL	Matrix Small Lateral Fixed-Left(3x5.5")	\$113
☐ PBSSR	Matrix Small Lateral Fixed-Right(3x5.5")	\$113
☐ PBLSL	Matrix Standard Lateral Fixed-Left(7.5x5")	\$116
☐ PBLSR	Matrix Standard Lateral Fixed-Right (7.5x5")	\$116
☐ PBOLL	Matrix Offset Lateral Fixed Left	\$125
☐ PBOLR	Matrix Offset Lateral Fixed Right	\$125

1. Only available on the E2 backrests.

MATRIX SWING AWAY LATERAL TRUNK SUPPORTS⁽¹⁾

Matrix Lateral Trunk Support Left

☐ EPADXL	Size XS (3.25" x 4.25")	\$115
☐ EPADSL	Size S (3.5" x 5.25")	\$115
☐ EPADML	Size M (3.75" x 5.75")	\$115
☐ EPADLL	Size L (4.25" x 6.5")	\$115

Must select one option:

☐ SAHEBL	Matrix Elite Swing Away Hardware Left	\$250
☐ OSAHL	Matrix Offset Elite Swing Away Hardware Left	\$250

Matrix Lateral Trunk Support Right

☐ EPADXR	Size XS (3.25" x 4.25")	\$115
☐ EPADSR	Size S (3.5" x 5.25")	\$115
☐ EPADMR	Size M (3.75" x 5.75")	\$115
☐ EPADLR	Size L (4.25" x 6.5")	\$115

Must select one option:

☐ SAHEBR	Matrix Elite Swing Away Hardware Right	\$250
☐ OSAHR	Matrix Offset Elite Swing Away Hardware Right	\$250

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

HIP AND KNEE SUPPORTS

Lateral Hip Support Left⁽¹⁾

☐ HSPXLH	Size XS (5.5" W x 4" H)	\$120
☐ HSPSMLH	Size S/M (7" W x 4" H)	\$120
☐ HSPSLH	Size S (7" W x 6" H)	\$120
☐ HSPMLH	Size M (9" W x 4" H)	\$120
☐ HSPLLH	Size L (13" W x 4" H)	\$120

Lateral Hip Support Right⁽¹⁾

☐ HSPXRH	Size XS (5.5" W x 4" H)	\$120
☐ HSPSMRH	Size S/M (7" W x 4" H)	\$120
☐ HSPSRH	Size S (7" W x 6" H)	\$120
☐ HSPMRH	Size M (9" W x 4" H)	\$120
☐ HSPLRH	Size L (13" W x 4" H)	\$120

Lateral Knee Support Left⁽¹⁾

☐ HSPXLK	Size XS (5.5" W x 4" H)	\$120
☐ HSPSMLK	Size S/M (7" W x 4" H)	\$120
☐ HSPSLK	Size S (7" W x 6" H)	\$120
☐ HSPMLK	Size M (9" W x 4" H)	\$120
☐ HSPLLK	Size L (13" W x 4" H)	\$120

Lateral Knee Support Right⁽¹⁾

☐ HSPXRK	Size XS (5.5" W x 4" H)	\$120
☐ HSPSMRK	Size S/M (7" W x 4" H)	\$120
☐ HSPSRK	Size S (7" W x 6" H)	\$120
☐ HSPMRK	Size M (9" W x 4" H)	\$120
☐ HSPLRK	Size L (13" W x 4" H)	\$120

Add Removable Hardware:

MAXX Style Quick Release, Multi-Axis Mounting Hardware for Hip Supports

☐ MHSHLH	Left mounted at Hip	\$215
☐ MHSHLK	Left mounted at Knee	\$215
☐ MHSHRH	Right mounted at Hip	\$215
☐ MHSHRK	Right mounted at Knee	\$215

Swing Away Hardware for Hip Supports

☐ SHSHLH	Left mounted at Hip	\$215
☐ SHSHLK	Left mounted at Knee	\$215
☐ SHSHRH	Right mounted at Hip	\$215
☐ SHSHRK	Right mounted at Knee	\$215

Extended Maxx Style Quick Release, Multi-Axis Mounting Hardware for Hip Supports⁽²⁾

☐ EMQHLH	Left mounted at Hip	\$215
☐ EMQHLLK	Left mounted at Knee	\$215
☐ EMQHRRH	Right mounted at Hip	\$215
☐ EMQHRRK	Right mounted at Knee	\$215

Lift Off Hardware for Hip Support (Hip Position Only)

☐ LOHL	Left	\$195
☐ LOHR	Right	\$195

1. Comes standard with Fixed Mount Hardware.

2. Recommended when finished seat cushion thickness will be greater than 3.5" high.

ELBOW BLOCKS

Left Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

☐ EBXL	Size XS (3.25" x 4.25")	\$195
☐ EBSL	Size S (3.5" x 5.25")	\$195
☐ EBML	Size M (3.75" x 5.75")	\$195
☐ EBLL	Size L (4.25" x 6.5")	\$195
☐ V3	Left Verticle Extension 3"	
☐ V5	Left Verticle Extension 5"	
☐ H5	Left Horizontal Extension 5"	
☐ H7	Left Horizontal Extension 7"	

Right Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

☐ EBX	Size XS (3.25" x 4.25")	\$195
☐ EBSR	Size S (3.5" x 5.25")	\$195
☐ EBMR	Size M (3.75" x 5.75")	\$195
☐ EBLR	Size L (4.25" x 6.5")	\$195
☐ V3	Right Verticle Extension 3"	
☐ V5	Right Verticle Extension 5"	
☐ H5	Right Horizontal Extension 5"	
☐ H7	Right Horizontal Extension 7"	

1. Prevents arm from slipping.

ADDITIONAL ARM ACCESSORIES

☐ UCH	Universal Cup Holder ⁽¹⁾	\$50
☐ MACESL	Multi Axis Upper Extremity Support ⁽²⁾ Left Side	\$225
☐ MACESR	Multi Axis Upper Extremity Support ⁽²⁾ Right Side	\$225
☐ OAS2L	1" Offset Armrest Spacer Left	\$125
☐ OAS2R	1" Offset Armrest Spacer Right	\$125
☐ OAS1L	1/2" Offset Armrest Spacer Left	\$125
☐ OAS1R	1/2" Offset Armrest Spacer Right	\$125
☐ MPAPA	Multi Purpose Arm Pad Adapter Block ⁽³⁾ (each) Quantity ()	\$35

1. Not available on the joystick side.

2. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.

3. Not available with Standard Arm Pads or Waterfall Arm Pads.

ARM ADDUCTOR PAD⁽¹⁾

☐ AADSL	Left Small (6" x 2.25")	\$135
☐ AADSR	Right Small (6" x 2.25")	\$135
☐ AADLL	Left Large (8.5" x 2.25")	\$135
☐ AADLR	Right Large (8.5" x 2.25")	\$135

1. Only available with Modular or Flat Multi Position Arm Pads.

ARTICULATING TRAYS⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾

☐ USAT	Universal Style Articulating Tray	\$1,350 (14"W x Adjustable 6"-10"D)
Add: ☐ TMA	Trilogy Mounting Adapter	\$50
☐ AVT	Articulating Compact LTV Vent Tray	\$1,350
☐ AVBH	Articulating Vent Bag Holder ⁽⁶⁾	\$1,350
☐ VBHNA	Non-Articulating Vent Bag Holder	\$990
Add: ☐ ABMC	Accessory Bag Mounting and Clips	\$100

1. The user weight limit is 225 lbs.

2. The maximum weight capacity of the vent tray is 35 lbs.

3. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs.

4. Not available with TRRO and TRRO-E.

5. Additional Battery Tray is recommended (RXBT).

6. Motion Concepts does not supply the Vent Bag.

GIMBAL VENT TRAYS⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾

☐ GUT	Gimbal Universal Tray	\$1,350 (14"W x Adjustable 6"-10"D)
Add: ☐ GTMA	Trilogy Mounting Adapter	\$50
☐ GVT	Gimbal LTV Vent Tray	\$1,350
☐ GVBH	Gimbal Vent Bag Holder ⁽⁶⁾	\$1,350
Add: ☐ GABMC	Accessory Bag Mounting and Clips	\$100

1. The Gimbal Vent Trays are not available with Tilt Only and/or Elevate Only Systems.

2. The user weight limit is 200 lbs.

3. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs.

4. The maximum weight capacity of the vent tray is 20 lbs.

5. Not available with TRRO and TRRO-E.

6. Motion Concepts does not supply the Vent Bag.

O2 HOLDERS⁽¹⁾

☐ NAO2H	Non-Articulating O ² Gas Holder	\$395
☐ AOHG	Articulating O ² Gas Holder Cylinder Style ⁽²⁾⁽³⁾	\$1,250
☐ AOHL	Articulating Liquid O ² Holder (6.5" x 4.5") ⁽²⁾⁽³⁾	\$1,250
☐ NAIPL	Non-Articulating IV Pole User's Left	\$395
☐ NAIPR	Non-Articulating IV Pole User's Right	\$395
☐ AIVP	Articulating IV Pole ⁽²⁾⁽³⁾	\$1,250

1. The user weight limit is 225 lbs. The maximum weight capacity of the O2 holder is 20 lbs.

2. Required for use on Recline Systems

3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems or with Single Post Backs.

Ultra Low Maxx Power Positioning on Invacare TDX-SP2 Base

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ADDITIONAL OPTIONS & ACCESSORIES

- | | | |
|--------------------------------|--|-------|
| ☐ THS | Transfer Handles 4" - 6.5" | \$350 |
| ☐ THL | Transfer Handles 9" - 11.5" | \$350 |
| ☐ FFLTL | Flip and Fold Away 1/2 Lap Tray
Left Mounting(14"W x 10"D) ⁽¹⁾⁽²⁾ | \$795 |
| ☐ FFLTR | Flip and Fold Away 1/2 Lap Tray
Right Mounting(14"W x 10"D) ⁽¹⁾⁽²⁾ | \$795 |

1. Not available on joystick side and only available with a flat arm pad or modular arm pad.

2. The maximum weight capacity of the tray is 10 lbs.

SUBMIT

PRINT FORM

SAVE

CLEAR FORM

NOTES

[illegible]