



ULTRA LOW MAXX REHAB SEAT ON  
STORM TOQUE RX-20 BASE  
CDN PRICE LIST AND ORDER FORM

Price Effective February 1st, 2026



PRINT FORM

SAVE

CLEAR FORM

REQUIRED INFORMATION


Quote  Order

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_  
 Dealer Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_ Tag \_\_\_\_\_  
 Purchasing Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 E-mail: \_\_\_\_\_  
 RTS/Therapist: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 City: \_\_\_\_\_ Province: \_\_\_\_\_  
 Postal Code: \_\_\_\_\_ Client Gender: M F  
 Special Client Conditions: \_\_\_\_\_

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all REQUIRED\* measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness	



Customer Service: 1.866.748-7943 | Fax: 1.905.695.0583 | www.motionconcepts.com

INVACARE POWER BASES

- SRX-20 Storm Torque RX-20 ..... \$7,575  
 Select STF(1):  SLOW 16.75"^(2) ..... \$250  
 LOW 17.25"  
 MED 18.25"  
 TALL 19.25"
- 22TRY Narrow Base 22.4"^(3) ..... N/C
- 24TRY Standard Base 24.4" ..... N/C
- UW250 Weight Capacity up to 250^(4) lbs ..... STD
- U300 Weight Capacity up to 300 lbs ..... STD

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size.
2. Must order Standard Frame battery tray type. Only available on 17" seat depths and greater.
3. User weight limited to 250 lbs (UW250).
4. Required when ordering Elevate or Narrow Frame.

PERFORMANCE MOTORS

- HSMTR Performance Motors - 7.5 mph^(1)(2)(3) ..... \$995

1. Not available on the NARROWFRM.
2. Extended high speed use on rugged terrain is not recommended.
3. Performance/High Speed Motors may have less torque and may have less range than standard speed motors

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets^(1) ..... STD

1. For unoccupied use only

CASTER OPTIONS

- AB1228-3 8" x 2.5" Black Casters w/Solid Inlay ..... N/C

TIRE OPTIONS

- AB1431-3 14" x 3" Black Tire with Gel Foam Inserts ..... STD
- AB1430-3 14" x 3" Black Tire - Pneumatic Tire ..... N/C

INSTALLATION OF BATTERIES

- M50-12 SLD M MK Group 50 Battery (On Chair) ..... \$650
- M24 SLD G FT 24 Gel Flat Top Battery (On Chair) ..... \$840

BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger ..... N/C
- OMIT Omit Battery Charger ..... N/C

LINX ELECTRONICS ACCESSORIE

- LAK LINX Access Key^(1) ..... \$50
- MUSB Motion Dual USB Charger ..... N/C

1. Required for programming

ADDITIONAL ACCESSORIES

- UCH Universal Cup Holder^(1) ..... \$50

1. Not available on the joystick side.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**FRAME FINISH**

- 170P Obsidian Black ..... N/C
- 171P Ruby Red ..... N/C
- 172P Peacock Blue ..... N/C
- 173P Pink Flamingo ..... N/C
- 174P Purple Passion ..... N/C
- 175P Denim Navy ..... N/C
- 176P Lime Twist ..... N/C
- 177P Arctic Blue ..... N/C

**RIM INSERT COLORS - MUST PICK ONE**

- 170PR Obsidian Black ..... N/C
- 171PR Ruby Red ..... N/C
- 172PR Peacock Blue ..... N/C
- 173PR Pink Flamingo ..... N/C
- 174PR Purple Passion ..... N/C
- 175PR Denim Navy ..... N/C
- 176PR Lime Twist ..... N/C
- 177PR Arctic Blue ..... N/C
- 178PR Banana Split ..... N/C
- 179PR Orange Marmalade ..... N/C
- 180PR Morning Mist ..... N/C

**CONTROLLERS**

- PWH Harness Required for Expandable System ..... STD
- GTRAC-LX LiNX G-TRAC Module ..... \$1,200

**NON-EXPANDABLE DRIVE CONTROL**

- REM215 LED Drive Only Remote/Joystick with Lights<sup>(1)</sup> ..... STD  
Three drive profiles, control of two seat functions

**EXPANDABLE JOYSTICKS<sup>(1)</sup>**

- REM216 LED Expandable Remote/Joystick with Lights<sup>(3)</sup> ..... \$875  
Three drive profiles, control of up to 6 actuators and light operation. Not for use with specialty controls.
- REM400 Color 3.5" Touch Screen Remote/Joystick<sup>(4)</sup> ..... \$1,000  
Four drive profiles, control of up to 6 actuators through the drive control and light operation.
- SW400TGL Toggle Kit for REM400<sup>(5)</sup> ..... \$200

1. Must choose EXPC Expandable Controller and PWH Harness.  
2. Not available with lights and indicators  
3. Must choose lighting package separately.  
4. REM400 or REM500 is required when using specialty controls, choose specialty .. controls purchased separately.  
5. Adds 2 toggle switches to REM400.

**ATTENDANT DRIVE CONTROL CHOICES**

- ACU Proportional Attendant Control<sup>(1)</sup> ..... \$1,024
- Select mounting position for ACU
  - RPBMR Right
  - RPBML Left
  - RPBMC Center

**MOUNTING FOR DRIVER CONTROLS**

- FM-ML Fixed Mounting Left ..... N/C
- FM-MR Fixed Mounting Right ..... N/C
- SAQR Motion Height Adjustable Swing-Away Quad Link (Right Hand) ..... \$344
- SAQL Motion Height Adjustable Swing-Away Quad Link (Left Hand) ..... \$344

**REHAB SEAT**

- D0195 Ultra Low Maxx Rehab Seat ..... \$2,995

**ELEVATING SEAT MODULE**

- MD12L 12" Elevating Lift for Rehab Seat<sup>(1)</sup> ..... \$1,995  
Weight capacity up to 250lbs
- SFCB Single Function Electronics<sup>(2)</sup>
  - Select Switch:
    - SPBO Single Push Button
    - STO Single Toggle
    - DPBO Dual Push Button
    - SOS Omit Switch
    - JOY Through Joystick
  - Select Switch Mounting:
    - ESML Left Armpad
    - ESMR Right Armpad
    - ESMJL Left Joystick (Mounts to Joystick Bracket)
    - ESMJR Right Joystick (Mounts to Joystick Bracket)
    - ESMD Do not mount
  - Optional:
    - RDE Reduced Drive for Elevating Seat<sup>(3)</sup> ..... \$400

1. Not available on seat sizes greater than 20" wide by 20" deep.  
2. Required when selecting 12" Elevating Lift. These electronics include drive lock-out.  
3. If reduced drive is not selected the system will be in drive lock-out when elevated.

**ADD: MODULES**

- PRETILT 5° Fixed Pre-Tilt<sup>(1)</sup> ..... N/C

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height. .

**Maxx Rehab Seat**  
Invacare Storm Torque RX-20 Base

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ADULT SEAT WIDTH**

- SW116 16" Wide (Adjustable from 16" to 20") ..... N/C
- SW117 17" Wide (Adjustable from 16" to 20") ..... N/C
- SW118 18" Wide (Adjustable from 16" to 20") ..... N/C
- SW119 19" Wide (Adjustable from 16" to 20") ..... N/C
- SW120 20" Wide (Adjustable from 16" to 20") ..... N/C
- W219 19" Wide (Adjustable from 19" to 22") ..... \$460
- W220 20" Wide (Adjustable from 19" to 22") ..... \$460
- W221 21" Wide (Adjustable from 19" to 22") ..... \$460
- W222 22" Wide (Adjustable from 19" to 22") ..... \$460

**ADULT ULTRA RAIL SEAT DEPTH**

- SD116 16" Deep (Adjustable from 16" to 21") ..... N/C
- SD117 17" Deep (Adjustable from 16" to 21") ..... N/C
- SD118 18" Deep (Adjustable from 16" to 21") ..... N/C
- SD119 19" Deep (Adjustable from 16" to 21") ..... N/C
- SD120 20" Deep (Adjustable from 16" to 21") ..... N/C
- D219 19" Deep (Adjustable from 19" to 22") ..... \$460
- D220 20" Deep (Adjustable from 19" to 22") ..... \$460
- D221 21" Deep (Adjustable from 19" to 22") ..... \$460
- D222 22" Deep (Adjustable from 19" to 22") ..... \$460
- D223 23" Deep (Adjustable from 19" to 22") ..... \$460

**LAP BELTS**

- BELT48 Push Button Style Lap Belt, Length 48" ..... STD
- BELT60 Push Button Style Lap Belt, Length 60" ..... N/C
- BELT71 Push Button Style Lap Belt, Length 71" ..... N/C
- PLB2 Bodypoint Padded Lap Belt 2 Point<sup>(1)</sup> ..... \$250
- PLB4 Bodypoint Padded Lap Belt 4 Point<sup>(1)</sup> ..... \$315

1. 67" length only option

**FINISHED BACK HEIGHT<sup>(1)(2)</sup>**

- 18" FBT18
- 19" FBT19
- 20" FBT20
- 21" FBT21
- 22" FBT22
- 23" FBT23
- 24" FBT24
- 25" FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.  
2. Finished back height is measured from seat pan to top of backrest.

**BACK CANES - NON RECLINE SYSTEMS**

- SBC18 18" Straight Back Canes ..... STD
- SBC20 20" Straight Back Canes ..... STD
- SBC22 22" Straight Back Canes ..... STD
- SBC24 24" Straight Back Canes ..... STD

Select Cane Angle

- 75° BA75
- 79° BA79
- 85° BA85
- 90° BA90
- 95° BA95
- 101° BA101
- 106° BA106
- 112° BA112
- 116° BA116

- ABC20 20" Angled 5° Back Canes ..... N/C
- ABC22 22" Angled 5° Back Canes ..... N/C
- ABC24 24" Angled 5° Back Canes ..... N/C

Select Cane Angle

- 81° BA81
- 84° BA84
- 91° BA91
- 96° BA96
- 101° BA101
- 106° BA106
- 111° BA111
- 117° BA117
- 121° BA121

**BACK OPTIONS**

- SRB Standard Rehab Back & Cushion ..... \$125
- ORBC Omit Rigid Back & Cushion<sup>(3)</sup> ..... N/C

**ARMRESTS - NON RECLINE SYSTEMS**

- DP1 Adult Dual Post 9.75" - 12.75" ..... STD
- DP2 Adult Dual Post 12.75" - 15.75" ..... STD
- DP3 Adult Dual Post 15.75" - 18.75" ..... STD

Select:

- FP Flat Pouch ..... STD
- GB Glove Boxes ..... N/C

Optional:

- OAL Outback Arm Option Left<sup>(1)</sup> each ..... \$195
- OAR Outback Arm Option Right<sup>(1)</sup> each ..... \$195

- RMCAN Ultra Rail Seat Mnt Flip Back Cantilever  
MAXX Tilt Arm 9" - 12.5"<sup>(2)</sup> ..... \$488

Optional:

- DPM Dual Post Module to RMCAN (pair) ..... \$250

- OCTA1 Outback Cantilever Tilt Armrest 9"-12"<sup>(3)</sup> ..... \$488
- OCTA2 Outback Cantilever Tilt Armrest 12"-15"<sup>(3)</sup> ..... \$488

1. Only available with Dual Post Armrests.  
2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.

**ARMREST TUBE LENGTH - PLEASE SPECIFY**

- ATFL Left Full ..... STD
- ATDL Left Desk ..... STD
- ATFR Right Full ..... STD
- ATDR Right Desk ..... STD

**ARMPADS**

Standard

- SAPFL Left Full (2.25" x 14") ..... STD
- SAPDL Left Desk (2.25" x 10") ..... STD
- SAPML Left Mini (2.25" x 7") ..... STD
- SAPFR Right Full (2.25" x 14") ..... STD
- SAPDR Right Desk (2.25" x 10") ..... STD
- SAPMR Right Mini (2.25" x 7") ..... STD

Waterfall Arm Pads<sup>(1)</sup>


- WAPFL Left Full (13") ..... N/C
- WAPDL Left Desk (9") ..... N/C
- WAPFR Right Full (13") ..... N/C
- WAPDR Right Desk (9") ..... N/C

Gel Pads

- GAPNL Narrow Left (2" x 12") ..... \$88
- GAPNR Narrow Right (2" x 12") ..... \$88
- GAPWL Wide Left (3.5" x 12") ..... \$96
- GAPWR Wide Right (3.5" x 12") ..... \$96
- GAPLL Long Left (3.5" x 14") ..... \$103
- GAPLR Long Right (3.5" x 14") ..... \$103


1. Full length pad requires full length arm tube if using Ultra Rail mounted Tilt Arms.

**MANUAL CENTER MOUNT FRONT RIGGING**

- MSFCMFP Maxx Style Fixed Center Mount Foot Platform<sup>(1)(2)</sup> ..... \$395
  - MIFPN Rubber Coated Footplate 10.5" W x 10" ..... N/C
  - ICP Individual calf pads ..... \$200
- Select seat pan to footplate:
- SPTF1 9" - 13" ..... N/C
  - SPTF14 13" - 19" ..... N/C
- 

1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
2. 3" ground clearance required.
3. Foot platform will be set at factory setting.

**FOOT PLATFORM OPTIONS FOR CENTER MOUNTS<sup>1)</sup>**



- SFFP Small Foot Platform (12"w x 8.5"d) ..... \$180
- MFFP Medium Platform (12"w x 10.75"d) ..... \$180
- LFFP Large Platform (14.5"w x 13"d) ..... \$250
- XLFFP Extra Large Platform (17"w x 13"d) ..... \$250
- XXLFFP Extra Extra Large Platform (21"w x 13"d) ..... \$250

1. May need to adjust leg position to avoid interference with front casters

**MANUAL ELEVATING CENTER MOUNT LEGREST<sup>(1)</sup>**

- MECML Manual Elevating Center Mount Legrest With Plantar Flex. .... \$595
- Select legrest set up position:
- 90SU Set at 90° (elevating range 90-40°) ..... STD
- Rubber Coated Footplate:
- MRCP Rubber coated footplate 10.5" W x 10" D ..... N/C
- Select seat pan to footplate for Rubber Coated Footplate
- CSPTF1 11.5" - 15" ..... N/C
  - CSPTF2 14" - 19" ..... N/C
- Individual Footplates:
- IFSL Small - (5" x 7.5" deep) - Left ..... \$165
  - IFSR Small - (5" x 7.5" deep) - Right ..... \$165
  - IFML Medium - (5.5" x 9.5" deep) - Left ..... \$165
  - IFMR Medium - (5.5" x 9.5" deep) - Right ..... \$165
  - IFLL Large - (6" x 11.5" deep) - Left ..... \$165
  - IFLR Large - (6" x 11.5" deep) - Right ..... \$165
- Select seat pan to footplate for Individual Footplates
- CSPTF3 9.5" - 14.5" ..... N/C
  - CSPTF4 12" - 19" ..... N/C
- Optional:
- CICP Individual calf pads ..... \$200

1. 3" ground clearance required.

**MAXX STYLE MANUAL ELEVATING LEGRESTS**

- MMPPSA MAXX Swing Away Manual Elevating<sup>(1)(2)</sup> ..... \$680
- Select Receivers:
- RS Straight Receivers. .... N/C
  - RF Flared 1" Outward Receivers ..... N/C
- Select Seat Pan to Footplate:
- SPTF5 Length: 10"-14" ..... N/C
  - SPTF2 Length: 13"-17" ..... N/C
  - SPTF7 Length: 16"-20" ..... N/C
- Options:
- FCPMS Flat Calf Pad Option ..... \$275
  - CCPMS Curved Calf Pad Option<sup>(3)</sup> ..... \$275

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

**SWING AWAY LEGREST**

- HD70 Heavy Duty 70° Swing Away Footrest<sup>(1)(2)</sup> ..... \$495
- Select Receivers:
- RS Straight Receivers. .... N/C
  - RF Flared 1" Outward Receivers ..... N/C
- Select Seat Pan to Footplate:
- SPTF11 Length: 6"-10" ..... N/C
  - SPTF8 Length: 11"-15" ..... N/C
  - SPTF2 Length: 13"-17" ..... N/C
  - SPTF7 Length: 16"-20" ..... N/C
- Options:
- HDFCP Flat Calf Pad Option ..... \$275
  - HD70CCP Curved Calf Pad Option<sup>(3)</sup> ..... \$275
  - DCS Dynamic Coil Springs ..... \$190

- I70SAF Invacare 70° Swing-away Footrests ..... \$325
- I70-CF Comes with composite foot plates. .... N/C
- 70HDAIH Invacare 70° Non-Tapered HD Swing-away Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards ..... \$457
- 70TAPAIH Invacare Tapered 70° Swing-away Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards ..... \$497
- ALPT Longer Pivot and Slide Tube for 70° Swingaways (Pair) ..... \$132
- RLSSL Residual Limb Support - Left<sup>(4)</sup> 6.75"W x 9.75"L x 2"H (each) ..... \$495
- RLSSR Residual Limb Support - Right<sup>(4)</sup> 6.75"W x 9.75"L x 2"H (each) ..... \$495

1. Must choose footplate on page 8
2. Heavy duty durable construction (4.5lbs ea.).
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.
4. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.



Maxx Rehab Seat  
 Invacare Storm Toque RX-20 Base

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRIX SEAT CUSHIONS**

Matrx Libra Seat Cushion ..... \$745  
 LC1616    LC1618    LC1620    LC1716    LC1718  
 LC1720    LC1816    LC1818    LC1820    LC1916  
 LC1918    LC1920    LC2016    LC2018    LC2020

Matrx Libra Seat Cushion ..... \$1,025  
 LC1622    LC1722    LC1822    LC1922    LC2022

Matrx Libra HD ..... \$1,025  
 LC2118    LC2120    LC2122    LC2218    LC2220  
 LC2320    LC2322    LC2418    LC2420    LC2422

Matrx PS Seat Cushion ..... \$569  
 PS1616    PS1618    PS1620    PS1716    PS1718  
 PS1720    PS1816    PS1818    PS1820    PS2016  
 PS2018    PS2020

Matrx PS HD Seat Cushion ..... \$975  
 PS2218    PS2220    PS2222    PS2418    PS2420  
 PS2422

Matrx PSP Seat Cushion ..... \$645  
 PSP1616    PSP1618    PSP1620    PSP1716    PSP1718  
 PSP1720    PSP1816    PSP1818    PSP1820    PSP2016  
 PSP2018    PSP2020

Matrx PSP HD Seat Cushion ..... \$995  
 PSP2218    PSP2220    PSP2222    PSP2418    PSP2420  
 PSP2422

Matrx Multi Cushion with reversable cover ..... \$745  
 MM1616    MM1618    MM1620    MM1716    MM1718  
 MM1720    MM1816    MM1818    MM1820    MM2016  
 MM2018    MM2020

Matrx Multi Cushion with reversable cover ..... \$995  
 MM2216    MM2218    MM2220    MM2222

**UPGRADABLE MATRIX BACK OPTIONS**

Matrx Elite Back (3" Contour) ..... \$766  
 PBE1614    PBE1616    PBE1618    PBE1620    PBE1714  
 PBE1716    PBE1718    PBE1720    PBE1814    PBE1816  
 PBE1818    PBE1820    PBE1916    PBE1918    PBE1920  
 PBE2014    PBE2016    PBE2018    PBE2020

Matrx Elite HD Back (3" Contour) ..... \$1,095  
 PBE2116-HD    PBE2120-HD    PBE2216-HD  
 PBE2220-HD    PBE2316-HD    PBE2320-HD  
 PBE2416-HD    PBE2420-HD

Matrx Elite Deep Back (6" Contour) ..... \$918  
 EDB1614    EDB1616    EDB1618    EDB1620    EDB1714  
 EDB1716    EDB1718    EDB1720    EDB1814    EDB1816  
 EDB1818    EDB1820    EDB1914    EDB1916    EDB1918  
 EDB1920    EDB2014    EDB2016    EDB2018    EDB2020

**UPGRADABLE MATRIX BACK OPTIONS CONT.**

Matrx Elite Deep HD Back (6" Contour) ..... \$1,095  
 EDB2116-HD    EDB2120-HD    EDB2216-HD  
 EDB2220-HD    EDB2316-HD    EDB2320-HD  
 EDB2416-HD    EDB2420-HD

Matrx Elite E2 Back (3" Contour) ..... \$965  
 E2S1614    E2S1616    E2S1714    E2S1716    E2S1814  
 E2S1816    E2S1914    E2S1916    E2S2014    E2S2016

Matrx Elite E2 Back (3" Contour) ..... \$995  
 E2S1618    E2S1620    E2S1718    E2S1720    E2S1818  
 E2S1820    E2S1918    E2S1920    E2S2018    E2S2020

Matrx Elite E2 Deep Back (6" Contour) ..... \$1,125  
 E2D1614    E2D1616    E2D1714    E2D1716    E2D1814  
 E2D1816    E2D1914    E2D1916    E2D2014    E2D2016

Matrx Elite E2 Deep Back (6" Contour) ..... \$1,155  
 E2D1618    E2D1620    E2D1718    E2D1720    E2D1818  
 E2D1820    E2D1918    E2D1920    E2D2018    E2D2020

Matrx Elite TR Back (3" Contour) ..... \$866  
 PBE1616-TR    PBE1618-TR    PBE1620-TR  
 PBE1716-TR    PBE1718-TR    PBE1720-TR  
 PBE1816-TR    PBE1818-TR    PBE1820-TR  
 PBE1916-TR    PBE1918-TR    PBE1920-TR  
 PBE2016-TR    PBE2018-TR    PBE2020-TR

Matrx Elite TR HD Back (3" Contour) ..... \$1,151  
 PBE2120-TRHD    PBE2220-TRHD  
 PBE2320-TRHD    PBE2420-TRHD

Matrx MAC Back<sup>(1)</sup> ..... \$1,495  
 MAC1616    MAC1618    MAC1620    MAC1622    MAC1716  
 MAC1718    MAC1720    MAC1722    MAC1816    MAC1818  
 MAC1820    MAC1822    MAC1916    MAC1918    MAC1920  
 MAC1922    MAC2016    MAC2018    MAC2020    MAC2022

1. Must select Pad Contour in Matrx MAC Back Pad Contour section below.

Matrx MAC Back w/ HUG<sup>(1)</sup> ..... \$1,695  
 MAC1616-HUG    MAC1618-HUG    MAC1620-HUG  
 MAC1622-HUG    MAC1716-HUG    MAC1718-HUG  
 MAC1720-HUG    MAC1722-HUG    MAC1816-HUG  
 MAC1818-HUG    MAC1820-HUG    MAC1822-HUG  
 MAC1916-HUG    MAC1918-HUG    MAC1920-HUG  
 MAC1922-HUG    MAC2016-HUG    MAC2018-HUG  
 MAC2020-HUG    MAC2022-HUG

1. Must select Pad Contour in Matrx MAC Back Pad Contour section below.

Matrx MAC Back Pad Contour  
 MLCSL   Left 3.5"    MLCXR   Right 3.5"  
 MLCDL   Left 6"    MLCXR   Right 6"  
 MLCXDL   Left 9"    MLCXDR   Right 9"

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRXL ELAN HEADREST PAD<sup>(1)</sup>**

<input type="checkbox"/> ESP6	Matrx Elan Standard Pad 6" .....	\$310
<input type="checkbox"/> ESP10	Matrx Elan Standard Pad 10" .....	\$310
<input type="checkbox"/> ESP14	Matrx Elan Standard Pad 14" .....	\$310
<input type="checkbox"/> ESP6-HUG	Matrx Elan Hug Standard Pad 6" .....	\$310
<input type="checkbox"/> ESP10-HUG	Matrx Elan Hug Standard Pad 10" .....	\$310
<input type="checkbox"/> ESP14-HUG	Matrx Elan Hug Standard Pad 14" .....	\$310

1. If mounting on PB Elite must order removable hardware

**HEADREST MOUNTING HARDWARE OPTIONS**

<input type="checkbox"/> LXSHW	Loxx Standard Hardware 15" .....	\$275
<input type="checkbox"/> MEHWRE	Elan Hardware 15" .....	\$285
<input type="checkbox"/> FPLXSHW	FLIPP Loxx Standard Hardware 15" .....	\$705
<input type="checkbox"/> FPLXSHW-EX	FLIPP Loxx Extended Hardware.....	\$705
<input type="checkbox"/> FPMEHW	FLIPP Elan Hardware 15" .....	\$715

**MATRXL FIXED LATERALS - FOR ELITE / ELITE TR BACKS**

<input type="checkbox"/> PBSSL	Matrx Small Lateral Fixed-Left(3x5.5").....	\$113
<input type="checkbox"/> PBSSR	Matrx Small Lateral Fixed-Right(3x5.5").....	\$113
<input type="checkbox"/> PBLSL	Matrx Standard Lateral Fixed-Left(7.5x5").....	\$116
<input type="checkbox"/> PBLSR	Matrx Standard Lateral Fixed-Right (7.5x5").....	\$116
<input type="checkbox"/> PBOLL	Matrx Offset Lateral Fixed Left .....	\$125
<input type="checkbox"/> PBOLR	Matrx Offset Lateral Fixed Right .....	\$125

1. Not available on the MAC backrests.

**MATRXL SWING AWAY LATERAL TRUNK SUPPORTS<sup>(1)</sup>**

**Matrx Lateral Trunk Support Left**

<input type="checkbox"/> EPADXL	Size XS (3.25" x 4.25") .....	\$115
<input type="checkbox"/> EPADSL	Size S (3.5" x 5.25") .....	\$115
<input type="checkbox"/> EPADM L	Size M (3.75" x 5.75") .....	\$115
<input type="checkbox"/> EPADLL	Size L (4.25" x 6.5") .....	\$115

Must select one option:

<input type="checkbox"/> SAHEBL	Matrx Elite Swing Away Hardware Left .....	\$250
<input type="checkbox"/> OSAHL	Matrx Offset Elite Swing Away Hardware Left.....	\$250

**Matrx Lateral Trunk Support Right**

<input type="checkbox"/> EPADXR	Size XS (3.25" x 4.25") .....	\$115
<input type="checkbox"/> EPADSR	Size S (3.5" x 5.25") .....	\$115
<input type="checkbox"/> EPADMR	Size M (3.75" x 5.75") .....	\$115
<input type="checkbox"/> EPADLR	Size L (4.25" x 6.5") .....	\$115

Must select one option:

<input type="checkbox"/> SAHEBR	Matrx Elite Swing Away Hardware Right .....	\$250
<input type="checkbox"/> OSAHR	Matrx Offset Elite Swing Away Hardware Right .....	\$250

1. Not available on the MAC backrests.

**HIP AND KNEE SUPPORTS**

**Lateral Hip Support Left<sup>(1)</sup>**

<input type="checkbox"/> HSPXLH	Size XS (5.5" W x 4" H).....	\$120
<input type="checkbox"/> HSPSMLH	Size S/M (7" W x 4" H) .....	\$120
<input type="checkbox"/> HSPSLH	Size S (7"W x 6"H) .....	\$120
<input type="checkbox"/> HSPMLH	Size M (9" W x 4"H) .....	\$120
<input type="checkbox"/> HSPLLH	Size L (13" W x 4"H) .....	\$120

**Lateral Hip Support Right<sup>(1)</sup>**

<input type="checkbox"/> HSPXRH	Size XS (5.5"W x 4"H) .....	\$120
<input type="checkbox"/> HSPSMRH	Size S/M (7"W x 4"H) .....	\$120
<input type="checkbox"/> HSPSRH	Size S (7"W x 6"H) .....	\$120
<input type="checkbox"/> HSPMRH	Size M (9"W x 4"H) .....	\$120
<input type="checkbox"/> HSPLRH	Size L (13"W x 4"H) .....	\$120

**Lateral Knee Support Left<sup>(1)</sup>**

<input type="checkbox"/> HSPXLK	Size XS (5.5"W x 4" H) .....	\$120
<input type="checkbox"/> HSPSMLK	Size S/M (7"W x 4" H) .....	\$120
<input type="checkbox"/> HSPSLK	Size S (7"W x 6" H).....	\$120
<input type="checkbox"/> HSPMLK	Size M (9"W x 4"H) .....	\$120
<input type="checkbox"/> HSPLLK	Size L (13"W x 4"H) .....	\$120

**Lateral Knee Support Right<sup>(1)</sup>**

<input type="checkbox"/> HSPXRK	Size XS (5.5"W x 4"H) .....	\$120
<input type="checkbox"/> HSPSMRK	Size S/M (7"W x 4"H) .....	\$120
<input type="checkbox"/> HSPSRK	Size S (7"W x 6"H) .....	\$120
<input type="checkbox"/> HSPMRK	Size M (9"W x 4"H) .....	\$120
<input type="checkbox"/> HSPLRK	Size L (13"W x 4"H) .....	\$120

Add Removable Hardware:

**MAXX Style Quick Release, Multi-Axis Mounting Hardware for Hip Supports**

<input type="checkbox"/> MHSHLH	Left mounted at Hip .....	\$215
<input type="checkbox"/> MHSHLK	Left mounted at Knee .....	\$215
<input type="checkbox"/> MHSHRH	Right mounted at Hip .....	\$215
<input type="checkbox"/> MHSHRK	Right mounted at Knee .....	\$215

**Extended Maxx Style Quick Release, Multi-Axis Mounting Hardware for Hip Supports<sup>(2)</sup>**

<input type="checkbox"/> EMQHLH	Left mounted at Hip .....	\$215
<input type="checkbox"/> EMQHLK	Left mounted at Knee .....	\$215
<input type="checkbox"/> EMQHRH	Right mounted at Hip .....	\$215
<input type="checkbox"/> EMQHRK	Right mounted at Knee .....	\$215

1. Comes standard with Fixed Mount Hardware.

2. Recommended when finished seat cushion thickness will be greater than 3.5" high.

