



Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

SAVE

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

SUBMIT

CLEAR FORM

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

905-695-0583

PRINT FORM

PRINT FORM **SAVE** **CLEAR FORM**

REQUIRED INFORMATION

Quote Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: M F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all REQUIRED* measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness	



Customer Service: 1.800.748.7943 | Fax: 1.905.695.0583 | www.motionconcepts.com

INVACARE POWER BASES

TDXSP2 TDX-SP2 Base \$7,575

Select STF⁽¹⁾:

SLOW	<input type="checkbox"/> 16.75" ⁽²⁾ \$250
LOW	<input type="checkbox"/> 17.25"
MED	<input type="checkbox"/> 18.25"
TALL	<input type="checkbox"/> 19.25"

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size.
2. Must order 24TRY battery tray type and Transport Brackets TRBKTS. Only available on 17" seat depths and greater.

USER WEIGHT LIMITS

U300 Weight Capacity up to 300 lbs STD

PERFORMANCE MOTORS

CLIENT COST

HSMTR Performance Motors - 7.5 mph⁽¹⁾⁽²⁾⁽³⁾ \$995

1. Must order 24TRY Battery Tray Type
2. Extended high speed use on rugged terrain is not recommended.
3. Performance/High Speed Motors may have less torque and may have less range than standard speed motors.

TRANSPORT TIE DOWN

CLIENT COST

TRBKTS Wheelchair Transport Brackets⁽¹⁾ STD

TRRO Transport Ready Option⁽²⁾⁽⁴⁾ \$499

TRRO-E Transport Ready Option for Elevate⁽²⁾⁽³⁾⁽⁴⁾ \$499

1. For unoccupied use only
2. TRRO meets ISO-7176 Part 19
3. Required for Elevate Systems to meet ISO-7176 Part 19, includes two additional tie-down loops on rail.
4. Must select a Motion Concepts or Matrix Elan headrest with hardware and Motion Concepts Standard Rehab, High Back or Matrix Upgradeable back. Comes with pelvic safety belt standard, please select TROBELT on page 4.

TIRE OPTIONS

B1431-3 14" x 3" Black Tire with Gel Foam Inserts STD

B1430-3T 14" x 3" Black Tire - Pneumatic Tire N/C

FORK OPTIONS

DSFK Double-Sided Fork Package STD

BATTERY TRAY TYPE

22TRY 22NF Style Tray - Narrow Base 24" N/C

24TRY Group 24 Style Tray - Standard Base 25.5" N/C

FRAME FINISH

<input type="checkbox"/> 170P		Obsidian Black	AADL COLOUR N/C
<input type="checkbox"/> 171P		Ruby Red	 N/C
<input type="checkbox"/> 172P		Peacock Blue	AADL COLOUR N/C
<input type="checkbox"/> 174P		Purple Passion	 N/C
<input type="checkbox"/> 175P		Denim Navy	AADL COLOUR N/C
<input type="checkbox"/> 176P		Lime Twist	 N/C
<input type="checkbox"/> 177P		Arctic Blue	 N/C

RIM INSERT COLORS - MUST PICK ONE

<input type="checkbox"/> 170PR		Obsidian Black N/C
<input type="checkbox"/> 171PR		Ruby Red N/C
<input type="checkbox"/> 172PR		Peacock Blue N/C
<input type="checkbox"/> 173PR		Pink Flamingo N/C
<input type="checkbox"/> 174PR		Purple Passion N/C
<input type="checkbox"/> 175PR		Denim Navy N/C
<input type="checkbox"/> 176PR		Lime Twist N/C
<input type="checkbox"/> 177PR		Arctic Blue N/C
<input type="checkbox"/> 178PR		Banana Split N/C
<input type="checkbox"/> 179PR		Orange Marmalade	N/C
<input type="checkbox"/> 180PR		Morning Mist N/C

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

LIGHTS AND INDICATORS

CLIENT COST

LIGHTS LED Lights and Indicators⁽¹⁾\$995

1. Must order REM216, REM400 or REM500

CONTROLLERS

EXPC Expandable ControllerSTD

PWH Harness Required for Expandable System N/C

GTRAC-LX LiNX G-TRAC Module..... \$1,200

EXPANDABLE DRIVE CONTROL

REM211 LED Expandable Remote/Joystick⁽¹⁾⁽²⁾ \$875
Three drive functions, control of up to 6 actuators through the drive control. Not for use with specialty controls.

REM216 LED Expandable Remote/Joystick with Lights⁽²⁾⁽³⁾\$875
Three drive functions, control of up to 6 actuators through the drive control and light operation. Not for use with specialty controls.

REM400 Colour 3.5" Touch Screen Remote/Joystick⁽²⁾⁽⁴⁾⁽⁵⁾ \$1,100

SW400TGL Toggle Kit for REM400⁽⁶⁾\$200

REM500 LiNX Colour 3.5" Touch Screen Display Only - No Driver Control⁽²⁾⁽³⁾⁽⁴⁾ \$1,100

Select mounting position

R5MR Right
 R5ML Left

1. Not available with lights and indicators
2. Must choose EXPC Expandable Controller and PWH Harness.
3. Must choose lighting package separately.
4. Can control lighting, must choose lighting package separately.
5. May also be used as display for specialty controls, chose specialty controls separately.
6. Adds 2 toggle switches to REM400.

SPECIALTY PROPORTIONAL CONTROLS

AADL CODE - W100

CREM LiNX Compact Remote\$1,097

CREM-LF LiNX Compact Remote Low Force\$1,097

Select mounting position

MR Right
 ML Left

1. Requires REM400 or REM500

ATTENDANT DRIVE CONTROL

PRIOR APPROVAL REQUIRED

ACU Proportional Attendant Control \$1,024

Select seating style for mounting type

TAM Tilt & Elevate Mount
 RAM Recline Mount

Select mounting position

RIGHT-ACU Right
 LEFT-ACU Left

EXPANDABLE MODULES

OUT500 Output Module \$645

IN200 LiNX Switch Input Module with Power ON/OFF for REM200 Series⁽¹⁾ \$645

IN500 LiNX 9 Pin Input Module with Sip-N-Puff Nozzle⁽²⁾⁽³⁾⁽⁴⁾\$1,800
AADL CODE - W100.....

PKG32666 Therafin Sip-N-Puff Breath Tube Kit **AADL CODE - W100**..... \$440

1. Not for use with Alternative Drive Controls.
2. Only one input module per chair, requires REM400 or REM500.
3. Input module also required for any 9-pin specialty control.
4. 9 pin specialty control and Sip-N-Puff cannot both be used together.

MOUNTING FOR DRIVER CONTROLS

FM-ML Fixed Mounting LeftN/C

FM-MR Fixed Mounting RightN/C

SAQR Motion Height Adjustable Swing-Away Quad Link (Right Hand)⁽¹⁾ \$344

SAQL Motion Height Adjustable Swing-Away Quad Link (Left Hand)⁽¹⁾ \$344

MRR Maxx Resolve Swing-Away Joystick Mount (Right Hand)⁽¹⁾ \$450

MRL Maxx Resolve Swing-Away Joystick Mount (Left Hand)⁽¹⁾ \$450

1. Not available with the REM500.

LINX ELECTRONIC ACCESSORIES

LAK LiNX Access Key⁽¹⁾ **NOT COVERED UNDER AADL**\$50

MUSB Motion Dual USB ChargerN/C

1.Required for programming

JOYSTICK TOPS

PC101A Bodypoint U Shaped Handle 3" \$120

PC102A Bodypoint U Shaped Handle 4" \$120

PC107A Bodypoint Rubber Dome..... \$120

1560 T Handle Flexible Joystick Extension \$120

1561 Straight Handle Flexible Joystick Extension \$120

1826 Chin Cup \$120

1. To remove/disengage the joystick knob from the REM400 Remote pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

INSTALLATION OF BATTERIES

22NFBATTERY 22 NF Battery (On Chair) **AADL PRICING**

24BATTERY 24 Gel Battery (On Chair) **AADL PRICING**

BATTERY CHARGER OPTION

110CHARGER 110 Volt Battery Charger.....N/C

L8900 Omit Battery ChargerN/C

Ultra Low Maxx Power Positioning
on Invacare TDX-SP2 Base

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

REHAB SEAT

D0195 Ultra Low Maxx Rehab Seat⁽¹⁾ \$3,995

1. Not available with 12" Elevating Seat Lift Module, Order Ultra Low Maxx Elevating Lift Seat Only.

POWER POSITIONING SYSTEMS

MDCGT-50 Ultra Low Maxx CG Tilt (50°) Only \$6,969

MDCGT-45 Ultra Low Maxx CG Tilt (45°) Only⁽¹⁾ \$6,969

MDCGTRESR-50 Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Extended Shear Reduction(ESR) \$11,569

MDCGTRESR-45⁽¹⁾ Ultra Low Maxx CG Tilt (45°) & Power Recline (170°) w/Extended Shear Reduction(ESR) \$11,569

MDCGTRPSB-50 Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Power Sliding Back(PSB) \$12,069

MDCGTRPSB-45⁽¹⁾ Ultra Low Maxx CG Tilt (45°) & Power Recline (170°) w/Power Sliding Back(PSB)⁽¹⁾ \$12,069

MDRESR Ultra Low Maxx Power Recline (170°) w/Extended Shear Reduction(ESR) \$10,350

MDRPSB Ultra Low Maxx Power Recline (170°) w/Power Sliding Back(PSB)⁽²⁾ \$10,850

MDES12 Ultra Low Maxx 12" Elevating Lift Seat Only⁽³⁾ \$4,995
Weight capacity up to 250lbs

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module
2. Only available on seat sized 17" wide and greater
3. Not available on seat sizes greater than 20" wide by 20" deep.

ADD: ELEVATING LIFT SEAT MODULE

CLIENT COST

ESM12 12" Elevating Lift Seat Module⁽¹⁾⁽²⁾⁽³⁾ \$2,495
Weight capacity up to 250lbs

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.
2. Not available on seat sizes greater than 20" wide by 20" deep
3. For weights greater than 175lbs tilt while elevated is not available.

MODULES MUST SELECT SYSTEM TO ADD MODULES

PRETILT 5° Fixed Pre-Tilt⁽¹⁾N/C

ERSSP Electronic Recline Start & Stop Program Function⁽²⁾ \$500

PRM8 8° Precline Module⁽³⁾ only available with ReclineN/C

PRM30 30° Precline Module⁽³⁾ only available with Recline \$500

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height
2. Requires Multi-Function M610 through Drive Control (MFCBD). Set range between 90° and 170°.
3. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available.

ADULT SEAT WIDTH SETTINGS

SW116 16" Wide (Adjustable from 16" to 20") N/C

SW117 17" Wide (Adjustable from 16" to 20") N/C

SW118 18" Wide (Adjustable from 16" to 20") N/C

SW119 19" Wide (Adjustable from 16" to 20") N/C

SW120 20" Wide (Adjustable from 16" to 20") N/C

W219 19" Wide (Adjustable from 19" to 22") \$460

W220 20" Wide (Adjustable from 19" to 22") \$460

W221 21" Wide (Adjustable from 19" to 22") \$460

W222 22" Wide (Adjustable from 19" to 22") \$460

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

SD116 16" Deep (Adjustable from 16" to 21") N/C

SD117 17" Deep (Adjustable from 16" to 21") N/C

SD118 18" Deep (Adjustable from 16" to 21") N/C

SD119 19" Deep (Adjustable from 16" to 21") N/C

SD120 20" Deep (Adjustable from 16" to 21") N/C

D219 19" Deep (Adjustable from 19" to 22") \$460

D220 20" Deep (Adjustable from 19" to 22") \$460

D221 21" Deep (Adjustable from 19" to 22") \$460

D222 22" Deep (Adjustable from 19" to 22") \$460

D223 23" Deep (Adjustable from 19" to 22") \$460

SMALL ADULT SEAT WIDTH SETTINGS⁽¹⁾

W314 14" Wide (Adjustable from 14" to 15") \$250

W315 15" Wide (Adjustable from 14" to 15") \$250

SMALL ADULT SEAT DEPTH SETTINGS⁽¹⁾

D314 14" Deep (Adjustable from 14" to 18") \$250

D315 15" Deep (Adjustable from 14" to 18") \$250

D316 16" Deep (Adjustable from 14" to 18") \$250

D317 17" Deep (Adjustable from 14" to 18") \$250

D318 18" Deep (Adjustable from 14" to 18") \$250

1. Small Adult Seat Size - weight limit is 175 lbs on these seat sizes and if recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

MATRIX SEATING OPTIONS

Matrix Cushions available on page 10

Matrix Backs available on page 11

Matrix Headrests available on page 12



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CHEST STRAPS

AADL SEATING BENEFIT

- MFUS BP Monoflex, Underarm RLS S(3" x 15.5") \$195
- MFUM BP Monoflex, Underarm RLS M(4" x 18.25") \$195
- MFUL BP Monoflex, Underarm RLS L(4.5" x 21.5") \$195
- MFCS BP Monoflex, Center RLS S(2.5" x 15.5") \$195
- MFCS BP Monoflex, Center RLS M(3.25" x 18.25") \$195
- MFCL BP Monoflex, Center RLS L(3.5" x 21.5") \$195
- PFLEX-M PostureFlex Abdominal Support
Height 6.5", length 16" to 23" \$295
- PFLEX-L PostureFlex Abdominal Support
Height 7", length 20" to 32" \$295
- QRTAB Quad release tab for easy independent
application and removal of PostureFlex \$35

LAP BELTS

- BELT60 Push Button Style Lap Belt, Length 60"⁽¹⁾ N/C
- BELT71 Push Button Style Lap Belt, Length 71" ⁽¹⁾ N/C
- PLB2 Bodypoint Padded Lap Belt 2 Point⁽²⁾ \$250
- PLB4 Bodypoint Padded Lap Belt 4 Point⁽²⁾ \$315
- TROBELT Safety Belt for use only with TRRO & TRRO-E N/C

1. Not available with TRRO and TRRO-E
2. 67" length only option

MOTION CONCEPTS HEADREST OPTIONS

- MCSH Motion Concepts Standard⁽¹⁾ N/C
- MCOH Motion Concepts Auto Style⁽¹⁾ \$235
- FM Fixed Mounting (non-removable)⁽²⁾ N/C
- OHR Omit: Headrest N/C

1. Must select Fixed Mounting (FM) or Matrix Headrest Hardware on page 12.
2. Only compatible with Standard Rehab, Contour or High Back Models. All other upgradeable backs require removeable hardware.

MATRIX HEADRESTS OPTIONS

- Matrix Elan Headrests and Hardware
available on page 12



FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ FOR TILT ONLY/ELEVATE ONLY

- 18" FBT18
- 19" FBT19
- 20" FBT20
- 21" FBT21
- 22" FBT22
- 23" FBT23
- 24" FBT24
- 25" FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ FOR RECLINE OR TILT/RECLINE

- 20" FBR20
- 21" FBR21
- 22" FBR22
- 23" FBR23
- 24" FBR24
- 25" FBR25
- 26" FBR26
- 27" FBR27

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK CANES (NON RECLINE SYSTEMS)

- SBC18 18" Straight Back Canes STD
- SBC20 20" Straight Back Canes STD
- SBC22 22" Straight Back Canes STD
- SBC24 24" Straight Back Canes STD

Select Cane Angle

- 75° BA75
- 79° BA79
- 85° BA85
- 90° BA90
- 95° BA95
- 101° BA101
- 106° BA106
- 112° BA112
- 116° BA116

- ABC20 20" Angled (5°) Back Canes N/C
- ABC22 22" Angled (5°) Back Canes N/C
- ABC24 24" Angled (5°) Back Canes N/C

Select Cane Angle

- 81° BA81
- 84° BA84
- 91° BA9
- 96° BA96
- 101° BA101
- 106° BA106
- 111° BA111
- 117° BA117
- 121° BA121

- EBC22 22" Endomorph Back Canes⁽¹⁾ \$250
- EBC24 24" Endomorph Back Canes⁽¹⁾ \$250
- EBCN1 1" Narrower N/C
- EBCN2 2" Narrower N/C
- EBCN4 4" Narrower N/C

1. Comes standard with Cantilever Armrest Cane Mounted.

CANE OPTIONS FOR RECLINE SYSTEMS

- BFS2 Recline Body Fit 2" Narrower⁽¹⁾ \$995
- BFS4 Recline Body Fit 4" Narrower⁽¹⁾ \$995
- AHH Adjustable Hooking Handles⁽²⁾ (pair) \$450
- RPMS Push Handles for Recline - Straight⁽³⁾ (pair) \$350
- RPMA Push Handles for Recline - Angled⁽³⁾ (pair) \$350

1. Minimum finished back widths of 16".
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

BACK OPTIONS

- SRB Standard Rehab Back & Cushion \$125
- PIPF Flat Planar Interface Plate⁽¹⁾ \$125
- PIPR Recessed Planar Interface Plate⁽¹⁾⁽²⁾ \$125
- ORBC Omit Rigid Back & Cushion⁽³⁾ N/C

1. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
2. Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
3. Not available with Recline/ESR or Recline/PSB

MATRIX BACKREST OPTIONS

AADL SEATING BENEFIT

- Matrix Backs available on page 11



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ARMRESTS (NON RECLINE SYSTEMS)

- DP1 Adult Dual Post Adj. 9.75"-12.75"STD
- DP2 Adult Dual Post Adj. 12.75"-15.75"STD
- DP3 Adult Dual Post Adj. 15.75"-18.75"STD
- Select: FP Flat PouchSTD
- GB Glove BoxesN/C
- OAL ADD: OutBack Arm Option Left⁽¹⁾ each\$195
- OAR ADD: OutBack Arm Option Right⁽¹⁾ each\$195
- RMCAN Ultra Rail Seat Mnt Flip Back Cantilever
MAXX Tilt Arm (9-12.5")⁽²⁾\$488
- DPM ADD: Dual Post Module to RMCAN (pair)\$250
- OCTA1 Outback Cantilever Tilt Armrest 9"-12"⁽³⁾\$488
- OCTA2 Outback Cantilever Titl Armrest 12"-15"⁽³⁾\$488
- CAFEB Cantilver Arms for Endomorph
Back Canes (pair) (see page 4)⁽³⁾STD

1. Only available with Dual Post Armrests.
2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.
3. Not available with Small Adult Sizes.

RECLINE ARMRESTS

- RC1 Reclining Adjustable Height. 9.5"-13"STD
- RC2 Reclining Adjustable Height. 13"-16"STD
- OAL ADD: OutBack Arm Option Left⁽¹⁾ each\$195
- OAR ADD: OutBack Arm Option Right⁽¹⁾ each\$195
- MCAN1 Maxx Style Cane Mounted Height, Width
and Angle Adjustable Cantilever Flip Back
Armrests for Recline Systems 9.5"-13"\$488
- MCAN2 Maxx Style Cane Mounted Height, Width
and Angle Adjustable Cantilever Flip Back
Armrests for Recline Systems 12.5"-16"\$488
- OCRA1 Outback Cantilever Recline Armrest 9"-12"⁽²⁾\$488
- OCRA2 Outback Cantilever Recline Armrest 12"-15"⁽²⁾\$488
- CANQ1 Cane Mounted Height & Angle Adjustable
Cantilver Flip Back Armrest 9.5"-12"⁽²⁾\$488
- CANQ2 Cane Mounted Height & Angle Adjustable
Cantilver Flip Back Armrest 10"-13"⁽²⁾\$488
- CANQ3 Cane Mounted Height & Angle Adjustable
Cantilever Flip Back Armrest 13"-16"⁽²⁾\$488

1. Only available with Reclining Armrests.
2. Not available with Small Adult Sizes.

ARMREST TUBE LENGTH - PLEASE SPECIFY

- ATFL Left Full STD
- ATDL Left Desk STD
- ATFR Right Full STD
- ATDR Right Desk STD

ARMPADS

- Standard**
- SAPFL Left Full (2.25" x 14") STD
 - SAPDL Left Desk (2.25" x 10") STD
 - SAPFR Right Full (2.25" x 14") STD
 - SAPDR Right Desk (2.25" x 10") STD

Waterfall Arm Pads⁽¹⁾

- WAPFL Left Full (13")N/C
- WAPDL Left Desk (9")N/C
- WAPFR Right Full (13")N/C
- WAPDR Right Desk (9")N/C

Modular Arm Pad - Left⁽²⁾

AADL SEATING BENEFIT

- MAPTFL Left Tray Full (14")\$95
- MAPTDL Left Tray Desk (10")\$95

Select Pad Insert below:

- ISSFL Iskin - Integrated soft skin foam
- SCVFL Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

AADL SEATING BENEFIT

- MAPTR Right Tray Full (14")\$95
- MAPTRD Right Tray Desk (10")\$95

Select Pad Insert below:

- ISSFR Iskin - Integrated soft skin foam
- SCVFR Startex Covered Visco foam

Flat Multi Position Arm Pads

CLIENT JUSTIFICATION

- MPAFL Left Full (14")\$150
- MPADL Left Desk (10")\$150
- MPAFR Right Full (14")\$150
- MPADR Right Desk (10")\$150

Gel Pads

AADL SEATING BENEFIT

- GAPNL Narrow Left (2" x 12")\$88
- GAPNR Narrow Right (2" x 12")\$88
- GAPWL Wide Left (3.5" x 12")\$96
- GAPWR Wide Right (3.5" x 12")\$96
- GAPLL Long Left (3.5" x 14")\$103
- GAPLR Long Right (3.5" x 14")\$103

Arm Troughs

AADL SEATING BENEFIT

- EATL Ergonomic Arm Trough - Left\$121
- EATR Ergonomic Arm Trough - Right\$121
- FHPL Flat Hand Pad - Left⁽³⁾\$91
- FHPR Flat Hand Pad - Right⁽³⁾\$91

Arm Pad Storage/Cup Holder⁽⁴⁾

CLIENT COST

- APSL Left\$295
- APSR Right\$295
- OMIT Omit Armpads(7/8" diameter tube)N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ELBOW BLOCKS

Elbow Blocks available on page 14

MANUAL CENTER MOUNT FRONT RIGGING

MSFCMFP Maxx Style Fixed Center Mount Foot Platform⁽¹⁾⁽²⁾\$395
With Quick Release Height Adjustment

IFP Comes with a rubber coated footplate (11.5" W x 10")N/C

ICP Individual calf pads\$200

Select seat pan to footplate:

SPTF1 9" - 13"N/C


SPTF14 13" -19"N/C

Set at:

SA70 Set at 70°N/C

SA90 Set at 90°N/C

SA97 Set at 97°N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
2. 3" ground clearance required.

POWER CENTER MOUNT FOOT PLATFORM

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression: (Add)	
Seat Pan to Footplate: (Total)	

LNX Power Center Mount Foot Platform

CLIENT COST

LNX 9" Articulation⁽¹⁾⁽²⁾ \$3,200

LSPTF1 10" - 14" Seat pan to footplate

LSPTF2 13" - 19" Seat pan to footplate

LNXA 6" Articulation⁽¹⁾⁽²⁾⁽⁶⁾ \$3,200

LASPTF2 10" - 14" Seat pan to footplate

Options:

LIFP Rubber Coated Footplate 11.5"W x 10" D STD

LICP Individual calf pads STD

OPCP One Piece Malleable Calf Panel⁽³⁾..... \$250

LNXPF Plantar Flex Option for LNX Foot platform⁽⁴⁾ \$250


Latitude Style Foot Platform

LSFP Latitude Style Foot Platform⁽⁵⁾ \$2,995
For finished knee to heel measurements of 6.5"-12"

LSPFP Latitude Style Padded Foot Plate \$250

- 3" ground clearance required.
- Available on Adult Systems 17" deep and greater and on Small Adult Sizes 15" deep and greater. For shorter seat depths call Customer Service for other options.
- Not available with Plantar Flex Option for LNX Foot platform, or with LNXA.
- Only available with Foot Platform Options
- Not available with TRRO and TRRO-E.
- For clients under 4' only.

FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP



SFFP Small Foot Platform⁽¹⁾(12"w x 8.5"d) \$180

MFFP Medium Platform⁽¹⁾(12"w x 10.75"d) \$180


LFFP Large Platform⁽¹⁾(14.5"w x 13"d) \$250

XLFFP Extra Large Platform⁽¹⁾(17"w x 13"d) \$250

XXLFFP Extra Extra Large Platform⁽¹⁾(21"w x 13"d) \$250

1. May need to adjust leg position to avoid interference with front castors

INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP



IFSL Small⁽¹⁾ - (5" x 7.5" deep) - Left \$165

IFSR Small⁽¹⁾ - (5" x 7.5" deep) - Right \$165

IFML Medium⁽¹⁾ - (5.5" x 9.5" deep) - Left \$165

IFMR Medium⁽¹⁾ - (5.5" x 9.5" deep) - Right \$165

IFLL Large⁽¹⁾ - (6" x 11.5" deep) - Left \$165

IFLR Large⁽¹⁾ - (6" x 11.5" deep) - Right \$165

MFOL Multi-Axis Optional Footplate Left (6"W x 10"D) \$225

MFOR Multi-Axis Optional Footplate Right (6"W x 10"D) \$225

1. Foot plates may interfere with casters in certain configurations

INVACARE CENTER MOUNT

Select seat pan to footplate:

CTMT-FJ 9.25" - 12" \$595

CTMT-F 12.25" -15" \$595

Select foot plates/footboard

SM Individual Foot Plates (5.25"W x 7.5"D) \$165

MDM Individual Foot Plates (5.25"W x 10"D) \$165

LRG Individual Foot Plates (5.9"W x 12"D) \$165

CTMTFB Footboard \$165

Optional:

CALF Invacare Calf Pads for use with CTMT⁽²⁾ \$268

- 3" ground clearance required.
- Not available with CTMT-FJ

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MAXX STYLE POWER ELEVATING LEGRESTS

CLIENT COST

- MSPPL MAXX Swing Away Power Elevating⁽¹⁾⁽²⁾ \$2,490
- Select Power Legrest Function:
 - INDA Individual Movement N/C
 - COMA Combined Movement..... N/C
 - SYNA Synchronized/Recline Movement N/C
- Select Receivers:
 - RS Straight Receivers..... N/C
 - RF Flared 1" Outward Receivers N/C
- Select Seat Pan to Footplate:
 - SPTF5 Length: 10"-14" N/C
 - SPTF2 Length: 13"-17" N/C
 - SPTF7 Length: 16"-20" N/C
 - FCPMS Flat Calf Pad Option (pair) \$275
 - CCPMS Curved Calf Pad Option⁽³⁾ (pair) \$275

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads.
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

MAXX STYLE MANUAL ELEVATING LEGRESTS

- MMPPSA MAXX Swing Away Manual Elevating⁽¹⁾⁽²⁾ \$680
- Select Receivers:
 - RS Straight Receivers.....N/C
 - RF Flared 1" Outward Receivers.....N/C
- Select Seat Pan to Footplate:
 - SPTF5 Length: 10"-14"N/C
 - SPTF2 Length: 13"-17"N/C
 - SPTF7 Length: 16"-20"N/C
 - FCPMS Flat Calf Pad Option (pair) \$275
 - CCPMS Curved Calf Pad Option (inside with 6")⁽³⁾ (pair) \$275

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

SWING AWAY LEGREST

- HD70 Heavy Duty 70° Swing Away Footrest^{(pair)(1)(2)}\$495
- Select Receivers:
 - RS Straight ReceiversN/C
 - RF Flared 1" Outward Receivers.....N/C
- Select Seat Pan to Footplate:
 - SPTF11 Length: 6"-10"N/C
 - SPTF8 Length: 11"-15".....N/C
 - SPTF2 Length: 13"-17".....N/C
 - SPTF7 Length: 16"-20".....N/C
 - HDFCP Flat Calf Pad Option (pair) \$275
 - HD70CCP Curved Calf Pad Option⁽³⁾ (inside with 6") (pair) \$275
 - DCS Dynamic Coil Springs (Pair) for Heavy Duty 70° Swingaway Footrests \$190
 - I70SAF Invacare 70° Swingaway Footrests (Pair)\$325
 - I70-CF Comes with composite foot plates..... N/C
 - 70HDAIH Invacare 70° Non-Tapered HD Swingaway Footrests with AT5543 Adjustable Angle Foot plates, Heel Loops and Impact Guards (Pair) \$457
 - 70HTAPAIH Invacare Tapered 70° Swingaway Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards(Pair) \$497
 - ALPT Longer Pivot and Slide Tube for 70° Swingaways (Pair) \$132
 - RLSL Residual Limb Support - Left⁽⁴⁾ (6.75"W x 9.75"L x 2"H) (each) \$495
 - RLSR Residual Limb Support - Right⁽⁴⁾ (6.75"W x 9.75"L x 2"H) (each) \$495

1. Must choose footplate on page 8
2. Heavy duty durable construction (4.5lbs ea).
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.
4. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.


Dealer Name: _____

Dealer Account Number: _____

PO #: _____


Tag: _____

FOOTPLATE FOR SWING AWAY/ELEVATING LEGS



1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

FOOTPLATE/CUT OUTS SIZE SELECTION



1. Certain footplate sizes may not fit all seat widths depending on legrests selected

ADJUSTABLE ANGLE FOOTPLATES

PRIOR APPROVAL

- AFSL Left Small Footplate (4.25"W x 8"D) \$165
- AFSR Right Small Footplate (4.25"W x 8"D) \$165
- AFML Left Medium Footplate (5.25"W x 8"D) \$165
- AFMR Right Medium Footplate (5.25"W x 8"D) \$165
- AFLL Left Large Footplate (6.25"W x 8"D) \$165
- AFLR Right Large Footplate (6.25"W x 8"D) \$165
- AFXL Left Extra Large Footplate (9"W x 11"D) \$165
- AFXR Right Extra Large Footplate (9"W x 11"D) \$165

MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES

PRIOR APPROVAL

- MFSL Left Small Footplate (4.25"W x 8"D) \$225
- MFSR Right Small Footplate (4.25"W x 8"D) \$225
- MFML Left Medium Footplate (5.25"W x 8"D) \$225
- MFMR Right Medium Footplate (5.25"W x 8"D) \$225
- MFLL Left Large Footplate (6.25"W x 8"D) \$225
- MFRL Right Large Footplate (6.25"W x 8"D) \$225
- MFXL Left Extra Large Footplate (9"W x 11"D) \$225
- MFXR Right Extra Large Footplate (9"W x 11"D) \$225

HEEL LOOPS, CUPS AND TOE STRAPS

- HLL Left Heel Loop (each) \$20
- HLR Right Heel Loop (each) \$20
- HCL Left Cup Loop (each) \$20
- HCR Right Cup Loop (each) \$20
- TSL Left Toe Strap (each) \$20
- TSR Right Toe Strap (each) \$20

ADJUSTABLE ANGLE FOOTPLATES/CUT OUTS

PRIOR APPROVAL

- AFCSL Left Small Footplate (7.25"W x 9"D) \$165
- AFCSR Right Small Footplate (7.25"W x 9"D) \$165
- AFCML Left Medium Footplate (7.5"W x 11"D) \$165
- AFCMR Right Medium Footplate (7.5"W x 11"D) \$165
- AFCLL Left Large Footplate (10"W x 14"D) \$165
- AFCLR Right Large Footplate (10"W x 14"D) \$165

MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES/ CUT OUTS

PRIOR APPROVAL

- MFCSL Left Small Footplate (7.25"W x 9"D) \$225
- MFCSR Right Small Footplate (7.25"W x 9"D) \$225
- MFCML Left Medium Footplate (7.5"W x 11"D) \$225
- MFCMR Right Medium Footplate (7.5"W x 11"D) \$225
- MFCLL Left Large Footplate (10"W x 14"D) \$225
- MFCLR Right Large Footplate (10"W x 14"D) \$225

CALF STRAPS AND CALF PANELS

- CSM Body Point Airmesh Calf Strap for Chairs
15"-18" wide with 3.5" H Pad⁽¹⁾ \$200
- CSL Body Point Airmesh Calf Strap for Chairs
18"-22" wide with 3.5" H Pad⁽¹⁾ \$200
- CPM Body Point Airmesh Calf Panel for Chairs
15"-18" wide with 35.5" W x 9" H Pad⁽¹⁾ \$295
- CPL Body Point Airmesh Calf Panel for Chairs
18"-22" wide with 38.5" W x 9" H Pad⁽¹⁾ \$315

THORACIC LATERALS AND HIP/KNEE SUPPORTS

- Matrx Laterals please see page 12
- Maxx Laterals and Hip/Knee supports
please see page 13

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ELECTRONIC SWITCH MOUNTING

- ESML Left Side (Armpad)
- ESMR Right Side (Armpad)
- ESMJL Joystick Left (Mounts to Joystick Bracket)
- ESMJR Joystick Right (Mounts to Joystick Bracket)
- ESMD Do not mount

SINGLE-DUAL FUNCTION ELECTRONICS

- SFCB Single/Dual Function Electronics⁽¹⁾..... STD

- Select Switch:
- SPBO Single Push Button
 - STO Single Toggle
 - DPBO Dual Push Button
 - SOS Omit Switch
 - JOY Through Joystick
 - SC Through Specialty Control

- Optional:
- RDE Reduced Drive for Elevating Seat⁽²⁾ CLIENT COST..... \$400

- SPBSA Additional Single Push Button..... \$150
- STA Additional Single Toggle Switch..... \$150
- DPBA Additional Dual Push Button..... \$150

1. These electronics include drive lock-out and tilt limit.
2. If reduced drive is not selected the system will be in drive lock-out when elevated.

MULTI-FUNCTION ELECTRONICS

- MFCB Multi-Function through Switch/Toggle⁽¹⁾ \$1,950

- Select Switch:
- QPBO Quad Push Button
 - FWTO Four-Way Toggle
 - EWRO 8-Way Rocker⁽²⁾

- MFCBD Multi-Function through Drive Control⁽¹⁾..... \$1,950

- Select Control:
- JOY Through Joystick
 - SC Through Specialty Control

- Optional:
- RDY Include Reduced Drive for Elevating Seat⁽³⁾
 - RDN Do not Induce Reduced Drive for Elevating Seat⁽³⁾

- QPBA Additional Quad Button..... \$150
- FWTA Additional Four-Way Toggle Switch \$150
- EWRA Additional 8-Way Rocker⁽²⁾ \$205
- 10WA Additional 10-Way Switch..... \$255

1. These electronics include seating attendant control, drive lock-out, tilt/recline limit and actuator speed control
2. Dual function is not available with 8-Way Rocker.
3. If reduced drive is not selected system will be in drive lock-out when elevated.

ELECTRONICS OPTIONS

- M616S M616 Converter⁽¹⁾..... \$125
- FTSQ Feather Touch Switch Single Function Quarter Size \$295
- FTMD Feather Touch Switch Multi-Function Dual Dime Size \$295
- FTMQ Feather Touch Switch Multi-Function Quad Dime Size \$295
- ASL320R ASL Pulse Mechanical Switch Red \$125
- ASL320G ASL Pulse Mechanical Switch Green \$125
- ASL320BK ASL Pulse Mechanical Switch Black..... \$125
- SESG Stealth Egg Switch - Single Function Green..... \$135
- SESB Stealth Egg Switch - Single Function Black \$135
- PSOS Power Shut Off Switch - Toggle⁽¹⁾..... \$175
- LJS Lazarus Jump Starter⁽¹⁾ \$150
- APPS Accessory Port Power Supply⁽¹⁾ (3 amp capacity) \$195
- TPMSMH 12 Piece Switch Mounting Hardware Kit (additional kits)..... \$39

1. Not available with Single-Dual Function Electronics.

VENT TRAYS, O2 HOLDERS, ACCESSORIES

- For Vent Trays, O2 Holders, IV Pole and Additional Accessories
Please see page 14

INSTALLATION

- MCI Motion Concepts Installation..... Included

PLEASE SPECIFY PAGES INCLUDED WITH ORDER

- Pages 1-9 Excludes Matrx and Maxx Seating Accessories Thoraic Supports, Hip Supports, Vent Trays and Accessories
- Pages 1-14 Complete Order Form

Total Retail Value \$ _____

SUBMIT

PRINT FORM **SAVE** **CLEAR FORM**



AADL SEATING BENEFIT

CUSHION COVER REFERENCE GUIDE

Cushion	Meshtex	Startex	Infection	Spacetex	Onyx
Libra/Adjust	N/A	STD	\$50	N/C	N/A
PS/PSP	N/A	STD-Reversible	STD-Reversible	N/A	N/A
Vi	STD	\$103	\$103	N/A	\$103

MATRIX SEAT CUSHION OPTIONS

Matrix Vi Seat Cushion 14" - 20" wide \$624

- MA1415-VI MA1417-VI MA1616-VI MA1618-VI
- MA1620-VI MA1716-VI MA1718-VI MA1720-VI
- MA1816-VI MA1818-VI MA1820-VI

Matrix Vi Seat Cushion 21" - 22" wide \$681

- MA2016-VI MA2018-VI MA2020-VI MA2218-VI
- MA2220-VI

Matrix Vi HD Seat Cushion 22" wide \$840

- MA2222-VI

Matrix Multi Cushion 14" - 20" wide/reversible cover

HCPCS code E2607 \$745

- MM1616 MM1618 MM1620 MM1716 MM1718
- MM1720 MM1816 MM1818 MM1820 MM2016
- MM2018 MM2020

MATRIX SEAT CUSHION OPTIONS

Matrix Libra Adjust Seat Cushion 16" -20" wide \$985

- LC1416-ADJ LC1516-ADJ LC1518-ADJ LC1520-ADJ
- LC1616-ADJ LC1618-ADJ LC1620-ADJ LC1716-ADJ
- LC1718-ADJ LC1720-ADJ LC1816-ADJ LC1818-ADJ
- LC1820-ADJ LC1916-ADJ LC1918-ADJ LC1920-ADJ
- LC2016-ADJ LC2018-ADJ LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" depth \$1,265

- LC1622-ADJ LC1722-ADJ LC1822-ADJ
- LC1922-ADJ LC2022-ADJ

Matrix Libra Adjust Seat Cushion 21" - 22" wide \$1,265

- LC2118-ADJ LC2120-ADJ LC2122-ADJ
- LC2218-ADJ LC2220-ADJ LC2222-ADJ

Matrix Libra Seat Cushion 16" - 20" wide \$745

- LC1416 LC1516 LC1518 LC1520 LC1616
- LC1618 LC1620 LC1716 LC1718 LC1720
- LC1816 LC1818 LC1820 LC1916 LC1918
- LC1920 LC2016 LC2018 LC2020

Matrix Libra Seat Cushion 22" deep \$1,025

- LC1622 LC1722 LC1822 LC1922 LC2022

Matrix Libra HD Seat Cushion 21" -22" wide \$1,025

- LC2118 LC2120 LC2122 LC2218 LC2220

Matrix PS Seat Cushion 16" - 20" wide \$569

- PS1416 PS1616 PS1618 PS1620 PS1716
- PS1718 PS1720 PS1816 PS1818 PS1820
- PS2016 PS2018 PS2020

Matrix PS HD Seat Cushion 22" - 22" wide \$975

- PS2218 PS2220 PS2222

Matrix PSP Seat Cushion 14" - 20" wide \$645

- PSP1416 PSP1616 PSP1618 PSP1620
- PSP1716 PSP1718 PSP1720 PSP1816
- PSP1818 PSP1820 PSP2016 PSP2018
- PSP2020

Matrix PSP HD Seat Cushion 22" - 22" wide \$995

- PSP2218 PSP2220 PSP2222

FABRIC OPTIONS FOR CUSHIONS

- SFLF Spacetex Fabric Option for
Libra Cushion N/C
- ICFL Infection Control Fabric Upgrade for
Libra Cushion \$50
- OCUMA Onyx Fabric Upgrade for
Matrix Vi Cushion \$103
- SFVS Smooth side Startex Fabric Upgrade for
Matrix Vi Cushion \$103
- SFVF Fabric side Startex Fabric Upgrade for
Matrix Vi Cushion \$103

BACK COVER REFERENCE GUIDE

Back	Meshtex	Infection	Onyx
High	STD	N/A	\$103
Elite/ Elite HD/Elite Deep	STD	\$103	N/A
Elite TR/ Elite TR HD	STD	\$103	N/A
PB/PB HD/PB Deep/PB HD Deep	STD	\$103	N/A

FABRIC OPTIONS FOR OPTIONS

- EBSS Super Soft Foam for Elite Back Upgrade \$103
- OCUHB Onyx Fabric Upgrade for High Backs \$103
- WSCU Infection Control Cover for Elite/PB Backs \$103
purchased with back
- WSC Infection Control Cover for Elite/PB Backs \$150
purchased with without back
- CHIEB Chest Harness Interface Kit for Elite Back⁽¹⁾⁽²⁾ \$103
- CHIPB Chest Harness Interface Kit for PB Back⁽¹⁾⁽²⁾ \$52

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harness



AADL SEATING BENEFIT

UPGRADABLE BACK OPTIONS

High Back 16" - 22" wide \$595

- | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> HB1616 | <input type="checkbox"/> HB1618 | <input type="checkbox"/> HB1620 | <input type="checkbox"/> HB1716 |
| <input type="checkbox"/> HB1718 | <input type="checkbox"/> HB1720 | <input type="checkbox"/> HB1816 | <input type="checkbox"/> HB1818 |
| <input type="checkbox"/> HB1820 | <input type="checkbox"/> HB1916 | <input type="checkbox"/> HB1918 | <input type="checkbox"/> HB1920 |
| <input type="checkbox"/> HB2016 | <input type="checkbox"/> HB2018 | <input type="checkbox"/> HB2020 | <input type="checkbox"/> HB2116 |
| <input type="checkbox"/> HB2118 | <input type="checkbox"/> HB2120 | <input type="checkbox"/> HB2216 | <input type="checkbox"/> HB2218 |
| <input type="checkbox"/> HB2220 | | | |

Matrix Elite Back 14" - 20" wide (3" Contour) \$766

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PBE1414 | <input type="checkbox"/> PBE1416 | <input type="checkbox"/> PBE1418 | <input type="checkbox"/> PBE1420 |
| <input type="checkbox"/> PBE1514 | <input type="checkbox"/> PBE1516 | <input type="checkbox"/> PBE1518 | <input type="checkbox"/> PBE1520 |
| <input type="checkbox"/> PBE1614 | <input type="checkbox"/> PBE1616 | <input type="checkbox"/> PBE1618 | <input type="checkbox"/> PBE1620 |
| <input type="checkbox"/> PBE1714 | <input type="checkbox"/> PBE1716 | <input type="checkbox"/> PBE1718 | <input type="checkbox"/> PBE1720 |
| <input type="checkbox"/> PBE1814 | <input type="checkbox"/> PBE1816 | <input type="checkbox"/> PBE1818 | <input type="checkbox"/> PBE1820 |
| <input type="checkbox"/> PBE1916 | <input type="checkbox"/> PBE1918 | <input type="checkbox"/> PBE1920 | <input type="checkbox"/> PBE2014 |
| <input type="checkbox"/> PBE2016 | <input type="checkbox"/> PBE2018 | <input type="checkbox"/> PBE2020 | |

Matrix Elite HD Back 21" - 22" wide (3" Contour) \$1,095

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> PBE2116-HD | <input type="checkbox"/> PBE2120-HD | <input type="checkbox"/> PBE2216-HD |
| <input type="checkbox"/> PBE2220-HD | | |

Matrix Elite Deep Back 14" - 20" wide (3" Contour) \$918

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> EDB1416 | <input type="checkbox"/> EDB1418 | <input type="checkbox"/> EDB1420 | <input type="checkbox"/> EDB1514 |
| <input type="checkbox"/> EDB1516 | <input type="checkbox"/> EDB1518 | <input type="checkbox"/> EDB1520 | <input type="checkbox"/> EDB1614 |
| <input type="checkbox"/> EDB1616 | <input type="checkbox"/> EDB1618 | <input type="checkbox"/> EDB1620 | <input type="checkbox"/> EDB1714 |
| <input type="checkbox"/> EDB1716 | <input type="checkbox"/> EDB1718 | <input type="checkbox"/> EDB1720 | <input type="checkbox"/> EDB1814 |
| <input type="checkbox"/> EDB1816 | <input type="checkbox"/> EDB1818 | <input type="checkbox"/> EDB1820 | <input type="checkbox"/> EDB1914 |
| <input type="checkbox"/> EDB1916 | <input type="checkbox"/> EDB1918 | <input type="checkbox"/> EDB1920 | <input type="checkbox"/> EDB2014 |
| <input type="checkbox"/> EDB2016 | <input type="checkbox"/> EDB2018 | <input type="checkbox"/> EDB2020 | |

Matrix Elite Deep HD Back 21" - 22" wide (6" Contour) \$1,095

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> EDB2116-HD | <input type="checkbox"/> EDB2120-HD | <input type="checkbox"/> EDB2216-HD |
| <input type="checkbox"/> EDB2220-HD | | |

Matrix Elite E2 Back 14" - 20" wide/14"/16"H(3" Contour) \$965

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2S1414 | <input type="checkbox"/> E2S1416 | <input type="checkbox"/> E2S1514 | <input type="checkbox"/> E2S1516 |
| <input type="checkbox"/> E2S1614 | <input type="checkbox"/> E2S1616 | <input type="checkbox"/> E2S1714 | <input type="checkbox"/> E2S1716 |
| <input type="checkbox"/> E2S1814 | <input type="checkbox"/> E2S1816 | <input type="checkbox"/> E2S1914 | <input type="checkbox"/> E2S1916 |
| <input type="checkbox"/> E2S2014 | <input type="checkbox"/> E2S2016 | | |

Matrix Elite E2 Back 14" - 20" wide/18"/20"H(3" Contour) \$995

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2S1418 | <input type="checkbox"/> E2S1420 | <input type="checkbox"/> E2S1518 | <input type="checkbox"/> E2S1520 |
| <input type="checkbox"/> E2S1618 | <input type="checkbox"/> E2S1620 | <input type="checkbox"/> E2S1718 | <input type="checkbox"/> E2S1720 |
| <input type="checkbox"/> E2S1818 | <input type="checkbox"/> E2S1820 | <input type="checkbox"/> E2S1918 | <input type="checkbox"/> E2S1920 |
| <input type="checkbox"/> E2S2018 | <input type="checkbox"/> E2S2020 | | |

Matrix Elite E2 Deep Back 14" - 20" wide/14"/16"H(6" Contour) \$1,125

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2D1414 | <input type="checkbox"/> E2D1416 | <input type="checkbox"/> E2D1514 | <input type="checkbox"/> E2D1516 |
| <input type="checkbox"/> E2D1614 | <input type="checkbox"/> E2D1616 | <input type="checkbox"/> E2D1714 | <input type="checkbox"/> E2D1716 |
| <input type="checkbox"/> E2D1814 | <input type="checkbox"/> E2D1816 | <input type="checkbox"/> E2D1914 | <input type="checkbox"/> E2D1916 |
| <input type="checkbox"/> E2D2014 | <input type="checkbox"/> E2D2016 | | |

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

UPGRADABLE BACK OPTION

Matrix Elite E2 Deep Back 14" - 20" wide/18"/20"H (6" Contour) \$1,155

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2D1418 | <input type="checkbox"/> E2D1420 | <input type="checkbox"/> E2D1518 | <input type="checkbox"/> E2D1520 |
| <input type="checkbox"/> E2D1618 | <input type="checkbox"/> E2D1620 | <input type="checkbox"/> E2D1718 | <input type="checkbox"/> E2D1720 |
| <input type="checkbox"/> E2D1818 | <input type="checkbox"/> E2D1820 | <input type="checkbox"/> E2D1918 | <input type="checkbox"/> E2D1920 |
| <input type="checkbox"/> E2D2018 | <input type="checkbox"/> E2D2020 | | |

Matrix Elite TR Back 15" - 20" wide (3" Contour) \$866

- | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> PBE1516-TR | <input type="checkbox"/> PBE1518-TR | <input type="checkbox"/> PBE1520-TR |
| <input type="checkbox"/> PBE1616-TR | <input type="checkbox"/> PBE1618-TR | <input type="checkbox"/> PBE1620-TR |
| <input type="checkbox"/> PBE1716-TR | <input type="checkbox"/> PBE1718-TR | <input type="checkbox"/> PBE1720-TR |
| <input type="checkbox"/> PBE1816-TR | <input type="checkbox"/> PBE1818-TR | <input type="checkbox"/> PBE1820-TR |
| <input type="checkbox"/> PBE1916-TR | <input type="checkbox"/> PBE1918-TR | <input type="checkbox"/> PBE1920-TR |
| <input type="checkbox"/> PBE2016-TR | <input type="checkbox"/> PBE2018-TR | <input type="checkbox"/> PBE2020-TR |

Matrix Elite TR HD Back 21" - 22" wide (3" Contour) \$1,151

- | | |
|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> PBE2120-TRHD | <input type="checkbox"/> PBE2220-TRHD |
|---------------------------------------|---------------------------------------|

Matrix PB Posture Back 14" - 20" wide (5" Contour) \$734

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> MPB1412 | <input type="checkbox"/> MPB1612 | <input type="checkbox"/> MPB1616 | <input type="checkbox"/> MPB1620 |
| <input type="checkbox"/> MPB1712 | <input type="checkbox"/> MPB1716 | <input type="checkbox"/> MPB1720 | <input type="checkbox"/> MPB1812 |
| <input type="checkbox"/> MPB1816 | <input type="checkbox"/> MPB1820 | <input type="checkbox"/> MPB2012 | <input type="checkbox"/> MPB2016 |
| <input type="checkbox"/> MPB2020 | | | |

Matrix PB Posture Back Deep 14" - 20" wide (7" Contour) \$892

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PBD1412 | <input type="checkbox"/> PBD1612 | <input type="checkbox"/> PBD1616 | <input type="checkbox"/> PBD1620 |
| <input type="checkbox"/> PBD1712 | <input type="checkbox"/> PBD1716 | <input type="checkbox"/> PBD1720 | <input type="checkbox"/> PBD1812 |
| <input type="checkbox"/> PBD1816 | <input type="checkbox"/> PBD1820 | <input type="checkbox"/> PBD2012 | <input type="checkbox"/> PBD2016 |
| <input type="checkbox"/> PBD2020 | | | |

Matrix PB Heavy Duty Back 20" - 22" wide (5" Contour).....\$1,051

- | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> HDB2016 | <input type="checkbox"/> HDB2020 | <input type="checkbox"/> HDB2216 | <input type="checkbox"/> HDB2220 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrix PB Heavy Duty Deep Back 20"-22" wide (5" Contour)\$1,209

- | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> DHDB2016 | <input type="checkbox"/> DHDB2020 | <input type="checkbox"/> DHDB2216 | <input type="checkbox"/> DHDB2220 |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|



MOTION CONCEPTS SINGLE POST BACKREST

- SPBMI Single Post Angle Adjustable Fold Down Back interface-Matrix Back⁽¹⁾⁽³⁾⁽⁴⁾
(Folds to 25° forward) \$750
 Matrix Backs available on page 11
- SPBCI Single Post Angle Adjustable Fold Down Back for Contour Back⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾
(Folds to 25° forward) \$995
Select Back: CB1617 CB1717 CB1817
 CB1917 CB2017
- OCUCB Onyx Fabric Upgrade for Contour Backs \$100

1. Single Post Back is an angle adjustable back between 73° and 109° in 6° increments from 17° forward to 19° back. The back folds down and is only available for 16"-20" seat widths & 21"-24" back height. A 4" gap from bottom of back to seat pan is required.
2. Includes a Contour Back with fabric cover.
3. This option is not available with a vent tray and has a 250lb weight capacity.
4. Cane mounted armrest are not available with this option.
5. Not available with TRRO and TRRO-E brackets.

HEADREST MOUNTING HARDWARE OPTIONS

- LXSHW Loxx Standard Hardware 15" \$275
- LXMMH Loxx Mini Hardware 11" \$275
- MEHWRE Elan Hardware 15" \$285
- MEMH Elan Mini Hardware 11" \$285
- MEHW-EXT Elan Extended Hardware \$320
- FPLXSHW FLIPP Loxx Standard Hardware 15" \$705
- FPLXSHW-EX FLIPP Loxx Extended Hardware \$705
- FPLXMMH FLIPP Loxx Mini Hardware 11" \$705
- FPLXMMH-EX FLIPP Loxx Mini Extended Hardware \$705
- FPMEHW FLIPP Elan Hardware 15" \$715
- FPMMH FLIPP Elan Mini Hardware 11" \$715
- FPMEHW-EXT FLIPP Elan Extended Hardware \$780

Dealer Name: _____
 Dealer Account Number: _____
 PO #: _____
 Tag: _____

MATRIX ELAN HEADREST PAD⁽¹⁾

Standard Pads

- ESP6 Matrix Elan Standard Pad 6" \$310
- ESP10 Matrix Elan Standard Pad 10" \$310
- ESP14 Matrix Elan Standard Pad 14" \$310
- ESP6-HUG Matrix Elan Hug Standard Pad 6" \$310
- ESP10-HUG Matrix Elan Hug Standard Pad 10" \$310
- ESP14-HUG Matrix Elan Hug Standard Pad 14" \$310
- SNUG14 Matrix SNUG Standard Pad 14" \$625
- SNUG14-A Matrix SNUG Asymmetrical Pad 14" \$645
- SNUG18 Matrix SNUG Standard Pad 18" \$625
- SNUG18-A Matrix SNUG Asymmetrical Pad 18" \$645

Standard Pads with Infection Control Cover

- ESP6-IC Matrix Elan Infection Control Pad 6" \$335
- ESP10-IC Matrix Elan Infection Control Pad 10" \$335
- ESP14-IC Matrix Elan Infection Control Pad 14" \$335
- ESP6-IC-HUG Matrix Elan Hug Infection Control Pad 6" \$335
- ESP10-IC-HUG Matrix Elan Hug Infection Control Pad 10" \$335
- ESP14-IC-HUG Matrix Elan Hug Infection Control Pad 14" \$335
- SNUG14-IC Matrix SNUG Standard Infection Control Pad 14" \$665
- SNUG14-A-IC Matrix SNUG Asymmetrical Infection Control Pad 14" \$685
- SNUG18-IC Matrix SNUG Standard Infection Control Pad 18" \$665
- SNUG18-A-IC Matrix SNUG Asymmetrical Infection Control Pad 18" \$685

Occipital /4-Point Pads

- EOP9 Matrix Elan Occipital Pad 9" \$341
- EOP12 Matrix Elan Occipital Pad 12" \$341
- E4PTSM Matrix Elan 4-Point Pad 9"W x 7"H \$448
- E4POINT Matrix Elan 4-Point Pad 11"W x 10"H \$448

Occipital / 4-Point Pads with Infection Control Cover

- EOP9-IC Matrix Elan Occipital Pad 9" with Infection Control \$366
- EOP12-IC Matrix Elan Occipital Pad 12" with Infection Control \$366
- E4PTSM-IC Matrix Elan 4-Point Pad 9"W x 7"H with Infection Control \$473
- E4POINT-IC Matrix Elan 4-Point Pad 11"W x 10"H with Infection Control \$473

1. If mounting on PB Elite must order removable hardware

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS

AADL SEATING BENEFIT

<input type="checkbox"/> E2FLL	Matrx E2 Fixed Laterals-Left(6.5x5") ⁽¹⁾	\$195
<input type="checkbox"/> E2FLR	Matrx E2 Fixed Laterals-Right(6.5x5") ⁽¹⁾	\$195
<input type="checkbox"/> PBSSL	Matrx Small Lateral Fixed-Left(3x5.5")	\$113
<input type="checkbox"/> PBSSR	Matrx Small Lateral Fixed-Right(3x5.5")	\$113
<input type="checkbox"/> PBLSL	Matrx Standard Lateral Fixed-Left(7.5x5")	\$116
<input type="checkbox"/> PBLSR	Matrx Standard Lateral Fixed-Right (7.5x5")\$	116
<input type="checkbox"/> PBOLL	Matrx Offset Lateral Fixed Left \$125	
<input type="checkbox"/> PBOLR	Matrx Offset Lateral Fixed Right	\$125

1. Only available on the E2 backrests.

MATRIX SWING AWAY LATERALS - FORELITE AND ELITE TR BACKS

AADL SEATING BENEFIT

Matrx Lateral Trunk Support Left⁽¹⁾

<input type="checkbox"/> EPADXL	Size XS (3.25x4.25")	\$115
<input type="checkbox"/> EPADSL	Size S (3.5x5.25")	\$115
<input type="checkbox"/> EPADM L	Size M (3.75x5.75")	\$115
<input type="checkbox"/> EPADL L	Size L (4.25x6.5")	\$115

Matrx Lateral Trunk Support Right⁽¹⁾

<input type="checkbox"/> EPADXR	Size XS (3.25x4.25")	\$115
<input type="checkbox"/> EPADSR	Size S (3.5x5.25")	\$115
<input type="checkbox"/> EPADMR	Size M (3.75x5.75")	\$115
<input type="checkbox"/> EPADLR	Size L (4.25x6.5")	\$115

1. Must choose swing away hardware below

MATRIX SWING AWAY HARDWARE - FORELITE AND ELITE TR BACKS

AADL SEATING BENEFIT

<input type="checkbox"/> SAEHBL	Matrx Elite Swing Away Hardware Left	\$250
<input type="checkbox"/> SAHEBR	Matrx Elite Swing Away Hardware Right	\$250
<input type="checkbox"/> OSAHL	Matrx Offset Elite Swing Away Hardware Left	\$250
<input type="checkbox"/> OSAHR	Matrx Offset Elite Swing Away Hardware Right	\$250

MAXX LATERALS TRUNK SUPPORTS⁽¹⁾

AADL SEATING BENEFIT

Left Fixed Mount Lateral Trunk Support⁽¹⁾

Select Left Trunk Support Size

<input type="checkbox"/> LATXL	Size XS (3.25" x 4.25")	\$100
<input type="checkbox"/> LATSL	Size S (3.5" x 5.25")	\$100
<input type="checkbox"/> LATML	Size M (3.75" x 5.75")	\$100
<input type="checkbox"/> LATLL	Size L (4.25" x 6.5")	\$100

Right Fixed Mount Lateral Trunk Support⁽¹⁾

Select Right Trunk Support Size

<input type="checkbox"/> LATXR	Size XS (3.25" x 4.25")	\$100
<input type="checkbox"/> LATSR	Size S (3.5" x 5.25")	\$100
<input type="checkbox"/> LATMR	Size M (3.75" x 5.75")	\$100
<input type="checkbox"/> LATLR	Size L (4.25" x 6.5")	\$100

1. These laterals are not available with Matrx PB, PB Elite Deep and PB Deep Backs.

MAXX STYLE SWING-AWAY TRUNK SUPPORT HARDWARE⁽¹⁾

AADL SEATING BENEFIT

<input type="checkbox"/> SALATL	Left MAXX Style Swing Away HDW	\$195
Select Left Trunk Support Size		
<input type="checkbox"/> LATXL	Size XS (3.25" x 4.25")	\$100
<input type="checkbox"/> LATSL	Size S (3.5" x 5.25")	\$100
<input type="checkbox"/> LATML	Size M (3.75" x 5.75")	\$100
<input type="checkbox"/> LATLL	Size L (4.25" x 6.5")	\$100
<input type="checkbox"/> SALATR	Right MAXX Style Swing Away HDW ⁽¹⁾	\$195
Select Right Trunk Support Size		
<input type="checkbox"/> LATXR	Size XS (3.25" x 4.25")	\$100
<input type="checkbox"/> LATSR	Size S (3.5" x 5.25")	\$100
<input type="checkbox"/> LATMR	Size M (3.75" x 5.75")	\$100
<input type="checkbox"/> LATLR	Size L (4.25" x 6.5")	\$100

ATLL Additional Telescoping Link - Left⁽²⁾

ATLR Additional Telescoping Link - Right⁽²⁾

1. These laterals are not available with Matrx PB, PB Elite Deep and PB Deep Backs
2. Up to 2.75" medial offset

FIXED HIP/KNEE SUPPORTS

AADL SEATING BENEFIT

Lateral Hip Support with Fixed Mounting Hardware Left⁽¹⁾

Select Pad Size and mounting position at Hip:

<input type="checkbox"/> HSPXLH	Size XS (5.5" W x 4" H) mounted at Hip	\$120
<input type="checkbox"/> HSPSMLH	Size S/M (7" W x 4" H) mounted at Hip	\$120
<input type="checkbox"/> HSPSLH	Size S (7"W x 6"H) mounted at Hip	\$120
<input type="checkbox"/> HSPMLH	Size M (9" W x 4"H) mounted at Hip	\$120
<input type="checkbox"/> HSPLLH	Size L (13" W x 4"H) mounted at Hip	\$120

Select Pad Size and mounting position at Knee:

<input type="checkbox"/> HSPXLK	Size XS (5.5"W x 4" H) mounted at Knee	\$120
<input type="checkbox"/> HSPSMLK	Size S/M (7"W x 4" H) mounted at Knee	\$120
<input type="checkbox"/> HSPSLK	Size S (7"W x 6" H) mounted at Knee	\$120
<input type="checkbox"/> HSPMLK	Size M (9"W x 4"H) mounted at Knee	\$120
<input type="checkbox"/> HSPLLK	Size L (13"W x 4"H) mounted at Knee	\$120

FIXED HIP/KNEE SUPPORTS-CONTINUED

AADL SEATING BENEFIT

Lateral Hip Support with Fixed Mounting Hardware Right⁽¹⁾

Select Pad Size and mounting position:

<input type="checkbox"/> HSPXRH	Size XS (5.5"W x 4"H) mounted at Hip	\$120
<input type="checkbox"/> HSPSMRH	Size S/M (7"W x 4"H) mounted at Hip	\$120
<input type="checkbox"/> HSPSRH	Size S (7"W x 6"H) mounted at Hip	\$120
<input type="checkbox"/> HSPMRH	Size M (9"W x 4"H) mounted at Hip	\$120
<input type="checkbox"/> HSPLRH	Size L (13"W x 4"H) mounted at Hip	\$120

Select Pad Size and mounting position at Knee:

<input type="checkbox"/> HSPXRK	Size XS (5.5"W x 4"H) mounted at Knee	\$120
<input type="checkbox"/> HSPSMRK	Size S/M (7"W x 4"H) mounted at Knee	\$120
<input type="checkbox"/> HSPSRK	Size S (7"W x 6"H) mounted at Knee	\$120
<input type="checkbox"/> HSPMRK	Size M (9"W x 4"H) mounted at Knee	\$120
<input type="checkbox"/> HSPLRK	Size L (13"W x 4"H) mounted at Knee	\$120

1. Must choose swing away hardware below if required.

REMOVABLE HIP/KNEE SUPPORT HARDWARE

AADL SEATING BENEFIT

Lift Off Removable Hardware for Hip Support

- LOHL Left \$195
- LOHR Right \$195

MAXX Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports

Select mounting position:

- MHSHLH Left mounted at Hip \$215
- MHSHLK Left mounted at Knee \$215
- MHSHRH Right mounted at Hip \$215
- MHSHRK Right mounted at Knee \$215

Swing Away Removable Hardware for Hip Supports

Select mounting position:

- SHSHLH Left mounted at Hip \$215
- SHSHLK Left mounted at Knee \$215
- SHSHRH Right mounted at Hip \$215
- SHSHRK Right mounted at Knee \$215

Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports⁽¹⁾

Select mounting position:

- EMQHLH Left mounted at Hip \$215
- EMQHLK Left mounted at Knee \$215
- EMQHRH Right mounted at Hip \$215
- EMQHRK Right mounted at Knee \$215

1. Recommended when finished seat cushion thickness will be greater than 3.5" high.

ELBOW BLOCKS⁽¹⁾

AADL SEATING BENEFIT

Left Elbow Block

Select Pad Size and Both Extensions:

- EBXL Size XS (3.25" x 4.25") \$195
- EBSL Size S (3.5" x 5.25") \$195
- EBML Size M (3.75" x 5.75") \$195
- EBLL Size L (4.25" x 6.5") \$195

- V3 Left Verticle Extension 3"
- V5 Left Verticle Extension 5"
- H5 Left Horizontal Extension 5"
- H7 Left Horizontal Extension 7"

Right Elbow Block

Select Pad Size and Both Extensions:

- EBXR Size XS (3.25" x 4.25") \$195
- EBSR Size S (3.5" x 5.25") \$195
- EBMR Size M (3.75" x 5.75") \$195
- EBLR Size L (4.25" x 6.5") \$195

- V3 Right Verticle Extension 3"
- V5 Right Verticle Extension 5"
- H5 Right Horizontal Extension 5"
- H7 Right Horizontal Extension 7"

1. Prevents arm from slipping.

ADDITIONAL ARM ACCESSORIES

- UCH Universal Cup Holder⁽¹⁾ CLIENT COST \$50
- MACESL Multi Axis Upper Extremity Support⁽²⁾(each)
Left Side \$225
- MACESR Multi Axis Upper Extremity Support⁽²⁾ (each)
Right Side \$225
- OAS2L 1" Offset Armrest Spacer Left(each) \$125
- OAS2R 1" Offset Armrest Spacer Right(each) \$125
- OAS1L 1/2" Offset Armrest Spacer Left(each) \$125
- OAS1R 1/2" Offset Armrest Spacer Right(each) \$125
- MPAPA Multi Purpose Arm Pad
Adapter Block⁽³⁾ (each) Quantity () \$35

1. Not available on the joystick side.
2. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
3. Not available with Standard Arm Pads or Waterfall Arm Pads.

ARM ADDUCTOR PAD⁽¹⁾

- AADSL Left Small (6" x 2.25") \$135
- AADSR Right Small (6" x 2.25") \$135
- AADLL Left Large (8.5" x 2.25") \$135
- AADLR Right Large (8.5" x 2.25") \$135

1. Only available with Modular or Flat Multi Position Arm Pads.

VENT TRAYS⁽¹⁾⁽²⁾⁽³⁾

- ARTVTH Vent Tray - Articulating Horizontal⁽⁴⁾ \$1,350
(14"W x 15"D x 10"H)
- ARTVTV Vent Tray - Articulating Vertical⁽⁴⁾ \$1,350
(14.625"W x 10.125"D x 13"H)
- TRILVT Vent Tray - Trilogy 100/200 \$1,350
(14"W x 6"D)

1. The user weight limit is 225 lbs. Not available with TRRO or TRRO-E brackets.
2. The maximum weight capacity of the vent tray is 35 lbs.
3. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. Required for use on Recline Systems

COMPACT-LAPTOP VENT TRAYS⁽¹⁾⁽²⁾⁽³⁾

- CVT Vent Tray \$1,350
- VBHA Vent Bag Holder Articulating⁽⁴⁾⁽⁵⁾ \$990
- VBHNA Vent Bag Holder Non-Articulating \$990

1. The user weight limit is 225 lbs. Not available with TRRO or TRRO-E brackets.
2. The maximum weight capacity of the vent tray is 20 lbs.
3. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. Required for use on Recline Systems.
5. Motion Concepts does not supply the Vent Bag

