



PRINT FORM

SAVE

CLEAR FORM

Contact Motion Concepts Customer Service: 1.888.433.6818 | quotes@motionconcepts.com

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

These Accessory Additions are only available with the purchase of a new power positioning system and completed order form with all Client Specifications completed. Subject to approval.

ACT4 electronics are not available for use with any power accessory.

WHEELCHAIR ACCESSORIES

BRAND NEW!

SELF CARE

☐ HK Hydration Kit⁽¹⁾ \$250

Vacuum insulated 24 oz stainless steel bottle keeps drinks cold for up to 18 hours. Features a fully adjustable 36" neck with drinking tube. The check valve prevents liquid from flowing back into the bottle. For those with limited lung function, keeping liquid primed and at-the-ready makes staying hydrated a breeze.

☐ CRVMR Rear View Mirror Right Mounted \$45

☐ CRVML Rear View Mirror Left Mounted \$45

☐ PH X-Grip Phone Holder \$75
Holds up to 3.25" wide Phone

The specialized mounting brackets allow the Phone to be seamlessly integrated via BALL MOUNT to the Arm Tube of your Motion Concepts Power Positioning System. The mounts can be reconfigured to accommodate a variety of mounting position

☐ MUB Medical Utility Bag w/ Mounting Hardware . \$250

1. Not available with TRRO, TRRO-E, XRO, XRO-E brackets

WHEELCHAIR COVER⁽¹⁾

☐ TAG571850 Wheelchair Cover \$105


Weather-Resistant black Nylon cover. An adjustable strap at the bottom of the cover provides a way to adjust the size.

1. Not available with TRRO, TRRO-E, XRO, XRO-E brackets

REAR STORAGE BASKET⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾

☐ RSB Removable Rear Storage Basket. . . . \$1,295

Exterior dimensions: 18"W x 6"H x 7"D
Interior dimensions: 17.8"W x 6"H x 6.8"D
Interior capacity: 3 gallon
Max Load: 20 lbs



1. Tall items extending beyond the basket's height may be damaged upon tilting.
2. Only available on TDX-SP2, ROVI X3/A3 bases.
3. Not available on Stander (MPSS) systems
4. Not available with TRRO, TRRO-E, XRO, XRO-E brackets

MISCELLANEOUS ACCESSORIES

☐ CSCC Cane and Crutch Holder. \$395

- ☐ Left Mount
☐ Right Mount

☐ CFH Flag Holder⁽¹⁾. \$65

- ☐ Left Mount
☐ Right Mount

☐ AMC Accessory Mount with Clip (QTY ____)⁽¹⁾. . . \$50

1. Not available with TRRO, TRRO-E, XRO, XRO-E brackets

PLEASE SPECIFY PAGES INCLUDED WITH ORDER

☐ Page 1 Excludes Power Accessories on pages 2-5

☐ Pages 1-5 Complete Order Form

Total Retail Value \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM

WHEELCHAIR POWER ACCESSORIES**LATERAL TILT MODULE⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾⁽⁷⁾⁽⁸⁾**☐ **LATM Power Lateral Tilt Module(up to 15°) . \$3,995**

All options must be selected on the system order form prior to selecting this option.

☐ **LATDY Reduced Drive while in Lateral Tilt - Yes**☐ **LATDN Reduced Drive while in Lateral Tilt - No**

1. Weight capacity is 225 lbs.
2. Only available on TDX SP2, Aviva FX, Aviva RX, 6000ZHD, Edge HD (must use 22 batteries), and Rovi A3 Bases.
3. When ordering a lateral tilt module on a Storm short base frame, standard fixed casters must be ordered on the base.
4. Only available on 16" to 20" seat widths and 16" to 22" seat depths.
5. Only available with Multi-function electronics Part # MFCBD.
6. Must order the Group 24 battery tray for the TDX-SP2, which will utilize Group 22 batteries to accommodate the lateral tilt module.
7. Vents, articulating vents and articulating O2 holders are not available with Lateral Tilt.
8. Minimum Seat to floor height for bases: TDX SP2 - 20.25", Aviva FX/RX - 18.75", 6000ZHD/Edge HD - 19.75", ROVI A3 - 19.75"

POWER ANTERIOR TILT MODULE⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾⁽⁷⁾⁽⁸⁾⁽⁹⁾☐ **PAT45 Power Anterior Tilt Module(up to 45°) \$3,995**

All options must be selected on the system order form prior to selecting this option.

☐ **SOCM Synchronized Option for Center Mount Front Rigging with Anterior Tilt⁽¹⁰⁾ \$1,350**
MANDATORY For LNX, Maxx Style Fixed and Power Elevating Calf Panel.

☐ **SOR Synchronized Option for Recline with Anterior Tilt N/C**
 For use on Tilt/Recline Systems, recommended with PAT45.

1. Strongly recommend selecting Recline with Power Anterior Tilt.
2. Weight capacity 250 lbs.
3. Only available on 16" to 20" seat widths and 16" to 22" seat depths.
4. Only available with Multi-Function electronics Part # MFCBD.
5. All front rigging will be set a 70° unless using the Synchronized Center Mount option. This module may limit front rigging options.
6. Minimum 18" seat to floor height at zero degrees.
7. Minimum finished knee to heel to touch floor is 13".
8. Not available with 12" Power Adjustable Seat Height Module.
9. Power Telescoping and Power Flip up foot platforms are not available.
10. This is an upgradeable option only, must be selected all 90° Front Rigging.

FIXED ANTERIOR ASSIST⁽¹⁾⁽²⁾

☐ **FA10DLO 10° Fixed Anterior Assist/Drive Lock Out⁽³⁾⁽⁴⁾ \$750**
 All options must be selected on the system order form prior to selecting this option.

☐ **FA5DLO 5° Fixed Anterior Assist/Drive Lock Out \$295**

1. This module may limit front rigging options. The 10° Fixed Anterior Assist adds approximately 3.6" to the seat to floor height at zero degrees(level). The 5° Fixed Anterior Assist adds approximately 2.2" to the seat to floor height at zero degrees(level).
2. Drive lock out is necessary when there is less than 3" ground clearance.
3. Only available with Multi-Function electronics Part # MFCBD.
4. Depending on seat depth and knee to heel, drive in pre-tilt may be required and/or restricted anterior may be necessary.

UPFRONT MODULE⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾⁽⁷⁾⁽⁸⁾⁽⁹⁾ **ROVI BASES ONLY**

☐ **UM Upfront Module (up to 45°) with Anterior Tilt and 7" of Seat Elevation \$5,995**

System moves Forward 5" while elevating forward and includes a synchronized leg mechanism to tuck when elevated forward and when extends forward in anterior tilt. All options must be selected on the system order form prior to selecting this option.

☐ **RDY Reduced Drive for Upfront Module - Included**

1. Weight capacity is 250 lbs.
2. Minimum Seat to Floor Height (STF) is 18-18.5".
3. Only available with LNX Power Center Mount Front Rigging or Maxx Style Manual Center Mount foot platform.
4. Recline is recommended to assist with forward transfer.
5. Only available on 16" to 20" seat widths and 16" to 22" seat depths.
6. Only available with Multi-Function electronics Part # MFCBD.
7. Articulating Vents and O2 Holders are not available.
8. Power Sliding Back (PSB) is not available with the upfront module.
9. Power Telescoping and Power Flip up foot platforms are not available.

POWER SWING AWAY JOYSTICKS

☐ **PSAJH Power Swing Away Joystick Horizontal⁽¹⁾⁽²⁾ (Quad Link Style) \$2,995**

☐ **PSAJV Power Swing Away Joystick Vertical Flip Up⁽¹⁾⁽²⁾ (Mid Line Style) \$2,995**

Select Mounting Position:

☐ **MLPS Mount on Left Side**
☐ **MRPS Mount on Right Side**

☒ **SESBS Stealth Egg Switch - Black**

☐ **JPOPS Operate Power Joystick/Chin Boom with Joystick Power "off" option⁽³⁾ \$750**

1. Comes with an egg switch.
2. Counts as an additional electronic function.
3. Allows Power Swing Away joystick to operate when chair is turned off.

POWER SWING AWAY CHIN BOOM

☐ **PSACB Power Swing Away Chin Boom⁽¹⁾⁽²⁾ Joystick Power "on" operation \$2,995**

Select Mounting Position:

☐ **MLCB Mount on Left Side**
☐ **MRCB Mount on Right Side**

☐ **JPOCB Operate Power Joystick/Chin Boom with Joystick Power "off" option⁽³⁾ \$750**

1. Comes with an egg switch
2. Counts as an additional electronic function.
3. Allows Power Swing Away Chin Boom to operate when chair is turned off.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

POWER RETRACTABLE HEADREST⁽¹⁾⁽²⁾⁽³⁾

- ☐ **PRH Power Retractable Headrest.....\$3,995**
The Headrest pad moves up 3.5" and forward 4".

1. Finished Back Height must be 22" high or greater.
2. Counts as an additional electronic function.
3. Not available with Power Sliding Back (PSB).
4. Not compatible with aftermarket backs, for use with Motion or Matrix backrests.
5. Seat to top of head measurement required

POWER CHEST RESTRAINT⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾

- ☐ **PCB Power Chest Restraint.....\$3,995**

1. Weight Capacity is 200 lbs
2. Counts as an additional electronic function.
3. Compatible with Dual Post Armrest and Reclining Adjustable Height Armrest.
4. Foam wrapped bar.
5. Not available with XRO/XROE, or TRRO/TRRO-E

TRANSFER ARM⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾⁽⁷⁾⁽⁸⁾

- ☐ **PTAL Power Transfer Armrest - Left Side.....\$1,995**
☐ **PTAR Power Transfer Armrest - Right Side.....\$1,995**
☐ **MTAL Manual Transfer Armrest - Left Side.....\$995**
☐ **MTAR Manual Transfer Armrest - Right Side.....\$995**
☐ **PFBAL Power Flip Back Armrest - Left Side.....\$1,995**
☐ **PFBAR Power Flip Back Armrest - Right Side.....\$1,995**

1. Weight capacity is 175 lbs.
2. Only with Adjust Dual Post Sytle Armrest for Tilt/Recline, arm height 10.5"-14.5"
3. Requires Multi-function electronics ordered on chair Part # MFCBD.
4. Minimum seat depth available is 16" and 16"-20" seat width.
5. Minimum seat to floor height 18"
6. Transfer only available in the chair home position.
7. Not available with Upfront Module.
8. Not available with XRO/XROE, or TRRO/TRRO-E

POWER COMMUNICATION DEVICE⁽¹⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾

- ☐ **PCDSL Power Communication Device Single Actuator - Left Joystick Power "on" operation.....\$2,995**
☐ **PCDSR Power Communication Device Single Actuator - Right Joystick Power "on" operation.....\$2,995**
☐ **PCDDL Power Communication Device Dual Actuator - Left Joystick Power "on" operation.....\$3,995**
☐ **PCDDR Power Communication Device Dual Actuator - Right Joystick Power "on" operation.....\$3,995**
☐ **JPOCD Operate Power Joystick/Chin Boom Communication Device Joystick Power "off" option⁽⁷⁾\$750**
☐ **AE Additional Electronics\$3,995**

1. Considered an additional electronic function and requires Multi-function electronics ordered on chair Part # MFCBD.
2. Consider as an additional two functions.
3. May require additional electronics if system has more than 6 actuators (functions).
4. Minimum seat depth available is 16".
5. Weight capacity 20 lbs.
6. Comes with an egg switch.
7. Allows Power Communication Device to operate when chair is turned off.

POWER BACK PACK⁽¹⁾

- ☐ **PBPHL Power Back Pack Holder with Single Back Pack mounted on Left Side\$2,995**
☐ **PBPHR Power Back Pack Holder with Single Back Pack mounted on Right Side\$2,995**

1. Weight limit 20lbs.

POWER FOOT PLATFORM OPTIONS FOR LNX

For use with LNX Foot platform ordered separately

- ☐ **LNXPTF Power Telescoping Foot Platform for LNX⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾⁽⁷⁾\$2,250**
Powers lower leg length vertically up to a maximum of 4" extension.
☐ **LNXPFUF Power Flip up Foot Platform for LNX⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾⁽⁷⁾\$1,995**

1. Available with one piece calf panel only and not available with a finished knee-to-heel of 13" or less. Not available with Split Footplates.
2. Comes standard with a 12" x 10.75" footplate made of durable aluminum 1/4" thick with heel loops.
3. The front rigging will be set at approximately 80° on the Storm 3G.
4. Minimum finished knee to heel is 13".
5. Minimum seat to floor height is 17".
6. Not available on Upfront Module.
7. Requires Multi-function electronics ordered on chair Part # MFCBD.

POWER FOOT PLATFORM OPTIONS⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾⁽⁷⁾

Fixed Center Mount footplatform is included

- ☐ **FCMFPTF Power Telescoping Foot Platform with Fixed Center Mount.....\$2,845**
Powers lower leg length vertically up to a maximum of 4" extension.
☐ **FCMFPUF Power Flip up Foot Platform with Fixed Center Mount.....\$2,250**

1. Includes the Fixed Center Mount Foot platform which does not elevate.
2. Comes standard with a 12" x 10.75" aluminum footplate with heel loops. Optional sizes available per standard order form.
3. Not available on Storm 3G unless Flared Front Casters are added to base.
4. Minimum finished knee to heel is 13".
5. Minimum seat to floor height is 17".
6. Not available on Upfront Module.
7. Requires Multi-function electronics ordered on chair Part # MFCBD.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____


Tag: _____

VERTICAL LIFTING PLATFORM⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾⁽⁷⁾

- ☐VLP Vertical Lifting Platform 11" of total travel \$4,995
- ☐OCUFP Onyx Pad for Foot Platform \$450

Specify preferred STOP and START positions

Area	Measurement in Inches
A. Stop	
Travel	11"
B. Start	



1. Standard Tilt-to-Drive, or Elevate-to-Drive (elevate only systems) lockout included
2. Comes standard with a 12" x 10.75" aluminum footplate. Optional sizes available per standard order form
3. Requires Multi-function electronics ordered on chair Part # MFCBD
4. 175 lb weight capacity
5. Minimum seat to floor height is 17.75"
6. Not available with Small Adult seat sizes
7. Not available on Rehab systems or Upfront Module

POWER CALF PANEL ONLY⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾

- ☐PECP Power Elevating Calf Panel Only (no articulation) \$2,995

1. No footplate available.
2. Calf panel measures up to 15" in length and measures 13" wide at the bottom and tapers up to within 2" narrower than the seat width.
3. Minimum seat depth 17" deep.
4. Not available for use with Dual Front Rigging adaptor.

POWER ELEVATING & ARTICULATING CALF PANEL⁽¹⁾⁽²⁾⁽³⁾⁽⁶⁾⁽⁷⁾

- ☐PEACP Power Elevating & Articulating Calf Panel⁽⁴⁾.....\$4,345
- ☐PEACPDF Power Elevating & Articulating Calf Panel with Power Down Footplate^{(12"W x 10.75")⁽⁴⁾⁽⁵⁾..\$4,995}

Options

- ☐MCP Meshtex Calf Panel. STD
- ☐VCP Vinyl Calf Panel. STD

1. Calf panel measures up to 15" in length and measures 13" wide at the bottom and tapers up to within 2" narrower than the seat width.
2. Minimum seat depth 17" deep.
3. Knee to heel (KTH) = 12"-16"
4. 3.5" of articulation
5. Power down footplate travels 4" maximum.
6. Not available for use with Dual Front Rigging adaptor.
7. Only available with Multi-Function electronics

CUSTOMIZED FOOT PLATFORM OPTIONS

- ☐CFP Custom Foot Platform⁽¹⁾\$550

	Area	Measurement in Inches
	1.	
	2.	
	3.	

- ☐OCUFP Onyx Pad for Any Foot Platform⁽²⁾⁽³⁾\$450

1. Available for Maxx/LNX Center Mount.
2. Available for all Maxx/LNX Center Mount Foot Platform options.
3. Pad thickness is 1"

CUSTOMIZED LEGREST OPTIONS⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾

- ☐PALSAMS Power Assisted Legrest Swing Away Mechanism controlled through Switch.....\$3,995
- ☐PALSAMJ Power Assisted Legrest Swing Away Mechanism controlled through Joystick\$3,995

1. Considered an additional electronic function and requires Multi-function electronics ordered on chair Part # MFCBD.
2. Compatible with Motion HD 70° Swing Away footrest, Maxx Style Swing Away Legrest Power or Manual.
3. Legrests are not removeable.
4. Minimum 16" seat width and minimum 16" seat depth.

MANUAL ELEVATING FOOT CAGE

- ☐MMEFC MAXX Manual Elevating Foot Cage⁽¹⁾ . \$1,995 25° of Elevation
- ☐FCFP Foot Plate 8.5" w x 7.5" d STD
- ☐FCHR Heel Rest 7.5" w N/C

1. Starting position will vary depending on Base, and Seat Depth of System.

OPTIONAL CALF PANEL⁽¹⁾⁽²⁾

- ☐OSCPL Optional Sized Flat One Piece Calf Panel for LNX\$450

Maximum pad size is seat width less 2" and 10" in length. With minimum 13" finished knee to heel.

Specify preferred size each Height: _____ Width: _____

1. Subject to clearance available.
2. Startex fabric over polyurethane foam is standard.

KNEE BOLSTER FOR LNX & MAXX STYLE FIXED CENTER MOUNT⁽¹⁾⁽²⁾

- ☐KBS Small - 8" wide Pads\$1,750
- ☐KBL Large - 10" wide Pads\$1,750

1. Not available on Fixed Center Mount.

Dealer Name: _____
Dealer Account Number: _____
PO #: _____
Tag: _____

SPECIALTY LEG PAD OPTIONS

Lateral Leg Support for LNX One Piece Calf Panel (pair)

- ☐ LLSLNXX Extra Small Lateral Pad (5.5"W x 4"H)\$650
- ☐ LLSLNXS Small Lateral Pad (7"W x 6"H)\$650
- ☐ LLSLNXM Medium Lateral Pad (9"W x 4"H)\$650
- ☐ LLSLNLX Large Lateral Pad (13"W x 4"H)\$650

Lateral Leg Supports for Power Elevating Calf Pane (pair)

- ☐ LLSPCPX Extra Small Lateral Pad (5.5"W x 4"H)\$650
- ☐ LLSPCPS Small Lateral Pad (7"W x 6"H)\$650
- ☐ LLSPCPM Medium Lateral Pad (9"W x 4"H)\$650
- ☐ LLSPCPL Large Lateral Pad (13"W x 4"H)\$650
- ☐ OSICP Optional Sized Individual Calf Pads for Motion Concepts Front Rigging\$650

Specify preferred size each Height: _____ Width: _____

SPECIALTY ARM PAD OPTIONS

- ☐ OSFAP Optional Sized Flat Arm Pads\$650

Specify preferred size each Length: _____ Width: _____

- ☐ OSLPEB Optional Sized Lateral Pad or Elbow Block ...\$650

Specify preferred size each Length: _____ Width: _____

- 1. Subject to clearance available.
- 2. Startex fabric over polyurethane foam is standard.

SPECIALTY FOOT PAD OPTIONS

- ☐ OPIF Onyx Pad for Individual Foot Plates (Pair)⁽¹⁾⁽²⁾⁽³⁾\$650

- 1. Pad thickness is 1"
- 2. Not available on Maxx Independent Adjustable Positioning Foot plates. Not available on Stander systems or with Power Anterior or Upfront modules.
- 3. Not recommended with Heel Loops or Toe Straps

ADDITIONAL BATTERY TRAYS

- ☐ SLBT40 Storm Long Base Additional Front Mounted Battery Tray for Group 40 Battery\$500
- ☐ SLBTU1 Storm Short Base Additional Front Mounted Battery Tray for U1 Battery.....\$500
- ☐ S322 Storm Base 3 x Group 22 Battery Tray Conversion.....\$950
- ☐ TDXU1 TDX-SP Base Rear Caster Mounted U1 Battery Tray⁽¹⁾\$500
- ☐ EDGE322 Edge HD 3 x Group 22 Battery Tray Conversion⁽²⁾\$500

- 1. Available on Wide Base TDX-SP Only and must have factory Transport brackets - TBKTS.
- 2. Deletes factory battery boxes.

ADDITIONAL ELECTRONIC OPTIONS⁽¹⁾

- ☐ SDC Step Down Converter\$500
Converts 24v to 12v
- ☐ RPVPC Richardson Powertech Vent Power Center . \$1,950
 - ☐ 90111 For use with TRILOGY vent
 - ☐ 90110 For use with OTHER vent

1. Should this product be used in life support applications Motion Concepts will not assume any risk and shall be held harmless against all damages.

MEMORY SEATING OPTIONS⁽¹⁾

- Single Function Memory Seating⁽²⁾ \$495
 - ☐ SFMSS Through Switch
 - ☐ SFMSJ Through LiNX Joystick⁽³⁾⁽⁴⁾
 - Optional:
 - ☐ SFMSAS 4-Way Splitter for Additional Phono Switches \$125

- Multiple Function Memory Seating \$695

- ☐ MFMSS Through Switch
- ☐ MFMSJ Through LiNX Joystick⁽³⁾⁽⁴⁾
- Optional:
 - ☐ MFMSAS 4-Way Splitter for Additional Phono Switches \$125

- 1. Four customizable and programmable positions. Comes with a Quad Push Button switch with LED indicator. Only available with Multi-Function Through Joy (MFCBD) electronics. Not available on systems with Lateral Tilt, Anterior Tilt, or Upfront Modules. Will not operate Power Adjustable Seat Height modules (elevating seat function) or Power Sliding Back (PSB).
- 2. For Single Actuator systems.
- 3. Only available on TDX-SP2, TDX-SP2 HD, AVIVA RX, and AVIVA FX bases.
- 4. Requires Expandable Joystick and Controller.

SAFETY LIGHTS

- ☐ FSLS Front Safety Lights - through Switch⁽¹⁾\$495
- ☐ FSLJ Front Safety Lights - through Joystick⁽²⁾\$995
- ☐ RSLS Rear Safety Lights - through Switch⁽¹⁾\$495
- ☐ RSLJ Rear Safety Lights - through Joystick⁽²⁾\$995
- ☐ LPC Front & Rear Light Combo Package Credit* (\$150)

*Purchase both Front and Rear Lights and receive a credit.

- ☐ GEL Ground Effects Lighting through Switch⁽¹⁾\$350

- 1. Operated through a separate power supply and single push button switch.
- 2. Only available with Multi-Function electronics
- 3. Not available with Power Sliding Back (PSB).

Total Retail Value \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM