



## Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



### Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

**SAVE**

### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

**SUBMIT**

**CLEAR FORM**

### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**888-433-6818**

**PRINT FORM**



# ULTRA LOW MAXX CG POWER POSITIONING SYSTEM ON TDX-SP2 BASE

US PRICE LIST AND ORDER FORM - Price Effective May 1, 2025



PRINT FORM

SAVE

CLEAR FORM

Contact Motion Concepts Customer Service: 1.888.433.6818 | quotes@motionconcepts.com

## REQUIRED INFORMATION

☐ Quote ☐ Order

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_ Tag: \_\_\_\_\_

Purchasing Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_

ATP/Therapist: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_


Zip: \_\_\_\_\_ Client Gender: ☐ M ☐ F

Special Client Conditions: \_\_\_\_\_

## REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all **REQUIRED\*** measurements

|                               |  |
|-------------------------------|--|
| <b>Client Height*</b>         |  |
| <b>Client Weight (lbs)*</b>   |  |
| <b>A. Seat to Shoulder</b>    |  |
| <b>B. Trunk Depth</b>         |  |
| <b>C. Chest Width</b>         |  |
| <b>D. Knee to Back</b>        |  |
| <b>E. Seat to Top of Head</b> |  |
| <b>F. Back to Joystick</b>    |  |
| <b>G. Seat to Elbow</b>       |  |
| <b>H. Hip Width</b>           |  |
| <b>I. Knee to Heel*</b>       |  |
| <b>Cushion Thickness*</b>     |  |



## REHAB SEAT BASE OPTION

- ☐ TDXSP2 TDX-SP2 Base for Rehab Seat System<sup>(1)</sup>  
HCPCS code K0848 ..... \$10,995

1. Must select Ultra Low Maxx Rehab Seat (ULMRSN) system.

## POWER BASE OPTIONS

- ☐ TDXSP2-CG TDX-SP2 Base for Single Actuator System  
HCPCS code K0856. .... \$7,595

- ☐ TDXSP2-MCG TDX-SP2 Base for Multiple Actuator System  
HCPCS code K0861 ..... \$7,595

## BASE WIDTH

- ☐ NARROWFRM Narrow Base 24" (MK 22 NF Batteries) ..... N/C  
☐ STDFRM Standard Base 25.5" (MK 24 NF Batteries) ..... N/C

## SEAT-TO-FLOOR HEIGHT<sup>(1)</sup>

- ☐ SLOW 16.75" Seat-To-Floor<sup>(2)</sup> ..... \$250  
☐ LOW 17.25" Seat-To-Floor ..... N/C  
☐ MED 18.25" Seat-To-Floor ..... N/C  
☐ TALL 19.25" Seat-To-Floor ..... N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"  
2. This seat to floor height is only available with the Group 24 Battery Tray (24TRY) and on 17" seat depths and greater.

## USER WEIGHT LIMIT

- ☒ U300 Weight Capacity up to 300lbs ..... STD

## WHEELCHAIR TRANSPORT BRACKETS

- ☐ TRBKTS Wheelchair Transport Brackets<sup>(1)</sup> ..... STD  
☐ TRRO Transport Ready Option<sup>(2)(4)</sup> ..... \$275  
☐ TRRO-E Transport Ready Option for Elevate<sup>(2)(3)(4)</sup> ..... \$275

1. For unoccupied use only

2. TRRO meets ISO-7176 Part 19

3. Required for Elevate Systems to meet ISO-7176 Part 19, includes two additional tie-down loops on rail.

4. Must select a Motion Concepts or Matrix Elan headrest with hardware and Motion Concepts Standard Rehab, High Back or Matrix Upgradeable back. Comes with pelvic safety belt standard, please select TROBELT on page 5.

## TIRE OPTIONS

- ☐ B1431-3 14" x 3" Black Tire with Gel Foam Inserts ..... STD  
☐ B1430-3 14" x 3" Black Tire - Pneumatic Tire ..... N/C

## FORK OPTIONS

- ☐ DSFK Double-Sided Fork Package ..... STD

## BATTERY CHARGER OPTION

- ☐ I10CHARGER 110 Volt Battery Charger ..... N/C  
☐ L8900 Omit Battery Charger ..... (\$80)

## HIGHSPEED MOTORS AND ENHANCED DRIVING

- ☐ HSMTR Performance Motors 7.5 mph<sup>(1)</sup> ..... \$1,400  
☐ GTRAC-LX LiNX G-TRAC Module<sup>(2)</sup>  
HCPCS code K0108 ..... \$995

1. Must order Standard Base 25.5"

2. Recommended for use with Specialty Drive Controls

## LIGHTS AND INDICATORS

- ☐ LIGHTS LED Lights and Indicators<sup>(1)</sup> ..... \$700

1. Must order REM215, REM216, REM400 or REM500

# Ultra Low Maxx CG Power Positioning System

## TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### FRAME FINISH

|                               |                                                                                   |                     |     |
|-------------------------------|-----------------------------------------------------------------------------------|---------------------|-----|
| <input type="checkbox"/> 170P |  | Obsidian Black..... | N/C |
| <input type="checkbox"/> 171P |  | Ruby Red.....       | N/C |
| <input type="checkbox"/> 172P |  | Peacock Blue.....   | N/C |
| <input type="checkbox"/> 174P |  | Purple Passion..... | N/C |
| <input type="checkbox"/> 175P |  | Denim Navy.....     | N/C |
| <input type="checkbox"/> 176P |  | Lime Twist.....     | N/C |
| <input type="checkbox"/> 177P |  | Arctic Blue.....    | N/C |

### RIM INSERT COLORS - MUST PICK ONE

|                                |                                                                                    |                       |     |
|--------------------------------|------------------------------------------------------------------------------------|-----------------------|-----|
| <input type="checkbox"/> 170PR |   | Obsidian Black.....   | N/C |
| <input type="checkbox"/> 171PR |   | Ruby Red.....         | N/C |
| <input type="checkbox"/> 172PR |   | Peacock Blue.....     | N/C |
| <input type="checkbox"/> 173PR |   | Pink Flamingo.....    | N/C |
| <input type="checkbox"/> 174PR |   | Purple Passion.....   | N/C |
| <input type="checkbox"/> 175PR |   | Denim Navy.....       | N/C |
| <input type="checkbox"/> 176PR |   | Lime Twist.....       | N/C |
| <input type="checkbox"/> 177PR |   | Arctic Blue.....      | N/C |
| <input type="checkbox"/> 178PR |   | Banana Split.....     | N/C |
| <input type="checkbox"/> 179PR |   | Orange Marmalade..... | N/C |
| <input type="checkbox"/> 180PR |  | Morning Mist.....     | N/C |

### NON-EXPANDABLE SINGLE/DUAL JOYSTICKS

|                                 |                                                          |       |                                                                         |
|---------------------------------|----------------------------------------------------------|-------|-------------------------------------------------------------------------|
| <input type="checkbox"/> REM110 | LED Drive Only Remote/Joystick <sup>(1)</sup> .....      | N/C   | Single drive profile                                                    |
| <input type="checkbox"/> REM210 | LED Non-Expandable Remote/Joystick <sup>(1)</sup> .....  | N/C   | Three drive profiles, control of two seat functions                     |
| <input type="checkbox"/> REM215 | LED Non-Expandable Remote with Lights <sup>(2)</sup> ... | \$350 | Three drive profiles, control of two seat functions and light operation |

1. Not available with lights and indicators.
2. Must choose lighting package separately

### EXPANDABLE JOYSTICKS<sup>(1)</sup>

|                                   |                                                                    |         |                                                                                                                |
|-----------------------------------|--------------------------------------------------------------------|---------|----------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> REM211   | LED Expandable Remote/Joystick <sup>(2)</sup> .....                | N/C     | Three drive profiles, control of up to 6 actuators<br>Not for use with specialty controls                      |
| <input type="checkbox"/> REM216   | LED Expandable Remote/Joystick<br>with Lights <sup>(3)</sup> ..... | \$350   | Three drive profiles, control of up to 6 actuators<br>and light operation. Not for use with specialty controls |
| <input type="checkbox"/> REM400   | Color 3.5" Touch Screen<br>Remote/Joystick <sup>(4)</sup> .....    | \$1,000 | Four drive profiles, control of up to 6 actuators through the<br>drive control and light operation.            |
| <input type="checkbox"/> SW400TGL | Toggle Kit for REM400 <sup>(5)</sup> .....                         | \$200   |                                                                                                                |

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Not available with lights and indicators
3. Must choose lighting package separately.
4. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.
5. Adds 2 toggle switches to REM400.

### EXPANDABLE CONTROLLER AND HARNESS

|                               |                                                               |       |
|-------------------------------|---------------------------------------------------------------|-------|
| <input type="checkbox"/> EXPC | Expandable Controller <sup>(1)</sup><br>HCPCS code E2377..... | \$700 |
| <input type="checkbox"/> PWH  | Harness for Expandable Controller<br>HCPCS code E2313.....    | \$450 |

1. This item is required when ordering power tilt, power recline, and power elevating legrest together OR when ordering an Alternative Drive Control.

### EXPANDABLE SPECIALTY DISPLAY

|                                 |                                                                    |         |
|---------------------------------|--------------------------------------------------------------------|---------|
| <input type="checkbox"/> REM500 | LiNX Color 3.5" Touch Screen<br>Display <sup>(1)(2)(3)</sup> ..... | \$1,500 |
|---------------------------------|--------------------------------------------------------------------|---------|

Select mounting position for REM500

|                                |                    |       |
|--------------------------------|--------------------|-------|
| <input type="checkbox"/> R5MR  | Right.....         | N/C   |
| <input type="checkbox"/> R5ML  | Left.....          | N/C   |
| <input type="checkbox"/> R5RPB | Rear Push Bar..... | \$140 |

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Driver Control not included.
3. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.

### SPECIALTY PROPORTIONAL CONTROLS<sup>(1)(2)</sup>

|                                  |                                                        |         |
|----------------------------------|--------------------------------------------------------|---------|
| <input type="checkbox"/> CREM    | LiNX Compact Remote HCPCS code E2373...                | \$1,500 |
| <input type="checkbox"/> CREM-LF | LiNX Compact Remote Low Force<br>HCPCS code E2373..... | \$1,500 |

1. Requires REM400 or REM500
2. Select Mount in Mounting For Driver Control section

### ATTENDANT DRIVE CONTROL CHOICES

|                              |                                                                        |       |
|------------------------------|------------------------------------------------------------------------|-------|
| <input type="checkbox"/> ACU | Proportional Attendant Control <sup>(1)</sup><br>HCPCS code E2331..... | \$900 |
|------------------------------|------------------------------------------------------------------------|-------|

Select mounting position for ACU

|                                |        |
|--------------------------------|--------|
| <input type="checkbox"/> RPBMR | Right  |
| <input type="checkbox"/> RPBML | Left   |
| <input type="checkbox"/> RPBMC | Center |

1. Requires Expandable Controller.

### THIRD PARTY SPECIALTY CONTROLS

Motion Concepts Install recommended

Installed at Motion Concepts: ☐ YES ☐ NO

☐ Type: \_\_\_\_\_

|                                                   |                                |                               |
|---------------------------------------------------|--------------------------------|-------------------------------|
| <input type="checkbox"/> Specialty Control Mount: | <input type="checkbox"/> Right | <input type="checkbox"/> Left |
| <input type="checkbox"/> Display Mount:           | <input type="checkbox"/> Right | <input type="checkbox"/> Left |
| <input type="checkbox"/> Mode Switch Mount:       | <input type="checkbox"/> Right | <input type="checkbox"/> Left |
| <input type="checkbox"/> On/Off Mount:            | <input type="checkbox"/> Right | <input type="checkbox"/> Left |

Notes

\_\_\_\_\_  
\_\_\_\_\_

# Ultra Low Maxx CG Power Positioning System

## TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### JOYSTICK MOUNT

#### Arm Mounted

- ☐ FR Fixed Joystick Mounting - Right ..... STD
- ☐ FL Fixed Joystick Mounting - Left ..... STD
- ☐ MRR Maxx Resolve Swing-Away Joystick Mount  
(Right Hand) HCPCS code E1028 .....\$395
- ☐ MRL Maxx Resolve Swing-Away Joystick Mount  
(Left Hand) HCPCS code E1028 .....\$395
- ☐ SAQR Motion Height Adjustable Swing-Away  
Quad Link (Right Hand)  
HCPCS code E1028 .....\$235
- ☐ SAQL Motion Height Adjustable Swing-Away  
Quad Link (Left Hand)  
HCPCS code E1028 .....\$235

#### Rear Push Bar Mounted

- ☐ RPBMR Right .....\$140
- ☐ RPBML Left .....\$140
- ☐ RPBMC Center .....\$140

### JOYSTICK KNOB - OPTIONAL

HCPCS code E2323

- ☐ PC101A-20 BodyPoint 3" U-Shaped .....\$120
- ☐ PC102A-20 BodyPoint 4" U-Shaped .....\$120
- ☐ PC104A-20 BodyPoint Straight Handle .....\$120
- ☐ PC107A-20 BodyPoint Domed Rubber .....\$120
- ☐ PC110A-20 BodyPoint Grooved Mushroom .....\$120

1. To remove the joystick knob from the REM400, pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

### ELECTRONICS - EXPANDABLE MODULES

- ☐ OUT500 Output Module ..... \$500
- ☐ IN200 LiNX Switch Input Module with  
Power ON/OFF for REM200 Series<sup>(1)(2)</sup> ..... \$500
- ☐ IN500H LiNX 9 Pin Input Module for  
Alternative Driver Controls<sup>(2)(3)(4)(6)</sup> ..... \$1,500
- ☐ IN500 LiNX 9 Pin Input Module with  
Sip-N-Puff Nozzle<sup>(2)(3)(4)(6)</sup> HCPCS code E2325. .... \$1,500
- ☐ PKG32666 Therafin Sip-N-Puff Breath Tube Kit ..... \$380  
HCPCS code E2326
- ☐ COVER-L Large Recline Shroud Cover<sup>(5)</sup> ..... N/C

1. Not for use with Alternative Drive Controls. (Uses switch ports).

2. Only one input module per chair.

3. Input module also required for any 9-pin specialty control.

4. 9 pin specialty control and Sip-N-Puff cannot both be used together.

5. May be required when choosing OUT500, IN200 or IN500 options.

6. Requires REM400 or REM500 for compatibility.

### LINX ELECTRONIC ACCESSORIES

- ☐ LAK LiNX Access Key<sup>(1)</sup> .....\$100

1. Required for programming

### NOTES

## Ultra Low Maxx CG Power Positioning System TDX-SP2

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

### REHAB SEAT SYSTEM<sup>(1)</sup>

☐ ULMRSN Ultra Low Maxx Rehab Seat.....N/C

1. Must order Rehab Seat Base (TDXSP2)

### POWER POSITIONING SYSTEMS<sup>(1)</sup>

☐ ULMT Ultra Low Maxx CG Tilt (50°) Only System<sup>(2)</sup>  
HCPCS code E1002.....\$6,595

☐ ULMTRE Ultra Low Maxx CG Tilt (50°) & Power Recline  
(170°) w/Extended Shear Reduction(ESR)<sup>(2)</sup>  
HCPCS code E1007.....\$10,795

☐ ULMTRP Ultra Low Maxx CG Tilt (50°) & Power Recline  
(170°) w/Power Sliding Back(PSB)<sup>(2)(3)</sup>  
HCPCS code E1008.....\$10,995

☐ ULMRE Ultra Low Maxx Power Recline (170°)  
w/Extended Shear Reduction(ESR)  
HCPCS code E1004.....\$5,450

☐ ULMRP Ultra Low Maxx Power Recline (170°)  
w/Power Sliding Back(PSB)<sup>(3)</sup>  
HCPCS code E1005.....\$6,250

☐ ULME Ultra Low Maxx 12" Elevating  
Lift Seat Only<sup>(4)</sup> HCPCS code E2298.....\$4,995  
Weight capacity up to 250lbs

☐ ULME10 Ultra Low Maxx 10" Elevating  
Lift Seat Only<sup>(5)</sup> HCPCS code E2298.....\$5,995  
Weight capacity up to 300lbs

1. Not available with Rehab Seat Base TDXSP2.
2. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module
3. Requires Multi-Function M610 through Drive Control (MFCBD). Not available on Multi Function Act4-LiNX Electronics.
4. Not available on seat sizes greater than 20" wide by 20" deep.
5. Not available on seat sizes greater than 22" wide by 22" deep. Not available on Multi Function Act4-LiNX Electronics.

### ADD: POWER ADJUSTABLE SEAT HEIGHT MODULE

☐ ESM12 12" Power Adjustable Seat Height Module<sup>(1)(2)(3)</sup>  
HCPCS code E2298.....\$2,495  
Weight capacity up to 250lbs

☐ ESM10 10" HD Power Adjustable Seat Height Module<sup>(4)(5)</sup>  
HCPCS code E2298.....\$2,995  
Weight capacity up to 300lbs

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.
2. Not available on seat sizes greater than 20" wide by 20" deep
3. For weights greater than 175lbs tilt while elevated is not available.
4. Only available with tilt, uses an integrated Tilt/Lift mechanism, tilt is 45°.
5. Not available on seat sizes greater than 22" wide by 22" deep. Only available on Standard Base (STDfrm). Not available on Multi Function Act4-LiNX Electronics.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### MODULES - MUST SELECT SYSTEM TO ADD MODULES

☐ PRETILT 5° Fixed Pre-Tilt<sup>(1)</sup>.....N/C  
☐ PRM8 8° Precline Module<sup>(2)</sup>.....N/C  
☐ PRM30 30° Precline Module<sup>(2)</sup>  
HCPCS code K0108.....\$500  
☐ ERSSP Programable Electronic Recline Function  
Start & Stop<sup>(3)</sup>.....\$400

Specify:

Start: \_\_\_\_\_degrees Stop: \_\_\_\_\_degrees

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height.
2. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available. Must order Recline System
3. Requires Multi-Function M610 through Drive Control (MFCBD). Set range between 90° and 170°. Will be set to factory settings unless specified.

### ADULT SEAT WIDTH SETTINGS

☐ SW116 16" Wide (Adjustable from 16" to 20").....N/C  
☐ SW117 17" Wide (Adjustable from 16" to 20").....N/C  
☐ SW118 18" Wide (Adjustable from 16" to 20").....N/C  
☐ SW119 19" Wide (Adjustable from 16" to 20").....N/C  
☐ SW120 20" Wide (Adjustable from 16" to 20").....N/C  
☐ SW219 19" Wide (Adjustable from 19" to 22").....N/C  
☐ SW220 20" Wide (Adjustable from 19" to 22").....N/C  
☐ SW221 21" Wide (Adjustable from 19" to 22").....N/C  
☐ SW222 22" Wide (Adjustable from 19" to 22").....N/C

### ADULT ULTRA RAIL SEAT DEPTH SETTINGS

☐ SD116 16" Deep (Adjustable from 16" to 21").....N/C  
☐ SD117 17" Deep (Adjustable from 16" to 21").....N/C  
☐ SD118 18" Deep (Adjustable from 16" to 21").....N/C  
☐ SD119 19" Deep (Adjustable from 16" to 21").....N/C  
☐ SD120 20" Deep (Adjustable from 16" to 21").....N/C  
☐ SD121 21" Deep (Adjustable from 16" to 21").....N/C  
☐ SD219 19" Deep (Adjustable from 19" to 23").....N/C  
☐ SD220 20" Deep (Adjustable from 19" to 23").....N/C  
☐ SD221 21" Deep (Adjustable from 19" to 23").....N/C  
☐ SD222 22" Deep (Adjustable from 19" to 23").....N/C  
☐ SD223 23" Deep (Adjustable from 19" to 23").....N/C

### SMALL ADULT SEAT WIDTH SETTINGS<sup>(1)</sup>

☐ SW314 14" Wide (Adjustable from 14" to 15").....N/C  
☐ SW315 15" Wide (Adjustable from 14" to 15").....N/C

1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

### SMALL ADULT SEAT DEPTH SETTINGS<sup>(1)</sup>

☐ SD314 14" Deep (Adjustable from 14" to 18").....N/C  
☐ SD315 15" Deep (Adjustable from 14" to 18").....N/C  
☐ SD316 16" Deep (Adjustable from 14" to 18").....N/C  
☐ SD317 17" Deep (Adjustable from 14" to 18").....N/C  
☐ SD318 18" Deep (Adjustable from 14" to 18").....N/C

1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

# Ultra Low Maxx CG Power Positioning System

## TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### MATRIX SEAT CUSHION OPTIONS

Matrix PS Seat Cushion 14" - 20" wide/reversable cover  
HCPCS code E2605.....\$345

☐ PS1414    ☐ PS1416    ☐ PS1616    ☐ PS1618    ☐ PS1620  
☐ PS1716    ☐ PS1718    ☐ PS1720    ☐ PS1816    ☐ PS1818  
☐ PS1820    ☐ PS2016    ☐ PS2018    ☐ PS2020

Matrix PS Seat Cushion 21"-22" wide  
HCPCS code E2606.....\$545  
☐ PS2218    ☐ PS2220    ☐ PS2222

Matrix PSP Seat Cushion 14" - 20" wide/reversable cover  
HCPCS code E2607.....\$395  
☐ PSP1414    ☐ PSP1416    ☐ PSP1616    ☐ PSP1618    ☐ PSP1620  
☐ PSP1716    ☐ PSP1718    ☐ PSP1720    ☐ PSP1816    ☐ PSP1818  
☐ PSP1820    ☐ PSP2016    ☐ PSP2018    ☐ PSP2020

Matrix PSP Seat Cushion 21"-22" wide  
HCPCS code E2608.....\$595  
☐ PSP2218    ☐ PSP2220    ☐ PSP2222

Matrix Vi Seat Cushion 14" - 20" wide  
HCPCS code E2607.....\$425  
☐ KA1414-VI    ☐ MA1415-VI    ☐ MA1417-VI    ☐ MA1616-VI    ☐ MA1618-VI  
☐ MA1620-VI    ☐ MA1716-VI    ☐ MA1718-VI    ☐ MA1720-VI    ☐ MA1816-VI  
☐ MA1818-VI    ☐ MA1820-VI    ☐ MA2016-VI    ☐ MA2018-VI  
☐ MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide  
HCPCS code E2608.....\$475  
☐ MA2218-VI    ☐ MA2220-VI

Matrix Vi HD Seat Cushion 22" wide  
HCPCS code E2608.....\$625  
☐ MA2222-VI

Matrix Libra Seat Cushion 14" - 20" wide  
HCPCS code E2624.....\$485  
☐ LC1414    ☐ LC1416    ☐ LC1516    ☐ LC1518    ☐ LC1520  
☐ LC1616    ☐ LC1618    ☐ LC1620    ☐ LC1716    ☐ LC1718  
☐ LC1720    ☐ LC1816    ☐ LC1818    ☐ LC1820    ☐ LC1916  
☐ LC1918    ☐ LC1920    ☐ LC2016    ☐ LC2018    ☐ LC2020

Matrix Libra HD Seat Cushion 21" wide  
HCPCS code E2624.....\$525  
☐ LC2118    ☐ LC2120

Matrix Libra Seat Cushion 22" deep  
HCPCS code E2624.....\$595  
☐ LC1622    ☐ LC1722    ☐ LC1822    ☐ LC1922    ☐ LC2022  
☐ LC2122

Matrix Libra HD Seat Cushion 22" wide  
HCPCS code E2625.....\$595  
☐ LC2218    ☐ LC2220    ☐ LC2222

### MATRIX SEAT CUSHION OPTIONS

Matrix Libra Adjust Seat Cushion 14" -20" wide  
HCPCS code E2624.....\$625

☐ LC1414-ADJ    ☐ LC1416-ADJ    ☐ LC1516-ADJ    ☐ LC1518-ADJ  
☐ LC1520-ADJ    ☐ LC1616-ADJ    ☐ LC1618-ADJ    ☐ LC1620-ADJ  
☐ LC1716-ADJ    ☐ LC1718-ADJ    ☐ LC1720-ADJ    ☐ LC1816-ADJ  
☐ LC1818-ADJ    ☐ LC1820-ADJ    ☐ LC1916-ADJ    ☐ LC1918-ADJ  
☐ LC1920-ADJ    ☐ LC2016-ADJ    ☐ LC2018-ADJ    ☐ LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" deep  
HCPCS code E2624.....\$735  
☐ LC1622-ADJ    ☐ LC1722-ADJ    ☐ LC1822-ADJ    ☐ LC1922-ADJ  
☐ LC2022-ADJ    ☐ LC2122-ADJ

Matrix Libra Adjust HD Seat Cushion 22" wide  
HCPCS code E2625.....\$735  
☐ LC2118-ADJ    ☐ LC2120-ADJ

Matrix Multi Cushion 14" - 20" wide w/reversable cover  
HCPCS code E2607.....\$445  
☐ MM1414    ☐ MM1416    ☐ MM1616    ☐ MM1618    ☐ MM1620  
☐ MM1716    ☐ MM1718    ☐ MM1720    ☐ MM1816    ☐ MM1818  
☐ MM1820    ☐ MM2016    ☐ MM2018    ☐ MM2020

Matrix Multi Cushion 22" wide w/reversable cover  
HCPCS code E2608.....\$595  
☐ MM2216    ☐ MM2218    ☐ MM2220    ☐ MM2222

☐ OCCC "ON-CHAIR" Cushion Credit<sup>(1)</sup> ..... (\$90)

1. Not available on Matrix PS or PSP Cushions.

### FABRIC OPTIONS FOR CUSHIONS

| Cushion      | Meshtex | Startex        | Infection Control | Spacetex | Onyx  | Fusion |
|--------------|---------|----------------|-------------------|----------|-------|--------|
| Libra/Adjust | N/A     | STD            | \$50              | N/A      | N/A   | N/A    |
| PS/PSP       | N/A     | STD Reversible | STD Reversible    | N/A      | N/A   | N/A    |
| Vi           | STD     | \$50           | \$50              | N/A      | \$100 | \$150  |

☐ LC-SP Spacetex Fabric Option for  
Libra Cushion .....N/C

☐ LC-IC Infection Control Fabric Upgrade for  
Libra Cushion .....\$50

☐ FSU-MA Fusion Fabric Upgrade for Matrix Vi Cushion....\$150

☐ OCU-MA Onyx Fabric Upgrade for Matrix Vi Cushion ....\$100

☐ STARTEX-S Smooth side Startex Fabric Upgrade for  
Matrix Vi Cushion .....\$50

☐ STARTEX-F Fabric side Startex Fabric Upgrade for  
Matrix Vi Cushion .....\$50

## Ultra Low Maxx CG Power Positioning System TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### LAP BELTS

|                                  |                                                                                     |
|----------------------------------|-------------------------------------------------------------------------------------|
| <input type="checkbox"/> BELT48  | Push Button Style Lap Belt, Length 48" <sup>(1)</sup> . . . . . STD                 |
| <input type="checkbox"/> BELT60  | Push Button Style Lap Belt, Length 60" <sup>(1)</sup> . . . . . STD                 |
| <input type="checkbox"/> BELT71  | Push Button Style Lap Belt, Length 71" <sup>(1)</sup> . . . . . STD                 |
| <input type="checkbox"/> PLB2    | Bodypoint Padded Lap Belt 2 Point <sup>(2)</sup><br>HCPCS code E0978. . . . . \$190 |
| <input type="checkbox"/> PLB4    | Bodypoint Padded Lap Belt 4 Point <sup>(2)</sup><br>HCPCS code E0978. . . . . \$235 |
| <input type="checkbox"/> TROBELT | Safety Belt for use only with TRRO & TRRO-E . . . N/C                               |

1. Not available with TRRO and TRRO-E.

2. 67" length only option

### MATRIX HEADREST PADS<sup>(1)</sup>

HCPCS code E0955

#### Standard Pads

|                                    |                                                  |
|------------------------------------|--------------------------------------------------|
| <input type="checkbox"/> ESP6      | Matrix Elan Standard Pad 6" . . . . . \$235      |
| <input type="checkbox"/> ESP10     | Matrix Elan Standard Pad 10" . . . . . \$235     |
| <input type="checkbox"/> ESP14     | Matrix Elan Standard Pad 14" . . . . . \$235     |
| <input type="checkbox"/> ESP6-HUG  | Matrix Elan Hug Standard Pad 6" . . . . . \$235  |
| <input type="checkbox"/> ESP10-HUG | Matrix Elan Hug Standard Pad 10" . . . . . \$235 |
| <input type="checkbox"/> ESP14-HUG | Matrix Elan Hug Standard Pad 14" . . . . . \$235 |
| <input type="checkbox"/> SNUG14    | Matrix SNUG Standard Pad 14" . . . . . \$415     |
| <input type="checkbox"/> SNUG14-A  | Matrix SNUG Asymmetrical Pad 14" . . . . . \$415 |
| <input type="checkbox"/> SNUG18    | Matrix SNUG Standard Pad 18" . . . . . \$415     |
| <input type="checkbox"/> SNUG18-A  | Matrix SNUG Asymmetrical Pad 18" . . . . . \$415 |

#### Standard Pads with Infection Control Cover

|                                       |                                                                       |
|---------------------------------------|-----------------------------------------------------------------------|
| <input type="checkbox"/> ESP6-IC      | Matrix Elan Infection Control Pad 6" . . . . . \$235                  |
| <input type="checkbox"/> ESP10-IC     | Matrix Elan Infection Control Pad 10" . . . . . \$235                 |
| <input type="checkbox"/> ESP14-IC     | Matrix Elan Infection Control Pad 14" . . . . . \$235                 |
| <input type="checkbox"/> ESP6-IC-HUG  | Matrix Elan Hug Infection Control Pad 6" . . . . . \$235              |
| <input type="checkbox"/> ESP10-IC-HUG | Matrix Elan Hug Infection Control Pad 10" . . . . . \$235             |
| <input type="checkbox"/> ESP14-IC-HUG | Matrix Elan Hug Infection Control Pad 14" . . . . . \$235             |
| <input type="checkbox"/> SNUG14-IC    | Matrix SNUG Standard<br>Infection Control Pad 14" . . . . . \$445     |
| <input type="checkbox"/> SNUG14-A-IC  | Matrix SNUG Asymmetrical<br>Infection Control Pad 14" . . . . . \$445 |
| <input type="checkbox"/> SNUG18-IC    | Matrix SNUG Standard<br>Infection Control Pad 18" . . . . . \$445     |
| <input type="checkbox"/> SNUG18-A-IC  | Matrix SNUG Asymmetrical<br>Infection Control Pad 18" . . . . . \$445 |

#### Occipital /4-Point Pads

|                                  |                                                     |
|----------------------------------|-----------------------------------------------------|
| <input type="checkbox"/> EOP9    | Matrix Elan Occipital Pad 9" . . . . . \$235        |
| <input type="checkbox"/> EOP12   | Matrix Elan Occipital Pad 12" . . . . . \$235       |
| <input type="checkbox"/> E4PTSM  | Matrix Elan 4-Point Pad 9"W x 7"H . . . . . \$325   |
| <input type="checkbox"/> E4POINT | Matrix Elan 4-Point Pad 11"W x 10"H . . . . . \$325 |

#### Occipital / 4-Point Pads with Infection Control Cover

|                                     |                                                                               |
|-------------------------------------|-------------------------------------------------------------------------------|
| <input type="checkbox"/> EOP9-IC    | Matrix Elan Occipital Pad 9"<br>with Infection Control. . . . . \$235         |
| <input type="checkbox"/> EOP12-IC   | Matrix Elan Occipital Pad 12"<br>with Infection Control. . . . . \$235        |
| <input type="checkbox"/> E4PTSM-IC  | Matrix Elan 4-Point Pad 9"W x 7"H<br>with Infection Control. . . . . \$325    |
| <input type="checkbox"/> E4POINT-IC | Matrix Elan 4-Point Pad 11"W x 10"H<br>with Infection Control . . . . . \$325 |

1. If mounting on PB Elite must order removable hardware

### MOTION CONCEPTS HEADREST PAD

HCPCS code E0955

|                               |                                                           |
|-------------------------------|-----------------------------------------------------------|
| <input type="checkbox"/> MCSH | Motion Concepts Standard . . . . . \$295                  |
| <input type="checkbox"/> MCOH | Motion Concepts Auto Style <sup>(1)</sup> . . . . . \$350 |

1. Comes Standard with Fixed Mounting hardware; for use on Matrix PB or Elite  
Backs you must select the removable hardware upgrade below.

### HEADREST MOUNTING HARDWARE

HCPCS code E1028


|                                     |                                                   |
|-------------------------------------|---------------------------------------------------|
| <input type="checkbox"/> FM         | Fixed Mounting <sup>(1)</sup> . . . . . N/C       |
| <input type="checkbox"/> LXSHW      | Loxx Standard Hardware 15" . . . . . \$215        |
| <input type="checkbox"/> LXMMH      | Loxx Mini Hardware 11" . . . . . \$215            |
| <input type="checkbox"/> MEHW       | Elan Hardware 15" . . . . . \$220                 |
| <input type="checkbox"/> MEMH       | Elan Mini Hardware 11" . . . . . \$225            |
| <input type="checkbox"/> MEHW-EXT   | Elan Extended Hardware . . . . . \$295            |
| <input type="checkbox"/> FPLXSHW    | FLIPP Loxx Standard Hardware 15" . . . . . \$515  |
| <input type="checkbox"/> FPLXSHW-EX | FLIPP Loxx Extended Hardware . . . . . \$525      |
| <input type="checkbox"/> FPLXMMH    | FLIPP Loxx Mini Hardware 11" . . . . . \$515      |
| <input type="checkbox"/> FPLXMMH-EX | FLIPP Loxx Mini Extended Hardware . . . . . \$525 |
| <input type="checkbox"/> FPMEHW     | FLIPP Elan Hardware 15" . . . . . \$525           |
| <input type="checkbox"/> FPMEH      | FLIPP Elan Mini Hardware 11" . . . . . \$525      |
| <input type="checkbox"/> FPMEHW-EXT | FLIPP Elan Extended Hardware . . . . . \$605      |

|                              |                              |
|------------------------------|------------------------------|
| <input type="checkbox"/> OHR | OMIT: Headrest . . . . . N/C |
|------------------------------|------------------------------|

1. HCPCS code does not apply.

### TILT ONLY BACK CANE/SPREADER BAR OPTION<sup>(1)</sup>

|                                |                                                                                                 |
|--------------------------------|-------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> MSBSB | Maxx Style Back Adjustable Spreader Bar . . . . . STD<br>Adjust height easily on back canes     |
| <input type="checkbox"/> CSBSB | Classic Style Back Adjustable Spreader Bar(2) . . \$100                                         |
| <input type="checkbox"/> OMR   | Optional Multi-Rail. . . . . \$150<br>Standard Rehab Back with Classic Style Spreader Bar Only. |



## Ultra Low Maxx CG Power Positioning System TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### STANDARD BACK CANES - FOR TILT SYSTEMS

|                                |                         |     |
|--------------------------------|-------------------------|-----|
| <input type="checkbox"/> SBC18 | 18" Straight Back Canes | STD |
| <input type="checkbox"/> SBC20 | 20" Straight Back Canes | STD |
| <input type="checkbox"/> SBC22 | 22" Straight Back Canes | STD |
| <input type="checkbox"/> SBC24 | 24" Straight Back Canes | STD |

#### Select Cane Angle

|                                     |                                     |                                     |                                     |                                   |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-----------------------------------|
| <input type="checkbox"/> BA75 75°   | <input type="checkbox"/> BA79 79°   | <input type="checkbox"/> BA85 85°   | <input type="checkbox"/> BA90 90°   | <input type="checkbox"/> BA95 95° |
| <input type="checkbox"/> BA101 101° | <input type="checkbox"/> BA106 106° | <input type="checkbox"/> BA112 112° | <input type="checkbox"/> BA116 116° |                                   |

|                                |                            |     |
|--------------------------------|----------------------------|-----|
| <input type="checkbox"/> ABC20 | 20" Angled (5°) Back Canes | N/C |
| <input type="checkbox"/> ABC22 | 22" Angled (5°) Back Canes | N/C |
| <input type="checkbox"/> ABC24 | 24" Angled (5°) Back Canes | N/C |

#### Select Cane Angle

|                                     |                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> BA81 81°   | <input type="checkbox"/> BA84 84°   | <input type="checkbox"/> BA91 91°   | <input type="checkbox"/> BA96 96°   | <input type="checkbox"/> BA101 101° |
| <input type="checkbox"/> BA106 106° | <input type="checkbox"/> BA111 111° | <input type="checkbox"/> BA117 117° | <input type="checkbox"/> BA121 121° |                                     |

|                                |                                         |       |
|--------------------------------|-----------------------------------------|-------|
| <input type="checkbox"/> EBC22 | 22" Endomorph Back Canes <sup>(1)</sup> | \$350 |
| <input type="checkbox"/> EBC24 | 24" Endomorph Back Canes <sup>(1)</sup> | \$350 |

#### Please Select

|                                |             |     |
|--------------------------------|-------------|-----|
| <input type="checkbox"/> EBCN1 | 1" Narrower | N/C |
| <input type="checkbox"/> EBCN2 | 2" Narrower | N/C |
| <input type="checkbox"/> EBCN4 | 4" Narrower | N/C |

|                                |                                                  |       |
|--------------------------------|--------------------------------------------------|-------|
| <input type="checkbox"/> TBC18 | 18" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup> | \$250 |
| <input type="checkbox"/> TBC20 | 20" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup> | \$250 |
| <input type="checkbox"/> TBC22 | 22" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup> | \$250 |
| <input type="checkbox"/> TBC24 | 24" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup> | \$250 |

|                               |                                                 |       |
|-------------------------------|-------------------------------------------------|-------|
| <input type="checkbox"/> FDBC | Fold Down Back Canes w/ Pull Pin <sup>(6)</sup> | \$250 |
|-------------------------------|-------------------------------------------------|-------|

1. Comes standard with Cantilever Armrest Cane Mounted.
2. Quick Height Adjustment on Backrest and Maxx Spreader Bars not available.
3. Minimum finished back widths of 16".
4. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
5. 4" of total vertical mounting adjustment.
6. 250 lbs weight capacity. Must select canes above.

### CANE OPTIONS - FOR RECLINE SYSTEMS

|                               |                                                           |       |
|-------------------------------|-----------------------------------------------------------|-------|
| <input type="checkbox"/> BFS2 | Recline Body Fit 2" Narrower <sup>(1)</sup>               | \$995 |
| <input type="checkbox"/> BFS4 | Recline Body Fit 4" Narrower <sup>(1)</sup>               | \$995 |
| <input type="checkbox"/> AHH  | Adjustable Hooking Handles <sup>(2)</sup> (pair)          | \$450 |
| <input type="checkbox"/> RPHS | Push Handles for Recline - Straight <sup>(3)</sup> (pair) | \$350 |
| <input type="checkbox"/> RPHA | Push Handles for Recline - Angled <sup>(3)</sup> (pair)   | \$350 |

1. Minimum finished back widths of 16".
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

### FINISHED BACK HEIGHT<sup>(1)(2)</sup> - FOR TILT ONLY/ELEVATE ONLY

|                                    |                                    |                                    |                                    |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| <input type="checkbox"/> FBT18 18" | <input type="checkbox"/> FBT19 19" | <input type="checkbox"/> FBT20 20" | <input type="checkbox"/> FBT21 21" |
| <input type="checkbox"/> FBT22 22" | <input type="checkbox"/> FBT23 23" | <input type="checkbox"/> FBT24 24" | <input type="checkbox"/> FBT25 25" |

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

### FINISHED BACK HEIGHT<sup>(1)(2)</sup> - FOR RECLINE OR TILT/RECLINE

|                                    |                                    |                                    |                                    |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| <input type="checkbox"/> FBR20 20" | <input type="checkbox"/> FBR21 21" | <input type="checkbox"/> FBR22 22" | <input type="checkbox"/> FBR23 23" |
| <input type="checkbox"/> FBR24 24" | <input type="checkbox"/> FBR25 25" | <input type="checkbox"/> FBR26 26" | <input type="checkbox"/> FBR27 27" |

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

### STANDARD BACK OPTIONS

|                               |                                                                |        |
|-------------------------------|----------------------------------------------------------------|--------|
| <input type="checkbox"/> SRB  | Standard Rehab Back & Cushion                                  | N/C    |
| <input type="checkbox"/> SS   | Additional Super Soft/HR Foam <sup>(1)</sup>                   | \$450  |
| <input type="checkbox"/> MFU  | Matrx Fabric Upcharge                                          | \$450  |
| <input type="checkbox"/> OBC  | OMIT: Back Cushion                                             | (\$25) |
| <input type="checkbox"/> PIPF | Flat Planar Interface/I-Back Plate <sup>(2)(3)(5)</sup>        | \$125  |
| <input type="checkbox"/> PIPR | Recessed Planar/I-Back Interface Plate <sup>(2)(3)(5)(6)</sup> | \$125  |
| <input type="checkbox"/> ORBC | OMIT: Rigid Back & Cushion <sup>(4)</sup>                      | (\$75) |

1. Additional foam added may impact seat depth
2. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
3. Required with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
4. Not available with Recline/ESR or Recline/PSB.
5. Not available with TRRO and TRRO-E.
6. Fits Comfort Company, Freedom, Pindot and Moulded Backs.

### MOTION CONCEPTS SINGLE POST BACKREST<sup>(1)(2)(3)</sup>

|                                |                                                                                     |                                 |
|--------------------------------|-------------------------------------------------------------------------------------|---------------------------------|
| <input type="checkbox"/> SPBMI | Single Post Angle Adjustable Fold Down<br>Back interface-Matrx Back.                | \$750<br>(Folds to 25° forward) |
| <input type="checkbox"/> SPBCI | Single Post Angle Adjustable Fold Down<br>Back for Contour/High Back <sup>(4)</sup> | \$750<br>(Folds to 25° forward) |

1. Single Post Back is an angle adjustable back between 73° and 109° in 6° increments from 17° forward to 19° back. The back folds down and is only available for 16"-20" seat widths & 21"-24" back height. A 4" gap from bottom of back to seat pan is required.
2. This option is not available with Recline and/or a vent tray and has a 250lb weight capacity.
3. Cane mounted armrest are not available with this option.
4. Not available with TRRO/TRRO-E

### MOTION CONCEPTS REHAB SEAT SLING BACK<sup>(1)</sup>

|                                    |                                   |         |
|------------------------------------|-----------------------------------|---------|
| <input type="checkbox"/> SLINGB-TT | Sling Back with Tension Fabric    | \$950   |
| <input type="checkbox"/> SLINGB-TA | Sling Back with Adjustable Fabric | \$1,250 |

1. Only available on Rehab Seat System. Not available with Elevate. Cane height must be 2" taller than Finished Back Height.

### FREEDOM BACK AND SEAT INTERFACE HARDWARE<sup>(1)</sup>

#### Freedom Back Interface Hardware for Tilt Only

|                                 |                |                                                      |
|---------------------------------|----------------|------------------------------------------------------|
| <input type="checkbox"/> FBIFEM | Flush Mount    | \$125<br>Back can be seat width or up to 4" narrower |
| <input type="checkbox"/> FBITRM | Recessed Mount | \$125<br>Back must be 4" less than seat width        |

#### Freedom Back Interface Hardware with Recline

|                                 |                |                                                      |
|---------------------------------|----------------|------------------------------------------------------|
| <input type="checkbox"/> FBIRFM | Flush Mount    | \$125<br>Back can be seat width or up to 5" narrower |
| <input type="checkbox"/> FBIRRM | Recessed Mount | \$125<br>Back must be 5" less than seat width        |

|                               |                                 |                                                      |
|-------------------------------|---------------------------------|------------------------------------------------------|
| <input type="checkbox"/> FSIP | Freedom Seat Interface Hardware | \$125<br>Seat can be up to 3" less than system width |
|-------------------------------|---------------------------------|------------------------------------------------------|

1. Not available with TRRO and TRRO-E.

### RIDE STYLE MULTI-AXIS INTERFACE<sup>(1)</sup>

|                               |                                                                        |       |
|-------------------------------|------------------------------------------------------------------------|-------|
| <input type="checkbox"/> RSMI | Ride Style Multi-Axis Interface for Recline<br>Systems 14" to 20" wide | \$250 |
|-------------------------------|------------------------------------------------------------------------|-------|

1. Only available with TRRO and TRRO-E when using with ELITE back.

**Ultra Low Maxx CG Power Positioning System**  
TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**UPGRADABLE BACK OPTIONS**

Contour Back 15" wide - 20" wide

HCPCS code E2613 ..... \$495

- |                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> CB1512 | <input type="checkbox"/> CB1515 | <input type="checkbox"/> CB1517 | <input type="checkbox"/> CB1612 |
| <input type="checkbox"/> CB1615 | <input type="checkbox"/> CB1617 | <input type="checkbox"/> CB1621 | <input type="checkbox"/> CB1712 |
| <input type="checkbox"/> CB1715 | <input type="checkbox"/> CB1717 | <input type="checkbox"/> CB1721 | <input type="checkbox"/> CB1812 |
| <input type="checkbox"/> CB1815 | <input type="checkbox"/> CB1817 | <input type="checkbox"/> CB1821 | <input type="checkbox"/> CB1912 |
| <input type="checkbox"/> CB1915 | <input type="checkbox"/> CB1917 | <input type="checkbox"/> CB1921 | <input type="checkbox"/> CB2012 |
| <input type="checkbox"/> CB2015 | <input type="checkbox"/> CB2017 | <input type="checkbox"/> CB2021 |                                 |

Contour Back 21"-22" wide

HCPCS code E2613/E2614 ..... \$595

- |                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> CB2115 | <input type="checkbox"/> CB2117 | <input type="checkbox"/> CB2121 | <input type="checkbox"/> CB2215 |
| <input type="checkbox"/> CB2217 | <input type="checkbox"/> CB2221 |                                 |                                 |

High Back 16" - 22" W

HCPCS code E2620 ..... \$595

- |                                 |                                 |                                 |                                 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> HB1616 | <input type="checkbox"/> HB1618 | <input type="checkbox"/> HB1620 | <input type="checkbox"/> HB1716 |
| <input type="checkbox"/> HB1718 | <input type="checkbox"/> HB1720 | <input type="checkbox"/> HB1816 | <input type="checkbox"/> HB1818 |
| <input type="checkbox"/> HB1820 | <input type="checkbox"/> HB1916 | <input type="checkbox"/> HB1918 | <input type="checkbox"/> HB1920 |
| <input type="checkbox"/> HB2016 | <input type="checkbox"/> HB2018 | <input type="checkbox"/> HB2020 | <input type="checkbox"/> HB2116 |
| <input type="checkbox"/> HB2118 | <input type="checkbox"/> HB2120 | <input type="checkbox"/> HB2216 | <input type="checkbox"/> HB2218 |
| <input type="checkbox"/> HB2220 |                                 |                                 |                                 |

Matrix Elite Back 14" - 20" W (3" Contour)

HCPCS code E2620 ..... \$625

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PBE1414 | <input type="checkbox"/> PBE1416 | <input type="checkbox"/> PBE1418 | <input type="checkbox"/> PBE1420 |
| <input type="checkbox"/> PBE1514 | <input type="checkbox"/> PBE1516 | <input type="checkbox"/> PBE1518 | <input type="checkbox"/> PBE1520 |
| <input type="checkbox"/> PBE1614 | <input type="checkbox"/> PBE1616 | <input type="checkbox"/> PBE1618 | <input type="checkbox"/> PBE1620 |
| <input type="checkbox"/> PBE1714 | <input type="checkbox"/> PBE1716 | <input type="checkbox"/> PBE1718 | <input type="checkbox"/> PBE1720 |
| <input type="checkbox"/> PBE1814 | <input type="checkbox"/> PBE1816 | <input type="checkbox"/> PBE1818 | <input type="checkbox"/> PBE1820 |
| <input type="checkbox"/> PBE1916 | <input type="checkbox"/> PBE1918 | <input type="checkbox"/> PBE1920 | <input type="checkbox"/> PBE2014 |
| <input type="checkbox"/> PBE2016 | <input type="checkbox"/> PBE2018 | <input type="checkbox"/> PBE2020 |                                  |

Matrix Elite HD Back 21" - 22" W (3" Contour)

HCPCS code E2620 for 21" E2621 for 22" ..... \$795

- |                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> PBE2116-HD | <input type="checkbox"/> PBE2120-HD | <input type="checkbox"/> PBE2216-HD | <input type="checkbox"/> PBE2220-HD |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|

Matrix Elite TR Back 15" - 20" W (3" Contour)

HCPCS code E2620 ..... \$650

- |                                     |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> PBE1516-TR | <input type="checkbox"/> PBE1518-TR | <input type="checkbox"/> PBE1520-TR | <input type="checkbox"/> PBE1616-TR |
| <input type="checkbox"/> PBE1618-TR | <input type="checkbox"/> PBE1620-TR | <input type="checkbox"/> PBE1716-TR | <input type="checkbox"/> PBE1718-TR |
| <input type="checkbox"/> PBE1720-TR | <input type="checkbox"/> PBE1816-TR | <input type="checkbox"/> PBE1818-TR | <input type="checkbox"/> PBE1820-TR |
| <input type="checkbox"/> PBE1916-TR | <input type="checkbox"/> PBE1918-TR | <input type="checkbox"/> PBE1920-TR | <input type="checkbox"/> PBE2016-TR |
| <input type="checkbox"/> PBE2018-TR | <input type="checkbox"/> PBE2020-TR |                                     |                                     |

Matrix Elite TR HD Back 21" - 22" W (3" Contour)

HCPCS code E2620 for 21" E2621 for 22" ..... \$845

- |                                       |                                       |
|---------------------------------------|---------------------------------------|
| <input type="checkbox"/> PBE2120-TRHD | <input type="checkbox"/> PBE2220-TRHD |
|---------------------------------------|---------------------------------------|

**UPGRADABLE BACK OPTIONS**

Matrix Elite Deep Back 14" - 20" W (6" Contour)

HCPCS code E2620 ..... \$760

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> EDB1416 | <input type="checkbox"/> EDB1418 | <input type="checkbox"/> EDB1420 | <input type="checkbox"/> EDB1514 |
| <input type="checkbox"/> EDB1516 | <input type="checkbox"/> EDB1518 | <input type="checkbox"/> EDB1520 | <input type="checkbox"/> EDB1614 |
| <input type="checkbox"/> EDB1616 | <input type="checkbox"/> EDB1618 | <input type="checkbox"/> EDB1620 | <input type="checkbox"/> EDB1714 |
| <input type="checkbox"/> EDB1716 | <input type="checkbox"/> EDB1718 | <input type="checkbox"/> EDB1720 | <input type="checkbox"/> EDB1814 |
| <input type="checkbox"/> EDB1816 | <input type="checkbox"/> EDB1818 | <input type="checkbox"/> EDB1820 | <input type="checkbox"/> EDB1914 |
| <input type="checkbox"/> EDB1916 | <input type="checkbox"/> EDB1918 | <input type="checkbox"/> EDB1920 | <input type="checkbox"/> EDB2014 |
| <input type="checkbox"/> EDB2016 | <input type="checkbox"/> EDB2018 | <input type="checkbox"/> EDB2020 |                                  |

Matrix Elite Deep HD Back 21" - 22" W (6" Contour)

HCPCS code E2620 for 21" E2621 for 22" ..... \$930

- |                                     |                                     |                                     |
|-------------------------------------|-------------------------------------|-------------------------------------|
| <input type="checkbox"/> EDB2116-HD | <input type="checkbox"/> EDB2118-HD | <input type="checkbox"/> EDB2120-HD |
| <input type="checkbox"/> EDB2216-HD | <input type="checkbox"/> EDB2218-HD | <input type="checkbox"/> EDB2220-HD |

Matrix Elite E2S Back 14" - 20" W x 14"/16"H (3" Contour)

HCPCS code E2615 ..... \$625

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2S1414 | <input type="checkbox"/> E2S1416 | <input type="checkbox"/> E2S1514 | <input type="checkbox"/> E2S1516 |
| <input type="checkbox"/> E2S1614 | <input type="checkbox"/> E2S1616 | <input type="checkbox"/> E2S1714 | <input type="checkbox"/> E2S1716 |
| <input type="checkbox"/> E2S1814 | <input type="checkbox"/> E2S1816 | <input type="checkbox"/> E2S1914 | <input type="checkbox"/> E2S1916 |
| <input type="checkbox"/> E2S2014 | <input type="checkbox"/> E2S2016 |                                  |                                  |

Matrix Elite E2S Back 14" - 20" W x 18"/20"H (3" Contour)

HCPCS code E2615 ..... \$645

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2S1418 | <input type="checkbox"/> E2S1420 | <input type="checkbox"/> E2S1518 | <input type="checkbox"/> E2S1520 |
| <input type="checkbox"/> E2S1618 | <input type="checkbox"/> E2S1620 | <input type="checkbox"/> E2S1718 | <input type="checkbox"/> E2S1720 |
| <input type="checkbox"/> E2S1818 | <input type="checkbox"/> E2S1820 | <input type="checkbox"/> E2S1918 | <input type="checkbox"/> E2S1920 |
| <input type="checkbox"/> E2S2018 | <input type="checkbox"/> E2S2020 |                                  |                                  |

Matrix Elite E2 Lateral Back

14" - 20" W x 14"/16"H (3" Contour) HCPCS code E2620 ..... \$725

- |                                   |                                   |                                   |                                   |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> E2LB1414 | <input type="checkbox"/> E2LB1416 | <input type="checkbox"/> E2LB1514 | <input type="checkbox"/> E2LB1516 |
| <input type="checkbox"/> E2LB1614 | <input type="checkbox"/> E2LB1616 | <input type="checkbox"/> E2LB1714 | <input type="checkbox"/> E2LB1716 |
| <input type="checkbox"/> E2LB1814 | <input type="checkbox"/> E2LB1816 | <input type="checkbox"/> E2LB1914 | <input type="checkbox"/> E2LB1916 |
| <input type="checkbox"/> E2LB2014 | <input type="checkbox"/> E2LB2016 |                                   |                                   |

Matrix Elite E2 Lateral Back

14" - 20" W x 18"/20"H (3" Contour) HCPCS code E2620 ..... \$745

- |                                   |                                   |                                   |                                   |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> E2LB1418 | <input type="checkbox"/> E2LB1420 | <input type="checkbox"/> E2LB1518 | <input type="checkbox"/> E2LB1520 |
| <input type="checkbox"/> E2LB1618 | <input type="checkbox"/> E2LB1620 | <input type="checkbox"/> E2LB1718 | <input type="checkbox"/> E2LB1720 |
| <input type="checkbox"/> E2LB1818 | <input type="checkbox"/> E2LB1820 | <input type="checkbox"/> E2LB1918 | <input type="checkbox"/> E2LB1920 |
| <input type="checkbox"/> E2LB2018 | <input type="checkbox"/> E2LB2020 |                                   |                                   |

Matrix Elite E2 Deep Back

14" - 20" W x 14"/16"H (6" Contour) HCPCS code E2620 ..... \$755

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2D1414 | <input type="checkbox"/> E2D1416 | <input type="checkbox"/> E2D1514 | <input type="checkbox"/> E2D1516 |
| <input type="checkbox"/> E2D1614 | <input type="checkbox"/> E2D1616 | <input type="checkbox"/> E2D1714 | <input type="checkbox"/> E2D1716 |
| <input type="checkbox"/> E2D1814 | <input type="checkbox"/> E2D1816 | <input type="checkbox"/> E2D1914 | <input type="checkbox"/> E2D1916 |
| <input type="checkbox"/> E2D2014 | <input type="checkbox"/> E2D2016 |                                  |                                  |

# Ultra Low Maxx CG Power Positioning System

## TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### UPGRADABLE BACK OPTIONS

Matrix Elite E2 Deep Back  
14" - 20" W x 18"/20"H (6" Contour) HCPCS code E2620. ....\$755

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> E2D1418 | <input type="checkbox"/> E2D1420 | <input type="checkbox"/> E2D1518 | <input type="checkbox"/> E2D1520 |
| <input type="checkbox"/> E2D1618 | <input type="checkbox"/> E2D1620 | <input type="checkbox"/> E2D1718 | <input type="checkbox"/> E2D1720 |
| <input type="checkbox"/> E2D1818 | <input type="checkbox"/> E2D1820 | <input type="checkbox"/> E2D1918 | <input type="checkbox"/> E2D1920 |
| <input type="checkbox"/> E2D2018 | <input type="checkbox"/> E2D2020 |                                  |                                  |

Matrix PB Posture Back  
14" - 20" W (5" Contour) HCPCS code E2615. ....\$625

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> MPB1412 | <input type="checkbox"/> MPB1612 | <input type="checkbox"/> MPB1616 | <input type="checkbox"/> MPB1620 |
| <input type="checkbox"/> MPB1712 | <input type="checkbox"/> MPB1716 | <input type="checkbox"/> MPB1720 | <input type="checkbox"/> MPB1812 |
| <input type="checkbox"/> MPB1816 | <input type="checkbox"/> MPB1820 | <input type="checkbox"/> MPB2012 | <input type="checkbox"/> MPB2016 |
| <input type="checkbox"/> MPB2020 |                                  |                                  |                                  |

Matrix PB Heavy Duty Back<sup>(1)</sup>  
20" - 22" W (5" Contour) HCPCS code E2615 for 21" E2616 for 22" .... \$835

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> HDB2016 | <input type="checkbox"/> HDB2020 | <input type="checkbox"/> HDB2216 | <input type="checkbox"/> HDB2220 |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|

Matrix PB Posture Back Deep<sup>(1)</sup>  
14" - 20" W (7" Contour) HCPCS code E2615 .... \$645

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> PBD1412 | <input type="checkbox"/> PBD1612 | <input type="checkbox"/> PBD1616 | <input type="checkbox"/> PBD1620 |
| <input type="checkbox"/> PBD1712 | <input type="checkbox"/> PBD1716 | <input type="checkbox"/> PBD1720 | <input type="checkbox"/> PBD1812 |
| <input type="checkbox"/> PBD1816 | <input type="checkbox"/> PBD1820 | <input type="checkbox"/> PBD2012 | <input type="checkbox"/> PBD2016 |
| <input type="checkbox"/> PBD2020 |                                  |                                  |                                  |

Matrix PB Heavy Duty Deep Back<sup>(1)</sup>  
20"-22" W (5" Contour) ....\$935

- |                                   |                                   |                                   |                                   |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> DHDB2016 | <input type="checkbox"/> DHDB2020 | <input type="checkbox"/> DHDB2216 | <input type="checkbox"/> DHDB2220 |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|

Matrix MAC Back<sup>(2)</sup> HCPCS code E2620. ....\$955

- |                                  |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| <input type="checkbox"/> MAC1616 | <input type="checkbox"/> MAC1618 | <input type="checkbox"/> MAC1620 | <input type="checkbox"/> MAC1622 |
| <input type="checkbox"/> MAC1716 | <input type="checkbox"/> MAC1718 | <input type="checkbox"/> MAC1720 | <input type="checkbox"/> MAC1722 |
| <input type="checkbox"/> MAC1816 | <input type="checkbox"/> MAC1818 | <input type="checkbox"/> MAC1820 | <input type="checkbox"/> MAC1822 |
| <input type="checkbox"/> MAC1916 | <input type="checkbox"/> MAC1918 | <input type="checkbox"/> MAC1920 | <input type="checkbox"/> MAC1922 |
| <input type="checkbox"/> MAC2016 | <input type="checkbox"/> MAC2018 | <input type="checkbox"/> MAC2020 | <input type="checkbox"/> MAC2022 |

Matrix MAC Back w/ HUG<sup>(2)</sup> HCPCS code E2620. .... \$1,055

- |                                      |                                      |                                      |
|--------------------------------------|--------------------------------------|--------------------------------------|
| <input type="checkbox"/> MAC1616-HUG | <input type="checkbox"/> MAC1618-HUG | <input type="checkbox"/> MAC1620-HUG |
| <input type="checkbox"/> MAC1622-HUG | <input type="checkbox"/> MAC1716-HUG | <input type="checkbox"/> MAC1718-HUG |
| <input type="checkbox"/> MAC1720-HUG | <input type="checkbox"/> MAC1722-HUG | <input type="checkbox"/> MAC1816-HUG |
| <input type="checkbox"/> MAC1818-HUG | <input type="checkbox"/> MAC1820-HUG | <input type="checkbox"/> MAC1822-HUG |
| <input type="checkbox"/> MAC1916-HUG | <input type="checkbox"/> MAC1918-HUG | <input type="checkbox"/> MAC1920-HUG |
| <input type="checkbox"/> MAC1922-HUG | <input type="checkbox"/> MAC2016-HUG | <input type="checkbox"/> MAC2018-HUG |
| <input type="checkbox"/> MAC2020-HUG | <input type="checkbox"/> MAC2022-HUG |                                      |

Matrix MAC Lateral Contour

- |                                 |           |                                 |            |
|---------------------------------|-----------|---------------------------------|------------|
| <input type="checkbox"/> MLCSL  | Left 3.5" | <input type="checkbox"/> MLCSR  | Right 3.5" |
| <input type="checkbox"/> MLCDL  | Left 6"   | <input type="checkbox"/> MLCDR  | Right 6"   |
| <input type="checkbox"/> MLCXDL | Left 9"   | <input type="checkbox"/> MLCXDR | Right 9"   |

☐ OCBC "ON-CHAIR" Matrix Back Credit .....(\$125)

1. The Heavy Duty PB Back and PB Deep Back are not available on Recline systems with seat widths of 21" and 22"

2. Must select Pad Contour in Matrix MAC Lateral Contour section.

### OTHER UPGRADABLE BACK OPTIONS

| Cushion                         | Meshtex | 4-Way Stretch Polyester | Infection Control | Onyx  | Fusion |
|---------------------------------|---------|-------------------------|-------------------|-------|--------|
| Contour                         | N/A     | STD                     | N/A               | \$100 | N/A    |
| High                            | STD     | N/A                     | N/A               | \$100 | \$100  |
| E2, Elite, Elite HD, Elite Deep | STD     | N/A                     | \$100             | N/A   | N/A    |
| E2, Elite TR, Elite TR HD       | STD     | N/A                     | \$100             | N/A   | N/A    |

- |                                  |                                                                                  |       |
|----------------------------------|----------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> EBSS    | Super Soft Foam for Elite Back Upgrade. ....                                     | \$100 |
| <input type="checkbox"/> FFUHB   | Fusion Fabric Upgrade for High Backs. ....                                       | \$100 |
| <input type="checkbox"/> OCUCB   | Onyx Fabric Upgrade for Contour Backs. ....                                      | \$100 |
| <input type="checkbox"/> OCUHB   | Onyx Fabric Upgrade for High Backs. ....                                         | \$100 |
| <input type="checkbox"/> CHIEB   | Chest Harness Interface Kit for Elite Back <sup>(1)(2)</sup> ....                | \$95  |
| <input type="checkbox"/> CHIPB   | Chest Harness Interface Kit for PB Back <sup>(1)(2)</sup> ....                   | \$45  |
| <input type="checkbox"/> E2CHI-S | Chest Harness Interface Kit for E2 Back for 14"-15" Backs <sup>(1)(2)</sup> .... | \$85  |
| <input type="checkbox"/> E2CHI-M | Chest Harness Interface Kit for E2 Back for 16"-20" Backs <sup>(1)(2)</sup> .... | \$85  |
| <input type="checkbox"/> E2BKB   | Contoured anatomical PSIS/Sacrolumbar Support. ....                              | \$55  |

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harness

### CHEST STRAPS - ABDOMINAL SUPPORT HCPCS code E0960

- |                                  |                                                                                   |       |
|----------------------------------|-----------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> MFUS    | BP Monoflex, Underarm RLS S(3" x 15.5") ....                                      | \$195 |
| <input type="checkbox"/> MFUM    | BP Monoflex, Underarm RLS M(4" x 18.25") ....                                     | \$195 |
| <input type="checkbox"/> MFUL    | BP Monoflex, Underarm RLS L(4.5" x 21.5") ....                                    | \$195 |
| <input type="checkbox"/> MFCS    | BP Monoflex, Center RLS S(2.5" x 15.5") ....                                      | \$195 |
| <input type="checkbox"/> MFCM    | BP Monoflex, Center RLS M(3.25" x 18.25") ....                                    | \$195 |
| <input type="checkbox"/> MFCL    | BP Monoflex, Center RLS L(3.5" x 21.5") ....                                      | \$195 |
| <input type="checkbox"/> PFLEX-M | PostureFlex Abdominal Support Height 6.5", length 16" to 23" <sup>(1)</sup> ....  | \$245 |
| <input type="checkbox"/> PFLEX-L | PostureFlex Abdominal Support Height 7", length 20" to 32" <sup>(1)</sup> ....    | \$245 |
| <input type="checkbox"/> QRTAB   | Quad release tab for easy independent application and removal of PostureFlex .... | \$35  |

1. Only available on the E2 and E2 Deep Backs.

## Ultra Low Maxx CG Power Positioning System TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### TILT ARMRESTS

#### Adult Dual Post

|                              |                          |       |
|------------------------------|--------------------------|-------|
| <input type="checkbox"/> DP1 | 9.75" - 12.75"           | STD   |
| <input type="checkbox"/> DP2 | 12.75" - 15.75"          | STD   |
| <input type="checkbox"/> DP3 | 15.75" - 18.75"          | STD   |
| <input type="checkbox"/> FP  | Flat Pouch               | STD   |
| <input type="checkbox"/> GB  | Glove Boxes              | N/C   |
| <input type="checkbox"/> OAL | Outback Arm Option Left  | \$195 |
| <input type="checkbox"/> OAR | Outback Arm Option Right | \$195 |

#### Ultra Rail Seat Mount Flip Back Cantilever MAXX Tilt Arm<sup>(1)</sup>

|                                 |                                  |       |
|---------------------------------|----------------------------------|-------|
| <input type="checkbox"/> RMCAN  | 9" - 12.5"                       | \$488 |
| <input type="checkbox"/> RMCAN2 | 11" - 14"                        | \$488 |
| <input type="checkbox"/> DPM    | Dual Post Module to RMCAN (pair) | \$250 |

#### Outback Cantilever Tilt Armrest

|                                |         |       |
|--------------------------------|---------|-------|
| <input type="checkbox"/> OCTA1 | 9"-12"  | \$488 |
| <input type="checkbox"/> OCTA2 | 12"-15" | \$488 |
| <input type="checkbox"/> OCTA2 | 6"-12"  | \$488 |

#### Cantilver Arms for Endomorph Back Canes<sup>(1)</sup>

|                                |            |     |
|--------------------------------|------------|-----|
| <input type="checkbox"/> CAFEB | See page 4 | STD |
|--------------------------------|------------|-----|

1. Not available with Small Adult Sizes.

### RECLINE ARMRESTS

#### Reclining Adjustable Height

|                              |                          |       |
|------------------------------|--------------------------|-------|
| <input type="checkbox"/> RC1 | 9.5"-13"                 | STD   |
| <input type="checkbox"/> RC2 | 13"-16"                  | STD   |
| <input type="checkbox"/> OAL | Outback Arm Option Left  | \$195 |
| <input type="checkbox"/> OAR | Outback Arm Option Right | \$195 |

#### Maxx Style Cane Mounted Height, Width, Angle Adj.

#### Cantilever Flip Back Armrests<sup>(1)</sup>

|                                |           |       |
|--------------------------------|-----------|-------|
| <input type="checkbox"/> MCAN1 | 9.5"-13"  | \$488 |
| <input type="checkbox"/> MCAN2 | 12.5"-16" | \$488 |

#### Outback Cantilever Recline Armrest

|                                |         |       |
|--------------------------------|---------|-------|
| <input type="checkbox"/> OCRA1 | 9"-12"  | \$488 |
| <input type="checkbox"/> OCRA2 | 12"-15" | \$488 |
| <input type="checkbox"/> OCRA3 | 7" - 9" | \$488 |

#### Cane Mounted Height and Angle Adj.

#### Cantilver Flip Back Armrest<sup>(1)</sup>

|                                |          |       |
|--------------------------------|----------|-------|
| <input type="checkbox"/> CANQ1 | 9.5"-12" | \$488 |
| <input type="checkbox"/> CANQ2 | 10"-13"  | \$488 |
| <input type="checkbox"/> CANQ3 | 13"-16"  | \$488 |

#### Adult Dual Post On Recline

|                                |                 |       |
|--------------------------------|-----------------|-------|
| <input type="checkbox"/> DPOR1 | 9.75" - 12.75"  | \$250 |
| <input type="checkbox"/> DPOR2 | 12.75" - 15.75" | \$250 |
| <input type="checkbox"/> DPOR3 | 15.75" - 18.75" | \$250 |
| <input type="checkbox"/> FP    | Flat Pouch      | STD   |
| <input type="checkbox"/> GB    | Glove Boxes     | N/C   |

1. Not available with Small Adult Sizes.

### ARMREST TUBE LENGTH - PLEASE SPECIFY

|                               |            |     |
|-------------------------------|------------|-----|
| <input type="checkbox"/> ATFL | Left Full  | STD |
| <input type="checkbox"/> ATDL | Left Desk  | STD |
| <input type="checkbox"/> ATFR | Right Full | STD |
| <input type="checkbox"/> ATDR | Right Desk | STD |

### ARMPADS<sup>(5)</sup>

#### Standard

|                                |                          |     |
|--------------------------------|--------------------------|-----|
| <input type="checkbox"/> SAPFL | Left Full (2.25" x 14")  | STD |
| <input type="checkbox"/> SAPDL | Left Desk (2.25" x 10")  | STD |
| <input type="checkbox"/> SAPML | Left Mini (2.25" x 7")   | STD |
| <input type="checkbox"/> SAPFR | Right Full (2.25" x 14") | STD |
| <input type="checkbox"/> SAPDR | Right Desk (2.25" x 10") | STD |
| <input type="checkbox"/> SAPMR | Right Mini (2.25" x 7")  | STD |

#### Waterfall Arm Pads<sup>(1)</sup>

|                                |                  |      |
|--------------------------------|------------------|------|
| <input type="checkbox"/> WAPFL | Left Full (13")  | \$50 |
| <input type="checkbox"/> WAPDL | Left Desk (9")   | \$50 |
| <input type="checkbox"/> WAPFR | Right Full (13") | \$50 |
| <input type="checkbox"/> WAPDR | Right Desk (9")  | \$50 |

#### Modular Arm Pad - Left<sup>(2)</sup>

|                                 |                      |      |
|---------------------------------|----------------------|------|
| <input type="checkbox"/> MAPTFL | Left Tray Full (14") | \$55 |
| <input type="checkbox"/> MAPTDL | Left Tray Desk (10") | \$55 |

Select Pad Insert below:

|                                |                                   |
|--------------------------------|-----------------------------------|
| <input type="checkbox"/> ISSFL | Iskin - Integrated soft skin foam |
| <input type="checkbox"/> SCVFL | Startex Covered Visco foam        |

#### Modular Arm Pad - Right<sup>(2)</sup>

|                                 |                       |      |
|---------------------------------|-----------------------|------|
| <input type="checkbox"/> MAPTFR | Right Tray Full (14") | \$55 |
| <input type="checkbox"/> MAPTDR | Right Tray Desk (10") | \$55 |

Select Pad Insert below:

|                                |                                   |
|--------------------------------|-----------------------------------|
| <input type="checkbox"/> ISSFR | Iskin - Integrated soft skin foam |
| <input type="checkbox"/> SCVFR | Startex Covered Visco foam        |

#### Flat Multi Position Arm Pads

|                                |                  |      |
|--------------------------------|------------------|------|
| <input type="checkbox"/> MPAFL | Left Full (14")  | \$55 |
| <input type="checkbox"/> MPADL | Left Desk (10")  | \$55 |
| <input type="checkbox"/> MPAML | Left Mini (7")   | \$55 |
| <input type="checkbox"/> MPAFR | Right Full (14") | \$55 |
| <input type="checkbox"/> MPADR | Right Desk (10") | \$55 |
| <input type="checkbox"/> MPAMR | Right Mini (7")  | \$55 |

#### Gel Pads

|                                |                              |      |
|--------------------------------|------------------------------|------|
| <input type="checkbox"/> GAPNL | Narrow Left Desk (2" x 12")  | \$75 |
| <input type="checkbox"/> GAPNR | Narrow Right Desk (2" x 12") | \$75 |
| <input type="checkbox"/> GAPWL | Wide Left Desk (3.5" x 12")  | \$75 |
| <input type="checkbox"/> GAPWR | Wide Right Desk (3.5" x 12") | \$75 |
| <input type="checkbox"/> GAPLL | Long Left Full (3.5" x 14")  | \$90 |
| <input type="checkbox"/> GAPLR | Long Right Full (3.5" x 14") | \$90 |

#### Arm Troughs

|                               |                                                       |      |
|-------------------------------|-------------------------------------------------------|------|
| <input type="checkbox"/> EATL | Ergonomic Arm Trough - Left<br>HCPCS code E2209       | \$85 |
| <input type="checkbox"/> EATR | Ergonomic Arm Trough - Right<br>HCPCS code E2209      | \$85 |
| <input type="checkbox"/> FHPL | Flat Hand Pad - Left <sup>(3)</sup> HCPCS code K0108  | \$55 |
| <input type="checkbox"/> FHPR | Flat Hand Pad - Right <sup>(3)</sup> HCPCS code K0108 | \$55 |

#### Arm Pad Storage/Cup Holder<sup>(4)</sup>

|                               |       |       |
|-------------------------------|-------|-------|
| <input type="checkbox"/> APSL | Left  | \$295 |
| <input type="checkbox"/> APSR | Right | \$295 |

☐ OAP Omit Armpads(7/8" diameter tube) .....(\$25)

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.

2. Requires interchangeable pad insert.

3. Must order with ergonomic arm trough.

4. Not available on the joystick side, or with Outback Cantilever arms.

5. Full length arm pads used on the same side as the Joystick result in a "Back to Joystick" measurement of 20" minimum (24" including switch mounting)

## Ultra Low Maxx CG Power Positioning System TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### ARM ADDUCTOR PAD<sup>(1)(2)</sup>

|                                |                            |       |
|--------------------------------|----------------------------|-------|
| <input type="checkbox"/> AADSL | Left Small (6" x 2.25")    | \$100 |
| <input type="checkbox"/> AADSR | Right Small (6" x 2.25")   | \$100 |
| <input type="checkbox"/> AADLL | Left Large (8.5" x 2.25")  | \$100 |
| <input type="checkbox"/> AADLR | Right Large (8.5" x 2.25") | \$100 |

1. Only available with Modular or Flat Multi Position Arm Pads.
2. 1" thick pad, unlimited fore/aft adjustment.

### ELBOW BLOCKS

#### Left Elbow Block<sup>(1)</sup>

Select Pad Size and Both Extensions:

|                               |                         |       |
|-------------------------------|-------------------------|-------|
| <input type="checkbox"/> EBXL | Size XS (3.25" x 4.25") | \$195 |
| <input type="checkbox"/> EBSL | Size S (3.5" x 5.25")   | \$195 |
| <input type="checkbox"/> EBML | Size M (3.75" x 5.75")  | \$195 |
| <input type="checkbox"/> EBL  | Size L (4.25" x 6.5")   | \$195 |

|                              |                              |
|------------------------------|------------------------------|
| <input type="checkbox"/> LV3 | Left Vertical Extension 3"   |
| <input type="checkbox"/> LV5 | Left Vertical Extension 5"   |
| <input type="checkbox"/> LH5 | Left Horizontal Extension 5" |
| <input type="checkbox"/> LH7 | Left Horizontal Extension 7" |

#### Right Elbow Block<sup>(1)</sup>

Select Pad Size and Both Extensions:

|                               |                         |       |
|-------------------------------|-------------------------|-------|
| <input type="checkbox"/> EBXR | Size XS (3.25" x 4.25") | \$195 |
| <input type="checkbox"/> EBSR | Size S (3.5" x 5.25")   | \$195 |
| <input type="checkbox"/> EBMR | Size M (3.75" x 5.75")  | \$195 |
| <input type="checkbox"/> EBLR | Size L (4.25" x 6.5")   | \$195 |

|                              |                               |
|------------------------------|-------------------------------|
| <input type="checkbox"/> RV3 | Right Vertical Extension 3"   |
| <input type="checkbox"/> RV5 | Right Vertical Extension 5"   |
| <input type="checkbox"/> RH5 | Right Horizontal Extension 5" |
| <input type="checkbox"/> RH7 | Right Horizontal Extension 7" |

1. Prevents arm from slipping.

### ADDITIONAL ARM ACCESSORIES

|                                 |                                                                        |       |
|---------------------------------|------------------------------------------------------------------------|-------|
| <input type="checkbox"/> MACESL | Multi Axis Upper Extremity Support <sup>(1)</sup> (each)<br>Left Side  | \$225 |
| <input type="checkbox"/> MACESR | Multi Axis Upper Extremity Support <sup>(1)</sup> (each)<br>Right Side | \$225 |
| <input type="checkbox"/> OAS2L  | 1" Offset Armrest Spacer Left(each)                                    | \$125 |
| <input type="checkbox"/> OAS2R  | 1" Offset Armrest Spacer Right(each)                                   | \$125 |
| <input type="checkbox"/> OAS1L  | 1/2" Offset Armrest Spacer Left(each)                                  | \$125 |
| <input type="checkbox"/> OAS1R  | 1/2" Offset Armrest Spacer Right(each)                                 | \$125 |
| <input type="checkbox"/> MPAPA  | Adapter Block <sup>(2)</sup> (each) Quantity ( )                       | \$35  |

1. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
2. Not available with Standard Arm Pads or Waterfall Arm Pads

### MAXX LATERAL TRUNK SUPPORTS<sup>(1)</sup>

#### Lateral Trunk Support Left HCPCS code E0956

|                                |                         |       |
|--------------------------------|-------------------------|-------|
| <input type="checkbox"/> LATXL | Size XS (3.25" x 4.25") | \$100 |
| <input type="checkbox"/> LATSL | Size S (3.5" x 5.25")   | \$100 |
| <input type="checkbox"/> LATML | Size M (3.75" x 5.75")  | \$100 |
| <input type="checkbox"/> LATLL | Size L (4.25" x 6.5")   | \$100 |

Add:

|                                 |                                                                          |       |
|---------------------------------|--------------------------------------------------------------------------|-------|
| <input type="checkbox"/> SALATL | MAXX Style Swing Away Hardware - Left <sup>(2)</sup><br>HCPCS code E1028 | \$195 |
| <input type="checkbox"/> ATLL   | Additional Telescoping Link - Left <sup>(3)</sup>                        | \$75  |

#### Lateral Trunk Support Right HCPCS code E0956

|                                |                         |       |
|--------------------------------|-------------------------|-------|
| <input type="checkbox"/> LATXR | Size XS (3.25" x 4.25") | \$100 |
| <input type="checkbox"/> LATSR | Size S (3.5" x 5.25")   | \$100 |
| <input type="checkbox"/> LATMR | Size M (3.75" x 5.75")  | \$100 |
| <input type="checkbox"/> LATLR | Size L (4.25" x 6.5")   | \$100 |

Add:

|                                 |                                                                           |       |
|---------------------------------|---------------------------------------------------------------------------|-------|
| <input type="checkbox"/> SALATR | MAXX Style Swing Away Hardware - Right <sup>(2)</sup><br>HCPCS code E1028 | \$195 |
| <input type="checkbox"/> ATLR   | Additional Telescoping Link - Right <sup>(3)</sup>                        | \$75  |

1. Comes standard with Fixed Mount Hardware. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Comes standard with 1" telescoping link which gives 0-2.75" medial offset.
3. The amount of offset available may be impacted by back thickness.

### MATRIX FIXED LATERAL TRUNK SUPPORTS<sup>(1)</sup>

HCPCS code E1028

|                                |                                                       |       |
|--------------------------------|-------------------------------------------------------|-------|
| <input type="checkbox"/> E2FLL | Matrix E2 Fixed Laterals-Left(6.5x5") <sup>(2)</sup>  | \$145 |
| <input type="checkbox"/> E2FLR | Matrix E2 Fixed Laterals-Right(6.5x5") <sup>(2)</sup> | \$145 |
| <input type="checkbox"/> PBSSL | Matrix Small Lateral Fixed-Left (3x5.5")              | \$100 |
| <input type="checkbox"/> PBSSR | Matrix Small Lateral Fixed-Right (3x5.5")             | \$100 |
| <input type="checkbox"/> PBLSL | Matrix Standard Lateral Fixed-Left (7.5x5")           | \$100 |
| <input type="checkbox"/> PBLSR | Matrix Standard Lateral Fixed-Right (7.5x5")          | \$100 |
| <input type="checkbox"/> PBOLL | Matrix Offset Lateral Fixed Left (7.5x5")             | \$110 |
| <input type="checkbox"/> PBOLR | Matrix Offset Lateral Fixed Right (7.5x5")            | \$110 |

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Only available with E2 backs.

### MATRIX SWING AWAY LATERAL TRUNK SUPPORTS<sup>(1)</sup>

#### Matrix Lateral Trunk Support Left HCPCS code E0956

|                                 |                         |       |
|---------------------------------|-------------------------|-------|
| <input type="checkbox"/> EPADXL | Size XS (3.25" x 4.25") | \$100 |
| <input type="checkbox"/> EPADSL | Size S (3.5" x 5.25")   | \$100 |
| <input type="checkbox"/> EPADML | Size M (3.75" x 5.75")  | \$100 |
| <input type="checkbox"/> EPADLL | Size L (4.25" x 6.5")   | \$100 |

Must select one option:

|                                 |                                                                  |       |
|---------------------------------|------------------------------------------------------------------|-------|
| <input type="checkbox"/> SAHEBL | Matrix Elite Swing Away Hardware Left<br>HCPCS code E1028        | \$195 |
| <input type="checkbox"/> OSAHL  | Matrix Offset Elite Swing Away Hardware Left<br>HCPCS code E1028 | \$195 |

#### Matrix Lateral Trunk Support Right HCPCS code E0956

|                                 |                         |       |
|---------------------------------|-------------------------|-------|
| <input type="checkbox"/> EPADXR | Size XS (3.25" x 4.25") | \$100 |
| <input type="checkbox"/> EPADSR | Size S (3.5" x 5.25")   | \$100 |
| <input type="checkbox"/> EPADMR | Size M (3.75" x 5.75")  | \$100 |
| <input type="checkbox"/> EPADLR | Size L (4.25" x 6.5")   | \$100 |

Must select one option:

|                                 |                                                                   |       |
|---------------------------------|-------------------------------------------------------------------|-------|
| <input type="checkbox"/> SAHEBR | Matrix Elite Swing Away Hardware Right<br>HCPCS code E1028        | \$195 |
| <input type="checkbox"/> OSAHR  | Matrix Offset Elite Swing Away Hardware Right<br>HCPCS code E1028 | \$195 |

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

## Ultra Low Maxx CG Power Positioning System TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### HIP AND KNEE SUPPORTS

#### Lateral Hip Support Left HCPCS code E0953

|                                  |                               |       |
|----------------------------------|-------------------------------|-------|
| <input type="checkbox"/> HSPXLH  | Size XS (5.5" W x 4" H) ..... | \$120 |
| <input type="checkbox"/> HSPSMLH | Size S/M (7" W x 4" H) .....  | \$120 |
| <input type="checkbox"/> HSPSLH  | Size S (7" W x 6" H) .....    | \$120 |
| <input type="checkbox"/> HSPMLH  | Size M (9" W x 4" H) .....    | \$120 |
| <input type="checkbox"/> HSPLLH  | Size L (13" W x 4" H) .....   | \$120 |

#### Lateral Hip Support Right<sup>(1)</sup> HCPCS code E0953

|                                  |                               |       |
|----------------------------------|-------------------------------|-------|
| <input type="checkbox"/> HSPXRH  | Size XS (5.5" W x 4" H) ..... | \$120 |
| <input type="checkbox"/> HSPSMRH | Size S/M (7" W x 4" H) .....  | \$120 |
| <input type="checkbox"/> HSPSRH  | Size S (7" W x 6" H) .....    | \$120 |
| <input type="checkbox"/> HSPMRH  | Size M (9" W x 4" H) .....    | \$120 |
| <input type="checkbox"/> HSPLRH  | Size L (13" W x 4" H) .....   | \$120 |

#### Lateral Knee Support Left<sup>(1)</sup> HCPCS code E0953

|                                  |                               |       |
|----------------------------------|-------------------------------|-------|
| <input type="checkbox"/> HSPXLK  | Size XS (5.5" W x 4" H) ..... | \$120 |
| <input type="checkbox"/> HSPSMLK | Size S/M (7" W x 4" H) .....  | \$120 |
| <input type="checkbox"/> HSPSLK  | Size S (7" W x 6" H) .....    | \$120 |
| <input type="checkbox"/> HSPMLK  | Size M (9" W x 4" H) .....    | \$120 |
| <input type="checkbox"/> HSPLLK  | Size L (13" W x 4" H) .....   | \$120 |

#### Lateral Knee Support Right<sup>(1)</sup> HCPCS code E0953

|                                  |                               |       |
|----------------------------------|-------------------------------|-------|
| <input type="checkbox"/> HSPXRK  | Size XS (5.5" W x 4" H) ..... | \$120 |
| <input type="checkbox"/> HSPSMRK | Size S/M (7" W x 4" H) .....  | \$120 |
| <input type="checkbox"/> HSPSRK  | Size S (7" W x 6" H) .....    | \$120 |
| <input type="checkbox"/> HSPMRK  | Size M (9" W x 4" H) .....    | \$120 |
| <input type="checkbox"/> HSPLRK  | Size L (13" W x 4" H) .....   | \$120 |

#### Add Removable Hardware:

##### MAXX Style Quick Release, Multi-Axis Mounting

##### Hardware for Hip Supports HCPCS code E1028

|                                 |                             |       |
|---------------------------------|-----------------------------|-------|
| <input type="checkbox"/> MHSHLH | Left mounted at Hip .....   | \$215 |
| <input type="checkbox"/> MHSHLK | Left mounted at Knee .....  | \$215 |
| <input type="checkbox"/> MHSHRH | Right mounted at Hip .....  | \$215 |
| <input type="checkbox"/> MHSHRK | Right mounted at Knee ..... | \$215 |

##### Swing Away Hardware for Hip Supports HCPCS code E1028

|                                 |                             |       |
|---------------------------------|-----------------------------|-------|
| <input type="checkbox"/> SHSHLH | Left mounted at Hip .....   | \$215 |
| <input type="checkbox"/> SHSHLK | Left mounted at Knee .....  | \$215 |
| <input type="checkbox"/> SHSHRH | Right mounted at Hip .....  | \$215 |
| <input type="checkbox"/> SHSHRK | Right mounted at Knee ..... | \$215 |

##### Extended Maxx Style Quick Release, Multi-Axis Mounting

##### Hardware for Hip Supports<sup>(2)</sup> HCPCS code E1028

|                                 |                             |       |
|---------------------------------|-----------------------------|-------|
| <input type="checkbox"/> EMQHLH | Left mounted at Hip .....   | \$215 |
| <input type="checkbox"/> EMQHLK | Left mounted at Knee .....  | \$215 |
| <input type="checkbox"/> EMQHRH | Right mounted at Hip .....  | \$215 |
| <input type="checkbox"/> EMQHRK | Right mounted at Knee ..... | \$215 |

##### Lift Off Hardware for Hip Support

##### (Hip Position Only) HCPCS code E1028

|                               |             |       |
|-------------------------------|-------------|-------|
| <input type="checkbox"/> LOHL | Left .....  | \$195 |
| <input type="checkbox"/> LOHR | Right ..... | \$195 |

1. Comes standard with Fixed Mount Hardware.

2. Recommended when finished seat cushion thickness will be greater than 3.5" high.

### POWER CENTER MOUNT FOOT PLATFORM

| MUST BE COMPLETED              | Measurement in Inches |
|--------------------------------|-----------------------|
| Knee to Heel:                  |                       |
| Cushion Thickness: (Subtract)  |                       |
| Cushion Compression: (Add)     |                       |
| Seat Pan to Footplate: (Total) |                       |

#### LNK Power Center Mount Foot Platform HCPCS code E1012

|                                             |                                            |         |
|---------------------------------------------|--------------------------------------------|---------|
| <input type="checkbox"/> LNX                | 9" Articulation <sup>(1)(2)</sup> .....    | \$2,995 |
| <input type="checkbox"/> LSPTF1             | 10" - 14" Seat pan to footplate            |         |
| <input type="checkbox"/> LSPTF2             | 13" - 19" Seat pan to footplate            |         |
| <input type="checkbox"/> LNXA               | 6" Articulation <sup>(1)(2)(6)</sup> ..... | \$2,995 |
| <input checked="" type="checkbox"/> LASPTF2 | 10" - 14" Seat pan to footplate            |         |

#### Options:

|                                |                                                                |       |
|--------------------------------|----------------------------------------------------------------|-------|
| <input type="checkbox"/> LIFP  | Rubber Coated Footplate 11.5" W x 10" D .....                  | STD   |
| <input type="checkbox"/> LICP  | Individual calf pads .....                                     | STD   |
| <input type="checkbox"/> OPCP  | One Piece Malleable Calf Panel <sup>(3)</sup> .....            | \$250 |
| <input type="checkbox"/> LNXPF | Plantar Flex Option for LNK Foot platform <sup>(4)</sup> ..... | \$250 |

#### Latitude Style Foot Platform

|                                |                                                    |         |
|--------------------------------|----------------------------------------------------|---------|
| <input type="checkbox"/> LSFP  | Latitude Style Foot Platform <sup>(5)</sup> .....  | \$2,995 |
|                                | For finished knee to heel measurements of 6.5"-12" |         |
| <input type="checkbox"/> LSPFP | Latitude Style Padded Foot Plate .....             | \$250   |

1. 3" ground clearance required.

2. Only available on adult range seat systems 16" deep or greater.

For shorter seat depths call Customer Service for other options.

3. Not available with Plantar Flex Option for LNK Foot platform, or with LNXA.

4. Only available with Foot Platform Options.

5. Not available with TRRO and TRRO-E.

6. For clients under 4' only.

### MANUAL ELEVATING CENTER MOUNT LEGREST

|                                |                                                                  |       |
|--------------------------------|------------------------------------------------------------------|-------|
| <input type="checkbox"/> MECML | Manual Elevating Center Mount Legrest<br>With Plantar Flex ..... | \$595 |
|--------------------------------|------------------------------------------------------------------|-------|

#### Select legrest set up position:

|                                |                                                               |     |
|--------------------------------|---------------------------------------------------------------|-----|
| <input type="checkbox"/> 90SU  | Set at 90° (elevating range 90-40°) .....                     | STD |
| <input type="checkbox"/> 110TU | Tucked at 110° (elevating range 110-60°) <sup>(1)</sup> ..... | N/C |

#### Rubber Coated Footplate:

|                               |                                               |       |
|-------------------------------|-----------------------------------------------|-------|
| <input type="checkbox"/> MRCP | Rubber coated footplate 11.5" W x 10" D ..... | \$165 |
|-------------------------------|-----------------------------------------------|-------|

#### Select seat pan to footplate for Rubber Coated Footplate

|                                 |                   |     |
|---------------------------------|-------------------|-----|
| <input type="checkbox"/> CSPTF1 | 11.5" - 15" ..... | N/C |
| <input type="checkbox"/> CSPTF2 | 14" - 19" .....   | N/C |

#### Individual Footplates:

|                               |                                           |       |
|-------------------------------|-------------------------------------------|-------|
| <input type="checkbox"/> IFSL | Small - (5" x 7.5" deep) - Left .....     | \$165 |
| <input type="checkbox"/> IFSR | Small - (5" x 7.5" deep) - Right .....    | \$165 |
| <input type="checkbox"/> IFML | Medium - (5.5" x 9.5" deep) - Left .....  | \$165 |
| <input type="checkbox"/> IFMR | Medium - (5.5" x 9.5" deep) - Right ..... | \$165 |
| <input type="checkbox"/> IFLL | Large - (6" x 11.5" deep) - Left .....    | \$165 |
| <input type="checkbox"/> IFLR | Large - (6" x 11.5" deep) - Right .....   | \$165 |


#### Select seat pan to footplate for Individual Footplates

|                                 |                    |     |
|---------------------------------|--------------------|-----|
| <input type="checkbox"/> CSPTF3 | 9.5" - 14.5" ..... | N/C |
| <input type="checkbox"/> CSPTF4 | 12" - 19" .....    | N/C |

#### Optional:

|                               |                            |       |
|-------------------------------|----------------------------|-------|
| <input type="checkbox"/> CICP | Individual calf pads ..... | \$200 |
|-------------------------------|----------------------------|-------|

1. Not available on seat depths 16", 17".



# Ultra Low Maxx CG Power Positioning System TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

## STANDARD FRONT RIGGING


- ☐ FCMP Fixed Center Mount Platform<sup>(1)</sup> ..... N/C  
☐ FIFP Comes with a rubber coated footplate  
 11.5" W x 10" D ..... N/C

Select seat pan to footplate:

- ☐ FSPTF1 9" - 13" ..... N/C  
☐ FSPTF14 13" -19" ..... N/C

Set at:

- ☐ FSA70 Set at 70° ..... N/C  
☐ FSA83 Set at 83° ..... N/C  
☐ FSA90 Set at 90° ..... N/C  
☐ FSA97 Set at 97° ..... N/C



- ☐ I70SAF Invacare 70° Swingaway Footrests (pair) ..... N/C  
☐ I70-CF Comes with composite foot plates. .... N/C

1. 3" ground clearance required

## MANUAL CENTER MOUNT FRONT RIGGING

- ☐ MSFCMFP Maxx Style Fixed Center Mount  
 Foot Platform<sup>(1)(2)</sup> ..... \$395


- ☐ MIFP Comes with a rubber coated footplate  
 11.5" W x 10" D ..... N/C  
☐ MICP Comes standard with individual calf pads ..... N/C

Select seat pan to footplate:

- ☐ MSPTF1 9" - 13" ..... N/C  
☐ MSPTF14 13" -19" ..... N/C

Set at:


- ☐ MSA70 Set at 70° ..... N/C  
☐ MSA83 Set at 83° ..... N/C  
☐ MSA90 Set at 90° ..... N/C  
☐ MSA97 Set at 97° ..... N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.

2. 3" ground clearance required.





## FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP



- ☐ SFFP Small Foot Platform (12"w x 8.5"d) ..... \$180  
☐ MFFP Medium Platform (12"w x 10.75"d) ..... \$180  
☐ LFFP Large Platform (14.5"w x 13"d) ..... \$250  
☐ XLFFP Extra Large Platform (17"w x 13"d) ..... \$250  
☐ XXLFFP Extra Extra Large Platform (21"w x 13"d) ..... \$250

1. May need to adjust leg position to avoid interference with front casters

## INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP



- ☐ IFSL Small (5" W x 7.5" D) - Left ..... \$165  
☐ IFSR Small (5" W x 7.5" D) - Right ..... \$165  
☐ IFML Medium (5.5" W x 9.5" D) - Left ..... \$165  
☐ IFMR Medium (5.5" W x 9.5" D) - Right ..... \$165  
☐ IFLL Large (6" W x 11.5" D) - Left ..... \$165  
☐ IFLR Large (6" W x 11.5" D) - Right ..... \$165

Add:

- ☐ TSSL Left Toe Strap Small 6.5" ..... \$125  
☐ TSSR Right Toe Strap Small 6.5" ..... \$125  
☐ TSML Left Toe Strap Medium 8.5" ..... \$125  
☐ TSMR Right Toe Strap Medium 8.5" ..... \$125  
☐ TSLL Left Toe Strap Large 10.5" ..... \$125  
☐ TSLR Right Toe Strap Large 10.5" ..... \$125

- ☐ MFOL Multi-Axis Footplate (6"W x 10"D) - Left ..... \$225  
☐ MFOR Multi-Axis Footplate (6"W x 10"D) - Right ..... \$225

1. Foot plates may interfere with casters in certain configurations

## MAXX INDEPENDENT ADJUSTABLE POSITIONING



- 5.5" - 6" WIDTH ADJUSTMENT
- 9.5" - 11.5" DEPTH ADJUSTMENT
- UP TO 9° ADDUCTION
- UP TO 45° ABDUCTION

- ☐ MPFPL Maxx Independent Adjustable Positioning  
 Foot plate<sup>(1)</sup> - Left ..... \$200  
☐ MPFPR Maxx Independent Adjustable Positioning  
 Foot plate<sup>(1)</sup> - Right ..... \$200  
☐ MTSL Left Toe Strap (each) HCPCS code E0952 ..... \$35  
☐ MTSR Right Toe Strap (each) HCPCS code E0952 ..... \$35

1. Foot plates may interfere with casters in certain configurations.

## Ultra Low Maxx CG Power Positioning System TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### INVACARE CENTER MOUNT<sup>(1)</sup>

#### Invacare Center Mount with Footboard<sup>(3)</sup>

Select seat pan to footplate:

- ☐ CTMT-FJ-FB Invacare Center Mount 9.25" - 12" ..... \$760  
☐ CTMT-F-FB Invacare Center Mount 12.25" -15" ..... \$760

#### Invacare Center Mount with Individual Footplates

Select seat pan to footplate:

- ☐ CTMT-FJ Invacare Center Mount 9.25" - 12" ..... \$595  
☐ CTMT-F Invacare Center Mount 12.25" -15" ..... \$595

Select Footplates:

- ☐ SM Individual Foot Plates (5.25"W x 7.5"D)..... \$165  
☐ MDM Individual Foot Plates (5.25"W x 10"D)..... \$165  
☐ LRG Individual Foot Plates (5.9"W x 12"D)..... \$165

Optional:

- ☐ CALF Invacare Calf Pads for use with CTMT<sup>(2)</sup> ..... \$268

1. 3" ground clearance required.  
2. Not available with CTMT-FJ, CTMT-FJ-FB  
3. Footboard Included

### MAXX STYLE POWER ELEVATING LEGRESTS

- ☐ MSPPL MAXX Swing Away  
Power Elevating<sup>(1)</sup> HCPCS code E1010 ..... \$2,490

Select Power Legrest Function:

- ☐ INDA Individual Movement ..... N/C  
☐ COMA Combined Movement ..... N/C  
☐ SYNA Synchronized/Recline Movement..... N/C

Select Receivers:

- ☐ MRS Straight Receivers ..... N/C  
☐ MRF Flared 1" Outward Receivers ..... N/C  
☐ MRF2 Flared 2" Outward Receivers..... \$250

Select Seat Pan to Footplate:

- ☐ PSPTF5 Length: 10"-14" ..... N/C  
☐ PSPTF2 Length: 13"-17"..... N/C  
☐ PSPTF7 Length: 16"-20"..... N/C

Optional:

- ☐ PMCP Moulded I-Skin Calf Pads ..... STD  
☐ PFCPMS Flat Calf Pad Option (pair) ..... \$275  
☐ PCCPMS Curved Calf Pad Option<sup>(2)</sup> (pair) ..... \$275

1. Must choose footplate on page 15.  
2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

### MAXX STYLE MANUAL ELEVATING LEGRESTS

- ☐ MMPPSA MAXX Swing Away  
Manual Elevating<sup>(1)</sup> HCPCS code E0990 ..... \$680

Select Receivers:

- ☐ ERS Straight Receivers ..... N/C  
☐ ERF Flared 1" Outward Receivers ..... N/C  
☐ ERF2 Flared 2" Outward Receivers..... \$250

Select Seat Pan to Footplate:

- ☐ ESPTF5 Length: 10"-14" ..... N/C  
☐ ESPTF2 Length: 13"-17"..... N/C  
☐ ESPTF7 Length: 16"-20"..... N/C

Optional:

- ☐ EMCP Moulded I-Skin Calf Pads ..... STD  
☐ EFCPMS Flat Calf Pad Option (pair) ..... \$275  
☐ ECCPMS Curved Calf Pad Option<sup>(2)</sup> (pair) ..... \$275

1. Must choose footplate on page 15.  
2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

### SWING AWAY LEGREST

- ☐ HD70 Heavy Duty 70° Swing Away Footrest(pair)<sup>(1)(2)</sup> . \$495

Select Receivers:

- ☐ RS Straight Receivers ..... N/C  
☐ RF Flared 1" Outward Receivers ..... N/C  
☐ RF2 Flared 2" Outward Receivers..... \$250

Select Seat Pan to Footplate:

- ☐ HSPTF11 Length: 6"-10" ..... N/C  
☐ HSPTF8 Length: 10"-14" ..... N/C  
☐ HSPTF2 Length: 12"-16" ..... N/C  
☐ HSPTF7 Length: 15"-19"..... N/C

Options:

- ☐ HD70FCP Flat Calf Pad Option (pair) ..... \$275  
☐ HD70CCP Curved Calf Pad Option<sup>(3)</sup> (pair) ..... \$275

- ☐ I70SAFT Invacare 70° Swingaway Footrests..... \$325  
☐ I70-CFT Comes with composite foot plates..... N/C

- ☐ RLRL Residual Limb Support - Left (8.5" x 6")<sup>(4)</sup>  
HCPCS code E1020 ..... \$495  
☐ RLRS Residual Limb Support - Right (8.5" x 6")<sup>(4)</sup>  
HCPCS code E1020 ..... \$495

1. Must choose footplate on page 15.  
2. Heavy duty durable construction (4.5lbs ea.).  
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.  
4. Swing away and removable, 50° of angle adjustment, 4.5" height adjustment, 3" horizontal and 2" of for/aft adjustment.

### DUAL FRONT RIGGING ADAPTOR<sup>(1)</sup>

- ☐ DFRA Dual Front Rigging Receivers..... \$450

1. Must purchase each set of front rigging separately. Not available on no charge front rigging options.

# Ultra Low Maxx CG Power Positioning System

## TDX-SP2

Dealer Name: \_\_\_\_\_





Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### FOOTPLATES - SWING AWAY/ELEVATING LEGRESTS<sup>(1)</sup>

HCPCS code K0040






#### Angle Adjustable Footplates

|                               |                                                |       |
|-------------------------------|------------------------------------------------|-------|
| <input type="checkbox"/> AFSL | Left Small Footplate (4.25"W x 8"D) .....      | \$165 |
| <input type="checkbox"/> AFSR | Right Small Footplate (4.25"W x 8"D) .....     | \$165 |
| <input type="checkbox"/> AFML | Left Medium Footplate (5.25"W x 8"D) .....     | \$165 |
| <input type="checkbox"/> AFMR | Right Medium Footplate (5.25"W x 8"D) .....    | \$165 |
| <input type="checkbox"/> AFLL | Left Large Footplate (6.25"W x 8"D) .....      | \$165 |
| <input type="checkbox"/> AFLR | Right Large Footplate (6.25"W x 8"D) .....     | \$165 |
| <input type="checkbox"/> AFXL | Left Extra Large Footplate (9"W x 11"D) .....  | \$165 |
| <input type="checkbox"/> AFXR | Right Extra Large Footplate (9"W x 11"D) ..... | \$165 |

#### Multi-Axis Adjustable Angle Footplates

|                               |                                                |       |
|-------------------------------|------------------------------------------------|-------|
| <input type="checkbox"/> MFSL | Left Small Footplate (4.25"W x 8"D) .....      | \$225 |
| <input type="checkbox"/> MFSR | Right Small Footplate (4.25"W x 8"D) .....     | \$225 |
| <input type="checkbox"/> MFML | Left Medium Footplate (5.25"W x 8"D) .....     | \$225 |
| <input type="checkbox"/> MFMR | Right Medium Footplate (5.25"W x 8"D) .....    | \$225 |
| <input type="checkbox"/> MFL  | Left Large Footplate (6.25"W x 8"D) .....      | \$225 |
| <input type="checkbox"/> MFLR | Right Large Footplate (6.25"W x 8"D) .....     | \$225 |
| <input type="checkbox"/> MFXL | Left Extra Large Footplate (9"W x 11"D) .....  | \$225 |
| <input type="checkbox"/> MFXR | Right Extra Large Footplate (9"W x 11"D) ..... | \$225 |



#### Cut Out Angle Adjustable Footplates

|                                |                                             |       |
|--------------------------------|---------------------------------------------|-------|
| <input type="checkbox"/> AFCSL | Left Small Footplate (7.25"W x 9"D) .....   | \$165 |
| <input type="checkbox"/> AFCSR | Right Small Footplate (7.25"W x 9"D) .....  | \$165 |
| <input type="checkbox"/> AFCML | Left Medium Footplate (7.5"W x 11"D) .....  | \$165 |
| <input type="checkbox"/> AFCMR | Right Medium Footplate (7.5"W x 11"D) ..... | \$165 |
| <input type="checkbox"/> AFCLL | Left Large Footplate (10"W x 14"D) .....    | \$165 |
| <input type="checkbox"/> AFCRL | Right Large Footplate (10"W x 14"D) .....   | \$165 |

#### Cut Out Multi-Axis Adjustable Angle Footplates

|                                |                                             |       |
|--------------------------------|---------------------------------------------|-------|
| <input type="checkbox"/> MFCSL | Left Small Footplate (7.25"W x 9"D) .....   | \$225 |
| <input type="checkbox"/> MFCSR | Right Small Footplate (7.25"W x 9"D) .....  | \$225 |
| <input type="checkbox"/> MFCML | Left Medium Footplate (7.5"W x 11"D) .....  | \$225 |
| <input type="checkbox"/> MFCMR | Right Medium Footplate (7.5"W x 11"D) ..... | \$225 |
| <input type="checkbox"/> MFCLL | Left Large Footplate (10"W x 14"D) .....    | \$225 |
| <input type="checkbox"/> MFCLR | Right Large Footplate (10"W x 14"D) .....   | \$225 |

1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

### HEEL LOOPS, CUPS AND TOE STRAPS

|                              |                                               |      |
|------------------------------|-----------------------------------------------|------|
| <input type="checkbox"/> HLL | Left Heel Loop (each) HCPCS code E0951 .....  | \$20 |
| <input type="checkbox"/> HLR | Right Heel Loop (each) HCPCS code E0951 ..... | \$20 |
| <input type="checkbox"/> HCL | Left Cup Loop (each) HCPCS code E0951 .....   | \$20 |
| <input type="checkbox"/> HCR | Right Cup Loop (each) HCPCS code E0951 .....  | \$20 |
| <input type="checkbox"/> TSL | Left Toe Strap (each) HCPCS code E0952 .....  | \$20 |
| <input type="checkbox"/> TSR | Right Toe Strap (each) HCPCS code E0952 ..... | \$20 |

### CALF STRAPS AND CALF PANELS<sup>(1)</sup>

|                              |                                                                                                       |       |
|------------------------------|-------------------------------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> CSM | Body Point Airmesh Calf Strap for Chairs<br>15"-18" wide with 3.5" H Pad <sup>(1)</sup> .....         | \$150 |
| <input type="checkbox"/> CSL | Body Point Airmesh Calf Strap for Chairs<br>18"-22" wide with 3.5" H Pad <sup>(1)</sup> .....         | \$150 |
| <input type="checkbox"/> CPM | Body Point Airmesh Calf Panel for Chairs<br>15"-18" wide with 35.5" W x 9" H Pad <sup>(1)</sup> ..... | \$220 |
| <input type="checkbox"/> CPL | Body Point Airmesh Calf Panel for Chairs<br>18"-22" wide with 38.5" W x 9" H Pad <sup>(1)</sup> ..... | \$235 |

1. Not available for Center Mounted Front Rigging.

### STANDARD FRONT RECEIVERS

|                                 |                                                |     |
|---------------------------------|------------------------------------------------|-----|
| <input type="checkbox"/> PHDR   | Pride HD Drop-in Receivers Only .....          | N/C |
| <input type="checkbox"/> PSAR   | Pride Style #8 Swing Away Receivers Only ..... | N/C |
| <input type="checkbox"/> IAR    | Invacare 70° Swing Away Receivers Only .....   | N/C |
| <input type="checkbox"/> Q70R   | Quickie 70° Swing Away Receivers Only .....    | N/C |
| <input type="checkbox"/> Q70HDR | Quickie 70HD Swing Away Receivers Only .....   | N/C |

### SINGLE-DUAL FUNCTION ELECTRONICS

|                                |                                                       |       |
|--------------------------------|-------------------------------------------------------|-------|
| <input type="checkbox"/> SFCB  | Single/Dual Function Electronics <sup>(1)</sup> ..... | STD   |
| <b>Select Switch:</b>          |                                                       |       |
| <input type="checkbox"/> SPBSO | Single Push Button                                    |       |
| <input type="checkbox"/> STO   | Single Toggle                                         |       |
| <input type="checkbox"/> DPBO  | Dual Push Button                                      |       |
| <input type="checkbox"/> OSO   | Omit Switch                                           |       |
| <b>Optional:</b>               |                                                       |       |
| <input type="checkbox"/> RDE   | Reduced Drive for Elevating Seat <sup>(2)</sup> ..... | \$400 |

|                                |                                                                             |         |
|--------------------------------|-----------------------------------------------------------------------------|---------|
| <input type="checkbox"/> SFCBA | Single/Dual Function Drive Control <sup>(1)</sup><br>HCPCS code E2310 ..... | \$1,795 |
| <b>Optional:</b>               |                                                                             |         |
| <input type="checkbox"/> RDES  | Reduced Drive for Elevating Seat <sup>(2)</sup> .....                       | \$400   |

|                                |                                                                              |         |
|--------------------------------|------------------------------------------------------------------------------|---------|
| <input type="checkbox"/> DFCBD | Dual Function Through Drive Control <sup>(1)</sup><br>HCPCS code E2311 ..... | \$2,495 |
| <b>Optional:</b>               |                                                                              |         |
| <input type="checkbox"/> RDED  | Reduced Drive for Elevating Seat <sup>(2)</sup> .....                        | \$400   |

1. These electronics include drive lock-out and tilt limit.

2. If reduced drive is not selected the system will be in drive lock-out when elevated.

# Ultra Low Maxx CG Power Positioning System

## TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### MULTI- FUNCTION ELECTRONICS

- ☐ MFCBS Multi-Function through Switch<sup>(1)</sup>.....N/C  
Select Switch:  
☐ QPBO Quad Push Button.....N/C  
☐ FWT Four-Way Toggle Switch .....\$150  
☐ EWR 8-Way Rocker .....\$205  
Optional:  
☐ RDEM Reduced Drive for Elevating Seat<sup>(2)</sup> .....\$400
- ☐ ACT4S Multi-Function Act4-LiNX Through Switch<sup>(1)(4)</sup> ...N/C  
Select Switch:  
☐ AQPBO Quad Push Button.....N/C  
☐ AFWTA Four-Way Toggle Switch .....\$150  
Optional:  
☐ ARDE Reduced Drive for Elevating Seat<sup>(2)</sup> .....\$400
- ☐ MFCBD Multi-Function M610 Through  
Drive Control<sup>(1)(3)</sup> HCPCS code E2311.....\$2,995  
- Up to 6 Seating Functions  
- On Board Programming Seating Functions  
- Specialty Products (i.e. PSB Power Sliding Back)  
- Includes reduced drive for Elevating Seat
- ☐ ACT4 Multi-Function Act4-LiNX Through  
Drive Control<sup>(1)(4)(5)</sup> HCPCS code E2311.....\$2,995  
- Up to 5 Seating Functions  
- Programming through LAK key/Cell Phone/Laptop  
- Includes reduced drive for Elevating Seat

1. Requires Expandable Controller (EXPC) & Harness for Expandable System (PWH), select on page 2.
2. If reduced drive is not selected system will be in drive lock-out when elevated.
3. These electronics include seating attendant control, drive lock-out, tilt/recline limit and actuator speed control.
4. Only available on Multi-function Seating Systems.
5. The LAK key is not included, please select on page 3.

### ELECTRONIC SWITCH MOUNTING

- ☐ ESML Left Side (Armpad)  
☐ ESMR Right Side (Armpad)  
☐ ESMJL Joystick Left (Mounts to Joystick Bracket)  
☐ ESMJR Joystick Right (Mounts to Joystick Bracket)  
☐ ESMO Do not mount

### ELECTRONICS OPTIONS

#### Additional Switches

- ☐ SPBSA Single Push Button .....\$150  
☐ STA Single Toggle Switch<sup>(4)</sup> .....\$150  
☐ DPBA Dual Push Button<sup>(4)</sup> .....\$150  
☐ QPBA Quad Push Button<sup>(4)</sup> .....\$150  
☐ FWT Four-Way Toggle Switch<sup>(4)</sup> .....\$150  
☐ EWRA 8-Way Rocker<sup>(1)(4)</sup> .....\$205  
☐ TWRA 10-Way Rocker<sup>(1)(3)</sup> .....\$205  
☐ EWPA 8-Way Push Button<sup>(1)(4)</sup> .....\$255  
☐ 10WA 10-Way Switch<sup>(1)(4)</sup> .....\$255  
☐ ASL320R ASL Pulse Mechanical Switch - Red.....\$125  
☐ ASL320G ASL Pulse Mechanical Switch - Green .....\$125  
☐ ASL320BK ASL Pulse Mechanical Switch - Black.....\$125  
☐ SESG Stealth Egg Switch - Green .....\$135  
☐ SESB Stealth Egg Switch - Black.....\$135  
☐ FTSQ Feather Touch - Single Quarter Size .....\$295  
☐ FTMD Feather Touch - Dual Dime Size<sup>(4)</sup> .....\$295  
☐ FTMQ Feather Touch - Quad Dime Size<sup>(4)</sup> .....\$295

#### Specify Use for Additional Switch:

- ☐ STG Seating  
☐ JM Joystick Mode<sup>(2)(3)</sup>  
☐ JOO Joystick On/Off<sup>(2)(3)</sup>

#### Specify Mounting for Additional Switch:

- ☐ ESMLA Left Side (Armpad)  
☐ ESMRA Right Side (Armpad)  
☐ ESMJLA Joystick Left (Mounts to Joystick Bracket)  
☐ ESMJRA Joystick Right (Mounts to Joystick Bracket)  
☐ ESMO Do not mount

#### Miscellaneous

- ☐ MJBS Motion Junction Box for Switch Control<sup>(1)(5)(6)</sup>  
For Mode and/or Emergency Stop.....\$175
- ☐ S46L Single to Dual Switch Adapter .....\$265  
Allows for dual use of a single switch
- ☐ M616S M616 Converter<sup>(1)</sup> HCPCS code E2399.....\$125  
Allows up to 4 mono phono switches into  
Multi-Function box
- ☐ LJS Lazarus Jump Starter HCPCS code A9270.....\$150
- ☐ APPS Accessory Port Power Supply<sup>(1)</sup>  
3 amp capacity.....\$195
- ☐ TSPMH 12 Piece Switch Mounting Hardware  
Kit (additional kits).....\$39

1. Not available with Single-Dual Function Electronics.
2. IN200 is required with REM210, REM211, REM216 joysticks.
3. Not available with REM110 joystick.
4. Not available for Joystick Mode or Joystick On/Off. Must Select Seating (STG)
5. Not available on Multi Function Act4-LiNX Electronics.
6. Cannot be used for power On/Off.

## Ultra Low Maxx CG Power Positioning System

### TDX-SP2

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

#### ARTICULATING TRAYS<sup>(1)</sup>

- ☐ USAT Universal Style Articulating Tray ..... \$1,350  
14"W x Adjustable 6"-10"D
- Add:
- ☐ TMA Trilogy Mounting Adapter ..... \$50
- ☐ AVT Articulating Vent (LTV) Tray ..... \$1,350
- ☐ AVBH Articulating Vent Bag Holder<sup>(2)</sup> ..... \$1,350
- Add:
- ☐ ABMC Accessory Bag Mounting and Clips ..... \$100

1. The user weight limit is 225 lbs. The maximum weight capacity of the vent tray is 35 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. Not available with TRRO and TRRO-E.
2. Motion Concepts does not supply the Vent Bag.

#### GIMBAL VENT TRAYS<sup>(1)</sup>

- ☐ GUT Gimbal Universal Tray ..... \$1,350  
14"W x Adjustable 6"-10"D
- Add:
- ☐ GTMA Trilogy Mounting Adapter ..... \$50
- ☐ GVT Gimbal Vent (LTV) Tray ..... \$1,350
- ☐ GVBH Gimbal Vent Bag Holder<sup>(2)</sup> ..... \$1,350
- Add:
- ☐ GABMC Accessory Bag Mounting and Clips ..... \$100

1. The Gimbal Vent Trays are not available with Tilt Only and/or Elevate Only Systems. The user weight limit is 200 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. The maximum weight capacity of the vent tray is 20 lbs. Not available with TRRO and TRRO-E.
2. Motion Concepts does not supply the Vent Bag.

#### O2 HOLDERS<sup>(1)</sup>

- ☐ NAO2H Non-Articulating O2 Gas Holder ..... \$395
- ☐ AOHG Articulating O2 Gas Holder Cylinder Style<sup>(2)(3)</sup> ..... \$1,250
- ☐ AOHL Articulating Liquid O2 Holder (6.5" x 4.5")<sup>(2)(3)</sup> ..... \$1,250
- ☐ NAIPL Non-Articulating IV Pole User's Left<sup>(4)</sup> ..... \$395
- ☐ NAIPR Non-Articulating IV Pole User's Right<sup>(4)</sup> ..... \$395
- ☐ AIVP Articulating IV Pole<sup>(2)(3)(4)</sup> ..... \$1,250

1. The user weight limit is 225 lbs. The maximum weight capacity of the O2 holder is 20 lbs. Not available with TRRO and TRRO-E.
2. Required for use on Recline Systems
3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. The maximum weight capacity of the IV pole is 5 lbs.

#### VENTILATOR INFORMATION

Required when ordering a Vent Tray

- ☐ BIB Built in battery
- ☐ DTSVPS Dealer to supply vent power source

#### REAR STORAGE BASKET<sup>(1)(2)</sup>

- ☐ RSB Removable Rear Storage Basket ..... \$1,295  
Exterior dimensions: 18"W x 6"H x 7"D  
Interior dimensions: 17.8"W x 6"H x 6.8"D  
Interior capacity: 3 gallon  
Max Load: 20 lbs

1. Tall items extending beyond the basket's height may be damaged upon tilting.
2. Not available with TRRO, TRRO-E brackets

#### ACCESSORIES

- ☐ MUSB Motion Dual USB Charger ..... N/C
- ☐ UCH Universal Cup Holder<sup>(4)</sup> ..... \$35
- ☐ PH X-Grip Phone Holder ..... \$75
- ☐ AMC Accessory Mount with Clip (QTY \_\_\_\_ ) <sup>(1)</sup> ..... \$50
- ☐ CRVMR Rear View Mirror Right Mounted ..... \$45
- ☐ CRVML Rear View Mirror Left Mounted ..... \$45
- ☐ TAG571850 Wheelchair Cover<sup>(1)</sup> ..... \$105
- ☐ HK Hydration Kit<sup>(2)</sup> ..... \$250
- ☐ MUB Medical Utility Bag w/ Mounting Hardware ..... \$250
- ☐ THS Transfer Handles 4" - 6.5" ..... \$350
- ☐ THL Transfer Handles 9" - 11.5" ..... \$350
- ☐ FFLTL Flip and Fold Away 1/2 Lap Tray  
Left Mounting(14"W x 10"D)<sup>(1)(2)(3)</sup> ..... \$795  
Flip and Fold Away 1/2 Lap Tray  
Right Mounting(14"W x 10"D)<sup>(1)(2)(3)</sup> ..... \$795
- ☐ FFLTR
- ☐ CFH Flag Holder <sup>(1)</sup> ..... \$65
- ☐ Left Mount
- ☐ Right Mount
- ☐ CSCC Cane and Crutch Holder ..... \$395
- ☐ Left Mount
- ☐ Right Mount

1. Not available with TRRO, TRRO-E brackets
2. Not available on joystick side and only available with a flat arm pad or modular arm pad. Not available with Outback Cantilever Arms.
3. The maximum weight capacity of the tray is 10 lbs.
4. Not available on the joystick side

#### BATTERIES

- ☐ M22NF SLD G MK Group 22 NF Gel Battery (ea) Quantity( )  
HCPCS code E2361 ..... \$265
- ☐ M24 SLD G MK Group 24 NF Gel Battery (ea) Quantity( )  
HCPCS code E2363 ..... \$365

#### INSTALLATION

- ☐ MBI Motion Concepts Installation<sup>(1)</sup> ..... N/C

1. Batteries may be required if items on Accessories order form are selected at time of installation, please call to confirm.

#### NOTES

---

---

---

---

---

---

RETAIL \$ \_\_\_\_\_

**SUBMIT**

**PRINT FORM**

**SAVE**

**CLEAR FORM**