





INVACARE® MATRX®

P B

BACK SPECIFICATIONS



MODEL	BACK SHELL DIMENSIONS								OVERALL DIMENSIONS (CONTACT AREA)									
	A WIDTH		B HEIGHT		C LATERAL CONTOUR HEIGHT		D HEIGHT		E INNER MAXIMUM WIDTH		F OUTER MAXIMUM WIDTH		G HEIGHT		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lb	kg	lb	kg
MPB1412	11¾	30	11¼	28.5	3½	9	2	5	10½	27	13	33	13	33	2.6	1.2	300	136
MPB1416	11¾	30	15¼	38.5	6½	16.5	2	5	10½	27	13	33	17	43	3.4	1.5	300	136
MPB1612	13¾	35	11	28	6	15	1¼	3	12½	32	15	38	13	33	3.0	1.4	300	136
MPB1616	13¾	35	15	38	6½	16.5	2¼	5.5	12½	32	15	38	17	43	3.8	1.7	300	136
MPB1620	13¾	35	19	48	10	26	3	7.5	12½	32	15	38	20	51	4.6	2.1	300	136
MPB1712	14¾	37.5	11	28	6	15	1¼	3	13½	34.5	16	41	13	33	3.2	1.5	300	136
MPB1716	14¾	37.5	15	38	8	20	2¼	5.5	13½	34.5	16	41	17	43	4.0	1.8	300	136
MPB1720	14¾	37.5	19	48	11	28	3	7.5	13½	34.5	16	41	20	51	4.8	2.2	300	136
MPB1812	15¾	40	11	28	6	15	1¼	3	14½	37	17	43	13	33	3.4	1.5	300	136
MPB1816	15¾	40	15	38	6½	16.5	2¼	5.5	14½	37	17	43	17	43	4.2	1.9	300	136
MPB1820	15¾	40	19	48	11	28	3	7.5	14½	37	17	43	20	51	5.0	2.3	300	136
MPB2012	17¾	45	11	28	5½	14	1¼	3	16½	42	19	48	13	33	3.8	1.7	300	136
MPB2016	17¾	45	15	38	7½	19	2¼	5.5	16½	42	19	48	17	43	4.6	2.1	300	136
MPB2020	17¾	45	19	48	10	26	3	7.5	16½	42	19	48	20	51	5.4	2.5	300	136

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd.
E. Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: 905-695-0138

© 2014 Invacare Corporation. All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form 12-564 Rev. 06/14 140437