





INVACARE® MATRX®

MX 1

BACK SPECIFICATIONS



MODEL	BACK SHELL DIMENSIONS						OVERALL DIMENSIONS (CONTACT AREA)									
	A		B		C		D		E		F		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY	
	WIDTH		HEIGHT		LATERAL CONTOUR HEIGHT		INNER MAXIMUM WIDTH		OUTER MAXIMUM WIDTH		HEIGHT		lb	kg	lb	kg
IMX1409	13	33	8	20	5½	14	11½	29	14	36	10	26	1.1	0.5	250	113.4
IMX1412	13	33	11	28	6½	16.5	11½	29	14	36	13	33	1.2	0.5	250	113.4
IMX1416	13	33	15	38	6½	16.5	11½	29	14	36	17	43	1.3	0.6	250	113.4
IMX1509	14	36	8	20	5½	14	12½	32	15	38	10	26	1.1	0.5	250	113.4
IMX1512	14	36	11	28	5½	14	12½	32	15	38	13	33	1.2	0.5	250	113.4
IMX1516	14	36	15	38	7	18	12½	32	15	38	17	43	1.4	0.6	250	113.4
IMX1609	15	38	8	20	5½	14	13½	34.5	16	41	10	26	1.2	0.5	250	113.4
IMX1612	15	38	11	28	5½	14	13½	34.5	16	41	13	33	1.4	0.6	250	113.4
IMX1616	15	38	15	38	7	18	13½	34.5	16	41	17	43	1.6	0.7	250	113.4
IMX1709	16	41	8	20	5½	14	14½	37	17	43	10	26	1.2	0.5	250	113.4
IMX1712	16	41	11	28	5½	14	14½	37	17	43	13	33	1.4	0.6	250	113.4
IMX1716	16	41	15	38	7	18	14½	37	17	43	17	43	1.7	0.8	250	113.4
IMX1809	17	43	8	20	5½	14	15½	39.5	18	46	10	26	1.3	0.6	250	113.4
IMX1812	17	43	11	28	5½	14	15½	39.5	18	46	13	33	1.6	0.7	250	113.4
IMX1816	17	43	15	38	7	18	15½	39.5	18	46	17	43	1.9	0.9	250	113.4

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd.
E. Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: 905-695-0138

© 2014 Invacare Corporation. All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form 12-575 Rev. 06/14 140445