





INVACARE® MATRX®

MINI

BACK SPECIFICATIONS



MODEL	BACK SHELL DIMENSIONS						OVERALL DIMENSIONS (CONTACT AREA)									
	A WIDTH		B HEIGHT		C LATERAL CONTOUR HEIGHT		D INNER MAXIMUM WIDTH		E OUTER MAXIMUM WIDTH		F HEIGHT		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lb.	kg.	lb.	kg.
MINI1007	8	20	6	15	2.5	6	7	18	10	25	7	18	0.9	0.41	150	68
MINI1010	8	20	9	23	2.5	6	7	18	10	25	10	25	1	0.45	150	68
MINI1012	8	20	11	28	2.5	6	7	18	10	25	12	30	1.3	0.59	150	68
MINI1107	9	23	6	15	2.5	6	8	20	11	28	7	18	1.1	0.50	150	68
MINI1110	9	23	9	23	2.5	6	8	20	11	28	10	25	1.2	0.54	150	68
MINI1112	9	23	11	28	2.5	6	8	20	11	28	12	30	1.4	0.64	150	68
MINI1207	10	25	6	15	2.5	6	9	23	12	30	7	18	1.1	0.50	150	68
MINI1210	10	25	9	23	2.5	6	9	23	12	30	10	25	1.2	0.54	150	68
MINI1212	10	25	11	28	2.5	6	9	23	12	30	12	30	1.6	0.73	150	68
MINI1214	10	25	13	33	2.5	6	9	23	12	30	14	36	1.8	0.82	150	68
MINI1307	11	28	6	15	2.5	6	10	25	13	33	7	18	1.2	0.54	150	68
MINI1310	11	28	9	23	2.5	6	10	25	13	33	10	25	1.4	0.64	150	68
MINI1312	11	28	11	28	2.5	6	10	25	13	33	12	30	1.7	0.77	150	68
MINI1314	11	28	13	33	2.5	6	10	25	13	33	14	36	2	0.91	150	68
MINI1407	12	30	6	15	2.5	6	11	28	14	36	7	18	1.3	0.59	150	68
MINI1410	12	30	9	23	2.5	6	11	28	14	36	10	25	1.6	0.73	150	68
MINI1412	12	30	11	28	2.5	6	11	28	14	36	12	30	1.8	0.82	150	68
MINI1414	12	30	13	33	2.5	6	11	28	14	36	14	36	2.1	0.95	150	68