


INVACARE® MATRX®

PS (Posture Seat)

CUSHION



The **Invacare Matrx PS** (Posture Seat) cushion is an exceptionally lightweight cushion designed to provide superior positioning, stability and comfort. The design features an anatomically contoured shape with a unique waffled ischial relief recess and an innovative reversible outer cover. It is a great choice for a very wide range of individual needs, including pediatric, geriatric, bariatric and highly active.

Call to order

1.800.333.6900
invacare.com

1.888.433.6818
motionconcepts.com



matrx®
SEATING SERIES



Features and Functional Benefits

- Startex™ outer cover – Reversible can be used smooth or fabric side out
- Moisture resistant and breathable to keep cool and dry, compatible with a wide range of disinfecting wipes
- Includes hook and loop strips
- Startex™ inner cover for extra moisture protection
- Molded high resilient polyurethane foam made with Ultra-fresh™ anti-microbial
- Rear cushion radius and waterfall front edge
- Gently contoured to provide neutral leg position
- Unique waffled relief zone allows immersion in ischial area
- Designed to provide off-loading to hips and thighs for support and stability
- Pre-ischial contour to help prevent sliding

HCPSC Code: Pediatric/Standard - E2605
Heavy-Duty - E2606



POK – Pelvic Obliquity Kit
(2 pieces) accommodates up to one inch of obliquity



Cushion Rigidizer – lightweight and low profile

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Motion Concepts
www.motionconcepts.com
info@motionconcepts.com

USA
700 Ensinger Rd., Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905)-695-0134
Fax: 905-695-0138

© 2023
Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Matrix is a registered trademark of Motion Concepts. Otto Bock is a trademark of Otto Bock HealthCare. Specifications are subject to change without notification.
Form No. 18-296 rev. 1/12
160959



Model #	Width	Depth	Height	Height at Leg Abductor/Adductor	Product Trough	Weight	Limited Capacity
Warranty							
PS (Pediatric)	10" to 14"	10", 12", 14", 16"	3.5"/3.5"	2.75"	1 lb. and up	150 lb.	2 Years
PS (Standard)	16", 17", 18", 20"	16", 18", 20"	3.5"/3.5"	2.75"	2 lb. and up	300 lb.	2 Years
PS (Heavy-Duty)	22", 24", 26"	18", 20", 22"	3.5"/3.5"	3.00"	4.2 lb. and up	600 lb.	2 Years

