

INVACARE® MATRX®

MX2

BACK



The **Invacare® Matrx® MX2 Back** — Designed for active users who demand the best. Carbon fiber specially engineered to provide unmatched stability for peak performance.

Call to order

1.800.333.6900
invacare.com

1.888.443.6818
motionconcepts.com



matrx®
SEATING SERIES



Back Shell

- Carbon fibers oriented to provide flex and suspension
- Shaped for perfect fit, maximum stability for function
- Cutouts for breathability
- Slotted for easy height adjustment
- Ultra-lightweight

3" contour depth



Models & Sizes

- 14", 15", 16", 17" & 18" wide x 6", 9", 12", 16" high
- Breathable Spacetex cover, 9" and higher models with leatherette side panels
- Includes sacral/lumbar positioning pad
- Optional privacy flap available



Hardware Options

- Height, depth, angle and rotation adjustable for precise positioning
- Width adjustable to fit chairs 1" wider or narrower than back width
- Fixed Mount Hardware – constructed with carbon fiber plate for added suspension, minimum weight
- Quick Release Hardware – allows easy removal for portability

Fixed Mount



Quick Release



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com
info@motionconcepts.com

USA
700 Ensminger Rd., Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905)-695-0134
Fax: 905-695-0138

© 2017
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed to
Invacare Corporation unless
otherwise noted. Matrx
is a registered trademark
of Motion Concepts.
Specifications are subject
to change without
notification.
Form No. 17-050
170217



Model #	Height	Width	Depth* Adjustment	Angle Adjustment	Product Weight	Weight Capacity	Limited Warranty
IMX2(WH)-CFH	6", 9", 12" & 16"	14" -18"	2"	30° range	1.5 lb. and up	250 lb.	2 Year
IMX2(WH)-QRH			4"	+/- 15°	(with hardware)		

*Limitation of depth/angle adjustment if used on chairs narrower than model width.

